



666 - Afinion Glycohemoglobin	AFF-3																								
Glycohemoglobin	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%
PG: Afinion Glycohemoglobin	157	6.36	0.76	4.10-8.70	11.88																				
747 - Auto Diff 7 (AcT 5 Diff & Horiba 5-Part Diff)	BCX-11					BCX-12					BCX-13					BCX-14					BCX-15				
White Blood Cell Count	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%
PG: Beckman Coulter AcT 5 diff	23	9.70	0.29	8.20-11.20	3.09	20	6.90	0.47	5.90-7.90	6.81	22	24.50	1.47	20.80-28.20	6.00	23	17.71	0.49	15.00-20.40	2.77	23	2.59	0.06	2.20-3.00	2.31
PG: Horiba Pentra 60	16	9.66	0.35	8.20-11.20	3.61	16	6.80	0.21	5.80-7.80	3.09	16	24.90	0.71	21.20-28.60	2.85	16	17.58	0.57	15.00-20.20	3.24	16	2.56	0.10	2.20-3.00	3.85
PG: Horiba Pentra 60C+	6	9.70	0.13	8.20-11.20	1.34	6	6.85	0.23	5.90-7.90	3.33	6	25.10	0.47	21.30-28.90	1.83	6	17.63	0.38	15.00-20.20	2.16	6	2.63	0.05	2.20-3.00	1.92
SG: Select Horiba Instruments (5-part Diff)	22	9.67	0.31	8.20-11.20	3.20	22	6.81	0.21	5.80-7.80	3.09	22	24.95	0.65	21.30-28.80	2.60	22	17.60	0.52	15.00-20.20	2.95	22	2.58	0.09	2.20-3.00	3.46
SG: International Horiba Instruments	25	9.70	0.31	8.20-11.20	3.20	25	6.84	0.22	5.80-7.80	3.24	25	24.92	0.65	21.20-28.60	2.61	25	17.60	0.51	15.00-20.20	2.90	25	2.60	0.10	2.20-3.00	3.85
Lymphocyte	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%
PG: Beckman Coulter AcT 5 diff	23	22.80	2.16	16.30-29.30	9.47	19	35.02	3.55	24.40-45.60	10.14	22	36.99	2.84	28.50-45.50	7.68	21	15.04	2.01	9.00-21.00	13.40	22	23.28	3.47	12.90-33.70	14.89
PG: Horiba Pentra 60	16	25.71	2.02	19.60-31.80	7.86	16	42.28	3.95	30.40-54.10	9.34	16	40.21	2.58	32.50-47.90	6.42	16	17.79	3.08	8.60-27.00	17.30	16	26.56	2.63	18.70-34.50	9.89
PG: Horiba Pentra 60C+	6	23.63	3.72	12.40-34.80	15.76	6	39.62	3.39	29.40-49.80	8.56	6	37.77	2.27	31.00-44.60	6.01	6	15.47	3.59	4.70-26.30	23.16	6	28.08	1.75	22.90-33.40	6.23
SG: Select Horiba Instruments (5-part Diff)	22	25.15	2.67	17.10-33.10	10.64	22	41.55	3.92	29.80-53.40	9.42	22	39.54	2.69	31.40-47.60	6.81	22	17.15	3.31	7.30-27.10	19.24	22	26.98	2.48	19.60-34.40	9.19
SG: International Horiba Instruments	25	25.80	3.20	16.20-35.40	12.40	25	42.38	4.90	27.70-57.10	11.56	25	40.43	3.96	28.50-52.30	9.80	25	18.33	4.66	4.30-32.30	25.46	25	27.33	3.06	18.10-36.50	11.21
Monocyte	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%
PG: Beckman Coulter AcT 5 diff	21	0.99	0.47	0.00-2.42	47.98	18	1.29	0.53	0.00-2.87	40.78	22	1.72	0.62	0.00-3.58	35.99	23	1.57	0.85	0.00-4.13	54.27	19	0.37	0.19	0.00-0.94	51.62
PG: Horiba Pentra 60	16	0.56	0.30	0.00-1.46	53.39	16	0.95	0.39	0.00-2.12	41.05	16	1.46	0.47	0.06-2.86	31.85	15	0.84	0.28	0.00-1.68	33.33	15	0.27	0.13	0.00-0.67	49.26
PG: Horiba Pentra 60C+	5	0.52	0.18	0.00-1.06	34.42	6	0.88	0.41	0.00-2.10	46.25	6	1.47	0.56	0.00-3.17	38.44	6	0.98	0.20	0.37-1.59	20.82	5	0.60	0.50	0.00-2.09	82.83
SG: Select Horiba Instruments (5-part Diff)	21	0.55	0.27	0.00-1.36	49.27	22	0.93	0.39	0.00-2.09	41.51	22	1.46	0.48	0.02-2.90	32.88	21	0.88	0.26	0.09-1.67	30.00	19	0.34	0.27	0.00-1.15	79.71
SG: International Horiba Instruments	24	0.53	0.26	0.00-1.31	49.06	25	0.97	0.38	0.00-2.12	39.38	25	1.42	0.48	0.00-2.85	33.52	24	0.86	0.25	0.10-1.62	29.53	22	0.32	0.26	0.00-1.09	80.00
Granulocyte	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%
PG: Beckman Coulter AcT 5 diff	21	61.62	2.17	55.10-68.10	3.52	17	52.78	3.69	41.70-63.90	6.99	20	49.83	3.44	39.50-60.10	6.91	21	67.77	3.42	57.50-78.10	5.04	21	52.17	3.37	42.10-62.30	6.46
PG: Horiba Pentra 60	16	65.73	3.13	56.30-75.10	4.76	16	50.45	4.21	37.90-63.10	8.32	16	52.79	2.97	43.90-61.70	5.63	16	75.95	4.58	62.20-89.60	6.03	16	59.75	2.10	53.50-66.10	3.49
PG: Horiba Pentra 60C+	6	68.10	2.92	59.30-76.90	4.29	6	52.80	2.82	44.30-61.30	5.34	6	55.18	2.14	48.80-61.60	3.88	6	79.35	3.32	69.30-89.30	4.19	6	58.90	1.65	54.00-63.80	2.80
SG: Select Horiba Instruments (5-part Diff)	22	66.38	3.19	56.80-76.00	4.80	22	51.09	3.96	39.20-63.00	7.75	22	53.44	2.93	44.60-62.20	5.49	22	76.88	4.47	63.50-90.30	5.81	22	59.52	1.98	53.60-65.40	3.33
SG: International Horiba Instruments	25	65.69	3.63	54.80-76.60	5.53	25	50.38	4.58	36.70-64.10	9.09	25	52.64	3.76	41.30-63.90	7.15	25	75.80	5.21	60.20-91.40	6.87	25	59.22	2.35	52.20-66.30	3.97
Eosinophil	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%
PG: Beckman Coulter AcT 5 diff	19	6.57	1.32	2.60-10.60	20.00	17	6.54	0.98	3.60-9.40	15.08	22	6.18	0.76	3.90-8.50	12.26	23	4.61	0.94	1.80-7.40	20.43	19	10.68	1.59	5.90-15.50	14.86
PG: Horiba Pentra 60	14	6.94	1.28	3.00-10.80	18.70	14	6.16	0.97	3.30-9.10	15.65	16	5.28	0.95	2.50-8.10	17.92	15	4.52	0.76	2.20-6.80	16.89	13	12.31	1.75	7.10-17.60	14.23
PG: Horiba Pentra 60C+	6	7.27	1.56	2.60-12.00	21.37	5	6.94	1.22	3.20-10.60	17.68	6	5.30	0.74	3.10-7.50	13.96	6	3.68	0.56	2.00-5.40	15.14	5	12.08	1.49	7.60-16.60	12.31
SG: Select Horiba Instruments (5-part Diff)	20	7.04	1.34	3.00-11.00	19.14	19	6.36	1.07	3.20-9.60	16.72	22	5.28	0.88	2.70-7.90	16.60	21	4.28	0.80	1.90-6.70	18.60	18	12.24	1.64	7.30-17	



SG: Select Horiba Instruments (5-part Diff)	22	12.48	0.50	11.00-14.00	4.00	22	11.75	0.52	10.20-13.40	4.41	22	12.36	0.54	10.80-14.00	4.35	22	11.83	0.46	10.40-13.20	3.90	22	14.39	0.63	12.50-16.30	4.38
SG: International Horiba Instruments	25	12.45	0.49	11.00-14.00	3.92	25	11.79	0.54	10.20-13.40	4.58	25	12.38	0.51	10.90-13.90	4.11	25	11.86	0.46	10.50-13.30	3.87	25	14.38	0.63	12.50-16.30	4.38
Platelet Count	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%
PG: Beckman Coulter AcT 5 diff	23	90.48	5.49	68.00-113.00	6.11	20	314.15	12.50	236.00-393.00	3.98	22	509.23	13.99	382.00-636.00	2.75	23	503.48	14.23	377.00-629.00	2.82	23	73.87	4.87	56.00-93.00	6.62
PG: Horiba Pentra 60	16	88.69	4.45	67.00-111.00	5.06	16	307.00	14.32	230.00-384.00	4.66	16	496.38	22.02	372.00-620.00	4.44	16	499.88	12.52	375.00-625.00	2.50	16	73.94	4.91	56.00-93.00	6.62
PG: Horiba Pentra 60C+	6	84.17	2.93	63.00-105.00	3.45	6	299.17	11.09	224.00-374.00	3.71	6	484.83	13.45	364.00-606.00	2.78	6	488.17	17.32	366.00-610.00	3.55	6	69.33	6.74	52.00-86.00	9.71
SG: Select Horiba Instruments (5-part Diff)	22	87.45	4.52	65.00-109.00	5.17	22	304.86	13.73	229.00-381.00	4.49	22	493.23	20.43	370.00-616.00	4.14	22	496.68	14.56	373.00-621.00	2.94	22	72.68	5.69	55.00-91.00	7.81
SG: International Horiba Instruments	25	87.44	4.25	65.00-109.00	4.94	25	305.28	13.13	229.00-381.00	4.30	25	493.76	19.24	371.00-618.00	3.89	25	497.44	13.81	373.00-621.00	2.78	25	72.00	5.71	54.00-90.00	7.92
702 - Blood Gases	BG-11					BG-12					BG-13					BG-14					BG-15				
Blood Gas - pCO2	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%
PG: Intl Siemens RAPIDPoint 400/405/500 System	5	77.05	9.83	71.00-83.00	12.73	5	50.60	6.22	46.00-56.00	12.16	5	50.55	3.46	46.00-56.00	6.86	5	46.20	2.55	41.00-51.00	5.43	5	72.65	6.15	67.00-79.00	8.49
PG: Radiometer ABL80 FLEX CO-OX	5	88.60	4.53	82.00-96.00	5.06	5	59.55	1.75	55.00-65.00	3.00	5	56.55	1.89	52.00-62.00	3.33	5	53.72	1.23	49.00-59.00	2.22	5	87.32	2.58	80.00-94.00	2.99
Blood Gas - pH	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%
PG: Radiometer ABL80 FLEX CO-OX	5	6.87	0.01	6.83-6.91	0.13	5	7.28	0.01	7.24-7.32	0.11	5	7.26	0.01	7.22-7.30	0.14	5	7.26	0.00	7.22-7.30	0.00	5	6.88	0.01	6.84-6.92	0.07
AG: All Group	6	6.88	0.01	6.84-6.92	0.22	6	7.29	0.02	7.25-7.33	0.33	6	7.27	0.02	7.23-7.31	0.32	6	7.26	0.01	7.22-7.30	0.15	6	6.88	0.01	6.84-6.92	0.22
Blood Gas - pO2	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%
PG: Intl Siemens RAPIDPoint 400/405/500 System	5	56.55	3.61	46.00-68.00	6.32	5	226.30	37.76	113.00-339.00	16.73	5	184.15	25.24	108.00-260.00	13.70	5	75.05	2.76	67.00-83.00	3.73	5	56.85	2.62	49.00-65.00	4.56
PG: Radiometer ABL80 FLEX CO-OX	5	34.25	3.59	23.00-45.00	10.59	5	220.75	10.15	191.00-251.00	4.57	5	181.75	7.80	159.00-205.00	4.29	5	53.50	3.11	45.00-63.00	5.74	5	35.00	7.53	13.00-58.00	21.43
714 - Complete Cardiac Markers;715 - Basic Cardiac Markers;711 - BNP;720 - Cardiac Markers I;740 - D-Dimer	CDM-11N					CDM-12N					CDM-13N					CDM-14N					CDM-15N				
BNP	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%
PG: Abbott i-Stat	23	2249.52	151.56	1794.80-2704.20	6.74	23	45.52	4.17	33.00-58.00	9.16	23	4064.65	219.07	3407.50-4721.90	5.39	23	1203.43	79.69	964.30-1442.50	6.62	23	5000.00	0.00	5000.00-5000.00	0.00
PG: AlereTriage	27	909.00	115.80	561.60-1256.40	12.74	31	55.85	7.46	33.50-78.30	13.35	30	1645.00	189.09	1077.70-2212.30	11.49	31	496.61	65.66	299.60-693.60	13.22	24	2667.50	363.91	1575.80-3759.20	13.64
AG: All Group	52	1516.31	688.73	0.00-3582.50	45.42	56	50.67	8.94	23.90-77.50	17.63	54	2715.01	1224.41	0.00-6388.20	45.10	55	797.68	356.57	0.00-1867.40	44.70	48	3833.75	1205.77	216.50-7451.10	31.45
CK-MB (Quantitative)	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%
PG: Abbott i-Stat	5	21.70	0.00	21.70-21.70	0.00	5	2.80	0.00	2.80-2.80	0.00	5	37.40	0.00	37.40-37.40	0.00	5	13.80	0.00	13.80-13.80	0.00	5	77.30	0.00	77.30-77.30	0.00
PG: AlereTriage	35	9.91	1.83	4.40-15.40	18.48	35	1.55	0.39	0.40-2.80	24.38	36	17.76	2.61	10.00-25.60	14.66	36	5.49	0.94	2.70-8.30	17.09	36	30.31	7.53	7.70-52.90	24.82
PG: TOSOH AIA	5	32.00	0.00	32.00-32.00	0.00	5	4.40	0.00	4.40-4.40	0.00	5	62.20	0.00	62.20-62.20	0.00	5	20.00	0.00	20.00-20.00	0.00	5	119.20	0.00	119.20-119.20	0.00
PG: Intl OCD Vitros 250,350,400,500,750,950	5	20.60	0.00	20.60-20.60	0.00	5	4.10	0.00	4.10-4.10	0.00	5	36.00	0.00	36.00-36.00	0.00	5	11.90	0.00	11.90-11.90	0.00	5	66.60	0.00	66.60-66.60	0.00
Myoglobin (Quantitative)	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%
PG: AlereTriage	30	131.89	22.99	62.90-200.90	17.43	30	29.80	5.30	13.90-45.70	17.79	30	209.20	32.32	112.20-306.20	15.45	30	83.32	12.28	46.50-120.10	14.74	30	325.53	44.64	191.60-459.40	13.71
AG: All Group	32	132.61	22.96	63.70-201.50	17.32	32	29.95	5.29	14.10-45.90	17.63	32	212.22	35.11	106.90-317.50	16.55	32	83.32	12.37	46.20-120.40	14.85	32	338.59	69.60	129.80-547.40	20.56
Troponin I (Quantitative)	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%
PG: Abbott i-Stat	35	14.46	1.82	9.02-19.91	12.55	35	0.30	0.03	0.20-0.40	10.89	27	28.83	3.82	17.36-40.30	13.26	27	6.84	0.96	3.96-9.72	14.02	27	43.83	5.61	27.01-60.65	12.79
PG: AlereTriage	36	2.41	0.59	0.65-4.17	24.35	36	0.05	0.00	0.04-0.06	3.40	36	6.68	0.82	4.23-9.13	12.22	34	0.76	0.14	0.35-						



SG: Abaxis Piccolo Instruments	50	122.84	1.58	119.00-127.00	1.30	12	142.92	1.73	139.00-147.00	1.19	12	157.17	1.80	153.00-161.00	1.15	12	132.42	2.23	128.00-136.00	1.67	12	149.25	2.80	145.00-153.00	1.88
SG: Select Abbott Architect Instruments	8	117.75	0.71	114.00-122.00	0.59	9	138.00	0.71	134.00-142.00	0.51	9	150.89	1.17	147.00-155.00	0.79	9	128.00	0.87	124.00-132.00	0.70	9	143.50	0.94	140.00-148.00	0.63
SG: International OCD Vitros Instruments	16	122.81	1.28	119.00-127.00	1.06	16	144.31	2.18	140.00-148.00	1.53	13	158.31	2.43	154.00-162.00	1.52	16	133.56	1.93	130.00-138.00	1.42	16	150.69	2.96	147.00-155.00	1.99
SG: International Siemens Instruments	18	118.72	3.23	115.00-123.00	2.69	18	137.67	2.17	134.00-142.00	1.59	18	149.11	2.08	145.00-153.00	1.41	18	128.06	2.46	124.00-132.00	1.95	18	142.44	2.23	138.00-146.00	1.55
SG: International Abaxis Piccolo Instruments	53	122.81	1.57	119.00-127.00	1.30	16	142.75	1.53	139.00-147.00	1.05	16	157.13	1.59	153.00-161.00	1.02	16	132.63	2.00	129.00-137.00	1.50	16	149.19	2.48	145.00-153.00	1.68
SG: Select Gen Chem Instruments	35	117.09	1.46	113.00-121.00	1.28	36	137.00	2.08	133.00-141.00	1.53	36	149.94	2.12	146.00-154.00	1.40	36	127.52	1.35	124.00-132.00	1.09	36	142.77	1.62	139.00-147.00	1.12
AG: All Group	136	119.19	3.33	115.00-123.00	2.77	109	138.50	3.97	135.00-143.00	2.88	95	151.65	3.34	148.00-156.00	2.17	109	128.88	3.75	125.00-133.00	2.87	109	144.56	4.00	141.00-149.00	2.76
Free T3	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%
PG: Beckman ACCESS	6	2.97	0.73	0.80-5.20	24.33	6	4.83	0.58	3.10-6.50	12.08	6	5.02	0.24	4.30-5.70	4.80	6	4.07	0.32	3.10-5.10	7.80	6	4.68	0.27	3.90-5.50	5.74
PG: TOSOH AIA	5	4.84	1.15	1.40-8.30	23.96	5	9.79	0.27	9.00-10.60	2.76	5	11.20	0.42	9.90-12.50	3.75	5	8.40	0.28	7.60-9.20	3.33	5	10.07	0.56	8.40-11.80	5.54
SG: Select Beckman Coulter Instruments (Chem)	9	2.85	0.69	0.70-4.90	24.64	9	4.70	0.57	3.00-6.40	12.13	9	5.01	0.28	4.20-5.80	5.60	9	4.02	0.41	2.80-5.20	10.25	9	4.65	0.23	4.00-5.40	4.89
SG: Select Abbott Architect Instruments	5	3.78	0.30	2.90-4.70	7.89	5	4.80	0.42	3.50-6.10	8.75	5	4.85	0.31	4.00-5.80	6.33	5	4.63	0.47	3.20-6.00	10.00	5	4.73	0.38	3.60-5.80	8.09
PG: Intl OCD Vitros ECI/ECiQ	5	6.63	1.66	1.70-11.50	25.00	5	16.30	1.61	11.50-21.10	9.88	5	17.70	1.76	12.40-23.00	9.94	5	14.97	2.06	8.80-21.20	13.73	5	16.73	1.72	11.50-21.90	10.30
PG: OCD Vitros 5600	5	7.80	0.00	7.80-7.80	0.00	5	15.00	0.00	15.00-15.00	0.00	5	16.70	0.00	16.70-16.70	0.00	5	13.30	0.00	13.30-13.30	0.00	5	15.40	0.00	15.40-15.40	0.00
AG: All Group	16	3.01	0.73	0.80-5.20	24.33	16	4.86	0.69	2.80-7.00	14.08	16	5.16	0.62	3.30-7.10	11.92	16	4.30	0.56	2.60-6.00	13.02	16	4.86	0.56	3.20-6.60	11.43
Potassium (K+)	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%
PG: Abaxis Piccolo (Non-Waived)	9	2.83	0.16	2.30-3.30	5.71	12	4.72	0.23	4.20-5.20	4.89	12	5.88	0.21	5.40-6.40	3.56	11	3.82	0.23	3.30-4.30	6.05	12	5.21	0.26	4.70-5.70	5.00
PG: Abaxis Piccolo (Waived)	37	2.80	0.14	2.30-3.30	5.00																				
PG: Abbott Architect	6	2.97	0.05	2.50-3.50	1.67	6	4.63	0.05	4.10-5.10	1.09	6	5.70	0.06	5.20-6.20	1.05	6	3.80	0.00	3.30-4.30	0.00	6	5.05	0.06	4.60-5.60	1.18
PG: Alfa Was. Ace/NexCT/Alera/Axcel / Alfa Wasserman Reagent	25	2.90	0.13	2.40-3.40	4.14	25	4.60	0.09	4.10-5.10	1.96	26	5.74	0.28	5.20-6.20	4.91	25	3.75	0.08	3.30-4.30	2.11	25	5.11	0.13	4.60-5.60	2.55
PG: Clinical Data Envoy 500	6	2.87	0.08	2.40-3.40	2.76	6	4.52	0.13	4.00-5.00	2.89	6	5.57	0.20	5.10-6.10	3.57	6	3.70	0.06	3.20-4.20	1.62	6	4.98	0.07	4.50-5.50	1.60
PG: Siemens Dimension/AR/ES/RxL/Xpand/EXL	16	2.89	0.05	2.40-3.40	1.72	16	4.56	0.06	4.10-5.10	1.30	16	5.58	0.05	5.10-6.10	0.89	16	3.71	0.06	3.20-4.20	1.62	16	4.99	0.06	4.50-5.50	1.20
PG: Beckman AU	15	2.91	0.03	2.40-3.40	1.03	15	4.51	0.06	4.00-5.00	1.33	15	5.53	0.06	5.00-6.00	1.09	15	3.72	0.04	3.20-4.20	1.08	15	4.94	0.05	4.40-5.40	1.02
PG: Roche Cobas Integra / Roche Reagent	6	2.93	0.05	2.40-3.40	1.72	6	4.60	0.09	4.10-5.10	1.96	6	5.67	0.12	5.20-6.20	2.11	6	3.75	0.06	3.30-4.30	1.32	6	5.05	0.11	4.60-5.60	1.96
PG: OCD Vitros 250,350,400,500,750,950	10	3.01	0.03	2.50-3.50	1.00	10	4.82	0.06	4.30-5.30	1.25	10	6.02	0.10	5.50-6.50	1.67	10	3.92	0.06	3.40-4.40	1.54	10	5.34	0.08	4.80-5.80	1.51
SG: Select Beckman Coulter Instruments (Chem)	18	2.90	0.04	2.40-3.40	1.38	18	4.51	0.07	4.00-5.00	1.33	18	5.55	0.07	5.00-6.00	1.27	18	3.71	0.06	3.20-4.20	1.35	18	4.95	0.05	4.40-5.40	1.02
SG: Select OCD Vitros Instruments	13	2.98	0.07	2.50-3.50	2.33	13	4.78	0.11	4.30-5.30	2.29	13	5.97	0.15	5.50-6.50	2.50	13	3.88	0.11	3.40-4.40	2.82	13	5.28	0.15	4.80-5.80	2.83
SG: Abaxis Piccolo Instruments	46	2.81	0.14	2.30-3.30	5.00	12	4.72	0.23	4.20-5.20	4.89	12	5.88	0.21	5.40-6.40	3.56	11	3.82	0.23	3.30-4.30	6.05	12	5.21	0.26	4.70-5.70	5.00
SG: Select Abbott Architect Instruments	9	2.97	0.05	2.50-3.50	1.67	9	4.63	0.05	4.10-5.10	1.09	9	5.70	0.05	5.20-6.20	0.88	9	3.80	0.00	3.30-4.30	0.00	9	5.07	0.05	4.60-5.60	0.98
SG: International OCD Vitros Instruments	16	2.98	0.07	2.50-3.50	2.33	16	4.78	0.10	4.30-5.30	2.08	16	5.98	0.14	5.50-6.50	2.33	16	3.89	0.10	3.40-4.40	2.56	16	5.30	0.14	4.80-5.80	2.64
SG: International Siemens Instruments	18	2.88	0.06	2.40-3.40	2.07	18	4.55	0.06	4.10-5.10	1.30	18	5.57	0.06	5.10-6.10	1.07	18	3.70	0.06	3.20-4.20	1.62	18	4.99	0.06	4.50-5.50	1.20
SG: International Abaxis Piccolo Instruments	50	2.81	0.14	2.30-3.30	5.00	16	4.73	0.22	4.20-5.20	4.68	16	5.89	0.20	5.40-6.40	3.39	15	3.80	0.23	3.30-4.30	6.05	16	5.18	0.23	4.70-5.70	4.42
SG: Select Gen Chem Instruments	35	2.91	0.06	2.40-3.40	2.07	35	4.54	0.09	4.00-5.00	2.00	35	5.59	0.12	5.10-6.10	2.14	35	3.73	0.06	3.20-4.20	1.62	35	4.98	0.07	4.50-5.50	1.40
AG: All Group	141	2.88	0.12	2.40-3.40	4.14	107	4.61	0.14	4.10-5.10	3.04	108	5.71	0.22	5.20-6.20	3.86					107	5.08	0.17	4.60-5.60	3.33	
Chloride (CL-)	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%
PG: Abaxis Piccolo (Non-Waived)	12	92.33	1.88	87.00-97.00	2.07	12	101.50	1.73	97.00-107.00	1.67	12	107.42	1.78	102.00-112.00	1.68	12	96.92	2.06	92.00-102.00	2.16	12	104.00	2.17	99.00-109.00	2.12
PG: Abaxis Piccolo (Waived)	38	91.50	2.04	87.00-97.00	2.17																				
PG: Abbott Architect	6	88.67	0.82	85.00-93.00	0.90	6	98.33	1.03	93.00-103.00	1.02	6	104.83	1.17	100.00-110.00	1.14	6	93.83	0.98	89.00-99.00	1.06	6	101.33	0.82	96.00-106.00	0.79
PG: Alfa Was. Ace/NexCT/Alera/Axcel / Alfa Wasserman Reagent	25	86.46	3.70	82.00-90.00	4.30	25	98.58	2.91	94.00-104.00	2.93	25	107.71	2.94	103.00-113.00	2.69	25	93.35	3.28	88.00-98.00	3.55	25	102.46	3.15	97.00-107.00	3.04
PG: Clinical Data Envoy 500	6	88.50	2.17	85.00-93.00	2.47	6	100.75	2.68	96.00-106.00	2.67	6	110.33	2.73	105.00-116.00	2.45	6	94.45	1.86	89.00-99.00	2.02	6	105.45	2.63	100.00-110.00	2.48
PG: Siemens Dimension/AR/ES/RxL/Xpand/EXL	16	84.81	0.98	81.00-89.00	1.18	16	99.19	0.98	94.00-104.00	1.01	16	107.63	1.20	103.00-113.00	1.11	16	92.13	1.02	87.00-97.00	1.09	16	102.63	1.20	98.00-108.00	1.17
PG: Beckman AU	14	86.93	1.77	83.00-91.00	2.07	14	97.79	1.53	93.00-103.00	1.53	14	104.29	1.82	99.00-109.00	1.73	14	93.14	1.29	88.00-98.00	1.40	14	100.86	1.56	96.00-106.00	1.58
PG: OCD Vitros 250,350,400,500,750,950	10	87.80	1.62	84.00-92.00	1.82	10	98.60	2.07	94.00-104.00	2.12	10	106.00	1.83	101.00-111.00	1.70	10	93.80	1.81	89.00-99.00	1.91	10	102.70	2.06	98.00-108.00	2.04
SG: Select Beckman Coulter Instruments (Chem)	17	87.12	1.69	83.00-91.00	1.95	17	97.82	1.43	93.00-103.00	1.43	17	104.29	1.65	99.00-109.00	1.54	17	93.18	1.24	88.00-98.00	1.29	17	100.88	1.45	96.00-106.00	1.49
SG: Select Roche Instruments / Roche Reagent	6	87.50	1.05	84.00-92.00	1.14	6	98.62	1.76	94.00-104.00	1.82	6	106.52	2.29	102.00-112.00	2.15	6	93.03	1.33	88.00-98.00	1.40	6	102.43	1.98	97.00-107.00	1.96
SG: Select OCD Vitros Instruments	13	88.08	1.66	84.00-92.00	1.93	13	98.77	1.92	94.00-104.00	1.92	13	106.62	2.57	102.00-112.00	2.43	13	94.00	1.68	89.00-99.00	1.81	13	102.92	2.22	98.00-108.00	2.14
SG: Abaxis Piccolo Instruments	50	91.70	2.01	87.00-97.00	2.17	12	101.50	1.73	97.00-107.00	1.67	12	107.42	1.78	102.00-112.00	1.68	12	96.92	2.06	92.00-102.00	2.16	12	104.00	2.17	99.00-109.00	2.12
SG: Select Abbott Architect Instruments	9	88.89	0.78	85.00-93.00	0.90	9	98.33	0.87	93.00-103.0																



PG: Abaxis Piccolo (Non-Waived)	11	12.36	1.21	8.00-16.00	10.00	11	22.18	1.33	18.00-26.00	5.91	11	27.55	0.69	26.00-30.00	2.50	11	17.18	1.25	13.00-21.00	7.65	11	24.55	1.37	21.00-29.00	5.60
PG: Abaxis Piccolo (Waived)	38	12.24	0.63	10.00-14.00	5.00																				
PG: Abbott Architect	6	12.67	1.21	9.00-17.00	9.23	6	20.33	1.51	16.00-25.00	7.50	6	24.83	1.47	21.00-30.00	6.00	6	16.50	1.38	13.00-21.00	8.24	6	21.83	1.84	17.00-27.00	8.18
PG: Alfa Was. Ace/NexCT/Alera/Axcel / Alfa Wasserman Reagent	24	13.45	1.70	8.00-18.00	13.08	25	21.56	2.31	15.00-29.00	10.45	25	26.68	2.71	19.00-35.00	10.00	24	17.12	2.15	10.00-24.00	12.94	25	23.31	2.72	15.00-31.00	11.74
PG: Clinical Data Envoy 500	6	14.62	2.00	9.00-21.00	13.33	6	24.52	2.84	17.00-33.00	11.20	6	30.77	2.59	23.00-39.00	8.39	6	18.68	1.89	13.00-25.00	10.00	6	25.95	2.42	19.00-33.00	9.23
PG: Siemens Dimension/AR/ES/RxL/Xpand/EXL	16	15.74	2.06	10.00-22.00	13.13	16	24.03	1.47	20.00-29.00	6.25	16	28.79	2.34	22.00-36.00	7.93	16	19.83	1.78	15.00-25.00	9.00	16	25.59	2.67	18.00-34.00	10.38
PG: Beckman AU	15	13.07	1.62	8.00-18.00	12.31	15	22.27	1.98	16.00-28.00	9.09	15	27.40	2.61	19.00-35.00	9.63	15	17.67	1.68	13.00-23.00	9.44	15	23.60	2.17	17.00-31.00	9.17
PG: Roche Cobas Integra / Roche Reagent	6	13.22	1.17	9.00-17.00	9.23	6	20.32	1.50	16.00-25.00	7.50	6	24.72	1.95	19.00-31.00	7.60	6	15.95	1.42	12.00-20.00	8.75	6	21.47	1.53	17.00-26.00	7.14
PG: OCD Vitros 250,350,400,500,750,950	10	11.70	0.82	10.00-14.00	6.67	10	18.40	1.78	13.00-23.00	10.00	10	22.20	2.04	16.00-28.00	9.09	10	14.10	0.88	11.00-17.00	6.43	10	18.50	1.78	14.00-24.00	9.47
SG: Select Beckman Coulter Instruments (Chem)	18	13.06	1.59	8.00-18.00	12.31	18	22.17	1.82	17.00-27.00	8.18	18	27.28	2.40	20.00-34.00	8.89	18	17.67	1.53	14.00-23.00	8.33	18	23.72	2.02	18.00-30.00	8.33
SG: Select OCD Vitros Instruments	13	11.69	0.75	10.00-14.00	6.67	13	18.54	1.61	14.00-24.00	8.42	13	22.46	2.02	16.00-28.00	9.09	13	14.38	1.04	11.00-17.00	7.14	13	19.00	1.83	14.00-24.00	9.47
SG: Abaxis Piccolo Instruments	49	12.27	0.79	10.00-14.00	6.67	11	22.18	1.33	18.00-26.00	5.91	11	27.55	0.69	26.00-30.00	2.50	11	17.18	1.25	13.00-21.00	7.65	11	24.55	1.37	21.00-29.00	5.60
SG: Select Abbott Architect Instruments	9	13.00	1.32	9.00-17.00	10.00	9	20.67	1.50	17.00-26.00	7.14	9	25.33	1.80	20.00-30.00	7.20	9	16.78	1.39	13.00-21.00	8.24	9	22.33	2.06	16.00-28.00	9.55
SG: International OCD Vitros Instruments	15	11.67	0.90	9.00-15.00	7.50	15	18.60	1.55	15.00-24.00	7.89	15	22.47	1.89	16.00-28.00	8.64	15	14.40	0.99	11.00-17.00	7.14	15	18.87	1.77	14.00-24.00	9.47
SG: International Abaxis Piccolo Instruments	53	12.19	0.81	10.00-14.00	6.67	15	21.67	1.50	18.00-27.00	6.82	15	27.13	0.99	24.00-30.00	3.70	15	16.87	1.25	13.00-21.00	7.06	15	24.07	1.53	20.00-29.00	6.25
SG: Select Gen Chem Instruments	35	13.39	1.75	8.00-18.00	13.85	35	22.17	2.24	15.00-29.00	10.00	35	27.36	2.80	19.00-35.00	10.37	35	17.63	1.67	13.00-23.00	9.44	34	23.76	2.34	17.00-31.00	9.58
AG: All Group	143	13.12	1.79	8.00-18.00	13.85	106	21.76	2.45	15.00-29.00	10.91	106	26.68	3.00	18.00-36.00	11.11	105	17.30	2.23	10.00-24.00	12.94	105	23.29	2.95	14.00-32.00	12.61
Calcium	Labs	Mean	SD	Range	CV %	Labs	Mean	SD	Range	CV %	Labs	Mean	SD	Range	CV %	Labs	Mean	SD	Range	CV %	Labs	Mean	SD	Range	CV %
PG: Abaxis Piccolo (Non-Waived)	11	7.82	0.15	6.80-8.80	1.92	11	10.10	0.20	9.10-11.10	1.98	11	11.42	0.30	10.40-12.40	2.63	11	8.96	0.18	8.00-10.00	2.00	11	10.65	0.30	9.70-11.70	2.80
PG: Abaxis Piccolo (Waived)	38	7.89	0.19	6.90-8.90	2.41																				
PG: Abbott Architect	6	7.65	0.06	6.70-8.70	0.65	6	10.22	0.07	9.20-11.20	0.78	6	11.87	0.14	10.90-12.90	1.18	6	8.93	0.12	7.90-9.90	1.35	6	10.92	0.12	9.90-11.90	1.10
PG: Alfa Was. Ace/NexCT/Alera/Axcel / Alfa Wasserman Reagent	26	7.95	0.59	6.90-8.90	7.47	26	10.15	0.49	9.20-11.20	4.80	26	11.57	0.41	10.60-12.60	3.53	26	9.22	0.56	8.20-10.20	5.98	26	10.78	0.47	9.80-11.80	4.35
PG: Clinical Data Envoy 500	6	7.42	0.17	6.40-8.40	2.30	6	9.83	0.27	8.80-10.80	2.76	6	11.42	0.28	10.40-12.40	2.46	6	8.54	0.29	7.50-9.50	3.41	6	10.66	0.21	9.70-11.70	1.96
PG: Siemens Dimension/AR/ES/RxL/Xpand/EXL	15	7.45	0.26	6.50-8.50	3.47	15	9.79	0.27	8.80-10.80	2.76	15	11.33	0.35	10.30-12.30	3.10	15	8.64	0.29	7.60-9.60	3.37	15	10.38	0.32	9.40-11.40	3.08
PG: Beckman AU	15	7.52	0.20	6.50-8.50	2.67	15	9.95	0.22	9.00-11.00	2.20	15	11.65	0.79	10.60-12.60	6.81	15	8.78	0.20	7.80-9.80	2.27	15	10.63	0.23	9.60-11.60	2.17
PG: Roche Cobas Integra / Roche Reagent	6	7.65	0.24	6.70-8.70	3.12	6	10.07	0.22	9.10-11.10	2.18	6	11.57	0.31	10.60-12.60	2.67	6	8.82	0.23	7.80-9.80	2.61	6	10.67	0.27	9.70-11.70	2.52
PG: OCD Vitros 250,350,400,500,750,950	10	7.43	0.10	6.40-8.40	1.22	10	10.00	0.11	9.00-11.00	1.10	10	11.44	0.14	10.40-12.40	1.23	10	8.86	0.13	7.90-9.90	1.46	10	10.67	0.14	9.70-11.70	1.31
SG: Select Beckman Coulter Instruments (Chem)	18	7.54	0.19	6.50-8.50	2.53	18	9.96	0.21	9.00-11.00	2.10	18	11.61	0.74	10.60-12.60	6.38	18	8.77	0.18	7.80-9.80	2.05	18	10.62	0.23	9.60-11.60	2.17
SG: Select OCD Vitros Instruments	13	7.45	0.10	6.40-8.40	1.35	13	9.99	0.11	9.00-11.00	1.10	13	11.39	0.19	10.40-12.40	1.67	13	8.82	0.15	7.80-9.80	1.70	13	10.62	0.18	9.60-11.60	1.70
SG: Abaxis Piccolo Instruments	49	7.88	0.18	6.90-8.90	2.28	11	10.10	0.20	9.10-11.10	1.98	11	11.42	0.30	10.40-12.40	2.63	11	8.96	0.18	8.00-10.00	2.00	11	10.65	0.30	9.70-11.70	2.80
SG: Select Abbott Architect Instruments	9	7.62	0.10	6.60-8.60	1.32	9	10.13	0.16	9.10-11.10	1.58	9	11.77	0.19	10.80-12.80	1.61	9	8.84	0.18	7.80-9.80	2.05	9	10.82	0.18	9.80-11.80	1.67
SG: International OCD Vitros Instruments	16	7.44	0.09	6.40-8.40	1.22	16	9.98	0.11	9.00-11.00	1.10	16	11.39	0.17	10.40-12.40	1.49	16	8.81	0.14	7.80-9.80	1.59	16	10.63	0.17	9.60-11.60	1.60
SG: International Abaxis Piccolo Instruments	50	7.88	0.18	6.90-8.90	2.28	12	10.09	0.19	9.10-11.10	1.88	12	11.40	0.29	10.40-12.40	2.54	12	8.94	0.19	7.90-9.90	2.13	12	10.63	0.30	9.60-11.60	2.83
SG: Select Gen Chem Instruments	35	7.54	0.18	6.50-8.50	2.40	35	9.99	0.23	9.00-11.00	2.30	35	11.62	0.55	10.60-12.60	4.74	35	8.76	0.22	7.80-9.80	2.50	35	10.69	0.24	9.70-11.70	2.24
AG: All Group	144	7.71	0.36	6.70-8.70	4.68	106	10.02	0.32	9.00-11.00	3.20	106	11.51	0.43	10.50-12.50	3.74	106	8.89	0.39	7.90-9.90	4.38	106	10.65	0.34	9.70-11.70	3.18
T3 Uptake	Labs	Mean	SD	Range	CV %	Labs	Mean	SD	Range	CV %	Labs	Mean	SD	Range	CV %	Labs	Mean	SD	Range	CV %	Labs	Mean	SD	Range	CV %
PG: OCD Vitros 5600	5	65.00	0.00	65.00-65.00	0.00	5	65.00	0.00	65.00-65.00	0.00	5	65.00	0.00	65.00-65.00	0.00	5	65.00	0.00	65.00-65.00	0.00	5	65.00	0.00	65.00-65.00	0.00
AG: All Group	7	0.53	0.07	0.33-0.73	12.56	7	0.47	0.06	0.28-0.66	13.29	7	0.43	0.05	0.28-0.58	11.52	7	0.49	0.07	0.29-0.68	13.28	7	0.46	0.05	0.30-0.61	11.43
Creatinine	Labs	Mean	SD	Range	CV %	Labs	Mean	SD	Range	CV %	Labs	Mean	SD	Range	CV %	Labs	Mean	SD	Range	CV %	Labs	Mean	SD	Range	CV %
PG: Abaxis Piccolo (Non-Waived)	10	1.13	0.12	0.83-1.43	10.91	12	2.89	0.09	2.50-3.30	3.11	12	4.05	0.16	3.44-4.66	3.90	12	2.05	0.08	1.74-2.36	3.90	12	3.35	0.14	2.85-3.85	4.12
PG: Abaxis Piccolo (Waived)	36	1.19	0.11	0.89-1.49	9.41																				
PG: Abbott Architect	6	0.92	0.04	0.62-1.22	4.57	6	3.01	0.02	2.50-3.50	0.73	6	4.32	0.08	3.67-4.97	1.86	6	1.98	0.03	1.68-2.28	1.50	6	3.57	0.03	3.03-4.11	0.95
PG: Alfa Was. Ace/NexCT/Alera/Axcel / Alfa Wasserman Reagent	25	1.21	0.07	0.91-1.51	5.37	25	2.94	0.09	2.50-3.30	2.89	25	3.89	0.23	3.31-4.47	5.90	25	2.12	0.08	1.80-2.44	3.81	26	3.30	0.26	2.80-3.80	7.94
PG: Clinical Data Envoy 500	6	1.10	0.08	0.80-1.40	7.27	6	2.73	0.11	2.30-3.10	4.07	6	3.78	0.13	3.21-4.35	3.42	6	1.95	0.07	1.65-2.25	3.68	6	3.20	0.12	2.72-3.68	3.75
PG: Siemens Dimension/AR/ES/RxL/Xpand/EXL	11	1.14	0.06	0.84-1.44	5.61	11	2.88	0.11	2.50-3.30	3.79	11	4.02	0.12	3.42-4.62	3.01	11	1.96	0.07	1.66-2.26	3.67	11	3.35	0.13	2.85-3.85	3.94
PG: Siemens Dimension/AR/ES/RxL/Xpand/EXL / IFCC	5	1.13	0.07	0.83-1.43	6.36	5	2.93	0.07	2.49-3.37	2.46	5	4.06	0.11	3.50-4.70	2.76	5	1.99	0.07	1.70-2.30	3.62	5	3.38	0.11	2.90-3.90	3.24
PG: Beckman AU	15	1.07	0.05	0.77-1.37	4.55	15	2.88	0.07	2.50-3.30	2.41	15	3.93	0.11	3.34-4.52	2.67	15	1.98	0.05	1.68-2.28	2.63	15	3.32	0.08	2.82-3.82	2.53
PG: Roche Cobas Integra / Roche Reagent	6	1.09	0.08	0.79-1.39	7.27	6	2.77	0.08	2.40-3.20	2.92	6	3.68	0.11	3.13-4.23	3.07	6	1.95	0.05	1.65-2.25	2.72	6	3.17	0.10	2.69-3.65	3.



PG: Siemens Dimension/AR/ES/RxL/Xpand/EXL	16	65.50	2.58	59.00-73.00	3.94	16	143.56	4.43	130.00-158.00	3.06	16	192.25	5.95	173.00-211.00	3.07	16	105.25	3.21	95.00-116.00	3.05	16	163.63	5.25	148.00-180.00	3.23
PG: Beckman AU	15	65.80	1.26	59.00-73.00	1.97	15	146.00	2.78	131.00-161.00	1.92	15	196.47	4.00	176.00-216.00	2.04	15	106.33	1.95	95.00-117.00	1.89	15	167.13	3.46	150.00-184.00	2.10
PG: Roche Cobas Integra / Roche Reagent	6	65.33	1.37	59.00-72.00	2.15	6	146.83	3.54	132.00-162.00	2.38	6	195.67	3.93	176.00-216.00	1.99	6	106.00	2.19	95.00-117.00	2.08	6	166.50	3.62	150.00-184.00	2.16
PG: OCD Vitros 250,350,400,500,750,950	10	60.60	2.27	55.00-67.00	3.77	10	138.70	2.36	125.00-153.00	1.73	9	187.33	2.96	168.00-206.00	1.60	10	101.10	1.29	91.00-111.00	1.29	10	159.30	2.41	143.00-175.00	1.51
SG: Select Beckman Coulter Instruments (Chem)	18	65.39	1.85	59.00-72.00	2.92	18	145.44	3.43	131.00-160.00	2.34	18	195.89	4.65	176.00-216.00	2.40	18	105.89	2.59	95.00-117.00	2.45	18	166.28	4.25	149.00-183.00	2.59
SG: Select OCD Vitros Instruments	13	60.54	2.33	55.00-67.00	3.77	13	139.62	4.89	126.00-154.00	3.50	12	188.42	5.88	169.00-207.00	3.14	13	101.62	3.15	92.00-112.00	3.14	13	159.85	5.32	144.00-176.00	3.31
SG: Abaxis Piccolo Instruments	49	66.16	1.06	59.00-73.00	1.67	12	141.83	1.70	128.00-156.00	1.20	12	189.25	1.87	170.00-208.00	1.01	12	103.50	1.24	94.00-114.00	1.15	12	161.75	1.54	146.00-178.00	0.93
SG: Select Abbott Architect Instruments	9	62.33	0.71	56.00-68.00	1.13	9	139.67	1.73	126.00-154.00	1.21	9	188.44	2.19	169.00-207.00	1.17	9	101.78	1.56	92.00-112.00	1.57	9	159.44	2.07	143.00-175.00	1.32
SG: International OCD Vitros Instruments	16	60.19	2.29	54.00-66.00	3.83	16	138.56	4.94	125.00-153.00	3.53	15	186.80	6.25	168.00-206.00	3.32	16	100.88	3.24	91.00-111.00	3.17	16	158.75	5.40	143.00-175.00	3.40
SG: International Abaxis Piccolo Instruments	53	66.22	1.05	59.00-73.00	1.67	16	141.89	1.53	128.00-156.00	1.06	16	189.18	1.68	170.00-208.00	0.90	16	103.85	1.28	94.00-114.00	1.25	16	161.85	1.34	146.00-178.00	0.80
SG: Select Gen Chem Instruments	35	64.83	3.41	59.00-72.00	5.23	35	144.86	5.32	131.00-160.00	3.66	35	195.43	7.43	176.00-215.00	3.79	35	105.54	4.25	95.00-117.00	3.96	35	165.60	6.17	149.00-183.00	3.73
AG: All Group	146	65.57	3.02	59.00-73.00	4.55	109	144.46	5.14	130.00-158.00	3.54	109	193.30	7.88	174.00-212.00	4.09	110	105.63	4.13	95.00-117.00	3.87	110	164.25	7.81	148.00-180.00	4.76
Triiodothyronine (T3)	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%
PG: TOSOH AIA	5	2.71	0.01	2.70-2.70	0.37	5	4.12	0.02	4.00-4.20	0.49	5	4.86	0.20	4.30-5.50	4.08	5	3.66	0.06	3.50-3.90	1.62	5	4.47	0.05	4.30-4.70	1.11
PG: OCD Vitros 5600	5	4.10	0.00	4.10-4.10	0.00	5	7.80	0.00	7.80-7.80	0.00	5	7.80	0.00	7.80-7.80	0.00	5	7.00	0.00	7.00-7.00	0.00	5	7.80	0.00	7.80-7.80	0.00
AG: All Group	8	1.08	0.30	0.20-2.00	27.27	8	1.58	0.32	0.60-2.60	20.00	8	1.81	0.38	0.70-2.90	21.11	8	1.32	0.26	0.50-2.10	20.00	8	1.68	0.34	0.70-2.70	20.00
TSH	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%
PG: Abbott Architect	6	1.19	0.08	0.95-1.43	6.72	6	4.24	0.20	3.63-4.85	4.76	6	5.99	0.39	4.83-7.15	6.46	6	2.71	0.14	2.27-3.15	5.35	6	4.98	0.33	4.00-5.96	6.59
PG: Beckman ACCESS	19	1.20	0.07	0.98-1.42	6.00	19	4.36	0.29	3.48-5.24	6.74	19	6.28	0.44	4.97-7.59	6.94	19	2.81	0.20	2.23-3.40	6.94	19	5.15	0.31	4.23-6.07	5.96
PG: Siemens Dimension/AR/ES/RxL/Xpand/EXL	13	1.25	0.12	0.89-1.61	9.52	13	4.46	0.42	3.20-5.72	9.39	11	6.21	0.54	4.59-7.83	8.71	13	2.92	0.29	2.05-3.79	9.97	13	5.29	0.56	3.62-6.96	10.55
PG: Qualigen FastPak Thyroid	19	1.10	0.12	0.74-1.46	10.91	20	4.27	0.39	3.10-5.44	9.13	19	6.08	0.51	4.55-7.61	8.39	20	2.84	0.26	2.07-3.61	9.01	19	5.12	0.43	3.83-6.41	8.40
PG: TOSOH AIA	19	1.62	0.13	1.22-2.02	8.15	19	5.64	0.50	4.13-7.15	8.92	19	7.93	0.72	5.77-10.09	9.07	19	3.70	0.35	2.64-4.76	9.54	19	6.56	0.58	4.80-8.31	8.92
SG: Select Beckman Coulter Instruments (Chem)	24	1.20	0.07	1.00-1.40	5.67	24	4.38	0.27	3.57-5.19	6.19	24	6.33	0.41	5.11-7.55	6.43	24	2.82	0.17	2.29-3.34	6.21	24	5.16	0.28	4.33-5.99	5.39
SG: Select Roche Instruments	6	1.45	0.11	1.13-1.77	7.45	6	4.66	0.30	3.75-5.58	6.55	6	6.37	0.34	5.34-7.40	5.37	6	3.16	0.22	2.50-3.82	6.96	6	5.42	0.33	4.42-6.42	6.13
SG: Select Siemens Immulite Instruments	5	1.50	0.11	1.19-1.82	7.00	5	5.31	0.37	4.20-6.42	6.99	5	7.31	0.43	6.02-8.60	5.90	5	3.45	0.13	3.07-3.83	3.65	5	6.05	0.42	4.80-7.30	6.89
PG: Roche Cobas e 411	6	1.45	0.11	1.13-1.77	7.45	6	4.66	0.30	3.75-5.58	6.55	6	6.37	0.34	5.34-7.40	5.37	6	3.16	0.22	2.50-3.82	6.96	6	5.42	0.33	4.42-6.42	6.13
PG: Abbott Architect i1000 SR	6	1.17	0.06	1.00-1.34	4.79	6	4.22	0.16	3.75-4.68	3.67	6	6.00	0.26	5.21-6.79	4.37	6	2.72	0.12	2.37-3.07	4.30	6	4.96	0.18	4.42-5.50	3.61
SG: Select Abbott Architect Instruments	14	1.18	0.06	0.99-1.37	5.34	14	4.25	0.17	3.75-4.75	3.93	14	6.01	0.29	5.14-6.88	4.84	14	2.73	0.12	2.36-3.10	4.54	14	4.97	0.23	4.27-5.67	4.67
PG: NanoEnTek FrenD System	7	0.61	0.06	0.43-0.79	10.00	6	2.34	0.27	1.54-3.14	11.41	5	3.89	0.32	2.92-4.86	8.33	7	1.53	0.35	0.48-2.58	22.94	6	3.07	0.14	2.66-3.48	4.43
PG: Intl OCD Vitros ECi/ECiQ	5	1.61	0.03	1.53-1.69	1.61	5	5.84	0.20	5.24-6.44	3.42	5	8.10	0.25	7.35-8.85	3.07	5	3.84	0.06	3.66-4.02	1.59	5	6.79	0.18	6.25-7.33	2.64
SG: Select Gen Chem Instruments	16	1.19	0.06	1.00-1.38	5.38	16	4.28	0.19	3.72-4.84	4.37	16	6.09	0.36	5.01-7.17	5.91	16	2.75	0.12	2.38-3.12	4.47	16	4.99	0.23	4.30-5.68	4.59
AG: All Group	101	1.29	0.20	0.68-1.90	15.81	101	4.63	0.59	2.87-6.39	12.68	97	6.53	0.77	4.21-8.85	11.82	103	3.03	0.41	1.79-4.27	13.63	102	5.48	0.71	3.36-7.60	12.92
Urea Nitrogen (BUN)	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%
PG: Abaxis Piccolo (Non-Waived)	12	8.92	0.52	7.00-11.00	5.56	12	22.67	0.65	21.00-25.00	3.04	12	31.75	1.05	29.00-35.00	3.44	12	15.83	0.58	14.00-18.00	3.75	12	26.50	0.91	25.00-29.00	3.33
PG: Abaxis Piccolo (Waived)	37	8.80	0.54	7.00-11.00	5.56																				
PG: Abbott Architect	6	10.00	0.00	8.00-12.00	0.00	6	25.50	0.55	24.00-28.00	1.92	6	34.83	0.75	32.00-38.00	2.29	6	17.50	0.55	16.00-20.00	2.78	6	29.83	1.17	27.00-33.00	4.00
PG: Alfa Was. Ace/NexCT/Alera/Axcel / Alfa Wasserman Reagent	25	9.72	0.61	8.00-12.00	6.00	25	24.28	0.98	22.00-26.00	4.17	26	33.00	2.13	30.00-36.00	6.36	25	17.24	0.66	15.00-19.00	4.12	26	27.50	2.42	25.00-31.00	8.57
PG: Clinical Data Envoy 500	5	9.98	0.04	8.00-12.00	0.00	5	24.42	0.58	22.00-26.00	2.50	5	33.42	1.18	30.00-36.00	3.64	5	17.28	0.55	15.00-19.00	2.94	5	27.88	1.00	25.00-31.00	3.57
PG: Siemens Dimension/AR/ES/RxL/Xpand/EXL	16	9.75	0.93	8.00-12.00	9.00	16	25.00	1.37	23.00-27.00	5.60	16	34.31	1.99	31.00-37.00	5.88	16	17.38	0.96	15.00-19.00	5.88	16	28.69	1.54	26.00-32.00	5.17
PG: Beckman AU	15	10.27	0.46	8.00-12.00	5.00	15	25.13	0.64	23.00-27.00	2.40	15	34.60	1.12	32.00-38.00	3.14	15	17.73	0.80	16.00-20.00	4.44	15	29.07	1.16	26.00-32.00	4.14
PG: Roche Cobas Integra / Roche Reagent	6	9.33	0.52	7.00-11.00	5.56	6	23.67	0.52	22.00-26.00	2.08	6	33.00	1.26	30.00-36.00	3.94	6	16.33	0.82	14.00-18.00	5.00	6	27.50	1.05	25.00-31.00	3.57
PG: OCD Vitros 250,350,400,500,750,950	10	8.10	0.32	6.00-10.00	3.75	10	20.60	0.52	19.00-23.00	2.38	10	28.30	0.68	25.00-31.00	2.50	10	14.50	0.53	13.00-17.00	3.33	10	23.80	0.63	22.00-26.00	2.50
SG: Select Beckman Coulter Instruments (Chem)	18	10.33	0.48	8.00-12.00	5.00	18	25.22	0.65	23.00-27.00	2.40	18	34.72	1.13	32.00-38.00	3.14	18	17.78	0.81	16.00-20.00	4.44	18	29.17	1.10	26.00-32.00	3.79
SG: Select OCD Vitros Instruments	13	8.00	0.41	6.00-10.00	5.00	13	20.54	0.52	19.00-23.00	2.38	13	28.38	0.77	25.00-31.00	2.86	13	14.46	0.52	12.00-16.00	3.57	13	23.85	0.69	22.00-26.00	2.92
SG: Abaxis Piccolo Instruments	49	8.83	0.53	7.00-11.00	5.56	12	22.67	0.65	21.00-25.00	3.04	12	31.75	1.05	29.00-35.00	3.44	12	15.83	0.58	14.00-18.00	3.75	12	26.50	0.91	25.00-29.00	3.33
SG: Select Abbott Architect Instruments	9	9.89	0.33	8.00-12.00	3.00	9	25.22	0.67	23.00-27.00	2.80	9	34.67	0.87	32.00-38.00	2.57	9	17.44	0.53	15.00-19.00	2.94	9	29.44	1.13	26.00-32.00	3.79
SG: International OCD Vitros Instruments	16	8.00	0.36	6.00-10.00	5.00	16	20.44	0.51	18.00-22.00	2.50	16	28.31	0.79	25.00-31.00	2.86	16	14.44	0.51	12.00-16.00	3.57	16				



SG: Select Gen Chem Instruments / BCP	11	2.01	0.05	1.80-2.20	2.50	11	2.95	0.05	2.60-3.20	1.72	11	3.53	0.08	3.10-3.90	2.29	11	2.48	0.04	2.30-2.80	1.60	11	3.19	0.05	2.90-3.50	1.56
MG: BCG	57	2.20	0.12	2.00-2.40	5.45	78	3.50	0.31	3.10-3.90	8.86	78	4.30	0.34	3.90-4.70	7.91	62	3.00	0.25	2.70-3.30	8.33	78	3.90	0.34	3.50-4.30	8.72
MG: BCG	57	2.20	0.12	2.00-2.40	5.45	78	3.50	0.31	3.10-3.90	8.86	78	4.30	0.34	3.90-4.70	7.91	62	3.00	0.25	2.70-3.30	8.33	78	3.90	0.34	3.50-4.30	8.72
Alkaline Phosphatase	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%
PG: Abaxis Piccolo (Non-Waived)	5	47.80	1.64	34.00-62.00	3.33	5	157.00	7.87	110.00-204.00	5.03	5	218.40	10.50	153.00-283.00	4.82	5	106.20	6.54	74.00-138.00	6.13	5	181.80	10.21	127.00-237.00	5.60
PG: Abaxis Piccolo (Waived)	39	48.15	2.80	34.00-62.00	5.83																				
PG: Abbott Architect	6	55.50	3.15	39.00-73.00	5.54	6	205.83	7.20	144.00-268.00	3.50	6	299.33	12.68	209.00-389.00	4.25	6	132.67	4.41	93.00-173.00	3.31	6	246.33	8.66	172.00-320.00	3.54
PG: Alfa Was. Ace/NexCT/Alera/Axcel / Alfa Wasserman Reagent	24	52.71	3.42	37.00-69.00	6.42	25	189.92	31.03	133.00-247.00	16.32	25	278.44	23.83	195.00-361.00	8.56	24	125.29	8.23	88.00-163.00	6.56	25	226.88	25.80	159.00-295.00	11.37
PG: Clinical Data Envoy 500	5	51.60	4.98	36.00-68.00	9.62	5	180.20	7.76	126.00-234.00	4.33	5	273.40	20.70	191.00-355.00	7.58	5	119.20	5.93	83.00-155.00	4.96	5	223.00	10.84	156.00-290.00	4.84
PG: Siemens Dimension/AR/ES/RxL/Xpand/EXL	10	58.90	2.92	41.00-77.00	4.92	10	221.10	6.66	155.00-287.00	3.03	10	317.80	8.39	223.00-413.00	2.64	10	141.40	5.04	99.00-183.00	3.55	10	262.20	7.33	183.00-341.00	2.79
PG: Siemens Dimension/AR/ES/RxL/Xpand/EXL / IFCC	6	57.00	2.97	40.00-74.00	5.26	6	213.67	7.12	150.00-278.00	3.32	6	310.33	10.15	217.00-403.00	3.29	6	135.50	3.33	95.00-177.00	2.43	6	253.67	7.06	178.00-330.00	2.80
PG: Beckman AU	15	48.13	5.10	34.00-62.00	10.63	15	177.87	15.69	125.00-231.00	8.82	15	257.80	22.48	181.00-335.00	8.72	15	114.93	9.12	81.00-150.00	7.91	15	211.67	17.54	148.00-276.00	8.25
PG: Roche Cobas Integra / Roche Reagent	6	54.33	2.58	38.00-70.00	4.81	6	202.97	11.91	142.00-264.00	5.86	6	293.22	18.36	205.00-381.00	6.28	6	131.05	8.20	92.00-170.00	6.26	6	240.03	13.07	168.00-312.00	5.46
PG: OCD Vitros 250,350,400,500,750,950	9	49.00	4.55	34.00-64.00	9.39	9	149.78	13.76	105.00-195.00	9.20	9	184.00	19.20	129.00-239.00	10.43	9	106.33	8.53	74.00-138.00	8.02	9	167.89	15.60	118.00-218.00	9.29
SG: Select Beckman Coulter Instruments (Chem)	18	48.39	4.79	34.00-62.00	10.00	18	179.17	15.28	125.00-233.00	8.55	18	259.50	21.72	182.00-338.00	8.35	18	115.50	9.09	81.00-151.00	7.84	18	212.61	17.09	149.00-277.00	8.03
SG: Select OCD Vitros Instruments	12	49.50	4.01	35.00-65.00	8.00	12	151.25	12.05	106.00-196.00	8.01	12	186.58	17.15	131.00-243.00	9.20	12	106.75	7.35	75.00-139.00	6.92	12	169.50	13.66	119.00-221.00	8.06
SG: Abaxis Piccolo Instruments	44	48.11	2.68	34.00-62.00	5.63	5	157.00	7.87	110.00-204.00	5.03	5	218.40	10.50	153.00-283.00	4.82	5	106.20	6.54	74.00-138.00	6.13	5	181.80	10.21	127.00-237.00	5.60
SG: Select Abbott Architect Instruments	9	54.56	3.13	39.00-72.00	5.64	9	202.44	9.77	141.00-263.00	4.85	9	293.33	15.54	205.00-381.00	5.29	9	130.00	6.84	91.00-169.00	5.23	9	241.44	12.38	169.00-313.00	5.15
SG: International OCD Vitros Instruments	15	49.20	3.78	34.00-64.00	7.76	15	151.13	10.86	106.00-196.00	7.22	15	186.27	15.48	130.00-242.00	8.33	15	105.93	6.90	74.00-138.00	6.51	15	169.87	12.16	119.00-221.00	7.18
SG: International Abaxis Piccolo Instruments	45	48.29	2.90	34.00-62.00	6.04	6	156.17	7.33	109.00-203.00	4.68	6	217.33	9.75	152.00-282.00	4.52	6	106.33	5.85	74.00-138.00	5.57	6	181.67	9.14	127.00-237.00	5.00
SG: Select Gen Chem Instruments	34	51.06	5.24	36.00-66.00	10.20	34	187.35	17.15	131.00-243.00	9.14	34	272.91	24.82	191.00-355.00	9.08	34	120.94	10.39	85.00-157.00	8.60	34	223.74	19.68	157.00-291.00	8.79
AG: All Group	136	51.24	4.85	36.00-66.00	9.61	98	188.05	27.51	132.00-244.00	14.63	98	269.08	42.49	188.00-350.00	15.80	97	123.14	12.94	86.00-160.00	10.49	98	222.51	31.91	156.00-290.00	14.30
ALT/SGPT	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%
PG: Abaxis Piccolo (Non-Waived)	8	33.25	2.31	26.00-40.00	6.97	8	119.38	2.20	95.00-143.00	1.85	8	170.25	4.03	136.00-204.00	2.35	8	77.00	3.12	62.00-92.00	4.03	8	140.50	2.00	113.00-169.00	1.42
PG: Abaxis Piccolo (Waived)	39	34.59	2.02	28.00-42.00	5.71																				
PG: Abbott Architect	6	32.83	1.33	26.00-40.00	3.94	6	134.17	3.87	107.00-161.00	2.91	6	197.17	4.49	158.00-236.00	2.28	6	83.67	2.50	67.00-101.00	2.98	6	160.83	3.76	129.00-193.00	2.36
PG: Alfa Was. Ace/NexCT/Alera/Axcel / Alfa Wasserman Reagent	25	25.56	1.80	21.00-31.00	6.92	25	106.64	4.12	86.00-128.00	3.83	26	157.12	10.30	126.00-188.00	6.56	24	66.46	2.86	53.00-79.00	4.39	26	127.88	12.46	102.00-154.00	9.77
PG: Clinical Data Envoy 500	5	34.40	1.14	27.00-41.00	3.24	5	134.00	4.64	107.00-161.00	3.43	5	195.20	7.79	156.00-234.00	4.00	5	85.60	2.88	69.00-103.00	3.37	5	162.20	5.80	130.00-194.00	3.58
PG: Siemens Dimension/AR/ES/RxL/Xpand/EXL	7	36.00	2.24	29.00-43.00	6.11	7	135.86	3.67	109.00-163.00	2.72	7	198.71	4.07	159.00-239.00	2.06	7	87.14	4.63	70.00-104.00	5.29	7	161.43	4.61	129.00-193.00	2.86
PG: Beckman AU	15	28.67	1.29	23.00-35.00	4.48	15	116.33	3.72	93.00-139.00	3.19	15	170.27	4.50	136.00-204.00	2.65	15	73.00	2.59	58.00-88.00	3.56	15	139.20	4.60	111.00-167.00	3.31
PG: Roche Cobas Integra / Roche Reagent	6	30.83	1.17	25.00-37.00	3.87	6	127.33	3.27	102.00-152.00	2.60	6	185.50	5.65	149.00-223.00	3.01	6	79.83	1.94	64.00-96.00	2.38	6	151.33	3.98	121.00-181.00	2.65
PG: OCD Vitros 250,350,400,500,750,950	8	51.50	3.25	42.00-62.00	6.35	9	135.11	8.34	108.00-162.00	6.15	9	187.00	9.84	150.00-224.00	5.24	9	92.44	7.94	74.00-110.00	8.59	9	158.11	8.62	126.00-190.00	5.44
SG: Select Beckman Coulter Instruments (Chem)	18	29.78	3.39	24.00-36.00	11.33	18	118.44	6.94	94.00-142.00	5.85	18	172.67	8.16	138.00-208.00	4.74	18	74.61	5.25	60.00-90.00	6.93	18	141.44	7.66	113.00-169.00	5.46
SG: Select OCD Vitros Instruments	10	51.20	3.08	41.00-61.00	6.08	12	134.25	9.77	107.00-161.00	7.31	12	186.08	9.88	149.00-223.00	5.32	12	90.83	9.09	73.00-109.00	10.00	12	156.25	9.95	125.00-187.00	6.35
SG: Abaxis Piccolo Instruments	47	34.36	2.11	27.00-41.00	6.18	8	119.38	2.20	95.00-143.00	1.85	8	170.25	4.03	136.00-204.00	2.35	8	77.00	3.12	62.00-92.00	4.03	8	140.50	2.00	113.00-169.00	1.42
SG: Select Abbott Architect Instruments	9	33.00	1.50	26.00-40.00	4.55	9	134.67	3.74	108.00-162.00	2.74	9	197.44	4.45	158.00-236.00	2.23	9	84.11	2.71	67.00-101.00	3.21	9	161.11	3.41	129.00-193.00	2.11
SG: International OCD Vitros Instruments	12	51.42	2.84	41.00-61.00	5.49	15	135.33	9.04	108.00-162.00	6.67	15	186.80	8.90	150.00-224.00	4.76	15	92.33	8.69	74.00-110.00	9.46	15	157.33	9.16	126.00-188.00	5.86
SG: International Abaxis Piccolo Instruments	48	34.40	2.10	27.00-41.00	6.18	9	119.56	2.13	96.00-144.00	1.75	9	170.00	3.84	136.00-204.00	2.24	9	77.22	2.99	62.00-92.00	3.90					
SG: Select Gen Chem Instruments	33	31.55	3.33	26.00-38.00	10.31	33	126.06	10.37	101.00-151.00	8.25	33	183.94	14.50	147.00-221.00	7.88	33	79.33	6.79	63.00-95.00	8.61	33	150.94	12.43	121.00-181.00	8.21
AG: All Group	128	31.98	4.25	26.00-38.00	13.44	101	123.59	13.13	99.00-149.00	10.56	102	178.86	18.03	143.00-215.00	10.06	100	78.89	9.87	63.00-95.00	12.53	102	146.98	16.20	118.00-176.00	11.02
AST/SGOT	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%
PG: Abaxis Piccolo (Non-Waived)	8	131.50	3.12	106.00-158.00	2.35	8	203.75	4.98	163.00-245.00	2.45	8	250.25	4.59	200.00-300.00	1.84	8	168.13	3.94	134.00-202.00	2.32	8	223.25	4.98	178.00-268.00	2.24
PG: Abaxis Piccolo (Waived)	39	127.18	3.05	102.00-152.00	2.44																				
PG: Abbott Architect	6	130.83	3.25	105.00-157.00	2.52	6	207.17	3.92	166.00-248.00	1.88	6	254.83	5.15	204.00-306.00	2.04	6	169.17	3.66	135.00-203.00	2.19	6	226.33	5.43	181.00-271.00	2.39
PG: Alfa Was. Ace/NexCT/Alera/Axcel / Alfa Wasserman Reagent	25	118.56	3.88	95.00-143.00	3.28	26	185.42	13.97	148.00-222.00	7.57	26	228.85	10.05	183.00-275.00	4.41	26	156.58	17.90	126.00-188.00	11.40	26	204.31	11.25	163.00-245.00	5.54
PG: Clinical Data Envoy 500	5	134.80	5.72	108.00-162.00	4.22	5	209.60	8.50	168.00-252.00	4.05	5	257.20	10.13	206.00-308.00	3.93	5	173.60	7.02	139.00-209.00	4.02	5				



SG: Select Beckman Coulter Instruments (Chem)	17	0.74	0.09	0.30-1.10	12.86	18	2.46	0.22	2.00-3.00	8.80	18	3.34	0.28	2.60-4.00	8.48	18	1.65	0.16	1.30-2.10	9.41	18	2.87	0.24	2.30-3.50	8.28
SG: Select OCD Vitros Instruments	12	0.62	0.12	0.20-1.00	20.00	12	2.39	0.19	1.90-2.90	7.92	11	3.27	0.32	2.60-4.00	9.70	12	1.53	0.14	1.10-1.90	9.33	12	2.76	0.22	2.20-3.40	7.86
SG: Abaxis Piccolo Instruments	42	0.73	0.10	0.30-1.10	14.29	5	2.24	0.22	1.80-2.60	10.00	5	3.02	0.39	2.40-3.60	13.00	5	1.52	0.15	1.10-1.90	10.00	5	2.62	0.29	2.10-3.10	11.15
SG: Select Abbott Architect Instruments	9	0.67	0.05	0.30-1.10	7.14	9	2.24	0.12	1.80-2.60	5.45	9	2.98	0.20	2.40-3.60	6.67	9	1.54	0.11	1.10-1.90	7.33	9	2.58	0.17	2.10-3.10	6.54
SG: International OCD Vitros Instruments	15	0.62	0.11	0.20-1.00	18.33	15	2.40	0.21	1.90-2.90	8.75	14	3.29	0.33	2.60-4.00	10.00	15	1.54	0.14	1.10-1.90	10.00	15	2.77	0.24	2.20-3.40	8.57
SG: International Abaxis Piccolo Instruments	43	0.73	0.10	0.30-1.10	14.29	6	2.18	0.24	1.80-2.60	10.91	6	2.93	0.41	2.30-3.50	14.14	6	1.48	0.16	1.10-1.90	10.67	6	2.55	0.32	2.10-3.10	11.92
SG: Select Gen Chem Instruments	33	0.70	0.11	0.30-1.10	15.71	34	2.34	0.26	1.80-2.80	11.30	34	3.18	0.34	2.60-3.80	10.63	34	1.58	0.18	1.20-2.00	11.25	34	2.73	0.30	2.20-3.20	11.11
AG: All Group	135	0.71	0.17	0.30-1.10	24.29	100	2.33	0.33	1.80-2.80	14.35	100	3.19	0.41	2.60-3.80	12.81	101	1.57	0.23	1.20-2.00	14.38	100	2.74	0.35	2.20-3.20	12.96
Total Protein	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%
PG: Abaxis Piccolo (Non-Waived)	5	3.26	0.09	3.00-3.60	2.73	5	5.66	0.13	5.10-6.30	2.28	5	7.14	0.15	6.40-7.80	2.11	5	4.48	0.08	4.00-5.00	1.78	5	6.30	0.10	5.70-6.90	1.59
PG: Abaxis Piccolo (Waived)	38	3.23	0.06	2.90-3.50	1.88																				
PG: Abbott Architect	6	3.03	0.05	2.70-3.30	1.67	6	5.62	0.12	5.00-6.20	2.14	6	7.28	0.13	6.60-8.00	1.78	6	4.28	0.10	3.90-4.70	2.33	6	6.29	0.11	5.70-6.90	1.75
PG: Alfa Was. Ace/NexCT/Alera/Axcel / Alfa Wasserman Reagent	25	3.16	0.09	2.90-3.50	2.81	26	5.60	0.50	5.00-6.20	8.93	26	7.20	0.33	6.50-7.90	4.58	25	4.46	0.10	4.00-5.00	2.22	26	6.28	0.38	5.70-6.90	6.03
PG: Clinical Data Envoy 500	5	3.20	0.16	2.90-3.50	5.00	5	5.60	0.19	5.00-6.20	3.39	5	7.40	0.22	6.70-8.10	2.97	5	4.46	0.21	4.00-5.00	4.67	5	6.38	0.28	5.80-7.00	4.38
PG: Siemens Dimension/AR/ES/RxL/Xpand/EXL	16	3.21	0.11	2.90-3.50	3.44	16	5.79	0.12	5.20-6.40	2.07	16	7.41	0.12	6.70-8.10	1.62	16	4.49	0.09	4.00-5.00	2.00	16	6.48	0.14	5.80-7.20	2.15
PG: Beckman AU	15	3.04	0.09	2.70-3.30	3.00	15	5.51	0.14	5.00-6.00	2.55	15	7.04	0.18	6.30-7.70	2.57	15	4.28	0.12	3.90-4.70	2.79	15	6.16	0.16	5.60-6.80	2.58
PG: Roche Cobas Integra / Roche Reagent	6	3.10	0.06	2.80-3.40	1.94	6	5.52	0.07	5.00-6.00	1.45	6	7.02	0.13	6.30-7.70	1.86	6	4.35	0.08	4.00-4.80	1.82	6	6.15	0.11	5.60-6.80	1.61
PG: OCD Vitros 250,350,400,500,750,950	9	3.11	0.09	2.80-3.40	2.90	9	5.47	0.20	5.00-6.00	3.64	9	6.96	0.19	6.30-7.70	2.71	9	4.32	0.13	3.90-4.70	3.02	9	6.18	0.22	5.60-6.80	3.55
SG: Select Beckman Coulter Instruments (Chem)	18	3.05	0.09	2.80-3.40	2.90	18	5.51	0.13	5.00-6.00	2.36	18	7.04	0.17	6.30-7.70	2.43	18	4.29	0.12	3.90-4.70	2.79	18	6.15	0.15	5.60-6.80	2.42
SG: Select OCD Vitros Instruments	12	3.09	0.11	2.80-3.40	3.55	12	5.43	0.19	4.90-5.90	3.52	12	6.95	0.19	6.30-7.70	2.71	12	4.31	0.12	3.90-4.70	2.79	12	6.14	0.21	5.50-6.70	3.44
SG: Abaxis Piccolo Instruments	43	3.23	0.06	2.90-3.50	1.88	5	5.66	0.13	5.10-6.30	2.28	5	7.14	0.15	6.40-7.80	2.11	5	4.48	0.08	4.00-5.00	1.78	5	6.30	0.10	5.70-6.90	1.59
SG: Select Abbott Architect Instruments	9	3.02	0.04	2.70-3.30	1.33	9	5.61	0.09	5.00-6.20	1.61	9	7.27	0.11	6.60-8.00	1.51	9	4.28	0.08	3.90-4.70	1.86	9	6.28	0.09	5.70-6.90	1.43
SG: International OCD Vitros Instruments	15	3.10	0.13	2.80-3.40	4.19	15	5.45	0.20	5.00-6.00	3.64	15	6.97	0.21	6.30-7.70	3.00	15	4.32	0.14	3.90-4.70	3.26	15	6.16	0.23	5.60-6.80	3.71
SG: International Abaxis Piccolo Instruments	44	3.23	0.06	2.90-3.50	1.88	6	5.67	0.12	5.10-6.30	2.11	6	7.17	0.15	6.50-7.90	2.08	6	4.48	0.07	4.00-5.00	1.78	6	6.32	0.10	5.70-6.90	1.59
SG: Select Gen Chem Instruments	34	3.05	0.13	2.80-3.40	4.19	34	5.55	0.14	5.00-6.00	2.55	34	7.15	0.22	6.40-7.80	3.10	34	4.31	0.14	3.90-4.70	3.26	34	6.21	0.18	5.60-6.80	2.90
AG: All Group	136	3.15	0.12	2.90-3.50	3.75	99	5.59	0.30	5.00-6.20	5.36	99	7.17	0.26	6.50-7.90	3.61	98	4.39	0.14	4.00-4.80	3.18	99	6.27	0.26	5.70-6.90	4.13
Phosphorous	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%
PG: Abaxis Piccolo (Non-Waived)	5	3.20	0.12	2.80-3.60	3.75	5	5.44	0.21	4.80-6.00	3.89	5	6.54	0.11	6.20-6.80	1.69	5	4.38	0.16	3.90-4.90	3.64	5	5.88	0.43	4.60-7.20	7.29
PG: Alfa Was. Ace/NexCT/Alera/Axcel / Alfa Wasserman Reagent	6	3.00	0.34	2.00-4.00	11.00	6	5.03	0.53	3.40-6.60	10.60	6	6.32	0.71	4.20-8.40	11.27	6	4.02	0.40	2.80-5.20	10.00	6	5.40	0.67	3.40-7.40	12.41
PG: Siemens Dimension/AR/ES/RxL/Xpand/EXL	5	2.68	0.08	2.50-2.90	2.96	5	4.66	0.15	4.30-5.20	3.19	5	5.84	0.17	5.30-6.30	2.93	5	3.68	0.13	3.30-4.10	3.51	5	5.10	0.14	4.70-5.50	2.75
PG: Beckman AU	9	2.67	0.34	1.70-3.80	12.96	8	4.55	0.30	3.60-5.40	6.67	9	5.90	0.54	4.30-7.50	9.15	9	3.66	0.33	2.70-4.70	8.92	9	5.10	0.43	3.80-6.40	8.43
SG: Select Beckman Coulter Instruments (Chem)	11	2.65	0.32	1.60-3.60	12.31	10	4.56	0.28	3.80-5.40	6.09	11	5.88	0.50	4.40-7.40	8.47	11	3.63	0.31	2.70-4.50	8.61	11	5.10	0.40	3.90-6.30	7.84
SG: Abaxis Piccolo Instruments	6	3.18	0.12	2.80-3.60	3.75	5	5.44	0.21	4.80-6.00	3.89	5	6.54	0.11	6.20-6.80	1.69	5	4.38	0.16	3.90-4.90	3.64	5	5.88	0.43	4.60-7.20	7.29
PG: Intl OCD Vitros 250,350,400,500,750,950	5	2.67	0.06	2.50-2.90	2.22	5	4.53	0.23	3.80-5.20	5.11	5	5.83	0.40	4.60-7.00	6.90	5	3.63	0.15	3.20-4.00	4.17	5	5.17	0.31	4.30-6.10	5.96
SG: International Abaxis Piccolo Instruments	7	3.26	0.22	2.60-4.00	6.67	6	5.53	0.29	4.60-6.40	5.27	6	6.68	0.37	5.60-7.80	5.52	6	4.47	0.26	3.70-5.30	5.78	6	5.98	0.46	4.60-7.40	7.67
SG: Select Gen Chem Instruments	17	2.67	0.32	1.70-3.70	11.85	16	4.64	0.37	3.50-5.70	8.04	17	5.94	0.51	4.40-7.40	8.47	17	3.65	0.30	2.70-4.50	8.33	17	5.13	0.41	3.90-6.30	8.04
AG: All Group	41	2.79	0.32	1.80-3.80	11.43	39	4.79	0.46	3.40-6.20	9.58	40	6.04	0.52	4.40-7.60	8.67	40	3.80	0.39	2.70-4.90	10.00	40	5.26	0.51	3.80-6.80	9.62
Cholesterol Total	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%
PG: Alfa Was. Ace/NexCT/Alera/Axcel / Alfa Wasserman Reagent	23	77.96	5.12	70.00-86.00	6.54	22	139.64	3.02	126.00-154.00	2.14	24	177.00	9.71	159.00-195.00	5.48	23	108.35	2.72	97.00-119.00	2.50	24	153.13	10.56	138.00-168.00	6.93
PG: Siemens Dimension/AR/ES/RxL/Xpand/EXL	14	72.00	3.31	65.00-79.00	4.58	14	134.07	4.55	121.00-147.00	3.36	14	172.29	5.63	155.00-189.00	3.26	14	103.14	3.65	93.00-113.00	3.59	14	149.86	4.45	135.00-165.00	3.00
PG: Beckman AU	15	72.00	2.75	65.00-79.00	3.89	15	136.00	10.42	122.00-150.00	7.65	15	172.93	5.27	156.00-190.00	3.06	15	104.67	7.41	95.00-116.00	7.05	15	150.73	4.71	136.00-166.00	3.11
PG: Roche Cobas Integra / Roche Reagent	6	73.17	2.23	66.00-80.00	3.01	6	133.17	4.88	120.00-146.00	3.68	6	173.00	6.16	156.00-190.00	3.58	6	103.00	3.10	93.00-113.00	3.01	6	149.50	5.01	135.00-165.00	3.33
PG: OCD Vitros 250,350,400,500,750,950	9	59.11	2.93	53.00-65.00	4.92	9	132.44	3.40	119.00-145.00	2.58	9	174.22	3.56	157.00-191.00	2.07	9	100.22	3.07	90.00-110.00	3.10	9	152.56	4.88	138.00-168.00	3.20
SG: Select Beckman Coulter Instruments (Chem)	18	72.22	2.69	65.00-79.00	3.75	18	136.11	9.49	122.00-150.00	6.99	18	173.61	5.11	157.00-191.00	2.93	18	104.89	6.81	95.00-116.00	6.48	18	151.22	4.45	136.00-166.00	2.98
SG: Select OCD Vitros Instruments	12	58.58	2.78	53.00-65.00	4.75	12	132.58	2.94	120.00-146.00	2.18	12	175.17	3.51	158.00-193.00	2.00	12	99.83	2.73	90.00-110.00	2.70	12	151.83	4.49	137.00-167.00	2.96
SG: Abaxis Piccolo Instruments	8	75.00	2.56	68.00-83.00	3.47																				
SG: Select Abbott Architect Instruments	6	75.50	1.05	68.00-84.00	1.32	6	137.67	1.03	124.00-152.00	0.72	6	177.50	1.23	160.00-196.00	0.67	6	106.67	0.82	96.00-118.00	0.75	6	153.67	1.03	139.00-169.00	0.65
SG: International OCD Vitros Instruments	15	58.87	2.77	53.00-65.00	4.75	15	132.73	2.87	120.00-146.00	2.18	15	174.73	3.39	158.00-193.00	1.94	15	100.07	2.66	90.00-110.00	2.70	15	151.87	4.00	137.00-167.00	2.63
SG: International Abaxis Piccolo Instruments	9	75.67																							



PG: Alfa Was. Ace/NexCT/Alera/Axcel / Alfa Wasserman Reagent	11	21.38	2.15	14.00-28.00	10.48	10	47.94	3.99	36.00-60.00	8.33	12	63.03	9.38	35.00-91.00	14.92	12	31.90	6.06	14.00-50.00	19.06	13	49.97	9.94	20.00-80.00	19.80
PG: Siemens Dimension/AR/ES/RxL/Xpand/EXL / Direct Method	6	16.83	2.48	10.00-25.00	14.71	6	33.17	2.32	26.00-40.00	6.97	6	44.33	3.20	34.00-54.00	7.27	6	25.50	2.43	19.00-33.00	9.23	6	38.00	2.37	31.00-45.00	6.32
PG: Beckman AU / Calc	9	15.33	2.87	6.00-24.00	19.33	6	35.67	2.25	29.00-43.00	6.39	7	50.29	2.21	43.00-57.00	4.40	6	26.33	2.58	18.00-34.00	10.00	7	42.14	3.29	32.00-52.00	7.86
PG: Beckman AU / Direct Method	5	10.50	1.00	8.00-14.00	9.09	5	20.00	1.41	16.00-24.00	7.00	5	26.00	0.82	24.00-28.00	3.08	5	15.75	0.50	15.00-18.00	3.13	5	22.75	0.50	22.00-25.00	2.17
PG: OCD Vitros 250,350,400,500,750,950 / Calc	5	1.00	3.46	0.00-12.00	350.00	6	16.50	2.43	10.00-24.00	14.12	6	18.17	6.18	0.00-37.00	34.44	6	10.83	1.60	6.00-16.00	14.55	6	17.67	4.46	5.00-32.00	25.00
PG: OCD Vitros 250,350,400,500,750,950 / Direct Method	5	-13.00	0.00	0.00--13.00	0.00	5	2.00	0.00	2.00-2.00	0.00	5	4.00	0.00	4.00-4.00	0.00					5	7.00	0.00	7.00-7.00	0.00	
SG: Select Beckman Coulter Instruments (Chem) / Calc	10	15.50	2.76	8.00-24.00	17.50	7	35.43	2.15	29.00-41.00	6.00	8	51.00	2.88	42.00-60.00	5.69	7	26.64	2.50	20.00-35.00	9.26	8	42.42	3.15	33.00-51.00	7.38
SG: Select Abbott Architect Instruments / Direct Method	5	13.00	1.00	10.00-16.00	7.69	5	26.00	2.00	20.00-32.00	7.69	5	33.33	2.52	26.00-41.00	7.58	5	19.67	1.53	16.00-25.00	7.50	5	29.00	2.00	23.00-35.00	6.90
PG: OCD Vitros 5600 / Calc	5	6.00	0.00	6.00-6.00	0.00	5	20.50	2.12	15.00-27.00	10.00	5	26.00	4.24	13.00-39.00	16.15	5	11.50	3.54	2.00-23.00	29.17	5	22.50	2.12	17.00-29.00	9.13
SG: Select Gen Chem Instruments / Calc	14	16.00	4.04	4.00-28.00	25.00	11	35.73	6.59	16.00-56.00	18.33	12	51.75	7.61	29.00-75.00	14.62	10	25.35	3.61	14.00-36.00	14.40	12	41.95	6.50	23.00-62.00	15.48
MG: Calc	34	18.00	4.70	4.00-32.00	26.11	33	36.00	11.20	2.00-70.00	31.11	31	53.00	14.30	10.00-96.00	26.98	39	23.00	9.80	0.00-52.00	42.61	36	40.00	14.20	0.00-83.00	35.50
MG: Direct Method	18	15.00	3.80	4.00-26.00	25.33	18	29.00	7.10	8.00-50.00	24.48	18	39.00	9.90	9.00-69.00	25.38	18	22.00	5.30	6.00-38.00	24.09	18	33.00	8.30	8.00-58.00	25.15
Triglycerides	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%
PG: Alfa Was. Ace/NexCT/Alera/Axcel / Alfa Wasserman Reagent	23	120.35	5.12	90.00-150.00	4.25	24	181.00	15.56	136.00-226.00	8.62	24	221.63	14.22	167.00-278.00	6.40	23	152.35	7.00	114.00-190.00	4.61	24	197.46	14.11	148.00-246.00	7.16
PG: Clinical Data Envoy 500	5	122.00	4.58	92.00-153.00	3.77	5	187.80	7.79	141.00-235.00	4.15	5	229.80	13.48	173.00-288.00	5.87	5	151.40	10.11	113.00-189.00	6.69	5	204.80	10.11	154.00-256.00	4.93
PG: Siemens Dimension/AR/ES/RxL/Xpand/EXL	14	107.71	6.79	81.00-135.00	6.30	14	170.21	6.48	128.00-213.00	3.82	14	211.07	7.89	158.00-264.00	3.74	14	138.29	5.66	104.00-173.00	4.13	14	187.07	8.32	140.00-234.00	4.44
PG: Beckman AU	14	117.07	5.12	88.00-146.00	4.36	14	181.21	7.89	136.00-226.00	4.36	14	222.36	11.19	167.00-278.00	5.05	14	148.36	5.81	111.00-185.00	3.92	14	197.36	9.35	148.00-246.00	4.72
PG: Roche Cobas Integra / Roche Reagent	6	117.17	5.64	88.00-146.00	4.79	6	178.33	9.65	134.00-223.00	5.39	6	214.50	12.34	161.00-269.00	5.72	6	147.67	6.62	111.00-185.00	4.46	6	193.00	10.32	145.00-241.00	5.34
PG: OCD Vitros 250,350,400,500,750,950	9	125.78	8.33	95.00-158.00	6.59	9	201.22	13.38	151.00-251.00	6.67	9	253.44	19.02	190.00-316.00	7.51	9	164.44	11.24	123.00-205.00	6.83	9	222.89	15.55	167.00-279.00	7.00
SG: Select Beckman Coulter Instruments (Chem)	17	117.47	5.09	88.00-146.00	4.36	17	182.53	8.14	137.00-229.00	4.43	17	224.65	11.82	169.00-281.00	5.24	17	149.47	6.26	112.00-186.00	4.23	17	199.00	9.66	149.00-249.00	4.87
SG: Select OCD Vitros Instruments	12	123.42	9.40	92.00-154.00	7.64	12	197.08	14.43	148.00-246.00	7.31	12	249.50	20.24	188.00-313.00	8.08	12	161.17	12.57	121.00-201.00	7.83	12	219.00	17.17	164.00-274.00	7.85
SG: Select Abbott Architect Instruments	6	111.00	3.03	83.00-139.00	2.70	6	175.00	6.54	131.00-219.00	3.71	6	217.50	10.69	164.00-273.00	4.91	6	143.50	4.97	108.00-180.00	3.47	6	193.17	7.99	145.00-241.00	4.15
SG: International OCD Vitros Instruments	15	122.93	8.46	92.00-154.00	6.91	15	195.93	13.28	147.00-245.00	6.79	15	247.73	18.77	186.00-310.00	7.58	15	160.00	11.61	120.00-200.00	7.25	15	218.73	15.49	164.00-274.00	7.08
SG: International Abaxis Piccolo Instruments	9	126.44	7.11	95.00-158.00	5.63																				
SG: Select Gen Chem Instruments	30	116.70	5.59	88.00-146.00	4.79	30	181.23	8.65	136.00-226.00	4.81	30	223.47	11.95	167.00-279.00	5.34	30	148.13	7.04	111.00-185.00	4.73	30	198.23	9.79	149.00-248.00	4.95
SG: Select Cholestech Instruments	5	114.00	4.36	86.00-143.00	3.86	5	171.33	7.23	128.00-214.00	4.21	5	218.00	15.40	164.00-273.00	7.06	5	148.33	7.23	111.00-185.00	4.86	5	195.33	8.14	146.00-244.00	4.15
AG: All Group	96	117.80	8.05	89.00-148.00	6.86	93	181.84	14.05	137.00-228.00	7.69	93	224.75	17.88	169.00-281.00	7.96	92	149.97	10.31	113.00-188.00	6.87	93	199.45	15.21	149.00-249.00	7.64
Amylase	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%
PG: Abaxis Piccolo (Waived)	30	27.03	1.52	19.00-35.00	5.56																				
PG: Alfa Was. Ace/NexCT/Alera/Axcel / Alfa Wasserman Reagent	6	28.17	1.84	20.00-36.00	6.43	6	96.67	3.33	68.00-126.00	3.40	6	140.67	2.88	99.00-183.00	2.06	6	62.83	1.94	44.00-82.00	3.02	6	115.17	3.19	81.00-150.00	2.78
PG: Siemens Dimension/AR/ES/RxL/Xpand/EXL	8	29.25	0.71	20.00-38.00	2.41	8	104.38	1.92	73.00-135.00	1.83	8	150.63	2.56	106.00-196.00	1.72	8	66.88	1.13	47.00-87.00	1.64	8	123.75	2.19	87.00-161.00	1.77
PG: Beckman AU	8	21.88	0.99	15.00-29.00	4.55	8	76.25	3.24	53.00-99.00	4.21	8	109.38	4.63	76.00-142.00	4.22	8	49.25	2.38	34.00-64.00	4.90	8	89.88	4.52	63.00-117.00	5.00
PG: OCD Vitros 250,350,400,500,750,950	5	30.00	0.00	21.00-39.00	0.00	5	57.00	1.41	40.00-74.00	2.46	5	89.50	0.71	63.00-117.00	0.78	5	37.50	0.71	27.00-49.00	1.84	5	67.00	1.41	47.00-87.00	2.09
SG: Select Beckman Coulter Instruments (Chem)	10	22.60	2.76	16.00-30.00	12.17	10	78.20	7.21	55.00-101.00	9.23	10	112.20	10.96	78.00-146.00	9.82	10	50.60	4.84	36.00-66.00	9.41	10	92.30	8.92	64.00-120.00	9.67
SG: Abaxis Piccolo Instruments	32	27.09	1.49	19.00-35.00	5.56																				
SG: International OCD Vitros Instruments	5	30.00	0.00	21.00-39.00	0.00	5	57.75	1.26	41.00-75.00	2.24	5	87.75	3.86	62.00-114.00	4.43	5	37.50	1.29	27.00-49.00	3.42	5	69.00	2.58	48.00-90.00	3.77
SG: International Abaxis Piccolo Instruments	33	27.03	1.51	19.00-35.00	5.56	5	90.33	1.53	63.00-117.00	1.67	5	131.00	2.00	92.00-170.00	1.53	5	60.00	1.73	42.00-78.00	2.83	5	107.33	3.21	75.00-139.00	2.99
SG: Select Gen Chem Instruments	15	24.73	3.96	18.00-33.00	16.00	15	85.73	13.07	60.00-112.00	15.23	15	123.47	19.78	86.00-160.00	16.10	15	55.93	9.18	39.00-73.00	16.43	15	101.93	16.44	71.00-133.00	16.08
AG: All Group	67	27.09	2.69	19.00-35.00	10.00	29	95.00	10.56	67.00-124.00	11.16	28	138.46	14.69	97.00-179.00	10.65	32	60.41	7.79	42.00-78.00	13.00	32	110.00	15.75	77.00-143.00	14.36
Creatine Kinase	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%
PG: Abaxis Piccolo (Waived)	22	65.00	8.58	46.00-85.00	13.23																				
PG: Alfa Was. Ace/NexCT/Alera/Axcel / Alfa Wasserman Reagent	6	54.00	2.83	38.00-70.00	5.19	6	170.33	7.76	119.00-221.00	4.59	7	228.00	24.06	160.00-296.00	10.57	6	111.17	4.58	78.00-144.00	4.14	6	200.33	9.22	140.00-260.00	4.60
PG: Beckman AU	7	45.29	3.15	32.00-59.00	6.89	7	163.86	9.84	115.00-213.00	5.98	7	235.00	13.66	165.00-306.00	5.83	7	104.71	6.29	74.00-137.00	6.00	7	194.86	10.12	137.00-254.00	5.18
SG: Select Beckman Coulter Instruments (Chem)	8	46.50	4.50	33.00-61.00	9.57	8	167.00	12.73	117.00-217.00	7.60	8	238.75	16.51	167.00-311.00	6.90	8	107.13	8.97	75.00-139.00	8.41	8	198.25	13.41	139.00-257.00	6.77
SG: Select OCD Vitros Instruments	5	40.75	3.95	29.00-53.00	9.51	5	149.00	15.56	104.00-194.00	10.47	5	206.25	18.17	144.00-268.00	8.83	5	99.75	9.36	70.00-130.00	9.40	5	175.75	16.52	123.00-229.00	9.38
SG: Abaxis Piccolo Instruments	25	65.12	8.05	46.00-85.00	12.46																				
SG: International OCD Vitros Instruments	7	41.57	4.61	29.00-55.00	10.95	7	147.86	11.45	104.00-192.00	7.77	7	203.00	14.66	142.00-264.00	7.24	7	99.14	7.08	69.00-129.00	7.17	7	173.14	12.98	121.00-	



Bilirubin Direct	Labs	Mean	SD	Range	CV %	Labs	Mean	SD	Range	CV %	Labs	Mean	SD	Range	CV %	Labs	Mean	SD	Range	CV %	Labs	Mean	SD	Range	CV %
PG: Abbott Architect	5	0.37	0.04	0.24-0.50	11.62	5	1.03	0.05	0.89-1.17	4.47	5	1.25	0.09	1.00-1.50	6.80	5	0.74	0.05	0.60-0.88	6.22	5	1.16	0.08	0.91-1.41	7.24
PG: Alfa Was. Ace/NexCT/Alera/Axcel / Alfa Wasserman Reagent	15	0.33	0.06	0.15-0.51	17.88	17	0.82	0.16	0.35-1.29	19.02	17	1.04	0.17	0.53-1.55	16.35	16	0.61	0.10	0.32-0.90	15.74	17	0.94	0.15	0.48-1.40	16.38
PG: Siemens Dimension/AR/ES/RxL/Xpand/EXL	5	0.20	0.01	0.17-0.23	4.50	5	0.50	0.07	0.29-0.71	14.20	5	0.60	0.07	0.39-0.81	11.83	5	0.33	0.04	0.20-0.47	13.64	5	0.55	0.05	0.40-0.70	9.09
PG: Siemens Dimension/AR/ES/RxL/Xpand/EXL / IFCC	5	0.20	0.01	0.17-0.23	4.50	5	0.56	0.05	0.41-0.71	8.93	5	0.72	0.04	0.59-0.86	6.25	5	0.41	0.01	0.37-0.45	3.17	5	0.66	0.06	0.49-0.83	8.79
PG: Beckman AU	14	0.30	0.06	0.12-0.48	19.67	14	0.77	0.10	0.46-1.08	13.51	14	0.92	0.14	0.51-1.33	14.89	14	0.58	0.09	0.32-0.84	14.83	14	0.86	0.11	0.53-1.19	12.79
PG: Roche Cobas Integra / Roche Reagent	5	0.17	0.04	0.04-0.30	25.29	5	0.54	0.06	0.38-0.71	10.19	5	0.70	0.11	0.38-1.02	15.14	5	0.35	0.06	0.18-0.52	16.57	5	0.64	0.07	0.44-0.84	10.47
SG: Select Beckman Coulter Instruments (Chem)	15	0.29	0.06	0.12-0.46	19.66	15	0.76	0.11	0.43-1.09	14.34	15	0.91	0.14	0.49-1.33	15.27	15	0.57	0.09	0.29-0.85	16.32	15	0.85	0.12	0.50-1.20	13.65
SG: Select OCD Vitros Instruments	7	0.39	0.14	0.00-0.80	34.62	7	0.97	0.18	0.43-1.51	18.56	7	1.33	0.18	0.79-1.87	13.53	7	0.64	0.15	0.19-1.09	23.59	7	1.09	0.16	0.62-1.56	14.40
SG: Select Abbott Architect Instruments	8	0.36	0.05	0.22-0.50	13.33	8	1.02	0.07	0.82-1.22	6.47	8	1.23	0.13	0.85-1.61	10.41	8	0.73	0.07	0.53-0.93	9.32	8	1.15	0.10	0.84-1.46	9.04
SG: International OCD Vitros Instruments	10	0.39	0.13	0.00-0.78	33.59	9	0.93	0.19	0.36-1.50	20.32	9	1.30	0.19	0.72-1.88	14.92	9	0.63	0.15	0.18-1.08	23.97	9	1.05	0.17	0.53-1.57	16.48
SG: Select Gen Chem Instruments	26	0.31	0.06	0.12-0.50	20.00	26	0.84	0.15	0.38-1.30	18.10	26	1.02	0.20	0.44-1.61	19.12	26	0.62	0.11	0.29-0.95	17.58	26	0.95	0.17	0.43-1.47	18.32
AG: All Group	62	0.28	0.08	0.05-0.51	27.86	66	0.77	0.20	0.18-1.36	25.45	65	0.95	0.24	0.23-1.67	25.37	60	0.53	0.14	0.13-0.94	25.47	67	0.87	0.22	0.22-1.52	24.83
Lactic Acid	Labs	Mean	SD	Range	CV %	Labs	Mean	SD	Range	CV %	Labs	Mean	SD	Range	CV %	Labs	Mean	SD	Range	CV %	Labs	Mean	SD	Range	CV %
PG: Abaxis Piccolo (Non-Waived)	5	0.82	0.04	0.70-0.90	5.00	5	3.68	0.08	3.50-3.90	2.16	5	5.60	0.19	5.00-6.20	3.39	5	2.37	0.04	2.30-2.50	1.67	5	4.40	0.07	4.20-4.60	1.59
AG: All Group	6	0.83	0.05	0.70-1.00	6.25	6	3.70	0.09	3.40-4.00	2.43	6	5.55	0.21	5.00-6.20	3.75	6	2.43	0.14	2.00-2.80	5.83	6	4.42	0.07	4.20-4.60	1.82
Magnesium	Labs	Mean	SD	Range	CV %	Labs	Mean	SD	Range	CV %	Labs	Mean	SD	Range	CV %	Labs	Mean	SD	Range	CV %	Labs	Mean	SD	Range	CV %
PG: Abaxis Piccolo (Non-Waived)	6	1.57	0.05	1.20-2.00	3.13	6	2.87	0.05	2.20-3.60	1.72	6	3.75	0.06	2.80-4.80	1.32	6	2.22	0.04	1.70-2.80	1.82	6	3.17	0.08	2.40-4.00	2.50
PG: Alfa Was. Ace/NexCT/Alera/Axcel / Alfa Wasserman Reagent	6	1.57	0.10	1.20-2.00	6.25	6	2.77	0.15	2.10-3.50	5.36	6	3.57	0.14	2.70-4.50	3.89	6	2.17	0.05	1.70-2.80	2.27	6	3.13	0.15	2.30-3.90	4.84
PG: Siemens Dimension/AR/ES/RxL/Xpand/EXL	7	1.53	0.05	1.10-1.90	3.33	7	2.94	0.08	2.20-3.60	2.76	7	3.73	0.05	2.80-4.60	1.35	7	2.24	0.08	1.70-2.80	3.64	7	3.27	0.08	2.50-4.10	2.42
PG: Beckman AU	6	1.57	0.05	1.20-2.00	3.13	6	2.90	0.06	2.20-3.60	2.07	6	3.72	0.10	2.80-4.60	2.70	6	2.27	0.05	1.70-2.90	2.17	6	3.22	0.07	2.40-4.00	2.50
SG: Select Beckman Coulter Instruments (Chem)	8	1.55	0.05	1.20-2.00	3.13	8	2.84	0.13	2.10-3.50	4.64	8	3.64	0.17	2.70-4.50	4.72	8	2.23	0.09	1.70-2.80	4.09	8	3.16	0.13	2.40-4.00	4.06
SG: Abaxis Piccolo Instruments	6	1.57	0.05	1.20-2.00	3.13	6	2.87	0.05	2.20-3.60	1.72	6	3.75	0.06	2.80-4.80	1.32	6	2.22	0.04	1.70-2.80	1.82	6	3.17	0.08	2.40-4.00	2.50
PG: Intl OCD Vitros 250,350,400,500,750,950	5	1.55	0.07	1.20-2.00	4.38	5	2.95	0.07	2.30-3.80	2.33	5	3.75	0.07	2.80-4.80	1.84	5	2.30	0.00	1.70-2.90	0.00	5	3.30	0.00	2.50-4.10	0.00
SG: International Abaxis Piccolo Instruments	7	1.57	0.05	1.20-2.00	3.13	7	2.86	0.05	2.20-3.60	1.72	7	3.73	0.08	2.80-4.60	2.16	7	2.21	0.04	1.70-2.80	1.82	7	3.17	0.08	2.40-4.00	2.50
SG: Select Gen Chem Instruments	15	1.49	0.12	1.10-1.90	8.00	15	2.85	0.23	2.10-3.50	8.21	15	3.63	0.23	2.70-4.50	6.39	15	2.17	0.20	1.70-2.80	9.55	15	3.11	0.27	2.30-3.90	8.71
AG: All Group	39	1.54	0.10	1.10-1.90	6.67	39	2.86	0.17	2.20-3.60	5.86	39	3.66	0.17	2.80-4.60	4.59	39	2.20	0.14	1.70-2.80	6.36	39	3.16	0.19	2.40-4.00	5.94
Uric Acid	Labs	Mean	SD	Range	CV %	Labs	Mean	SD	Range	CV %	Labs	Mean	SD	Range	CV %	Labs	Mean	SD	Range	CV %	Labs	Mean	SD	Range	CV %
PG: Abaxis Piccolo (Waived)	7	1.85	0.16	1.60-2.20	8.42																				
PG: Alfa Was. Ace/NexCT/Alera/Axcel / Alfa Wasserman Reagent	12	3.52	0.24	2.90-4.10	6.86	12	6.30	0.30	5.20-7.40	4.76	12	8.34	0.24	6.90-9.70	2.89	12	4.73	0.18	3.90-5.50	3.83	12	7.06	0.19	5.90-8.30	2.68
PG: Siemens Dimension/AR/ES/RxL/Xpand/EXL	12	2.60	0.10	2.20-3.00	3.85	12	5.88	0.24	4.90-6.90	4.07	12	7.86	0.34	6.60-9.20	4.30	12	4.22	0.19	3.50-4.90	4.52	12	6.67	0.20	5.60-7.80	2.99
PG: Beckman AU	9	2.76	0.07	2.30-3.30	2.50	9	6.22	0.13	5.10-7.30	2.10	9	8.34	0.12	6.90-9.70	1.45	9	4.47	0.09	3.70-5.30	2.00	9	7.11	0.16	5.90-8.30	2.25
PG: Roche Cobas Integra / Roche Reagent	5	2.60	0.08	2.20-3.00	3.08	5	6.13	0.22	5.10-7.10	3.61	5	8.25	0.24	6.90-9.70	2.89	5	4.35	0.17	3.70-5.10	3.86	5	7.03	0.30	5.80-8.20	4.29
PG: OCD Vitros 250,350,400,500,750,950	7	2.50	0.12	2.10-2.90	4.80	7	5.91	0.12	4.90-6.90	2.03	7	7.91	0.16	6.60-9.20	2.03	7	4.29	0.15	3.60-5.00	3.49	7	6.76	0.10	5.60-8.00	1.47
SG: Select Beckman Coulter Instruments (Chem)	12	2.73	0.09	2.20-3.20	3.33	12	6.16	0.18	5.10-7.30	2.90	12	8.26	0.24	6.90-9.70	2.89	12	4.44	0.10	3.70-5.10	2.27	12	7.03	0.21	5.80-8.20	3.00
SG: Select OCD Vitros Instruments	9	2.51	0.11	2.10-2.90	4.40	9	5.96	0.14	5.00-7.00	2.33	9	7.93	0.15	6.60-9.20	1.90	9	4.30	0.13	3.60-5.00	3.02					
SG: Select Abbott Architect Instruments	5	2.82	0.08	2.30-3.30	2.86	5	6.60	0.10	5.50-7.70	1.52	5	9.02	0.13	7.50-10.50	1.44	5	4.72	0.08	3.90-5.50	1.70	5	7.60	0.10	6.30-8.90	1.32
SG: International OCD Vitros Instruments	12	2.50	0.10	2.10-2.90	4.00	12	5.92	0.14	4.90-6.90	2.37	12	7.90	0.14	6.60-9.20	1.77	12	4.29	0.12	3.60-5.00	2.79	12	6.78	0.12	5.60-8.00	1.76
SG: Select Gen Chem Instruments	22	2.71	0.14	2.20-3.20	5.19	22	6.23	0.28	5.10-7.30	4.52	22	8.41	0.42	7.00-9.80	5.00	22	4.48	0.17	3.70-5.30	3.78	22	7.15	0.32	5.90-8.30	4.51
AG: All Group	68	2.71	0.47	2.20-3.20	17.41	61	6.13	0.30	5.10-7.10	4.92	61	8.20	0.39	6.80-9.60	4.76	61	4.44	0.25	3.70-5.10	5.45	61	6.97	0.30	5.80-8.20	4.29
Iron	Labs	Mean	SD	Range	CV %	Labs	Mean	SD	Range	CV %	Labs	Mean	SD	Range	CV %	Labs	Mean	SD	Range	CV %	Labs	Mean	SD	Range	CV %
PG: Siemens Dimension/AR/ES/RxL/Xpand/EXL	6	47.17	0.75	38.00-56.00	1.70	6	118.67	1.03	95.00-143.00	0.84	6	163.83	0.98	131.00-197.00	0.61	6	83.50	0.55	67.00-101.00	0.60	6	137.50	0.55	110.00-166.00	0.36
PG: Beckman AU	6	51.67	1.97	42.00-62.00	3.85	6	131.67	6.15	106.00-158.00	4.70	6	183.67	8.34	147.00-221.00	4.51	6	94.17	4.49	75.00-113.00	4.79	6	152.67	6.12	122.00-184.00	3.99
SG: Select Beckman Coulter Instruments (Chem)	8	51.13	2.36	41.00-61.00	4.71	8	131.00	5.98	105.00-157.00	4.58	8	182.63	8.09	146.00-220.00	4.43	8	92.75	4.98	74.00-112.00	5.38	8	152.13	6.03	122.00-182.00	3.95
SG: Select Abbott Architect Instruments	6	50.67	1.37	41.00-61.00	2.75	6	124.67	2.07	100.00-150.00	1.68	6	171.17	3.76	137.00-205.00	2.22	6	87.50	1.64	70.00-106.00	1.82	6	144.50	3.27	116.00-174.00	2.28
PG: Intl OCD Vitros 250,350,400,500,750,950	5	31.00	0.00	25.00-37.00	0.00	5	130.00	18.25	104.00-156.00	14.00	5	189.00	14.80	151.00-227.00	7.83	5	81.33	3.21	65.00-97.00	3.95	5	152.67	14.22	122.00-184.00	9.28
SG: Select Gen Chem Instruments	18	49.78	3.10	40.00-60.00	6.20	18	125.28	8.23	100.00-150.00	6.56	18	173.83	11.42	139.00-209.00	6.55	18	88.28	6.23	70.00-106.00	7.05	18	145.61	8.96	117.00-175.00	6.16
AG: All Group	30	48.37	3.75	38.00-58.00	7.92	30	124.10	8.58	99.00-149.00	6.94	30	171.93	11.66	138.00-206.00	6.80	30	86.13	5.93	69.00-103.00	6.86	30	143.50	8.65	115.00-173.00	5.97
TIBC	Labs	Mean	SD	Range	CV %	Labs	Mean	SD	Range	CV %	Labs	Mean	SD	Range	CV %	Labs	Mean	SD	Range	CV %	Labs	Mean	SD	Range	CV %
PG: Siemens Dimension/AR/ES/RxL/Xpand/EXL	5	95.60																							



PG: Abbott Cell-Dyn Ruby	7	53.03	0.77	50.70-55.30	1.45	7	78.77	0.33	77.80-79.80	0.42	7	67.93	0.75	65.70-70.20	1.10	7	53.29	0.80	50.90-55.70	1.50	7	78.57	0.43	77.30-79.90	0.55
Eosinophil	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%
SG: Select Abbott Cell Dyn Instruments (5-part Diff)	8	2.66	0.37	1.60-3.80	13.70	8	4.02	0.28	3.20-4.80	7.00	8	4.48	0.26	3.80-5.30	5.56	8	2.56	0.43	1.30-3.90	16.54	8	4.02	0.18	3.50-4.50	4.50
PG: Abbott Cell-Dyn Ruby	7	2.75	0.28	2.00-3.60	10.00	7	4.09	0.21	3.50-4.70	5.12	7	4.45	0.26	3.60-5.20	5.91	7	2.58	0.46	1.20-4.00	17.69	7	4.01	0.20	3.40-4.60	5.00
Basophil	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%
SG: Select Abbott Cell Dyn Instruments (5-part Diff)	6	0.12	0.06	0.00-0.29	47.50	8	0.24	0.05	0.08-0.40	21.67	7	0.34	0.14	0.00-0.76	41.18	7	0.26	0.27	0.00-1.06	102.69	7	0.22	0.10	0.00-0.52	45.45
PG: Abbott Cell-Dyn Ruby	6	0.18	0.20	0.00-0.77	108.89	7	0.25	0.04	0.13-0.37	16.00	7	0.34	0.14	0.00-0.76	41.18	7	0.26	0.27	0.00-1.06	102.69	7	0.26	0.19	0.00-0.83	73.46
Red Blood Cell Count	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%
SG: Select Abbott Cell Dyn Instruments (5-part Diff)	8	2.00	0.02	1.88-2.12	1.20	8	5.22	0.15	4.91-5.53	2.89	8	4.58	0.16	4.31-4.85	3.43	8	2.00	0.04	1.88-2.12	1.90	8	5.20	0.16	4.89-5.51	3.08
PG: Abbott Cell-Dyn Ruby	7	2.00	0.02	1.88-2.12	1.20	7	5.26	0.08	4.94-5.58	1.52	7	4.62	0.11	4.34-4.90	2.32	7	2.01	0.03	1.89-2.13	1.59	7	5.25	0.08	4.93-5.57	1.60
Hemoglobin	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%
SG: Select Abbott Cell Dyn Instruments (5-part Diff)	8	5.58	0.16	5.20-6.00	2.86	8	16.41	0.46	15.30-17.50	2.80	8	13.39	0.30	12.50-14.30	2.24	8	5.61	0.12	5.20-6.00	2.14	8	16.34	0.47	15.20-17.40	2.88
PG: Abbott Cell-Dyn Ruby	7	5.55	0.15	5.20-6.00	2.68	7	16.39	0.49	15.30-17.50	2.99	7	13.43	0.30	12.50-14.30	2.24	7	5.60	0.13	5.20-6.00	2.32	7	16.37	0.50	15.30-17.50	3.05
Hematocrit	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%
SG: Select Abbott Cell Dyn Instruments (5-part Diff)	8	14.68	0.38	13.80-15.60	2.59	6	42.45	1.23	40.00-45.00	2.89	6	36.70	1.47	34.50-38.90	4.01	8	14.60	0.47	13.70-15.50	3.22	8	42.19	1.71	39.70-44.70	4.05
PG: Abbott Cell-Dyn Ruby	7	14.71	0.39	13.80-15.60	2.65	7	42.84	1.53	40.20-45.40	3.57	6	36.70	1.47	34.50-38.90	4.01	7	14.69	0.43	13.80-15.60	2.93	7	42.63	1.27	40.00-45.20	2.98
MCV	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%
SG: Select Abbott Cell Dyn Instruments (5-part Diff)	8	73.14	1.52	68.50-77.70	2.08	8	81.17	2.17	74.70-87.70	2.67	8	80.06	2.00	74.10-86.10	2.50	8	72.99	1.39	68.80-77.20	1.90	8	81.15	1.74	76.00-86.40	2.14
PG: Abbott Cell-Dyn Ruby	7	73.24	1.61	68.40-78.00	2.20	7	81.33	2.29	74.40-88.20	2.82	7	80.19	2.13	73.80-86.60	2.64	7	73.09	1.47	68.70-77.50	2.01	7	81.23	1.86	75.60-86.80	2.29
RDW	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%
SG: Select Abbott Cell Dyn Instruments (5-part Diff)	8	13.03	0.44	11.70-14.30	3.38	8	10.45	0.23	9.80-11.20	2.19	8	11.81	0.31	10.90-12.70	2.63	8	13.11	0.43	11.80-14.40	3.28	8	10.44	0.33	9.40-11.40	3.17
PG: Abbott Cell-Dyn Ruby	7	12.97	0.44	11.70-14.30	3.38	7	10.41	0.23	9.70-11.10	2.21	7	11.77	0.32	10.90-12.70	2.63	7	13.09	0.46	11.70-14.50	3.51	7	10.34	0.22	9.70-10.90	2.04
Platelet Count	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%
SG: Select Abbott Cell Dyn Instruments (5-part Diff)	8	82.31	3.24	62.00-103.00	3.90	8	494.38	27.43	371.00-618.00	5.55	8	267.88	17.73	201.00-335.00	6.60	8	81.97	6.17	62.00-103.00	7.56	8	500.75	25.24	376.00-626.00	5.03
PG: Abbott Cell-Dyn Ruby	7	81.64	2.84	62.00-103.00	3.41	7	494.29	29.63	371.00-618.00	5.99	7	269.14	18.76	202.00-336.00	6.99	7	80.69	5.38	61.00-101.00	6.67	7	503.43	26.00	377.00-629.00	5.17
683 - Antistreptolysin O (ASO);684 - C-Reactive Protein	CR-3N					CR-6N																			
CRP Quantitative	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%
600 - Cholestech LDX/CardioChek Pkg	CT-3																								
CholesterolTotal (Waived)	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%
PG: Cholestech LDX (Waived)	78	170.26	22.35	103.00-237.00	13.12																				
PG: Polymer CardioChek	81	147.57	15.35	102.00-194.00	10.34																				
Glucose (Waived)	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%
PG: Cholestech LDX (Waived)	11	50.73	2.41	44.00-58.00	4.71																				
PG: Polymer CardioChek	5	22.50	3.54	13.00-34.00	15.22																				
Triglycerides (Waived)	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%
PG: Cholestech LDX (Waived)	68	108.65	16.25	60.00-158.00	14.86																				
PG: Polymer CardioChek	77	105.17	12.17	68.00-142.00	11.62																				
Cholesterol HDL (Waived)	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%
PG: Cholestech LDX (Waived)	71	20.32	3.54	10.00-31.00	17.50																				
PG: Polymer CardioChek / Direct Method	75	84.25	13.87	42.00-126.00	16.55																				
Cholesterol LDL (Waived)	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%	Labs									



SG: Select Beckman Coulter Instruments (5-part Diff)	5	2.56	0.06	2.41-2.71	2.34	5	5.50	0.07	5.17-5.83	1.22	5	4.14	0.08	3.89-4.39	2.00	5	2.53	0.04	2.38-2.68	1.42	5	5.50	0.10	5.17-5.83	1.73
Hemoglobin	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%
SG: Select Beckman Coulter Instruments (5-part Diff)	5	6.30	0.08	5.90-6.70	1.27	5	17.08	0.13	15.90-18.30	0.76	5	11.85	0.06	11.10-12.70	0.50	5	6.38	0.05	6.00-6.80	0.78	5	17.05	0.06	15.80-18.20	0.35
Hematocrit	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%
SG: Select Beckman Coulter Instruments (5-part Diff)	5	19.42	0.46	18.20-20.60	2.37	5	53.97	0.55	50.80-57.20	1.02	5	35.53	0.67	33.40-37.60	1.89	5	19.22	0.29	18.00-20.40	1.51	5	54.05	0.90	50.80-57.20	1.67
Platelet Count	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%
SG: Select Beckman Coulter Instruments (5-part Diff)	5	87.75	1.50	66.00-110.00	1.70	5	463.50	18.45	348.00-580.00	3.97	5	261.75	12.45	197.00-328.00	4.73	5	87.25	1.50	65.00-109.00	1.72	5	480.50	18.50	361.00-601.00	3.85
MCV	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%
SG: Select Beckman Coulter Instruments (5-part Diff)	5	75.93	0.74	73.70-78.10	0.97	5	98.20	0.40	97.00-99.40	0.41	5	86.03	0.29	85.10-86.90	0.34	5	75.97	0.21	75.40-76.60	0.28	5	98.15	0.17	97.70-98.70	0.17
RDW	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%
SG: Select Beckman Coulter Instruments (5-part Diff)	5	16.28	0.10	16.00-16.60	0.61	5	13.03	0.19	12.40-13.60	1.46	5	14.10	0.18	13.60-14.60	1.28	5	16.13	0.13	15.70-16.50	0.81	5	13.08	0.24	12.40-13.80	1.83
730 - Sed Rate (ESR)	ES-3																								
Sedimentation Rate	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%
PG: Clinical Data (Vital Diag) Excyte M/10	15	7.13	1.25	3.00-11.00	17.14																				
PG: Modified Westergren - Plastic	14	8.29	1.38	4.00-12.00	17.50																				
PG: Streck ESR-Auto Plus	6	11.33	1.75	6.00-16.00	16.36																				
PG: Westergren - Glass	10	8.70	2.16	2.00-16.00	24.44																				
PG: Westergren - Plastic	17	7.71	2.11	2.00-14.00	26.25																				
PG: Polymedco Sediplast Westergren ESR	10	7.00	2.06	1.00-13.00	30.00																				
PG: Hypoguard (ARKRAY) Dispette 2 (Labsco)	11	7.36	1.86	1.00-13.00	27.14																				
PG: Sedi-Rate (Globe Scientific)	5	9.00	2.00	3.00-15.00	22.22																				
PG: Streck ESR-10 Manual	14	9.00	1.47	5.00-14.00	16.67																				
AG: All Group	111	8.05	2.22	1.00-15.00	27.50																				
672 - Glycohemoglobin	GH-5																								
Glycohemoglobin	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%
PG: Bayer DCA 2000 series	17	9.60	0.45	8.30-10.90	4.69																				
PG: Siemens Dimension/AR/ES/RxL/Xpand/EXL	13	10.96	0.31	10.10-11.90	2.82																				
PG: Chek Diagnostics (Bayer) A1C Now	36	6.61	0.70	4.50-8.70	10.61																				
SG: Select Beckman Coulter Instruments (Chem)	7	11.09	0.69	9.00-13.20	6.22																				
SG: Select Roche Instruments	5	10.98	0.29	10.10-11.90	2.64																				
PG: Siemens DCA Vantage	111	9.72	0.49	8.20-11.20	5.05																				
SG: Select Siemens (Bayer) DCA Instruments	120	9.72	0.33	8.70-10.70	3.40																				
PG: Bayer A1CNow+	35	6.70	0.97	3.80-9.60	14.48																				
SG: Select A1C Now Instruments	37	6.61	0.69	4.50-8.70	10.45																				
PG: TOSOH G8	6	9.68	0.32	8.70-10.70	3.30																				
SG: Select Abbott Architect Instruments	5	11.23	1.32	7.20-15.20	11.79																				
SG: Select Gen Chem Instruments	12	11.11	0.86	8.50-13.70	7.84																				
PG: PTS A1CNow	18	7.15	0.87	4.60-9.80	12.08																				
AG: All Group	233	8.99	1.54	4.30-13.70	17.22																				
761 - High Sensitivity CRP	HCR-3																								
hsCRP	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%
AG: All Group	5	1.10	0.11	0.77-1.43	9.91																				
602 - Hemoglobin and Hematocrit - Waived;614 - Auto Diff II (All Other 3-part);610 - Regular Hematology;601 - Waived/PPM Tests II;609 - Waived/PPM Tests I	HD-11					HD-12					HD-13					HD-14					HD-15				
Granulocyte	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%
PG: Abbott Cell-Dyn 1700	13	43.37	2.14	37.00-49.80	4.93	13	83.29	1.34	79.30-87.30	1.61	13	30.81	2.40	23.60-38.00	7.82	13	37.40	2.77	29.10-45.70	7.41	13	78.86	1.24	75.20-82.60	1.56
PG: Abbott Cell-Dyn 1800	19	44.02	3.19	34.40-53.60	7.25	19	83.53	2.15	77.00-90.00	2.57	19	31.45	2.26	24.60-38.20	7.20	19	37.41	2.52	29.80-45.00	6.74	19	78.40	2.11	72.10-84.70	2.69
PG: CDS Medonic 620	9	37.86	2.16	31.40-44.40	5.70	9	79.98	2.15	73.50-86.50	2.69	9	25.04	1.39	20.80-29.20	5.56	9	29.10	1.66	24.10-34.10	5.70	9	73.92	2.79	65.50-82.30	3.78
PG: Beckman Coulter AcT Diff/Diff 2	184	39.66	3.36	29.60-49.80	8.46	184	81.71	4.65	67.80-95.70	5.69	184	28.32	1.44	24.00-32.60	5.09	184	33.08	2.97	24.20-42.00	8.97	184	76.66	0.97	73.80-79.60	



AG: All Group	480	18.59	0.77	17.50-19.70	4.14	558	41.26	1.66	38.80-43.80	4.02	558	38.77	1.53	36.50-41.10	3.94	558	28.53	1.69	26.80-30.20	5.93	558	30.89	1.37	29.00-32.80	4.43
Hemoglobin	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%
PG: Abbott Cell-Dyn 1700	14	6.76	0.12	6.30-7.30	1.76	14	14.84	0.29	13.80-15.80	1.96	14	13.59	0.32	12.60-14.60	2.35	14	10.10	0.20	9.40-10.80	1.98	14	11.07	0.20	10.30-11.90	1.80
PG: Abbott Cell-Dyn 1800	20	6.81	0.20	6.30-7.30	2.94	20	14.90	0.87	13.90-15.90	5.84	20	13.72	0.28	12.70-14.70	2.04	20	10.17	0.27	9.50-10.90	2.65	20	11.20	0.22	10.40-12.00	1.96
PG: CDS Medonic 620	10	6.60	0.08	6.10-7.10	1.21	10	14.75	0.27	13.80-15.80	1.82	10	13.46	0.16	12.60-14.40	1.19	10	9.88	0.13	9.20-10.60	1.31	10	10.94	0.11	10.10-11.70	1.01
PG: Beckman Coulter AcT Diff/Diff 2	195	6.56	0.20	6.10-7.10	3.03	195	14.67	0.35	13.70-15.70	2.38	195	13.47	0.27	12.60-14.40	2.00	195	9.93	0.22	9.20-10.60	2.12	195	10.89	0.27	10.10-11.70	2.48
PG: Horiba Micros/45/60	58	6.59	0.14	6.10-7.10	2.12	58	14.58	0.27	13.60-15.60	1.85	57	13.45	0.26	12.50-14.30	1.94	58	9.91	0.15	9.20-10.60	1.52	58	10.89	0.18	10.10-11.70	1.65
SG: Select Abbott Cell Dyn Instruments (3-part Diff)	134	6.57	0.22	6.10-7.10	3.18	134	14.62	0.50	13.60-15.60	3.36	134	13.46	0.28	12.60-14.40	2.07	134	9.92	0.25	9.20-10.60	2.53	134	10.88	0.27	10.10-11.70	2.48
PG: Drew Scientific D3	15	6.55	0.24	6.00-7.00	3.69	15	14.50	0.60	13.50-15.50	4.14	15	13.55	0.40	12.60-14.40	2.96	15	9.99	0.31	9.30-10.70	3.10	15	10.92	0.33	10.10-11.70	3.03
PG: Abbott Cell-Dyn Emerald	100	6.49	0.17	6.00-7.00	2.62	100	14.53	0.38	13.50-15.50	2.62	100	13.40	0.23	12.50-14.30	1.72	100	9.84	0.21	9.10-10.50	2.14	100	10.79	0.22	10.00-11.60	2.04
PG: CDS M Series	144	6.68	0.17	6.20-7.20	2.54	144	14.97	0.52	13.90-16.10	3.47	143	13.51	0.38	12.60-14.40	2.81	144	9.90	0.67	9.20-10.60	6.77	144	11.02	0.31	10.20-11.80	2.82
SG: Select CDS Instruments	154	6.68	0.17	6.20-7.20	2.54	154	14.96	0.51	13.90-16.10	3.40	153	13.51	0.37	12.60-14.40	2.74	154	9.90	0.65	9.20-10.60	6.57	154	11.01	0.30	10.20-11.80	2.73
AG: All Group	557	6.60	0.20	6.10-7.10	3.03	557	14.72	0.46	13.70-15.70	3.13	555	13.48	0.30	12.60-14.40	2.30	557	9.92	0.39	9.20-10.60	3.94	557	10.92	0.28	10.10-11.70	2.57
Lymphocyte	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%
PG: Abbott Cell-Dyn 1700	12	38.22	2.66	30.23-46.21	6.97	12	8.43	0.83	5.95-10.91	9.82	13	49.62	3.22	39.97-59.27	6.48	13	41.57	3.60	30.77-52.37	8.66	12	12.38	0.71	10.26-14.50	5.71
PG: Abbott Cell-Dyn 1800	19	34.28	5.08	19.05-49.51	14.80	18	7.61	0.70	5.50-9.72	9.24	19	45.09	5.39	28.93-61.25	11.95	19	38.14	4.87	23.54-52.74	12.76	19	11.38	1.58	6.64-16.12	13.89
PG: CDS Medonic 620	9	53.02	2.33	46.03-60.01	4.39	9	14.54	3.00	5.55-23.52	20.60	9	65.74	1.27	61.92-69.56	1.94	9	60.08	0.88	57.45-62.71	1.46	9	17.73	1.47	13.33-22.13	8.27
PG: Beckman Coulter AcT Diff/Diff 2	184	44.59	1.67	39.59-49.59	3.74	179	8.36	0.66	6.39-10.33	7.85	184	58.00	1.46	53.63-62.37	2.51	184	48.54	1.69	43.47-53.61	3.48	184	14.48	0.63	12.60-16.36	4.32
PG: Horiba Micros/45/60	55	28.77	2.49	21.30-36.24	8.65	57	9.05	0.87	6.44-11.66	9.61	57	39.79	3.25	30.03-49.55	8.18	56	32.72	2.93	23.92-41.52	8.97	57	10.58	0.84	8.07-13.09	7.92
SG: Select Abbott Cell Dyn Instruments (3-part Diff)	130	37.48	3.17	27.97-46.99	8.46	112	10.48	1.69	5.40-15.56	16.16	131	49.74	3.62	38.87-60.61	7.28	130	41.70	3.14	32.28-51.12	7.53	129	13.05	1.21	9.42-16.68	9.26
PG: Drew Scientific D3	15	41.73	2.75	33.49-49.96	6.58	14	10.20	0.81	7.77-12.63	7.94	15	53.43	1.61	48.61-58.25	3.01	15	44.95	1.64	40.03-49.87	3.65	15	13.51	0.58	11.76-15.26	4.32
PG: Abbott Cell-Dyn Emerald	99	38.00	2.32	31.05-44.95	6.10	85	11.32	1.03	8.24-14.40	9.06	99	50.64	2.41	43.42-57.86	4.76	98	42.41	2.03	36.31-48.51	4.79	98	13.46	0.80	11.06-15.86	5.95
PG: CDS M Series	143	52.66	1.96	46.78-58.54	3.72	140	14.47	3.05	5.33-23.61	21.06	143	65.04	4.37	51.94-78.14	6.71	143	58.94	2.48	51.50-66.38	4.21	142	17.85	0.81	15.43-20.27	4.53
SG: Select CDS Instruments	152	52.68	1.98	46.75-58.61	3.75	149	14.47	3.03	5.37-23.57	20.97	152	65.08	4.25	52.34-77.82	6.53	152	59.01	2.43	51.73-66.29	4.11	151	17.84	0.85	15.28-20.40	4.79
AG: All Group	516	44.16	7.08	22.91-65.41	16.04	454	9.84	2.02	3.77-15.91	20.55	516	56.87	7.42	34.59-79.14	13.06	525	48.56	8.20	23.97-73.15	16.88	524	14.76	2.34	7.73-21.79	15.87
MCV	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%
PG: Abbott Cell-Dyn 1700	14	89.73	2.20	83.10-96.30	2.45	14	83.26	2.02	77.20-89.40	2.42	14	83.03	1.87	77.40-88.60	2.25	14	72.05	1.81	66.70-77.50	2.51	14	70.71	1.42	66.40-75.00	2.01
PG: Abbott Cell-Dyn 1800	20	92.65	1.15	89.30-96.20	1.24	20	85.67	1.22	82.00-89.40	1.42	20	85.14	1.17	81.60-88.60	1.37	20	74.09	0.94	71.30-76.90	1.27	20	73.26	4.36	60.20-86.40	5.95
PG: CDS Medonic 620	9	85.66	1.85	80.20-91.30	2.16	9	81.37	1.63	76.50-86.30	2.00	9	79.66	1.66	74.70-84.70	2.08	9	69.32	1.37	65.20-73.40	1.98	9	69.20	1.28	65.40-73.00	1.85
PG: Beckman Coulter AcT Diff/Diff 2	195	89.28	1.18	85.80-92.80	1.32	195	82.97	1.16	79.50-86.50	1.40	195	82.15	1.12	78.70-85.50	1.36	195	72.40	1.02	69.30-75.50	1.41	195	70.58	4.32	57.60-83.60	6.12
PG: Horiba Micros/45/60	50	83.54	1.03	80.40-86.60	1.23	50	79.28	0.83	76.80-81.80	1.05	50	79.12	0.82	76.60-81.60	1.04	50	69.30	0.84	66.80-71.80	1.21	50	67.18	0.77	64.90-69.50	1.15
SG: Select Abbott Cell Dyn Instruments (3-part Diff)	133	94.12	2.21	87.40-100.80	2.36	133	86.46	1.78	81.20-91.80	2.06	133	85.70	1.61	80.90-90.50	1.88	133	75.04	1.67	70.00-80.00	2.23	133	73.16	2.17	66.70-79.70	2.96
PG: Drew Scientific D3	15	93.09	1.33	89.10-97.10	1.43	15	85.14	1.16	81.60-88.60	1.36	15	84.49	1.08	81.30-87.70	1.28	15	74.31	1.04	71.20-77.40	1.40	15	72.51	1.01	69.50-75.50	1.39
PG: Abbott Cell-Dyn Emerald	99	95.04	1.29	91.10-98.90	1.36	99	87.07	1.22	83.40-90.80	1.40	99	86.20	1.20	82.60-89.80	1.39	99	75.66	1.13	72.30-79.10	1.49	99	73.48	1.19	69.90-77.10	1.62
PG: CDS M Series	143	85.50	1.40	81.30-89.70	1.64	143	81.33	1.33	77.30-85.30	1.64	143	79.19	1.44	74.90-83.50	1.82	143	68.87	2.14	62.50-75.30	3.11	143	69.13	2.23	62.40-75.80	3.23
SG: Select CDS Instruments	152	85.51	1.42	81.20-89.80	1.66	152	81.33	1.34	77.30-85.30	1.65	152	79.22	1.46	74.80-83.60	1.84	152	68.90	2.10	62.60-75.20	3.05	152	69.13	2.19	62.60-75.60	3.15
AG: All Group	546	89.01	3.94	77.20-100.80	4.43	546	83.10	2.65	75.10-91.00	3.19	546	82.00	2.88	73.40-90.60	3.51	546	71.85	2.86	63.20-80.40	3.98	546	70.55	3.54	59.90-81.30	5.03
Monocyte/Mixed	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%
PG: Abbott Cell-Dyn 1700	13	18.97	2.13	12.60-25.40	11.21	13	8.08	1.01	5.10-11.10	12.47	13	19.58	2.05	13.50-25.80	10.46	13	21.03	1.60	16.20-25.80	7.62	13	8.95	0.87	6.40-11.60	9.67
PG: Abbott Cell-Dyn 1800	18	22.00	1.98	16.10-27.90	9.00	17	8.39	0.75	6.20-10.70	8.93	18	23.92	2.89	15.20-32.60	12.09	18	24.78	2.15	18.40-31.30	8.67	18	10.20	0.80	7.80-12.60	7.84
PG: CDS Medonic 620	9	9.12	1.58	4.40-13.80	17.36	9	5.48	1.21	1.90-9.10	22.00	9	9.21	1.85	3.60-14.80	20.11	9	10.82	1.95	5.00-16.60	18.06	9	8.34	1.47	3.90-12.70	17.71
PG: Beckman Coulter AcT Diff/Diff 2	184	15.55	1.57	10.90-20.30	10.06	184	9.50	0.96	6.60-12.40	10.11	184	13.68	1.20	10.10-17.30	8.76	184	18.21	1.50	13.70-22.70	8.24	184	8.82	0.88	6.20-11.40	10.00
PG: Horiba Micros/45/60	58	18.54	1.68	13.50-23.50	9.08	58	3.42	0.55	1.70-5.00	16.18	58	21.48	3.36	11.40-31.60	15.63	58	19.68	1.34	15.70-23.70	6.80	58	5.54	0.39	4.40-6.60	6.91
SG: Select Abbott Cell Dyn Instruments (3-part Diff)	121	13.93	3.06	4.70-23.10	22.01	120	3.91	1.85	0.00-9.50	47.44	119	15.40	2.92	6.60-24.20	18.96	118	15.39	3.24	5.70-25.10	21.04	110	4.96	1.29	1.10-8.90	26.00
PG: Abbott Cell-Dyn Emerald	99	12.67	1.42	8.40-17.00	11.18	99	3.08	0.26	2.30-3.90	8.39	99	14.43	1.87	8.80-20.00	12.99	99	14.20	1.69	9.20-19.20	11.83	99	4.55	0.42	3.30-5.90	9.13
PG: CDS M Series	139	10.50	1.74	5.30-15.70	16.57	138	5.92	1.09	2.60-9.20	18.47	136	8.90	1.46	4.50-13.30	16.40	139	10.62	1.81	5.10-16.10	17.17	143	9.76	1.08	6.60-13.00	11.02
AG: All Group	527</																								



SG: Select CDS Instruments	152	12.06	0.77	9.80-14.40	6.36	152	14.06	1.62	9.20-19.00	11.49	151	14.43	0.80	12.00-16.80	5.56	152	14.34	1.20	10.70-17.90	8.39	152	14.75	1.35	10.60-18.80	9.18
AG: All Group	532	12.48	0.83	10.00-15.00	6.64	532	14.28	1.26	10.50-18.10	8.81	531	14.65	0.97	11.80-17.60	6.60	532	15.04	1.33	11.00-19.00	8.87	532	15.45	1.41	11.20-19.60	9.16
Red Blood Cell Count	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%
PG: Abbott Cell-Dyn 1700	14	2.16	0.07	2.03-2.29	3.01	14	4.99	0.11	4.69-5.29	2.18	14	4.76	0.11	4.47-5.05	2.25	14	4.04	0.06	3.80-4.28	1.51	14	4.41	0.10	4.15-4.67	2.20
PG: Abbott Cell-Dyn 1800	20	2.19	0.07	2.06-2.32	3.11	19	4.87	0.29	4.58-5.16	5.91	20	4.68	0.09	4.40-4.96	2.01	20	4.00	0.10	3.76-4.24	2.58	20	4.38	0.11	4.12-4.64	2.47
PG: CDS Medonic 620	10	2.10	0.03	1.97-2.23	1.19	10	5.00	0.04	4.70-5.30	0.84	10	4.78	0.04	4.49-5.07	0.84	10	3.96	0.06	3.72-4.20	1.39	10	4.35	0.05	4.09-4.61	1.17
PG: Beckman Coulter AcT Diff/Diff 2	195	2.13	0.07	2.00-2.26	3.33	195	5.03	0.15	4.73-5.33	2.90	195	4.80	0.14	4.51-5.09	2.83	195	4.05	0.13	3.81-4.29	3.26	195	4.45	0.16	4.18-4.72	3.64
PG: Horiba Micros/45/60	59	2.07	0.05	1.95-2.19	2.22	59	4.96	0.11	4.66-5.26	2.24	59	4.70	0.10	4.42-4.98	2.21	59	3.94	0.07	3.70-4.18	1.75	59	4.36	0.09	4.10-4.62	2.04
SG: Select Abbott Cell Dyn Instruments (3-part Diff)	134	2.09	0.09	1.96-2.22	4.40	133	4.84	0.20	4.55-5.13	4.07	134	4.62	0.14	4.34-4.90	3.03	134	3.92	0.13	3.68-4.16	3.42	134	4.30	0.13	4.04-4.56	3.05
PG: Drew Scientific D3	15	2.15	0.06	2.02-2.28	2.56	15	4.94	0.17	4.64-5.24	3.46	15	4.72	0.11	4.44-5.00	2.29	15	4.02	0.07	3.78-4.26	1.72	15	4.40	0.13	4.14-4.66	2.98
PG: Abbott Cell-Dyn Emerald	100	2.06	0.08	1.94-2.18	3.74	100	4.81	0.18	4.52-5.10	3.66	100	4.59	0.14	4.31-4.87	2.98	100	3.89	0.13	3.66-4.12	3.37	100	4.28	0.13	4.02-4.54	2.97
PG: CDS M Series	144	2.11	0.06	1.98-2.24	2.65	144	5.01	0.18	4.71-5.31	3.63	144	4.76	0.13	4.47-5.05	2.69	144	3.92	0.27	3.68-4.16	6.79	144	4.34	0.09	4.08-4.60	2.17
SG: Select CDS Instruments	154	2.11	0.05	1.98-2.24	2.56	154	5.01	0.18	4.71-5.31	3.51	154	4.76	0.12	4.47-5.05	2.61	154	3.93	0.26	3.69-4.17	6.56	154	4.34	0.09	4.08-4.60	2.12
AG: All Group	558	2.11	0.07	1.98-2.24	3.51	557	4.97	0.18	4.67-5.27	3.66	558	4.73	0.15	4.45-5.01	3.11	558	3.97	0.18	3.73-4.21	4.56	558	4.38	0.14	4.12-4.64	3.26
White Blood Cell Count	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%
PG: Abbott Cell-Dyn 1700	14	2.65	0.18	2.30-3.10	6.67	14	26.34	1.21	22.40-30.20	4.60	14	5.41	0.21	4.60-6.20	3.89	14	4.43	0.21	3.70-5.10	4.77	14	16.54	0.69	14.00-19.00	4.18
PG: Abbott Cell-Dyn 1800	20	2.38	0.08	2.00-2.80	3.33	20	24.90	1.67	21.20-28.60	6.71	20	4.84	0.28	4.10-5.50	5.83	20	3.95	0.22	3.40-4.60	5.50	20	15.43	0.66	13.10-17.70	4.22
PG: CDS Medonic 620	10	2.31	0.06	2.00-2.60	2.61	10	25.71	0.68	21.80-29.60	2.65	10	4.78	0.24	4.10-5.50	5.00	10	3.91	0.21	3.30-4.50	5.38	10	15.75	0.41	13.40-18.20	2.59
PG: Beckman Coulter AcT Diff/Diff 2	195	2.55	0.13	2.10-2.90	5.20	195	26.31	0.94	22.40-30.20	3.57	195	5.24	0.19	4.40-6.00	3.65	195	4.25	0.21	3.70-4.90	4.88	195	16.00	0.56	13.60-18.40	3.50
PG: Horiba Micros/45/60	59	2.43	0.09	2.00-2.80	3.75	59	25.72	0.65	21.80-29.60	2.53	59	5.08	0.18	4.30-5.90	3.53	59	4.06	0.12	3.50-4.70	2.68	59	15.58	0.49	13.30-17.90	3.14
SG: Select Abbott Cell Dyn Instruments (3-part Diff)	134	2.50	0.15	2.10-2.90	6.00	134	24.91	1.57	21.20-28.60	6.31	134	5.14	0.28	4.30-5.90	5.49	133	4.15	0.21	3.60-4.80	5.00	134	15.30	0.73	13.00-17.60	4.77
PG: Drew Scientific D3	15	2.58	0.12	2.20-3.00	4.23	15	24.77	1.14	21.10-28.50	4.60	15	5.18	0.16	4.40-6.00	3.08	15	4.22	0.10	3.60-4.80	2.38	15	15.41	0.39	13.10-17.70	2.53
PG: Abbott Cell-Dyn Emerald	100	2.51	0.14	2.10-2.90	5.60	100	24.71	1.50	21.00-28.40	6.07	100	5.16	0.23	4.40-6.00	4.42	99	4.16	0.16	3.60-4.80	3.81	100	15.10	0.55	12.80-17.40	3.64
PG: CDS M Series	144	2.34	0.16	2.00-2.60	6.52	144	25.94	1.83	22.00-29.80	7.07	142	4.78	0.17	4.10-5.50	3.54	144	3.88	0.30	3.30-4.50	7.69	144	15.78	1.25	13.40-18.20	7.91
SG: Select CDS Instruments	154	2.34	0.15	2.00-2.60	6.52	154	25.92	1.77	22.00-29.80	6.83	152	4.78	0.17	4.10-5.50	3.54	154	3.88	0.29	3.30-4.50	7.44	154	15.78	1.22	13.40-18.20	7.72
AG: All Group	558	2.47	0.16	2.10-2.90	6.40	558	25.76	1.47	21.90-29.70	5.70	556	5.07	0.28	4.30-5.90	5.49	557	4.11	0.27	3.50-4.70	6.59	558	15.71	0.86	13.30-18.10	5.48
MCH	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%
PG: Horiba Micros/45/60	5	31.65	0.07	31.50-31.90	0.22	5	29.65	0.78	27.40-32.00	2.63	5	29.05	0.07	28.80-29.20	0.24	5	25.30	0.14	24.90-25.70	0.55	5	25.15	0.21	24.60-25.80	0.83
MCHC	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%
PG: Horiba Micros/45/60	5	37.80	0.42	36.50-39.10	1.11	5	37.40	0.71	35.30-39.50	1.90	5	36.75	0.21	36.20-37.40	0.57	5	36.55	0.07	36.30-36.70	0.19	5	37.50	0.85	35.00-40.00	2.27
Hematocrit (Waived)	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%
PG: Clay Adam Auto,Readacrit,Triac	5	18.36	2.04	12.30-24.50	11.09																				
PG: Microhematocrit, Spun	19	16.74	1.56	12.00-212																					



PG: Abbott i-Stat (Waived) / Thermal	26	17.77	0.65	15.90-19.80	3.65																				
Glucose	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%
PG: Abbott i-Stat / Thermal	11	37.18	1.17	31.20-43.20	3.24	9	87.44	1.59	78.00-96.00	1.84	9	195.78	2.86	176.00-216.00	1.48						9	37.56	1.33	32.00-44.00	3.54
SG: International Abbott i-Stat Instruments / Thermal	41	37.03	1.01	31.00-43.00	2.73	13	87.85	1.51	79.00-97.00	1.72	13	195.56	2.52	176.00-216.00	1.28	13	51.26	0.86	45.30-57.30	1.76	13	37.64	1.10	32.00-44.00	2.93
PG: Abbott i-Stat (Waived) / Thermal	26	36.77	0.81	30.80-42.80	2.23																				
Hematocrit	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%
PG: Abbott i-Stat / Thermal	6	19.17	0.75	18.00-20.40	3.91	6	36.50	0.84	34.30-38.70	2.30	6	32.67	1.03	30.70-34.70	3.15	6	15.50	0.55	14.60-16.40	3.55	6	18.83	0.75	17.70-19.90	3.99
SG: International Abbott i-Stat Instruments / Thermal	20	18.75	0.72	17.70-19.90	3.83	10	36.30	0.68	34.10-38.50	1.85	10	32.60	0.84	30.60-34.60	2.58	10	15.60	0.52	14.70-16.50	3.33	10	18.80	0.63	17.70-19.90	3.35
PG: Abbott i-Stat (Waived) / Thermal	10	18.40	0.70	17.30-19.50	3.80																				
Hemoglobin	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%
PG: Abbott i-Stat / Thermal	6	6.57	0.31	6.10-7.10	4.70	6	12.45	0.39	11.60-13.40	3.12	6	11.20	0.48	10.40-12.00	4.29	6	5.18	0.17	4.80-5.60	3.27	6	6.33	0.31	5.90-6.70	4.92
SG: International Abbott i-Stat Instruments / Thermal	20	6.42	0.27	6.00-6.80	4.22	10	12.35	0.32	11.40-13.20	2.60	10	11.14	0.38	10.30-11.90	3.42	10	5.24	0.17	4.80-5.60	3.27	10	6.36	0.26	6.00-6.80	4.06
PG: Abbott i-Stat (Waived) / Thermal	10	6.29	0.26	5.90-6.70	4.13																				
Ionized Calcium	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%
PG: Abbott i-Stat / Thermal	6	2.49	0.01	2.46-2.52	0.40	6	1.30	0.02	1.25-1.35	1.23	6	0.97	0.02	0.91-1.03	2.06	6	0.60	0.01	0.58-0.62	1.33	6	2.50	0.00	2.50-2.50	0.00
SG: International Abbott i-Stat Instruments / Thermal	33	2.50	0.01	2.47-2.53	0.36	10	1.31	0.01	1.27-1.35	1.07	10	0.97	0.02	0.92-1.02	1.65	10	0.60	0.01	0.58-0.62	1.00	10	2.50	0.00	2.50-2.50	0.00
PG: Abbott i-Stat (Waived) / Thermal	23	2.50	0.01	2.47-2.53	0.36																				
Potassium (K+)	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%
PG: Abbott i-Stat / Thermal	14	3.41	0.03	2.90-3.90	0.88	12	4.40	0.00	3.90-4.90	0.00	12	4.60	0.00	4.10-5.10	0.00	12	2.25	0.05	1.80-2.80	2.17	12	3.40	0.00	2.90-3.90	0.00
SG: International Abbott i-Stat Instruments / Thermal	43	3.41	0.03	2.90-3.90	0.88	16	4.40	0.00	3.90-4.90	0.00	16	4.60	0.00	4.10-5.10	0.00	16	2.26	0.05	1.80-2.80	2.17	16	3.40	0.00	2.90-3.90	0.00
PG: Abbott i-Stat (Waived) / Thermal	25	3.41	0.03	2.90-3.90	0.88																				
Sodium (NA+)	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%
PG: Abbott i-Stat / Thermal	11	131.73	0.91	128.00-136.00	0.68	9	138.11	1.17	134.10-142.10	0.87	9	156.00	1.00	152.00-160.00	0.64	9	126.56	0.73	123.00-131.00	0.58	9	131.56	0.73	128.00-136.00	0.55
SG: International Abbott i-Stat Instruments / Thermal	41	131.90	0.63	128.00-136.00	0.47	13	138.23	1.01	134.20-142.20	0.73	13	156.00	0.82	152.00-160.00	0.53	13	126.69	0.75	123.00-131.00	0.63	13	131.77	0.72	128.00-136.00	0.55
PG: Abbott i-Stat (Waived) / Thermal	26	131.92	0.48	128.00-136.00	0.38																				
Urea Nitrogen (BUN)	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%
PG: Abbott i-Stat / Thermal	8	5.38	0.52	3.40-7.40	9.63	6	35.50	0.55	32.30-38.70	1.55	6	55.67	0.52	50.70-60.70	0.93	6	9.33	0.52	7.30-11.30	5.59	6	5.50	0.55	3.50-7.50	10.00
SG: International Abbott i-Stat Instruments / Thermal	38	5.71	0.44	3.70-7.70	7.72	10	35.47	0.42	32.30-38.70	1.18	10	55.56	0.51	50.60-60.60	0.92	10	9.41	0.40	7.40-11.40	4.26	10	5.60	0.44	3.60-7.60	7.86
PG: Abbott i-Stat (Waived) / Thermal	26	5.81	0.40	3.80-7.80	6.90																				
HCO3	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%
PG: Abbott i-Stat	5	13.93	0.15	13.50-14.30	1.08	5	26.00	0.36	24.90-27.10	1.38	5	23.23	0.40	22.00-24.40	1.72	5	21.20	0.85	18.60-23.80	4.01	5	13.93	0.15	13.50-14.30	1.08
SG: International Abbott i-Stat Instruments	7	14.26	0.35	13.30-15.40	2.45	7	26.36	0.62	24.50-28.30	2.35	7	23.61	0.49	22.10-25.10	2.08	7	21.46	0.63	19.60-23.40	2.93	7	14.10	0.25	13.30-14.80	1.77
Lactate	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%
SG: International Abbott i-Stat Instruments	5	5.02	0.14	4.59-5.45	2.85	5	0.76	0.03	0.68-0.84	3.55	5	2.26	0.04	2.14-2.38	1.73	5	6.47	0.12	6.12-6.82	1.82	5	5.03	0.10	4.74-5.32	1.91
PG: Intl Abbott i-Stat	5	5.02	0.17	4.52-5.51	3.29	5	0.77	0.02	0.71-0.83	2.47	5	2.25	0.04	2.14-2.36	1.60	5	6.47	0.14	6.06-6.88	2.10	5	5.02	0.11	4.70-5.34	2.13
Anion Gap	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%
SG: International Abbott i-Stat Instruments / Thermal	8	6.13	1.55	1.00-11.00	26.67																				
Creatinine	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%	Labs	Mean</			

Print Date: 11/21/2018



PG: Beckman ACCESS	11	740.73	50.91	588.00-894.00	6.87																				
PG: OCD Vitros ECI/ECIQ	5	851.20	46.81	711.00-991.00	5.50																				
SG: Select Beckman Coulter Instruments (Chem)	15	739.67	51.20	586.00-894.00	6.92																				
SG: Select Roche Instruments	5	928.36	70.15	717.00-1139.00	7.56																				
SG: Select OCD Vitros Instruments	6	821.83	83.23	572.00-1072.00	10.12																				
SG: Select Abbott Architect Instruments	11	905.82	115.21	560.00-1252.00	12.72																				
AG: All Group	43	834.53	106.67	515.00-1155.00	12.78																				
685 - Special Chemistry II	SPC-5					SPC-6																			
Insulin	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%
PG: Siemens Immulite 1000	5	60.90	0.00	60.90-60.90	0.00	5	24.00	0.00	24.00-24.00	0.00															
PG: Siemens Immulite 2000	5	42.50	0.00	42.50-42.50	0.00	5	17.20	0.00	17.20-17.20	0.00															
AG: All Group	5	107.62	5.51	91.10-124.10	5.12	5	37.07	2.03	31.00-43.20	5.47															
PTH	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%
AG: All Group	8	147.43	21.68	82.40-212.40	14.71	8	17.86	7.45	0.00-40.30	41.62															
Vitamin D	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%
PG: Beckman ACCESS	17	103.60	11.43	69.31-137.89	11.03	17	29.88	4.17	17.37-42.39	13.95															
PG: TOSOH AIA	7	97.14	4.48	83.69-110.59	4.61	7	41.96	2.25	35.22-48.70	5.36															
SG: Select Beckman Coulter Instruments (Chem)	21	105.07	11.54	70.45-139.69	10.98	21	29.39	4.56	15.70-43.08	15.52															
PG: Roche Cobas e 411	6	65.71	13.99	23.72-107.69	21.30	6	19.54	1.58	14.79-24.29	8.10															
PG: Abbott Architect i1000 SR	9	99.86	2.92	91.09-108.63	2.93	9	27.41	0.92	24.64-30.18	3.37															
SG: Select Abbott Architect Instruments	15	98.51	3.69	87.43-109.59	3.75	15	27.19	1.35	23.13-31.25	4.97															
PG: Qualigen FastPack Vitamin D	6	61.05	4.67	47.04-75.06	7.65	6	28.68	6.28	9.85-47.51	21.89															
AG: All Group	60	93.82	18.23	39.13-148.51	19.43	61	28.74	6.15	10.28-47.20	21.41															
752 - Auto Diff 8 (Sysmex 5-part Diff)	SXT-11					SXT-12					SXT-13					SXT-14					SXT-15				
White Blood Cell Count	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%
PG: Sysmex XS-1000i	11	3.80	0.11	3.20-4.40	2.89	11	20.62	0.39	17.50-23.70	1.89	11	7.80	0.15	6.60-9.00	1.92	11	3.80	0.07	3.20-4.40	1.84	11	20.69	0.59	17.60-23.80	2.85
SG: Select Sysmex Instruments (5-part Diff)	13	3.75	0.20	3.10-4.30	5.41	13	20.64	0.36	17.50-23.70	1.75	13	7.72	0.29	6.50-8.90	3.77	13	3.75	0.15	3.20-4.40	3.95	13	20.72	0.55	17.60-23.80	2.66
Lymphocyte	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%
PG: Sysmex XS-1000i	11	16.48	0.97	13.60-19.40	5.94	11	19.23	1.35	15.10-23.30	7.03	11	29.53	0.59	27.70-31.30	2.00	11	17.17	0.79	14.90-19.50	4.53	11	19.10	1.43	14.80-23.40	7.43
SG: Select Sysmex Instruments (5-part Diff)	13	16.58	1.25	12.90-20.40	7.53	13	19.55	1.50	15.00-24.00	7.74	13	29.38	1.83	23.90-34.90	6.22	13	17.13	1.03	14.00-20.20	6.02	13	19.44	1.55	14.70-24.00	7.99
Neutrophil	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%
PG: Sysmex XS-1000i	11	58.94	1.02	55.80-62.00	1.73	11	56.19	1.80	50.80-61.60	3.20	11	48.28	0.60	46.50-50.10	1.24	11	58.23	0.93	55.40-61.00	1.60	11	56.54	1.42	52.20-60.80	2.51
SG: Select Sysmex Instruments (5-part Diff)	13	59.32	2.84	50.80-67.80	4.79	13	56.50	2.21	49.90-63.10	3.89	13	48.75	2.72	40.50-56.90	5.59	13	58.85	2.85	50.30-67.50	4.84	13	56.94	1.79	51.50-62.30	3.15
Mono/Mixed	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%
PG: Sysmex XS-1000i	11	0.69	0.15	0.30-1.10	21.43	11	1.07	0.19	0.60-1.60	16.36	11	1.66	0.35	0.70-2.80	20.59	11	0.64	0.16	0.10-1.10	26.67	11	1.06	0.16	0.60-1.60	14.55
SG: Select Sysmex Instruments (5-part Diff)	12	0.71	0.16	0.20-1.20	22.86	12	1.04	0.21	0.40-1.60	21.00	12	1.63	0.35	0.60-2.60	21.88	12	0.63	0.15	0.20-1.00	25.00	12	1.03	0.19	0.40-1.60	19.00
Eosinophil	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%
PG: Sysmex XS-1000i	11	15.17	0.54	13.60-16.80	3.55	11	14.37	0.67	12.40-16.40	4.65	11	13.00	0.31	12.10-13.90	2.38	11	14.97	0.52	13.40-16.60	3.47	11	14.11	0.58	12.40-15.80	4.11
SG: Select Sysmex Instruments (5-part Diff)	13	15.14	0.50	13.60-16.60	3.31	13	14.31	0.64	12.40-16.20	4.48	13	13.01	0.30	12.10-13.90	2.31	13	15.00	0.49	13.50-16.50	3.27	13	14.06	0.58	12.40-15.80	4.11
Basophil	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%
PG: Sysmex XS-1000i	11	8.45	1.02	5.39-11.51	12.06	11	9.14	0.87	6.52-11.76	9.54	11	7.53	0.76	5.24-9.82	10.13	11	8.99	0							



Neutrophil	Labs	Mean	SD	Range	CV %	Labs	Mean	SD	Range	CV %	Labs	Mean	SD	Range	CV %	Labs	Mean	SD	Range	CV %	Labs	Mean	SD	Range	CV %
PG: Sysmex K-800,1000,4500,KX-21	17	69.65	5.01	54.60-84.80	7.20	17	23.81	0.84	21.30-26.30	3.53	17	56.32	4.18	43.80-68.80	7.42	16	71.11	1.55	66.40-75.80	2.18	16	23.66	0.84	21.20-26.20	3.54
PG: Sysmex pocH-100i	22	73.32	1.28	69.50-77.10	1.75	22	24.81	0.73	22.60-27.00	2.98	22	57.90	1.05	54.80-61.00	1.81	22	73.75	1.47	69.40-78.20	1.99	22	24.65	0.93	21.90-27.50	3.77
SG: Select Sysmex Instruments (3-part Diff)	117	71.26	2.69	63.20-79.40	3.77	124	24.04	1.30	20.10-27.90	5.42	119	55.77	2.31	48.90-62.70	4.14	117	70.79	4.86	56.20-85.40	6.86	122	23.97	1.31	20.10-27.90	5.46
PG: Sysmex XP-300	78	71.03	1.84	65.50-76.50	2.59	78	23.69	0.71	21.60-25.80	3.00	78	55.00	1.46	50.60-59.40	2.65	78	70.48	2.02	64.40-76.60	2.87	78	23.67	0.82	21.20-26.20	3.46
SG: International Sysmex Instruments	26	73.42	1.25	69.70-77.10	1.69	26	24.89	0.80	22.50-27.30	3.21	26	57.83	1.01	54.80-60.80	1.75	26	73.82	1.39	69.70-77.90	1.87	26	24.60	0.87	22.00-27.20	3.54
Red Blood Cell Count	Labs	Mean	SD	Range	CV %	Labs	Mean	SD	Range	CV %	Labs	Mean	SD	Range	CV %	Labs	Mean	SD	Range	CV %	Labs	Mean	SD	Range	CV %
PG: Sysmex K-800,1000,4500,KX-21	16	2.40	0.03	2.26-2.54	1.33	17	6.00	0.10	5.64-6.36	1.65	16	4.21	0.05	3.96-4.46	1.09	16	2.40	0.04	2.26-2.54	1.75	17	5.87	0.41	5.52-6.22	6.97
PG: Sysmex pocH-100i	23	2.42	0.05	2.27-2.57	2.15	23	6.06	0.15	5.70-6.42	2.44	23	4.26	0.09	4.00-4.52	2.16	23	2.44	0.06	2.29-2.59	2.30	23	6.05	0.14	5.69-6.41	2.23
SG: Select Sysmex Instruments (3-part Diff)	119	2.42	0.13	2.27-2.57	5.17	127	5.98	0.13	5.62-6.34	2.14	126	4.24	0.22	3.99-4.49	5.26	119	2.41	0.09	2.27-2.55	3.73	125	5.96	0.18	5.60-6.32	3.04
PG: Sysmex XP-300	80	2.42	0.15	2.27-2.57	6.16	80	5.95	0.11	5.59-6.31	1.82	80	4.23	0.27	3.98-4.48	6.45	80	2.40	0.10	2.26-2.54	4.25	79	5.94	0.07	5.58-6.30	1.21
SG: International Sysmex Instruments	27	2.43	0.05	2.28-2.58	2.10	27	6.06	0.14	5.70-6.42	2.33	27	4.25	0.09	4.00-4.50	2.12	27	2.43	0.06	2.28-2.58	2.35	27	6.04	0.13	5.68-6.40	2.19
Hemoglobin	Labs	Mean	SD	Range	CV %	Labs	Mean	SD	Range	CV %	Labs	Mean	SD	Range	CV %	Labs	Mean	SD	Range	CV %	Labs	Mean	SD	Range	CV %
PG: Sysmex K-800,1000,4500,KX-21	16	6.09	0.14	5.70-6.50	2.30	17	18.46	0.49	17.20-19.80	2.65	16	12.02	0.27	11.20-12.80	2.25	16	6.13	0.17	5.70-6.50	2.79	16	18.31	0.40	17.00-19.60	2.19
PG: Sysmex pocH-100i	23	6.20	0.12	5.80-6.60	1.94	23	18.47	0.45	17.20-19.80	2.43	23	12.03	0.20	11.20-12.80	1.67	23	6.20	0.13	5.80-6.60	2.10	23	18.38	0.35	17.10-19.70	1.96
SG: Select Sysmex Instruments (3-part Diff)	119	6.13	0.32	5.70-6.50	5.25	127	18.43	0.39	17.10-19.70	2.12	126	12.13	0.70	11.30-12.90	5.79	119	6.16	0.24	5.80-6.60	3.87	125	18.39	0.32	17.10-19.70	1.74
PG: Sysmex XP-300	80	6.11	0.38	5.70-6.50	6.23	80	18.40	0.36	17.10-19.70	1.96	80	12.19	0.86	11.30-13.10	7.05	80	6.15	0.27	5.80-6.60	4.35	80	18.39	0.28	17.10-19.70	1.52
SG: International Sysmex Instruments	27	6.21	0.13	5.80-6.60	2.10	27	18.47	0.42	17.20-19.80	2.27	27	12.04	0.19	11.20-12.80	1.58	27	6.20	0.12	5.80-6.60	1.94	27	18.38	0.34	17.10-19.70	1.85
Lymphocyte	Labs	Mean	SD	Range	CV %	Labs	Mean	SD	Range	CV %	Labs	Mean	SD	Range	CV %	Labs	Mean	SD	Range	CV %	Labs	Mean	SD	Range	CV %
PG: Sysmex K-800,1000,4500,KX-21	16	11.26	1.25	7.60-15.10	11.06	17	60.86	0.63	59.00-62.80	1.03	16	29.21	0.82	26.70-31.70	2.81	14	10.86	0.68	8.90-12.90	6.24	16	60.91	0.80	58.50-63.30	1.31
PG: Sysmex pocH-100i	22	9.93	0.92	7.10-12.70	9.29	22	61.50	0.92	58.70-64.30	1.50	22	28.31	0.94	25.50-31.10	3.32	22	10.18	1.07	7.00-13.40	10.49	22	61.33	0.63	59.40-63.20	1.03
SG: Select Sysmex Instruments (3-part Diff)	114	10.80	1.26	7.00-14.60	11.67	125	61.33	1.98	55.40-67.20	3.23	119	29.30	1.07	26.10-32.50	3.65	112	10.94	1.17	7.40-14.40	10.73	123	61.24	0.74	59.00-63.40	1.21
PG: Sysmex XP-300	76	10.96	1.25	7.30-14.80	11.36	79	61.11	0.67	59.10-63.10	1.10	79	29.62	0.94	26.80-32.40	3.18	76	11.18	1.18	7.70-14.70	10.54	79	61.21	0.69	59.10-63.30	1.13
SG: International Sysmex Instruments	26	9.81	0.92	7.00-12.60	9.39	26	61.41	0.87	58.80-64.00	1.42	26	28.34	0.88	25.70-30.90	3.11	26	10.16	1.03	7.10-13.30	10.10	26	61.33	0.62	59.40-63.20	1.01
Hematocrit	Labs	Mean	SD	Range	CV %	Labs	Mean	SD	Range	CV %	Labs	Mean	SD	Range	CV %	Labs	Mean	SD	Range	CV %	Labs	Mean	SD	Range	CV %
PG: Sysmex K-800,1000,4500,KX-21	16	16.98	0.34	16.00-18.00	2.00	17	50.38	1.29	47.40-53.40	2.56	16	32.56	0.58	30.60-34.60	1.78	16	17.01	0.35	16.00-18.00	2.06	17	49.07	3.58	46.20-52.00	7.29
PG: Sysmex pocH-100i	23	17.88	0.53	16.80-19.00	2.96	23	52.08	1.69	49.00-55.20	3.24	23	34.04	1.08	32.00-36.00	3.18	23	17.96	0.60	16.90-19.10	3.33	22	51.73	1.62	48.60-54.80	3.13
SG: Select Sysmex Instruments (3-part Diff)	119	17.38	0.91	16.40-18.40	5.23	127	50.90	2.10	47.80-54.00	4.13	126	33.25	1.87	31.20-35.20	5.63	119	17.35	0.90	16.30-18.30	5.14	124	50.41	1.85	47.40-53.40	3.67
PG: Sysmex XP-300	80	17.32	1.02	16.30-18.30	5.90	80	50.20	1.00	47.20-53.20	1.99	80	32.99	2.11	31.00-35.00	6.39	80	17.24	0.96	16.20-18.20	5.58	79	50.13	0.75	47.10-53.10	1.50
SG: International Sysmex Instruments	27	17.91	0.53	16.80-19.00	2.96	27	52.09	1.60	49.00-55.20	3.07	27	34.05	1.02	32.10-36.10	2.99	27	17.93	0.59	16.80-19.00	3.30	26	51.73	1.54	48.60-54.80	2.98
Mono/Mixed	Labs	Mean	SD	Range	CV %	Labs	Mean	SD	Range	CV %	Labs	Mean	SD	Range	CV %	Labs	Mean	SD	Range	CV %	Labs	Mean	SD	Range	CV %
PG: Sysmex K-800,1000,4500,KX-21	17	17.83	1.43	13.50-22.10	8.03	17	15.33	1.05	12.20-18.50	6.86	17	15.59	1.10	12.30-18.90	7.05	16	18.06	1.22	14.40-21.80	6.74	17	15.47	0.73	13.30-17.70	4.71
PG: Sysmex pocH-100i	22	16.75	1.49	12.20-21.20	8.92	22	13.70	0.90	11.00-16.40	6.57	22	13.79	0.71	11.70-15.90	5.14	22	16.07	1.12	12.70-19.50	6.96	22	14.02	0.72	11.80-16.20	5.14
SG: Select Sysmex Instruments (3-part Diff)	118	17.76	1.60	13.00-22.60	8.93	118	14.95	1.01	11.90-17.90	6.78	116	15.08	1.06	11.90-18.30	7.02	115	17.85	1.79	12.50-23.30	10.00	118	14.97	0.86	12.40-17.60	5.73
PG: Sysmex XP-300	79	18.03	1.56	13.30-22.70	8.67	79	15.21	0.73	13.00-17.40	4.80	77	15.33	0.84	12.80-17.80	5.49	77	18.32	1.73	13.10-23.50	9.45	79	15.12	0.72	12.90-17.30	4.77
SG: International Sysmex Instruments	26	16.77	1.37	12.70-20.90	8.15	26	13.70	0.91	11.00-16.40	6.64	26	13.83	0.68	11.80-15.80	4.86	26	16.02	1.10	12.70-19.30	6.88	26	14.06	0.68	12.10-16.10	4.82
MCV	Labs	Mean	SD	Range	CV %	Labs	Mean	SD	Range	CV %	Labs	Mean	SD	Range	CV %	Labs	Mean	SD	Range	CV %	Labs	Mean	SD	Range	CV %
PG: Sysmex K-800,1000,4500,KX-21	17	71.48	3.14	62.10-80.90	4.39	17	84.02	1.83	78.50-89.50	2.18	17	77.51	1.47	73.10-81.90	1.90	17	71.95	4.66	57.90-85.90	6.48	17	83.64	1.43	79.30-87.90	1.71
PG: Sysmex pocH-100i	18	73.92	1.45	69.60-78.20	1.95	18	86.10	1.72	80.90-91.30	2.00	18	80.14	1.57	75.40-84.80	1.96	18	73.87	1.47	69.50-78.30	1.99	18	85.87	1.53	81.30-90.50	1.78
SG: Select Sysmex Instruments (3-part Diff)	114	71.99	1.71	66.90-77.10	2.36	121	84.77	1.45	80.50-89.10	1.71	121	78.41	1.48	74.00-82.80	1.89	115	72.29	2.92	63.50-81.10	4.04	120	84.63	1.41	80.40-88.80	1.67
PG: Sysmex XP-300	79	71.66	0.84	69.20-74.20	1.17	79	84.45	0.91	81.70-87.10	1.08	79	77.99	0.88	75.40-80.60	1.13	79	71.82	2.00	65.80-77.80	2.79	79	84.37	0.91	81.70-87.10	1.08
SG: International Sysmex Instruments	22	74.00	1.36	69.90-78.10	1.84	22	86.16	1.59	81.40-91.00	1.84	22	80.23	1.45	75.90-84.50	1.81	22	73.90	1.36	69.80-78.00	1.84	22	86.00	1.42	81.70-90.30	1.65
Platelet Count	Labs	Mean	SD	Range	CV %	Labs	Mean	SD	Range	CV %	Labs	Mean	SD	Range	CV %	Labs	Mean	SD	Range	CV %	Labs	Mean	SD	Range	CV %
PG: Sysmex K-800,1000,4500,KX-21	16	60.38	4.24	45.00-75.00	7.00	17	401.29	13.05	301.00-501.00	3.27	17	186.94	32.12	140.00-234.00	17.17	16	62.50	4.69	47.00-79.00	7.46	17	389.18	52.25	292.00-486.00	13.44
PG: Sysmex pocH-100i	23	59.70	4.94	45.00-75.00	8.17	23	379.61	45.06	285.00-475.00	11.87	23	189.09	8.48	142.00-236.00	4.50	22	61.59	4.99	47.00-78.00	8.06	23	392.30	96.06	294.00-490.00	24.52
SG: Select Sysmex Instruments (3-part Diff)	117	60.21	4.02	45.00-75.00	6.67	124	394.54	25.25	296.00-494.00	6.38	125	192.95	14.70	145.00-241.00	7.62	115	61.03	4.60	46.00-76.00	7.54	123	394.20	46.76	296.00-493.00	11.88
PG: Sysmex XP-300	78	60.32	3.71	45.00-75.00	6.17	77	399.08	15.17	299.00-499.00	3.81	78	195.03	9.38	146.00-244.00	4.82	78	61.31	7.86	46.00-76.00	12.95	77	397.17	14.79	298.00-496.00	3.73
SG: International Sysmex Instruments	27</																								



PG: Beckman AU	7	64.50	5.56	48.00-82.00	8.62																				
PG: Siemens DCA Vantage	18	70.74	4.96	56.00-86.00	7.04																				
SG: Select Siemens (Bayer) DCA Instruments	20	71.01	4.76	57.00-85.00	6.76																				
PG: Afinion ACR	10	67.72	3.51	58.00-79.00	5.15																				
SG: Select Gen Chem Instruments	11	66.68	5.94	49.00-85.00	8.81																				
SG: Select Siemens Clinitek/Visual	33	133.33	47.87	0.00-277.00	36.02																				
AG: All Group	77	77.28	15.71	30.00-124.00	20.39																				
Urine Microalbumin	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%
PG: Siemens Clinitek 50	6	103.33	36.15	0.00-211.00	35.05																				
PG: Siemens Clinitek Status	27	90.37	25.34	14.00-166.00	28.11																				
PG: Siemens Dimension/AR/ES/RxL/Xpand/EXL	8	115.78	13.56	75.00-157.00	11.72																				
PG: Beckman AU	7	116.10	7.34	94.00-138.00	6.29																				
SG: Select Beckman Coulter Instruments (Chem)	9	117.97	7.69	95.00-141.00	6.53																				
PG: Siemens DCA Vantage	18	113.39	7.88	89.00-137.00	6.99																				
SG: Select Siemens (Bayer) DCA Instruments	20	114.04	7.74	91.00-137.00	6.75																				
PG: Afinion ACR	10	122.95	10.52	92.00-155.00	8.54																				
SG: Select Gen Chem Instruments	14	112.79	11.06	80.00-146.00	9.82																				
SG: Select Siemens Clinitek/Visual	36	96.39	44.92	0.00-231.00	46.77																				
AG: All Group	93	103.41	26.70	23.00-183.00	25.92																				
728 - Urine Chemistry	URC-5					URC-6																			
Creatinine - Urine	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%
AG: All Group	5	190.40	15.53	143.80-237.00	8.16	5	118.13	12.48	80.70-155.50	10.57															
Total Protein - Urine	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%
AG: All Group	5	109.42	0.92	106.60-112.20	0.84	5	52.60	2.91	43.90-61.30	5.53															
601 - Waived/PPM Tests I;609 - Waived/PPM Tests I;664 - Whole Blood Glucose;669 - Whole Blood Glucose	WBG-11					WBG-12					WBG-13					WBG-14					WBG-15				
Glucose Wh Blood (Waived)	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%
PG: Abbott Freestyle Lite	11	278.82	13.41	239.00-319.00	4.80																				
PG: Bayer Ascensia Contour	137	234.76	41.35	111.00-359.00	17.62	15	340.80	32.41	244.00-438.00	9.50															
PG: Bayer Ascensia Elite XL	10	236.80	30.40	146.00-328.00	12.83																				
PG: Lifescan One Touch Ultra	35	294.97	51.15	141.00-449.00	17.36																				
PG: Lifescan Sure Step, others	5	305.20	9.18	277.00-333.00	3.02																				
PG: McKesson/SunMark True Track	5	574.60	26.72	495.00-655.00	4.64	5	584.00	0.00	584.00-584.00	0.00															
PG: Roche Accu-Chek Advantage/Comp/GTS	23	278.65	37.19	167.00-391.00	13.33																				
PG: Roche AccuChek Aviva	22	294.64	9.61	266.00-324.00	3.25																				
PG: Roche Accu-Chek, others	17	290.18	16.14	242.00-338.00	5.55																				
SG: Select Lifescan One Touch Instruments	46	290.85	50.15	140.00-442.00	17.25	5	400.00	55.76	233.00-567.00	13.95															
SG: Select Lifescan SureStep Instruments	6	297.67	20.20	237.00-359.00	6.78																				
SG: Select Abbott Freestyle Instruments	14	276.29	13.50	236.00-317.00	4.89																				
PG: Roche AccuChek Inform	66	298.00	12.74	260.00-336.00	4.26	35	389.40	22.30	322.00-456.00	5.73															
PG: Consult Diag Quintet Blood Glucose	22	297.05	24.14	225.00-369.00	8.11	7	457.00	58.44	282.00-632.00	12.78															
PG: Medisense Precision Xceed Pro	26	274.88	29.47	187.00-364.00	10.73																				
SG: Select Medisense Instruments	32	276.63	28.70	191.00-363.00	10.36																				
SG: Select Roche WBG Meters	278	293.85	22.22	227.00-361.00	7.55	51	385.51	24.27	313.00-459.00	6.30	6	60.83	2.23	54.00-68.00	3.61	6	100.00	1.90	94.00-106.00	1.90	6	294.33	9.79	265.00-323.00	3.33
SG: Select Bayer WBG Meters	167	235.17	40.13	115.00-355.00	17.06	21	341.81	29.56	253.00-431.00	8.65	5	39.50	1.00	37.00-43.00	2.50	5	78.00	12.33	41.00-115.00	15.77	5	243.00	34.52	140.00-347.00	14.20
SG: Select True Track WBG Instruments	7	575.86	22.82	508.00-644.00	3.96																				
PG: Nova Biomedical StatStrip Glucose meter	6	251.67	15.11	207.00-297.00	5.99																				
SG: Select Nova Biomedical StatStrip Instruments	9	250.33	12.18	213.00-287.00	4.88																				
PG: Bayer Contour Next	13	236.54	26.01	159.00-315.00	10.97	5	337.80	18.94	281.00-395.00	5.59															
PG: Roche AccuChek Performa	145	295.15	22.93	226.00-364.00	7.76	7	372.29	19.65	313.00-431.00	5.30															
SG: Select TrueResult Meters	7	278.14	21.68	213.00-343.00	7.81																				
PG: True Metrix Blood Glucose Meter	10	278.60	19.32	221.00-337.00	6.92	5	598.00	0.00	598.00-598.00	0.00															
PG: McKesson True Metrix Pro Blood Glucose Meter	30	275.53	33.49	176.00-377.00	12.14	5	559.00	0.00	559.00-559.00	0.00															
SG: Select ARKRAY Whole Blood Glucose	5	312.00	56.30	143.00-481.00	18.04																				
PG: EvenCare G3 Blood Glucose System	5	278.80	10.83	247.00-311.00	3.87																				
SG: Select EvenCare Whole Blood Glucose	7	284.86	13.66	244.00-326.00	4.81																				
AG: All Group	610	277.11	30.48	186.00-369.00	11.01	80	377.05	27.93	293.00-461.00	7.40	22	48.45	8.54	23.00-74.00	17.71	23	87.74	11.09	55.00-121.00	12.61	23	277.52	31.75	183.00-373.00	11.40



789 - Bacterial Antigen Detection	Labs	Response	Number	Percent
		Strep B Ag Detect		
AD-11		* Denotes Correct Response		
AG: All Methods	2			
		*Negative	2	100.00%
Cepheid GeneXpert	1	*Negative	1	100.00%
BD Max	1	*Negative	1	100.00%
AD-12		* Denotes Correct Response		
AG: All Methods	2			
		*Negative	2	100.00%
Cepheid GeneXpert	1	*Negative	1	100.00%
BD Max	1	*Negative	1	100.00%
AD-13		* Denotes Correct Response		
AG: All Methods	2			
		*Positive	2	100.00%
Cepheid GeneXpert	1	*Positive	1	100.00%
BD Max	1	*Positive	1	100.00%
AD-14		* Denotes Correct Response		
AG: All Methods	1			
		*Negative	1	100.00%
Cepheid GeneXpert	1	*Negative	1	100.00%
AD-15		* Denotes Correct Response		
AG: All Methods	1			
		*Positive	1	100.00%
Cepheid GeneXpert	1	*Positive	1	100.00%
		C. difficile Antigen Detection		
AD-11		* Denotes Correct Response		
AG: All Methods	1			
		*Positive	1	100.00%
Meridian illumigene C. diff	1	*Positive	1	100.00%
AD-12		* Denotes Correct Response		
AG: All Methods	1			
		*Negative	1	100.00%
Meridian illumigene C. diff	1	*Negative	1	100.00%
AD-13		* Denotes Correct Response		
AG: All Methods	1			
		*Negative	1	100.00%
Meridian illumigene C. diff	1	*Negative	1	100.00%
AD-14		* Denotes Correct Response		
AG: All Methods	1			
		*Positive	1	100.00%
Meridian illumigene C. diff	1	*Positive	1	100.00%
AD-15		* Denotes Correct Response		
AG: All Methods	1			
		*Negative	1	100.00%
Meridian illumigene C. diff	1	*Negative	1	100.00%
765 - Allergen Testing	Labs	Response	Number	Percent
		Ash - White Alrqn		
AL-3		* Denotes Correct Response		
AG: All Methods	3			
		*1	1	33.33%
		*2	2	66.67%
Hitachi CLA-1	3	*1	1	33.33%
Hitachi CLA-1		*2	2	66.67%
		Birch/Alder Mix Alrqn		
AL-3		* Denotes Correct Response		
AG: All Methods	4			
		*3	3	75.00%
		*4	1	25.00%
Hitachi CLA-1	4	*3	3	75.00%
Hitachi CLA-1		*4	1	25.00%
		Cedar - Mountain Alrqn		
AL-3		* Denotes Correct Response		
AG: All Methods	4			
		*2	1	25.00%
		*3	3	75.00%
Hitachi CLA-1	4	*2	1	25.00%



Hitachi CLA-1		*3	3	75.00%
		Cottonwood - Eastern Alrgn		
AL-3		* Denotes Correct Response		
AG: All Methods	4			
		*1	1	25.00%
		*2	3	75.00%
Hitachi CLA-1	4	*1	1	25.00%
Hitachi CLA-1		*2	3	75.00%
		Elm - White Alrgn		
AL-3		* Denotes Correct Response		
AG: All Methods	4			
		*2	4	100.00%
Hitachi CLA-1	4	*2	4	100.00%
		Mulberry Mix Alrgn		
AL-3		* Denotes Correct Response		
AG: All Methods	3			
		*2	3	100.00%
Hitachi CLA-1	3	*2	3	100.00%
		Pine Mix Alrgn		
AL-3		* Denotes Correct Response		
AG: All Methods	2			
		*1	1	50.00%
		*2	1	50.00%
Hitachi CLA-1	2	*1	1	50.00%
Hitachi CLA-1		*2	1	50.00%
		Privet Alrgn		
AL-3		* Denotes Correct Response		
AG: All Methods	2			
		*1	1	50.00%
		*2	1	50.00%
Hitachi CLA-1	2	*1	1	50.00%
Hitachi CLA-1		*2	1	50.00%
		Sycamore - American Alrgn		
AL-3		* Denotes Correct Response		
AG: All Methods	2			
		*2	2	100.00%
Hitachi CLA-1	2	*2	2	100.00%
		Cocklebur Alrgn		
AL-3		* Denotes Correct Response		
AG: All Methods	2			
		*1	1	50.00%
		*2	1	50.00%
Hitachi CLA-1	2	*1	1	50.00%
Hitachi CLA-1		*2	1	50.00%
		Pigweed Alrgn		
AL-3		* Denotes Correct Response		
AG: All Methods	4			
		*2	2	50.00%
		*3	2	50.00%
Hitachi CLA-1	4	*2	2	50.00%
Hitachi CLA-1		*3	2	50.00%
		Ragweed - Short Alrgn		
AL-3		* Denotes Correct Response		
AG: All Methods	5			
		*2	1	20.00%
		*3	1	20.00%
		*4	3	60.00%
Hitachi CLA-1	5	*2	1	20.00%
Hitachi CLA-1		*3	1	20.00%
Hitachi CLA-1		*4	3	60.00%
		Sheep Sorrel Alrgn		
AL-3		* Denotes Correct Response		
AG: All Methods	5			
		*2	3	60.00%
		*3	1	20.00%
		*4	1	20.00%
Hitachi CLA-1	5	*2	3	60.00%



Hitachi CLA-1		*3	1	20.00%
Hitachi CLA-1		*4	1	20.00%
		Mite D Farinae Alrgn		
AL-3		* Denotes Correct Response		
AG: All Methods	5			
		*4	5	100.00%
Hitachi CLA-1	5	*4	5	100.00%
		Tomato Alrgn		
AL-3		* Denotes Correct Response		
AG: All Methods	3			
		*2	3	100.00%
Hitachi CLA-1	3	*2	3	100.00%
		Wheat Alrgn		
AL-3		* Denotes Correct Response		
AG: All Methods	4			
		*3	4	100.00%
Hitachi CLA-1	4	*3	4	100.00%
		Tuna Alrgn		
AL-3		* Denotes Correct Response		
AG: All Methods	2			
		*2	2	100.00%
Hitachi CLA-1	2	*2	2	100.00%
		Shellfish Mix Alrgn		
AL-3		* Denotes Correct Response		
AG: All Methods	3			
		*4	3	100.00%
Hitachi CLA-1	3	*4	3	100.00%
		Rice Alrgn		
AL-3		* Denotes Correct Response		
AG: All Methods	4			
		*3	4	100.00%
Hitachi CLA-1	4	*3	4	100.00%
		Potato Alrgn		
AL-3		* Denotes Correct Response		
AG: All Methods	2			
		*3	1	50.00%
		*4	1	50.00%
Hitachi CLA-1	2	*3	1	50.00%
Hitachi CLA-1		*4	1	50.00%
		Orange Alrgn		
AL-3		* Denotes Correct Response		
AG: All Methods	3			
		*3	3	100.00%
Hitachi CLA-1	3	*3	3	100.00%
		Onion Mix Alrgn		
AL-3		* Denotes Correct Response		
AG: All Methods	3			
		*2	2	66.67%
		*3	1	33.33%
Hitachi CLA-1	3	*2	2	66.67%
Hitachi CLA-1		*3	1	33.33%
		Oat Alrgn		
AL-3		* Denotes Correct Response		
AG: All Methods	2			
		*3	1	50.00%
		*4	1	50.00%
Hitachi CLA-1	2	*3	1	50.00%
Hitachi CLA-1		*4	1	50.00%
		Milk Alrgn		
AL-3		* Denotes Correct Response		
AG: All Methods	4			
		*4	4	100.00%
Hitachi CLA-1	4	*4	4	100.00%
		Garlic Alrgn		
AL-3		* Denotes Correct Response		
AG: All Methods	2			
		*3	1	50.00%



		*4	1	50.00%
Hitachi CLA-1	2	*3	1	50.00%
Hitachi CLA-1		*4	1	50.00%
		Egg - whole Alrgn		
AL-3		* Denotes Correct Response		
AG: All Methods	4			
		*3	1	25.00%
		*4	3	75.00%
Hitachi CLA-1	4	*3	1	25.00%
Hitachi CLA-1		*4	3	75.00%
		Corn Alrgn		
AL-3		* Denotes Correct Response		
AG: All Methods	3			
		*1	1	33.33%
		*2	1	33.33%
		*3	1	33.33%
Hitachi CLA-1	3	*1	1	33.33%
Hitachi CLA-1		*2	1	33.33%
Hitachi CLA-1		*3	1	33.33%
		Chocolate Alrgn		
AL-3		* Denotes Correct Response		
AG: All Methods	2			
		*2	1	50.00%
		*3	1	50.00%
Hitachi CLA-1	2	*2	1	50.00%
Hitachi CLA-1		*3	1	50.00%
		Beef Alrgn		
AL-3		* Denotes Correct Response		
AG: All Methods	2			
		*3	2	100.00%
Hitachi CLA-1	2	*3	2	100.00%
		Bean - White Alrgn		
AL-3		* Denotes Correct Response		
AG: All Methods	2			
		*4	2	100.00%
Hitachi CLA-1	2	*4	2	100.00%
		Barley Alrgn		
AL-3		* Denotes Correct Response		
AG: All Methods	2			
		*3	2	100.00%
Hitachi CLA-1	2	*3	2	100.00%
		Apple Alrgn		
AL-3		* Denotes Correct Response		
AG: All Methods	2			
		*2	1	50.00%
		*3	1	50.00%
Hitachi CLA-1	2	*2	1	50.00%
Hitachi CLA-1		*3	1	50.00%
		Almond Alrgn		
AL-3		* Denotes Correct Response		
AG: All Methods	2			
		*3	1	50.00%
		*4	1	50.00%
Hitachi CLA-1	2	*3	1	50.00%
Hitachi CLA-1		*4	1	50.00%
		Alternaria alternata Alrgn		
AL-3		* Denotes Correct Response		
AG: All Methods	5			
		*2	4	80.00%
		*4	1	20.00%
Hitachi CLA-1	5	*2	4	80.00%
Hitachi CLA-1		*4	1	20.00%
		Bermuda Grass Alrgn		
AL-3		* Denotes Correct Response		
AG: All Methods	4			
		*3	2	50.00%
		*4	2	50.00%



Hitachi CLA-1	4	*3	2	50.00%
Hitachi CLA-1		*4	2	50.00%
		Cat Epithelium Alrgn		
AL-3		* Denotes Correct Response		
AG: All Methods	2			
		*4	2	100.00%
Hitachi CLA-1	2	*4	2	100.00%
		Dog Epithelium Alrgn		
AL-3		* Denotes Correct Response		
AG: All Methods	2			
		*4	2	100.00%
Hitachi CLA-1	2	*4	2	100.00%
		Dog Dander Alrgn		
AL-3		* Denotes Correct Response		
AG: All Methods	3			
		*3	1	33.33%
		*4	2	66.67%
Hitachi CLA-1	3	*3	1	33.33%
Hitachi CLA-1		*4	2	66.67%
		House Dust Mite Alrgn		
AL-3		* Denotes Correct Response		
AG: All Methods	4			
		*2	1	25.00%
		*3	2	50.00%
		*4	1	25.00%
Hitachi CLA-1	4	*2	1	25.00%
Hitachi CLA-1		*3	2	50.00%
Hitachi CLA-1		*4	1	25.00%
		Timothy Grass Alrgn		
AL-3		* Denotes Correct Response		
AG: All Methods	4			
		*4	4	100.00%
Hitachi CLA-1	4	*4	4	100.00%
		Peanut (Arachis hyogaea) Alrgn		
AL-3		* Denotes Correct Response		
AG: All Methods	4			
		*4	4	100.00%
Hitachi CLA-1	4	*4	4	100.00%
		Cat Dander Alrgn		
AL-3		* Denotes Correct Response		
AG: All Methods	3			
		*4	3	100.00%
Hitachi CLA-1	3	*4	3	100.00%
		Soybean Alrgn		
AL-3		* Denotes Correct Response		
AG: All Methods	4			
		*4	4	100.00%
Hitachi CLA-1	4	*4	4	100.00%
		Cockroach Mix Alrgn		
AL-3		* Denotes Correct Response		
AG: All Methods	5			
		*2	3	60.00%
		*3	1	20.00%
		*4	1	20.00%
Hitachi CLA-1	5	*2	3	60.00%
Hitachi CLA-1		*3	1	20.00%
Hitachi CLA-1		*4	1	20.00%
		Oak, White Alrgn		
AL-3		* Denotes Correct Response		
AG: All Methods	2			
		*3	1	50.00%
		*4	1	50.00%
Hitachi CLA-1	2	*3	1	50.00%
Hitachi CLA-1		*4	1	50.00%
		Walnut/Hickory/Pecan Mix Alrgn		
AL-3		* Denotes Correct Response		
AG: All Methods	4			



		*2	3	75.00%
		*3	1	25.00%
Hitachi CLA-1	4	*2	3	75.00%
Hitachi CLA-1		*3	1	25.00%
683 - Antistreptolysin O (ASO)	Labs	Response	Number	Percent
		ASO (Qual) (International)		
AS-11		* Denotes Correct Response		
AG: All Methods	1			
		*Positive	1	100.00%
Wampole ColorCard ASO	1	*Positive	1	100.00%
AS-12		* Denotes Correct Response		
AG: All Methods	1			
		*Negative	1	100.00%
Wampole ColorCard ASO	1	*Negative	1	100.00%
AS-13		* Denotes Correct Response		
AG: All Methods	1			
		*Positive	1	100.00%
Wampole ColorCard ASO	1	*Positive	1	100.00%
AS-14		* Denotes Correct Response		
AG: All Methods	1			
		*Negative	1	100.00%
Wampole ColorCard ASO	1	*Negative	1	100.00%
AS-15		* Denotes Correct Response		
AG: All Methods	1			
		*Positive	1	100.00%
Wampole ColorCard ASO	1	*Positive	1	100.00%
975 - Bacterial Vaginosis	Labs	Response	Number	Percent
		Bacterial Vaginosa (Sialidase)		
BV-5		* Denotes Correct Response		
AG: All Methods	102			
		*Negative	94	92.16%
		Positive	8	7.84%
Sekisui Diagnostics OSOM BVBlue	102	*Negative	94	92.16%
Sekisui Diagnostics OSOM BVBlue		Positive	8	7.84%
BV-6		* Denotes Correct Response		
AG: All Methods	102			
		*Positive	90	88.24%
		Negative	12	11.76%
Sekisui Diagnostics OSOM BVBlue	102	*Positive	90	88.24%
Sekisui Diagnostics OSOM BVBlue		Negative	12	11.76%
661 - Urine Colony Count;784 - Presump ID/Colony Count - Urine;785 - Presump ID/Colony Count-Urine & Strep Combo;775 - Urine Culture;778 - Urine & Throat Cultures;779 - Urine & Genital Cultures;781 - Urine, Throat & Genital Cultures;782 - Complete Bacteri	Labs	Response	Number	Percent
		Colony Count		
CC-11		* Denotes Correct Response		
AG: All Methods	30			
		*10,000 - 100,000	4	13.33%
		*Greater than 100,000	24	80.00%
		No growth	2	6.67%
Calibrated Loop - Manual	9	*10,000 - 100,000	1	11.11%
Uri-Chek (Troy Bio)	6	*10,000 - 100,000	2	33.33%
Uricult (Orion)	13	*10,000 - 100,000	1	7.69%
Calibrated Loop - Manual	9	*Greater than 100,000	8	88.89%
Uri-Chek (Troy Bio)	6	*Greater than 100,000	4	66.67%
Uricult (Orion)	13	*Greater than 100,000	10	76.92%
Uri-Three (Culture Kits	1	*Greater than 100,000	1	100.00%
Uricult (Lifesign)	1	*Greater than 100,000	1	100.00%
Uricult (Orion)	13	No growth	2	15.38%
784 - Presump ID/Colony Count - Urine;785 - Presump ID/Colony Count-Urine & Strep Combo	Labs	Response	Number	Percent
		Colony Count		
CC-11		* Denotes Correct Response		
AG: All Methods	30			
		*10,000 - 100,000	4	13.33%
		*Greater than 100,000	24	80.00%
		No growth	2	6.67%
Calibrated Loop - Manual	9	*10,000 - 100,000	1	11.11%



Uri-Chek (Troy Bio)	6	*10,000 - 100,000	2	33.33%
Uricult (Orion)	13	*10,000 - 100,000	1	7.69%
Calibrated Loop - Manual	9	*Greater than 100,000	8	88.89%
Uri-Chek (Troy Bio)	6	*Greater than 100,000	4	66.67%
Uricult (Orion)	13	*Greater than 100,000	10	76.92%
Uri-Three (Culture Kits)	1	*Greater than 100,000	1	100.00%
Uricult (Lifesign)	1	*Greater than 100,000	1	100.00%
Uricult (Orion)	13	No growth	2	15.38%
CC-12		* Denotes Correct Response		
AG: All Methods	5			
		*Greater than 100,000	5	100.00%
Calibrated Loop - Manual	4	*Greater than 100,000	4	100.00%
Uri-Chek (Troy Bio)	1	*Greater than 100,000	1	100.00%
CC-13		* Denotes Correct Response		
AG: All Methods	5			
		*Greater than 100,000	5	100.00%
Calibrated Loop - Manual	4	*Greater than 100,000	4	100.00%
Uri-Chek (Troy Bio)	1	*Greater than 100,000	1	100.00%
CC-14		* Denotes Correct Response		
AG: All Methods	5			
		*Greater than 100,000	5	100.00%
Calibrated Loop - Manual	4	*Greater than 100,000	4	100.00%
Uri-Chek (Troy Bio)	1	*Greater than 100,000	1	100.00%
CC-15		* Denotes Correct Response		
AG: All Methods	5			
		*10,000 - 100,000	4	80.00%
		*Greater than 100,000	1	20.00%
Calibrated Loop - Manual	4	*10,000 - 100,000	3	75.00%
Uri-Chek (Troy Bio)	1	*10,000 - 100,000	1	100.00%
Calibrated Loop - Manual	4	*Greater than 100,000	1	25.00%
		Urine Presump Id		
CC-11		* Denotes Correct Response		
AG: All Methods	2			
		*Gram Negative Organsim	1	50.00%
		*Presumptive Proteus sp.	1	50.00%
Calibrated Loop - Manual	2	*Gram Negative Organsim	1	50.00%
Calibrated Loop - Manual		*Presumptive Proteus sp.	1	50.00%
CC-12		* Denotes Correct Response		
AG: All Methods	2			
		*Gram Negative Organsim	1	50.00%
		*Presumptive Klebsiella sp.	1	50.00%
Calibrated Loop - Manual	2	*Gram Negative Organsim	1	50.00%
Calibrated Loop - Manual		*Presumptive Klebsiella sp.	1	50.00%
CC-13		* Denotes Correct Response		
AG: All Methods	2			
		*Gram Positive Organism	1	50.00%
		Presumptive Proteus sp.	1	50.00%
Calibrated Loop - Manual	2	*Gram Positive Organism	1	50.00%
Calibrated Loop - Manual		Presumptive Proteus sp.	1	50.00%
CC-14		* Denotes Correct Response		
AG: All Methods	2			
		*Yeast	1	50.00%
		Gram Positive Organism	1	50.00%
Calibrated Loop - Manual	2	*Yeast	1	50.00%
Calibrated Loop - Manual		Gram Positive Organism	1	50.00%
CC-15		* Denotes Correct Response		
AG: All Methods	2			
		*Gram Negative Organsim	1	50.00%
		*Presumptive E. Coli	1	50.00%
Calibrated Loop - Manual	2	*Gram Negative Organsim	1	50.00%
Calibrated Loop - Manual		*Presumptive E. Coli	1	50.00%
		Susceptibility Amikacin		
CC-11		* Denotes Correct Response		
AG: All Methods	1			
		*Sensitive	1	100.00%
Susceptibility	1	*Sensitive	1	100.00%
		Susceptibility Amox/Clavulanate		
CC-11		* Denotes Correct Response		



AG: All Methods	4			
		*Sensitive	4	100.00%
Susceptibility	4	*Sensitive	4	100.00%
		Susceptibility Ampicillin		
CC-11		* Denotes Correct Response		
AG: All Methods	5			
		*Resistant	5	100.00%
Susceptibility	5	*Resistant	5	100.00%
		Susceptibility Ampicillin/Sulbactam		
CC-11		* Denotes Correct Response		
AG: All Methods	2			
		*Intermediate	2	100.00%
Susceptibility	2	*Intermediate	2	100.00%
		Susceptibility Cefazolin		
CC-11		* Denotes Correct Response		
AG: All Methods	2			
		*Resistant	2	100.00%
Susceptibility	2	*Resistant	2	100.00%
		Susceptibility Cefepime		
CC-11		* Denotes Correct Response		
AG: All Methods	2			
		*Sensitive	2	100.00%
Susceptibility	2	*Sensitive	2	100.00%
		Susceptibility Cefpodoxime		
CC-11		* Denotes Correct Response		
AG: All Methods	1			
		*Sensitive	1	100.00%
Susceptibility	1	*Sensitive	1	100.00%
		Susceptibility Cefazidime		
CC-11		* Denotes Correct Response		
AG: All Methods	2			
		*Sensitive	2	100.00%
Susceptibility	2	*Sensitive	2	100.00%
		Susceptibility Cefuroxime		
CC-11		* Denotes Correct Response		
AG: All Methods	3			
		*Resistant	3	100.00%
Susceptibility	3	*Resistant	3	100.00%
		Susceptibility Cephalothin		
CC-11		* Denotes Correct Response		
AG: All Methods	4			
		*Resistant	4	100.00%
Susceptibility	4	*Resistant	4	100.00%
		Susceptibility Ciprofloxacin		
CC-11		* Denotes Correct Response		
AG: All Methods	7			
		*Sensitive	7	100.00%
Susceptibility	7	*Sensitive	7	100.00%
		Susceptibility Gentamicin		
CC-11		* Denotes Correct Response		
AG: All Methods	5			
		*Sensitive	5	100.00%
Susceptibility	5	*Sensitive	5	100.00%
		Susceptibility Levofloxacin		
CC-11		* Denotes Correct Response		
AG: All Methods	4			
		*Sensitive	4	100.00%
Susceptibility	4	*Sensitive	4	100.00%
		Susceptibility Nitrofurantoin		
CC-11		* Denotes Correct Response		
AG: All Methods	6			
		*Intermediate	2	33.33%
		*Resistant	3	50.00%
		Sensitive	1	16.67%
Susceptibility	6	*Intermediate	2	33.33%
Susceptibility		*Resistant	3	50.00%
Susceptibility		Sensitive	1	16.67%



		Susceptibility Norfloxacin/Noroxin		
CC-11		* Denotes Correct Response		
AG: All Methods	1			
		*Sensitive	1	100.00%
Susceptibility	1	*Sensitive	1	100.00%
		Susceptibility Piperacillin		
CC-11		* Denotes Correct Response		
AG: All Methods	1			
		*Sensitive	1	100.00%
Susceptibility	1	*Sensitive	1	100.00%
		Susceptibility Piperacillin/Tazobactam		
CC-11		* Denotes Correct Response		
AG: All Methods	2			
		*Sensitive	2	100.00%
Susceptibility	2	*Sensitive	2	100.00%
		Susceptibility Sulfamethoxazole		
CC-11		* Denotes Correct Response		
AG: All Methods	1			
		*Sensitive	1	100.00%
Susceptibility	1	*Sensitive	1	100.00%
		Susceptibility Tobramycin		
CC-11		* Denotes Correct Response		
AG: All Methods	3			
		*Sensitive	3	100.00%
Susceptibility	3	*Sensitive	3	100.00%
		Susceptibility Trimethoprim		
CC-11		* Denotes Correct Response		
AG: All Methods	1			
		*Sensitive	1	100.00%
Susceptibility	1	*Sensitive	1	100.00%
		Susceptibility Trimeth./Sulfameth.		
CC-11		* Denotes Correct Response		
AG: All Methods	5			
		*Sensitive	5	100.00%
Susceptibility	5	*Sensitive	5	100.00%
		Susceptibility Ertapenem		
CC-11		* Denotes Correct Response		
AG: All Methods	1			
		*Sensitive	1	100.00%
Susceptibility	1	*Sensitive	1	100.00%
963 - C. difficile Toxin/Ag (5 specimens);964 - C. difficile Toxin/Ag (2 specimens)	Labs	Response	Number	Percent
		C. difficile Antigen Detection		
CDF-11		* Denotes Correct Response		
AG: All Methods	1			
		*Negative	1	100.00%
Cepheid GeneXpert - Xpert C. diff	1	*Negative	1	100.00%
CDF-12		* Denotes Correct Response		
AG: All Methods	1			
		*Positive	1	100.00%
Cepheid GeneXpert - Xpert C. diff	1	*Positive	1	100.00%
794 - Vaginal Pathogens II (5 specimens);795 - Vaginal Pathogens II (2 specimens);977 - Trichomonas Ag Det (Wavied)	Labs	Response	Number	Percent
		Chlamydia - Non Culture Tech		
CGT-11		* Denotes Correct Response		
AG: All Methods	20			
		*Negative	20	100.00%
BD Probe Tec ET	1	*Negative	1	100.00%
Quidel QuickVue Chlamydia	2	*Negative	2	100.00%
Cepheid GeneXpert	15	*Negative	15	100.00%
BD Viper XTR	1	*Negative	1	100.00%
BD Max	1	*Negative	1	100.00%
CGT-12		* Denotes Correct Response		
AG: All Methods	20			
		*Positive	19	95.00%
		Negative	1	5.00%
BD Probe Tec ET	1	*Positive	1	100.00%



Quidel QuickVue Chlamydia	2	*Positive	1	50.00%
Cepheid GeneXpert	15	*Positive	15	100.00%
BD Viper XTR	1	*Positive	1	100.00%
BD Max	1	*Positive	1	100.00%
Quidel QuickVue Chlamydia	2	Negative	1	50.00%
CGT-13		* Denotes Correct Response		
AG: All Methods	20			
		*Positive	19	95.00%
		Negative	1	5.00%
BD Probe Tec ET	1	*Positive	1	100.00%
Quidel QuickVue Chlamydia	2	*Positive	2	100.00%
Cepheid GeneXpert	15	*Positive	14	93.33%
BD Viper XTR	1	*Positive	1	100.00%
BD Max	1	*Positive	1	100.00%
Cepheid GeneXpert	15	Negative	1	6.67%
CGT-14		* Denotes Correct Response		
AG: All Methods	20			
		*Positive	19	95.00%
		Negative	1	5.00%
BD Probe Tec ET	1	*Positive	1	100.00%
Quidel QuickVue Chlamydia	2	*Positive	1	50.00%
Cepheid GeneXpert	15	*Positive	15	100.00%
BD Viper XTR	1	*Positive	1	100.00%
BD Max	1	*Positive	1	100.00%
Quidel QuickVue Chlamydia	2	Negative	1	50.00%
CGT-15		* Denotes Correct Response		
AG: All Methods	20			
		*Negative	20	100.00%
BD Probe Tec ET	1	*Negative	1	100.00%
Quidel QuickVue Chlamydia	2	*Negative	2	100.00%
Cepheid GeneXpert	15	*Negative	15	100.00%
BD Viper XTR	1	*Negative	1	100.00%
BD Max	1	*Negative	1	100.00%
		GC - Non Culture Tech		
CGT-11		* Denotes Correct Response		
AG: All Methods	18			
		*Positive	18	100.00%
BD Probe Tec ET	1	*Positive	1	100.00%
Cepheid GeneXpert	15	*Positive	15	100.00%
BD Viper XTR	1	*Positive	1	100.00%
BD Max	1	*Positive	1	100.00%
CGT-12		* Denotes Correct Response		
AG: All Methods	18			
		*Negative	18	100.00%
BD Probe Tec ET	1	*Negative	1	100.00%
Cepheid GeneXpert	15	*Negative	15	100.00%
BD Viper XTR	1	*Negative	1	100.00%
BD Max	1	*Negative	1	100.00%
CGT-13		* Denotes Correct Response		
AG: All Methods	18			
		*Positive	17	94.44%
		Negative	1	5.56%
BD Probe Tec ET	1	*Positive	1	100.00%
Cepheid GeneXpert	15	*Positive	14	93.33%
BD Viper XTR	1	*Positive	1	100.00%
BD Max	1	*Positive	1	100.00%
Cepheid GeneXpert	15	Negative	1	6.67%
CGT-14		* Denotes Correct Response		
AG: All Methods	18			
		*Negative	18	100.00%
BD Probe Tec ET	1	*Negative	1	100.00%
Cepheid GeneXpert	15	*Negative	15	100.00%
BD Viper XTR	1	*Negative	1	100.00%
BD Max	1	*Negative	1	100.00%
CGT-15		* Denotes Correct Response		
AG: All Methods	18			
		*Positive	18	100.00%
BD Probe Tec ET	1	*Positive	1	100.00%



Cepheid GeneXpert	15	*Positive	15	100.00%
BD Viper XTR	1	*Positive	1	100.00%
BD Max	1	*Positive	1	100.00%
		T. vaginalis - Non Culture Tech		
CGT-11		* Denotes Correct Response		
AG: All Methods	3			
		*Positive	3	100.00%
Cepheid GeneXpert	2	*Positive	2	100.00%
BD Max	1	*Positive	1	100.00%
CGT-12		* Denotes Correct Response		
AG: All Methods	3			
		*Negative	3	100.00%
Cepheid GeneXpert	2	*Negative	2	100.00%
BD Max	1	*Negative	1	100.00%
CGT-13		* Denotes Correct Response		
AG: All Methods	3			
		*Positive	3	100.00%
Cepheid GeneXpert	2	*Positive	2	100.00%
BD Max	1	*Positive	1	100.00%
CGT-14		* Denotes Correct Response		
AG: All Methods	3			
		*Negative	3	100.00%
Cepheid GeneXpert	2	*Negative	2	100.00%
BD Max	1	*Negative	1	100.00%
CGT-15		* Denotes Correct Response		
AG: All Methods	3			
		*Positive	3	100.00%
Cepheid GeneXpert	2	*Positive	2	100.00%
BD Max	1	*Positive	1	100.00%
		Trichomonas Aq Det (Waived)		
CGT-11		* Denotes Correct Response		
AG: All Methods	100			
		*Positive	93	93.00%
		Negative	7	7.00%
Sekisui Diagnostics OSOM Trichomonas	100	*Positive	93	93.00%
Sekisui Diagnostics OSOM Trichomonas		Negative	7	7.00%
CGT-12		* Denotes Correct Response		
AG: All Methods	100			
		*Negative	94	94.00%
		Positive	6	6.00%
Sekisui Diagnostics OSOM Trichomonas	100	*Negative	94	94.00%
Sekisui Diagnostics OSOM Trichomonas		Positive	6	6.00%
639 - Complete Urinalysis Package;638 - Urine Sediment;641 - Clinical Microscopy;601 - Waived/PPM Tests II;609 - Waived/PPM Tests I	Labs	Response	Number	Percent
		Urine Sediment Id		
CM-19		* Denotes Correct Response		
AG: All Methods	645			
		*Calcium oxalate crystal	634	98.29%
		Ammonium magnesium phosphate (Triple phosphate)	6	0.93%
		Calcium oxalate crystal	634	98.29%
		Calcium phosphate crystal	1	0.16%
		Cystine crystal	1	0.16%
		Hippuric acid	1	0.16%
		Mucous strands	1	0.16%
		Uric acid crystal	1	0.16%
Microscopy	645	*Calcium oxalate crystal	634	98.29%
Microscopy		Ammonium magnesium phosphate (Triple phosphate)	6	0.93%
Microscopy		Calcium oxalate crystal	634	98.29%
Microscopy		Calcium phosphate crystal	1	0.16%
Microscopy		Cystine crystal	1	0.16%
Microscopy		Hippuric acid	1	0.16%
Microscopy		Mucous strands	1	0.16%
Microscopy		Uric acid crystal	1	0.16%
CM-20		* Denotes Correct Response		
AG: All Methods	649			
		*Ident Unknown - Referred	3	0.46%
		*Mucous strands	602	92.76%



		Calcium oxalate crystal	1	0.15%
		Cholesterol crystal	1	0.15%
		Fiber (exogenous)/fecal contamination	21	3.24%
		Granular cast	1	0.15%
		Hyaline cast	11	1.69%
		Waxy cast	2	0.31%
		Yeast/fungi	7	1.08%
Microscopy	649	*Ident Unknown - Referred	3	0.46%
Microscopy		*Mucous strands	602	92.76%
Microscopy		Calcium oxalate crystal	1	0.15%
Microscopy		Cholesterol crystal	1	0.15%
Microscopy		Fiber (exogenous)/fecal contamination	21	3.24%
Microscopy		Granular cast	1	0.15%
Microscopy		Hyaline cast	11	1.69%
Microscopy		Waxy cast	2	0.31%
Microscopy		Yeast/fungi	7	1.08%
CM-21		* Denotes Correct Response		
AG: All Methods	650			
		*Trichomonas vaginalis	639	98.31%
		Ammonium biurate crystal	2	0.31%
		Red blood cell (RBC) (Erythrocyte)	2	0.31%
		Spermatozoa	1	0.15%
		Squamous epithelial cell	2	0.31%
		Starch granules	1	0.15%
		Trichomonas vaginalis	639	98.31%
		White blood cell (WBC) (Leukocyte)	3	0.46%
Microscopy	650	*Trichomonas vaginalis	639	98.31%
Microscopy		Ammonium biurate crystal	2	0.31%
Microscopy		Red blood cell (RBC) (Erythrocyte)	2	0.31%
Microscopy		Spermatozoa	1	0.15%
Microscopy		Squamous epithelial cell	2	0.31%
Microscopy		Starch granules	1	0.15%
Microscopy		Trichomonas vaginalis	639	98.31%
Microscopy		White blood cell (WBC) (Leukocyte)	3	0.46%
CM-22		* Denotes Correct Response		
AG: All Methods	646			
		*Squamous epithelial cell	625	96.75%
		Renal tubular epithelial (RTE)	3	0.46%
		Spermatozoa	1	0.15%
		Squamous epithelial cell	625	96.75%
		Transitional epithelial cell (urothelial)	12	1.86%
		Trichomonas vaginalis	2	0.31%
		White blood cell (WBC) (Leukocyte)	2	0.31%
		Yeast/fungi	1	0.15%
Microscopy	646	*Squamous epithelial cell	625	96.75%
Microscopy		Renal tubular epithelial (RTE)	3	0.46%
Microscopy		Spermatozoa	1	0.15%
Microscopy		Squamous epithelial cell	625	96.75%
Microscopy		Transitional epithelial cell (urothelial)	12	1.86%
Microscopy		Trichomonas vaginalis	2	0.31%
Microscopy		White blood cell (WBC) (Leukocyte)	2	0.31%
Microscopy		Yeast/fungi	1	0.15%
		Vaginal Wet Prep		
CM-23		* Denotes Correct Response		
AG: All Methods	897			
		*Spermatozoa present	885	98.66%
		Clue cell	3	0.33%
		Spermatozoa absent	2	0.22%
		Squamous epithelial cell	2	0.22%
		Trichomonas vaginalis	2	0.22%
		Yeast/fungi present	3	0.33%
Microscopy	897	*Spermatozoa present	885	98.66%
Microscopy		Clue cell	3	0.33%
Microscopy		Spermatozoa absent	2	0.22%
Microscopy		Squamous epithelial cell	2	0.22%
Microscopy		Trichomonas vaginalis	2	0.22%
Microscopy		Yeast/fungi present	3	0.33%
CM-24		* Denotes Correct Response		
AG: All Methods	897			



		*Bacteria present	823	91.75%
		Bacteria absent	4	0.45%
		Fiber (exogenous)/fecal contamination	1	0.11%
		Mucous strand	1	0.11%
		Pollen grain	1	0.11%
		Squamous epithelial cell	8	0.89%
		Starch granule	3	0.33%
		Trichomonas vaginalis	1	0.11%
		Uric acid crystal	1	0.11%
		Waxy cast	1	0.11%
		White Blood Cell (WBC) (Leukocyte)	2	0.22%
		Yeast/fungi absent	1	0.11%
		Yeast/fungi present	50	5.57%
Microscopy	897	*Bacteria present	823	91.75%
Microscopy		Bacteria absent	4	0.45%
Microscopy		Fiber (exogenous)/fecal contamination	1	0.11%
Microscopy		Mucous strand	1	0.11%
Microscopy		Pollen grain	1	0.11%
Microscopy		Squamous epithelial cell	8	0.89%
Microscopy		Starch granule	3	0.33%
Microscopy		Trichomonas vaginalis	1	0.11%
Microscopy		Uric acid crystal	1	0.11%
Microscopy		Waxy cast	1	0.11%
Microscopy		White Blood Cell (WBC) (Leukocyte)	2	0.22%
Microscopy		Yeast/fungi absent	1	0.11%
Microscopy		Yeast/fungi present	50	5.57%
		Skin Scraping		
CM-25		* Denotes Correct Response		
AG: All Methods	588			
		*No scabies present	4	0.68%
		*Yeast/fungi present	582	98.98%
		Yeast/fungi absent	2	0.34%
Microscopy	588	*No scabies present	4	0.68%
Microscopy		*Yeast/fungi present	582	98.98%
Microscopy		Yeast/fungi absent	2	0.34%
		Fern Test		
CM-26		* Denotes Correct Response		
AG: All Methods	384			
		*Negative Fern Test	383	99.74%
		Positive Fern Test	1	0.26%
Microscopy	384	*Negative Fern Test	383	99.74%
Microscopy		Positive Fern Test	1	0.26%
772 - Dermatophyte Screen		Response	Number	Percent
		Dermatophyte Screen		
DM-11		* Denotes Correct Response		
AG: All Methods	15			
		*Positive, Dermatophytes Present	12	80.00%
		Negative, No Dermatophytes Present	3	20.00%
Dermatophytes	15	*Positive, Dermatophytes Present	12	80.00%
Dermatophytes		Negative, No Dermatophytes Present	3	20.00%
DM-12		* Denotes Correct Response		
AG: All Methods	14			
		*Negative, No Dermatophytes Present	14	100.00%
Dermatophytes	14	*Negative, No Dermatophytes Present	14	100.00%
DM-13		* Denotes Correct Response		
AG: All Methods	15			
		*Positive, Dermatophytes Present	12	80.00%
		Negative, No Dermatophytes Present	3	20.00%
Dermatophytes	15	*Positive, Dermatophytes Present	12	80.00%
Dermatophytes		Negative, No Dermatophytes Present	3	20.00%
DM-14		* Denotes Correct Response		
AG: All Methods	14			
		*Negative, No Dermatophytes Present	13	92.86%
		Positive, Dermatophytes Present	1	7.14%
Dermatophytes	14	*Negative, No Dermatophytes Present	13	92.86%
Dermatophytes		Positive, Dermatophytes Present	1	7.14%
DM-15		* Denotes Correct Response		
AG: All Methods	14			



		*Positive, Dermatophytes Present	13	92.86%
		Negative, No Dermatophytes Present	1	7.14%
Dermatophytes	14	*Positive, Dermatophytes Present	13	92.86%
Dermatophytes		Negative, No Dermatophytes Present	1	7.14%
777 - Genital Culture;796 - Genital Culture (2 specimens)	Labs	Response	Number	Percent
		Genital Culture		
GC-11		* Denotes Correct Response		
AG: All Methods	3			
		*Negative for GC	1	33.33%
		*No growth, GC selective media	2	66.67%
Cultures	3	*Negative for GC	1	33.33%
Cultures		*No growth, GC selective media	2	66.67%
GC-12		* Denotes Correct Response		
AG: All Methods	3			
		*Presumptive positive for GC, GR(-) diplococci, oxidase positive organism	1	33.33%
		*Presumptive positive for GC, growth referred for ID/confirmation	2	66.67%
Cultures	3	*Presumptive positive for GC, GR(-) diplococci, oxidase positive organism	1	33.33%
Cultures		*Presumptive positive for GC, growth referred for ID/confirmation	2	66.67%
GC-13		* Denotes Correct Response		
AG: All Methods	3			
		*No growth, GC selective media	3	100.00%
Cultures	3	*No growth, GC selective media	3	100.00%
GC-14		* Denotes Correct Response		
AG: All Methods	5			
		*Group B Streptococcus not present	1	20.00%
		*Negative for GC	2	40.00%
		*No growth, GC selective media	2	40.00%
Cultures	5	*Group B Streptococcus not present	1	20.00%
Cultures		*Negative for GC	2	40.00%
Cultures		*No growth, GC selective media	2	40.00%
GC-15		* Denotes Correct Response		
AG: All Methods	4			
		*Positive for GC	1	25.00%
		*Presumptive positive for GC, GR(-) diplococci, oxidase positive organism	1	25.00%
		*Presumptive positive for GC, growth referred for ID/confirmation	2	50.00%
Cultures	4	*Positive for GC	1	25.00%
Cultures		*Presumptive positive for GC, GR(-) diplococci, oxidase positive organism	1	25.00%
Cultures		*Presumptive positive for GC, growth referred for ID/confirmation	2	50.00%
668 - Gram Stain	Labs	Response	Number	Percent
		Gram Stain		
GS-11		* Denotes Correct Response		
AG: All Methods	8			
		*Gram Positive	8	100.00%
Gram Stain	8	*Gram Positive	8	100.00%
GS-12		* Denotes Correct Response		
AG: All Methods	8			
		*Gram Negative	8	100.00%
Gram Stain	8	*Gram Negative	8	100.00%
GS-13		* Denotes Correct Response		
AG: All Methods	8			
		*Gram Positive	8	100.00%
Gram Stain	8	*Gram Positive	8	100.00%
GS-14		* Denotes Correct Response		
AG: All Methods	8			
		*Gram Positive	8	100.00%
Gram Stain	8	*Gram Positive	8	100.00%
GS-15		* Denotes Correct Response		
AG: All Methods	8			
		*Gram Negative	8	100.00%
Gram Stain	8	*Gram Negative	8	100.00%
		Gram Stain Morph		
GS-11		* Denotes Correct Response		



AG: All Methods	7			
		*Cocci	5	71.43%
		*Cocci in Chains	1	14.29%
		*Cocco-Bacilli	1	14.29%
Gram Stain	7	*Cocci	5	71.43%
Gram Stain		*Cocci in Chains	1	14.29%
Gram Stain		*Cocco-Bacilli	1	14.29%
GS-12		* Denotes Correct Response		
AG: All Methods	7			
		*Rods/Bacilli	7	100.00%
Gram Stain	7	*Rods/Bacilli	7	100.00%
GS-13		* Denotes Correct Response		
AG: All Methods	7			
		*Cocci	5	71.43%
		*Cocci in Pairs	2	28.57%
Gram Stain	7	*Cocci	5	71.43%
Gram Stain		*Cocci in Pairs	2	28.57%
GS-14		* Denotes Correct Response		
AG: All Methods	7			
		*Rods/Bacilli	7	100.00%
Gram Stain	7	*Rods/Bacilli	7	100.00%
GS-15		* Denotes Correct Response		
AG: All Methods	7			
		*Cocci	3	42.86%
		*Diplococci	4	57.14%
Gram Stain	7	*Cocci	3	42.86%
Gram Stain		*Diplococci	4	57.14%
678 - Serum hCG	Labs	Response	Number	Percent
		HCG - Qual Serum		
HCG-11		* Denotes Correct Response		
AG: All Methods	86			
		*Positive	86	100.00%
Beck Coulter ICON 25	7	*Positive	7	100.00%
Fisher HC Sure-Vue hCG STAT	2	*Positive	2	100.00%
Sekisui Diagnostics OSOM	2	*Positive	2	100.00%
Henry Schein One Step + (hCG)	5	*Positive	5	100.00%
Immunostics Immuno/hCG	1	*Positive	1	100.00%
Inverness Signify HCG	2	*Positive	2	100.00%
LifeSign Status hCG	1	*Positive	1	100.00%
McKesson Medi-Lab hCG - Urine	1	*Positive	1	100.00%
Quidel QuickVue + hCG	2	*Positive	2	100.00%
Quidel QuickVue One-Step hCG Combo (U/S)	30	*Positive	30	100.00%
Select Diag hCG Combo	3	*Positive	3	100.00%
SP Brand (Cardinal) HCG Combo	14	*Positive	14	100.00%
Wampole Clearview	1	*Positive	1	100.00%
Polymedco Poly stat hCG Serum	1	*Positive	1	100.00%
Stanbio True 20 Plus	1	*Positive	1	100.00%
Consult Diag hCG Combo Cassette	8	*Positive	8	100.00%
Alere hCG Combo Cassette (25 mIU/mL)	1	*Positive	1	100.00%
McKesson Medi-Lab hCG - Serum	4	*Positive	4	100.00%
HCG-12		* Denotes Correct Response		
AG: All Methods	86			
		*Positive	86	100.00%
Beck Coulter ICON 25	7	*Positive	7	100.00%
Fisher HC Sure-Vue hCG STAT	2	*Positive	2	100.00%
Sekisui Diagnostics OSOM	2	*Positive	2	100.00%
Henry Schein One Step + (hCG)	5	*Positive	5	100.00%
Immunostics Immuno/hCG	1	*Positive	1	100.00%
Inverness Signify HCG	2	*Positive	2	100.00%
LifeSign Status hCG	1	*Positive	1	100.00%
McKesson Medi-Lab hCG - Urine	1	*Positive	1	100.00%
Quidel QuickVue + hCG	2	*Positive	2	100.00%
Quidel QuickVue One-Step hCG Combo (U/S)	30	*Positive	30	100.00%
Select Diag hCG Combo	3	*Positive	3	100.00%
SP Brand (Cardinal) HCG Combo	14	*Positive	14	100.00%
Wampole Clearview	1	*Positive	1	100.00%
Polymedco Poly stat hCG Serum	1	*Positive	1	100.00%
Stanbio True 20 Plus	1	*Positive	1	100.00%



Consult Diag hCG Combo Cassette	8	*Positive	8	100.00%
Alere hCG Combo Cassette (25 mIU/mL)	1	*Positive	1	100.00%
McKesson Medi-Lab hCG - Serum	4	*Positive	4	100.00%
HCG-13		* Denotes Correct Response		
AG: All Methods	86			
		*Negative	86	100.00%
Beck Coulter ICON 25	7	*Negative	7	100.00%
Fisher HC Sure-Vue hCG STAT	2	*Negative	2	100.00%
Sekisui Diagnostics OSOM	2	*Negative	2	100.00%
Henry Schein One Step + (hCG)	5	*Negative	5	100.00%
Immunostics Immuno/hCG	1	*Negative	1	100.00%
Inverness Signify HCG	2	*Negative	2	100.00%
LifeSign Status hCG	1	*Negative	1	100.00%
McKesson Medi-Lab hCG - Urine	1	*Negative	1	100.00%
Quidel QuickVue + hCG	2	*Negative	2	100.00%
Quidel QuickVue One-Step hCG Combo (U/S)	30	*Negative	30	100.00%
Select Diag hCG Combo	3	*Negative	3	100.00%
SP Brand (Cardinal) HCG Combo	14	*Negative	14	100.00%
Wampole Clearview	1	*Negative	1	100.00%
Polymedco Poly stat hCG Serum	1	*Negative	1	100.00%
Stanbio True 20 Plus	1	*Negative	1	100.00%
Consult Diag hCG Combo Cassette	8	*Negative	8	100.00%
Alere hCG Combo Cassette (25 mIU/mL)	1	*Negative	1	100.00%
McKesson Medi-Lab hCG - Serum	4	*Negative	4	100.00%
HCG-14		* Denotes Correct Response		
AG: All Methods	86			
		*Negative	86	100.00%
Beck Coulter ICON 25	7	*Negative	7	100.00%
Fisher HC Sure-Vue hCG STAT	2	*Negative	2	100.00%
Sekisui Diagnostics OSOM	2	*Negative	2	100.00%
Henry Schein One Step + (hCG)	5	*Negative	5	100.00%
Immunostics Immuno/hCG	1	*Negative	1	100.00%
Inverness Signify HCG	2	*Negative	2	100.00%
LifeSign Status hCG	1	*Negative	1	100.00%
McKesson Medi-Lab hCG - Urine	1	*Negative	1	100.00%
Quidel QuickVue + hCG	2	*Negative	2	100.00%
Quidel QuickVue One-Step hCG Combo (U/S)	30	*Negative	30	100.00%
Select Diag hCG Combo	3	*Negative	3	100.00%
SP Brand (Cardinal) HCG Combo	14	*Negative	14	100.00%
Wampole Clearview	1	*Negative	1	100.00%
Polymedco Poly stat hCG Serum	1	*Negative	1	100.00%
Stanbio True 20 Plus	1	*Negative	1	100.00%
Consult Diag hCG Combo Cassette	8	*Negative	8	100.00%
Alere hCG Combo Cassette (25 mIU/mL)	1	*Negative	1	100.00%
McKesson Medi-Lab hCG - Serum	4	*Negative	4	100.00%
HCG-15		* Denotes Correct Response		
AG: All Methods	86			
		*Positive	86	100.00%
Beck Coulter ICON 25	7	*Positive	7	100.00%
Fisher HC Sure-Vue hCG STAT	2	*Positive	2	100.00%
Sekisui Diagnostics OSOM	2	*Positive	2	100.00%
Henry Schein One Step + (hCG)	5	*Positive	5	100.00%
Immunostics Immuno/hCG	1	*Positive	1	100.00%
Inverness Signify HCG	2	*Positive	2	100.00%
LifeSign Status hCG	1	*Positive	1	100.00%
McKesson Medi-Lab hCG - Urine	1	*Positive	1	100.00%
Quidel QuickVue + hCG	2	*Positive	2	100.00%
Quidel QuickVue One-Step hCG Combo (U/S)	30	*Positive	30	100.00%
Select Diag hCG Combo	3	*Positive	3	100.00%
SP Brand (Cardinal) HCG Combo	14	*Positive	14	100.00%
Wampole Clearview	1	*Positive	1	100.00%
Polymedco Poly stat hCG Serum	1	*Positive	1	100.00%
Stanbio True 20 Plus	1	*Positive	1	100.00%
Consult Diag hCG Combo Cassette	8	*Positive	8	100.00%
Alere hCG Combo Cassette (25 mIU/mL)	1	*Positive	1	100.00%
McKesson Medi-Lab hCG - Serum	4	*Positive	4	100.00%
		HCG - Qual Serum (International)		
HCG-11		* Denotes Correct Response		
AG: All Methods	3			



		*Positive	3	100.00%
Henry Schein One Step + (hCG)	1	*Positive	1	100.00%
SAS Ultra HCG	2	*Positive	2	100.00%
HCG-12		* Denotes Correct Response		
AG: All Methods	3			
		*Positive	3	100.00%
Henry Schein One Step + (hCG)	1	*Positive	1	100.00%
SAS Ultra HCG	2	*Positive	2	100.00%
HCG-13		* Denotes Correct Response		
AG: All Methods	3			
		*Negative	3	100.00%
Henry Schein One Step + (hCG)	1	*Negative	1	100.00%
SAS Ultra HCG	2	*Negative	2	100.00%
HCG-14		* Denotes Correct Response		
AG: All Methods	3			
		*Negative	3	100.00%
Henry Schein One Step + (hCG)	1	*Negative	1	100.00%
SAS Ultra HCG	2	*Negative	2	100.00%
HCG-15		* Denotes Correct Response		
AG: All Methods	3			
		*Positive	3	100.00%
Henry Schein One Step + (hCG)	1	*Positive	1	100.00%
SAS Ultra HCG	2	*Positive	2	100.00%
933 - HIV Antibody/Antigen Combo (2 specimens);934 - HIV Antibody/Antigen Combo (5 specimens)	Labs	Response	Number	Percent
		HIV Antibody (Waived)		
HIV-11		* Denotes Correct Response		
AG: All Methods	107			
		*Negative	101	94.39%
		Positive	6	5.61%
Alere Clearview HIV 1/2 STAT PAK (WB)	47	*Negative	47	100.00%
Orasure OraQuick Anti-HIV	9	*Negative	8	88.89%
Trinity Biotech Uni-Gold Recombigen HIV	2	*Negative	2	100.00%
Alere Determine HIV Ag/Ab	33	*Negative	29	87.88%
Chembio HIV 1/2 STAT-PAK Dipstick	4	*Negative	4	100.00%
Chembio HIV 1/2 STAT-PAK Cassette	5	*Negative	5	100.00%
bioLytical INSTI HIV Antibody	6	*Negative	6	100.00%
Alere Clearview Complete HIV 1/2 (mod)	1	Positive	1	100.00%
Orasure OraQuick Anti-HIV	9	Positive	1	11.11%
Alere Determine HIV Ag/Ab	33	Positive	4	12.12%
HIV-12		* Denotes Correct Response		
AG: All Methods	107			
		*Negative	107	100.00%
Alere Clearview Complete HIV 1/2 (mod)	1	*Negative	1	100.00%
Alere Clearview HIV 1/2 STAT PAK (WB)	47	*Negative	47	100.00%
Orasure OraQuick Anti-HIV	9	*Negative	9	100.00%
Trinity Biotech Uni-Gold Recombigen HIV	2	*Negative	2	100.00%
Alere Determine HIV Ag/Ab	33	*Negative	33	100.00%
Chembio HIV 1/2 STAT-PAK Dipstick	4	*Negative	4	100.00%
Chembio HIV 1/2 STAT-PAK Cassette	5	*Negative	5	100.00%
bioLytical INSTI HIV Antibody	6	*Negative	6	100.00%
		HIV Antibody (International)		
HIV-11		* Denotes Correct Response		
AG: All Methods	5			
		*Negative	5	100.00%
Alere Clearview HIV 1/2 STAT PAK (WB)	1	*Negative	1	100.00%
Intl OCD Vitros ECI/ECiQ	3	*Negative	3	100.00%
Intl Chembio HIV 1/2 STAT-PAK Dipstick	1	*Negative	1	100.00%
HIV-12		* Denotes Correct Response		
AG: All Methods	5			
		*Negative	5	100.00%
Alere Clearview HIV 1/2 STAT PAK (WB)	1	*Negative	1	100.00%
Intl OCD Vitros ECI/ECiQ	3	*Negative	3	100.00%
Intl Chembio HIV 1/2 STAT-PAK Dipstick	1	*Negative	1	100.00%
		HIV Antigen (p24) (Waived)		
HIV-11		* Denotes Correct Response		
AG: All Methods	33			
		*Positive	33	100.00%



Alere Determine HIV Ag/Ab	33	*Positive	33	100.00%
HIV-12		* Denotes Correct Response		
AG: All Methods	33			
		*Negative	33	100.00%
Alere Determine HIV Ag/Ab	33	*Negative	33	100.00%
934 - HIV Antibody/Antigen Combo (5 specimens)	Labs	Response	Number	Percent
		HIV Antibody (Waived)		
HIV-11		* Denotes Correct Response		
AG: All Methods	107			
		*Negative	101	94.39%
		Positive	6	5.61%
Alere Clearview HIV 1/2 STAT PAK (WB)	47	*Negative	47	100.00%
Orasure OraQuick Anti-HIV	9	*Negative	8	88.89%
Trinity Biotech Uni-Gold Recombigen HIV	2	*Negative	2	100.00%
Alere Determine HIV Ag/Ab	33	*Negative	29	87.88%
Chembio HIV 1/2 STAT-PAK Dipstick	4	*Negative	4	100.00%
Chembio HIV 1/2 STAT-PAK Cassette	5	*Negative	5	100.00%
bioLytical INSTI HIV Antibody	6	*Negative	6	100.00%
Alere Clearview Complete HIV 1/2 (mod)	1	Positive	1	100.00%
Orasure OraQuick Anti-HIV	9	Positive	1	11.11%
Alere Determine HIV Ag/Ab	33	Positive	4	12.12%
HIV-12		* Denotes Correct Response		
AG: All Methods	107			
		*Negative	107	100.00%
Alere Clearview Complete HIV 1/2 (mod)	1	*Negative	1	100.00%
Alere Clearview HIV 1/2 STAT PAK (WB)	47	*Negative	47	100.00%
Orasure OraQuick Anti-HIV	9	*Negative	9	100.00%
Trinity Biotech Uni-Gold Recombigen HIV	2	*Negative	2	100.00%
Alere Determine HIV Ag/Ab	33	*Negative	33	100.00%
Chembio HIV 1/2 STAT-PAK Dipstick	4	*Negative	4	100.00%
Chembio HIV 1/2 STAT-PAK Cassette	5	*Negative	5	100.00%
bioLytical INSTI HIV Antibody	6	*Negative	6	100.00%
HIV-13		* Denotes Correct Response		
AG: All Methods	9			
		*Positive	9	100.00%
Alere Clearview Complete HIV 1/2 (mod)	1	*Positive	1	100.00%
Alere Clearview HIV 1/2 STAT PAK (WB)	4	*Positive	4	100.00%
Orasure OraQuick Anti-HIV	1	*Positive	1	100.00%
Trinity Biotech Uni-Gold Recombigen HIV	1	*Positive	1	100.00%
Alere Determine HIV Ag/Ab	1	*Positive	1	100.00%
bioLytical INSTI HIV Antibody	1	*Positive	1	100.00%
HIV-14		* Denotes Correct Response		
AG: All Methods	9			
		*Negative	9	100.00%
Alere Clearview Complete HIV 1/2 (mod)	1	*Negative	1	100.00%
Alere Clearview HIV 1/2 STAT PAK (WB)	4	*Negative	4	100.00%
Orasure OraQuick Anti-HIV	1	*Negative	1	100.00%
Trinity Biotech Uni-Gold Recombigen HIV	1	*Negative	1	100.00%
Alere Determine HIV Ag/Ab	1	*Negative	1	100.00%
bioLytical INSTI HIV Antibody	1	*Negative	1	100.00%
HIV-15		* Denotes Correct Response		
AG: All Methods	9			
		*Positive	9	100.00%
Alere Clearview Complete HIV 1/2 (mod)	1	*Positive	1	100.00%
Alere Clearview HIV 1/2 STAT PAK (WB)	4	*Positive	4	100.00%
Orasure OraQuick Anti-HIV	1	*Positive	1	100.00%
Trinity Biotech Uni-Gold Recombigen HIV	1	*Positive	1	100.00%
Alere Determine HIV Ag/Ab	1	*Positive	1	100.00%
bioLytical INSTI HIV Antibody	1	*Positive	1	100.00%
		HIV Antibody (International)		
HIV-11		* Denotes Correct Response		
AG: All Methods	5			
		*Negative	5	100.00%
Alere Clearview HIV 1/2 STAT PAK (WB)	1	*Negative	1	100.00%
Intl OCD Vitros ECI/ECiQ	3	*Negative	3	100.00%
Intl Chembio HIV 1/2 STAT-PAK Dipstick	1	*Negative	1	100.00%
HIV-12		* Denotes Correct Response		
AG: All Methods	5			



		*Negative	5	100.00%
Alere Clearview HIV 1/2 STAT PAK (WB)	1	*Negative	1	100.00%
Intl OCD Vitros ECi/ECiQ	3	*Negative	3	100.00%
Intl Chembio HIV 1/2 STAT-PAK Dipstick	1	*Negative	1	100.00%
HIV-13		* Denotes Correct Response		
AG: All Methods	3			
		*Positive	3	100.00%
Intl OCD Vitros ECi/ECiQ	3	*Positive	3	100.00%
HIV-14		* Denotes Correct Response		
AG: All Methods	3			
		*Negative	3	100.00%
Intl OCD Vitros ECi/ECiQ	3	*Negative	3	100.00%
HIV-15		* Denotes Correct Response		
AG: All Methods	3			
		*Positive	3	100.00%
Intl OCD Vitros ECi/ECiQ	3	*Positive	3	100.00%
		HIV Antigen (p24) (Waived)		
HIV-11		* Denotes Correct Response		
AG: All Methods	33			
		*Positive	33	100.00%
Alere Determine HIV Ag/Ab	33	*Positive	33	100.00%
HIV-12		* Denotes Correct Response		
AG: All Methods	33			
		*Negative	33	100.00%
Alere Determine HIV Ag/Ab	33	*Negative	33	100.00%
HIV-13		* Denotes Correct Response		
AG: All Methods	1			
		*Negative	1	100.00%
Alere Determine HIV Ag/Ab	1	*Negative	1	100.00%
HIV-14		* Denotes Correct Response		
AG: All Methods	1			
		*Negative	1	100.00%
Alere Determine HIV Ag/Ab	1	*Negative	1	100.00%
HIV-15		* Denotes Correct Response		
AG: All Methods	1			
		*Negative	1	100.00%
Alere Determine HIV Ag/Ab	1	*Negative	1	100.00%
644 - H. pylori	Labs	Response	Number	Percent
		H. pylori Ab Detect		
HP-3		* Denotes Correct Response		
AG: All Methods	66			
		*Negative	65	98.48%
		Positive	1	1.52%
Beck Coulter FlexSure	3	*Negative	3	100.00%
Beck Coulter ICON HP	2	*Negative	2	100.00%
McKesson Medi-Lab Perf H. pylori	6	*Negative	6	100.00%
Polymedco Poly stat (WB - Waived)	1	*Negative	1	100.00%
Henry Schein One Step+ (H. pylori)	8	*Negative	8	100.00%
Select Diag H. pylori Cassette	2	*Negative	2	100.00%
Quidel QuickVue H. pylori gII (S/P)	3	*Negative	3	100.00%
Quidel QuickVue H. pylori gII (WB)	9	*Negative	8	88.89%
SP Brand (Cardinal) HP Rapid Test	9	*Negative	9	100.00%
Alere Clearview II S/P (mod)	3	*Negative	3	100.00%
Alere Clearview II WB (waived)	3	*Negative	3	100.00%
LifeSign Status H. pylori	2	*Negative	2	100.00%
Instant-View H. pylori	1	*Negative	1	100.00%
Consult Diag H pylori	13	*Negative	13	100.00%
Jant Accutest H pylori	1	*Negative	1	100.00%
Quidel QuickVue H. pylori gII (WB)	9	Positive	1	11.11%
603 - Blood Cell Identification;610 - Regular Hematology	Labs	Response	Number	Percent
		Blood Cell Identification		
ID-13		* Denotes Correct Response		
AG: All Methods	61			
		*Lymphocyte	59	96.72%
		Polymphocyte	2	3.28%
Blood Cell Microscopy	61	*Lymphocyte	59	96.72%
Blood Cell Microscopy		Polymphocyte	2	3.28%
ID-14		* Denotes Correct Response		



AG: All Methods	61			
		*Segmented or band neutrophil	61	100.00%
Blood Cell Microscopy	61	*Segmented or band neutrophil	61	100.00%
ID-15		* Denotes Correct Response		
AG: All Methods	61			
		*Immature/abnormal cell - Referred	2	3.28%
		*Nucleated red cell	59	96.72%
Blood Cell Microscopy	61	*Immature/abnormal cell - Referred	2	3.28%
Blood Cell Microscopy		*Nucleated red cell	59	96.72%
ID-16		* Denotes Correct Response		
AG: All Methods	60			
		*Monocyte	54	90.00%
		Immature/abnormal cell - Referred	1	1.67%
		Metamyelocyte (juvenile)	5	8.33%
Blood Cell Microscopy	60	*Monocyte	54	90.00%
Blood Cell Microscopy		Immature/abnormal cell - Referred	1	1.67%
Blood Cell Microscopy		Metamyelocyte (juvenile)	5	8.33%
ID-17		* Denotes Correct Response		
AG: All Methods	61			
		*Hypochromasia	44	72.13%
		*Immature/abnormal cell - Referred	1	1.64%
		Anisocytosis	2	3.28%
		Red Blood Cell, normal (RBC) (Erythrocyte)	14	22.95%
Blood Cell Microscopy	61	*Hypochromasia	44	72.13%
Blood Cell Microscopy		*Immature/abnormal cell - Referred	1	1.64%
Blood Cell Microscopy		Anisocytosis	2	3.28%
Blood Cell Microscopy		Red Blood Cell, normal (RBC) (Erythrocyte)	14	22.95%
ID-18		* Denotes Correct Response		
AG: All Methods	61			
		*Blast cell	18	29.51%
		*Identification Unknown - Referred	2	3.28%
		*Immature/abnormal cell - Referred	16	26.23%
		*Lymphoblast	13	21.31%
		Lymphocyte, immature	3	4.92%
		Lymphocyte, reactive (atypical, variant, Downey)	2	3.28%
		Monoblast	1	1.64%
		Myelocyte cell	1	1.64%
		Plasmoblast	1	1.64%
		Prolymphocyte	3	4.92%
		Promonocyte	1	1.64%
Blood Cell Microscopy	61	*Blast cell	18	29.51%
Blood Cell Microscopy		*Identification Unknown - Referred	2	3.28%
Blood Cell Microscopy		*Immature/abnormal cell - Referred	16	26.23%
Blood Cell Microscopy		*Lymphoblast	13	21.31%
Blood Cell Microscopy		Lymphocyte, immature	3	4.92%
Blood Cell Microscopy		Lymphocyte, reactive (atypical, variant, Downey)	2	3.28%
Blood Cell Microscopy		Monoblast	1	1.64%
Blood Cell Microscopy		Myelocyte cell	1	1.64%
Blood Cell Microscopy		Plasmoblast	1	1.64%
Blood Cell Microscopy		Prolymphocyte	3	4.92%
Blood Cell Microscopy		Promonocyte	1	1.64%
		Blood Cell Identification (International)		
ID-13		* Denotes Correct Response		
AG: All Methods	1			
		*Lymphocyte	1	100.00%
Intl Blood Cell Microscopy	1	*Lymphocyte	1	100.00%
ID-14		* Denotes Correct Response		
AG: All Methods	1			
		*Segmented or band neutrophil	1	100.00%
Intl Blood Cell Microscopy	1	*Segmented or band neutrophil	1	100.00%
ID-15		* Denotes Correct Response		
AG: All Methods	1			
		*Nucleated red cell	1	100.00%
Intl Blood Cell Microscopy	1	*Nucleated red cell	1	100.00%
ID-16		* Denotes Correct Response		
AG: All Methods	1			
		*Monocyte	1	100.00%
Intl Blood Cell Microscopy	1	*Monocyte	1	100.00%



ID-17		* Denotes Correct Response		
AG: All Methods	1			
		*Hypochromasia	1	100.00%
Intl Blood Cell Microscopy	1	*Hypochromasia	1	100.00%
ID-18		* Denotes Correct Response		
AG: All Methods	1			
		*Lymphoblast	1	100.00%
Intl Blood Cell Microscopy	1	*Lymphoblast	1	100.00%
609 - Waived/PPM Tests 675 - Infectious Mononucleosis; 753 - Infectious Mononucleosis (Waived) (1 specimens); 754 - Infectious Mono/RF Combo	Labs	Response	Number	Percent
		Infectious Mono (Waived)		
IM-11		* Denotes Correct Response		
AG: All Methods	500			
		*Positive	493	98.60%
		Negative	7	1.40%
Alere Acceava (S)	1	*Positive	1	100.00%
Alere Acceava (WB)	55	*Positive	55	100.00%
Sekisui Diagnostics OSOM Mono (WB)	205	*Positive	202	98.54%
Henry Schein One-Step+ Mono	47	*Positive	46	97.87%
Instant Technologies i-Screen (WB)	1	*Positive	1	100.00%
Inverness Signify Mono	4	*Positive	4	100.00%
LifeSign Status Inf Mono	32	*Positive	32	100.00%
McKesson Medi-Lab Perf - Inf Mono	11	*Positive	11	100.00%
Polymedco Poly stat Mono - (WB)	3	*Positive	3	100.00%
Quidel QuickVue+ (WB)	11	*Positive	11	100.00%
SP Brand (Cardinal) Rapid Test	15	*Positive	14	93.33%
Alere Clearview Mono Cassette	20	*Positive	20	100.00%
Pulse Scientific	2	*Positive	2	100.00%
Diag Test Group Clarity Mono	3	*Positive	3	100.00%
Fisher HC Sure-View Mono	6	*Positive	6	100.00%
Immunostics Inc. ImmunoDetector Mono (WB)	2	*Positive	2	100.00%
Meridian Immunocard STAT! Mono (WB)	6	*Positive	6	100.00%
NDC ProAdvantage Mono	2	*Positive	2	100.00%
Consult Diag Mono Cassette	59	*Positive	57	96.61%
Alere Acceava (WB) Mono Cassette	14	*Positive	14	100.00%
Medline Mono Cassette (WB)	1	*Positive	1	100.00%
Sekisui Diagnostics OSOM Mono (WB)	205	Negative	3	1.46%
Henry Schein One-Step+ Mono	47	Negative	1	2.13%
SP Brand (Cardinal) Rapid Test	15	Negative	1	6.67%
Consult Diag Mono Cassette	59	Negative	2	3.39%
		Infectious Mononucleosis		
IM-11		* Denotes Correct Response		
AG: All Methods	45			
		*Positive	45	100.00%
Beck Coulter ICON (WB)	1	*Positive	1	100.00%
Alere Acceava (WB)	1	*Positive	1	100.00%
Sekisui Diagnostics OSOM Mono (WB)	9	*Positive	9	100.00%
LifeSign Status Inf Mono	2	*Positive	2	100.00%
Lifesign Unistep	1	*Positive	1	100.00%
McKesson Medi-Lab Perf - Inf Mono	2	*Positive	2	100.00%
Polymedco Poly stat Mono - (WB)	1	*Positive	1	100.00%
Polymedco Poly stat Mono (S/P)	2	*Positive	2	100.00%
Quidel QuickVue+ (S/P)	2	*Positive	2	100.00%
Quidel QuickVue+ (WB)	7	*Positive	7	100.00%
Seradyn Color Slide II	3	*Positive	3	100.00%
SP Brand (Cardinal) Rapid Test	1	*Positive	1	100.00%
Alere Clearview Mono Cassette	8	*Positive	8	100.00%
Sekisui Diagnostics OSOM Mono Non-Waived	1	*Positive	1	100.00%
Consult Diag Mono Cassette	4	*Positive	4	100.00%
675 - Infectious Mononucleosis; 754 - Infectious Mono/RF Combo	Labs	Response	Number	Percent
		Infectious Mononucleosis		
IM-11		* Denotes Correct Response		
AG: All Methods	45			
		*Positive	45	100.00%
Beck Coulter ICON (WB)	1	*Positive	1	100.00%
Alere Acceava (WB)	1	*Positive	1	100.00%
Sekisui Diagnostics OSOM Mono (WB)	9	*Positive	9	100.00%



LifeSign Status Inf Mono	2	*Positive	2	100.00%
Lifesign Unistep	1	*Positive	1	100.00%
McKesson Medi-Lab Perf - Inf Mono	2	*Positive	2	100.00%
Polymedco Poly stat Mono - (WB)	1	*Positive	1	100.00%
Polymedco Poly stat Mono (S/P)	2	*Positive	2	100.00%
Quidel QuickVue+ (S/P)	2	*Positive	2	100.00%
Quidel QuickVue+ (WB)	7	*Positive	7	100.00%
Seradyn Color Slide II	3	*Positive	3	100.00%
SP Brand (Cardinal) Rapid Test	1	*Positive	1	100.00%
Alere Clearview Mono Cassette	8	*Positive	8	100.00%
Sekisui Diagnostics OSOM Mono Non-Waived	1	*Positive	1	100.00%
Consult Diag Mono Cassette	4	*Positive	4	100.00%
IM-12		* Denotes Correct Response		
AG: All Methods	45			
		*Negative	45	100.00%
Beck Coulter ICON (WB)	1	*Negative	1	100.00%
Alere Acceava (WB)	1	*Negative	1	100.00%
Sekisui Diagnostics OSOM Mono (WB)	9	*Negative	9	100.00%
LifeSign Status Inf Mono	2	*Negative	2	100.00%
Lifesign Unistep	1	*Negative	1	100.00%
McKesson Medi-Lab Perf - Inf Mono	2	*Negative	2	100.00%
Polymedco Poly stat Mono - (WB)	1	*Negative	1	100.00%
Polymedco Poly stat Mono (S/P)	2	*Negative	2	100.00%
Quidel QuickVue+ (S/P)	2	*Negative	2	100.00%
Quidel QuickVue+ (WB)	7	*Negative	7	100.00%
Seradyn Color Slide II	3	*Negative	3	100.00%
SP Brand (Cardinal) Rapid Test	1	*Negative	1	100.00%
Alere Clearview Mono Cassette	8	*Negative	8	100.00%
Sekisui Diagnostics OSOM Mono Non-Waived	1	*Negative	1	100.00%
Consult Diag Mono Cassette	4	*Negative	4	100.00%
IM-13		* Denotes Correct Response		
AG: All Methods	45			
		*Positive	45	100.00%
Beck Coulter ICON (WB)	1	*Positive	1	100.00%
Alere Acceava (WB)	1	*Positive	1	100.00%
Sekisui Diagnostics OSOM Mono (WB)	9	*Positive	9	100.00%
LifeSign Status Inf Mono	2	*Positive	2	100.00%
Lifesign Unistep	1	*Positive	1	100.00%
McKesson Medi-Lab Perf - Inf Mono	2	*Positive	2	100.00%
Polymedco Poly stat Mono - (WB)	1	*Positive	1	100.00%
Polymedco Poly stat Mono (S/P)	2	*Positive	2	100.00%
Quidel QuickVue+ (S/P)	2	*Positive	2	100.00%
Quidel QuickVue+ (WB)	7	*Positive	7	100.00%
Seradyn Color Slide II	3	*Positive	3	100.00%
SP Brand (Cardinal) Rapid Test	1	*Positive	1	100.00%
Alere Clearview Mono Cassette	8	*Positive	8	100.00%
Sekisui Diagnostics OSOM Mono Non-Waived	1	*Positive	1	100.00%
Consult Diag Mono Cassette	4	*Positive	4	100.00%
IM-14		* Denotes Correct Response		
AG: All Methods	45			
		*Negative	45	100.00%
Beck Coulter ICON (WB)	1	*Negative	1	100.00%
Alere Acceava (WB)	1	*Negative	1	100.00%
Sekisui Diagnostics OSOM Mono (WB)	9	*Negative	9	100.00%
LifeSign Status Inf Mono	2	*Negative	2	100.00%
Lifesign Unistep	1	*Negative	1	100.00%
McKesson Medi-Lab Perf - Inf Mono	2	*Negative	2	100.00%
Polymedco Poly stat Mono - (WB)	1	*Negative	1	100.00%
Polymedco Poly stat Mono (S/P)	2	*Negative	2	100.00%
Quidel QuickVue+ (S/P)	2	*Negative	2	100.00%
Quidel QuickVue+ (WB)	7	*Negative	7	100.00%
Seradyn Color Slide II	3	*Negative	3	100.00%
SP Brand (Cardinal) Rapid Test	1	*Negative	1	100.00%
Alere Clearview Mono Cassette	8	*Negative	8	100.00%
Sekisui Diagnostics OSOM Mono Non-Waived	1	*Negative	1	100.00%
Consult Diag Mono Cassette	4	*Negative	4	100.00%
IM-15		* Denotes Correct Response		
AG: All Methods	45			
		*Positive	45	100.00%



Beck Coulter ICON (WB)	1	*Positive	1	100.00%
Alere Acceava (WB)	1	*Positive	1	100.00%
Sekisui Diagnostics OSOM Mono (WB)	9	*Positive	9	100.00%
LifeSign Status Inf Mono	2	*Positive	2	100.00%
Lifesign Unistep	1	*Positive	1	100.00%
McKesson Medi-Lab Perf - Inf Mono	2	*Positive	2	100.00%
Polymedco Poly stat Mono - (WB)	1	*Positive	1	100.00%
Polymedco Poly stat Mono (S/P)	2	*Positive	2	100.00%
Quidel QuickVue+ (S/P)	2	*Positive	2	100.00%
Quidel QuickVue+ (WB)	7	*Positive	7	100.00%
Seradyn Color Slide II	3	*Positive	3	100.00%
SP Brand (Cardinal) Rapid Test	1	*Positive	1	100.00%
Alere Clearview Mono Cassette	8	*Positive	8	100.00%
Sekisui Diagnostics OSOM Mono Non-Waived	1	*Positive	1	100.00%
Consult Diag Mono Cassette	4	*Positive	4	100.00%
		Infectious Mono (Waived)		
IM-11		* Denotes Correct Response		
AG: All Methods	500			
		*Positive	493	98.60%
		Negative	7	1.40%
Alere Acceava (S)	1	*Positive	1	100.00%
Alere Acceava (WB)	55	*Positive	55	100.00%
Sekisui Diagnostics OSOM Mono (WB)	205	*Positive	202	98.54%
Henry Schein One-Step+ Mono	47	*Positive	46	97.87%
Instant Technologies i-Screen (WB)	1	*Positive	1	100.00%
Inverness Signify Mono	4	*Positive	4	100.00%
LifeSign Status Inf Mono	32	*Positive	32	100.00%
McKesson Medi-Lab Perf - Inf Mono	11	*Positive	11	100.00%
Polymedco Poly stat Mono - (WB)	3	*Positive	3	100.00%
Quidel QuickVue+ (WB)	11	*Positive	11	100.00%
SP Brand (Cardinal) Rapid Test	15	*Positive	14	93.33%
Alere Clearview Mono Cassette	20	*Positive	20	100.00%
Pulse Scientific	2	*Positive	2	100.00%
Diag Test Group Clarity Mono	3	*Positive	3	100.00%
Fisher HC Sure-Vue Mono	6	*Positive	6	100.00%
Immunostics Inc. ImmunoDetector Mono (WB)	2	*Positive	2	100.00%
Meridian Immunocard STAT! Mono (WB)	6	*Positive	6	100.00%
NDC ProAdvantage Mono	2	*Positive	2	100.00%
Consult Diag Mono Cassette	59	*Positive	57	96.61%
Alere Acceava (WB) Mono Cassette	14	*Positive	14	100.00%
Medline Mono Cassette (WB)	1	*Positive	1	100.00%
Sekisui Diagnostics OSOM Mono (WB)	205	Negative	3	1.46%
Henry Schein One-Step+ Mono	47	Negative	1	2.13%
SP Brand (Cardinal) Rapid Test	15	Negative	1	6.67%
Consult Diag Mono Cassette	59	Negative	2	3.39%
		Infectious Mononucleosis (International)		
IM-11		* Denotes Correct Response		
AG: All Methods	2			
		*Positive	2	100.00%
Alere Clearview Mono Cassette	1	*Positive	1	100.00%
Fisher HC Sure-Vue Mono	1	*Positive	1	100.00%
IM-12		* Denotes Correct Response		
AG: All Methods	2			
		*Negative	2	100.00%
Alere Clearview Mono Cassette	1	*Negative	1	100.00%
Fisher HC Sure-Vue Mono	1	*Negative	1	100.00%
IM-13		* Denotes Correct Response		
AG: All Methods	2			
		*Positive	2	100.00%
Alere Clearview Mono Cassette	1	*Positive	1	100.00%
Fisher HC Sure-Vue Mono	1	*Positive	1	100.00%
IM-14		* Denotes Correct Response		
AG: All Methods	2			
		*Negative	2	100.00%
Alere Clearview Mono Cassette	1	*Negative	1	100.00%
Fisher HC Sure-Vue Mono	1	*Negative	1	100.00%
IM-15		* Denotes Correct Response		
AG: All Methods	2			



		*Positive	2	100.00%
Alere Clearview Mono Cassette	1	*Positive	1	100.00%
Fisher HC Sure-View Mono	1	*Positive	1	100.00%
637 - KOH - Glass Slides (Skin Hair Nail)	Labs	Response	Number	Percent
		KOH Preparation		
KOH-5		* Denotes Correct Response		
AG: All Methods	20			
		*Fungal elements/yeast observed	20	100.00%
KOH Glass Slide	20	*Fungal elements/yeast observed	20	100.00%
KOH-6		* Denotes Correct Response		
AG: All Methods	20			
		*Fungal elements/yeast observed	20	100.00%
KOH Glass Slide	20	*Fungal elements/yeast observed	20	100.00%
758 - Mycoplasma Antibody	Labs	Response	Number	Percent
		Mycoplasma Antibody		
MY-3		* Denotes Correct Response		
AG: All Methods	4			
		*Positive	4	100.00%
Meridian ImmunoCard	4	*Positive	4	100.00%
791 - Fecal Occult Blood Add-On;601 - Waived/PPM Tests II;609 - Waived/PPM Tests I;647 - Mini-Waived Package	Labs	Response	Number	Percent
		Fecal Occult Blood		
OC-3		* Denotes Correct Response		
MG: Guaiac methods	722			
		*Positive	679	94.04%
		Negative	43	5.96%
OC-3		* Denotes Correct Response		
MG: Immunochemical methods	99			
		*Positive	91	91.92%
		Negative	8	8.08%
950 - Post Vasectomy Sperm Analysis	Labs	Response	Number	Percent
		Sperm Count (Post Vasectomy)		
PVAS-5		* Denotes Correct Response		
AG: All Methods	13			
		*Positive for Sperm	13	100.00%
Hemocytometer	4	*Positive for Sperm	4	100.00%
Post Vasectomy Microscopy	9	*Positive for Sperm	9	100.00%
PVAS-6		* Denotes Correct Response		
AG: All Methods	13			
		*Negative for Sperm	13	100.00%
Hemocytometer	4	*Negative for Sperm	4	100.00%
Post Vasectomy Microscopy	9	*Negative for Sperm	9	100.00%
754 - Infectious Mono/RF Combo;653 - Rheumatoid Factor	Labs	Response	Number	Percent
		Rheumatoid Factor		
RF-11		* Denotes Correct Response		
AG: All Methods	9			
		*N/A (Negative)	9	100.00%
Siemens RapiTex	2	*N/A (Negative)	2	100.00%
Seradyn Color Slide	2	*N/A (Negative)	2	100.00%
Wampole ColorCard RF	2	*N/A (Negative)	2	100.00%
Wampole Rheumatex	2	*N/A (Negative)	2	100.00%
Remel SeraTest RF-Latex	1	*N/A (Negative)	1	100.00%
RF-12		* Denotes Correct Response		
AG: All Methods	9			
		*N/A (Negative)	9	100.00%
Siemens RapiTex	2	*N/A (Negative)	2	100.00%
Seradyn Color Slide	2	*N/A (Negative)	2	100.00%
Wampole ColorCard RF	2	*N/A (Negative)	2	100.00%
Wampole Rheumatex	2	*N/A (Negative)	2	100.00%
Remel SeraTest RF-Latex	1	*N/A (Negative)	1	100.00%
RF-13		* Denotes Correct Response		
AG: All Methods	9			
		*Positive	9	100.00%
Siemens RapiTex	2	*Positive	2	100.00%
Seradyn Color Slide	2	*Positive	2	100.00%
Wampole ColorCard RF	2	*Positive	2	100.00%



Wampole Rheumatex	2	*Positive	2	100.00%
Remel SeraTest RF-Latex	1	*Positive	1	100.00%
RF-14		* Denotes Correct Response		
AG: All Methods	9			
		*N/A (Negative)	9	100.00%
Siemens RapiTex	2	*N/A (Negative)	2	100.00%
Seradyn Color Slide	2	*N/A (Negative)	2	100.00%
Wampole ColorCard RF	2	*N/A (Negative)	2	100.00%
Wampole Rheumatex	2	*N/A (Negative)	2	100.00%
Remel SeraTest RF-Latex	1	*N/A (Negative)	1	100.00%
RF-15		* Denotes Correct Response		
AG: All Methods	9			
		*Positive	9	100.00%
Siemens RapiTex	2	*Positive	2	100.00%
Seradyn Color Slide	2	*Positive	2	100.00%
Wampole ColorCard RF	2	*Positive	2	100.00%
Wampole Rheumatex	2	*Positive	2	100.00%
Remel SeraTest RF-Latex	1	*Positive	1	100.00%
		Rheumatoid Factor (International)		
RF-11		* Denotes Correct Response		
AG: All Methods	3			
		*Negative	3	100.00%
Wampole ColorCard RF	1	*Negative	1	100.00%
Fisher HC Sure-Vue RF	1	*Negative	1	100.00%
Immunostics Rhuemacol/LEX	1	*Negative	1	100.00%
RF-12		* Denotes Correct Response		
AG: All Methods	3			
		*Negative	3	100.00%
Wampole ColorCard RF	1	*Negative	1	100.00%
Fisher HC Sure-Vue RF	1	*Negative	1	100.00%
Immunostics Rhuemacol/LEX	1	*Negative	1	100.00%
RF-13		* Denotes Correct Response		
AG: All Methods	3			
		*Positive	3	100.00%
Wampole ColorCard RF	1	*Positive	1	100.00%
Fisher HC Sure-Vue RF	1	*Positive	1	100.00%
Immunostics Rhuemacol/LEX	1	*Positive	1	100.00%
RF-14		* Denotes Correct Response		
AG: All Methods	3			
		*Negative	3	100.00%
Wampole ColorCard RF	1	*Negative	1	100.00%
Fisher HC Sure-Vue RF	1	*Negative	1	100.00%
Immunostics Rhuemacol/LEX	1	*Negative	1	100.00%
RF-15		* Denotes Correct Response		
AG: All Methods	3			
		*Positive	3	100.00%
Wampole ColorCard RF	1	*Positive	1	100.00%
Fisher HC Sure-Vue RF	1	*Positive	1	100.00%
Immunostics Rhuemacol/LEX	1	*Positive	1	100.00%
		Rheumatoid Factor (Semi Quant)		
RF-11		* Denotes Correct Response		
AG: All Methods	1			
		*Negative	1	100.00%
Fisher HC Sure-Vue RF (Semi Quant)	1	*Negative	1	100.00%
RF-12		* Denotes Correct Response		
AG: All Methods	1			
		*Negative	1	100.00%
Fisher HC Sure-Vue RF (Semi Quant)	1	*Negative	1	100.00%
RF-13		* Denotes Correct Response		
AG: All Methods	1			
		*80 IU/mL	1	100.00%
Fisher HC Sure-Vue RF (Semi Quant)	1	*80 IU/mL	1	100.00%
RF-14		* Denotes Correct Response		
AG: All Methods	1			
		*Negative	1	100.00%
Fisher HC Sure-Vue RF (Semi Quant)	1	*Negative	1	100.00%
RF-15		* Denotes Correct Response		
AG: All Methods	1			



		*160 IU/mL	1	100.00%
Fisher HC Sure-Vue RF (Semi Quant)	1	*160 IU/mL	1	100.00%
645 - Fetal Membrane Rupture	Labs	Response	Number	Percent
		Fetal Membranes Rupture PAMG-1		
ROM-5		* Denotes Correct Response		
AG: All Methods	18			
		*Negative	18	100.00%
AmniSure ROM Test	18	*Negative	18	100.00%
ROM-6		* Denotes Correct Response		
AG: All Methods	18			
		*Positive	18	100.00%
AmniSure ROM Test	18	*Positive	18	100.00%
601 - Waived/PPM Tests II;609 - Waived/PPM Tests I;785 - Presump ID/Colony Count-Urine & Strep Combo;691 - Group A Strep Antigen (5 specimens);699 - Group A Strep Antigen (1 specimen);783 - Group A Strep Antigen/Add-on (2 specimens);665 - Group A Strep An	Labs	Response	Number	Percent
		Strep A Ag Det. (Waived)		
RS-11		* Denotes Correct Response		
AG: All Methods	879			
		*Negative	869	98.86%
		Positive	10	1.14%
Abbott Signify Strep A - Waived	6	*Negative	6	100.00%
Beck Coulter ICON DS Strep	1	*Negative	1	100.00%
Alere Acceava Strep A Dipstick	115	*Negative	113	98.26%
Fisher HC Sure-Vue Strep A (Waived)	9	*Negative	9	100.00%
Sekisui Diagnostics OSOM	38	*Negative	38	100.00%
Sekisui Diagnostics OSOM Ultra Strep	45	*Negative	45	100.00%
Henry Schein OneStep + Rapid Strep Cass	3	*Negative	3	100.00%
Henry Schein OneStep+ Rapid Strep Dip	73	*Negative	71	97.26%
Innovacon One Step Strep Ag	1	*Negative	1	100.00%
Instant Technologies i Strep	1	*Negative	1	100.00%
Inverness Signify Strep A Dipstick	7	*Negative	7	100.00%
McKesson Strep A Cas	4	*Negative	4	100.00%
McKesson Strep A Dip	31	*Negative	31	100.00%
McKesson Strep A Test - Twist	2	*Negative	2	100.00%
Permaxim Rediscreen Strep A Rapid	1	*Negative	1	100.00%
Polymedco Poly Stat Strep A (W)	1	*Negative	1	100.00%
Quidel QuickVue Dipstick Strep	76	*Negative	76	100.00%
Quidel QuickVue In-Line - 25 Test	63	*Negative	62	98.41%
Quidel QuickVue+ Strep A	33	*Negative	32	96.97%
Select Diag Strep A Dip	4	*Negative	4	100.00%
SP (Cardinal) Rapid Strep Cassette	5	*Negative	5	100.00%
SP Brand (Cardinal) Rapid Strep A Dipstick	83	*Negative	83	100.00%
Alere Clearview Exact II Dip	35	*Negative	34	97.14%
Sekisui Diagnostics OSOM Ultra	7	*Negative	7	100.00%
Consult Diag Strep A Dipstick	121	*Negative	120	99.17%
Consult Diag Strep A Twist	1	*Negative	1	100.00%
MooreBrand Strep A Dipstick	6	*Negative	5	83.33%
NDC ProAdvantage Strep Strips	1	*Negative	1	100.00%
Meridian ImmunoCard STAT Strep A	1	*Negative	1	100.00%
LifeSign Bio Strep A Test STRip	1	*Negative	1	100.00%
Quidel Sofia/Sofia 2 Strep A FIA	31	*Negative	30	96.77%
BD Veritor Strep A	38	*Negative	38	100.00%
Alere i Strep A	11	*Negative	11	100.00%
Medline Strep A Dipstick	3	*Negative	3	100.00%
Roche Cobas Liat Strep A	21	*Negative	21	100.00%
Alere Acceava Strep A Dipstick	115	Positive	2	1.74%
Henry Schein OneStep+ Rapid Strep Dip	73	Positive	2	2.74%
Quidel QuickVue In-Line - 25 Test	63	Positive	1	1.59%
Quidel QuickVue+ Strep A	33	Positive	1	3.03%
Alere Clearview Exact II Dip	35	Positive	1	2.86%
Consult Diag Strep A Dipstick	121	Positive	1	0.83%
MooreBrand Strep A Dipstick	6	Positive	1	16.67%
Quidel Sofia/Sofia 2 Strep A FIA	31	Positive	1	3.23%
RS-12		* Denotes Correct Response		
AG: All Methods	55			
		*Positive	54	98.18%
		Negative	1	1.82%



Alere Acceava Strep A Dipstick	2	*Positive	2	100.00%
Sekisui Diagnostics OSOM	8	*Positive	8	100.00%
Sekisui Diagnostics OSOM Ultra Strep	4	*Positive	4	100.00%
Henry Schein OneStep+ Rapid Strep Dip	1	*Positive	1	100.00%
Innovacon One Step Strep Ag	1	*Positive	1	100.00%
McKesson Strep A Dip	1	*Positive	1	100.00%
Quidel QuickVue Dipstick Strep	1	*Positive	1	100.00%
Quidel QuickVue In-Line - 25 Test	3	*Positive	3	100.00%
Quidel QuickVue+ Strep A	1	*Positive	1	100.00%
SP Brand (Cardinal) Rapid Strep A Dipstick	4	*Positive	4	100.00%
Consult Diag Strep A Dipstick	7	*Positive	7	100.00%
BD Veritor Strep A	1	*Positive	1	100.00%
Roche Cobas Liat Strep A	20	*Positive	20	100.00%
MooreBrand Strep A Dipstick	1	Negative	1	100.00%
		Strep A Ag Det.		
RS-11		* Denotes Correct Response		
AG: All Methods	160			
		*Negative	160	100.00%
Abbott Signify Strep A - Waived	2	*Negative	2	100.00%
BD Check Strep A (Dip)	2	*Negative	2	100.00%
Beck Coulter ICON DS Strep	1	*Negative	1	100.00%
Alere Acceava Strep A Dipstick	3	*Negative	3	100.00%
Fisher HC Sure-Vue Select Strep Ag	1	*Negative	1	100.00%
Sekisui Diagnostics OSOM	16	*Negative	16	100.00%
Sekisui Diagnostics OSOM Ultra Strep	13	*Negative	13	100.00%
Henry Schein OneStep+ Rapid Strep Dip	10	*Negative	10	100.00%
Immunostics Immuno/Strep A Det	2	*Negative	2	100.00%
Inverness Signify Strep A Dipstick	1	*Negative	1	100.00%
McKesson Strep A Cas	2	*Negative	2	100.00%
McKesson Strep A Dip	5	*Negative	5	100.00%
Permaxim Rediscreen Strep A Rapid	1	*Negative	1	100.00%
PSS Rapid Test Strip Strep A	1	*Negative	1	100.00%
Quidel QuickVue Dipstick Strep	10	*Negative	10	100.00%
Quidel QuickVue In-Line - 25 Test	10	*Negative	10	100.00%
Quidel QuickVue+ Strep A	17	*Negative	17	100.00%
SP (Cardinal) Rapid Strep Cassette	1	*Negative	1	100.00%
SP Brand (Cardinal) Rapid Strep A Dipstick	2	*Negative	2	100.00%
Stanbio QuStick Strep A	2	*Negative	2	100.00%
Consult Diag Strep A Dipstick	18	*Negative	18	100.00%
Meridian ImmunoCard STAT Strep A	1	*Negative	1	100.00%
Quidel Sofia/Sofia 2 Strep A FIA	8	*Negative	8	100.00%
BD Veritor Strep A	6	*Negative	6	100.00%
Quidel Solana	6	*Negative	6	100.00%
Cepheid GeneXpert - Xpert Strep A	19	*Negative	19	100.00%
RS-12		* Denotes Correct Response		
AG: All Methods	160			
		*Positive	160	100.00%
Abbott Signify Strep A - Waived	2	*Positive	2	100.00%
BD Check Strep A (Dip)	2	*Positive	2	100.00%
Beck Coulter ICON DS Strep	1	*Positive	1	100.00%
Alere Acceava Strep A Dipstick	3	*Positive	3	100.00%
Fisher HC Sure-Vue Select Strep Ag	1	*Positive	1	100.00%
Sekisui Diagnostics OSOM	16	*Positive	16	100.00%
Sekisui Diagnostics OSOM Ultra Strep	13	*Positive	13	100.00%
Henry Schein OneStep+ Rapid Strep Dip	10	*Positive	10	100.00%
Immunostics Immuno/Strep A Det	2	*Positive	2	100.00%
Inverness Signify Strep A Dipstick	1	*Positive	1	100.00%
McKesson Strep A Cas	2	*Positive	2	100.00%
McKesson Strep A Dip	5	*Positive	5	100.00%
Permaxim Rediscreen Strep A Rapid	1	*Positive	1	100.00%
PSS Rapid Test Strip Strep A	1	*Positive	1	100.00%
Quidel QuickVue Dipstick Strep	10	*Positive	10	100.00%
Quidel QuickVue In-Line - 25 Test	10	*Positive	10	100.00%
Quidel QuickVue+ Strep A	17	*Positive	17	100.00%
SP (Cardinal) Rapid Strep Cassette	1	*Positive	1	100.00%
SP Brand (Cardinal) Rapid Strep A Dipstick	2	*Positive	2	100.00%
Stanbio QuStick Strep A	2	*Positive	2	100.00%
Consult Diag Strep A Dipstick	18	*Positive	18	100.00%
Meridian ImmunoCard STAT Strep A	1	*Positive	1	100.00%



Quidel Sofia/Sofia 2 Strep A FIA	8	*Positive	8	100.00%
BD Veritor Strep A	6	*Positive	6	100.00%
Quidel Solana	6	*Positive	6	100.00%
Cepheid GeneXpert - Xpert Strep A	19	*Positive	19	100.00%
		Strep A Ag Det. (International)		
RS-11		* Denotes Correct Response		
AG: All Methods	1			
		*Negative	1	100.00%
Alere Clearview Exact II Cas	1	*Negative	1	100.00%
RS-12		* Denotes Correct Response		
AG: All Methods	1			
		*Positive	1	100.00%
Alere Clearview Exact II Cas	1	*Positive	1	100.00%
691 - Group A Strep Antigen (5 specimens)	Labs	Response	Number	Percent
		Strep A Ag Det.		
RS-11		* Denotes Correct Response		
AG: All Methods	160			
		*Negative	160	100.00%
Abbott Signify Strep A - Waived	2	*Negative	2	100.00%
BD Check Strep A (Dip)	2	*Negative	2	100.00%
Beck Coulter ICON DS Strep	1	*Negative	1	100.00%
Alere Acceava Strep A Dipstick	3	*Negative	3	100.00%
Fisher HC Sure-Vue Select Strep Ag	1	*Negative	1	100.00%
Sekisui Diagnostics OSOM	16	*Negative	16	100.00%
Sekisui Diagnostics OSOM Ultra Strep	13	*Negative	13	100.00%
Henry Schein OneStep+ Rapid Strep Dip	10	*Negative	10	100.00%
Immunostics Immuno/Strep A Det	2	*Negative	2	100.00%
Inverness Signify Strep A Dipstick	1	*Negative	1	100.00%
McKesson Strep A Cas	2	*Negative	2	100.00%
McKesson Strep A Dip	5	*Negative	5	100.00%
Permaxim Rediscreen Strep A Rapid	1	*Negative	1	100.00%
PSS Rapid Test Strip Strep A	1	*Negative	1	100.00%
Quidel QuickVue Dipstick Strep	10	*Negative	10	100.00%
Quidel QuickVue In-Line - 25 Test	10	*Negative	10	100.00%
Quidel QuickVue+ Strep A	17	*Negative	17	100.00%
SP (Cardinal) Rapid Strep Cassette	1	*Negative	1	100.00%
SP Brand (Cardinal) Rapid Strep A Dipstick	2	*Negative	2	100.00%
Stanbio QuStick Strep A	2	*Negative	2	100.00%
Consult Diag Strep A Dipstick	18	*Negative	18	100.00%
Meridian ImmunoCard STAT Strep A	1	*Negative	1	100.00%
Quidel Sofia/Sofia 2 Strep A FIA	8	*Negative	8	100.00%
BD Veritor Strep A	6	*Negative	6	100.00%
Quidel Solana	6	*Negative	6	100.00%
Cepheid GeneXpert - Xpert Strep A	19	*Negative	19	100.00%
RS-12		* Denotes Correct Response		
AG: All Methods	160			
		*Positive	160	100.00%
Abbott Signify Strep A - Waived	2	*Positive	2	100.00%
BD Check Strep A (Dip)	2	*Positive	2	100.00%
Beck Coulter ICON DS Strep	1	*Positive	1	100.00%
Alere Acceava Strep A Dipstick	3	*Positive	3	100.00%
Fisher HC Sure-Vue Select Strep Ag	1	*Positive	1	100.00%
Sekisui Diagnostics OSOM	16	*Positive	16	100.00%
Sekisui Diagnostics OSOM Ultra Strep	13	*Positive	13	100.00%
Henry Schein OneStep+ Rapid Strep Dip	10	*Positive	10	100.00%
Immunostics Immuno/Strep A Det	2	*Positive	2	100.00%
Inverness Signify Strep A Dipstick	1	*Positive	1	100.00%
McKesson Strep A Cas	2	*Positive	2	100.00%
McKesson Strep A Dip	5	*Positive	5	100.00%
Permaxim Rediscreen Strep A Rapid	1	*Positive	1	100.00%
PSS Rapid Test Strip Strep A	1	*Positive	1	100.00%
Quidel QuickVue Dipstick Strep	10	*Positive	10	100.00%
Quidel QuickVue In-Line - 25 Test	10	*Positive	10	100.00%
Quidel QuickVue+ Strep A	17	*Positive	17	100.00%
SP (Cardinal) Rapid Strep Cassette	1	*Positive	1	100.00%
SP Brand (Cardinal) Rapid Strep A Dipstick	2	*Positive	2	100.00%
Stanbio QuStick Strep A	2	*Positive	2	100.00%
Consult Diag Strep A Dipstick	18	*Positive	18	100.00%
Meridian ImmunoCard STAT Strep A	1	*Positive	1	100.00%



Quidel Sofia/Sofia 2 Strep A FIA	8	*Positive	8	100.00%
BD Veritor Strep A	6	*Positive	6	100.00%
Quidel Solana	6	*Positive	6	100.00%
Cepheid GeneXpert - Xpert Strep A	19	*Positive	19	100.00%
RS-13		* Denotes Correct Response		
AG: All Methods	136			
		*Negative	135	99.26%
		Positive	1	0.74%
Abbott Signify Strep A - Waived	2	*Negative	2	100.00%
BD Check Strep A (Dip)	2	*Negative	2	100.00%
Beck Coulter ICON DS Strep	1	*Negative	1	100.00%
Alere Acceava Strep A Dipstick	3	*Negative	3	100.00%
Fisher HC Sure-Vue Select Strep Ag	1	*Negative	1	100.00%
Sekisui Diagnostics OSOM	16	*Negative	16	100.00%
Sekisui Diagnostics OSOM Ultra Strep	12	*Negative	12	100.00%
Henry Schein OneStep+ Rapid Strep Dip	10	*Negative	10	100.00%
Immunostics Immuno/Strep A Det	2	*Negative	2	100.00%
Inverness Signify Strep A Dipstick	1	*Negative	1	100.00%
McKesson Strep A Cas	1	*Negative	1	100.00%
McKesson Strep A Dip	5	*Negative	5	100.00%
Permaxim Rediscreen Strep A Rapid	1	*Negative	1	100.00%
PSS Rapid Test Strip Strep A	1	*Negative	1	100.00%
Quidel QuickVue Dipstick Strep	9	*Negative	9	100.00%
Quidel QuickVue In-Line - 25 Test	9	*Negative	8	88.89%
Quidel QuickVue+ Strep A	16	*Negative	16	100.00%
SP (Cardinal) Rapid Strep Cassette	1	*Negative	1	100.00%
SP Brand (Cardinal) Rapid Strep A Dipstick	2	*Negative	2	100.00%
Consult Diag Strep A Dipstick	13	*Negative	13	100.00%
Meridian ImmunoCard STAT Strep A	1	*Negative	1	100.00%
Quidel Sofia/Sofia 2 Strep A FIA	5	*Negative	5	100.00%
BD Veritor Strep A	1	*Negative	1	100.00%
Quidel Solana	4	*Negative	4	100.00%
Cepheid GeneXpert - Xpert Strep A	17	*Negative	17	100.00%
Quidel QuickVue In-Line - 25 Test	9	Positive	1	11.11%
RS-14		* Denotes Correct Response		
AG: All Methods	136			
		*Positive	136	100.00%
Abbott Signify Strep A - Waived	2	*Positive	2	100.00%
BD Check Strep A (Dip)	2	*Positive	2	100.00%
Beck Coulter ICON DS Strep	1	*Positive	1	100.00%
Alere Acceava Strep A Dipstick	3	*Positive	3	100.00%
Fisher HC Sure-Vue Select Strep Ag	1	*Positive	1	100.00%
Sekisui Diagnostics OSOM	16	*Positive	16	100.00%
Sekisui Diagnostics OSOM Ultra Strep	12	*Positive	12	100.00%
Henry Schein OneStep+ Rapid Strep Dip	10	*Positive	10	100.00%
Immunostics Immuno/Strep A Det	2	*Positive	2	100.00%
Inverness Signify Strep A Dipstick	1	*Positive	1	100.00%
McKesson Strep A Cas	1	*Positive	1	100.00%
McKesson Strep A Dip	5	*Positive	5	100.00%
Permaxim Rediscreen Strep A Rapid	1	*Positive	1	100.00%
PSS Rapid Test Strip Strep A	1	*Positive	1	100.00%
Quidel QuickVue Dipstick Strep	9	*Positive	9	100.00%
Quidel QuickVue In-Line - 25 Test	9	*Positive	9	100.00%
Quidel QuickVue+ Strep A	16	*Positive	16	100.00%
SP (Cardinal) Rapid Strep Cassette	1	*Positive	1	100.00%
SP Brand (Cardinal) Rapid Strep A Dipstick	2	*Positive	2	100.00%
Consult Diag Strep A Dipstick	13	*Positive	13	100.00%
Meridian ImmunoCard STAT Strep A	1	*Positive	1	100.00%
Quidel Sofia/Sofia 2 Strep A FIA	5	*Positive	5	100.00%
BD Veritor Strep A	1	*Positive	1	100.00%
Quidel Solana	4	*Positive	4	100.00%
Cepheid GeneXpert - Xpert Strep A	17	*Positive	17	100.00%
RS-15		* Denotes Correct Response		
AG: All Methods	136			
		*Positive	134	98.53%
		Negative	2	1.47%
Abbott Signify Strep A - Waived	2	*Positive	1	50.00%
BD Check Strep A (Dip)	2	*Positive	2	100.00%
Beck Coulter ICON DS Strep	1	*Positive	1	100.00%



Alere Acceava Strep A Dipstick	3	*Positive	3	100.00%
Fisher HC Sure-Vue Select Strep Ag	1	*Positive	1	100.00%
Sekisui Diagnostics OSOM	16	*Positive	16	100.00%
Sekisui Diagnostics OSOM Ultra Strep	12	*Positive	11	91.67%
Henry Schein OneStep+ Rapid Strep Dip	10	*Positive	10	100.00%
Immunostics Immuno/Strep A Det	2	*Positive	2	100.00%
Inverness Signify Strep A Dipstick	1	*Positive	1	100.00%
McKesson Strep A Cas	1	*Positive	1	100.00%
McKesson Strep A Dip	5	*Positive	5	100.00%
Permaxim Rediscreen Strep A Rapid	1	*Positive	1	100.00%
PSS Rapid Test Strip Strep A	1	*Positive	1	100.00%
Quidel QuickVue Dipstick Strep	9	*Positive	9	100.00%
Quidel QuickVue In-Line - 25 Test	9	*Positive	9	100.00%
Quidel QuickVue+ Strep A	16	*Positive	16	100.00%
SP (Cardinal) Rapid Strep Cassette	1	*Positive	1	100.00%
SP Brand (Cardinal) Rapid Strep A Dipstick	2	*Positive	2	100.00%
Consult Diag Strep A Dipstick	13	*Positive	13	100.00%
Meridian ImmunoCard STAT Strep A	1	*Positive	1	100.00%
Quidel Sofia/Sofia 2 Strep A FIA	5	*Positive	5	100.00%
BD Veritor Strep A	1	*Positive	1	100.00%
Quidel Solana	4	*Positive	4	100.00%
Cepheid GeneXpert - Xpert Strep A	17	*Positive	17	100.00%
Abbott Signify Strep A - Waived	2	Negative	1	50.00%
Sekisui Diagnostics OSOM Ultra Strep	12	Negative	1	8.33%
		Strep A Ag Det. (Waived)		
RS-11		* Denotes Correct Response		
AG: All Methods	879			
		*Negative	869	98.86%
		Positive	10	1.14%
Abbott Signify Strep A - Waived	6	*Negative	6	100.00%
Beck Coulter ICON DS Strep	1	*Negative	1	100.00%
Alere Acceava Strep A Dipstick	115	*Negative	113	98.26%
Fisher HC Sure-Vue Strep A (Waived)	9	*Negative	9	100.00%
Sekisui Diagnostics OSOM	38	*Negative	38	100.00%
Sekisui Diagnostics OSOM Ultra Strep	45	*Negative	45	100.00%
Henry Schein OneStep + Rapid Strep Cass	3	*Negative	3	100.00%
Henry Schein OneStep+ Rapid Strep Dip	73	*Negative	71	97.26%
Innovacon One Step Strep Ag	1	*Negative	1	100.00%
Instant Technologies i Strep	1	*Negative	1	100.00%
Inverness Signify Strep A Dipstick	7	*Negative	7	100.00%
McKesson Strep A Cas	4	*Negative	4	100.00%
McKesson Strep A Dip	31	*Negative	31	100.00%
McKesson Strep A Test - Twist	2	*Negative	2	100.00%
Permaxim Rediscreen Strep A Rapid	1	*Negative	1	100.00%
Polymedco Poly Stat Strep A (W)	1	*Negative	1	100.00%
Quidel QuickVue Dipstick Strep	76	*Negative	76	100.00%
Quidel QuickVue In-Line - 25 Test	63	*Negative	62	98.41%
Quidel QuickVue+ Strep A	33	*Negative	32	96.97%
Select Diag Strep A Dip	4	*Negative	4	100.00%
SP (Cardinal) Rapid Strep Cassette	5	*Negative	5	100.00%
SP Brand (Cardinal) Rapid Strep A Dipstick	83	*Negative	83	100.00%
Alere Clearview Exact II Dip	35	*Negative	34	97.14%
Sekisui Diagnostics OSOM Ultra	7	*Negative	7	100.00%
Consult Diag Strep A Dipstick	121	*Negative	120	99.17%
Consult Diag Strep A Twist	1	*Negative	1	100.00%
MooreBrand Strep A Dipstick	6	*Negative	5	83.33%
NDC ProAdvantage Strep Strips	1	*Negative	1	100.00%
Meridian ImmunoCard STAT Strep A	1	*Negative	1	100.00%
LifeSign Bio Strep A Test STrip	1	*Negative	1	100.00%
Quidel Sofia/Sofia 2 Strep A FIA	31	*Negative	30	96.77%
BD Veritor Strep A	38	*Negative	38	100.00%
Alere i Strep A	11	*Negative	11	100.00%
Medline Strep A Dipstick	3	*Negative	3	100.00%
Roche Cobas Liat Strep A	21	*Negative	21	100.00%
Alere Acceava Strep A Dipstick	115	Positive	2	1.74%
Henry Schein OneStep+ Rapid Strep Dip	73	Positive	2	2.74%
Quidel QuickVue In-Line - 25 Test	63	Positive	1	1.59%
Quidel QuickVue+ Strep A	33	Positive	1	3.03%
Alere Clearview Exact II Dip	35	Positive	1	2.86%



Consult Diag Strep A Dipstick	121	Positive	1	0.83%
MooreBrand Strep A Dipstick	6	Positive	1	16.67%
Quidel Sofia/Sofia 2 Strep A FIA	31	Positive	1	3.23%
RS-12		* Denotes Correct Response		
AG: All Methods	55			
		*Positive	54	98.18%
		Negative	1	1.82%
Alere Aceava Strep A Dipstick	2	*Positive	2	100.00%
Sekisui Diagnostics OSOM	8	*Positive	8	100.00%
Sekisui Diagnostics OSOM Ultra Strep	4	*Positive	4	100.00%
Henry Schein OneStep+ Rapid Strep Dip	1	*Positive	1	100.00%
Innovacon One Step Strep Ag	1	*Positive	1	100.00%
McKesson Strep A Dip	1	*Positive	1	100.00%
Quidel QuickVue Dipstick Strep	1	*Positive	1	100.00%
Quidel QuickVue In-Line - 25 Test	3	*Positive	3	100.00%
Quidel QuickVue+ Strep A	1	*Positive	1	100.00%
SP Brand (Cardinal) Rapid Strep A Dipstick	4	*Positive	4	100.00%
Consult Diag Strep A Dipstick	7	*Positive	7	100.00%
BD Veritor Strep A	1	*Positive	1	100.00%
Roche Cobas Liat Strep A	20	*Positive	20	100.00%
MooreBrand Strep A Dipstick	1	Negative	1	100.00%
RS-13		* Denotes Correct Response		
AG: All Methods	1			
		*Negative	1	100.00%
Quidel QuickVue In-Line - 25 Test	1	*Negative	1	100.00%
RS-14		* Denotes Correct Response		
AG: All Methods	1			
		*Positive	1	100.00%
Quidel QuickVue In-Line - 25 Test	1	*Positive	1	100.00%
RS-15		* Denotes Correct Response		
AG: All Methods	1			
		*Positive	1	100.00%
Quidel QuickVue In-Line - 25 Test	1	*Positive	1	100.00%
		Strep A Ag Det. (International)		
RS-11		* Denotes Correct Response		
AG: All Methods	1			
		*Negative	1	100.00%
Alere Clearview Exact II Cas	1	*Negative	1	100.00%
RS-12		* Denotes Correct Response		
AG: All Methods	1			
		*Positive	1	100.00%
Alere Clearview Exact II Cas	1	*Positive	1	100.00%
RS-13		* Denotes Correct Response		
AG: All Methods	1			
		*Negative	1	100.00%
Alere Clearview Exact II Cas	1	*Negative	1	100.00%
RS-14		* Denotes Correct Response		
AG: All Methods	1			
		*Positive	1	100.00%
Alere Clearview Exact II Cas	1	*Positive	1	100.00%
RS-15		* Denotes Correct Response		
AG: All Methods	1			
		*Positive	1	100.00%
Alere Clearview Exact II Cas	1	*Positive	1	100.00%
739 - Sickle Cell Screen	Labs	Response	Number	Percent
		Sickle Cell Screen		
SS-3		* Denotes Correct Response		
AG: All Methods	13			
		*Positive	13	100.00%
Pacific Hemostasis Sickie Screen	3	*Positive	3	100.00%
Streck SickieDex	10	*Positive	10	100.00%
679 - Syphilis Serology	Labs	Response	Number	Percent
		RPR		
SY-11		* Denotes Correct Response		
AG: All Methods	9			
		*Nonreactive	9	100.00%
ASI RPR	2	*Nonreactive	2	100.00%
Becton Dickinson BBL	3	*Nonreactive	3	100.00%



New Horizons	1	*Nonreactive	1	100.00%
Stanbio Quicktest	1	*Nonreactive	1	100.00%
Wampole Impact	1	*Nonreactive	1	100.00%
Intl ASI RPR	1	*Nonreactive	1	100.00%
SY-12		* Denotes Correct Response		
AG: All Methods	9			
		*Nonreactive	9	100.00%
ASI RPR	2	*Nonreactive	2	100.00%
Becton Dickinson BBL	3	*Nonreactive	3	100.00%
New Horizons	1	*Nonreactive	1	100.00%
Stanbio Quicktest	1	*Nonreactive	1	100.00%
Wampole Impact	1	*Nonreactive	1	100.00%
Intl ASI RPR	1	*Nonreactive	1	100.00%
SY-13		* Denotes Correct Response		
AG: All Methods	9			
		*Reactive	9	100.00%
ASI RPR	2	*Reactive	2	100.00%
Becton Dickinson BBL	3	*Reactive	3	100.00%
New Horizons	1	*Reactive	1	100.00%
Stanbio Quicktest	1	*Reactive	1	100.00%
Wampole Impact	1	*Reactive	1	100.00%
Intl ASI RPR	1	*Reactive	1	100.00%
SY-14		* Denotes Correct Response		
AG: All Methods	9			
		*Reactive	8	88.89%
		*Weakly Reactive	1	11.11%
ASI RPR	2	*Reactive	2	100.00%
Becton Dickinson BBL	3	*Reactive	3	100.00%
New Horizons	1	*Reactive	1	100.00%
Stanbio Quicktest	1	*Reactive	1	100.00%
Intl ASI RPR	1	*Reactive	1	100.00%
Wampole Impact	1	*Weakly Reactive	1	100.00%
SY-15		* Denotes Correct Response		
AG: All Methods	9			
		*Reactive	9	100.00%
ASI RPR	2	*Reactive	2	100.00%
Becton Dickinson BBL	3	*Reactive	3	100.00%
New Horizons	1	*Reactive	1	100.00%
Stanbio Quicktest	1	*Reactive	1	100.00%
Wampole Impact	1	*Reactive	1	100.00%
Intl ASI RPR	1	*Reactive	1	100.00%
		RPR/VDRL (titer)		
SY-11		* Denotes Correct Response		
AG: All Methods	3			
		*Nonreactive	3	100.00%
ASI RPR	1	*Nonreactive	1	100.00%
Becton Dickinson BBL	1	*Nonreactive	1	100.00%
New Horizons	1	*Nonreactive	1	100.00%
SY-12		* Denotes Correct Response		
AG: All Methods	3			
		*Nonreactive	3	100.00%
ASI RPR	1	*Nonreactive	1	100.00%
Becton Dickinson BBL	1	*Nonreactive	1	100.00%
New Horizons	1	*Nonreactive	1	100.00%
SY-13		* Denotes Correct Response		
AG: All Methods	3			
		*1:2	2	66.67%
		*1:4	1	33.33%
ASI RPR	1	*1:2	1	100.00%
New Horizons	1	*1:2	1	100.00%
Becton Dickinson BBL	1	*1:4	1	100.00%
SY-14		* Denotes Correct Response		
AG: All Methods	3			
		*1:2	1	33.33%
		*1:4	1	33.33%
		*1:8	1	33.33%
New Horizons	1	*1:2	1	100.00%
ASI RPR	1	*1:4	1	100.00%



Becton Dickinson BBL	1	*1:8	1	100.00%
SY-15		* Denotes Correct Response		
AG: All Methods	3			
		*1:2	2	66.67%
		*1:4	1	33.33%
ASI RPR	1	*1:2	1	100.00%
New Horizons	1	*1:2	1	100.00%
Becton Dickinson BBL	1	*1:4	1	100.00%
		RPR (International)		
SY-11		* Denotes Correct Response		
AG: All Methods	3			
		*Nonreactive	3	100.00%
BD Macro-Vue RPR	1	*Nonreactive	1	100.00%
Fisher HC Sure-Vue RPR	1	*Nonreactive	1	100.00%
Wampole Impact	1	*Nonreactive	1	100.00%
SY-12		* Denotes Correct Response		
AG: All Methods	3			
		*Nonreactive	3	100.00%
BD Macro-Vue RPR	1	*Nonreactive	1	100.00%
Fisher HC Sure-Vue RPR	1	*Nonreactive	1	100.00%
Wampole Impact	1	*Nonreactive	1	100.00%
SY-13		* Denotes Correct Response		
AG: All Methods	3			
		*Reactive	3	100.00%
BD Macro-Vue RPR	1	*Reactive	1	100.00%
Fisher HC Sure-Vue RPR	1	*Reactive	1	100.00%
Wampole Impact	1	*Reactive	1	100.00%
SY-14		* Denotes Correct Response		
AG: All Methods	3			
		*Reactive	3	100.00%
BD Macro-Vue RPR	1	*Reactive	1	100.00%
Fisher HC Sure-Vue RPR	1	*Reactive	1	100.00%
Wampole Impact	1	*Reactive	1	100.00%
SY-15		* Denotes Correct Response		
AG: All Methods	3			
		*Reactive	3	100.00%
BD Macro-Vue RPR	1	*Reactive	1	100.00%
Fisher HC Sure-Vue RPR	1	*Reactive	1	100.00%
Wampole Impact	1	*Reactive	1	100.00%
776 - Throat Culture;780 - Throat & Genital Cultures	Labs	Response	Number	Percent
		Throat Culture		
TC-11		* Denotes Correct Response		
AG: All Methods	64			
		*Beta hemolysis, referred for ID/confirmation	4	6.25%
		*Positive for Group A Streptococcus	29	45.31%
		*Presumptive Positive for Group A Streptococcus	23	35.94%
		*Presumptive Streptococcus species	2	3.13%
		Negative for Group A Streptococcus	5	7.81%
		Streptococcus, beta-hemolytic, not Group A	1	1.56%
Cultures	64	*Beta hemolysis, referred for ID/confirmation	4	6.25%
Cultures		*Positive for Group A Streptococcus	29	45.31%
Cultures		*Presumptive Positive for Group A Streptococcus	23	35.94%
Cultures		*Presumptive Streptococcus species	2	3.13%
Cultures		Negative for Group A Streptococcus	5	7.81%
Cultures		Streptococcus, beta-hemolytic, not Group A	1	1.56%
TC-12		* Denotes Correct Response		
AG: All Methods	64			
		*Growth, referred for identification	1	1.56%
		*Negative for Group A Streptococcus	46	71.88%
		*No beta hemolysis	5	7.81%
		*No Growth - Strep Selective medium	2	3.13%
		No Growth	1	1.56%
		Positive for Group A Streptococcus	3	4.69%
		Presumptive Positive for Group A Streptococcus	6	9.38%
Cultures	64	*Growth, referred for identification	1	1.56%
Cultures		*Negative for Group A Streptococcus	46	71.88%
Cultures		*No beta hemolysis	5	7.81%
Cultures		*No Growth - Strep Selective medium	2	3.13%



Cultures		No Growth	1	1.56%
Cultures		Positive for Group A Streptococcus	3	4.69%
Cultures		Presumptive Positive for Group A Streptococcus	6	9.38%
TC-13		* Denotes Correct Response		
AG: All Methods	65			
		*Beta hemolysis, referred for ID/confirmation	5	7.69%
		*Growth, referred for identification	1	1.54%
		*Negative for Group A Streptococcus	48	73.85%
		*Streptococcus, beta-hemolytic, not Group A	3	4.62%
		No beta hemolysis	1	1.54%
		Positive for Group A Streptococcus	3	4.62%
		Presumptive Positive for Group A Streptococcus	4	6.15%
Cultures	65	*Beta hemolysis, referred for ID/confirmation	5	7.69%
Cultures		*Growth, referred for identification	1	1.54%
Cultures		*Negative for Group A Streptococcus	48	73.85%
Cultures		*Streptococcus, beta-hemolytic, not Group A	3	4.62%
Cultures		No beta hemolysis	1	1.54%
Cultures		Positive for Group A Streptococcus	3	4.62%
Cultures		Presumptive Positive for Group A Streptococcus	4	6.15%
TC-14		* Denotes Correct Response		
AG: All Methods	81			
		*Beta hemolysis, referred for ID/confirmation	6	7.41%
		*Positive for Group A Streptococcus	42	51.85%
		*Presumptive Positive for Group A Streptococcus	25	30.86%
		*Presumptive Streptococcus species	2	2.47%
		*Streptococcus, beta-hemolytic, Group A	1	1.23%
		Negative for Group A Streptococcus	5	6.17%
Cultures	81	*Beta hemolysis, referred for ID/confirmation	6	7.41%
Cultures		*Positive for Group A Streptococcus	42	51.85%
Cultures		*Presumptive Positive for Group A Streptococcus	25	30.86%
Cultures		*Presumptive Streptococcus species	2	2.47%
Cultures		*Streptococcus, beta-hemolytic, Group A	1	1.23%
Cultures		Negative for Group A Streptococcus	5	6.17%
TC-15		* Denotes Correct Response		
AG: All Methods	81			
		*Negative for Group A Streptococcus	57	70.37%
		*No beta hemolysis	8	9.88%
		*No Growth - Strep Selective medium	12	14.81%
		Beta hemolysis, referred for ID/confirmation	2	2.47%
		No Growth	1	1.23%
		Presumptive Positive for Group A Streptococcus	1	1.23%
Cultures	81	*Negative for Group A Streptococcus	57	70.37%
Cultures		*No beta hemolysis	8	9.88%
Cultures		*No Growth - Strep Selective medium	12	14.81%
Cultures		Beta hemolysis, referred for ID/confirmation	2	2.47%
Cultures		No Growth	1	1.23%
Cultures		Presumptive Positive for Group A Streptococcus	1	1.23%
640 - Urinalysis;642 - Urinalysis/hCG;639 - Complete Urinalysis Package;601 - Waived/PPM Tests II;609 - Waived/PPM Tests I;647 - Mini-Waived Package	Labs	Response	Number	Percent
		Urine Sp Gravity		
UA-3		* Denotes Correct Response		
PG: Siemens Clinitek 10 / 100	12			
		*1.010	2	16.67%
		*1.015	3	25.00%
		*1.020	5	41.67%
		1.005	1	8.33%
		1.025	1	8.33%
UA-3		* Denotes Correct Response		
PG: Siemens Clinitek 50	80			
		*1.015	21	26.25%
		*1.020	45	56.25%
		*1.025	13	16.25%
		1.010	1	1.25%
UA-3		* Denotes Correct Response		
PG: Siemens Clinitek 500	6			
		*1.020	5	83.33%
		*1.025	1	16.67%
UA-3		* Denotes Correct Response		



PG: Siemens Clinitek Status	613			
		*1.020	488	79.61%
		*1.025	115	18.76%
		1.010	3	0.49%
		1.015	3	0.49%
		1.030	4	0.65%
UA-3		* Denotes Correct Response		
PG: Siemens Clinitek Status / Siemens - Plus Sys / Siemens Urine Dipstick	17			
		*1.015	2	11.76%
		*1.020	11	64.71%
		*1.025	4	23.53%
UA-3		* Denotes Correct Response		
PG: Siemens Multistix - Visual	214			
		*1.010	65	30.37%
		*1.015	78	36.45%
		*1.020	52	24.30%
		1.000	1	0.47%
		1.005	6	2.80%
		1.025	11	5.14%
		1.030	1	0.47%
UA-3		* Denotes Correct Response		
PG: Henry Schein Urispec	8			
		*1.010	3	37.50%
		*1.015	4	50.00%
		*1.020	1	12.50%
UA-3		* Denotes Correct Response		
PG: Consult Diagnostics Urine Reagent Strips	8			
		*1.010	1	12.50%
		*1.015	7	87.50%
UA-3		* Denotes Correct Response		
PG: Diag Test Group Clarity Urochek	30			
		*1.015	23	76.67%
		*1.020	2	6.67%
		1.000	1	3.33%
		1.005	2	6.67%
		1.025	2	6.67%
UA-3		* Denotes Correct Response		
PG: Roche Chemstrips - Visual	116			
		*1.010	81	69.83%
		*1.015	19	16.38%
		*1.020	6	5.17%
		1.000	2	1.72%
		1.005	6	5.17%
		1.025	2	1.72%
UA-3		* Denotes Correct Response		
PG: Roche Urisys	10			
		*1.010	6	60.00%
		*1.015	1	10.00%
		*1.020	2	20.00%
		1.025	1	10.00%
UA-3		* Denotes Correct Response		
PG: Roche Urisys / Roche - Conventional / Roche Urine Dipstick	1			
		*1.015	1	100.00%
UA-3		* Denotes Correct Response		
PG: Roche Urisys / Roche - Arbitrary / Roche Urine Arbitrary Result	2			
		*1.010	2	100.00%
UA-3		* Denotes Correct Response		
PG: Medi-Lab Performance Urine Strip - Visual	7			
		*1.010	2	28.57%
		*1.015	5	71.43%
UA-3		* Denotes Correct Response		
PG: Siemens Clinitek Advantus	21			
		*1.015	6	28.57%
		*1.020	14	66.67%
		1.025	1	4.76%
UA-3		* Denotes Correct Response		
PG: Consult Diag Urinalysis Test Strips	7			



		*1.015	5	71.43%
		*1.020	1	14.29%
		*1.025	1	14.29%
UA-3		* Denotes Correct Response		
PG: Consult Diag Urine Analyzer	15			
		*1.015	12	80.00%
		*1.020	3	20.00%
UA-3		* Denotes Correct Response		
PG: Consult Diagnostics 120 Urine Analyzer	44			
		*1.015	42	95.45%
		*1.020	2	4.55%
UA-3		* Denotes Correct Response		
SG: Select Siemens Clinitek/Visual	947			
		*1.010	71	7.50%
		*1.015	111	11.72%
		*1.020	609	64.31%
		*1.025	143	15.10%
		1.000	1	0.11%
		1.005	7	0.74%
		1.010	71	7.50%
		1.015	111	11.72%
		1.025	143	15.10%
		1.030	5	0.53%
UA-3		* Denotes Correct Response		
SG: Select Siemens Clinitek/Visual / Siemens - Plus Sys / Siemens Urine Dipstick	18			
		*1.015	2	11.11%
		*1.020	12	66.67%
		*1.025	4	22.22%
UA-3		* Denotes Correct Response		
SG: Select Roche Urinalysis Inst/Visual	132			
		*1.010	90	68.18%
		*1.015	21	15.91%
		*1.020	9	6.82%
		1.000	2	1.52%
		1.005	7	5.30%
		1.025	3	2.27%
UA-3		* Denotes Correct Response		
SG: Select Consult Diag Urinalysis Inst/Visual	30			
		*1.010	1	3.33%
		*1.015	24	80.00%
		*1.020	4	13.33%
		*1.025	1	3.33%
UA-3		* Denotes Correct Response		
SG: Select McKesson Urinalysis Inst/Visual	47			
		*1.015	45	95.74%
		*1.020	2	4.26%
UA-3		* Denotes Correct Response		
SG: Select Henry Schein Urinalysis Inst/Visual	11			
		*1.010	3	27.27%
		*1.015	5	45.45%
		*1.020	3	27.27%
UA-3		* Denotes Correct Response		
SG: Select Medi-Lab Performance UA Inst/Visual	8			
		*1.010	2	25.00%
		*1.015	6	75.00%
UA-3		* Denotes Correct Response		
MG: Method Group	1229			
		*1.010	173	14.08%
		*1.015	247	20.10%
		*1.020	634	51.59%
		*1.025	149	12.12%
		1.000	4	0.33%
		1.005	16	1.30%
		1.010	173	14.08%
		1.015	247	20.10%
		1.025	149	12.12%
		1.030	6	0.49%
UA-3		* Denotes Correct Response		



MG: Method Group	1			
		*1.015	1	100.00%
UA-3		* Denotes Correct Response		
MG: Method Group	18			
		*1.015	2	11.11%
		*1.020	12	66.67%
		*1.025	4	22.22%
UA-3		* Denotes Correct Response		
MG: Method Group	2			
		*1.010	2	100.00%
		Urine Bilirubin		
UA-3		* Denotes Correct Response		
PG: Siemens Clinitek 10 / 100	12			
		*Large (3+)	7	58.33%
		*Moderate (2+)	5	41.67%
UA-3		* Denotes Correct Response		
PG: Siemens Clinitek 50	80			
		*Large (3+)	29	36.25%
		*Moderate (2+)	48	60.00%
		*Trace/Small (1+)	2	2.50%
		Negative	1	1.25%
UA-3		* Denotes Correct Response		
PG: Siemens Clinitek 500	5			
		*Large (3+)	1	20.00%
		*Moderate (2+)	4	80.00%
UA-3		* Denotes Correct Response		
PG: Siemens Clinitek Status	606			
		*Large (3+)	197	32.51%
		*Moderate (2+)	400	66.01%
		Negative	1	0.17%
		Trace/Small (1+)	8	1.32%
UA-3		* Denotes Correct Response		
PG: Siemens Clinitek Status / Siemens - Plus Sys / Siemens Urine Dipstick	17			
		*1+	1	5.88%
		*2+	12	70.59%
		*3+	4	23.53%
UA-3		* Denotes Correct Response		
PG: Siemens Multistix - Visual	207			
		*Large (3+)	113	54.59%
		*Moderate (2+)	62	29.95%
		*Trace/Small (1+)	23	11.11%
		Negative	9	4.35%
UA-3		* Denotes Correct Response		
PG: Henry Schein Urispec	8			
		*Large (3+)	5	62.50%
		*Moderate (2+)	2	25.00%
		*Trace/Small (1+)	1	12.50%
UA-3		* Denotes Correct Response		
PG: Consult Diagnostics Urine Reagent Strips	8			
		*Large (3+)	2	25.00%
		*Moderate (2+)	2	25.00%
		*Trace/Small (1+)	4	50.00%
UA-3		* Denotes Correct Response		
PG: Diag Test Group Clarity Urochek	29			
		*Large (3+)	4	13.79%
		*Moderate (2+)	17	58.62%
		*Trace/Small (1+)	8	27.59%
UA-3		* Denotes Correct Response		
PG: Roche Chemstrips - Visual	117			
		*Large (3+)	58	49.57%
		*Moderate (2+)	52	44.44%
		Negative	1	0.85%
		Trace/Small (1+)	6	5.13%
UA-3		* Denotes Correct Response		
PG: Roche Urisys	10			
		*Large (3+)	3	30.00%
		*Moderate (2+)	7	70.00%



UA-3		* Denotes Correct Response		
PG: Roche Urisys / Roche - Conventional / Roche Urine Dipstick	1			
		*1 mg/dL	1	100.00%
UA-3		* Denotes Correct Response		
PG: Roche Urisys / Roche - Arbitrary / Roche Urine Arbitrary Result	2			
		*++	2	100.00%
UA-3		* Denotes Correct Response		
PG: UriScan Reagent Strips	5			
		*Large (3+)	3	60.00%
		*Moderate (2+)	1	20.00%
		*Trace/Small (1+)	1	20.00%
UA-3		* Denotes Correct Response		
PG: Medi-Lab Performance Urine Strip - Visual	7			
		*Moderate (2+)	5	71.43%
		*Trace/Small (1+)	2	28.57%
UA-3		* Denotes Correct Response		
PG: Siemens Clinitek Advantus	21			
		*Moderate (2+)	21	100.00%
UA-3		* Denotes Correct Response		
PG: Consult Diag Urinalysis Test Strips	7			
		*Moderate (2+)	2	28.57%
		*Trace/Small (1+)	5	71.43%
UA-3		* Denotes Correct Response		
PG: Consult Diag Urine Analyzer	15			
		*Moderate (2+)	9	60.00%
		*Trace/Small (1+)	5	33.33%
		Large (3+)	1	6.67%
UA-3		* Denotes Correct Response		
PG: Consult Diagnostics 120 Urine Analyzer	44			
		*Moderate (2+)	11	25.00%
		*Trace/Small (1+)	33	75.00%
UA-3		* Denotes Correct Response		
SG: Select Siemens Clinitek/Visual	932			
		*Large (3+)	347	37.23%
		*Moderate (2+)	541	58.05%
		*Trace/Small (1+)	33	3.54%
		Negative	11	1.18%
		Trace/Small (1+)	33	3.54%
UA-3		* Denotes Correct Response		
SG: Select Siemens Clinitek/Visual / Siemens - Plus Sys / Siemens Urine Dipstick	18			
		*1+	1	5.56%
		*2+	13	72.22%
		*3+	4	22.22%
UA-3		* Denotes Correct Response		
SG: Select Roche Urinalysis Inst/Visual	133			
		*Large (3+)	63	47.37%
		*Moderate (2+)	61	45.86%
		*Trace/Small (1+)	7	5.26%
		Negative	2	1.50%
		Trace/Small (1+)	7	5.26%
UA-3		* Denotes Correct Response		
SG: Select Consult Diag Urinalysis Inst/Visual	30			
		*Large (3+)	3	10.00%
		*Moderate (2+)	13	43.33%
		*Trace/Small (1+)	14	46.67%
		Large (3+)	3	10.00%
UA-3		* Denotes Correct Response		
SG: Select McKesson Urinalysis Inst/Visual	47			
		*Large (3+)	1	2.13%
		*Moderate (2+)	12	25.53%
		*Trace/Small (1+)	33	70.21%
		Negative	1	2.13%
UA-3		* Denotes Correct Response		
SG: Select Henry Schein Urinalysis Inst/Visual	11			
		*Large (3+)	6	54.55%
		*Moderate (2+)	4	36.36%
		*Trace/Small (1+)	1	9.09%



UA-3		* Denotes Correct Response		
SG: Select Medi-Lab Performance UA Inst/Visual	8			
		*Moderate (2+)	5	62.50%
		*Trace/Small (1+)	3	37.50%
UA-3		* Denotes Correct Response		
MG: Method Group	1214			
		*Large (3+)	431	35.50%
		*Moderate (2+)	658	54.20%
		*Trace/Small (1+)	110	9.06%
		Large (3+)	431	35.50%
		Negative	15	1.24%
		Trace/Small (1+)	110	9.06%
UA-3		* Denotes Correct Response		
MG: Method Group	1			
		*1 mg/dL	1	100.00%
UA-3		* Denotes Correct Response		
MG: Method Group	18			
		*1+	1	5.56%
		*2+	13	72.22%
		*3+	4	22.22%
UA-3		* Denotes Correct Response		
MG: Method Group	2			
		*++	2	100.00%
		Urine Ketones		
UA-3		* Denotes Correct Response		
PG: Siemens Clinitek 10 / 100	12			
		*Large (3+) / 80 mg/dL	9	75.00%
		*Moderate (2+) / 40 mg/dL	2	16.67%
		Trace/5 mg/dL	1	8.33%
UA-3		* Denotes Correct Response		
PG: Siemens Clinitek 50	80			
		*>=160 mg/dL	22	27.50%
		*Large (3+) / 80 mg/dL	55	68.75%
		Moderate (2+) / 40 mg/dL	2	2.50%
		Negative	1	1.25%
UA-3		* Denotes Correct Response		
PG: Siemens Clinitek 500	6			
		*Large (3+) / 80 mg/dL	6	100.00%
UA-3		* Denotes Correct Response		
PG: Siemens Clinitek Status	612			
		*>=160 mg/dL	314	51.31%
		*Large (3+) / 80 mg/dL	274	44.77%
		Moderate (2+) / 40 mg/dL	21	3.43%
		Negative	2	0.33%
		Small (1+) / 15mg/dL	1	0.16%
UA-3		* Denotes Correct Response		
PG: Siemens Clinitek Status / Siemens - Plus Sys / Siemens Urine Dipstick	17			
		*3+	8	47.06%
		*4+	9	52.94%
UA-3		* Denotes Correct Response		
PG: Siemens Multistix - Visual	216			
		*Large (3+) / 80 mg/dL	208	96.30%
		Moderate (2+) / 40 mg/dL	8	3.70%
UA-3		* Denotes Correct Response		
PG: Henry Schein Urispec	8			
		*Large (3+) / 80 mg/dL	6	75.00%
		*Moderate (2+) / 40 mg/dL	2	25.00%
UA-3		* Denotes Correct Response		
PG: Consult Diagnostics Urine Reagent Strips	8			
		*Large (3+) / 80 mg/dL	6	75.00%
		*Moderate (2+) / 40 mg/dL	2	25.00%
UA-3		* Denotes Correct Response		
PG: Diag Test Group Clarity Urochek	30			
		*Large (3+) / 80 mg/dL	28	93.33%
		Moderate (2+) / 40 mg/dL	2	6.67%
UA-3		* Denotes Correct Response		
PG: Roche Chemstrips - Visual	118			



		*Large (3+) / 80 mg/dL	95	80.51%
		*Moderate (2+) / 40 mg/dL	22	18.64%
		Negative	1	0.85%
UA-3		* Denotes Correct Response		
PG: Roche Urisys	10			
		*Large (3+) / 80 mg/dL	9	90.00%
		Moderate (2+) / 40 mg/dL	1	10.00%
UA-3		* Denotes Correct Response		
PG: Roche Urisys / Roche - Conventional / Roche Urine Dipstick	1			
		*50 mg/dL	1	100.00%
UA-3		* Denotes Correct Response		
PG: Roche Urisys / Roche - Arbitrary / Roche Urine Arbitrary Result	2			
		*++	2	100.00%
UA-3		* Denotes Correct Response		
PG: UriScan Reagent Strips	5			
		*Large (3+) / 80 mg/dL	4	80.00%
		Moderate (2+) / 40 mg/dL	1	20.00%
UA-3		* Denotes Correct Response		
PG: Medi-Lab Performance Urine Strip - Visual	7			
		*Large (3+) / 80 mg/dL	7	100.00%
UA-3		* Denotes Correct Response		
PG: Siemens Clinitek Advantus	21			
		*Large (3+) / 80 mg/dL	17	80.95%
		Moderate (2+) / 40 mg/dL	4	19.05%
UA-3		* Denotes Correct Response		
PG: Consult Diag Urinalysis Test Strips	7			
		*Large (3+) / 80 mg/dL	7	100.00%
UA-3		* Denotes Correct Response		
PG: Consult Diag Urine Analyzer	15			
		*Large (3+) / 80 mg/dL	14	93.33%
		Moderate (2+) / 40 mg/dL	1	6.67%
UA-3		* Denotes Correct Response		
PG: Consult Diagnostics 120 Urine Analyzer	44			
		*Large (3+) / 80 mg/dL	44	100.00%
UA-3		* Denotes Correct Response		
SG: Select Siemens Clinitek/Visual	948			
		*>=160 mg/dL	336	35.44%
		*Large (3+) / 80 mg/dL	570	60.13%
		*Moderate (2+) / 40 mg/dL	37	3.90%
		Moderate (2+) / 40 mg/dL	37	3.90%
		Negative	3	0.32%
		Small (1+) / 15mg/dL	1	0.11%
		Trace/5 mg/dL	1	0.11%
UA-3		* Denotes Correct Response		
SG: Select Siemens Clinitek/Visual / Siemens - Plus Sys / Siemens Urine Dipstick	18			
		*3+	9	50.00%
		*4+	9	50.00%
		3+	9	50.00%
UA-3		* Denotes Correct Response		
SG: Select Roche Urinalysis Inst/Visual	134			
		*Large (3+) / 80 mg/dL	108	80.60%
		*Moderate (2+) / 40 mg/dL	25	18.66%
		Moderate (2+) / 40 mg/dL	25	18.66%
		Negative	1	0.75%
UA-3		* Denotes Correct Response		
SG: Select Consult Diag Urinalysis Inst/Visual	30			
		*Large (3+) / 80 mg/dL	27	90.00%
		*Moderate (2+) / 40 mg/dL	3	10.00%
		Moderate (2+) / 40 mg/dL	3	10.00%
UA-3		* Denotes Correct Response		
SG: Select McKesson Urinalysis Inst/Visual	47			
		*Large (3+) / 80 mg/dL	46	97.87%
		*Moderate (2+) / 40 mg/dL	1	2.13%
UA-3		* Denotes Correct Response		
SG: Select Henry Schein Urinalysis Inst/Visual	11			
		*Large (3+) / 80 mg/dL	9	81.82%



		*Moderate (2+) / 40 mg/dL	2	18.18%
UA-3		* Denotes Correct Response		
SG: Select Medi-Lab Performance UA Inst/Visual	8			
		*Large (3+) / 80 mg/dL	8	100.00%
UA-3		* Denotes Correct Response		
MG: Method Group	1233			
		*>=160 mg/dL	336	27.25%
		*Large (3+) / 80 mg/dL	817	66.26%
		*Moderate (2+) / 40 mg/dL	74	6.00%
		Moderate (2+) / 40 mg/dL	74	6.00%
		Negative	4	0.32%
		Small (1+) / 15mg/dL	1	0.08%
		Trace/5 mg/dL	1	0.08%
UA-3		* Denotes Correct Response		
MG: Method Group	1			
		*50 mg/dL	1	100.00%
UA-3		* Denotes Correct Response		
MG: Method Group	18			
		*3+	9	50.00%
		*4+	9	50.00%
		3+	9	50.00%
UA-3		* Denotes Correct Response		
MG: Method Group	2			
		*++	2	100.00%
		Urine pH		
UA-3		* Denotes Correct Response		
PG: Siemens Clinitek 10 / 100	12			
		*7	3	25.00%
		*7.5	7	58.33%
		*8 or greater	2	16.67%
UA-3		* Denotes Correct Response		
PG: Siemens Clinitek 50	80			
		*7	4	5.00%
		*7.5	76	95.00%
UA-3		* Denotes Correct Response		
PG: Siemens Clinitek 500	6			
		*7.5	6	100.00%
UA-3		* Denotes Correct Response		
PG: Siemens Clinitek Status	613			
		*7	15	2.45%
		*7.5	595	97.06%
		*8 or greater	3	0.49%
UA-3		* Denotes Correct Response		
PG: Siemens Clinitek Status / Siemens - Plus Sys / Siemens Urine Dipstick	17			
		*7.5	16	94.12%
		7	1	5.88%
UA-3		* Denotes Correct Response		
PG: Siemens Multistix - Visual	215			
		*7	9	4.19%
		*7.5	184	85.58%
		*8 or greater	16	7.44%
		6	3	1.40%
		6.5	3	1.40%
UA-3		* Denotes Correct Response		
PG: Henry Schein Urispec	8			
		*7	8	100.00%
UA-3		* Denotes Correct Response		
PG: Consult Diagnostics Urine Reagent Strips	8			
		*7	2	25.00%
		*7.5	6	75.00%
UA-3		* Denotes Correct Response		
PG: Diag Test Group Clarity Urochek	30			
		*7	10	33.33%
		*7.5	20	66.67%
UA-3		* Denotes Correct Response		
PG: Roche Chemstrips - Visual	117			
		*7	81	69.23%



		*7.5	15	12.82%
		*8 or greater	21	17.95%
UA-3		* Denotes Correct Response		
PG: Roche Urisys / Roche - Conventional / Roche Urine Dipstick	1			
		*7	1	100.00%
UA-3		* Denotes Correct Response		
PG: Roche Urisys / Roche - Arbitrary / Roche Urine Arbitrary Result	2			
		*7	2	100.00%
UA-3		* Denotes Correct Response		
SG: Select Siemens Clinitek/Visual	949			
		*7	31	3.27%
		*7.5	890	93.78%
		*8 or greater	22	2.32%
		6	3	0.32%
		6.5	3	0.32%
UA-3		* Denotes Correct Response		
SG: Select Siemens Clinitek/Visual / Siemens - Plus Sys / Siemens Urine Dipstick	18			
		*7.5	17	94.44%
		7	1	5.56%
UA-3		* Denotes Correct Response		
SG: Select Roche Urinalysis Inst/Visual	133			
		*7	90	67.67%
		*7.5	20	15.04%
		*8 or greater	23	17.29%
		8 or greater	23	17.29%
UA-3		* Denotes Correct Response		
SG: Select Consult Diag Urinalysis Inst/Visual	30			
		*7	5	16.67%
		*7.5	23	76.67%
		*8 or greater	1	3.33%
		6.5	1	3.33%
UA-3		* Denotes Correct Response		
SG: Select McKesson Urinalysis Inst/Visual	47			
		*7	13	27.66%
		*7.5	34	72.34%
UA-3		* Denotes Correct Response		
SG: Select Henry Schein Urinalysis Inst/Visual	11			
		*7	9	81.82%
		*7.5	2	18.18%
UA-3		* Denotes Correct Response		
SG: Select Medi-Lab Performance UA Inst/Visual	8			
		*7	3	37.50%
		*7.5	5	62.50%
UA-3		* Denotes Correct Response		
MG: Method Group	1234			
		*7	170	13.78%
		*7.5	1009	81.77%
		*8 or greater	48	3.89%
		6	3	0.24%
		6.5	4	0.32%
		8 or greater	48	3.89%
UA-3		* Denotes Correct Response		
MG: Method Group	1			
		*7	1	100.00%
UA-3		* Denotes Correct Response		
MG: Method Group	18			
		*7.5	17	94.44%
		7	1	5.56%
UA-3		* Denotes Correct Response		
MG: Method Group	2			
		*7	2	100.00%
		Urine Protein - Qual		
UA-3		* Denotes Correct Response		
PG: Siemens Clinitek 10 / 100	12			
		*100 (2+)	10	83.33%
		*300-500 (3+)	2	16.67%



UA-3		* Denotes Correct Response		
PG: Siemens Clinitek 50	80			
		*100 (2+)	13	16.25%
		*300-500 (3+)	66	82.50%
		>/= 1000 (4+)	1	1.25%
UA-3		* Denotes Correct Response		
PG: Siemens Clinitek 500	6			
		*100 (2+)	2	33.33%
		*300-500 (3+)	4	66.67%
UA-3		* Denotes Correct Response		
PG: Siemens Clinitek Status	614			
		*100 (2+)	103	16.78%
		*300-500 (3+)	507	82.57%
		>/= 1000 (4+)	1	0.16%
		30 (1+)	3	0.49%
UA-3		* Denotes Correct Response		
PG: Siemens Clinitek Status / Siemens - Plus Sys / Siemens Urine Dipstick	17			
		*2+	5	29.41%
		*3+	12	70.59%
UA-3		* Denotes Correct Response		
PG: Siemens Multistix - Visual	217			
		*100 (2+)	100	46.08%
		*300-500 (3+)	103	47.47%
		>/= 1000 (4+)	9	4.15%
		30 (1+)	3	1.38%
		Negative	2	0.92%
UA-3		* Denotes Correct Response		
PG: Henry Schein Urispec	8			
		*100 (2+)	6	75.00%
		*30 (1+)	1	12.50%
		*300-500 (3+)	1	12.50%
UA-3		* Denotes Correct Response		
PG: Consult Diagnostics Urine Reagent Strips	8			
		*100 (2+)	6	75.00%
		*300-500 (3+)	2	25.00%
UA-3		* Denotes Correct Response		
PG: Diag Test Group Clarity Urochek	30			
		*100 (2+)	24	80.00%
		*300-500 (3+)	5	16.67%
		30 (1+)	1	3.33%
UA-3		* Denotes Correct Response		
PG: Roche Chemstrips - Visual	123			
		*100 (2+)	85	69.11%
		*30 (1+)	16	13.01%
		*300-500 (3+)	20	16.26%
		>/= 1000 (4+)	1	0.81%
		Negative	1	0.81%
UA-3		* Denotes Correct Response		
PG: Roche Urisys	10			
		*100 (2+)	7	70.00%
		*30 (1+)	1	10.00%
		*300-500 (3+)	2	20.00%
UA-3		* Denotes Correct Response		
PG: Roche Urisys / Roche - Conventional / Roche Urine Dipstick	1			
		*100 mg/dL	1	100.00%
UA-3		* Denotes Correct Response		
PG: Roche Urisys / Roche - Arbitrary / Roche Urine Arbitrary Result	2			
		*++	2	100.00%
UA-3		* Denotes Correct Response		
PG: UriScan Reagent Strips	6			
		*100 (2+)	2	33.33%
		*30 (1+)	2	33.33%
		*300-500 (3+)	1	16.67%
		>/= 1000 (4+)	1	16.67%
UA-3		* Denotes Correct Response		
PG: Medi-Lab Performance Urine Strip - Visual	7			
		*100 (2+)	7	100.00%



UA-3		* Denotes Correct Response		
PG: Siemens Clinitek Advantus	21			
		*100 (2+)	21	100.00%
UA-3		* Denotes Correct Response		
PG: Consult Diag Urinalysis Test Strips	7			
		*100 (2+)	4	57.14%
		*300-500 (3+)	2	28.57%
		30 (1+)	1	14.29%
UA-3		* Denotes Correct Response		
PG: Consult Diag Urine Analyzer	15			
		*100 (2+)	11	73.33%
		*300-500 (3+)	3	20.00%
		30 (1+)	1	6.67%
UA-3		* Denotes Correct Response		
PG: Consult Diagnostics 120 Urine Analyzer	44			
		*100 (2+)	41	93.18%
		30 (1+)	3	6.82%
UA-3		* Denotes Correct Response		
SG: Select Siemens Clinitek/Visual	956			
		*100 (2+)	252	26.36%
		*300-500 (3+)	684	71.55%
		>/= 1000 (4+)	11	1.15%
		30 (1+)	6	0.63%
		Negative	2	0.21%
		Trace	1	0.10%
UA-3		* Denotes Correct Response		
SG: Select Siemens Clinitek/Visual / Siemens - Plus Sys / Siemens Urine Dipstick	18			
		*2+	6	33.33%
		*3+	12	66.67%
UA-3		* Denotes Correct Response		
SG: Select Roche Urinalysis Inst/Visual	139			
		*100 (2+)	94	67.63%
		*30 (1+)	18	12.95%
		*300-500 (3+)	25	17.99%
		>/= 1000 (4+)	1	0.72%
		Negative	1	0.72%
UA-3		* Denotes Correct Response		
SG: Select Consult Diag Urinalysis Inst/Visual	30			
		*100 (2+)	21	70.00%
		*300-500 (3+)	7	23.33%
		30 (1+)	2	6.67%
UA-3		* Denotes Correct Response		
SG: Select McKesson Urinalysis Inst/Visual	47			
		*100 (2+)	44	93.62%
		30 (1+)	3	6.38%
UA-3		* Denotes Correct Response		
SG: Select Henry Schein Urinalysis Inst/Visual	11			
		*100 (2+)	9	81.82%
		*30 (1+)	1	9.09%
		*300-500 (3+)	1	9.09%
UA-3		* Denotes Correct Response		
SG: Select Medi-Lab Performance UA Inst/Visual	8			
		*100 (2+)	8	100.00%
UA-3		* Denotes Correct Response		
MG: Method Group	1247			
		*100 (2+)	471	37.77%
		*30 (1+)	33	2.65%
		*300-500 (3+)	726	58.22%
		>/= 1000 (4+)	13	1.04%
		30 (1+)	33	2.65%
		Negative	3	0.24%
		Trace	1	0.08%
UA-3		* Denotes Correct Response		
MG: Method Group	1			
		*100 mg/dL	1	100.00%
UA-3		* Denotes Correct Response		
MG: Method Group	18			



		*2+	6	33.33%
		*3+	12	66.67%
UA-3		* Denotes Correct Response		
MG: Method Group	2			
		*++	2	100.00%
		Urine Blood/Hemoglobin		
UA-3		* Denotes Correct Response		
PG: Siemens Clinitek 10 / 100	12			
		*Large (3+)	11	91.67%
		*Moderate (2+)	1	8.33%
UA-3		* Denotes Correct Response		
PG: Siemens Clinitek 50	80			
		*Large (3+)	77	96.25%
		*Moderate (2+)	3	3.75%
UA-3		* Denotes Correct Response		
PG: Siemens Clinitek 500	6			
		*Large (3+)	6	100.00%
UA-3		* Denotes Correct Response		
PG: Siemens Clinitek Status	613			
		*Large (3+)	550	89.72%
		*Moderate (2+)	59	9.62%
		Negative	2	0.33%
		Small (1+)	2	0.33%
UA-3		* Denotes Correct Response		
PG: Siemens Clinitek Status / Siemens - Plus Sys / Siemens Urine Dipstick	17			
		*2+	4	23.53%
		*3+	13	76.47%
UA-3		* Denotes Correct Response		
PG: Siemens Multistix - Visual	215			
		*Large (3+)	207	96.28%
		*Moderate (2+)	4	1.86%
		Negative	3	1.40%
		Small (1+)	1	0.47%
UA-3		* Denotes Correct Response		
PG: Henry Schein Urispec	8			
		*Large (3+)	7	87.50%
		*Moderate (2+)	1	12.50%
UA-3		* Denotes Correct Response		
PG: Diag Test Group Clarity Urochek	30			
		*Large (3+)	29	96.67%
		*Moderate (2+)	1	3.33%
UA-3		* Denotes Correct Response		
PG: Roche Chemstrips - Visual	122			
		*Large (3+)	106	86.89%
		*Moderate (2+)	14	11.48%
		Negative	2	1.64%
UA-3		* Denotes Correct Response		
PG: Roche Urisys	10			
		*Large (3+)	9	90.00%
		*Moderate (2+)	1	10.00%
UA-3		* Denotes Correct Response		
PG: Roche Urisys / Roche - Conventional / Roche Urine Dipstick	1			
		*250 Ery/uL	1	100.00%
UA-3		* Denotes Correct Response		
PG: Roche Urisys / Roche - Arbitrary / Roche Urine Arbitrary Result	2			
		*++	2	100.00%
UA-3		* Denotes Correct Response		
PG: UriScan Reagent Strips	6			
		*Large (3+)	5	83.33%
		Small (1+)	1	16.67%
UA-3		* Denotes Correct Response		
PG: Medi-Lab Performance Urine Strip - Visual	7			
		*Large (3+)	7	100.00%
UA-3		* Denotes Correct Response		
PG: Siemens Clinitek Advantus	21			
		*Large (3+)	21	100.00%



UA-3		* Denotes Correct Response		
PG: Consult Diag Urinalysis Test Strips	7			
		*Large (3+)	5	71.43%
		*Moderate (2+)	1	14.29%
		Small (1+)	1	14.29%
UA-3		* Denotes Correct Response		
PG: Consult Diag Urine Analyzer	15			
		*Large (3+)	12	80.00%
		*Moderate (2+)	2	13.33%
		Negative	1	6.67%
UA-3		* Denotes Correct Response		
PG: Consult Diagnostics 120 Urine Analyzer	44			
		*Large (3+)	44	100.00%
UA-3		* Denotes Correct Response		
SG: Select Siemens Clinitek/Visual	949			
		*Large (3+)	874	92.10%
		*Moderate (2+)	67	7.06%
		Negative	5	0.53%
		Small (1+)	3	0.32%
UA-3		* Denotes Correct Response		
SG: Select Siemens Clinitek/Visual / Siemens - Plus Sys / Siemens Urine Dipstick	18			
		*2+	4	22.22%
		*3+	14	77.78%
UA-3		* Denotes Correct Response		
SG: Select Roche Urinalysis Inst/Visual	138			
		*Large (3+)	120	86.96%
		*Moderate (2+)	16	11.59%
		Negative	2	1.45%
UA-3		* Denotes Correct Response		
SG: Select Consult Diag Urinalysis Inst/Visual	30			
		*Large (3+)	25	83.33%
		*Moderate (2+)	3	10.00%
		Negative	1	3.33%
		Small (1+)	1	3.33%
UA-3		* Denotes Correct Response		
SG: Select McKesson Urinalysis Inst/Visual	47			
		*Large (3+)	47	100.00%
UA-3		* Denotes Correct Response		
SG: Select Henry Schein Urinalysis Inst/Visual	11			
		*Large (3+)	10	90.91%
		*Moderate (2+)	1	9.09%
UA-3		* Denotes Correct Response		
SG: Select Medi-Lab Performance UA Inst/Visual	8			
		*Large (3+)	8	100.00%
UA-3		* Denotes Correct Response		
MG: Method Group	1238			
		*Large (3+)	1137	91.84%
		*Moderate (2+)	88	7.11%
		Negative	8	0.65%
		Small (1+)	5	0.40%
UA-3		* Denotes Correct Response		
MG: Method Group	1			
		*250 Ery/uL	1	100.00%
UA-3		* Denotes Correct Response		
MG: Method Group	18			
		*2+	4	22.22%
		*3+	14	77.78%
UA-3		* Denotes Correct Response		
MG: Method Group	2			
		*++	2	100.00%
		Urine Glucose (Reducing Sub)		
UA-3		* Denotes Correct Response		
PG: Siemens Clinitek 10 / 100	12			
		*Negative	12	100.00%
UA-3		* Denotes Correct Response		
PG: Siemens Clinitek 50	80			
		*Negative	80	100.00%



UA-3		* Denotes Correct Response		
PG: Siemens Clinitek 500	6			
		*Negative	6	100.00%
UA-3		* Denotes Correct Response		
PG: Siemens Clinitek Status	613			
		*Negative	613	100.00%
UA-3		* Denotes Correct Response		
PG: Siemens Clinitek Status / Siemens - Plus Sys / Siemens Urine Dipstick	17			
		*Negative	17	100.00%
UA-3		* Denotes Correct Response		
PG: Siemens Multistix - Visual	216			
		*Negative	214	99.07%
		<100 mg/dL	1	0.46%
		250 mg/dL	1	0.46%
UA-3		* Denotes Correct Response		
PG: Henry Schein Urispec	8			
		*Negative	8	100.00%
UA-3		* Denotes Correct Response		
PG: Consult Diagnostics Urine Reagent Strips	8			
		*Negative	8	100.00%
UA-3		* Denotes Correct Response		
PG: Diag Test Group Clarity Urochek	30			
		*Negative	30	100.00%
UA-3		* Denotes Correct Response		
PG: Roche Chemstrips - Visual	123			
		*Negative	121	98.37%
		250 mg/dL	2	1.63%
UA-3		* Denotes Correct Response		
PG: Roche Urisys	10			
		*Negative	10	100.00%
UA-3		* Denotes Correct Response		
PG: Roche Urisys / Roche - Conventional / Roche Urine Dipstick	1			
		*Normal	1	100.00%
UA-3		* Denotes Correct Response		
PG: Roche Urisys / Roche - Arbitrary / Roche Urine Arbitrary Result	2			
		*Normal	2	100.00%
UA-3		* Denotes Correct Response		
PG: UriScan Reagent Strips	6			
		*Negative	6	100.00%
UA-3		* Denotes Correct Response		
PG: Medi-Lab Performance Urine Strip - Visual	7			
		*Negative	7	100.00%
UA-3		* Denotes Correct Response		
PG: Siemens Clinitek Advantus	21			
		*Negative	21	100.00%
UA-3		* Denotes Correct Response		
PG: Siemens Clinitek Advantus / Siemens - Plus Sys / Siemens Urine Dipstick	1			
		*Negative	1	100.00%
UA-3		* Denotes Correct Response		
PG: Consult Diag Urinalysis Test Strips	7			
		*Negative	7	100.00%
UA-3		* Denotes Correct Response		
PG: Consult Diag Urine Analyzer	15			
		*Negative	15	100.00%
UA-3		* Denotes Correct Response		
PG: Consult Diagnostics 120 Urine Analyzer	44			
		*Negative	44	100.00%
UA-3		* Denotes Correct Response		
SG: Select Siemens Clinitek/Visual	954			
		*Negative	952	99.79%
		<100 mg/dL	1	0.10%
		250 mg/dL	1	0.10%
UA-3		* Denotes Correct Response		
SG: Select Siemens Clinitek/Visual / Siemens - Plus Sys / Siemens Urine Dipstick	18			
		*Negative	18	100.00%



UA-3		* Denotes Correct Response		
SG: Select Roche Urinalysis Inst/Visual	139			
		*Negative	137	98.56%
		250 mg/dL	2	1.44%
UA-3		* Denotes Correct Response		
SG: Select Consult Diag Urinalysis Inst/Visual	30			
		*Negative	30	100.00%
UA-3		* Denotes Correct Response		
SG: Select McKesson Urinalysis Inst/Visual	47			
		*Negative	47	100.00%
UA-3		* Denotes Correct Response		
SG: Select Henry Schein Urinalysis Inst/Visual	11			
		*Negative	11	100.00%
UA-3		* Denotes Correct Response		
SG: Select Medi-Lab Performance UA Inst/Visual	8			
		*Negative	8	100.00%
UA-3		* Denotes Correct Response		
MG: Method Group	1245			
		*Negative	1241	99.68%
		<100 mg/dL	1	0.08%
		250 mg/dL	3	0.24%
UA-3		* Denotes Correct Response		
MG: Method Group	1			
		*Normal	1	100.00%
UA-3		* Denotes Correct Response		
MG: Method Group	18			
		*Negative	18	100.00%
UA-3		* Denotes Correct Response		
MG: Method Group	2			
		*Normal	2	100.00%
		Urine Leuk Esterase		
UA-3		* Denotes Correct Response		
PG: Siemens Clinitek 10 / 100	12			
		*Negative	12	100.00%
UA-3		* Denotes Correct Response		
PG: Siemens Clinitek 50	80			
		*Negative	80	100.00%
UA-3		* Denotes Correct Response		
PG: Siemens Clinitek 500	6			
		*Negative	6	100.00%
UA-3		* Denotes Correct Response		
PG: Siemens Clinitek Status	611			
		*Negative	609	99.67%
		Large (3+)	1	0.16%
		Moderate (2+)	1	0.16%
UA-3		* Denotes Correct Response		
PG: Siemens Clinitek Status / Siemens - Plus Sys / Siemens Urine Dipstick	17			
		*Negative	17	100.00%
UA-3		* Denotes Correct Response		
PG: Siemens Multistix - Visual	216			
		*Negative	214	99.07%
		Trace	2	0.93%
UA-3		* Denotes Correct Response		
PG: Henry Schein Urispec	8			
		*Negative	8	100.00%
UA-3		* Denotes Correct Response		
PG: Consult Diagnostics Urine Reagent Strips	8			
		*Negative	8	100.00%
UA-3		* Denotes Correct Response		
PG: Diag Test Group Clarity Urochek	30			
		*Negative	30	100.00%
UA-3		* Denotes Correct Response		
PG: Roche Chemstrips - Visual	123			
		*Negative	121	98.37%
		Large (3+)	1	0.81%
		Trace	1	0.81%
UA-3		* Denotes Correct Response		



PG: Roche Urisys	10			
		*Negative	10	100.00%
UA-3		* Denotes Correct Response		
PG: Roche Urisys / Roche - Conventional / Roche Urine Dipstick	1			
		*Negative	1	100.00%
UA-3		* Denotes Correct Response		
PG: Roche Urisys / Roche - Arbitrary / Roche Urine Arbitrary Result	2			
		*Negative	2	100.00%
UA-3		* Denotes Correct Response		
PG: UriScan Reagent Strips	6			
		*Negative	5	83.33%
		Small (1+)	1	16.67%
UA-3		* Denotes Correct Response		
PG: Medi-Lab Performance Urine Strip - Visual	7			
		*Negative	7	100.00%
UA-3		* Denotes Correct Response		
PG: Siemens Clinitek Advantus	21			
		*Negative	21	100.00%
UA-3		* Denotes Correct Response		
PG: Siemens Clinitek Advantus / Siemens - Plus Sys / Siemens Urine Dipstick	1			
		*Negative	1	100.00%
UA-3		* Denotes Correct Response		
PG: Consult Diag Urinalysis Test Strips	7			
		*Negative	6	85.71%
		Moderate (2+)	1	14.29%
UA-3		* Denotes Correct Response		
PG: Consult Diag Urine Analyzer	15			
		*Negative	15	100.00%
UA-3		* Denotes Correct Response		
PG: Consult Diagnostics 120 Urine Analyzer	44			
		*Negative	43	97.73%
		Trace	1	2.27%
UA-3		* Denotes Correct Response		
SG: Select Siemens Clinitek/Visual	948			
		*Negative	944	99.58%
		Large (3+)	1	0.11%
		Moderate (2+)	1	0.11%
		Trace	2	0.21%
UA-3		* Denotes Correct Response		
SG: Select Siemens Clinitek/Visual / Siemens - Plus Sys / Siemens Urine Dipstick	18			
		*Negative	18	100.00%
UA-3		* Denotes Correct Response		
SG: Select Roche Urinalysis Inst/Visual	139			
		*Negative	137	98.56%
		Large (3+)	1	0.72%
		Trace	1	0.72%
UA-3		* Denotes Correct Response		
SG: Select Consult Diag Urinalysis Inst/Visual	30			
		*Negative	29	96.67%
		Moderate (2+)	1	3.33%
UA-3		* Denotes Correct Response		
SG: Select McKesson Urinalysis Inst/Visual	47			
		*Negative	46	97.87%
		Trace	1	2.13%
UA-3		* Denotes Correct Response		
SG: Select Henry Schein Urinalysis Inst/Visual	11			
		*Negative	11	100.00%
UA-3		* Denotes Correct Response		
MG: Method Group	1238			
		*Negative	1228	99.19%
		Large (3+)	2	0.16%
		Moderate (2+)	2	0.16%
		Small (1+)	1	0.08%
		Trace	5	0.40%
UA-3		* Denotes Correct Response		
MG: Method Group	1			



		*Negative	1	100.00%
UA-3		* Denotes Correct Response		
MG: Method Group	18			
		*Negative	18	100.00%
UA-3		* Denotes Correct Response		
MG: Method Group	2			
		*Negative	2	100.00%
		Urine Nitrite		
UA-3		* Denotes Correct Response		
PG: Siemens Clinitek 10 / 100	12			
		*Negative	12	100.00%
UA-3		* Denotes Correct Response		
PG: Siemens Clinitek 50	80			
		*Negative	80	100.00%
UA-3		* Denotes Correct Response		
PG: Siemens Clinitek 500	6			
		*Negative	6	100.00%
UA-3		* Denotes Correct Response		
PG: Siemens Clinitek Status	613			
		*Negative	612	99.84%
		Positive	1	0.16%
UA-3		* Denotes Correct Response		
PG: Siemens Clinitek Status / Siemens - Plus Sys / Siemens Urine Dipstick	17			
		*Negative	17	100.00%
UA-3		* Denotes Correct Response		
PG: Siemens Multistix - Visual	215			
		*Negative	214	99.53%
		Positive	1	0.47%
UA-3		* Denotes Correct Response		
PG: Henry Schein Urispec	8			
		*Negative	8	100.00%
UA-3		* Denotes Correct Response		
PG: Consult Diagnostics Urine Reagent Strips	8			
		*Negative	8	100.00%
UA-3		* Denotes Correct Response		
PG: Diag Test Group Clarity Urochek	29			
		*Negative	29	100.00%
UA-3		* Denotes Correct Response		
PG: Roche Chemstrips - Visual	123			
		*Negative	123	100.00%
UA-3		* Denotes Correct Response		
PG: Roche Urisys	10			
		*Negative	10	100.00%
UA-3		* Denotes Correct Response		
PG: Roche Urisys / Roche - Conventional / Roche Urine Dipstick	1			
		*Negative	1	100.00%
UA-3		* Denotes Correct Response		
PG: Roche Urisys / Roche - Arbitrary / Roche Urine Arbitrary Result	2			
		*Negative	2	100.00%
UA-3		* Denotes Correct Response		
PG: UriScan Reagent Strips	6			
		*Negative	6	100.00%
UA-3		* Denotes Correct Response		
PG: Medi-Lab Performance Urine Strip - Visual	7			
		*Negative	7	100.00%
UA-3		* Denotes Correct Response		
PG: Siemens Clinitek Advantus	21			
		*Negative	21	100.00%
UA-3		* Denotes Correct Response		
PG: Consult Diag Urinalysis Test Strips	7			
		*Negative	7	100.00%
UA-3		* Denotes Correct Response		
PG: Consult Diag Urine Analyzer	15			
		*Negative	15	100.00%
UA-3		* Denotes Correct Response		
PG: Consult Diagnostics 120 Urine Analyzer	44			



		*Negative	44	100.00%
UA-3		* Denotes Correct Response		
SG: Select Siemens Clinitek/Visual	949			
		*Negative	947	99.79%
		Positive	2	0.21%
UA-3		* Denotes Correct Response		
SG: Select Siemens Clinitek/Visual / Siemens - Plus Sys / Siemens Urine Dipstick	18			
		*Negative	18	100.00%
UA-3		* Denotes Correct Response		
SG: Select Roche Urinalysis Inst/Visual	139			
		*Negative	139	100.00%
UA-3		* Denotes Correct Response		
SG: Select Consult Diag Urinalysis Inst/Visual	30			
		*Negative	30	100.00%
UA-3		* Denotes Correct Response		
SG: Select McKesson Urinalysis Inst/Visual	47			
		*Negative	47	100.00%
UA-3		* Denotes Correct Response		
SG: Select Henry Schein Urinalysis Inst/Visual	11			
		*Negative	11	100.00%
UA-3		* Denotes Correct Response		
SG: Select Medi-Lab Performance UA Inst/Visual	8			
		*Negative	8	100.00%
UA-3		* Denotes Correct Response		
MG: Method Group	1238			
		*Negative	1236	99.84%
		Positive	2	0.16%
UA-3		* Denotes Correct Response		
MG: Method Group	1			
		*Negative	1	100.00%
UA-3		* Denotes Correct Response		
MG: Siemens - Plus Sys	18			
		*Negative	18	100.00%
UA-3		* Denotes Correct Response		
MG: Siemens - Plus Sys	2			
		*Negative	2	100.00%
		Urine Urobilogen		
UA-3		* Denotes Correct Response		
PG: Siemens Clinitek 10 / 100	12			
		*1 mg/dL	4	33.33%
		*2.0 mg/dL	4	33.33%
		*4 mg/dL	3	25.00%
		0.2/Normal	1	8.33%
UA-3		* Denotes Correct Response		
PG: Siemens Clinitek 50	80			
		*>= 8 mg/dL	2	2.50%
		*2.0 mg/dL	20	25.00%
		*4 mg/dL	51	63.75%
		0.2/Normal	5	6.25%
		1 mg/dL	2	2.50%
UA-3		* Denotes Correct Response		
PG: Siemens Clinitek 500	5			
		*2.0 mg/dL	2	40.00%
		*4 mg/dL	3	60.00%
UA-3		* Denotes Correct Response		
PG: Siemens Clinitek Status	605			
		*>= 8 mg/dL	20	3.31%
		*2.0 mg/dL	23	3.80%
		*4 mg/dL	557	92.07%
		0.2/Normal	5	0.83%
UA-3		* Denotes Correct Response		
PG: Siemens Clinitek Status / Siemens - Plus Sys / Siemens Urine Dipstick	7			
		*33 umol/L	1	14.29%
		*66 umol/L	5	71.43%
		16 umol/L	1	14.29%
UA-3		* Denotes Correct Response		



PG: Siemens Clinitek Status / Siemens - Plus Sys / Siemens Urine Urobilinogen	9			
		*4.0 EU/dL	9	100.00%
UA-3		* Denotes Correct Response		
PG: Siemens Multistix - Visual	206			
		*>= 8 mg/dL	28	13.59%
		*2.0 mg/dL	75	36.41%
		*4 mg/dL	78	37.86%
		0.2/Normal	15	7.28%
		1 mg/dL	10	4.85%
UA-3		* Denotes Correct Response		
PG: Henry Schein Urispec	7			
		*0.2/Normal	2	28.57%
		*2.0 mg/dL	5	71.43%
UA-3		* Denotes Correct Response		
PG: Diag Test Group Clarity Urochek	30			
		*0.2/Normal	7	23.33%
		*1 mg/dL	11	36.67%
		*2.0 mg/dL	7	23.33%
		*4 mg/dL	4	13.33%
		>= 8 mg/dL	1	3.33%
UA-3		* Denotes Correct Response		
PG: Roche Chemstrips - Visual	117			
		*1 mg/dL	52	44.44%
		*2.0 mg/dL	4	3.42%
		*4 mg/dL	48	41.03%
		>= 8 mg/dL	3	2.56%
		0.2/Normal	10	8.55%
UA-3		* Denotes Correct Response		
PG: Roche Urisys	10			
		*1 mg/dL	3	30.00%
		*2.0 mg/dL	2	20.00%
		*4 mg/dL	5	50.00%
UA-3		* Denotes Correct Response		
PG: Roche Urisys / Roche - Conventional / Roche Urine Dipstick	1			
		*1 mg/dL	1	100.00%
UA-3		* Denotes Correct Response		
PG: Roche Urisys / Roche - Arbitrary / Roche Urine Arbitrary Result	2			
		*+	2	100.00%
UA-3		* Denotes Correct Response		
PG: UriScan Reagent Strips	5			
		*0.2/Normal	2	40.00%
		*2.0 mg/dL	1	20.00%
		*4 mg/dL	1	20.00%
		>= 8 mg/dL	1	20.00%
UA-3		* Denotes Correct Response		
PG: Medi-Lab Performance Urine Strip - Visual	7			
		*0.2/Normal	1	14.29%
		*2.0 mg/dL	5	71.43%
		*4 mg/dL	1	14.29%
UA-3		* Denotes Correct Response		
PG: Siemens Clinitek Advantus	21			
		*1 mg/dL	2	9.52%
		*2.0 mg/dL	12	57.14%
		*4 mg/dL	7	33.33%
UA-3		* Denotes Correct Response		
PG: Consult Diag Urinalysis Test Strips	7			
		*0.2/Normal	2	28.57%
		*1 mg/dL	1	14.29%
		*2.0 mg/dL	1	14.29%
		*4 mg/dL	3	42.86%
UA-3		* Denotes Correct Response		
PG: Consult Diag Urine Analyzer	15			
		*0.2/Normal	9	60.00%
		*1 mg/dL	4	26.67%
		*4 mg/dL	2	13.33%
UA-3		* Denotes Correct Response		
PG: Consult Diagnostics 120 Urine Analyzer	43			



		*0.2/Normal	33	76.74%
		*1 mg/dL	10	23.26%
UA-3		* Denotes Correct Response		
SG: Select Siemens Clinitek/Visual	930			
		*>= 8 mg/dL	50	5.38%
		*1 mg/dL	18	1.94%
		*2.0 mg/dL	136	14.62%
		*4 mg/dL	700	75.27%
		0.2/Normal	26	2.80%
		1 mg/dL	18	1.94%
UA-3		* Denotes Correct Response		
SG: Select Siemens Clinitek/Visual / Siemens - Plus Sys / Siemens Urine Dipstick	8			
		*33 umol/L	2	25.00%
		*66 umol/L	5	62.50%
		16 umol/L	1	12.50%
UA-3		* Denotes Correct Response		
SG: Select Roche Urinalysis Inst/Visual	133			
		*0.2/Normal	12	9.02%
		*1 mg/dL	58	43.61%
		*2.0 mg/dL	6	4.51%
		*4 mg/dL	54	40.60%
		>/= 8 mg/dL	3	2.26%
		0.2/Normal	12	9.02%
UA-3		* Denotes Correct Response		
SG: Select Consult Diag Urinalysis Inst/Visual	30			
		*0.2/Normal	13	43.33%
		*1 mg/dL	7	23.33%
		*2.0 mg/dL	2	6.67%
		*4 mg/dL	8	26.67%
		0.2/Normal	13	43.33%
UA-3		* Denotes Correct Response		
SG: Select McKesson Urinalysis Inst/Visual	46			
		*0.2/Normal	34	73.91%
		*1 mg/dL	10	21.74%
		*2.0 mg/dL	2	4.35%
UA-3		* Denotes Correct Response		
SG: Select Henry Schein Urinalysis Inst/Visual	10			
		*0.2/Normal	2	20.00%
		*2.0 mg/dL	7	70.00%
		*4 mg/dL	1	10.00%
UA-3		* Denotes Correct Response		
SG: Select Medi-Lab Performance UA Inst/Visual	8			
		*0.2/Normal	2	25.00%
		*2.0 mg/dL	5	62.50%
		*4 mg/dL	1	12.50%
UA-3		* Denotes Correct Response		
MG: Method Group	1211			
		*>= 8 mg/dL	57	4.71%
		*0.2/Normal	105	8.67%
		*1 mg/dL	109	9.00%
		*2.0 mg/dL	169	13.96%
		*4 mg/dL	771	63.67%
		>/= 8 mg/dL	57	4.71%
		0.2/Normal	105	8.67%
		1 mg/dL	109	9.00%
UA-3		* Denotes Correct Response		
MG: Method Group	1			
		*1 mg/dL	1	100.00%
UA-3		* Denotes Correct Response		
MG: Method Group	17			
		*33 umol/L	2	11.76%
		*4.0 EU/dL	9	52.94%
		*66 umol/L	5	29.41%
		16 umol/L	1	5.88%
UA-3		* Denotes Correct Response		
MG: Method Group	2			
		*+	2	100.00%



		Urine Microalbumin		
UA-3		* Denotes Correct Response		
PG: Siemens Clinitek 50	17			
		*150 mg/L	16	94.12%
		Negative	1	5.88%
UA-3		* Denotes Correct Response		
PG: Siemens Clinitek Status	37			
		*100 mg/L (3+)	1	2.70%
		*150 mg/L	35	94.59%
		80 mg/L	1	2.70%
UA-3		* Denotes Correct Response		
PG: Roche Chemstrip Micral	28			
		*100 mg/L (3+)	23	82.14%
		*150 mg/L	3	10.71%
		50 mg/L (2+)	2	7.14%
UA-3		* Denotes Correct Response		
SG: Select Siemens Clinitek/Visual	60			
		*100 mg/L (3+)	4	6.67%
		*150 mg/L	53	88.33%
		20/30 mg/L (Small	1	1.67%
		80 mg/L	1	1.67%
		Negative	1	1.67%
UA-3		* Denotes Correct Response		
AG: All Methods	90			
		*100 mg/L (3+)	27	30.00%
		*150 mg/L	57	63.33%
		20/30 mg/L (Small	1	1.11%
		50 mg/L (2+)	2	2.22%
		80 mg/L	1	1.11%
		Negative	2	2.22%
Siemens Clinitek Status	37	*100 mg/L (3+)	1	2.70%
Siemens Multistix - Visual	4	*100 mg/L (3+)	3	75.00%
Roche Chemstrip Micral	28	*100 mg/L (3+)	23	82.14%
Siemens Clinitek 10 / 100	2	*150 mg/L	2	100.00%
Siemens Clinitek 50	17	*150 mg/L	16	94.12%
Siemens Clinitek Status	37	*150 mg/L	35	94.59%
Roche Chemstrip Micral	28	*150 mg/L	3	10.71%
Hemocue Albumin 201	1	*150 mg/L	1	100.00%
Siemens Multistix - Visual	4	20/30 mg/L (Small	1	25.00%
Roche Chemstrip Micral	28	50 mg/L (2+)	2	7.14%
Siemens Clinitek Status	37	80 mg/L	1	2.70%
Siemens Clinitek 50	17	Negative	1	5.88%
Roche Chemstrips - Visual	1	Negative	1	100.00%
		Fluid pH		
UA-3		* Denotes Correct Response		
PG: Nitrazine Paper	13			
		*7	2	15.38%
		*7.5	10	76.92%
		*8 or greater	1	7.69%
		Urine hCG		
UA-3		* Denotes Correct Response		
AG: All Methods	1070			
		*Negative	1062	99.25%
		Positive	8	0.75%
ASI ProPhase Plus	1	*Negative	1	100.00%
Siemens Clinitek Status	12	*Negative	12	100.00%
Beck Coulter ICON 20 HCG	6	*Negative	6	100.00%
Beck Coulter ICON 25 hCG	66	*Negative	66	100.00%
Alere Acceava HCG Basic II	24	*Negative	24	100.00%
Alere Acceava hCG Combo	5	*Negative	5	100.00%
Alere Acceava hCG-Urine	21	*Negative	21	100.00%
Conception Tech One Step hCG	2	*Negative	2	100.00%
Diag Test Group Clarity hCG Plus Urine Cassettes	5	*Negative	5	100.00%
Fisher HC Sure-Vue hCG STAT	4	*Negative	4	100.00%
Fisher HC Sure-Vue Urine hCG	12	*Negative	12	100.00%
Fisher HC Sure-Vure U/S hCG	6	*Negative	5	83.33%
Sekisui Diagnostics OSOM Urine	11	*Negative	11	100.00%
Sekisui Diagnostics OSOM Card Pregnancy	29	*Negative	29	100.00%



Sekisui Diagnostics OSOM hCG Combo	42	*Negative	42	100.00%
Henry Schein One Step + (hCG)	55	*Negative	55	100.00%
Immunostics Detector Combo	3	*Negative	3	100.00%
Immunostics hCG Detector-urine	3	*Negative	3	100.00%
Instant Tech. i Pregnancy	3	*Negative	2	66.67%
Inverness Clearview Combo II	2	*Negative	2	100.00%
Inv Clearview Urine - Cassette	3	*Negative	3	100.00%
Inverness Signify hCG combo	3	*Negative	3	100.00%
Inverness Signify hCG urine	4	*Negative	4	100.00%
LifeSign Status hCG	1	*Negative	1	100.00%
McKesson Medi-Lab hCG - Urine	33	*Negative	33	100.00%
Meridian Immunocard Stat!	1	*Negative	1	100.00%
Permaxim RediScreen	3	*Negative	3	100.00%
Polymedco Poly stat hCG Urine	1	*Negative	1	100.00%
Quidel QuickVue + hCG	25	*Negative	25	100.00%
Quidel QuickVue One-Step hCG Combo U/S	66	*Negative	66	100.00%
Quidel QuickVue One-Step Urine	143	*Negative	142	99.30%
Quidel RapidVue hCG - Urine	8	*Negative	8	100.00%
SAS Ultra HCG	2	*Negative	2	100.00%
Select Diag hCG Cassette (urine)	20	*Negative	20	100.00%
Select Diag hCG Combo	10	*Negative	9	90.00%
Select Diag hCG Dipstick	5	*Negative	5	100.00%
SP (Cardinal) Brand combo	16	*Negative	16	100.00%
SP (Cardinal) Brand-cassette	28	*Negative	28	100.00%
SP Brand (Cardinal) HCG Combo	46	*Negative	46	100.00%
Stanbio Qupid	4	*Negative	4	100.00%
Stanbio QuStick	1	*Negative	1	100.00%
Stanbio TRUE hCG	5	*Negative	5	100.00%
MooreBrand hCG Cassette	3	*Negative	3	100.00%
Nova Scientific Century S/U hCG	2	*Negative	2	100.00%
Diag Test Group Clarity hCG Urine Strips	1	*Negative	1	100.00%
Jant AccuTest	5	*Negative	5	100.00%
Instant-View Urine Pregnancy	1	*Negative	1	100.00%
Becton Dickinson Precise One-Step	1	*Negative	1	100.00%
Consult Diag hCG Urine Cassette	170	*Negative	167	98.24%
CLIAwaived Inc Pregnancy Urine Dip-Strip	4	*Negative	4	100.00%
Germaine AIM Step Combo Pregnancy	1	*Negative	1	100.00%
NDC ProAdvantage Urine hCG Cassette	5	*Negative	5	100.00%
Consult Diag hCG Urine Dipstick	35	*Negative	35	100.00%
SP (Cardinal) Brand Dipstick	10	*Negative	10	100.00%
Consult Diag hCG Combo Cassette (Urine)	12	*Negative	12	100.00%
McKesson Medi-Lab Performance hCG Urine Cassette	8	*Negative	8	100.00%
Alere hCG Combo Cassette (25 mIU/mL)	2	*Negative	2	100.00%
Alere hCG Cassette (25 mIU/mL)	10	*Negative	10	100.00%
Alere hCG Dipstick	7	*Negative	7	100.00%
Medline Pregnancy Urine Dip-Strip	2	*Negative	2	100.00%
Medline Pregnancy Urine Cassette Test	4	*Negative	4	100.00%
Quidel QuickVue hCG Urine Test	17	*Negative	17	100.00%
MooreBrand hCG Dipstick	2	*Negative	2	100.00%
Siemens Clinitest hCG	1	*Negative	1	100.00%
Henry Schein One Step + HCG Combo	2	*Negative	2	100.00%
Henry Schein One Step hCG Urine Cassette	22	*Negative	21	95.45%
Medline hCG Test Strip 25 mIU/mL	1	*Negative	1	100.00%
NCS hCG One Step Pregnancy Test	1	*Negative	1	100.00%
MediChoice Urine hCG Cassette	1	*Negative	1	100.00%
Fisher HC Sure-Vure U/S hCG	6	Positive	1	16.67%
Instant Tech. i Pregnancy	3	Positive	1	33.33%
Quidel QuickVue One-Step Urine	143	Positive	1	0.70%
Select Diag hCG Combo	10	Positive	1	10.00%
Consult Diag hCG Urine Cassette	170	Positive	3	1.76%
Henry Schein One Step hCG Urine Cassette	22	Positive	1	4.55%
661 - Urine Colony Count;784 - Presump ID/Colony Count - Urine;785 - Presump ID/Colony Count-Urine & Strep Combo;775 - Urine Culture;778 - Urine & Throat Cultures;779 - Urine & Genital Cultures;781 - Urine, Throat & Genital Cultures;782 - Complete Bacteri	Labs	Response	Number	Percent
		Susceptibility Cefazolin		
UC-11		* Denotes Correct Response		
AG: All Methods	2			
		*Resistant	2	100.00%



Susceptibility	2	*Resistant	2	100.00%
		Susceptibility Cefepime		
UC-11		* Denotes Correct Response		
AG: All Methods	2			
		*Sensitive	2	100.00%
Susceptibility	2	*Sensitive	2	100.00%
		Urine Culture		
UC-11		* Denotes Correct Response		
AG: All Methods	30			
		*Gram Negative organism	4	13.33%
		*Growth, referred for identification	14	46.67%
		*Presumptive Proteus species	2	6.67%
		*Proteus species NOS	3	10.00%
		*Uricult - Growth (Positive)	7	23.33%
Cultures	30	*Gram Negative organism	4	13.33%
Cultures		*Growth, referred for identification	14	46.67%
Cultures		*Presumptive Proteus species	2	6.67%
Cultures		*Proteus species NOS	3	10.00%
Cultures		*Uricult - Growth (Positive)	7	23.33%
UC-12		* Denotes Correct Response		
AG: All Methods	30			
		*Gram Negative organism	4	13.33%
		*Growth, referred for identification	15	50.00%
		*Klebsiella pneumoniae	3	10.00%
		*Presumptive Klebsiella species	1	3.33%
		*Uricult - Growth (Positive)	7	23.33%
Cultures	30	*Gram Negative organism	4	13.33%
Cultures		*Growth, referred for identification	15	50.00%
Cultures		*Klebsiella pneumoniae	3	10.00%
Cultures		*Presumptive Klebsiella species	1	3.33%
Cultures		*Uricult - Growth (Positive)	7	23.33%
UC-13		* Denotes Correct Response		
AG: All Methods	29			
		*Enterococcus faecium	2	6.90%
		*Gram Positive organism	2	6.90%
		*Growth, referred for identification	17	58.62%
		*Presumptive Enterococcus species	1	3.45%
		*Uricult - Growth (Positive)	7	24.14%
Cultures	29	*Enterococcus faecium	2	6.90%
Cultures		*Gram Positive organism	2	6.90%
Cultures		*Growth, referred for identification	17	58.62%
Cultures		*Presumptive Enterococcus species	1	3.45%
Cultures		*Uricult - Growth (Positive)	7	24.14%
UC-14		* Denotes Correct Response		
AG: All Methods	11			
		*Growth, referred for identification	3	27.27%
		*Uricult - Growth (Positive)	4	36.36%
		*Yeast, NOS	1	9.09%
		Gram Positive organism	2	18.18%
		Presumptive Staphylococcus species	1	9.09%
Cultures	11	*Growth, referred for identification	3	27.27%
Cultures		*Uricult - Growth (Positive)	4	36.36%
Cultures		*Yeast, NOS	1	9.09%
Cultures		Gram Positive organism	2	18.18%
Cultures		Presumptive Staphylococcus species	1	9.09%
UC-15		* Denotes Correct Response		
AG: All Methods	11			
		*Escherichia coli	2	18.18%
		*Gram Negative organism	2	18.18%
		*Growth, referred for identification	2	18.18%
		*Presumptive E. coli	1	9.09%
		*Uricult - Growth (Positive)	4	36.36%
Cultures	11	*Escherichia coli	2	18.18%
Cultures		*Gram Negative organism	2	18.18%
Cultures		*Growth, referred for identification	2	18.18%
Cultures		*Presumptive E. coli	1	9.09%
Cultures		*Uricult - Growth (Positive)	4	36.36%
		Susceptibility Cefpodoxime		



UC-11		* Denotes Correct Response		
AG: All Methods	1			
		*Sensitive	1	100.00%
Susceptibility	1	*Sensitive	1	100.00%
		Susceptibility Cefazidime		
UC-11		* Denotes Correct Response		
AG: All Methods	2			
		*Sensitive	2	100.00%
Susceptibility	2	*Sensitive	2	100.00%
		Susceptibility Cefuroxime		
UC-11		* Denotes Correct Response		
AG: All Methods	3			
		*Resistant	3	100.00%
Susceptibility	3	*Resistant	3	100.00%
		Susceptibility Cephalothin		
UC-11		* Denotes Correct Response		
AG: All Methods	4			
		*Resistant	4	100.00%
Susceptibility	4	*Resistant	4	100.00%
		Susceptibility Ciprofloxacin		
UC-11		* Denotes Correct Response		
AG: All Methods	7			
		*Sensitive	7	100.00%
Susceptibility	7	*Sensitive	7	100.00%
		Susceptibility Gentamicin		
UC-11		* Denotes Correct Response		
AG: All Methods	5			
		*Sensitive	5	100.00%
Susceptibility	5	*Sensitive	5	100.00%
		Susceptibility Levofloxacin		
UC-11		* Denotes Correct Response		
AG: All Methods	4			
		*Sensitive	4	100.00%
Susceptibility	4	*Sensitive	4	100.00%
		Susceptibility Nitrofurantoin		
UC-11		* Denotes Correct Response		
AG: All Methods	6			
		*Intermediate	2	33.33%
		*Resistant	3	50.00%
		Sensitive	1	16.67%
Susceptibility	6	*Intermediate	2	33.33%
Susceptibility		*Resistant	3	50.00%
Susceptibility		Sensitive	1	16.67%
		Susceptibility Norfloxacin/Noroxin		
UC-11		* Denotes Correct Response		
AG: All Methods	1			
		*Sensitive	1	100.00%
Susceptibility	1	*Sensitive	1	100.00%
		Susceptibility Piperacillin		
UC-11		* Denotes Correct Response		
AG: All Methods	1			
		*Sensitive	1	100.00%
Susceptibility	1	*Sensitive	1	100.00%
		Susceptibility Piperacillin/Tazobactam		
UC-11		* Denotes Correct Response		
AG: All Methods	2			
		*Sensitive	2	100.00%
Susceptibility	2	*Sensitive	2	100.00%
		Susceptibility Sulfamethoxazole		
UC-11		* Denotes Correct Response		
AG: All Methods	1			
		*Sensitive	1	100.00%
Susceptibility	1	*Sensitive	1	100.00%
		Susceptibility Tobramycin		
UC-11		* Denotes Correct Response		
AG: All Methods	3			
		*Sensitive	3	100.00%



Susceptibility	3	*Sensitive	3	100.00%
		Susceptibility Trimethoprim		
UC-11		* Denotes Correct Response		
AG: All Methods	1			
		*Sensitive	1	100.00%
Susceptibility	1	*Sensitive	1	100.00%
		Susceptibility Trimeth./Sulfameth.		
UC-11		* Denotes Correct Response		
AG: All Methods	5			
		*Sensitive	5	100.00%
Susceptibility	5	*Sensitive	5	100.00%
		Susceptibility Amikacin		
UC-11		* Denotes Correct Response		
AG: All Methods	1			
		*Sensitive	1	100.00%
Susceptibility	1	*Sensitive	1	100.00%
		Susceptibility Amox/Clavulanate		
UC-11		* Denotes Correct Response		
AG: All Methods	4			
		*Sensitive	4	100.00%
Susceptibility	4	*Sensitive	4	100.00%
		Susceptibility Ampicillin		
UC-11		* Denotes Correct Response		
AG: All Methods	5			
		*Resistant	5	100.00%
Susceptibility	5	*Resistant	5	100.00%
		Susceptibility Ampicillin/Sulbactam		
UC-11		* Denotes Correct Response		
AG: All Methods	2			
		*Intermediate	2	100.00%
Susceptibility	2	*Intermediate	2	100.00%
		Susceptibility Ertapenem		
UC-11		* Denotes Correct Response		
AG: All Methods	1			
		*Sensitive	1	100.00%
Susceptibility	1	*Sensitive	1	100.00%
703 - Urine Drug Screen;713 - Urine Drug Screen (50 ml sample volume)	Labs	Response	Number	Percent
		Barbiturates		
UDS-5		* Denotes Correct Response		
AG: All Methods	137			
		*Negative	137	100.00%
AlereTriage	3	*Negative	3	100.00%
E-Screen/Alere	26	*Negative	26	100.00%
Express Diag. Express Check	9	*Negative	9	100.00%
Exp Results Integ Multi-Drug Screen Cup	1	*Negative	1	100.00%
Alere iCup Drugs of Abuse	12	*Negative	12	100.00%
Instant Technologies i Screen	1	*Negative	1	100.00%
MedTox Profile II	1	*Negative	1	100.00%
Beckman AU	1	*Negative	1	100.00%
Innovacon Multi-Clin	10	*Negative	10	100.00%
MedTox Profile III ER	2	*Negative	2	100.00%
CLIAwaived Inc Rapid Drug Cup	3	*Negative	3	100.00%
ASC One Screen Drug Test Card	1	*Negative	1	100.00%
MedTox Profile V	1	*Negative	1	100.00%
Jant Accutest Multi Drug Screen	3	*Negative	3	100.00%
Amedica Biotech AmediCheck Instant Test Cup	4	*Negative	4	100.00%
Nova Century Scientific Nova Screen DOA Cassette	1	*Negative	1	100.00%
American Screening Corp Reveal Multi-Drug Test Cup	4	*Negative	4	100.00%
Mindray BS-200	2	*Negative	2	100.00%
Medica Easy RA	2	*Negative	2	100.00%
Alere iCassette	2	*Negative	2	100.00%
Roche Hitachi 717	1	*Negative	1	100.00%
US Diagnostics UScreen Cup	6	*Negative	6	100.00%
Bio-Rad Tox/See Rapid Urine Drug Screen	1	*Negative	1	100.00%
Diag Test Group Clarity DOA Screen	1	*Negative	1	100.00%
Noble Medical Insta-Screen Test Cup	3	*Negative	3	100.00%
Hemosure First Sign UDS Dipstick	22	*Negative	22	100.00%



Premier Biotech Bio-Cup	11	*Negative	11	100.00%
UCP Biosciences Home Drug Screening Test Cup	1	*Negative	1	100.00%
Abbott Architect c and ci Series	2	*Negative	2	100.00%
UDS-6		* Denotes Correct Response		
AG: All Methods	137			
		*Negative	137	100.00%
AlereTriage	3	*Negative	3	100.00%
E-Screen/Alere	26	*Negative	26	100.00%
Express Diag. Express Check	9	*Negative	9	100.00%
Exp Results Integ Multi-Drug Screen Cup	1	*Negative	1	100.00%
Alere iCup Drugs of Abuse	12	*Negative	12	100.00%
Instant Technologies i Screen	1	*Negative	1	100.00%
MedTox Profile II	1	*Negative	1	100.00%
Beckman AU	1	*Negative	1	100.00%
Innovacon Multi-Clin	10	*Negative	10	100.00%
MedTox Profile III ER	2	*Negative	2	100.00%
CLIAwaived Inc Rapid Drug Cup	3	*Negative	3	100.00%
ASC One Screen Drug Test Card	1	*Negative	1	100.00%
MedTox Profile V	1	*Negative	1	100.00%
Jant Accutest Multi Drug Screen	3	*Negative	3	100.00%
Amedica Biotech AmediCheck Instant Test Cup	4	*Negative	4	100.00%
Nova Century Scientific Nova Screen DOA Cassette	1	*Negative	1	100.00%
American Screening Corp Reveal Multi-Drug Test Cup	4	*Negative	4	100.00%
Mindray BS-200	2	*Negative	2	100.00%
Medica Easy RA	2	*Negative	2	100.00%
Alere iCassette	2	*Negative	2	100.00%
Roche Hitachi 717	1	*Negative	1	100.00%
US Diagnostics UScreen Cup	6	*Negative	6	100.00%
Bio-Rad Tox/See Rapid Urine Drug Screen	1	*Negative	1	100.00%
Diag Test Group Clarity DOA Screen	1	*Negative	1	100.00%
Noble Medical Insta-Screen Test Cup	3	*Negative	3	100.00%
Hemosure First Sign UDS Dipstick	22	*Negative	22	100.00%
Premier Biotech Bio-Cup	11	*Negative	11	100.00%
UCP Biosciences Home Drug Screening Test Cup	1	*Negative	1	100.00%
Abbott Architect c and ci Series	2	*Negative	2	100.00%
		Benzodiazepines		
UDS-5		* Denotes Correct Response		
AG: All Methods	152			
		*Negative	152	100.00%
AlereTriage	3	*Negative	3	100.00%
Siemens Dim/AR/ES/RxL/Xpand/EXL	1	*Negative	1	100.00%
E-Screen/Alere	26	*Negative	26	100.00%
Express Diag. Express Check	9	*Negative	9	100.00%
Exp Results Integ Multi-Drug Screen Cup	1	*Negative	1	100.00%
Alere iCup Drugs of Abuse	19	*Negative	19	100.00%
Instant Technologies i Screen	1	*Negative	1	100.00%
MedTox Profile II	1	*Negative	1	100.00%
Noble Med Split Specimen Cup	1	*Negative	1	100.00%
Beckman AU	1	*Negative	1	100.00%
Innovacon Multi-Clin	10	*Negative	10	100.00%
MedTox Profile III ER	2	*Negative	2	100.00%
CLIAwaived Inc Rapid Drug Cup	4	*Negative	4	100.00%
ASC One Screen Drug Test Card	1	*Negative	1	100.00%
MedTox Profile V	1	*Negative	1	100.00%
Jant Accutest Multi Drug Screen	3	*Negative	3	100.00%
Amedica Biotech AmediCheck Instant Test Cup	4	*Negative	4	100.00%
Nova Century Scientific Nova Screen DOA Cassette	1	*Negative	1	100.00%
American Screening Corp Reveal Multi-Drug Test Cup	4	*Negative	4	100.00%
Mindray BS-200	2	*Negative	2	100.00%
Medica Easy RA	2	*Negative	2	100.00%
Alere iCassette	2	*Negative	2	100.00%
Roche Hitachi 717	1	*Negative	1	100.00%
US Diagnostics UScreen Cup	7	*Negative	7	100.00%
Bio-Rad Tox/See Rapid Urine Drug Screen	1	*Negative	1	100.00%
Diag Test Group Clarity DOA Screen	1	*Negative	1	100.00%
ThermoScientific Indiko Plus	1	*Negative	1	100.00%
Noble Medical Insta-Screen Test Cup	3	*Negative	3	100.00%
Hemosure First Sign UDS Dipstick	22	*Negative	22	100.00%
Siemens Viva/ProE System	2	*Negative	2	100.00%



TFS DRI Assay	1	*Negative	1	100.00%
Premier Biotech Bio-Cup	11	*Negative	11	100.00%
UCP Biosciences Home Drug Screening Test Cup	1	*Negative	1	100.00%
Abbott Architect c and ci Series	2	*Negative	2	100.00%
UDS-6		* Denotes Correct Response		
AG: All Methods	152			
		*Negative	152	100.00%
AlereTriage	3	*Negative	3	100.00%
Siemens Dim/AR/ES/RxL/Xpand/EXL	1	*Negative	1	100.00%
E-Screen/Alere	26	*Negative	26	100.00%
Express Diag. Express Check	9	*Negative	9	100.00%
Exp Results Integ Multi-Drug Screen Cup	1	*Negative	1	100.00%
Alere iCup Drugs of Abuse	19	*Negative	19	100.00%
Instant Technologies i Screen	1	*Negative	1	100.00%
MedTox Profile II	1	*Negative	1	100.00%
Noble Med Split Specimen Cup	1	*Negative	1	100.00%
Beckman AU	1	*Negative	1	100.00%
Innovacon Multi-Clin	10	*Negative	10	100.00%
MedTox Profile III ER	2	*Negative	2	100.00%
CLIAwaived Inc Rapid Drug Cup	4	*Negative	4	100.00%
ASC One Screen Drug Test Card	1	*Negative	1	100.00%
MedTox Profile V	1	*Negative	1	100.00%
Jant Accutest Multi Drug Screen	3	*Negative	3	100.00%
Amedica Biotech AmediCheck Instant Test Cup	4	*Negative	4	100.00%
Nova Century Scientific Nova Screen DOA Cassette	1	*Negative	1	100.00%
American Screening Corp Reveal Multi-Drug Test Cup	4	*Negative	4	100.00%
Mindray BS-200	2	*Negative	2	100.00%
Medica Easy RA	2	*Negative	2	100.00%
Alere iCassette	2	*Negative	2	100.00%
Roche Hitachi 717	1	*Negative	1	100.00%
US Diagnostics UScreen Cup	7	*Negative	7	100.00%
Bio-Rad Tox/See Rapid Urine Drug Screen	1	*Negative	1	100.00%
Diag Test Group Clarity DOA Screen	1	*Negative	1	100.00%
ThermoScientific Indiko Plus	1	*Negative	1	100.00%
Noble Medical Insta-Screen Test Cup	3	*Negative	3	100.00%
Hemosure First Sign UDS Dipstick	22	*Negative	22	100.00%
Siemens Viva/ProE System	2	*Negative	2	100.00%
TFS DRI Assay	1	*Negative	1	100.00%
Premier Biotech Bio-Cup	11	*Negative	11	100.00%
UCP Biosciences Home Drug Screening Test Cup	1	*Negative	1	100.00%
Abbott Architect c and ci Series	2	*Negative	2	100.00%
		Amphetamines		
UDS-5		* Denotes Correct Response		
AG: All Methods	157			
		*Negative	157	100.00%
AlereTriage	3	*Negative	3	100.00%
Siemens Dim/AR/ES/RxL/Xpand/EXL	2	*Negative	2	100.00%
E-Screen/Alere	23	*Negative	23	100.00%
Express Diag. Express Check	9	*Negative	9	100.00%
Exp Results Integ Multi-Drug Screen Cup	1	*Negative	1	100.00%
Alere iCup Drugs of Abuse	22	*Negative	22	100.00%
Instant Technologies i Screen	1	*Negative	1	100.00%
MedTox Profile II	2	*Negative	2	100.00%
Noble Med Split Specimen Cup	1	*Negative	1	100.00%
Beckman AU	1	*Negative	1	100.00%
Innovacon Multi-Clin	9	*Negative	9	100.00%
MedTox Profile III ER	2	*Negative	2	100.00%
CLIAwaived Inc Rapid Drug Cup	4	*Negative	4	100.00%
ASC One Screen Drug Test Card	1	*Negative	1	100.00%
MedTox Profile V	1	*Negative	1	100.00%
Jant Accutest Multi Drug Screen	3	*Negative	3	100.00%
Amedica Biotech AmediCheck Instant Test Cup	4	*Negative	4	100.00%
Nova Century Scientific Nova Screen DOA Cassette	1	*Negative	1	100.00%
American Screening Corp Reveal Multi-Drug Test Cup	4	*Negative	4	100.00%
Mindray BS-200	2	*Negative	2	100.00%
Medica Easy RA	2	*Negative	2	100.00%
Alere iCassette	5	*Negative	5	100.00%
Roche Hitachi 717	1	*Negative	1	100.00%
US Diagnostics UScreen Cup	7	*Negative	7	100.00%



Bio-Rad Tox/See Rapid Urine Drug Screen	1	*Negative	1	100.00%
Diag Test Group Clarity DOA Screen	1	*Negative	1	100.00%
ThermoScientific Indiko Plus	1	*Negative	1	100.00%
Noble Medical Insta-Screen Test Cup	3	*Negative	3	100.00%
Hemosure First Sign UDS Dipstick	22	*Negative	22	100.00%
Siemens Viva/ProE System	3	*Negative	3	100.00%
TFS DRI Assay	1	*Negative	1	100.00%
Premier Biotech Bio-Cup	11	*Negative	11	100.00%
UCP Biosciences Home Drug Screening Test Cup	1	*Negative	1	100.00%
Abbott Architect c and ci Series	2	*Negative	2	100.00%
UDS-6		* Denotes Correct Response		
AG: All Methods	157			
		*Negative	155	98.73%
		Positive	2	1.27%
AlereTriage	3	*Negative	3	100.00%
Siemens Dim/AR/ES/RxL/Xpand/EXL	2	*Negative	2	100.00%
E-Screen/Alere	23	*Negative	23	100.00%
Express Diag. Express Check	9	*Negative	9	100.00%
Exp Results Integ Multi-Drug Screen Cup	1	*Negative	1	100.00%
Alere iCup Drugs of Abuse	22	*Negative	22	100.00%
Instant Technologies i Screen	1	*Negative	1	100.00%
MedTox Profile II	2	*Negative	2	100.00%
Noble Med Split Specimen Cup	1	*Negative	1	100.00%
Beckman AU	1	*Negative	1	100.00%
Innovacon Multi-Clin	9	*Negative	9	100.00%
MedTox Profile III ER	2	*Negative	2	100.00%
CLIAwaived Inc Rapid Drug Cup	4	*Negative	4	100.00%
ASC One Screen Drug Test Card	1	*Negative	1	100.00%
Jant Accutest Multi Drug Screen	3	*Negative	3	100.00%
Amedica Biotech AmediCheck Instant Test Cup	4	*Negative	4	100.00%
Nova Century Scientific Nova Screen DOA Cassette	1	*Negative	1	100.00%
American Screening Corp Reveal Multi-Drug Test Cup	4	*Negative	4	100.00%
Mindray BS-200	2	*Negative	2	100.00%
Medica Easy RA	2	*Negative	2	100.00%
Alere iCassette	5	*Negative	5	100.00%
Roche Hitachi 717	1	*Negative	1	100.00%
US Diagnostics UScreen Cup	7	*Negative	6	85.71%
Bio-Rad Tox/See Rapid Urine Drug Screen	1	*Negative	1	100.00%
Diag Test Group Clarity DOA Screen	1	*Negative	1	100.00%
ThermoScientific Indiko Plus	1	*Negative	1	100.00%
Noble Medical Insta-Screen Test Cup	3	*Negative	3	100.00%
Hemosure First Sign UDS Dipstick	22	*Negative	22	100.00%
Siemens Viva/ProE System	3	*Negative	3	100.00%
TFS DRI Assay	1	*Negative	1	100.00%
Premier Biotech Bio-Cup	11	*Negative	11	100.00%
UCP Biosciences Home Drug Screening Test Cup	1	*Negative	1	100.00%
Abbott Architect c and ci Series	2	*Negative	2	100.00%
MedTox Profile V	1	Positive	1	100.00%
US Diagnostics UScreen Cup	7	Positive	1	14.29%
		Cannabinoids (THC)		
UDS-5		* Denotes Correct Response		
AG: All Methods	188			
		*Negative	185	98.40%
		Positive	3	1.60%
AlereTriage	3	*Negative	3	100.00%
Siemens Dim/AR/ES/RxL/Xpand/EXL	2	*Negative	2	100.00%
E-Screen/Alere	51	*Negative	51	100.00%
Express Diag. Express Check	9	*Negative	9	100.00%
Exp Results Integ Multi-Drug Screen Cup	1	*Negative	1	100.00%
Alere iCup Drugs of Abuse	21	*Negative	20	95.24%
Instant Technologies i Screen	2	*Negative	2	100.00%
Integrated E-Z Split Key Cup	1	*Negative	1	100.00%
MedTox Profile II	2	*Negative	2	100.00%
Noble Med Split Specimen Cup	1	*Negative	1	100.00%
Beckman AU	1	*Negative	1	100.00%
Innovacon Multi-Clin	10	*Negative	8	80.00%
MedTox Profile III ER	2	*Negative	2	100.00%
CLIAwaived Inc Rapid Drug Cup	4	*Negative	4	100.00%
ASC One Screen Drug Test Card	1	*Negative	1	100.00%



MedTox Profile V	1	*Negative	1	100.00%
Jant Accutest Multi Drug Screen	3	*Negative	3	100.00%
Amedica Biotech AmediCheck Instant Test Cup	4	*Negative	4	100.00%
Nova Century Scientific Nova Screen DOA Cassette	1	*Negative	1	100.00%
American Screening Corp Reveal Multi-Drug Test Cup	4	*Negative	4	100.00%
Mindray BS-200	2	*Negative	2	100.00%
Medica Easy RA	2	*Negative	2	100.00%
Alere iCassette	5	*Negative	5	100.00%
Roche Hitachi 717	1	*Negative	1	100.00%
US Diagnostics UScreen Cup	7	*Negative	7	100.00%
Bio-Rad Tox/See Rapid Urine Drug Screen	1	*Negative	1	100.00%
Diag Test Group Clarity DOA Screen	1	*Negative	1	100.00%
ThermoScientific Indiko Plus	1	*Negative	1	100.00%
Noble Medical Insta-Screen Test Cup	3	*Negative	3	100.00%
CRL Stat Cup	1	*Negative	1	100.00%
Hemosure First Sign UDS Dipstick	22	*Negative	22	100.00%
Siemens Viva/ProE System	3	*Negative	3	100.00%
TFS DRI Assay	1	*Negative	1	100.00%
Premier Biotech Bio-Cup	11	*Negative	11	100.00%
UCP Biosciences Home Drug Screening Test Cup	1	*Negative	1	100.00%
Abbott Architect c and ci Series	2	*Negative	2	100.00%
Alere iCup Drugs of Abuse	21	Positive	1	4.76%
Innovacon Multi-Clin	10	Positive	2	20.00%
UDS-6		* Denotes Correct Response		
AG: All Methods	188			
		*Negative	187	99.47%
		Positive	1	0.53%
AlereTriage	3	*Negative	3	100.00%
Siemens Dim/AR/ES/RxL/Xpand/EXL	2	*Negative	2	100.00%
E-Screen/Alere	51	*Negative	51	100.00%
Express Diag. Express Check	9	*Negative	9	100.00%
Exp Results Integ Multi-Drug Screen Cup	1	*Negative	1	100.00%
Alere iCup Drugs of Abuse	21	*Negative	21	100.00%
Instant Technologies i Screen	2	*Negative	2	100.00%
Integrated E-Z Split Key Cup	1	*Negative	1	100.00%
MedTox Profile II	2	*Negative	2	100.00%
Noble Med Split Specimen Cup	1	*Negative	1	100.00%
Beckman AU	1	*Negative	1	100.00%
Innovacon Multi-Clin	10	*Negative	10	100.00%
MedTox Profile III ER	2	*Negative	2	100.00%
CLIAwaived Inc Rapid Drug Cup	4	*Negative	4	100.00%
ASC One Screen Drug Test Card	1	*Negative	1	100.00%
Jant Accutest Multi Drug Screen	3	*Negative	3	100.00%
Amedica Biotech AmediCheck Instant Test Cup	4	*Negative	4	100.00%
Nova Century Scientific Nova Screen DOA Cassette	1	*Negative	1	100.00%
American Screening Corp Reveal Multi-Drug Test Cup	4	*Negative	4	100.00%
Mindray BS-200	2	*Negative	2	100.00%
Medica Easy RA	2	*Negative	2	100.00%
Alere iCassette	5	*Negative	5	100.00%
Roche Hitachi 717	1	*Negative	1	100.00%
US Diagnostics UScreen Cup	7	*Negative	7	100.00%
Bio-Rad Tox/See Rapid Urine Drug Screen	1	*Negative	1	100.00%
Diag Test Group Clarity DOA Screen	1	*Negative	1	100.00%
ThermoScientific Indiko Plus	1	*Negative	1	100.00%
Noble Medical Insta-Screen Test Cup	3	*Negative	3	100.00%
CRL Stat Cup	1	*Negative	1	100.00%
Hemosure First Sign UDS Dipstick	22	*Negative	22	100.00%
Siemens Viva/ProE System	3	*Negative	3	100.00%
TFS DRI Assay	1	*Negative	1	100.00%
Premier Biotech Bio-Cup	11	*Negative	11	100.00%
UCP Biosciences Home Drug Screening Test Cup	1	*Negative	1	100.00%
Abbott Architect c and ci Series	2	*Negative	2	100.00%
MedTox Profile V	1	Positive	1	100.00%
		Cocaine Metabolites		
UDS-5		* Denotes Correct Response		
AG: All Methods	186			
		*Negative	184	98.92%
		Positive	2	1.08%
AlereTriage	3	*Negative	3	100.00%



Siemens Dim/AR/ES/RxL/Xpand/EXL	2	*Negative	2	100.00%
E-Screen/Alere	51	*Negative	51	100.00%
Express Diag. Express Check	9	*Negative	9	100.00%
Exp Results Integ Multi-Drug Screen Cup	1	*Negative	1	100.00%
Alere iCup Drugs of Abuse	22	*Negative	21	95.45%
Instant Technologies i Screen	2	*Negative	2	100.00%
Integrated E-Z Split Key Cup	1	*Negative	1	100.00%
MedTox Profile II	2	*Negative	2	100.00%
Noble Med Split Specimen Cup	1	*Negative	1	100.00%
Beckman AU	1	*Negative	1	100.00%
Innovacon Multi-Clin	9	*Negative	9	100.00%
MedTox Profile III ER	1	*Negative	1	100.00%
CLIAwaived Inc Rapid Drug Cup	4	*Negative	3	75.00%
ASC One Screen Drug Test Card	1	*Negative	1	100.00%
MedTox Profile V	1	*Negative	1	100.00%
Jant Accutest Multi Drug Screen	3	*Negative	3	100.00%
Amedica Biotech AmediCheck Instant Test Cup	4	*Negative	4	100.00%
Nova Century Scientific Nova Screen DOA Cassette	1	*Negative	1	100.00%
American Screening Corp Reveal Multi-Drug Test Cup	4	*Negative	4	100.00%
Mindray BS-200	2	*Negative	2	100.00%
Medica Easy RA	2	*Negative	2	100.00%
Alere iCassette	5	*Negative	5	100.00%
Roche Hitachi 717	1	*Negative	1	100.00%
US Diagnostics UScreen Cup	7	*Negative	7	100.00%
Bio-Rad Tox/See Rapid Urine Drug Screen	1	*Negative	1	100.00%
Diag Test Group Clarity DOA Screen	1	*Negative	1	100.00%
Noble Medical Insta-Screen Test Cup	3	*Negative	3	100.00%
CRL Stat Cup	1	*Negative	1	100.00%
Hemosure First Sign UDS Dipstick	22	*Negative	22	100.00%
Siemens Viva/ProE System	3	*Negative	3	100.00%
TFS DRI Assay	1	*Negative	1	100.00%
Premier Biotech Bio-Cup	11	*Negative	11	100.00%
UCP Biosciences Home Drug Screening Test Cup	1	*Negative	1	100.00%
Abbott Architect c and ci Series	2	*Negative	2	100.00%
Alere iCup Drugs of Abuse	22	Positive	1	4.55%
CLIAwaived Inc Rapid Drug Cup	4	Positive	1	25.00%
UDS-6		* Denotes Correct Response		
AG: All Methods	186			
		*Positive	184	98.92%
		Negative	2	1.08%
AlereTriage	3	*Positive	3	100.00%
Siemens Dim/AR/ES/RxL/Xpand/EXL	2	*Positive	2	100.00%
E-Screen/Alere	51	*Positive	51	100.00%
Express Diag. Express Check	9	*Positive	9	100.00%
Exp Results Integ Multi-Drug Screen Cup	1	*Positive	1	100.00%
Alere iCup Drugs of Abuse	22	*Positive	22	100.00%
Instant Technologies i Screen	2	*Positive	2	100.00%
Integrated E-Z Split Key Cup	1	*Positive	1	100.00%
MedTox Profile II	2	*Positive	2	100.00%
Noble Med Split Specimen Cup	1	*Positive	1	100.00%
Beckman AU	1	*Positive	1	100.00%
Innovacon Multi-Clin	9	*Positive	9	100.00%
MedTox Profile III ER	1	*Positive	1	100.00%
CLIAwaived Inc Rapid Drug Cup	4	*Positive	3	75.00%
ASC One Screen Drug Test Card	1	*Positive	1	100.00%
Jant Accutest Multi Drug Screen	3	*Positive	3	100.00%
Amedica Biotech AmediCheck Instant Test Cup	4	*Positive	4	100.00%
Nova Century Scientific Nova Screen DOA Cassette	1	*Positive	1	100.00%
American Screening Corp Reveal Multi-Drug Test Cup	4	*Positive	4	100.00%
Mindray BS-200	2	*Positive	2	100.00%
Medica Easy RA	2	*Positive	2	100.00%
Alere iCassette	5	*Positive	5	100.00%
Roche Hitachi 717	1	*Positive	1	100.00%
US Diagnostics UScreen Cup	7	*Positive	7	100.00%
Bio-Rad Tox/See Rapid Urine Drug Screen	1	*Positive	1	100.00%
Diag Test Group Clarity DOA Screen	1	*Positive	1	100.00%
Noble Medical Insta-Screen Test Cup	3	*Positive	3	100.00%
CRL Stat Cup	1	*Positive	1	100.00%
Hemosure First Sign UDS Dipstick	22	*Positive	22	100.00%



Siemens Viva/ProE System	3	*Positive	3	100.00%
TFS DRI Assay	1	*Positive	1	100.00%
Premier Biotech Bio-Cup	11	*Positive	11	100.00%
UCP Biosciences Home Drug Screening Test Cup	1	*Positive	1	100.00%
Abbott Architect c and ci Series	2	*Positive	2	100.00%
CLIAwaived Inc Rapid Drug Cup	4	Negative	1	25.00%
MedTox Profile V	1	Negative	1	100.00%
		Ethanol		
UDS-5		* Denotes Correct Response		
AG: All Methods	7			
		*Negative	7	100.00%
Mindray BS-200	2	*Negative	2	100.00%
Medica Easy RA	1	*Negative	1	100.00%
Roche Hitachi 717	1	*Negative	1	100.00%
Siemens Viva/ProE System	1	*Negative	1	100.00%
TFS DRI Assay	1	*Negative	1	100.00%
Abbott Architect c and ci Series	1	*Negative	1	100.00%
UDS-6		* Denotes Correct Response		
AG: All Methods	7			
		*Negative	7	100.00%
Mindray BS-200	2	*Negative	2	100.00%
Medica Easy RA	1	*Negative	1	100.00%
Roche Hitachi 717	1	*Negative	1	100.00%
Siemens Viva/ProE System	1	*Negative	1	100.00%
TFS DRI Assay	1	*Negative	1	100.00%
Abbott Architect c and ci Series	1	*Negative	1	100.00%
		Acetaminophen		
UDS-5		* Denotes Correct Response		
AG: All Methods	6			
		*Negative	6	100.00%
AlereTriage	3	*Negative	3	100.00%
E-Screen/Alere	1	*Negative	1	100.00%
Noble Med Split Specimen Cup	1	*Negative	1	100.00%
Jant Accutest Multi Drug Screen	1	*Negative	1	100.00%
UDS-6		* Denotes Correct Response		
AG: All Methods	6			
		*Negative	6	100.00%
AlereTriage	3	*Negative	3	100.00%
E-Screen/Alere	1	*Negative	1	100.00%
Noble Med Split Specimen Cup	1	*Negative	1	100.00%
Jant Accutest Multi Drug Screen	1	*Negative	1	100.00%
		MDMA (Ecstasy)		
UDS-5		* Denotes Correct Response		
AG: All Methods	96			
		*Negative	96	100.00%
E-Screen/Alere	4	*Negative	4	100.00%
Express Diag. Express Check	8	*Negative	8	100.00%
Alere iCup Drugs of Abuse	12	*Negative	12	100.00%
Innovacon Multi-Clin	10	*Negative	10	100.00%
CLIAwaived Inc Rapid Drug Cup	3	*Negative	3	100.00%
Jant Accutest Multi Drug Screen	3	*Negative	3	100.00%
Amedica Biotech AmediCheck Instant Test Cup	4	*Negative	4	100.00%
Nova Century Scientific Nova Screen DOA Cassette	1	*Negative	1	100.00%
American Screening Corp Reveal Multi-Drug Test Cup	2	*Negative	2	100.00%
Medica Easy RA	1	*Negative	1	100.00%
Alere iCassette	2	*Negative	2	100.00%
US Diagnostics UScreen Cup	7	*Negative	7	100.00%
Bio-Rad Tox/See Rapid Urine Drug Screen	1	*Negative	1	100.00%
Diag Test Group Clarity DOA Screen	1	*Negative	1	100.00%
Noble Medical Insta-Screen Test Cup	3	*Negative	3	100.00%
Hemosure First Sign UDS Dipstick	22	*Negative	22	100.00%
Premier Biotech Bio-Cup	11	*Negative	11	100.00%
UCP Biosciences Home Drug Screening Test Cup	1	*Negative	1	100.00%
UDS-6		* Denotes Correct Response		
AG: All Methods	96			
		*Negative	96	100.00%
E-Screen/Alere	4	*Negative	4	100.00%
Express Diag. Express Check	8	*Negative	8	100.00%



Alere iCup Drugs of Abuse	12	*Negative	12	100.00%
Innovacon Multi-Clin	10	*Negative	10	100.00%
CLIAwaived Inc Rapid Drug Cup	3	*Negative	3	100.00%
Jant Accutest Multi Drug Screen	3	*Negative	3	100.00%
Amedica Biotech AmediCheck Instant Test Cup	4	*Negative	4	100.00%
Nova Century Scientific Nova Screen DOA Cassette	1	*Negative	1	100.00%
American Screening Corp Reveal Multi-Drug Test Cup	2	*Negative	2	100.00%
Medica Easy RA	1	*Negative	1	100.00%
Alere iCassette	2	*Negative	2	100.00%
US Diagnostics UScreen Cup	7	*Negative	7	100.00%
Bio-Rad Tox/See Rapid Urine Drug Screen	1	*Negative	1	100.00%
Diag Test Group Clarity DOA Screen	1	*Negative	1	100.00%
Noble Medical Insta-Screen Test Cup	3	*Negative	3	100.00%
Hemosure First Sign UDS Dipstick	22	*Negative	22	100.00%
Premier Biotech Bio-Cup	11	*Negative	11	100.00%
UCP Biosciences Home Drug Screening Test Cup	1	*Negative	1	100.00%
		Methadone		
UDS-5		* Denotes Correct Response		
AG: All Methods	118			
		*Negative	116	98.31%
		Positive	2	1.69%
AlereTriage	3	*Negative	3	100.00%
Siemens Dim/AR/ES/RxL/Xpand/EXL	1	*Negative	1	100.00%
E-Screen/Alere	26	*Negative	26	100.00%
Express Diag. Express Check	9	*Negative	9	100.00%
Exp Results Integ Multi-Drug Screen Cup	1	*Negative	1	100.00%
Alere iCup Drugs of Abuse	14	*Negative	13	92.86%
Beckman AU	1	*Negative	1	100.00%
MedTox Profile III ER	2	*Negative	2	100.00%
CLIAwaived Inc Rapid Drug Cup	4	*Negative	3	75.00%
ASC One Screen Drug Test Card	1	*Negative	1	100.00%
MedTox Profile V	1	*Negative	1	100.00%
Jant Accutest Multi Drug Screen	3	*Negative	3	100.00%
Amedica Biotech AmediCheck Instant Test Cup	4	*Negative	4	100.00%
Nova Century Scientific Nova Screen DOA Cassette	1	*Negative	1	100.00%
American Screening Corp Reveal Multi-Drug Test Cup	4	*Negative	4	100.00%
Mindray BS-200	2	*Negative	2	100.00%
Medica Easy RA	2	*Negative	2	100.00%
Roche Hitachi 717	1	*Negative	1	100.00%
US Diagnostics UScreen Cup	6	*Negative	6	100.00%
Bio-Rad Tox/See Rapid Urine Drug Screen	1	*Negative	1	100.00%
Diag Test Group Clarity DOA Screen	1	*Negative	1	100.00%
Noble Medical Insta-Screen Test Cup	3	*Negative	3	100.00%
Hemosure First Sign UDS Dipstick	22	*Negative	22	100.00%
Siemens Viva/ProE System	1	*Negative	1	100.00%
TFS DRI Assay	1	*Negative	1	100.00%
UCP Biosciences Home Drug Screening Test Cup	1	*Negative	1	100.00%
Abbott Architect c and ci Series	2	*Negative	2	100.00%
Alere iCup Drugs of Abuse	14	Positive	1	7.14%
CLIAwaived Inc Rapid Drug Cup	4	Positive	1	25.00%
UDS-6		* Denotes Correct Response		
AG: All Methods	117			
		*Positive	107	91.45%
		Negative	10	8.55%
AlereTriage	3	*Positive	3	100.00%
Siemens Dim/AR/ES/RxL/Xpand/EXL	1	*Positive	1	100.00%
E-Screen/Alere	26	*Positive	26	100.00%
Express Diag. Express Check	9	*Positive	9	100.00%
Exp Results Integ Multi-Drug Screen Cup	1	*Positive	1	100.00%
Alere iCup Drugs of Abuse	14	*Positive	14	100.00%
Beckman AU	1	*Positive	1	100.00%
CLIAwaived Inc Rapid Drug Cup	4	*Positive	2	50.00%
ASC One Screen Drug Test Card	1	*Positive	1	100.00%
Jant Accutest Multi Drug Screen	3	*Positive	2	66.67%
Amedica Biotech AmediCheck Instant Test Cup	4	*Positive	4	100.00%
Nova Century Scientific Nova Screen DOA Cassette	1	*Positive	1	100.00%
American Screening Corp Reveal Multi-Drug Test Cup	4	*Positive	4	100.00%
Mindray BS-200	2	*Positive	2	100.00%
Medica Easy RA	2	*Positive	2	100.00%



US Diagnostics UScreen Cup	5	*Positive	4	80.00%
Bio-Rad Tox/See Rapid Urine Drug Screen	1	*Positive	1	100.00%
Diag Test Group Clarity DOA Screen	1	*Positive	1	100.00%
Noble Medical Insta-Screen Test Cup	3	*Positive	3	100.00%
Hemosure First Sign UDS Dipstick	22	*Positive	20	90.91%
Siemens Viva/ProE System	1	*Positive	1	100.00%
TFS DRI Assay	1	*Positive	1	100.00%
UCP Biosciences Home Drug Screening Test Cup	1	*Positive	1	100.00%
Abbott Architect c and ci Series	2	*Positive	2	100.00%
MedTox Profile III ER	2	Negative	2	100.00%
CLIAwaived Inc Rapid Drug Cup	4	Negative	2	50.00%
MedTox Profile V	1	Negative	1	100.00%
Jant Accutest Multi Drug Screen	3	Negative	1	33.33%
Roche Hitachi 717	1	Negative	1	100.00%
US Diagnostics UScreen Cup	5	Negative	1	20.00%
Hemosure First Sign UDS Dipstick	22	Negative	2	9.09%
		Methamphetamine		
UDS-5		* Denotes Correct Response		
AG: All Methods	156			
		*Negative	156	100.00%
AlereTriage	3	*Negative	3	100.00%
E-Screen/Alere	50	*Negative	50	100.00%
Express Diag. Express Check	9	*Negative	9	100.00%
Exp Results Integ Multi-Drug Screen Cup	1	*Negative	1	100.00%
Alere iCup Drugs of Abuse	21	*Negative	21	100.00%
Instant Technologies i Screen	1	*Negative	1	100.00%
Integrated E-Z Split Key Cup	1	*Negative	1	100.00%
MedTox Profile II	1	*Negative	1	100.00%
Noble Med Split Specimen Cup	1	*Negative	1	100.00%
MedTox Profile III ER	2	*Negative	2	100.00%
CLIAwaived Inc Rapid Drug Cup	4	*Negative	4	100.00%
ASC One Screen Drug Test Card	1	*Negative	1	100.00%
MedTox Profile V	1	*Negative	1	100.00%
Jant Accutest Multi Drug Screen	3	*Negative	3	100.00%
Amedica Biotech AmediCheck Instant Test Cup	4	*Negative	4	100.00%
Nova Century Scientific Nova Screen DOA Cassette	1	*Negative	1	100.00%
American Screening Corp Reveal Multi-Drug Test Cup	4	*Negative	4	100.00%
Alere iCassette	1	*Negative	1	100.00%
US Diagnostics UScreen Cup	7	*Negative	7	100.00%
Bio-Rad Tox/See Rapid Urine Drug Screen	1	*Negative	1	100.00%
Diag Test Group Clarity DOA Screen	1	*Negative	1	100.00%
Noble Medical Insta-Screen Test Cup	3	*Negative	3	100.00%
CRL Stat Cup	1	*Negative	1	100.00%
Hemosure First Sign UDS Dipstick	22	*Negative	22	100.00%
Premier Biotech Bio-Cup	11	*Negative	11	100.00%
UCP Biosciences Home Drug Screening Test Cup	1	*Negative	1	100.00%
UDS-6		* Denotes Correct Response		
AG: All Methods	155			
		*Negative	152	98.06%
		Positive	3	1.94%
AlereTriage	3	*Negative	3	100.00%
E-Screen/Alere	50	*Negative	50	100.00%
Express Diag. Express Check	9	*Negative	9	100.00%
Exp Results Integ Multi-Drug Screen Cup	1	*Negative	1	100.00%
Alere iCup Drugs of Abuse	21	*Negative	20	95.24%
Instant Technologies i Screen	1	*Negative	1	100.00%
Integrated E-Z Split Key Cup	1	*Negative	1	100.00%
MedTox Profile II	1	*Negative	1	100.00%
Noble Med Split Specimen Cup	1	*Negative	1	100.00%
MedTox Profile III ER	2	*Negative	2	100.00%
CLIAwaived Inc Rapid Drug Cup	4	*Negative	4	100.00%
ASC One Screen Drug Test Card	1	*Negative	1	100.00%
Jant Accutest Multi Drug Screen	3	*Negative	2	66.67%
Amedica Biotech AmediCheck Instant Test Cup	4	*Negative	4	100.00%
Nova Century Scientific Nova Screen DOA Cassette	1	*Negative	1	100.00%
American Screening Corp Reveal Multi-Drug Test Cup	4	*Negative	4	100.00%
Alere iCassette	1	*Negative	1	100.00%
US Diagnostics UScreen Cup	6	*Negative	6	100.00%
Bio-Rad Tox/See Rapid Urine Drug Screen	1	*Negative	1	100.00%



Diag Test Group Clarity DOA Screen	1	*Negative	1	100.00%
Noble Medical Insta-Screen Test Cup	3	*Negative	3	100.00%
CRL Stat Cup	1	*Negative	1	100.00%
Hemosure First Sign UDS Dipstick	22	*Negative	22	100.00%
Premier Biotech Bio-Cup	11	*Negative	11	100.00%
UCP Biosciences Home Drug Screening Test Cup	1	*Negative	1	100.00%
Alere iCup Drugs of Abuse	21	Positive	1	4.76%
MedTox Profile V	1	Positive	1	100.00%
Jant Accutest Multi Drug Screen	3	Positive	1	33.33%
		Methaqualone		
UDS-5		* Denotes Correct Response		
AG: All Methods	1			
		*Negative	1	100.00%
Beckman AU	1	*Negative	1	100.00%
UDS-6		* Denotes Correct Response		
AG: All Methods	1			
		*Negative	1	100.00%
Beckman AU	1	*Negative	1	100.00%
		Opiates		
UDS-5		* Denotes Correct Response		
AG: All Methods	188			
		*Negative	184	97.87%
		Positive	4	2.13%
AlereTriage	3	*Negative	3	100.00%
Siemens Dim/AR/ES/RxL/Xpand/EXL	2	*Negative	2	100.00%
E-Screen/Alere	51	*Negative	51	100.00%
Express Diag. Express Check	9	*Negative	9	100.00%
Exp Results Integ Multi-Drug Screen Cup	1	*Negative	1	100.00%
Alere iCup Drugs of Abuse	22	*Negative	20	90.91%
Instant Technologies i Screen	2	*Negative	2	100.00%
Integrated E-Z Split Key Cup	1	*Negative	1	100.00%
MedTox Profile II	1	*Negative	1	100.00%
Noble Med Split Specimen Cup	1	*Negative	1	100.00%
Beckman AU	1	*Negative	1	100.00%
Innovacon Multi-Clin	10	*Negative	10	100.00%
MedTox Profile III ER	2	*Negative	2	100.00%
CLIAwaived Inc Rapid Drug Cup	4	*Negative	3	75.00%
ASC One Screen Drug Test Card	1	*Negative	1	100.00%
MedTox Profile V	1	*Negative	1	100.00%
Jant Accutest Multi Drug Screen	3	*Negative	3	100.00%
Amedica Biotech AmediCheck Instant Test Cup	4	*Negative	4	100.00%
Nova Century Scientific Nova Screen DOA Cassette	1	*Negative	1	100.00%
American Screening Corp Reveal Multi-Drug Test Cup	4	*Negative	4	100.00%
Mindray BS-200	2	*Negative	2	100.00%
Medica Easy RA	2	*Negative	2	100.00%
Alere iCassette	5	*Negative	5	100.00%
Roche Hitachi 717	1	*Negative	1	100.00%
US Diagnostics UScreen Cup	7	*Negative	6	85.71%
Bio-Rad Tox/See Rapid Urine Drug Screen	1	*Negative	1	100.00%
Diag Test Group Clarity DOA Screen	1	*Negative	1	100.00%
ThermoScientific Indiko Plus	1	*Negative	1	100.00%
Noble Medical Insta-Screen Test Cup	3	*Negative	3	100.00%
CRL Stat Cup	1	*Negative	1	100.00%
Hemosure First Sign UDS Dipstick	22	*Negative	22	100.00%
Siemens Viva/ProE System	3	*Negative	3	100.00%
TFS DRI Assay	1	*Negative	1	100.00%
Premier Biotech Bio-Cup	11	*Negative	11	100.00%
UCP Biosciences Home Drug Screening Test Cup	1	*Negative	1	100.00%
Abbott Architect c and ci Series	2	*Negative	2	100.00%
Alere iCup Drugs of Abuse	22	Positive	2	9.09%
CLIAwaived Inc Rapid Drug Cup	4	Positive	1	25.00%
US Diagnostics UScreen Cup	7	Positive	1	14.29%
UDS-6		* Denotes Correct Response		
AG: All Methods	188			
		*Positive	168	89.36%
		Negative	20	10.64%
AlereTriage	3	*Positive	3	100.00%
Siemens Dim/AR/ES/RxL/Xpand/EXL	2	*Positive	2	100.00%
E-Screen/Alere	51	*Positive	34	66.67%



Express Diag. Express Check	9	*Positive	9	100.00%
Exp Results Integ Multi-Drug Screen Cup	1	*Positive	1	100.00%
Alere iCup Drugs of Abuse	22	*Positive	22	100.00%
Instant Technologies i Screen	2	*Positive	2	100.00%
Integrated E-Z Split Key Cup	1	*Positive	1	100.00%
MedTox Profile II	1	*Positive	1	100.00%
Noble Med Split Specimen Cup	1	*Positive	1	100.00%
Beckman AU	1	*Positive	1	100.00%
Innovacon Multi-Clin	10	*Positive	10	100.00%
MedTox Profile III ER	2	*Positive	2	100.00%
CLIAwaived Inc Rapid Drug Cup	4	*Positive	3	75.00%
ASC One Screen Drug Test Card	1	*Positive	1	100.00%
Jant Accutest Multi Drug Screen	3	*Positive	3	100.00%
Amedica Biotech AmediCheck Instant Test Cup	4	*Positive	4	100.00%
Nova Century Scientific Nova Screen DOA Cassette	1	*Positive	1	100.00%
American Screening Corp Reveal Multi-Drug Test Cup	4	*Positive	4	100.00%
Mindray BS-200	2	*Positive	2	100.00%
Medica Easy RA	2	*Positive	2	100.00%
Alere iCassette	5	*Positive	5	100.00%
Roche Hitachi 717	1	*Positive	1	100.00%
US Diagnostics UScreen Cup	7	*Positive	7	100.00%
Bio-Rad Tox/See Rapid Urine Drug Screen	1	*Positive	1	100.00%
Diag Test Group Clarity DOA Screen	1	*Positive	1	100.00%
ThermoScientific Indiko Plus	1	*Positive	1	100.00%
Noble Medical Insta-Screen Test Cup	3	*Positive	3	100.00%
CRL Stat Cup	1	*Positive	1	100.00%
Hemosure First Sign UDS Dipstick	22	*Positive	21	95.45%
Siemens Viva/ProE System	3	*Positive	3	100.00%
TFS DRI Assay	1	*Positive	1	100.00%
Premier Biotech Bio-Cup	11	*Positive	11	100.00%
UCP Biosciences Home Drug Screening Test Cup	1	*Positive	1	100.00%
Abbott Architect c and ci Series	2	*Positive	2	100.00%
E-Screen/Alere	51	Negative	17	33.33%
CLIAwaived Inc Rapid Drug Cup	4	Negative	1	25.00%
MedTox Profile V	1	Negative	1	100.00%
Hemosure First Sign UDS Dipstick	22	Negative	1	4.55%
		Oxycodone		
UDS-5		* Denotes Correct Response		
AG: All Methods	109			
		*Negative	109	100.00%
E-Screen/Alere	18	*Negative	18	100.00%
Express Diag. Express Check	8	*Negative	8	100.00%
Exp Results Integ Multi-Drug Screen Cup	1	*Negative	1	100.00%
Alere iCup Drugs of Abuse	17	*Negative	17	100.00%
Beckman AU	1	*Negative	1	100.00%
Innovacon Multi-Clin	10	*Negative	10	100.00%
MedTox Profile III ER	2	*Negative	2	100.00%
CLIAwaived Inc Rapid Drug Cup	4	*Negative	4	100.00%
ASC One Screen Drug Test Card	1	*Negative	1	100.00%
MedTox Profile V	1	*Negative	1	100.00%
Jant Accutest Multi Drug Screen	3	*Negative	3	100.00%
Amedica Biotech AmediCheck Instant Test Cup	4	*Negative	4	100.00%
Nova Century Scientific Nova Screen DOA Cassette	1	*Negative	1	100.00%
American Screening Corp Reveal Multi-Drug Test Cup	4	*Negative	4	100.00%
Mindray BS-200	2	*Negative	2	100.00%
Medica Easy RA	2	*Negative	2	100.00%
Alere iCassette	2	*Negative	2	100.00%
Roche Hitachi 717	1	*Negative	1	100.00%
US Diagnostics UScreen Cup	7	*Negative	7	100.00%
Bio-Rad Tox/See Rapid Urine Drug Screen	2	*Negative	2	100.00%
Diag Test Group Clarity DOA Screen	1	*Negative	1	100.00%
Noble Medical Insta-Screen Test Cup	3	*Negative	3	100.00%
Hemosure First Sign UDS Cup	1	*Negative	1	100.00%
TFS DRI Assay	1	*Negative	1	100.00%
Premier Biotech Bio-Cup	11	*Negative	11	100.00%
UCP Biosciences Home Drug Screening Test Cup	1	*Negative	1	100.00%
UDS-6		* Denotes Correct Response		
AG: All Methods	109			
		*Negative	104	95.41%



		Positive	5	4.59%
E-Screen/Alere	18	*Negative	18	100.00%
Express Diag. Express Check	8	*Negative	8	100.00%
Exp Results Integ Multi-Drug Screen Cup	1	*Negative	1	100.00%
Alere iCup Drugs of Abuse	17	*Negative	16	94.12%
Beckman AU	1	*Negative	1	100.00%
Innovacon Multi-Clin	10	*Negative	10	100.00%
MedTox Profile III ER	2	*Negative	2	100.00%
CLIAwaived Inc Rapid Drug Cup	4	*Negative	4	100.00%
ASC One Screen Drug Test Card	1	*Negative	1	100.00%
MedTox Profile V	1	*Negative	1	100.00%
Jant Accutest Multi Drug Screen	3	*Negative	3	100.00%
Amedica Biotech AmediCheck Instant Test Cup	4	*Negative	1	25.00%
Nova Century Scientific Nova Screen DOA Cassette	1	*Negative	1	100.00%
American Screening Corp Reveal Multi-Drug Test Cup	4	*Negative	4	100.00%
Mindray BS-200	2	*Negative	2	100.00%
Medica Easy RA	2	*Negative	2	100.00%
Alere iCassette	2	*Negative	2	100.00%
Roche Hitachi 717	1	*Negative	1	100.00%
US Diagnostics UScreen Cup	7	*Negative	6	85.71%
Bio-Rad Tox/See Rapid Urine Drug Screen	2	*Negative	2	100.00%
Diag Test Group Clarity DOA Screen	1	*Negative	1	100.00%
Noble Medical Insta-Screen Test Cup	3	*Negative	3	100.00%
Hemosure First Sign UDS Cup	1	*Negative	1	100.00%
TFS DRI Assay	1	*Negative	1	100.00%
Premier Biotech Bio-Cup	11	*Negative	11	100.00%
UCP Biosciences Home Drug Screening Test Cup	1	*Negative	1	100.00%
Alere iCup Drugs of Abuse	17	Positive	1	5.88%
Amedica Biotech AmediCheck Instant Test Cup	4	Positive	3	75.00%
US Diagnostics UScreen Cup	7	Positive	1	14.29%
		PCP (Phencyclidine)		
UDS-5		* Denotes Correct Response		
AG: All Methods	163			
		*Negative	163	100.00%
AlereTriage	3	*Negative	3	100.00%
Siemens Dim/AR/ES/RxL/Xpand/EXL	1	*Negative	1	100.00%
E-Screen/Alere	51	*Negative	51	100.00%
Express Diag. Express Check	9	*Negative	9	100.00%
Exp Results Integ Multi-Drug Screen Cup	1	*Negative	1	100.00%
Alere iCup Drugs of Abuse	12	*Negative	12	100.00%
Instant Technologies i Screen	2	*Negative	2	100.00%
Integrated E-Z Split Key Cup	1	*Negative	1	100.00%
MedTox Profile II	2	*Negative	2	100.00%
Beckman AU	1	*Negative	1	100.00%
Innovacon Multi-Clin	10	*Negative	10	100.00%
MedTox Profile III ER	2	*Negative	2	100.00%
CLIAwaived Inc Rapid Drug Cup	4	*Negative	4	100.00%
ASC One Screen Drug Test Card	1	*Negative	1	100.00%
MedTox Profile V	1	*Negative	1	100.00%
Jant Accutest Multi Drug Screen	3	*Negative	3	100.00%
Amedica Biotech AmediCheck Instant Test Cup	1	*Negative	1	100.00%
Nova Century Scientific Nova Screen DOA Cassette	1	*Negative	1	100.00%
American Screening Corp Reveal Multi-Drug Test Cup	4	*Negative	4	100.00%
Mindray BS-200	1	*Negative	1	100.00%
Medica Easy RA	1	*Negative	1	100.00%
Alere iCassette	4	*Negative	4	100.00%
Roche Hitachi 717	1	*Negative	1	100.00%
US Diagnostics UScreen Cup	3	*Negative	3	100.00%
Bio-Rad Tox/See Rapid Urine Drug Screen	1	*Negative	1	100.00%
Diag Test Group Clarity DOA Screen	1	*Negative	1	100.00%
Noble Medical Insta-Screen Test Cup	3	*Negative	3	100.00%
CRL Stat Cup	1	*Negative	1	100.00%
Hemosure First Sign UDS Dipstick	22	*Negative	22	100.00%
TFS DRI Assay	1	*Negative	1	100.00%
Premier Biotech Bio-Cup	11	*Negative	11	100.00%
UCP Biosciences Home Drug Screening Test Cup	1	*Negative	1	100.00%
Abbott Architect c and ci Series	2	*Negative	2	100.00%
UDS-6		* Denotes Correct Response		
AG: All Methods	162			



		*Negative	162	100.00%
AlereTriage	3	*Negative	3	100.00%
Siemens Dim/AR/ES/RxL/Xpand/EXL	1	*Negative	1	100.00%
E-Screen/Alere	51	*Negative	51	100.00%
Express Diag. Express Check	9	*Negative	9	100.00%
Exp Results Integ Multi-Drug Screen Cup	1	*Negative	1	100.00%
Alere iCup Drugs of Abuse	12	*Negative	12	100.00%
Instant Technologies i Screen	2	*Negative	2	100.00%
Integrated E-Z Split Key Cup	1	*Negative	1	100.00%
MedTox Profile II	2	*Negative	2	100.00%
Beckman AU	1	*Negative	1	100.00%
Innovacon Multi-Clin	10	*Negative	10	100.00%
MedTox Profile III ER	2	*Negative	2	100.00%
CLIAwaived Inc Rapid Drug Cup	4	*Negative	4	100.00%
ASC One Screen Drug Test Card	1	*Negative	1	100.00%
MedTox Profile V	1	*Negative	1	100.00%
Jant Accutest Multi Drug Screen	3	*Negative	3	100.00%
Amedica Biotech AmediCheck Instant Test Cup	1	*Negative	1	100.00%
Nova Century Scientific Nova Screen DOA Cassette	1	*Negative	1	100.00%
American Screening Corp Reveal Multi-Drug Test Cup	4	*Negative	4	100.00%
Mindray BS-200	1	*Negative	1	100.00%
Medica Easy RA	1	*Negative	1	100.00%
Alere iCassette	4	*Negative	4	100.00%
Roche Hitachi 717	1	*Negative	1	100.00%
US Diagnostics UScreen Cup	3	*Negative	3	100.00%
Bio-Rad Tox/See Rapid Urine Drug Screen	1	*Negative	1	100.00%
Diag Test Group Clarity DOA Screen	1	*Negative	1	100.00%
Noble Medical Insta-Screen Test Cup	3	*Negative	3	100.00%
CRL Stat Cup	1	*Negative	1	100.00%
Hemosure First Sign UDS Dipstick	21	*Negative	21	100.00%
TFS DRI Assay	1	*Negative	1	100.00%
Premier Biotech Bio-Cup	11	*Negative	11	100.00%
UCP Biosciences Home Drug Screening Test Cup	1	*Negative	1	100.00%
Abbott Architect c and ci Series	2	*Negative	2	100.00%
		Propoxyphene		
UDS-5		* Denotes Correct Response		
AG: All Methods	28			
		*Negative	28	100.00%
E-Screen/Alere	4	*Negative	4	100.00%
Exp Results Integ Multi-Drug Screen Cup	1	*Negative	1	100.00%
Alere iCup Drugs of Abuse	9	*Negative	9	100.00%
MedTox Profile II	1	*Negative	1	100.00%
Beckman AU	1	*Negative	1	100.00%
Innovacon Multi-Clin	7	*Negative	7	100.00%
CLIAwaived Inc Rapid Drug Cup	1	*Negative	1	100.00%
MedTox Profile V	1	*Negative	1	100.00%
Alere iCassette	2	*Negative	2	100.00%
Abbott Architect c and ci Series	1	*Negative	1	100.00%
UDS-6		* Denotes Correct Response		
AG: All Methods	28			
		*Negative	28	100.00%
E-Screen/Alere	4	*Negative	4	100.00%
Exp Results Integ Multi-Drug Screen Cup	1	*Negative	1	100.00%
Alere iCup Drugs of Abuse	9	*Negative	9	100.00%
MedTox Profile II	1	*Negative	1	100.00%
Beckman AU	1	*Negative	1	100.00%
Innovacon Multi-Clin	7	*Negative	7	100.00%
CLIAwaived Inc Rapid Drug Cup	1	*Negative	1	100.00%
MedTox Profile V	1	*Negative	1	100.00%
Alere iCassette	2	*Negative	2	100.00%
Abbott Architect c and ci Series	1	*Negative	1	100.00%
		Tricyclic Antidepressants		
UDS-5		* Denotes Correct Response		
AG: All Methods	51			
		*Positive	40	78.43%
		Negative	11	21.57%
AlereTriage	3	*Positive	3	100.00%
Express Diag. Express Check	8	*Positive	5	62.50%
Alere iCup Drugs of Abuse	11	*Positive	7	63.64%



Innovacon Multi-Clin	9	*Positive	7	77.78%
CLIAwaived Inc Rapid Drug Cup	1	*Positive	1	100.00%
MedTox Profile V	1	*Positive	1	100.00%
Jant Accutest Multi Drug Screen	3	*Positive	3	100.00%
Nova Century Scientific Nova Screen DOA Cassette	1	*Positive	1	100.00%
Alere iCassette	2	*Positive	2	100.00%
US Diagnostics UScreen Cup	4	*Positive	4	100.00%
Bio-Rad Tox/See Rapid Urine Drug Screen	2	*Positive	2	100.00%
Noble Medical Insta-Screen Test Cup	3	*Positive	2	66.67%
Mindray 480	1	*Positive	1	100.00%
UCP Biosciences Home Drug Screening Test Cup	1	*Positive	1	100.00%
Express Diag. Express Check	8	Negative	3	37.50%
Alere iCup Drugs of Abuse	11	Negative	4	36.36%
Innovacon Multi-Clin	9	Negative	2	22.22%
Amedica Biotech AmediCheck Instant Test Cup	1	Negative	1	100.00%
Noble Medical Insta-Screen Test Cup	3	Negative	1	33.33%
UDS-6		* Denotes Correct Response		
AG: All Methods	51			
		*Negative	50	98.04%
		Positive	1	1.96%
AlereTriage	3	*Negative	3	100.00%
Express Diag. Express Check	8	*Negative	8	100.00%
Alere iCup Drugs of Abuse	11	*Negative	10	90.91%
Innovacon Multi-Clin	9	*Negative	9	100.00%
CLIAwaived Inc Rapid Drug Cup	1	*Negative	1	100.00%
MedTox Profile V	1	*Negative	1	100.00%
Jant Accutest Multi Drug Screen	3	*Negative	3	100.00%
Amedica Biotech AmediCheck Instant Test Cup	1	*Negative	1	100.00%
Nova Century Scientific Nova Screen DOA Cassette	1	*Negative	1	100.00%
Alere iCassette	2	*Negative	2	100.00%
US Diagnostics UScreen Cup	4	*Negative	4	100.00%
Bio-Rad Tox/See Rapid Urine Drug Screen	2	*Negative	2	100.00%
Noble Medical Insta-Screen Test Cup	3	*Negative	3	100.00%
Mindray 480	1	*Negative	1	100.00%
UCP Biosciences Home Drug Screening Test Cup	1	*Negative	1	100.00%
Alere iCup Drugs of Abuse	11	Positive	1	9.09%
Buprenorphine				
UDS-5		* Denotes Correct Response		
AG: All Methods	46			
		*Positive	40	86.96%
		Negative	6	13.04%
Exp Results Integ Multi-Drug Screen Cup	1	*Positive	1	100.00%
Alere iCup Drugs of Abuse	13	*Positive	8	61.54%
CLIAwaived Inc Rapid Drug Cup	3	*Positive	2	66.67%
ASC One Screen Drug Test Card	1	*Positive	1	100.00%
Jant Accutest Multi Drug Screen	2	*Positive	2	100.00%
Mindray BS-200	2	*Positive	2	100.00%
Medica Easy RA	2	*Positive	2	100.00%
US Diagnostics UScreen Cup	2	*Positive	2	100.00%
Diag Test Group Clarity DOA Screen	1	*Positive	1	100.00%
Noble Medical Insta-Screen Test Cup	3	*Positive	3	100.00%
Hemosure First Sign UDS Cup	1	*Positive	1	100.00%
Siemens Viva/ProE System	3	*Positive	3	100.00%
TFS DRI Assay	1	*Positive	1	100.00%
Premier Biotech Bio-Cup	11	*Positive	11	100.00%
Alere iCup Drugs of Abuse	13	Negative	5	38.46%
CLIAwaived Inc Rapid Drug Cup	3	Negative	1	33.33%
UDS-6		* Denotes Correct Response		
AG: All Methods	46			
		*Negative	45	97.83%
		Positive	1	2.17%
Exp Results Integ Multi-Drug Screen Cup	1	*Negative	1	100.00%
Alere iCup Drugs of Abuse	13	*Negative	12	92.31%
CLIAwaived Inc Rapid Drug Cup	3	*Negative	3	100.00%
ASC One Screen Drug Test Card	1	*Negative	1	100.00%
Jant Accutest Multi Drug Screen	2	*Negative	2	100.00%
Mindray BS-200	2	*Negative	2	100.00%
Medica Easy RA	2	*Negative	2	100.00%
US Diagnostics UScreen Cup	2	*Negative	2	100.00%



Diag Test Group Clarity DOA Screen	1	*Negative	1	100.00%
Noble Medical Insta-Screen Test Cup	3	*Negative	3	100.00%
Hemosure First Sign UDS Cup	1	*Negative	1	100.00%
Siemens Viva/ProE System	3	*Negative	3	100.00%
TFS DRI Assay	1	*Negative	1	100.00%
Premier Biotech Bio-Cup	11	*Negative	11	100.00%
Alere iCup Drugs of Abuse	13	Positive	1	7.69%
		6-Monoacetylmorphine (6 AM)		
UDS-5		* Denotes Correct Response		
AG: All Methods	3			
		*Positive	3	100.00%
Siemens Viva/ProE System	3	*Positive	3	100.00%
UDS-6		* Denotes Correct Response		
AG: All Methods	3			
		*Negative	3	100.00%
Siemens Viva/ProE System	3	*Negative	3	100.00%
		EDDP		
UDS-5		* Denotes Correct Response		
AG: All Methods	13			
		*Negative	13	100.00%
Siemens Viva/ProE System	2	*Negative	2	100.00%
Premier Biotech Bio-Cup	11	*Negative	11	100.00%
UDS-6		* Denotes Correct Response		
AG: All Methods	13			
		*Negative	13	100.00%
Siemens Viva/ProE System	2	*Negative	2	100.00%
Premier Biotech Bio-Cup	11	*Negative	11	100.00%
		Hydrocodone		
UDS-5		* Denotes Correct Response		
AG: All Methods	1			
		*Negative	1	100.00%
Abbott Architect c and ci Series	1	*Negative	1	100.00%
UDS-6		* Denotes Correct Response		
AG: All Methods	1			
		*Negative	1	100.00%
Abbott Architect c and ci Series	1	*Negative	1	100.00%
770 - Viral Antigen Detection;792 - Viral Antigen Detection - Waived Methods;797 - Viral Antigen Detection II	Labs	Response	Number	Percent
		Rotavirus - Aq Detect		
VAD-11		* Denotes Correct Response		
AG: All Methods	9			
		*Negative	9	100.00%
Meridian Immunocard Stat! Rotavirus	5	*Negative	5	100.00%
SAS Rota-Test	4	*Negative	4	100.00%
VAD-12		* Denotes Correct Response		
AG: All Methods	9			
		*Negative	9	100.00%
Meridian Immunocard Stat! Rotavirus	5	*Negative	5	100.00%
SAS Rota-Test	4	*Negative	4	100.00%
VAD-13		* Denotes Correct Response		
AG: All Methods	9			
		*Positive	9	100.00%
Meridian Immunocard Stat! Rotavirus	5	*Positive	5	100.00%
SAS Rota-Test	4	*Positive	4	100.00%
VAD-14		* Denotes Correct Response		
AG: All Methods	9			
		*Positive	9	100.00%
Meridian Immunocard Stat! Rotavirus	5	*Positive	5	100.00%
SAS Rota-Test	4	*Positive	4	100.00%
VAD-15		* Denotes Correct Response		
AG: All Methods	9			
		*Negative	9	100.00%
Meridian Immunocard Stat! Rotavirus	5	*Negative	5	100.00%
SAS Rota-Test	4	*Negative	4	100.00%
		Influenza A - Aq Detect		
VAD-11		* Denotes Correct Response		
AG: All Methods	553			
		*Negative	542	98.01%



		Positive	11	1.99%
Binax Now Flu A	10	*Negative	10	100.00%
Alere BinaxNOW Influenza A+B Card	62	*Negative	60	96.77%
Sekisui Diagnostics OSOM Influenza A & B	19	*Negative	19	100.00%
Quidel QuickVue A+B	35	*Negative	35	100.00%
Quidel Sofia/Sofia 2 Influenza A+B	159	*Negative	155	97.48%
Alere Influenza A & B	15	*Negative	14	93.33%
BD Veritor Flu A+B	68	*Negative	68	100.00%
Henry Schein OneStep + Influenza A&B (distinguish)	9	*Negative	9	100.00%
Medline Influenza A and B	1	*Negative	1	100.00%
Cepheid GeneXpert - Xpert Flu	20	*Negative	20	100.00%
Alere i Influenza A & B	23	*Negative	23	100.00%
Roche Cobas Liat Influenza A/B	7	*Negative	7	100.00%
Sekisui Diagnostics OSOM Ultra Influenza A & B	88	*Negative	84	95.45%
Consult Diag Influenza A & B	22	*Negative	22	100.00%
Lifesign Status Influenza A & B	15	*Negative	15	100.00%
Alere BinaxNOW Influenza A+B Card	62	Positive	2	3.23%
Quidel Sofia/Sofia 2 Influenza A+B	159	Positive	4	2.52%
Alere Influenza A & B	15	Positive	1	6.67%
Sekisui Diagnostics OSOM Ultra Influenza A & B	88	Positive	4	4.55%
VAD-12		* Denotes Correct Response		
AG: All Methods	553			
		*Positive	537	97.11%
		Negative	16	2.89%
Binax Now Flu A	10	*Positive	10	100.00%
Alere BinaxNOW Influenza A+B Card	61	*Positive	58	95.08%
Sekisui Diagnostics OSOM Influenza A & B	19	*Positive	19	100.00%
Quidel QuickVue A+B	35	*Positive	35	100.00%
Quidel Sofia/Sofia 2 Influenza A+B	159	*Positive	150	94.34%
Alere Influenza A & B	15	*Positive	15	100.00%
BD Veritor Flu A+B	68	*Positive	68	100.00%
Henry Schein OneStep + Influenza A&B (distinguish)	9	*Positive	9	100.00%
Cepheid GeneXpert - Xpert Flu	21	*Positive	21	100.00%
Alere i Influenza A & B	23	*Positive	23	100.00%
Roche Cobas Liat Influenza A/B	7	*Positive	7	100.00%
Sekisui Diagnostics OSOM Ultra Influenza A & B	88	*Positive	85	96.59%
Consult Diag Influenza A & B	22	*Positive	22	100.00%
Lifesign Status Influenza A & B	15	*Positive	15	100.00%
Alere BinaxNOW Influenza A+B Card	61	Negative	3	4.92%
Quidel Sofia/Sofia 2 Influenza A+B	159	Negative	9	5.66%
Medline Influenza A and B	1	Negative	1	100.00%
Sekisui Diagnostics OSOM Ultra Influenza A & B	88	Negative	3	3.41%
VAD-13		* Denotes Correct Response		
AG: All Methods	63			
		*Negative	63	100.00%
Alere BinaxNOW Influenza A+B Card	5	*Negative	5	100.00%
Sekisui Diagnostics OSOM Influenza A & B	10	*Negative	10	100.00%
Quidel QuickVue A+B	13	*Negative	13	100.00%
Quidel Sofia/Sofia 2 Influenza A+B	18	*Negative	18	100.00%
Alere Influenza A & B	1	*Negative	1	100.00%
BD Veritor Flu A+B	5	*Negative	5	100.00%
Henry Schein OneStep + Influenza A&B (distinguish)	1	*Negative	1	100.00%
Cepheid GeneXpert - Xpert Flu	5	*Negative	5	100.00%
Alere i Influenza A & B	1	*Negative	1	100.00%
Sekisui Diagnostics OSOM Ultra Influenza A & B	1	*Negative	1	100.00%
Consult Diag Influenza A & B	3	*Negative	3	100.00%
VAD-14		* Denotes Correct Response		
AG: All Methods	63			
		*Positive	63	100.00%
Alere BinaxNOW Influenza A+B Card	5	*Positive	5	100.00%
Sekisui Diagnostics OSOM Influenza A & B	10	*Positive	10	100.00%
Quidel QuickVue A+B	13	*Positive	13	100.00%
Quidel Sofia/Sofia 2 Influenza A+B	18	*Positive	18	100.00%
Alere Influenza A & B	1	*Positive	1	100.00%
BD Veritor Flu A+B	5	*Positive	5	100.00%
Henry Schein OneStep + Influenza A&B (distinguish)	1	*Positive	1	100.00%
Cepheid GeneXpert - Xpert Flu	5	*Positive	5	100.00%
Alere i Influenza A & B	1	*Positive	1	100.00%
Sekisui Diagnostics OSOM Ultra Influenza A & B	1	*Positive	1	100.00%



Consult Diag Influenza A & B	3	*Positive	3	100.00%
VAD-15		* Denotes Correct Response		
AG: All Methods	63			
		*Negative	61	96.83%
		Positive	2	3.17%
Alere BinaxNOW Influenza A+B Card	5	*Negative	5	100.00%
Sekisui Diagnostics OSOM Influenza A & B	10	*Negative	9	90.00%
Quidel QuickVue A+B	13	*Negative	13	100.00%
Quidel Sofia/Sofia 2 Influenza A+B	18	*Negative	17	94.44%
Alere Influenza A & B	1	*Negative	1	100.00%
BD Veritor Flu A+B	5	*Negative	5	100.00%
Henry Schein OneStep + Influenza A&B (distinguish)	1	*Negative	1	100.00%
Cepheid GeneXpert - Xpert Flu	5	*Negative	5	100.00%
Alere i Influenza A & B	1	*Negative	1	100.00%
Sekisui Diagnostics OSOM Ultra Influenza A & B	1	*Negative	1	100.00%
Consult Diag Influenza A & B	3	*Negative	3	100.00%
Sekisui Diagnostics OSOM Influenza A & B	10	Positive	1	10.00%
Quidel Sofia/Sofia 2 Influenza A+B	18	Positive	1	5.56%
		Influenza A/B - Aq Detect		
VAD-11		* Denotes Correct Response		
AG: All Methods	23			
		*Negative	23	100.00%
Alere BinaxNOW Influenza A+B Card	1	*Negative	1	100.00%
Quidel QuickVue	13	*Negative	13	100.00%
Henry Schein OneStep + Influenza A/B	9	*Negative	9	100.00%
VAD-12		* Denotes Correct Response		
AG: All Methods	23			
		*Positive	22	95.65%
		Negative	1	4.35%
Alere BinaxNOW Influenza A+B Card	1	*Positive	1	100.00%
Quidel QuickVue	13	*Positive	13	100.00%
Henry Schein OneStep + Influenza A/B	9	*Positive	8	88.89%
Henry Schein OneStep + Influenza A/B		Negative	1	11.11%
VAD-13		* Denotes Correct Response		
AG: All Methods	2			
		Negative	1	50.00%
		Positive	1	50.00%
Quidel QuickVue	2	Negative	1	50.00%
Quidel QuickVue		Positive	1	50.00%
VAD-14		* Denotes Correct Response		
AG: All Methods	2			
		Positive	2	100.00%
Quidel QuickVue	2	Positive	2	100.00%
VAD-15		* Denotes Correct Response		
AG: All Methods	2			
		Positive	2	100.00%
Quidel QuickVue	2	Positive	2	100.00%
		Influenza B - Aq Detect		
VAD-11		* Denotes Correct Response		
AG: All Methods	552			
		*Negative	544	98.55%
		Positive	8	1.45%
Binax Now Flu B	10	*Negative	10	100.00%
Alere BinaxNOW Influenza A+B Card	62	*Negative	61	98.39%
Sekisui Diagnostics OSOM Influenza A & B	19	*Negative	19	100.00%
Quidel QuickVue A+B	35	*Negative	35	100.00%
Quidel Sofia/Sofia 2 Influenza A+B	159	*Negative	154	96.86%
Alere Influenza A & B	15	*Negative	15	100.00%
BD Veritor Flu A+B	68	*Negative	68	100.00%
Henry Schein OneStep + Influenza A&B (distinguish)	9	*Negative	9	100.00%
Cepheid GeneXpert - Xpert Flu	20	*Negative	20	100.00%
Alere i Influenza A & B	23	*Negative	23	100.00%
Roche Cobas Liat Influenza A/B	7	*Negative	7	100.00%
Sekisui Diagnostics OSOM Ultra Influenza A & B	87	*Negative	86	98.85%
Consult Diag Influenza A & B	22	*Negative	22	100.00%
Lifesign Status Influenza A & B	15	*Negative	15	100.00%
Alere BinaxNOW Influenza A+B Card	62	Positive	1	1.61%
Quidel Sofia/Sofia 2 Influenza A+B	159	Positive	5	3.14%



Medline Influenza A and B	1	Positive	1	100.00%
Sekisui Diagnostics OSOM Ultra Influenza A & B	87	Positive	1	1.15%
VAD-12		* Denotes Correct Response		
AG: All Methods	552			
		*Negative	541	98.01%
		Positive	11	1.99%
Binax Now Flu B	10	*Negative	10	100.00%
Alere BinaxNOW Influenza A+B Card	61	*Negative	60	98.36%
Sekisui Diagnostics OSOM Influenza A & B	19	*Negative	19	100.00%
Quidel QuickVue A+B	35	*Negative	35	100.00%
Quidel Sofia/Sofia 2 Influenza A+B	159	*Negative	157	98.74%
Alere Influenza A & B	15	*Negative	14	93.33%
BD Veritor Flu A+B	68	*Negative	67	98.53%
Henry Schein OneStep + Influenza A&B (distinguish)	9	*Negative	9	100.00%
Medline Influenza A and B	1	*Negative	1	100.00%
Cepheid GeneXpert - Xpert Flu	21	*Negative	21	100.00%
Alere i Influenza A & B	23	*Negative	22	95.65%
Roche Cobas Liat Influenza A/B	7	*Negative	7	100.00%
Sekisui Diagnostics OSOM Ultra Influenza A & B	87	*Negative	85	97.70%
Consult Diag Influenza A & B	22	*Negative	22	100.00%
Lifesign Status Influenza A & B	15	*Negative	12	80.00%
Alere BinaxNOW Influenza A+B Card	61	Positive	1	1.64%
Quidel Sofia/Sofia 2 Influenza A+B	159	Positive	2	1.26%
Alere Influenza A & B	15	Positive	1	6.67%
BD Veritor Flu A+B	68	Positive	1	1.47%
Alere i Influenza A & B	23	Positive	1	4.35%
Sekisui Diagnostics OSOM Ultra Influenza A & B	87	Positive	2	2.30%
Lifesign Status Influenza A & B	15	Positive	3	20.00%
VAD-13		* Denotes Correct Response		
AG: All Methods	63			
		*Positive	62	98.41%
		Negative	1	1.59%
Alere BinaxNOW Influenza A+B Card	5	*Positive	5	100.00%
Sekisui Diagnostics OSOM Influenza A & B	10	*Positive	10	100.00%
Quidel QuickVue A+B	13	*Positive	12	92.31%
Quidel Sofia/Sofia 2 Influenza A+B	18	*Positive	18	100.00%
Alere Influenza A & B	1	*Positive	1	100.00%
BD Veritor Flu A+B	5	*Positive	5	100.00%
Henry Schein OneStep + Influenza A&B (distinguish)	1	*Positive	1	100.00%
Cepheid GeneXpert - Xpert Flu	5	*Positive	5	100.00%
Alere i Influenza A & B	1	*Positive	1	100.00%
Sekisui Diagnostics OSOM Ultra Influenza A & B	1	*Positive	1	100.00%
Consult Diag Influenza A & B	3	*Positive	3	100.00%
Quidel QuickVue A+B	13	Negative	1	7.69%
VAD-14		* Denotes Correct Response		
AG: All Methods	63			
		*Negative	63	100.00%
Alere BinaxNOW Influenza A+B Card	5	*Negative	5	100.00%
Sekisui Diagnostics OSOM Influenza A & B	10	*Negative	10	100.00%
Quidel QuickVue A+B	13	*Negative	13	100.00%
Quidel Sofia/Sofia 2 Influenza A+B	18	*Negative	18	100.00%
Alere Influenza A & B	1	*Negative	1	100.00%
BD Veritor Flu A+B	5	*Negative	5	100.00%
Henry Schein OneStep + Influenza A&B (distinguish)	1	*Negative	1	100.00%
Cepheid GeneXpert - Xpert Flu	5	*Negative	5	100.00%
Alere i Influenza A & B	1	*Negative	1	100.00%
Sekisui Diagnostics OSOM Ultra Influenza A & B	1	*Negative	1	100.00%
Consult Diag Influenza A & B	3	*Negative	3	100.00%
VAD-15		* Denotes Correct Response		
AG: All Methods	63			
		*Positive	61	96.83%
		Negative	2	3.17%
Alere BinaxNOW Influenza A+B Card	5	*Positive	5	100.00%
Sekisui Diagnostics OSOM Influenza A & B	10	*Positive	9	90.00%
Quidel QuickVue A+B	13	*Positive	13	100.00%
Quidel Sofia/Sofia 2 Influenza A+B	18	*Positive	17	94.44%
Alere Influenza A & B	1	*Positive	1	100.00%
BD Veritor Flu A+B	5	*Positive	5	100.00%
Henry Schein OneStep + Influenza A&B (distinguish)	1	*Positive	1	100.00%



Cepheid GeneXpert - Xpert Flu	5	*Positive	5	100.00%
Alere i Influenza A & B	1	*Positive	1	100.00%
Sekisui Diagnostics OSOM Ultra Influenza A & B	1	*Positive	1	100.00%
Consult Diag Influenza A & B	3	*Positive	3	100.00%
Sekisui Diagnostics OSOM Influenza A & B	10	Negative	1	10.00%
Quidel Sofia/Sofia 2 Influenza A+B	18	Negative	1	5.56%
		RSV - Aq Detect		
VAD-11		* Denotes Correct Response		
AG: All Methods	284			
		*Positive	273	96.13%
		Negative	11	3.87%
BD Directigen RSV	1	*Positive	1	100.00%
BD EZ RSV	1	*Positive	1	100.00%
Alere BinaxNOW RSV	172	*Positive	162	94.19%
Quidel QuickVue RSV	26	*Positive	26	100.00%
Wampole Clearview RSV	5	*Positive	5	100.00%
SAS RSV Alert	1	*Positive	1	100.00%
Quidel QuickVue RSV 10	2	*Positive	2	100.00%
Quidel Sofia/Sofia 2 RSV	33	*Positive	33	100.00%
BD Veritor RSV	39	*Positive	38	97.44%
Cepheid GeneXpert - Xpert RSV	1	*Positive	1	100.00%
Alere i RSV	3	*Positive	3	100.00%
Alere BinaxNOW RSV	172	Negative	10	5.81%
BD Veritor RSV	39	Negative	1	2.56%
VAD-12		* Denotes Correct Response		
AG: All Methods	284			
		*Positive	272	95.77%
		Negative	12	4.23%
BD Directigen RSV	1	*Positive	1	100.00%
BD EZ RSV	1	*Positive	1	100.00%
Alere BinaxNOW RSV	172	*Positive	160	93.02%
Quidel QuickVue RSV	26	*Positive	26	100.00%
Wampole Clearview RSV	5	*Positive	5	100.00%
SAS RSV Alert	1	*Positive	1	100.00%
Quidel QuickVue RSV 10	2	*Positive	2	100.00%
Quidel Sofia/Sofia 2 RSV	33	*Positive	33	100.00%
BD Veritor RSV	39	*Positive	39	100.00%
Cepheid GeneXpert - Xpert RSV	1	*Positive	1	100.00%
Alere i RSV	3	*Positive	3	100.00%
Alere BinaxNOW RSV	172	Negative	12	6.98%
VAD-13		* Denotes Correct Response		
AG: All Methods	23			
		*Negative	23	100.00%
BD EZ RSV	1	*Negative	1	100.00%
Alere BinaxNOW RSV	8	*Negative	8	100.00%
Quidel QuickVue RSV	3	*Negative	3	100.00%
Wampole Clearview RSV	4	*Negative	4	100.00%
SAS RSV Alert	1	*Negative	1	100.00%
Quidel QuickVue RSV 10	2	*Negative	2	100.00%
Quidel Sofia/Sofia 2 RSV	3	*Negative	3	100.00%
Cepheid GeneXpert - Xpert RSV	1	*Negative	1	100.00%
VAD-14		* Denotes Correct Response		
AG: All Methods	23			
		*Negative	23	100.00%
BD EZ RSV	1	*Negative	1	100.00%
Alere BinaxNOW RSV	8	*Negative	8	100.00%
Quidel QuickVue RSV	3	*Negative	3	100.00%
Wampole Clearview RSV	4	*Negative	4	100.00%
SAS RSV Alert	1	*Negative	1	100.00%
Quidel QuickVue RSV 10	2	*Negative	2	100.00%
Quidel Sofia/Sofia 2 RSV	3	*Negative	3	100.00%
Cepheid GeneXpert - Xpert RSV	1	*Negative	1	100.00%
VAD-15		* Denotes Correct Response		
AG: All Methods	23			
		*Negative	23	100.00%
BD EZ RSV	1	*Negative	1	100.00%
Alere BinaxNOW RSV	8	*Negative	8	100.00%
Quidel QuickVue RSV	3	*Negative	3	100.00%
Wampole Clearview RSV	4	*Negative	4	100.00%



SAS RSV Alert	1	*Negative	1	100.00%
Quidel QuickVue RSV 10	2	*Negative	2	100.00%
Quidel Sofia/Sofia 2 RSV	3	*Negative	3	100.00%
Cepheid GeneXpert - Xpert RSV	1	*Negative	1	100.00%
		Adenovirus		
VAD-11		* Denotes Correct Response		
AG: All Methods	10			
		*Positive	10	100.00%
SAS Adeno Test	9	*Positive	9	100.00%
Quidel AdenoPlus	1	*Positive	1	100.00%
VAD-12		* Denotes Correct Response		
AG: All Methods	10			
		*Negative	10	100.00%
SAS Adeno Test	9	*Negative	9	100.00%
Quidel AdenoPlus	1	*Negative	1	100.00%
VAD-13		* Denotes Correct Response		
AG: All Methods	10			
		*Negative	10	100.00%
SAS Adeno Test	9	*Negative	9	100.00%
Quidel AdenoPlus	1	*Negative	1	100.00%
VAD-14		* Denotes Correct Response		
AG: All Methods	10			
		*Negative	10	100.00%
SAS Adeno Test	9	*Negative	9	100.00%
Quidel AdenoPlus	1	*Negative	1	100.00%
VAD-15		* Denotes Correct Response		
AG: All Methods	10			
		*Positive	10	100.00%
SAS Adeno Test	9	*Positive	9	100.00%
Quidel AdenoPlus	1	*Positive	1	100.00%
932 - Hepatitis C Antibody (Waived);767 - Hepatitis Markers	Labs	Response	Number	Percent
		HCV (Waived)		
VM-11		* Denotes Correct Response		
AG: All Methods	9			
		*Positive	8	88.89%
		Negative	1	11.11%
OraSure QraQuick HCV Rapid Antibody Test	9	*Positive	8	88.89%
OraSure QraQuick HCV Rapid Antibody Test		Negative	1	11.11%
VM-12		* Denotes Correct Response		
AG: All Methods	9			
		*Negative	8	88.89%
		Positive	1	11.11%
OraSure QraQuick HCV Rapid Antibody Test	9	*Negative	8	88.89%
OraSure QraQuick HCV Rapid Antibody Test		Positive	1	11.11%
767 - Hepatitis Markers	Labs	Response	Number	Percent
		HBsAg		
VM-14		* Denotes Correct Response		
AG: All Methods	1			
		*Negative	1	100.00%
Abbott Architect ci4100	1	*Negative	1	100.00%
VM-15		* Denotes Correct Response		
AG: All Methods	1			
		*Negative	1	100.00%
Abbott Architect ci4100	1	*Negative	1	100.00%
		HCV		
VM-11		* Denotes Correct Response		
AG: All Methods	1			
		*Positive	1	100.00%
Abbott Architect ci4100	1	*Positive	1	100.00%
VM-12		* Denotes Correct Response		
AG: All Methods	1			
		*Negative	1	100.00%
Abbott Architect ci4100	1	*Negative	1	100.00%
VM-13		* Denotes Correct Response		
AG: All Methods	1			
		*Positive	1	100.00%
Abbott Architect ci4100	1	*Positive	1	100.00%



VM-14		* Denotes Correct Response		
AG: All Methods	1			
		*Negative	1	100.00%
Abbott Architect ci4100	1	*Negative	1	100.00%
VM-15		* Denotes Correct Response		
AG: All Methods	1			
		*Negative	1	100.00%
Abbott Architect ci4100	1	*Negative	1	100.00%
		Anti-HBc		
VM-11		* Denotes Correct Response		
AG: All Methods	1			
		*Negative	1	100.00%
Abbott Architect ci4100	1	*Negative	1	100.00%
VM-12		* Denotes Correct Response		
AG: All Methods	1			
		*Positive	1	100.00%
Abbott Architect ci4100	1	*Positive	1	100.00%
VM-13		* Denotes Correct Response		
AG: All Methods	1			
		*Negative	1	100.00%
Abbott Architect ci4100	1	*Negative	1	100.00%
VM-14		* Denotes Correct Response		
AG: All Methods	1			
		*Negative	1	100.00%
Abbott Architect ci4100	1	*Negative	1	100.00%
VM-15		* Denotes Correct Response		
AG: All Methods	1			
		*Negative	1	100.00%
Abbott Architect ci4100	1	*Negative	1	100.00%
		HAV		
VM-11		* Denotes Correct Response		
AG: All Methods	1			
		*Negative	1	100.00%
Abbott Architect ci4100	1	*Negative	1	100.00%
VM-12		* Denotes Correct Response		
AG: All Methods	1			
		*Negative	1	100.00%
Abbott Architect ci4100	1	*Negative	1	100.00%
VM-13		* Denotes Correct Response		
AG: All Methods	1			
		*Negative	1	100.00%
Abbott Architect ci4100	1	*Negative	1	100.00%
VM-14		* Denotes Correct Response		
AG: All Methods	1			
		*Negative	1	100.00%
Abbott Architect ci4100	1	*Negative	1	100.00%
VM-15		* Denotes Correct Response		
AG: All Methods	1			
		*Negative	1	100.00%
Abbott Architect ci4100	1	*Negative	1	100.00%
		Anti-HBc (International)		
VM-11		* Denotes Correct Response		
AG: All Methods	3			
		*Negative	3	100.00%
Intl OCD Vitros ECi/ECiQ	3	*Negative	3	100.00%
VM-12		* Denotes Correct Response		
AG: All Methods	3			
		*Positive	3	100.00%
Intl OCD Vitros ECi/ECiQ	3	*Positive	3	100.00%
VM-13		* Denotes Correct Response		
AG: All Methods	3			
		*Positive	3	100.00%
Intl OCD Vitros ECi/ECiQ	3	*Positive	3	100.00%
VM-14		* Denotes Correct Response		
AG: All Methods	3			
		*Negative	3	100.00%
Intl OCD Vitros ECi/ECiQ	3	*Negative	3	100.00%



VM-15		* Denotes Correct Response		
AG: All Methods	3			
		*Negative	3	100.00%
Intl OCD Vitros ECi/ECiQ	3	*Negative	3	100.00%
		HBsAb (International)		
VM-11		* Denotes Correct Response		
AG: All Methods	3			
		*Negative	3	100.00%
Intl OCD Vitros ECi/ECiQ	3	*Negative	3	100.00%
VM-12		* Denotes Correct Response		
AG: All Methods	3			
		*Negative	3	100.00%
Intl OCD Vitros ECi/ECiQ	3	*Negative	3	100.00%
VM-13		* Denotes Correct Response		
AG: All Methods	3			
		*Negative	3	100.00%
Intl OCD Vitros ECi/ECiQ	3	*Negative	3	100.00%
VM-14		* Denotes Correct Response		
AG: All Methods	3			
		*Negative	3	100.00%
Intl OCD Vitros ECi/ECiQ	3	*Negative	3	100.00%
VM-15		* Denotes Correct Response		
AG: All Methods	3			
		*Positive	3	100.00%
Intl OCD Vitros ECi/ECiQ	3	*Positive	3	100.00%
		HBsAg (International)		
VM-11		* Denotes Correct Response		
AG: All Methods	3			
		*Positive	1	33.33%
		Negative	2	66.67%
Intl OCD Vitros ECi/ECiQ	3	*Positive	1	33.33%
Intl OCD Vitros ECi/ECiQ		Negative	2	66.67%
VM-12		* Denotes Correct Response		
AG: All Methods	3			
		*Positive	3	100.00%
Intl OCD Vitros ECi/ECiQ	3	*Positive	3	100.00%
VM-13		* Denotes Correct Response		
AG: All Methods	3			
		*Positive	3	100.00%
Intl OCD Vitros ECi/ECiQ	3	*Positive	3	100.00%
VM-14		* Denotes Correct Response		
AG: All Methods	3			
		*Negative	3	100.00%
Intl OCD Vitros ECi/ECiQ	3	*Negative	3	100.00%
VM-15		* Denotes Correct Response		
AG: All Methods	3			
		*Negative	3	100.00%
Intl OCD Vitros ECi/ECiQ	3	*Negative	3	100.00%
		HCV (International)		
VM-11		* Denotes Correct Response		
AG: All Methods	3			
		*Positive	3	100.00%
Intl OCD Vitros ECi/ECiQ	3	*Positive	3	100.00%
VM-12		* Denotes Correct Response		
AG: All Methods	3			
		*Negative	3	100.00%
Intl OCD Vitros ECi/ECiQ	3	*Negative	3	100.00%
VM-13		* Denotes Correct Response		
AG: All Methods	3			
		*Positive	3	100.00%
Intl OCD Vitros ECi/ECiQ	3	*Positive	3	100.00%
VM-14		* Denotes Correct Response		
AG: All Methods	3			
		*Negative	3	100.00%
Intl OCD Vitros ECi/ECiQ	3	*Negative	3	100.00%
VM-15		* Denotes Correct Response		
AG: All Methods	3			



		*Negative	3	100.00%
Intl OCD Vitros ECI/ECIQ	3	*Negative	3	100.00%
771 - Vaginal Pathogens I - Affirm VP	Labs	Response	Number	Percent
		Candida sp. - Non Culture		
VP-11		* Denotes Correct Response		
AG: All Methods	50			
		*Positive	49	98.00%
		Negative	1	2.00%
Affirm VP	49	*Positive	48	97.96%
BD Max	1	*Positive	1	100.00%
Affirm VP	49	Negative	1	2.04%
VP-12		* Denotes Correct Response		
AG: All Methods	50			
		*Negative	50	100.00%
Affirm VP	49	*Negative	49	100.00%
BD Max	1	*Negative	1	100.00%
VP-13		* Denotes Correct Response		
AG: All Methods	50			
		*Negative	49	98.00%
		Positive	1	2.00%
Affirm VP	49	*Negative	49	100.00%
BD Max	1	Positive	1	100.00%
VP-14		* Denotes Correct Response		
AG: All Methods	50			
		*Negative	50	100.00%
Affirm VP	49	*Negative	49	100.00%
BD Max	1	*Negative	1	100.00%
VP-15		* Denotes Correct Response		
AG: All Methods	50			
		*Positive	49	98.00%
		Negative	1	2.00%
Affirm VP	49	*Positive	49	100.00%
BD Max	1	Negative	1	100.00%
		G. vaginalis - Non Culture		
VP-11		* Denotes Correct Response		
AG: All Methods	50			
		*Positive	50	100.00%
Affirm VP	49	*Positive	49	100.00%
BD Max	1	*Positive	1	100.00%
VP-12		* Denotes Correct Response		
AG: All Methods	50			
		*Negative	50	100.00%
Affirm VP	49	*Negative	49	100.00%
BD Max	1	*Negative	1	100.00%
VP-13		* Denotes Correct Response		
AG: All Methods	50			
		*Positive	49	98.00%
		Negative	1	2.00%
Affirm VP	49	*Positive	49	100.00%
BD Max	1	Negative	1	100.00%
VP-14		* Denotes Correct Response		
AG: All Methods	50			
		*Negative	50	100.00%
Affirm VP	49	*Negative	49	100.00%
BD Max	1	*Negative	1	100.00%
VP-15		* Denotes Correct Response		
AG: All Methods	50			
		*Negative	49	98.00%
		Positive	1	2.00%
Affirm VP	49	*Negative	49	100.00%
BD Max	1	Positive	1	100.00%
		T. vaginalis - Non Culture Tech		
VP-11		* Denotes Correct Response		
AG: All Methods	50			
		*Negative	50	100.00%
Affirm VP	49	*Negative	49	100.00%
BD Max	1	*Negative	1	100.00%
VP-12		* Denotes Correct Response		



American Academy of Family Physicians
11400 Tomahawk Creek Parkway
Leawood, KS 66211-2672
1-800-274-7911
www.aafp.org/pt pt@aafp.org

Proficiency Testing Event Summary

AAFP PT 2018-C

AG: All Methods	50			
		*Positive	50	100.00%
Affirm VP	49	*Positive	49	100.00%
BD Max	1	*Positive	1	100.00%
VP-13		* Denotes Correct Response		
AG: All Methods	50			
		*Negative	50	100.00%
Affirm VP	49	*Negative	49	100.00%
BD Max	1	*Negative	1	100.00%
VP-14		* Denotes Correct Response		
AG: All Methods	50			
		*Negative	50	100.00%
Affirm VP	49	*Negative	49	100.00%
BD Max	1	*Negative	1	100.00%
VP-15		* Denotes Correct Response		
AG: All Methods	50			
		*Negative	50	100.00%
Affirm VP	49	*Negative	49	100.00%
BD Max	1	*Negative	1	100.00%