



666 - Afinion Glycohemoglobin	AFF-1																								
Glycohemoglobin	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%
PG: Afinion Glycohemoglobin	158	6.33	0.45	4.90-7.70	7.14																				
747 - Auto Diff 7 (AcT 5 Diff & Horiba 5-Part Diff)	BCX-1					BCX-2					BCX-3					BCX-4					BCX-5				
White Blood Cell Count	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%
PG: Beckman Coulter AcT 5 diff	24	3.96	0.13	3.40-4.60	3.25	24	24.74	0.83	21.00-28.40	3.36	24	5.03	0.19	4.30-5.80	3.60	24	17.35	0.63	14.70-19.90	3.64	24	7.60	0.25	6.50-8.70	3.29
PG: Horiba Pentra 60	16	3.94	0.08	3.30-4.50	2.05	16	24.87	0.58	21.20-28.60	2.33	16	4.93	0.13	4.20-5.60	2.65	16	17.40	0.41	14.80-20.00	2.36	16	7.49	0.16	6.40-8.60	2.13
SG: Select Horiba Instruments (5-part Diff)	21	3.94	0.09	3.30-4.50	2.31	21	24.86	0.56	21.20-28.60	2.25	21	4.95	0.13	4.30-5.80	2.40	21	16.88	2.19	14.40-19.40	12.96	20	7.50	0.16	6.40-8.60	2.13
SG: International Horiba Instruments	24	3.95	0.09	3.40-4.60	2.25	24	24.70	0.90	21.00-28.40	3.64	24	4.96	0.12	4.30-5.80	2.40	24	16.77	2.19	14.30-19.30	13.10	23	7.52	0.17	6.40-8.60	2.27
Lymphocyte	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%
PG: Beckman Coulter AcT 5 diff	24	28.94	3.56	18.20-39.60	12.32	23	35.67	3.87	24.10-47.30	10.84	24	39.68	3.31	29.80-49.60	8.34	23	21.93	2.23	15.20-28.60	10.18	24	21.08	2.04	15.00-27.20	9.67
PG: Horiba Pentra 60	16	30.17	3.90	18.50-41.90	12.91	15	38.99	3.64	28.10-49.90	9.33	15	42.55	3.79	31.10-53.90	8.92	16	24.54	1.55	19.90-29.10	6.33	16	24.71	2.31	17.80-31.60	9.35
SG: Select Horiba Instruments (5-part Diff)	21	30.60	3.95	18.80-42.50	12.91	20	37.92	3.90	26.20-49.60	10.32	20	41.76	3.74	30.60-53.00	8.95	21	24.04	1.69	18.90-29.10	7.04	21	23.91	2.75	15.60-32.10	11.51
SG: International Horiba Instruments	24	31.57	4.82	17.10-46.10	15.25	23	39.57	6.00	21.60-57.60	15.15	23	42.17	3.83	30.70-53.70	9.08	24	25.40	4.36	12.30-38.50	17.17	24	24.75	3.49	14.30-35.30	14.07
Monocyte	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%
PG: Beckman Coulter AcT 5 diff	21	1.44	0.69	0.00-3.52	48.13	24	4.71	1.78	0.00-10.06	37.88	24	1.93	1.07	0.00-5.14	55.39	24	2.85	1.02	0.00-5.92	35.89	24	2.19	1.23	0.00-5.88	56.21
PG: Horiba Pentra 60	15	1.23	0.43	0.00-2.52	34.96	16	3.91	1.01	0.89-6.93	25.73	16	1.29	0.69	0.00-3.35	53.10	15	1.78	0.54	0.17-3.39	30.11	16	1.23	0.62	0.00-3.08	50.08
SG: Select Horiba Instruments (5-part Diff)	20	1.17	0.51	0.00-2.69	43.16	21	3.70	1.13	0.31-7.09	30.51	21	1.24	0.67	0.00-3.24	53.87	20	1.74	0.60	0.00-3.55	34.71	21	1.11	0.60	0.00-2.91	53.96
SG: International Horiba Instruments	23	1.14	0.48	0.00-2.59	42.28	24	3.60	1.27	0.00-7.42	35.39	24	1.23	0.65	0.00-3.18	52.85	23	1.77	0.61	0.00-3.59	34.29	24	1.06	0.59	0.00-2.84	56.04
Granulocyte	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%
PG: Beckman Coulter AcT 5 diff	23	56.26	3.92	44.50-68.10	6.96	23	49.67	4.09	37.40-62.00	8.23	23	48.05	3.22	38.30-57.70	6.71	23	59.45	2.80	51.10-67.90	4.71	23	58.71	3.02	49.60-67.80	5.14
PG: Horiba Pentra 60	16	60.19	4.76	45.90-74.50	7.91	16	53.79	4.23	41.10-66.50	7.86	16	52.33	4.62	38.40-66.20	8.83	16	70.10	1.83	64.60-75.60	2.61	16	68.64	8.40	43.40-93.80	12.24
SG: Select Horiba Instruments (5-part Diff)	21	59.76	4.58	46.10-73.50	7.66	21	54.76	4.22	42.10-67.50	7.70	21	52.99	4.41	39.80-66.20	8.32	21	70.65	1.89	65.00-76.40	2.67	21	70.01	7.80	46.60-93.40	11.14
SG: International Horiba Instruments	24	58.61	5.33	42.60-74.60	9.10	24	53.28	5.84	35.80-70.80	10.94	24	52.51	4.42	39.20-65.80	8.42	24	68.84	5.22	53.10-84.50	7.59	24	69.36	7.55	46.80-92.10	10.88
Eosinophil	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%
PG: Beckman Coulter AcT 5 diff	22	7.30	1.25	3.50-11.10	17.12	24	5.53	0.81	3.10-7.90	14.73	22	4.36	0.72	2.20-6.60	16.36	24	3.57	0.52	2.00-5.20	14.44	21	4.00	0.46	2.60-5.40	11.50
PG: Horiba Pentra 60	16	8.04	2.02	1.90-14.10	25.25	16	3.77	0.81	1.40-6.20	21.32	15	4.08	0.72	1.90-6.30	17.56	16	2.88	0.38	1.80-4.00	13.10	16	2.86	0.64	1.00-4.80	22.07
SG: Select Horiba Instruments (5-part Diff)	21	8.14	1.87	2.50-13.70	23.09	21	3.89	0.81	1.50-6.30	20.77	20	4.18	0.71	2.10-6.30	16.90	21	2.90	0.35	1.90-3.90	12.07	21	2.84	0.58	1.10-4.50	20.71
SG: International Horiba Instruments	24	8.35	2.20	1.80-15.00	26.19	24	3.82	0.82	1.30-6.30	21.58	23	4.23	0.68	2.20-6.20	16.19	24	2.90	0.36	1.80-4.00	12.41	24	2.89	0.63	1.00-4.80	21.72
Basophil	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%
PG: Beckman Coulter AcT 5 diff	23	4.97	0.49	3.49-6.44	9.88	23	3.86	0.25	3.12-4.61	6.43	23	5.68	0.47	4.26-7.11	8.37	23	11.69	0.53	10.09-13.29	4.56	23	13.83	0.70	11.74-15.91	5.03
PG: Horiba Pentra 60	16	0.27	0.06	0.09-0.45	22.38	16	0.20	0.00	0.20-0.20	0.00	16	0.30	0.00	0.30-0.30	0.00	16	0.59	0.03	0.52-0.67	4.21	15	0.69	0.04	0.58-0.79	5.12
SG: Select Horiba Instruments (5-part Diff)	21	0.26	0.06	0.09-0.44	22.52	21	0.20	0.00	0.20-0.20	0.00	21	0.30	0.00	0.30-0.30	0.00	21	0.60	0.04	0.48-0.71	6.45	20	0.69	0.04	0.57-0.80	5.34
SG: International Horiba Instruments	24	0.25	0.06	0.08-0.43	23.15																				



Blood Gas - pH	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%
PG: Radiometer ABL80 FLEX CO-OX	5	7.26	0.01	7.22-7.30	0.14	5	7.26	0.01	7.22-7.30	0.14	5	7.28	0.01	7.24-7.32	0.08	5	6.88	0.01	6.84-6.92	0.09	5	7.26	0.01	7.22-7.30	0.14
AG: All Group	5	7.26	0.02	7.22-7.30	0.23	5	7.26	0.01	7.22-7.30	0.14	5	7.29	0.03	7.25-7.33	0.34	5	6.88	0.02	6.84-6.92	0.29	5	7.26	0.01	7.22-7.30	0.21
Blood Gas - pO2	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%
PG: Radiometer ABL80 FLEX CO-OX	5	56.50	4.66	43.00-71.00	8.25	5	188.00	4.32	175.00-201.00	2.29	5	225.25	6.95	204.00-246.00	3.07	5	35.75	5.19	20.00-52.00	14.44	5	56.50	3.70	46.00-68.00	6.49
AG: All Group	5	59.96	8.72	34.00-86.00	14.50	5	188.00	4.32	175.00-201.00	2.29	5	221.74	9.89	192.00-252.00	4.46	5	40.42	11.37	6.00-74.00	28.50	5	59.90	8.25	35.00-85.00	13.67
714 - Complete Cardiac Markers;715 - Basic Cardiac Markers;720 - Cardiac Markers I;740 - D-Dimer	CDM-1N					CDM-2N					CDM-3N					CDM-4N					CDM-5N				
	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%
PG: Abbott i-Stat	14	892.71	51.01	739.70-1045.70	5.71	14	39.00	3.82	27.50-50.50	9.79	14	1922.21	148.28	1477.40-2367.00	7.71	14	3494.86	269.03	2687.80-4302.00	7.70	14	434.43	45.13	299.00-569.80	10.39
PG: AlereTriage	29	673.17	101.08	369.90-976.50	15.02	29	40.52	5.75	23.30-57.80	14.20	29	1252.31	162.62	764.40-1740.20	12.99	29	2167.93	288.13	1303.50-3032.30	13.29	29	344.76	49.39	196.60-493.00	14.32
AG: All Group	46	738.66	135.25	333.00-1144.50	18.31	46	39.62	6.90	18.90-60.30	17.42	45	1468.80	346.01	430.80-2506.80	23.56	45	2608.06	674.19	585.50-4630.70	25.85	45	371.80	67.28	170.00-573.60	18.10
CK-MB (Quantitative)	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%
PG: Abbott i-Stat	5	19.70	0.00	19.70-19.70	0.00	5	4.50	0.00	4.50-4.50	0.00	5	37.60	0.00	37.60-37.60	0.00	5	85.20	0.00	85.20-85.20	0.00	5	11.20	0.00	11.20-11.20	0.00
PG: AlereTriage	35	12.78	2.25	6.10-19.60	17.58	35	2.50	0.48	1.10-3.90	19.20	35	22.09	3.71	11.00-33.20	16.79	35	40.03	7.39	17.80-62.20	18.48	35	7.02	0.95	4.20-9.80	13.57
PG: TOSOH AIA	5	35.00	0.00	35.00-35.00	0.00	5	8.10	0.00	8.10-8.10	0.00	5	62.10	0.00	62.10-62.10	0.00	5	111.80	0.00	111.80-111.80	0.00	5	21.40	0.00	21.40-21.40	0.00
PG: Beckman Coulter Dxl Series	5	35.20	0.00	35.20-35.20	0.00	5	7.60	0.00	7.60-7.60	0.00	5	62.10	0.00	62.10-62.10	0.00	5	112.00	0.00	112.00-112.00	0.00	5	19.10	0.00	19.10-19.10	0.00
PG: Intl OCD Vitros 250,350,400,500,750,950	5	20.40	0.00	20.40-20.40	0.00	5	4.80	0.00	4.80-4.80	0.00	5	35.60	0.00	35.60-35.60	0.00	5	66.20	0.00	66.20-66.20	0.00	5	11.80	0.00	11.80-11.80	0.00
Myoglobin (Quantitative)	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%
PG: AlereTriage	30	150.20	27.68	67.20-233.20	18.43	30	36.17	6.65	16.30-56.20	18.37	30	225.77	27.77	142.50-309.10	12.30	30	357.93	53.30	198.00-517.80	14.89	30	92.14	14.70	48.00-136.20	15.96
AG: All Group	33	148.97	27.68	66.00-232.00	18.58	33	35.73	6.63	15.80-55.60	18.57	33	226.54	30.14	136.10-316.90	13.31	33	366.81	65.58	170.10-563.50	17.88	33	91.18	14.82	46.70-135.70	16.25
Troponin I (Quantitative)	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%
PG: Abbott i-Stat	20	10.52	1.34	6.49-14.54	12.75	19	0.38	0.04	0.25-0.50	11.11	19	21.95	2.90	13.24-30.66	13.22	19	39.12	6.82	18.68-59.58	17.42	19	4.89	0.56	3.22-6.57	11.40
PG: AlereTriage	33	1.70	0.34	0.70-2.71	19.64	36	0.05	0.00	0.05-0.05	0.00	35	5.54	0.96	2.66-8.41	17.31	35	12.77	2.13	6.39-19.16	16.66	35	0.53	0.11	0.19-0.87	21.47
PG: TOSOH AIA	5	28.03	0.00	28.03-28.03	0.00	5	1.22	0.00	1.22-1.22	0.00	5	52.35	0.00	52.35-52.35	0.00	5	97.47	0.00	97.47-97.47	0.00	5	13.86	0.00	13.86-13.86	0.00
PG: Intl Abbott i-Stat	5	9785.00	159.26	9307.20-10262.8	1.63	5	377.50	28.72	291.33-463.67	7.61															
AG: All Group	35	1.67	0.35	0.62-2.73	20.96	38	0.05	0.00	0.04-0.06	6.53	37	5.38	1.16	1.91-8.85	21.51	37	12.40	2.62	4.53-20.27	21.15	37	0.53	0.11	0.20-0.86	20.93
D-dimer (Quantitative)	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%
PG: AlereTriage	48	933.46	87.34	671.00-1195.00	9.36	48	512.00	65.93	314.20-709.80	12.88	46	1332.83	136.44	923.50-1742.10	10.24	46	2061.52	208.63	1436.00-2688.00	10.12	46	727.65	57.81	554.30-901.10	7.94
PG: Beckman AU	5	220.00	0.00	220.00-220.00	0.00	5	75.00	0.00	75.00-75.00	0.00	5	370.00	0.00	370.00-370.00	0.00	5	755.00	0.00	755.00-755.00	0.00	5	135.00	0.00	135.00-135.00	0.00
PG: Roche H 232	5	500.00	50.00	350.00-650.00	10.00	5	283.33	28.87	196.70-369.90	10.21															
AG: All Group	53	896.72	162.26	410.00-1384.00	18.10	53	494.36	104.94	179.60-809.20	21.23	47	1312.34	194.77	728.00-1896.60	14.85	47	2033.72	280.89	1191.00-2877.00	13.81	47	715.04	103.65	404.00-1026.00	14.50
680 - Coagulation - Plasma	CG-1					CG-2					CG-3					CG-4					CG-5				
INR	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%
PG: Diagnostica Stago STA Satellite	16	2.66	0.14	2.30-3.20	5.56	16	1.11	0.05	1.00-1.30	4.55	16	1.35	0.12	0.90-1.70	9.23	16	4.								



AG: All Group	16	2.51	0.55	0.80-4.20	22.00	16	4.17	0.48	2.80-5.60	11.43	16	4.68	0.61	2.90-6.50	12.98	15	2.40	0.49	0.90-3.90	20.42	16	3.75	0.52	2.10-5.30	14.05
Potassium (K+)	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%
PG: Abaxis Piccolo (Non-Waived)	9	2.51	0.16	2.00-3.00	6.40	9	5.32	0.22	4.80-5.80	4.15	9	7.09	0.28	6.60-7.60	3.94	9	2.38	0.13	1.90-2.90	5.42	9	3.99	0.16	3.50-4.50	4.00
PG: Abaxis Piccolo (Waived)	30	2.45	0.18	2.00-3.00	7.20	30	5.25	0.23	4.80-5.80	4.34															
PG: Abbott Architect	6	2.67	0.05	2.20-3.20	1.85	6	5.38	0.04	4.90-5.90	0.74	6	7.22	0.07	6.70-7.70	1.11	6	2.60	0.00	2.10-3.10	0.00	6	3.97	0.05	3.50-4.50	1.25
PG: Alfa Wasserman Alera/Axcel	23	2.80	0.11	2.30-3.30	3.93	24	5.44	0.12	4.90-5.90	2.22	24	7.67	0.16	7.20-8.20	2.08	24	2.60	0.10	2.10-3.10	3.85	24	3.93	0.06	3.40-4.40	1.54
PG: Clinical Data Envoy 500	6	2.57	0.05	2.10-3.10	1.92	6	5.30	0.13	4.80-5.80	2.45	6	6.95	0.15	6.50-7.50	2.14	6	2.58	0.10	2.10-3.10	3.85	6	3.95	0.14	3.50-4.50	3.50
PG: Siemens Dimension/AR/ES/RxL/Xpand/EXL	15	2.54	0.05	2.00-3.00	2.00	15	5.35	0.06	4.90-5.90	1.11	15	7.21	0.07	6.70-7.70	0.97	15	2.55	0.05	2.10-3.10	1.92	15	3.95	0.05	3.50-4.50	1.25
PG: Beckman AU	17	2.58	0.05	2.10-3.10	1.92	17	5.25	0.09	4.70-5.70	1.73	17	7.04	0.12	6.50-7.50	1.71	16	2.60	0.05	2.10-3.10	1.92	16	3.90	0.05	3.40-4.40	1.28
PG: Roche Cobas Integra	6	2.60	0.00	2.10-3.10	0.00	6	5.37	0.05	4.90-5.90	0.93	6	7.17	0.05	6.70-7.70	0.69	6	2.60	0.00	2.10-3.10	0.00	6	4.00	0.00	3.50-4.50	0.00
PG: OCD Vitros 250,350,400,500,750,950	8	2.88	0.05	2.40-3.40	1.72	8	5.68	0.07	5.20-6.20	1.23	8	7.69	0.08	7.20-8.20	1.04	8	2.71	0.04	2.20-3.20	1.48	8	4.18	0.09	3.70-4.70	2.14
SG: Select Beckman Coulter Instruments (Chem)	21	2.57	0.06	2.10-3.10	2.31	21	5.25	0.09	4.70-5.70	1.73	21	7.07	0.13	6.60-7.60	1.83	20	2.59	0.06	2.10-3.10	2.31	20	3.90	0.05	3.40-4.40	1.28
SG: Select Roche Instruments	7	2.61	0.04	2.10-3.10	1.54	7	5.37	0.05	4.90-5.90	0.93	7	7.19	0.07	6.70-7.70	0.97	7	2.61	0.04	2.10-3.10	1.54	7	4.01	0.04	3.50-4.50	1.00
SG: Select OCD Vitros Instruments	10	2.84	0.08	2.30-3.30	2.86	10	5.61	0.15	5.10-6.10	2.68	10	7.61	0.19	7.10-8.10	2.50	10	2.67	0.10	2.20-3.20	3.33	10	4.12	0.14	3.60-4.60	3.41
SG: Abaxis Piccolo Instruments	39	2.46	0.17	2.00-3.00	7.20	39	5.27	0.23	4.80-5.80	4.34	9	7.09	0.28	6.60-7.60	3.94	9	2.38	0.13	1.90-2.90	5.42	9	3.99	0.16	3.50-4.50	4.00
SG: Select Abbott Architect Instruments	9	2.68	0.04	2.20-3.20	1.48	9	5.37	0.05	4.90-5.90	0.93	9	7.22	0.07	6.70-7.70	0.97	9	2.60	0.00	2.10-3.10	0.00	9	3.98	0.04	3.50-4.50	1.00
SG: International OCD Vitros Instruments	13	2.85	0.09	2.40-3.40	3.10	13	5.62	0.14	5.10-6.10	2.50	13	7.64	0.19	7.10-8.10	2.50	12	2.68	0.09	2.20-3.20	3.33	12	4.13	0.13	3.60-4.60	3.17
SG: International Siemens Instruments	16	2.56	0.08	2.10-3.10	3.08	16	5.35	0.06	4.90-5.90	1.11	16	7.21	0.07	6.70-7.70	0.97	16	2.55	0.05	2.10-3.10	1.92	16	3.95	0.05	3.50-4.50	1.25
SG: International Abaxis Piccolo Instruments	42	2.47	0.17	2.00-3.00	6.80	43	5.25	0.24	4.80-5.80	4.53	13	7.02	0.31	6.50-7.50	4.43	13	2.32	0.20	1.80-2.80	8.70	12	3.97	0.15	3.50-4.50	3.75
SG: Select Gen Chem Instruments	59	2.68	0.11	2.20-3.20	4.07	60	5.34	0.11	4.80-5.80	2.08	61	7.30	0.31	6.80-7.80	4.25	59	2.58	0.06	2.10-3.10	2.31	59	3.92	0.07	3.40-4.40	1.54
Chloride (CL-)	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%
PG: Abaxis Piccolo (Non-Waived)	9	96.89	1.45	92.00-102.00	1.55	9	107.11	1.76	102.00-112.00	1.68	9	117.00	2.06	111.00-123.00	1.79	9	90.89	1.27	86.00-96.00	1.43	9	97.89	2.37	93.00-103.00	2.45
PG: Abaxis Piccolo (Waived)	31	96.42	1.57	91.00-101.00	1.67	31	106.00	2.25	101.00-111.00	2.17															
PG: Abbott Architect	6	90.67	1.21	86.00-96.00	1.32	6	102.00	1.26	97.00-107.00	1.27	6	112.00	1.10	106.00-118.00	0.98	6	86.33	1.03	82.00-90.00	1.16	6	94.17	0.75	89.00-99.00	0.85
PG: Alfa Wasserman Alera/Axcel	23	89.30	2.04	85.00-93.00	2.25	23	104.03	2.09	99.00-109.00	2.02	23	119.51	1.76	114.00-126.00	1.50	23	82.88	1.43	79.00-87.00	1.69	23	93.40	1.79	88.00-98.00	1.94
PG: Clinical Data Envoy 500	6	91.08	2.15	86.00-96.00	2.42	5	107.20	4.09	102.00-112.00	3.83	5	119.25	2.75	113.00-125.00	2.35	5	84.28	2.26	80.00-88.00	2.74	5	95.64	3.74	91.00-101.00	3.85
PG: Siemens Dimension/AR/ES/RxL/Xpand/EXL	15	82.20	1.08	78.00-86.00	1.34	15	102.67	1.05	98.00-108.00	0.97	15	117.27	1.44	111.00-123.00	1.20	15	80.93	0.46	77.00-85.00	0.62	15	91.60	0.63	87.00-97.00	0.65
PG: Beckman AU	16	87.19	1.51	83.00-91.00	1.72	16	100.31	1.78	95.00-105.00	1.80	16	111.38	1.89	105.00-117.00	1.71	16	84.81	2.37	81.00-89.00	2.82	16	91.88	2.58	87.00-97.00	2.83
PG: Roche Cobas Integra	6	84.83	0.98	81.00-89.00	1.18	6	105.00	3.46	100.00-110.00	3.33	6	118.83	5.19	113.00-125.00	4.37	6	85.33	2.34	81.00-89.00	2.71	6	95.00	2.97	90.00-100.00	3.16
PG: OCD Vitros 250,350,400,500,750,950	8	89.50	1.07	86.00-95.00	1.22	8	102.00	1.31	97.00-107.00	1.27	8	114.25	1.67	108.00-120.00	1.49	8	85.38	1.30	81.00-89.00	1.53	8	93.50	1.41	89.00-99.00	1.49
SG: Select Beckman Coulter Instruments (Chem)	20	87.70	1.75	84.00-92.00	2.05	20	100.45	1.70	95.00-105.00	1.70	20	111.40	1.82	105.00-117.00	1.62	20	85.05	2.19	81.00-89.00	2.59	20	92.20	2.40	87.00-97.00	2.61
SG: Select Roche Instruments	7	85.43	1.81	81.00-89.00	2.12	7	104.57	3.36	100.00-110.00	3.24	7	117.86	5.40	112.00-124.00	4.58	7	85.43	2.15	81.00-89.00	2.47	7	94.86	2.73	90.00-100.00	2.84
SG: Select OCD Vitros Instruments	10	89.20	1.48	85.00-93.00	1.69	10	101.70	2.11	97.00-107.00	2.06	10	114.00	2.11	108.00-120.00	1.84	10	85.00	1.83	81.00-89.00	2.12	10	93.20	1.99	88.00-98.00	2.15
SG: Abaxis Piccolo Instruments	40	96.53	1.54	92.00-102.00	1.55	40	106.25	2.18	101.00-111.00	2.08	9	117.00	2.06	111.00-123.00	1.79	9	90.89	1.27	86.00-96.00	1.43	9	97.89	2.37	93.00-103.00	2.45
SG: Select Abbott Architect Instruments	9	91.00	1.12	86.00-96.00	1.21	9	102.11	1.05	97.00-107.00	1.08	9	112.22	0.97	106.00-118.00	0.89	9	86.44	0.88	82.00-90.00	1.05	9	94.44	0.73	89.00-99.00	0.74
SG: International OCD Vitros Instruments	13	89.15	1.68	85.00-93.00	1.91	13	101.62	2.18	97.00-107.00	2.16	13	113.92	2.25	108.00-120.00	2.02	13	86.23	3.92	82.00-90.00	4.53	13	92.92	2.78	88.00-98.00	3.01
SG: International Siemens Instruments	16	82.69	2.21	79.00-87.00	2.65	16	102.75	1.06	98.00-108.00	1.07	16	117.38	1.46	111.00-123.00	1.28	16	81.13	0.89	77.00-85.00	1.11	16	91.75	0.86	87.00-97.00	0.98
SG: International Abaxis Piccolo Instruments	44	96.52	1.62	92.00-102.00	1.65	44	106.45	2.28	101.00-111.00	2.17	13	116.77	2.28	111.00-123.00	1.97	13	90.46	1.39	86.00-95.00	1.56	13	98.23	2.20	93.00-103.00	2.24
SG: Select Gen Chem Instruments	60	89.24	2.23	85.00-93.00	2.47	59	103.02	3.30	98.00-108.00	3.20	58	115.96	5.28	110.00-122.00	4.57	59	84.41	2.38	80.00-88.00	2.86	59	93.50	2.61	89.00-99.00	2.77
Free T4	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%
PG: Abbott Architect	6	1.53	0.12	1.16-1.90	8.11	6	2.63	0.16	2.17-3.10	5.90	6	3.60	0.25	2.85-4.35	6.93	6	0.78	0.05	0.63-0.93	6.34	6	1.74	0.11	1.41-2.06	6.17
PG: Beckman ACCESS	15	1.02	0.11	0.70-1.33	10.45	15	2.18	0.17	1.68-2.68	7.67	15	2.66	0.14	2.23-3.09	5.37	15	0.84	0.09	0.56-1.12	11.05	15	1.70	0.14	1.28-2.11	8.14
PG: Siemens Dimension/AR/ES/RxL/Xpand/EXL	11	1.24	0.17	0.74-1.73	13.32	11	3.16	0.18	2.61-3.71	5.81	11	4.05	0.33	3.06-5.05	8.20	11	0.87	0.12	0.50-1.24	14.15	11	2.26	0.12	1.89-2.64	5.49
PG: Qualigen FastPak Thyroid	5	0.86	0.08	0.62-1.10	9.26	6	3.42	0.35	2.35-4.48	10.37	6	4.47	0.36	3.38-5.55	8.09	5	0.99	0.07	0.77-1.21	7.28	6	2.53	0.34	1.50-3.57	13.60
PG: TOSOH AIA	12	1.56	0.10	1.26-1.85	6.29	12	3.43	0.25	2.68-4.18	7.28	12	4.28	0.30	3.39-5.17	6.92	12	0.97	0.10	0.66-1.28	10.55	12	2.48	0.15	2.03-2.93	6.07
SG: Select Beckman Coulter Instruments (Chem)	20	1.02	0.09	0.74-1.29	9.16	20	2.19	0.15	1.75-2.62	6.68	20	2.70	0.15	2.25-3.14	5.47	20	0.85	0.08	0.60-1.10	9.87	20	1.70	0.12	1.34-2.06	7.08
SG: Select Siemens Instruments	14	1.26	0.15	0.80-1.72	12.18	14	3.04	0.29	2.15-3.92	9.71	14	3.89	0.44	2.59-5.20	11.18	14	0.91	0.14	0.49-1.33	15.24	14	2.22	0.15	1.77-2.66	6.66
PG: Roche Cobas e Series	6	1.60	0.07	1.41-1.80	4.09	6	2.93	0.14	2.50-3.35	4.82	6	3.95	0.24	3.23-4.67	6.08	6	0.98	0							



Calcium	Labs	Mean	SD	Range	CV %	Labs	Mean	SD	Range	CV %	Labs	Mean	SD	Range	CV %	Labs	Mean	SD	Range	CV %	Labs	Mean	SD	Range	CV %
PG: Abaxis Piccolo (Non-Waived)	8	8.03	0.18	7.00-9.00	2.25	8	10.58	0.17	9.60-11.60	1.60	8	12.15	0.28	11.20-13.20	2.30	8	7.63	0.34	6.60-8.60	4.61	8	9.15	0.27	8.20-10.20	2.93
PG: Abaxis Piccolo (Waived)	32	7.90	0.25	6.90-8.90	3.16	32	10.60	0.53	9.60-11.60	5.00															
PG: Abbott Architect	6	8.27	0.15	7.30-9.30	1.81	6	10.70	0.18	9.70-11.70	1.68	6	13.02	0.17	12.00-14.00	1.31	6	7.30	0.11	6.30-8.30	1.51	6	9.00	0.18	8.00-10.00	2.00
PG: Alfa Wasserman Alera/Axcel	24	8.10	0.15	7.10-9.10	1.85	24	10.70	0.18	9.70-11.70	1.68	24	12.68	0.22	11.70-13.70	1.73	24	7.55	0.13	6.50-8.50	1.73	24	9.23	0.17	8.20-10.20	1.85
PG: Clinical Data Envoy 500	6	8.22	0.10	7.20-9.20	1.22	6	10.18	0.22	9.20-11.20	2.16	6	12.44	0.30	11.40-13.40	2.42	6	6.60	0.19	5.60-7.60	2.88	6	8.52	0.20	7.50-9.50	2.35
PG: Siemens Dimension/AR/ES/RxL/Xpand/EXL	12	7.49	0.28	6.50-8.50	3.73	13	10.46	0.28	9.50-11.50	2.67	13	12.67	0.30	11.70-13.70	2.36	13	7.19	0.37	6.20-8.20	5.14	13	8.82	0.30	7.80-9.80	3.41
PG: Beckman AU	17	8.26	0.26	7.30-9.30	3.13	17	10.54	0.24	9.50-11.50	2.29	17	12.72	0.32	11.70-13.70	2.52	17	7.23	0.47	6.20-8.20	6.53	17	8.79	0.46	7.80-9.80	5.23
PG: Roche Cobas Integra	6	7.90	0.24	6.90-8.90	3.04	6	10.62	0.30	9.60-11.60	2.83	6	12.78	0.39	11.80-13.80	3.05	6	7.27	0.15	6.30-8.30	2.05	6	8.98	0.20	8.00-10.00	2.22
PG: OCD Vitros 250,350,400,500,750,950	8	8.06	0.23	7.10-9.10	2.84	8	10.79	0.22	9.80-11.80	2.04	8	12.96	0.23	12.00-14.00	1.77	8	7.19	0.16	6.20-8.20	2.22	8	9.09	0.17	8.10-10.10	1.87
SG: Select Beckman Coulter Instruments (Chem)	20	8.20	0.28	7.20-9.20	3.41	20	10.52	0.23	9.50-11.50	2.19	20	12.70	0.30	11.70-13.70	2.36	20	7.24	0.43	6.20-8.20	5.97	20	8.81	0.42	7.80-9.80	4.77
SG: Select Roche Instruments	7	7.91	0.22	6.90-8.90	2.78	7	10.63	0.28	9.60-11.60	2.64	7	12.83	0.38	11.80-13.80	2.97	7	7.29	0.15	6.30-8.30	2.05	7	9.00	0.19	8.00-10.00	2.11
SG: Select OCD Vitros Instruments	10	8.05	0.21	7.10-9.10	2.59	10	10.75	0.23	9.80-11.80	2.13	10	12.90	0.27	11.90-13.90	2.09	10	7.19	0.14	6.20-8.20	1.94	10	9.04	0.18	8.00-10.00	2.00
SG: Abaxis Piccolo Instruments	40	7.93	0.24	6.90-8.90	3.04	40	10.59	0.48	9.60-11.60	4.53	8	12.15	0.28	11.20-13.20	2.30	8	7.63	0.34	6.60-8.60	4.61	8	9.15	0.27	8.20-10.20	2.93
SG: Select Abbott Architect Instruments	9	8.27	0.12	7.30-9.30	1.45	9	10.72	0.18	9.70-11.70	1.68	9	13.08	0.17	12.10-14.10	1.30	9	7.28	0.12	6.30-8.30	1.64	9	9.00	0.15	8.00-10.00	1.67
SG: International OCD Vitros Instruments	13	8.02	0.20	7.00-9.00	2.50	13	10.72	0.23	9.70-11.70	2.15	13	12.86	0.26	11.90-13.90	2.02	12	7.18	0.15	6.20-8.20	2.08	13	8.95	0.37	8.00-10.00	4.11
SG: International Abaxis Piccolo Instruments	44	7.94	0.24	6.90-8.90	3.04	44	10.59	0.46	9.60-11.60	4.34	12	12.15	0.30	11.10-13.10	2.48	12	7.61	0.28	6.60-8.60	3.68	12	9.13	0.22	8.10-10.10	2.42
SG: Select Gen Chem Instruments	60	8.17	0.20	7.20-9.20	2.44	60	10.60	0.26	9.60-11.60	2.45	60	12.73	0.30	11.70-13.70	2.36	60	7.30	0.39	6.30-8.30	5.21	60	8.98	0.37	8.00-10.00	4.11
T3 Uptake	Labs	Mean	SD	Range	CV %	Labs	Mean	SD	Range	CV %	Labs	Mean	SD	Range	CV %	Labs	Mean	SD	Range	CV %	Labs	Mean	SD	Range	CV %
AG: All Group	6	0.30	0.06	0.12-0.49	20.24	7	0.46	0.05	0.32-0.60	10.22	7	0.42	0.05	0.28-0.56	11.03	7	0.54	0.09	0.26-0.82	17.34	7	0.50	0.07	0.28-0.72	14.56
Creatinine	Labs	Mean	SD	Range	CV %	Labs	Mean	SD	Range	CV %	Labs	Mean	SD	Range	CV %	Labs	Mean	SD	Range	CV %	Labs	Mean	SD	Range	CV %
PG: Abaxis Piccolo (Non-Waived)	7	0.63	0.05	0.33-0.93	8.33	9	3.11	0.09	2.64-3.58	2.99	9	4.83	0.17	4.11-5.55	3.54	9	0.56	0.13	0.30-0.90	23.75	9	1.87	0.13	1.57-2.17	7.06
PG: Abaxis Piccolo (Waived)	30	0.59	0.13	0.30-0.90	21.69	31	3.11	0.17	2.64-3.58	5.50															
PG: Abbott Architect	6	0.84	0.03	0.50-1.10	3.75	6	3.34	0.05	2.84-3.84	1.52	6	5.33	0.08	4.53-6.13	1.51	6	0.52	0.02	0.20-0.80	4.00	6	1.83	0.04	1.53-2.13	2.22
PG: Alfa Wasserman Alera/Axcel	24	0.68	0.05	0.40-1.00	7.35	24	3.11	0.15	2.64-3.58	4.89	24	4.62	0.21	3.93-5.31	4.57	24	0.66	0.05	0.40-1.00	7.14	24	1.95	0.07	1.65-2.25	3.50
PG: Clinical Data Envoy 500	6	0.58	0.02	0.30-0.90	3.33	6	3.12	0.11	2.65-3.59	3.65	6	4.71	0.16	4.00-5.42	3.40	6	0.58	0.05	0.30-0.90	8.33	6	1.89	0.01	1.59-2.19	0.79
PG: Siemens Dimension/AR/ES/RxL/Xpand/EXL	10	0.65	0.07	0.40-1.00	10.62	9	3.09	0.13	2.63-3.55	4.27	10	4.86	0.14	4.13-5.59	2.86	10	0.47	0.04	0.20-0.80	8.30	10	1.77	0.07	1.47-2.07	4.44
PG: Siemens Dimension/AR/ES/RxL/Xpand/EXL / IFCC	5	0.65	0.03	0.35-0.95	5.23	5	3.18	0.08	2.70-3.70	2.50	5	4.93	0.10	4.19-5.67	2.07	5	0.44	0.03	0.14-0.74	7.50	5	1.83	0.08	1.53-2.13	4.44
PG: Beckman AU	17	0.65	0.06	0.40-1.00	8.62	17	3.07	0.09	2.61-3.53	3.00	17	4.69	0.13	3.99-5.39	2.81	16	0.47	0.05	0.20-0.80	10.00	16	1.81	0.06	1.51-2.11	3.48
PG: Roche Cobas Integra	6	0.61	0.02	0.30-0.90	3.93	6	2.95	0.08	2.51-3.39	2.81	6	4.26	0.19	3.62-4.90	4.51	6	0.43	0.05	0.10-0.70	12.50	6	1.76	0.05	1.46-2.06	2.84
PG: OCD Vitros 250,350,400,500,750,950	8	0.69	0.01	0.40-1.00	2.03	8	3.15	0.10	2.68-3.62	3.23	8	4.67	0.10	3.97-5.37	2.18	8	0.50	0.01	0.20-0.80	2.80	7	1.84	0.07	1.54-2.14	3.89
SG: Select Beckman Coulter Instruments (Chem)	21	0.65	0.06	0.40-1.00	8.77	21	3.08	0.09	2.62-3.54	2.90	21	4.71	0.12	4.00-5.42	2.61	20	0.47	0.06	0.20-0.80	12.77	20	1.81	0.06	1.51-2.11	3.43
SG: Select Roche Instruments	7	0.59	0.05	0.30-0.90	8.33	7	2.91	0.12	2.47-3.35	4.14	7	4.24	0.19	3.60-4.88	4.52	7	0.41	0.07	0.10-0.70	17.50	7	1.74	0.08	1.44-2.04	4.71
SG: Select OCD Vitros Instruments	10	0.69	0.03	0.40-1.00	4.64	10	3.13	0.10	2.66-3.60	3.23	10	4.64	0.11	3.94-5.34	2.26	10	0.49	0.03	0.20-0.80	6.73	9	1.83	0.06	1.53-2.13	3.39
SG: Abaxis Piccolo Instruments	37	0.60	0.12	0.30-0.90	20.00	40	3.11	0.16	2.64-3.58	5.16	9	4.83	0.17	4.11-5.55	3.54	9	0.56	0.13	0.30-0.90	23.75	9	1.87	0.13	1.57-2.17	7.06
SG: Select Abbott Architect Instruments	9	0.84	0.03	0.50-1.10	3.75	9	3.32	0.07	2.82-3.82	2.12	9	5.30	0.11	4.50-6.09	2.08	9	0.51	0.01	0.20-0.80	4.00	9	1.82	0.03	1.52-2.12	1.67
SG: International OCD Vitros Instruments	13	0.68	0.04	0.40-1.00	5.71	13	3.14	0.09	2.67-3.61	2.96	13	4.65	0.10	3.95-5.35	2.22	12	0.49	0.03	0.19-0.79	6.12	11	1.83	0.06	1.50-2.10	3.33
SG: International Abaxis Piccolo Instruments	41	0.60	0.12	0.30-0.90	20.00	44	3.12	0.16	2.65-3.59	5.22	13	4.84	0.17	4.11-5.57	3.54	13	0.56	0.11	0.26-0.86	20.18	13	1.89	0.13	1.60-2.20	6.84
SG: Select Gen Chem Instruments	61	0.68	0.09	0.40-1.00	12.50	61	3.14	0.14	2.67-3.61	4.52	61	4.77	0.28	4.05-5.49	5.87	60	0.56	0.10	0.26-0.86	16.96	60	1.88	0.09	1.60-2.20	4.74
Thyroxine (T4)	Labs	Mean	SD	Range	CV %	Labs	Mean	SD	Range	CV %	Labs	Mean	SD	Range	CV %	Labs	Mean	SD	Range	CV %	Labs	Mean	SD	Range	CV %
PG: Beckman ACCESS	7	8.22	0.29	6.60-9.80	3.54	7	9.97	0.60	8.00-12.00	6.00	7	13.31	0.63	10.60-16.00	4.74	7	1.83	0.26	0.80-2.80	14.44	7	6.42	0.31	5.10-7.70	4.84
PG: TOSOH AIA	5	9.72	0.50	7.80-11.60	5.15	5	9.08	0.95	7.30-10.90	10.44	5	13.00	0.91	10.40-15.60	7.00	5	2.08	0.19	1.10-3.10	9.05	5	6.08	0.33	4.90-7.30	5.41
SG: Select Beckman Coulter Instruments (Chem)	8	8.15	0.34	6.50-9.70	4.20	8	9.99	0.56	8.00-12.00	5.60	8	13.40	0.63	10.70-16.10	4.70	8	1.89	0.29	0.90-2.90	15.26	8	6.38	0.31	5.10-7.70	4.84
SG: Select Siemens Instruments	5	12.00	0.91	9.60-14.40	7.58	5	9.78	0.50	7.80-11.80	5.10	5	14.35	0.84	11.50-17.30	5.83	5	2.05	0.51	1.10-3.10	24.29	5	5.98	0.50	4.80-7.20	8.33
PG: Intl OCD Vitros Eci/ECiQ	5	8.06	0.57	6.50-9.70	7.16	5	9.73	0.72	7.80-11.60	7.42	5	14.60	1.04	11.70-17.50	7.12	5	1.98	0.16	1.00-3.00	8.00	5	5.74	0.41	4.60-6.80	7.19
AG: All Group	17	9.39	1.44	7.50-11.30	15.32	17	9.76	0.63	7.80-11.80	6.43	17	13.43	0.74	10.70-16.10	5.52	17	2.05	0.29	1.10-3.10	13.81	17	6.25	0.31	5.00-7.60	4.92
Glucose	Labs	Mean	SD	Range	CV %	Labs	Mean	SD	Range	CV %	Labs	Mean	SD	Range	CV %	Labs	Mean	SD	Range	CV %	Labs	Mean	SD	Range	CV %
PG: Abaxis Piccolo (Non-Waived)	9	61.78	1.79	56.00-68.00	2.90	9	162.00	1.66	146.00-178.00	1.05	9	233.00	3.00	210.00-256.00	1.29	9	54.89	1.17	49.00-61.00	2.18	9	107.89	1.36	97.00-119.00	1.30
PG: Abaxis Piccolo (Waived)	31	61.34	1.05	55.00-67.00	1.80	30	162.04	1.65	146.00-178.00	0.99															
PG: Abbott Architect	6	56.83	0.75	51.00-63.00	1.40	6	159.83	1.84	144.00-176.00	1.13	6	234.50	3.45	212.00-259											



AG: All Group	10	1.30	0.40	0.10-2.50	30.77	8	1.52	0.33	0.50-2.50	22.00	8	1.98	0.50	0.50-3.50	25.00	8	0.80	0.23	0.10-1.50	28.75	8	1.20	0.25	0.50-2.00	20.83
TSH	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%
PG: Abbott Architect	6	2.46	0.14	2.03-2.89	5.77	6	4.03	0.27	3.22-4.84	6.67	6	6.36	0.38	5.21-7.51	6.01	6	0.37	0.03	0.29-0.45	7.30	6	2.24	0.15	1.80-2.68	6.56
PG: Beckman ACCESS	18	2.70	0.21	2.08-3.32	7.63	18	4.35	0.39	3.18-5.52	8.99	18	6.54	0.86	3.97-9.11	13.12	18	0.40	0.02	0.33-0.47	5.50	18	2.45	0.17	1.95-2.95	6.82
PG: Siemens Dimension/AR/ES/RxL/Xpand/EXL	13	2.73	0.11	2.40-3.06	3.99	12	4.32	0.37	3.21-5.43	8.54	12	6.63	0.77	4.31-8.95	11.66	13	0.41	0.04	0.29-0.53	9.51	13	2.48	0.21	1.84-3.12	8.55
PG: Qualigen FastPak Thyroid	13	2.28	0.20	1.68-2.88	8.82	13	4.44	0.28	3.60-5.28	6.33	12	6.43	0.71	4.30-8.56	11.06	12	0.27	0.04	0.15-0.39	14.81	12	2.43	0.17	1.91-2.95	7.08
PG: TOSOH AIA	16	3.01	0.27	2.20-3.82	9.00	15	4.97	0.52	3.40-6.54	10.50	16	7.59	0.70	5.48-9.70	9.25	16	0.49	0.05	0.35-0.63	9.39	15	2.94	0.26	2.16-3.72	8.81
SG: Select Beckman Coulter Instruments (Chem)	24	2.70	0.20	2.12-3.29	7.22	24	4.31	0.36	3.23-5.39	8.35	24	6.59	0.76	4.31-8.87	11.52	24	0.39	0.02	0.32-0.46	5.90	24	2.44	0.16	1.95-2.93	6.72
SG: Select Siemens Instruments	17	2.76	0.13	2.36-3.16	4.82	16	4.36	0.34	3.33-5.39	7.89	16	6.69	0.67	4.67-8.71	10.04	17	0.42	0.04	0.30-0.54	9.52	17	2.51	0.19	1.93-3.09	7.73
PG: Roche Cobas e Series	8	2.78	0.13	2.40-3.16	4.60	8	4.30	0.16	3.83-4.77	3.65	8	5.99	0.22	5.33-6.65	3.69	8	0.54	0.02	0.49-0.59	2.96	8	2.63	0.08	2.38-2.88	3.16
SG: Select Abbott Architect Instruments	13	2.44	0.17	1.94-2.94	6.89	13	4.00	0.28	3.17-4.83	6.88	13	6.28	0.35	5.22-7.34	5.64	13	0.36	0.03	0.26-0.46	9.17	13	2.21	0.15	1.76-2.66	6.83
PG: NanoEnTek Frend System	10	2.57	0.35	1.51-3.63	13.77	10	1.98	0.34	0.96-3.00	17.17	10	3.53	0.60	1.74-5.32	16.94	10	0.22	0.05	0.07-0.37	23.18	10	1.19	0.21	0.57-1.81	17.39
PG: Intl OCD Vitros ECI/ECIQ	5	2.68	0.21	2.06-3.30	7.76	5	5.88	0.34	4.86-6.90	5.80	5	8.72	0.48	7.27-10.17	5.55	5	0.40	0.03	0.33-0.48	6.25	5	3.35	0.07	3.15-3.55	1.97
AG: All Group	97	2.65	0.26	1.87-3.43	9.85	91	4.39	0.50	2.90-5.88	11.28	91	6.61	0.81	4.18-9.04	12.24	96	0.40	0.09	0.14-0.66	21.50	91	2.49	0.31	1.56-3.42	12.49
Urea Nitrogen (BUN)	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%
PG: Abaxis Piccolo (Non-Waived)	9	5.67	0.50	4.00-8.00	8.33	9	25.67	0.71	24.00-28.00	2.69	9	40.00	0.50	36.00-44.00	1.25	9	5.44	0.53	3.00-7.00	10.00	9	15.11	0.60	13.00-17.00	4.00
PG: Abaxis Piccolo (Waived)	32	5.90	0.62	4.00-8.00	10.00	31	25.62	0.51	24.00-28.00	1.92															
PG: Abbott Architect	6	6.17	0.41	4.00-8.00	6.67	6	28.50	1.05	26.00-32.00	3.45	6	44.50	0.84	41.00-49.00	1.78	6	6.00	0.00	4.00-8.00	0.00	6	17.28	0.45	15.00-19.00	2.35
PG: Alfa Wasserman Alera/Axcel	24	5.75	0.53	4.00-8.00	8.33	24	27.29	0.95	25.00-29.00	3.70	24	41.67	1.37	38.00-46.00	3.33	24	6.17	0.38	4.00-8.00	6.67	24	17.08	0.78	15.00-19.00	4.71
PG: Siemens Dimension/AR/ES/RxL/Xpand/EXL	15	6.47	0.64	4.00-8.00	10.00	15	28.60	1.06	26.00-32.00	3.79	15	43.87	1.69	40.00-48.00	3.86	15	6.20	0.68	4.00-8.00	11.67	15	17.40	0.83	15.00-19.00	4.71
PG: Beckman AU	17	6.47	0.87	4.00-8.00	15.00	17	28.24	1.30	25.00-31.00	4.64	17	42.53	1.88	39.00-47.00	4.42	16	6.06	0.57	4.00-8.00	10.00	16	17.25	1.00	15.00-19.00	5.88
PG: OCD Vitros 250,350,400,500,750,950	8	4.88	0.35	3.00-7.00	8.00	8	23.13	0.35	21.00-25.00	1.74	8	36.13	0.35	33.00-39.00	1.11	8	4.75	0.46	3.00-7.00	10.00	8	14.00	0.00	12.00-16.00	0.00
SG: Select Beckman Coulter Instruments (Chem)	21	6.57	0.81	5.00-9.00	11.43	21	28.38	1.24	25.00-31.00	4.29	21	43.00	1.98	39.00-47.00	4.65	20	6.10	0.55	4.00-8.00	10.00	20	17.30	0.92	15.00-19.00	5.29
SG: Select Roche Instruments	7	5.86	0.38	4.00-8.00	6.67	7	27.43	0.79	25.00-29.00	2.96	7	42.14	1.35	38.00-46.00	3.10	7	6.00	0.00	4.00-8.00	0.00	7	17.00	0.82	15.00-19.00	4.71
SG: Select OCD Vitros Instruments	10	4.80	0.42	3.00-7.00	8.00	10	22.80	0.79	21.00-25.00	3.48	10	35.80	0.79	33.00-39.00	2.22	10	4.70	0.48	3.00-7.00	10.00	10	13.80	0.42	12.00-16.00	2.86
SG: Abaxis Piccolo Instruments	41	5.85	0.60	4.00-8.00	10.00	40	25.63	0.55	24.00-28.00	1.92	9	40.00	0.50	36.00-44.00	1.25	9	5.44	0.53	3.00-7.00	10.00	9	15.11	0.60	13.00-17.00	4.00
SG: Select Abbott Architect Instruments	9	6.11	0.33	4.00-8.00	5.00	9	28.22	0.97	25.00-31.00	3.57	9	43.89	1.17	40.00-48.00	2.73	9	6.00	0.00	4.00-8.00	0.00	9	17.19	0.38	15.00-19.00	2.35
SG: International OCD Vitros Instruments	13	4.77	0.44	3.00-7.00	8.00	13	22.85	0.80	21.00-25.00	3.48	13	35.77	0.93	33.00-39.00	2.50	12	4.75	0.45	3.00-7.00	10.00	12	13.92	0.52	12.00-16.00	3.57
SG: International Abaxis Piccolo Instruments	45	5.82	0.58	4.00-8.00	10.00	44	25.67	0.57	24.00-28.00	2.31	13	39.94	0.55	36.00-44.00	1.25	13	5.52	0.47	4.00-8.00	8.33	13	15.15	0.55	13.00-17.00	4.00
SG: Select Gen Chem Instruments	60	5.97	0.86	4.00-8.00	15.00	60	27.77	1.20	25.00-31.00	4.29	60	42.47	1.75	38.00-46.00	4.29	59	5.99	0.56	4.00-8.00	10.00	59	17.07	0.88	15.00-19.00	5.29
Albumin	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%
PG: Abaxis Piccolo (Non-Waived) / BCP	7	5.96	0.20	5.40-6.60	3.33	7	3.49	0.11	3.10-3.90	3.14	7	4.31	0.17	3.90-4.70	3.95	7	1.91	0.09	1.70-2.10	4.74	7	2.74	0.11	2.40-3.00	4.07
PG: Abaxis Piccolo (Waived) / BCP	34	5.92	0.12	5.30-6.50	2.03	34	3.51	0.42	3.10-3.90	12.00															
PG: Abbott Architect / BCP	6	6.22	0.35	5.60-6.80	5.65	6	3.60	0.11	3.20-4.00	3.06	6	4.68	0.17	4.20-5.20	3.62	6	1.87	0.08	1.70-2.10	4.21	6	2.77	0.10	2.50-3.10	3.57
PG: Alfa Wasserman Alera/Axcel / BCG	24	6.78	0.16	6.10-7.50	2.35	24	4.26	0.09	3.90-4.70	2.09	24	5.58	0.12	5.00-6.20	2.14	24	2.09	0.03	1.90-2.30	1.43	24	3.20	0.10	2.90-3.50	3.13
PG: Clinical Data Envoy 500 / BCG	5	6.66	0.20	6.00-7.40	2.84	5	4.32	0.19	3.90-4.70	4.42	5	5.60	0.14	5.00-6.20	2.50	5	2.16	0.06	2.00-2.40	2.27	5	3.30	0.14	3.00-3.60	4.24
PG: Siemens Dimension/AR/ES/RxL/Xpand/EXL / BCG	10	6.65	0.14	6.00-7.40	2.09	10	3.61	0.09	3.20-4.00	2.50	10	4.83	0.15	4.30-5.30	3.13	10	1.74	0.07	1.50-1.90	4.12	10	2.68	0.06	2.40-3.00	2.22
PG: Siemens Dimension/AR/ES/RxL/Xpand/EXL / BCP	5	6.62	0.24	5.90-7.30	3.64	5	3.62	0.13	3.20-4.00	3.61	5	4.84	0.15	4.30-5.30	3.13	5	1.76	0.06	1.60-2.00	2.78	5	2.70	0.07	2.40-3.00	2.59
PG: Beckman AU / BCG	17	6.58	0.52	5.90-7.30	7.88	17	4.18	0.10	3.80-4.60	2.38	17	5.63	0.33	5.00-6.20	5.89	16	1.91	0.06	1.70-2.10	3.16	17	3.02	0.27	2.70-3.30	9.00
PG: Roche Cobas Integra / BCG	6	6.32	0.33	5.70-6.90	5.24	6	4.25	0.14	3.90-4.70	3.26	6	5.57	0.21	5.00-6.20	3.75	6	1.90	0.06	1.70-2.10	3.16	6	3.15	0.08	2.90-3.50	2.50
PG: OCD Vitros 250,350,400,500,750,950 / BCG	8	6.39	0.19	5.80-7.00	2.97	8	3.94	0.12	3.50-4.30	3.08	8	5.41	0.16	4.90-5.90	2.96	8	1.61	0.04	1.40-1.80	2.50	8	2.78	0.09	2.50-3.10	3.21
SG: Select Beckman Coulter Instruments (Chem) / BCG	16	6.63	0.16	5.90-7.30	2.42	18	4.19	0.10	3.80-4.60	2.38	18	5.53	0.18	5.00-6.00	3.27	18	1.91	0.06	1.70-2.10	3.16	18	3.08	0.10	2.80-3.40	3.23
SG: Select Roche Instruments / BCG	7	6.36	0.32	5.80-7.00	5.00	7	4.26	0.13	3.90-4.70	3.02	7	5.59	0.20	5.00-6.20	3.57	7	1.91	0.07	1.70-2.10	3.68	7	3.16	0.08	2.90-3.50	2.50
SG: Select OCD Vitros Instruments / BCG	9	6.34	0.22	5.70-6.90	3.49	10	3.95	0.11	3.60-4.40	2.75	10	5.47	0.20	5.00-6.00	3.45	10	1.60	0.05	1.40-1.80	3.13	10	2.76	0.08	2.50-3.10	2.86
SG: Abaxis Piccolo Instruments / BCP	41	5.93	0.13	5.30-6.50	2.20	41	3.50	0.39	3.10-3.90	10.86	7	4.31	0.17	3.90-4.70	3.95	7	1.91	0.09	1.70-2.10	4.74	7	2.74	0.11	2.40-3.00	4.07
SG: Select Abbott Architect Instruments / BCP	9	6.34	0.35	5.70-6.90	5.56	9	3.62	0.10	3.20-4.00	2.78	9	4.69	0.14	4.20-5.20	3.19	9	1.87	0.07	1.70-2.10	3.68	9	2.78	0.08	2.50-3.10	2.86
SG: International OCD Vitros Instruments / BCG	12	6.33	0.22	5.70-6.90	3.49	13	3.95	0.10	3.50-4.30	2.56	13	5.45	0.19	5.00-6.00	3.45	12	1.61	0.05	1.40-1.80	3.13	12	2.76	0.08	2.50-3.10	2.86
SG: International Abaxis Piccolo Instruments / BCP	45	5.94	0.14	5.30-6.50	2.37	45	3.50	0.37	3.10-3.90	10.57	11	4.33	0.14	3.90-4.70	3.26	11	1.93	0.08	1.70-2.10	4.21	11	2.75	0.10	2.50-3.10	3.57
SG: Select Gen Chem Instruments / BCG	46	6.72	0.19	6.00-7.40	2.84	48																			



MG: Method Group	123	100.00	11.30	70.00-130.00	11.30	123	184.00	23.50	129.00-239.00	12.77	88	285.00	42.70	200.00-371.00	14.98	86	33.00	4.40	23.00-43.00	13.33	87	115.00	12.10	81.00-150.00	10.52
MG: IFCC	7	110.00	6.70	77.00-143.00	6.09	7	213.00	9.00	149.00-277.00	4.23	7	328.00	15.30	230.00-426.00	4.66	7	34.00	2.10	24.00-44.00	6.18	7	127.00	6.00	89.00-165.00	4.72
ALT/SGPT	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%
PG: Abaxis Piccolo (Non-Waived)	6	28.50	3.27	23.00-35.00	11.38	6	141.67	4.27	114.00-170.00	3.03	6	220.17	3.25	176.00-264.00	1.50	6	23.00	3.03	18.00-28.00	13.04	6	83.17	2.93	66.00-100.00	3.49
PG: Abaxis Piccolo (Waived)	35	28.29	2.31	22.00-34.00	8.21	35	137.29	19.38	110.00-164.00	14.16															
PG: Abbott Architect	6	22.50	1.64	18.00-28.00	6.96	6	161.17	3.31	129.00-193.00	2.05	6	258.50	5.72	207.00-311.00	2.20	6	16.83	1.94	14.00-20.00	11.18	6	89.33	3.39	71.00-107.00	3.82
PG: Alfa Wasserman Alera/Axcel	24	7.75	3.26	6.00-10.00	41.25	24	124.33	17.58	99.00-149.00	14.19	24	204.21	6.19	163.00-245.00	3.04	24	12.79	1.59	10.00-16.00	12.31	24	68.67	2.85	55.00-83.00	4.20
PG: Beckman Coulter Dxc Series	5	34.00	1.00	27.00-41.00	2.94	5	159.00	2.65	127.00-191.00	1.64	5	248.00	5.20	198.00-298.00	2.10	5	24.00	1.00	19.00-29.00	4.17	5	92.33	2.31	74.00-110.00	2.50
PG: Clinical Data Envoy 500	5	23.82	3.91	19.00-29.00	16.25	5	156.66	2.39	126.00-188.00	1.53	5	242.72	3.25	194.00-292.00	1.36	5	17.52	0.48	14.00-22.00	2.78	5	102.80	34.23	82.00-124.00	33.20
PG: Siemens Dimension/AR/ES/RxL/Xpand/EXL	6	29.67	3.88	24.00-36.00	13.00	6	162.33	3.67	130.00-194.00	2.28	6	255.50	7.01	205.00-307.00	2.73	6	20.00	2.28	16.00-24.00	11.50	6	92.50	3.39	74.00-112.00	3.66
PG: Siemens Dimension/AR/ES/RxL/Xpand/EXL / IFCC	9	28.11	2.80	22.00-34.00	10.00	9	164.00	3.12	131.00-197.00	1.89	9	256.00	3.43	205.00-307.00	1.33	9	20.22	1.48	16.00-24.00	7.50	9	92.22	2.22	74.00-110.00	2.39
PG: Beckman AU	17	21.35	1.37	17.00-25.00	6.67	17	139.76	4.42	112.00-168.00	3.14	17	225.53	10.05	181.00-271.00	4.47	16	15.06	0.93	12.00-18.00	6.00	16	77.56	2.66	62.00-94.00	3.46
PG: Roche Cobas Integra	6	20.50	1.05	17.00-25.00	4.76	6	150.83	4.88	121.00-181.00	3.25	6	238.50	7.56	191.00-287.00	3.18	6	16.50	0.55	14.00-20.00	2.94	6	85.00	2.28	68.00-102.00	2.71
PG: OCD Vitros 250,350,400,500,750,950	6	30.17	3.97	24.00-36.00	13.33	7	159.57	5.29	128.00-192.00	3.31	7	237.71	7.11	190.00-286.00	2.98	7	37.43	2.88	30.00-44.00	7.84	7	100.00	4.24	80.00-120.00	4.20
SG: Select Roche Instruments	7	20.57	0.98	17.00-25.00	4.76	7	151.14	4.53	121.00-181.00	2.98	7	239.00	7.02	191.00-287.00	2.93	7	16.57	0.54	14.00-20.00	2.94	7	85.14	2.12	68.00-102.00	2.47
SG: Select OCD Vitros Instruments	8	29.63	4.75	24.00-36.00	15.67	9	157.00	8.20	126.00-188.00	5.22	9	236.78	7.45	190.00-284.00	3.12	8	37.50	2.67	30.00-46.00	7.11	9	97.44	8.85	78.00-116.00	9.07
SG: Abaxis Piccolo Instruments	41	28.32	2.42	22.00-34.00	8.57	41	137.93	18.00	110.00-166.00	13.04	6	220.17	3.25	176.00-264.00	1.50	6	23.00	3.03	18.00-28.00	13.04	6	83.17	2.93	66.00-100.00	3.49
SG: Select Abbott Architect Instruments	9	22.78	1.86	18.00-28.00	8.26	9	161.11	3.18	129.00-193.00	1.99	9	257.56	5.48	206.00-310.00	2.13	9	16.67	2.18	14.00-20.00	12.94	9	89.22	3.07	71.00-107.00	3.48
SG: International OCD Vitros Instruments	11	27.73	5.22	22.00-34.00	18.57	12	158.33	7.60	126.00-190.00	4.81	12	236.92	6.42	190.00-284.00	2.70	10	39.90	5.63	32.00-48.00	14.00	11	98.91	8.70	79.00-119.00	8.79
SG: International Abaxis Piccolo Instruments	42	28.33	2.40	22.00-34.00	8.57	42	138.10	17.81	110.00-166.00	12.90	7	221.00	3.70	177.00-265.00	1.67	7	23.43	2.99	18.00-28.00	13.04	7	84.00	3.46	67.00-101.00	4.17
SG: Select Gen Chem Instruments	33	21.76	2.69	18.00-26.00	12.27	58	139.14	18.67	111.00-167.00	13.45	58	224.32	21.31	179.00-269.00	9.51	54	14.55	2.27	12.00-18.00	15.33	56	77.13	9.09	62.00-92.00	11.82
MG: IFCC	8	26.00	2.40	21.00-31.00	9.23	11	162.00	8.50	130.00-194.00	5.25	11	255.00	12.50	204.00-306.00	4.90	11	20.00	2.10	16.00-24.00	10.50	11	91.00	4.90	73.00-109.00	5.38
AST/SGOT	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%
PG: Abaxis Piccolo (Non-Waived)	6	41.67	1.63	34.00-50.00	3.81	6	210.50	4.85	169.00-253.00	2.27	6	276.67	6.86	222.00-332.00	2.49	6	106.83	2.79	86.00-128.00	2.62	6	157.17	4.96	126.00-188.00	3.18
PG: Abaxis Piccolo (Waived)	35	40.37	1.91	32.00-48.00	4.75	35	201.57	28.72	162.00-242.00	14.21															
PG: Abbott Architect	6	37.00	0.63	30.00-44.00	1.62	6	211.67	3.88	170.00-254.00	1.84	6	282.33	5.28	226.00-338.00	1.88	6	104.17	1.17	83.00-125.00	1.15	6	156.67	3.08	126.00-188.00	1.97
PG: Alfa Wasserman Alera/Axcel	24	25.96	2.24	21.00-31.00	8.46	24	193.38	6.30	154.00-232.00	3.26	24	255.92	7.47	205.00-307.00	2.93	24	97.92	3.05	78.00-118.00	3.06	24	143.38	5.69	114.00-172.00	3.99
PG: Clinical Data Envoy 500	5	34.90	3.58	28.00-42.00	10.29	5	211.74	6.79	170.00-254.00	3.21	5	278.14	10.34	222.00-334.00	3.71	5	105.14	5.45	84.00-126.00	5.24	5	158.58	6.39	127.00-191.00	4.03
PG: Siemens Dimension/AR/ES/RxL/Xpand/EXL	15	32.40	2.29	26.00-38.00	7.19	15	218.87	6.73	175.00-263.00	3.06	15	293.00	8.60	234.00-352.00	2.94	15	106.13	5.00	85.00-127.00	4.72	14	162.57	5.91	130.00-196.00	3.62
PG: Beckman AU	17	33.00	1.62	26.00-40.00	4.85	17	187.76	7.37	150.00-226.00	3.94	17	249.65	8.93	200.00-300.00	3.56	17	96.35	13.31	77.00-115.00	13.85	17	136.41	10.74	109.00-163.00	7.87
PG: Roche Cobas Integra	6	33.50	0.55	27.00-41.00	1.47	6	217.00	7.29	174.00-260.00	3.36	6	285.17	10.74	228.00-342.00	3.75	6	108.33	4.03	86.00-130.00	3.70	6	162.00	5.76	130.00-194.00	3.58
PG: OCD Vitros 250,350,400,500,750,950	7	39.57	2.23	32.00-48.00	5.50	7	228.71	9.57	183.00-275.00	4.19	7	341.00	16.54	273.00-409.00	4.84	7	101.57	5.22	82.00-122.00	5.10	7	158.71	6.55	127.00-191.00	4.15
SG: Select Roche Instruments	7	33.57	0.54	27.00-41.00	1.47	7	216.43	6.83	173.00-259.00	3.15	7	284.86	9.84	228.00-342.00	3.44	7	108.14	3.72	86.00-130.00	3.43	7	161.57	5.38	130.00-194.00	3.33
SG: Select OCD Vitros Instruments	9	39.11	2.15	31.00-47.00	5.38	9	227.44	9.34	182.00-272.00	4.10	9	336.44	18.23	269.00-403.00	5.42	9	100.56	5.05	81.00-121.00	5.05	9	157.44	6.52	126.00-188.00	4.14
SG: Abaxis Piccolo Instruments	41	40.56	1.91	33.00-49.00	4.63	41	202.88	26.73	162.00-244.00	13.15	6	276.67	6.86	222.00-332.00	2.49	6	106.83	2.79	86.00-128.00	2.62	6	157.17	4.96	126.00-188.00	3.18
SG: Select Abbott Architect Instruments	9	37.00	0.71	30.00-44.00	1.89	9	210.11	7.18	168.00-252.00	3.43	9	282.00	5.03	226.00-338.00	1.77	9	104.44	1.67	83.00-125.00	1.63	9	156.67	2.92	126.00-188.00	1.85
SG: International OCD Vitros Instruments	12	39.00	1.95	31.00-47.00	5.13	12	227.17	9.17	182.00-272.00	4.05	12	337.75	18.00	270.00-406.00	5.33	11	100.27	4.76	80.00-120.00	4.80	11	157.18	6.45	126.00-188.00	4.08
SG: International Abaxis Piccolo Instruments	42	40.57	1.89	33.00-49.00	4.63	42	203.02	26.42	162.00-244.00	13.00	7	276.29	6.34	221.00-331.00	2.28	7	106.43	2.76	85.00-127.00	2.64	7	157.14	4.53	126.00-188.00	2.87
SG: Select Gen Chem Instruments	57	31.13	5.19	25.00-37.00	16.77	58	196.77	11.33	158.00-236.00	5.74	58	260.96	14.51	209.00-313.00	5.56	57	98.77	5.70	79.00-119.00	5.76	57	146.10	8.89	117.00-175.00	6.10
AG: All Group	131	34.99	5.67	28.00-42.00	16.29	132	204.46	19.77	163.00-245.00	9.71	97	275.87	26.88	221.00-331.00	9.75	96	101.33	6.52	81.00-121.00	6.44	95	151.57	10.70	122.00-182.00	7.04
Bilirubin Total	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%
PG: Abaxis Piccolo (Non-Waived)	5	0.53	0.10	0.10-0.90	20.00	5	2.70	0.08	2.20-3.20	2.96	5	3.97	0.12	3.20-4.80	3.00	5	0.53	0.05	0.10-0.90	10.00	5	1.55	0.13	1.20-2.00	8.13
PG: Abaxis Piccolo (Waived)	35	0.52	0.07	0.10-0.90	14.00	35	2.58	0.41	2.10-3.10	15.77															
PG: Abbott Architect	6	0.48	0.04	0.10-0.90	8.00	6	2.67	0.17	2.20-3.20	6.67	6	4.07	0.20	3.30-4.90	4.88	6	0.38	0.04	0.00-0.80	10.00	6	1.53	0.14	1.10-1.90	9.33
PG: Alfa Wasserman Alera/Axcel	24	0.65	0.06	0.30-1.10	8.57	24	3.46	0.17	2.80-4.20	4.86	24	5.23	0.23	4.20-6.20	4.42	24	0.51	0.06	0.10-0.90	12.00	24	1.92	0.11	1.50-2.30	5.79
PG: Clinical Data Envoy 500	5	0.40	0.07	0.00-0.80	17.50	5	2.62	0.32	2.10-3.10	12.31	5	4.00	0.49	3.20-4.80	12.25	5	0.34	0.09	0.00-0.70	30.00	5	1.38	0.18	1.	



SG: Abaxis Piccolo Instruments	38	11.01	0.25	9.90-12.10	2.27	37	6.89	0.10	6.20-7.60	1.45															
SG: Select Abbott Architect Instruments	9	11.23	0.29	10.10-12.30	2.59	9	6.98	0.22	6.30-7.70	3.14	9	9.87	0.23	8.90-10.90	2.32	9	2.62	0.12	2.30-2.90	4.62	9	4.77	0.16	4.30-5.30	3.33
SG: International OCD Vitros Instruments	11	11.35	0.44	10.20-12.40	3.89	12	6.63	0.14	5.90-7.30	2.12	12	9.22	0.20	8.30-10.10	2.17	11	2.75	0.10	2.40-3.00	3.70	11	4.76	0.12	4.30-5.30	2.50
SG: International Abaxis Piccolo Instruments	39	11.00	0.25	9.90-12.10	2.27	38	6.89	0.10	6.20-7.60	1.45															
SG: Select Gen Chem Instruments	60	10.81	0.46	9.70-11.90	4.26	60	6.87	0.23	6.20-7.60	3.33	60	9.58	0.34	8.60-10.60	3.54	58	2.71	0.11	2.40-3.00	4.07	59	4.80	0.16	4.30-5.30	3.33
AG: All Group	129	10.94	0.44	9.80-12.00	4.04	129	6.88	0.21	6.20-7.60	3.04	94	9.58	0.36	8.60-10.60	3.75	93	2.73	0.11	2.40-3.00	4.07	93	4.82	0.17	4.30-5.30	3.54
Phosphorous	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%
PG: Alfa Wasserman Alera/Axcel	6	1.38	0.17	0.90-1.90	12.14	5	3.36	0.09	3.10-3.70	2.65	6	4.52	0.10	4.20-4.80	2.22	6	1.57	0.08	1.40-1.80	5.00	6	2.42	0.07	2.20-2.60	3.33
PG: Beckman AU	8	1.53	0.07	1.30-1.70	4.67	9	3.63	0.49	2.10-5.10	13.61	9	4.98	0.53	3.40-6.60	10.60	7	1.43	0.05	1.30-1.50	3.57	8	2.35	0.36	1.30-3.50	15.00
SG: Select Beckman Coulter Instruments (Chem)	11	1.56	0.13	1.20-2.00	8.13	12	3.63	0.43	2.30-4.90	11.94	12	4.95	0.45	3.60-6.30	9.00	10	1.44	0.05	1.30-1.50	3.57	11	2.39	0.31	1.50-3.30	12.92
PG: Intl OCD Vitros 250,350,400,500,750,950	5	2.10	0.17	1.60-2.60	8.10	5	3.77	0.06	3.60-4.00	1.58	5	5.27	0.06	5.10-5.50	1.13	5	1.85	0.07	1.70-2.10	3.68	5	2.60	0.52	1.00-4.20	20.00
PG: Intl Abaxis Piccolo	5	2.22	0.17	1.70-2.70	7.73	5	3.94	0.14	3.50-4.30	3.59	5	5.41	0.11	5.10-5.70	2.04	5	1.97	0.16	1.50-2.50	8.00	5	3.04	0.09	2.70-3.30	3.00
SG: Select Gen Chem Instruments	23	1.55	0.20	0.90-2.10	13.33	23	3.53	0.34	2.50-4.50	9.71	24	4.83	0.40	3.60-6.00	8.33	22	1.48	0.09	1.20-1.80	6.00	23	2.43	0.24	1.70-3.10	10.00
AG: All Group	37	1.64	0.26	0.80-2.40	16.25	37	3.56	0.31	2.70-4.50	8.61	38	4.89	0.37	3.80-6.00	7.55	36	1.54	0.25	0.80-2.30	16.67	37	2.51	0.26	1.70-3.30	10.40
Cholesterol Total	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%
PG: Alfa Wasserman Alera/Axcel	21	259.05	4.65	233.00-285.00	1.81	21	154.48	2.25	139.00-169.00	1.43	21	213.33	3.67	192.00-234.00	1.74	21	62.24	1.87	56.00-68.00	3.06	21	109.33	2.58	98.00-120.00	2.39
PG: Siemens Dimension/AR/ES/RxL/Xpand/EXL	14	255.21	6.24	230.00-281.00	2.43	13	147.08	3.55	132.00-162.00	2.38	14	207.43	5.06	186.00-228.00	2.46	14	55.71	2.40	50.00-62.00	4.29	14	101.36	2.53	91.00-111.00	2.48
PG: Beckman AU	17	254.29	6.91	229.00-279.00	2.72	17	149.53	3.86	135.00-165.00	2.60	17	210.71	5.59	190.00-232.00	2.65	16	57.75	1.98	52.00-64.00	3.45	16	105.00	6.73	95.00-116.00	6.38
PG: OCD Vitros 250,350,400,500,750,950	7	279.14	4.78	251.00-307.00	1.72	7	150.86	4.74	136.00-166.00	3.11	7	216.00	2.58	194.00-238.00	1.20	7	47.14	3.93	42.00-52.00	8.30	7	101.71	2.63	92.00-112.00	2.55
SG: Select Beckman Coulter Instruments (Chem)	21	255.10	6.47	230.00-281.00	2.55	21	150.67	4.33	136.00-166.00	2.85	21	212.43	6.52	191.00-233.00	3.07	20	58.20	2.02	52.00-64.00	3.45	20	105.55	6.11	95.00-117.00	5.75
SG: Select Roche Instruments	7	249.29	5.94	224.00-274.00	2.37	7	149.29	3.77	134.00-164.00	2.55	7	208.57	5.56	188.00-230.00	2.68	7	58.57	1.99	53.00-65.00	3.39	7	104.14	2.41	94.00-114.00	2.31
SG: Select OCD Vitros Instruments	9	279.56	6.17	252.00-308.00	2.21	9	150.78	4.47	136.00-166.00	2.98	9	215.33	3.00	194.00-237.00	1.40	8	47.50	3.78	43.00-53.00	7.92	9	101.00	2.96	91.00-111.00	2.97
SG: Abaxis Piccolo Instruments	7	253.29	4.23	228.00-278.00	1.66	7	152.29	1.80	137.00-167.00	1.18															
PG: Polymer CardioChek	53	246.89	17.22	222.00-272.00	6.96	59	157.63	10.15	142.00-174.00	6.39															
SG: Select Abbott Architect Instruments	6	261.67	3.50	236.00-288.00	1.34	6	153.33	1.51	138.00-168.00	0.98	6	215.33	2.42	194.00-237.00	1.12	6	61.17	0.75	55.00-67.00	1.31	6	107.83	1.72	97.00-119.00	1.57
SG: International OCD Vitros Instruments	12	277.17	9.75	249.00-305.00	3.50	12	151.50	4.89	137.00-167.00	3.22	12	215.58	4.46	194.00-238.00	2.08	10	47.90	4.04	43.00-53.00	8.33	11	102.18	4.04	92.00-112.00	3.92
SG: International Abaxis Piccolo Instruments	8	253.25	3.92	228.00-278.00	1.54	8	152.63	1.92	138.00-168.00	1.24															
SG: Select Gen Chem Instruments	53	257.40	6.19	231.00-283.00	2.41	53	152.94	3.79	138.00-168.00	2.48	53	213.69	5.39	193.00-235.00	2.52	52	60.43	2.78	54.00-66.00	4.67	52	107.78	4.64	97.00-119.00	4.26
SG: Select Cholestech Instruments	5	286.60	9.86	258.00-316.00	3.45	5	162.40	6.31	146.00-178.00	3.89	5	226.80	11.78	204.00-250.00	5.20	5	100.00	0.00	90.00-110.00	0.00	5	110.60	2.97	100.00-122.00	2.70
Cholesterol HDL	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%
PG: Alfa Wasserman Alera/Axcel / Direct Method	21	84.48	3.72	59.00-109.00	4.40	21	65.62	3.07	46.00-86.00	4.70	21	91.24	4.39	64.00-118.00	4.84	21	26.62	0.86	19.00-35.00	3.33	21	45.81	2.11	32.00-60.00	4.57
PG: Siemens Dimension/AR/ES/RxL/Xpand/EXL / Direct Method	14	56.64	1.69	40.00-74.00	2.98	14	77.50	2.68	55.00-101.00	3.46	13	109.62	3.55	77.00-143.00	3.18	14	28.43	1.09	20.00-36.00	3.93	14	53.57	1.70	38.00-70.00	3.15
PG: Beckman AU / Direct Method	17	87.24	3.60	61.00-113.00	4.14	17	68.06	4.37	48.00-88.00	6.47	17	99.14	7.29	69.00-129.00	7.37	16	24.13	1.86	17.00-31.00	7.92	17	44.21	6.16	31.00-57.00	14.09
PG: Roche Cobas Integra / Direct Method	6	64.00	4.00	45.00-83.00	6.25	6	79.67	3.33	56.00-104.00	4.13	6	114.17	4.62	80.00-148.00	4.04	6	27.17	1.47	19.00-35.00	5.56	6	52.83	1.72	37.00-69.00	3.21
PG: OCD Vitros 250,350,400,500,750,950 / Direct Method	7	87.43	8.30	61.00-113.00	9.54	7	78.14	3.13	55.00-101.00	3.97	7	118.14	9.96	83.00-153.00	8.47	7	25.86	1.07	18.00-34.00	4.23	6	49.17	1.72	34.00-64.00	3.47
SG: Select Beckman Coulter Instruments (Chem) / Direct Method	21	87.10	3.51	61.00-113.00	4.02	21	68.05	4.26	48.00-88.00	6.32	21	99.11	6.77	69.00-129.00	6.87	20	24.00	1.75	17.00-31.00	7.08	21	44.60	5.66	32.00-59.00	12.67
SG: Select Roche Instruments / Direct Method	7	64.43	3.82	45.00-83.00	5.94	7	80.29	3.45	56.00-104.00	4.38	7	114.86	4.60	81.00-150.00	4.00	7	27.29	1.38	19.00-35.00	5.19	7	53.14	1.77	37.00-69.00	3.40
SG: Select OCD Vitros Instruments / Direct Method	9	86.11	7.69	60.00-112.00	8.95	9	76.89	3.76	54.00-100.00	4.94	8	116.63	10.17	82.00-152.00	8.72	9	25.67	1.12	18.00-34.00	4.23	8	48.63	2.07	34.00-64.00	4.29
SG: Abaxis Piccolo Instruments / Direct Method	7	50.86	6.31	36.00-66.00	12.35	7	64.29	5.25	45.00-83.00	8.28															
PG: Polymer CardioChek / Direct Method	67	107.72	9.23	76.00-140.00	8.52	68	101.53	13.55	71.00-133.00	13.33															
SG: Select Abbott Architect Instruments / Direct Method	6	85.50	2.51	60.00-112.00	2.91	6	68.00	3.29	48.00-88.00	4.85	6	97.00	4.38	68.00-126.00	4.54	5	25.00	2.12	18.00-33.00	8.40	6	46.00	2.00	32.00-60.00	4.35
SG: International OCD Vitros Instruments / Direct Method	12	86.25	7.25	60.00-112.00	8.49	12	76.58	3.31	54.00-100.00	4.29	11	116.36	9.75	81.00-151.00	8.36	11	25.64	1.03	18.00-34.00	3.85	10	48.60	1.96	34.00-64.00	4.08
SG: Select Gen Chem Instruments / Direct Method	53	85.45	4.32	60.00-111.00	5.06	53	67.13	3.85	47.00-87.00	5.67	53	95.66	6.57	67.00-125.00	6.88	51	25.10	1.98	18.00-33.00	8.00	53	45.41	3.95	32.00-59.00	8.67
MG: Direct Method	172	88.00	35.50	62.00-114.00	40.34	162	84.00	17.80	59.00-109.00	21.19	85	102.00	10.40	71.00-133.00	10.20	85	25.00	3.20	18.00-33.00	12.80	86	47.00	5.10	33.00-61.00	10.85
MG: Direct Method	172	88.00	35.50	62.00-114.00	40.34	162	84.00	17.80	59.00-109.00	21.19	85	102.00	10.40	71.00-133.00	10.20	85	25.00	3.20	18.00-33.00	12.80	86	47.00	5.10	33.00-61.00	10.85
AG: All Group	139	84.91	22.59	60.00-111.00	26.59	152	77.89	26.02	55.00-101.00	33.33	88	103.06	16.63	72.00-134.00	16.12	88	25.70	5.68	18.00-34.00	21.92	89	47.82	9.22	34.00-62.00	19.17
Cholesterol LDL	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%
PG: Alfa Wasserman Alera/Axcel / Calc	12	144.12	6.58	124.00-164.00	4.58	12	57.48	3.67	46.00-68.00	6.49	12	82.98	4.85	69.00-97.00	5.78	12	18.12	3.04	9.00-27.00	16.67	12	39.03	3.47	29.00-50.00	8.97
PG: Siemens Dimension/AR/ES/RxL/Xpand/EXL / Calc	7	148.71	24.38	76.00-222.00	16.38																				



PG: OCD Vitros 250,350,400,500,750,950	7	168.86	6.26	127.00-211.00	3.73	7	177.29	6.40	133.00-221.00	3.62	7	238.29	7.76	179.00-298.00	3.28	7	93.86	2.67	71.00-118.00	2.87	7	135.14	4.56	101.00-169.00	3.41
SG: Select Beckman Coulter Instruments (Chem)	20	150.65	8.26	113.00-189.00	5.50	20	161.90	5.77	122.00-203.00	3.58	20	209.20	6.33	157.00-261.00	3.01	19	88.79	4.71	67.00-111.00	5.28	20	124.30	10.97	93.00-155.00	8.87
SG: Select Roche Instruments	7	146.14	4.98	110.00-183.00	3.42	7	154.43	4.93	116.00-193.00	3.18	7	193.00	5.42	145.00-241.00	2.80	7	90.00	2.52	68.00-113.00	2.78	7	124.00	3.65	93.00-155.00	2.98
SG: Select OCD Vitros Instruments	9	166.00	8.59	125.00-208.00	5.18	9	174.67	8.08	131.00-219.00	4.63	9	234.67	10.44	176.00-294.00	4.43	9	91.67	5.77	69.00-115.00	6.30	9	132.67	7.21	100.00-166.00	5.41
SG: Abaxis Piccolo Instruments	7	160.43	3.31	120.00-200.00	2.06	7	175.00	2.16	131.00-219.00	1.26															
PG: Polymer CardioChek	63	158.98	12.24	119.00-199.00	7.67	61	179.18	15.61	134.00-224.00	8.72															
SG: Select Abbott Architect Instruments	6	144.67	3.14	109.00-181.00	2.14	6	153.00	3.74	115.00-191.00	2.42	6	197.83	4.79	149.00-248.00	2.42	6	84.33	3.50	63.00-105.00	4.17	6	119.33	3.62	89.00-149.00	3.03
SG: International OCD Vitros Instruments	12	167.33	7.71	125.00-209.00	4.61	12	176.17	7.41	132.00-220.00	4.20	12	235.08	10.03	176.00-294.00	4.26	11	91.91	5.19	69.00-115.00	5.65	12	136.75	12.83	103.00-171.00	9.34
SG: International Abaxis Piccolo Instruments	8	160.50	3.07	121.00-201.00	1.93	8	175.00	2.00	131.00-219.00	1.14															
SG: Select Gen Chem Instruments	53	150.62	6.89	113.00-189.00	4.57	53	160.28	6.22	120.00-200.00	3.88	53	204.30	7.75	153.00-255.00	3.82	52	90.52	5.40	68.00-114.00	5.93	53	125.33	7.79	94.00-156.00	6.24
AG: All Group	160	154.65	11.34	116.00-194.00	7.29	158	168.58	14.78	127.00-211.00	8.76	93	207.84	14.08	156.00-260.00	6.78	91	89.00	6.22	67.00-111.00	6.97	93	125.61	9.51	95.00-158.00	7.54
Amylase	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%
PG: Abaxis Piccolo (Waived)	22	75.50	2.52	53.00-99.00	3.29	21	100.57	2.66	71.00-131.00	2.67															
PG: Alfa Wasserman Alera/Axcel	9	85.78	2.33	60.00-112.00	2.67	9	105.22	4.12	74.00-137.00	3.90	9	168.67	5.05	118.00-220.00	2.96	9	13.89	1.45	10.00-18.00	10.71	9	60.56	2.30	43.00-79.00	3.77
PG: Siemens Dimension/AR/ES/RxL/Xpand/EXL	7	85.00	2.08	60.00-111.00	2.47	7	116.43	2.51	81.00-151.00	2.16	7	183.43	3.26	128.00-238.00	1.80	7	13.00	2.08	9.00-17.00	16.15	7	65.00	1.41	46.00-85.00	2.15
PG: Beckman AU	8	63.25	3.37	44.00-82.00	5.40	8	86.38	4.14	60.00-112.00	4.77	8	135.63	6.93	95.00-177.00	5.07	7	10.14	0.38	7.00-13.00	4.00	7	47.71	2.63	34.00-62.00	5.42
PG: OCD Vitros 250,350,400,500,750,950	5	101.00	7.07	71.00-131.00	7.03	5	73.00	8.48	51.00-95.00	11.64	5	121.50	0.71	85.00-159.00	0.57	5	30.00	0.00	21.00-39.00	0.00	5	41.50	4.95	29.00-55.00	11.67
SG: Abaxis Piccolo Instruments	24	75.33	2.51	53.00-98.00	3.33	23	100.22	3.25	70.00-130.00	3.20															
PG: Intl OCD Vitros 250,350,400,500,750,950	5	99.50	6.36	70.00-130.00	6.40	5	72.00	1.41	50.00-94.00	1.94	5	111.50	9.19	78.00-146.00	8.21	5	51.50	30.41	36.00-68.00	58.46	5	71.00	46.67	50.00-92.00	65.77
SG: International Abaxis Piccolo Instruments	25	75.24	2.50	53.00-98.00	3.33	24	99.88	3.59	70.00-130.00	3.60															
SG: Select Gen Chem Instruments	25	79.42	12.72	55.00-103.00	16.08	25	99.50	10.93	69.00-129.00	11.01	25	157.98	18.31	111.00-205.00	11.58	24	13.03	2.46	9.00-17.00	19.23	24	56.75	6.96	40.00-74.00	12.28
AG: All Group	60	79.61	10.10	56.00-104.00	12.63	59	100.96	10.60	71.00-131.00	10.50						37	14.07	4.44	10.00-18.00	31.43	37	57.43	7.66	40.00-74.00	13.51
Creatine Kinase	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%
PG: Abaxis Piccolo (Waived)	11	144.45	6.58	101.00-187.00	4.58	12	223.17	6.85	156.00-290.00	3.05															
PG: Alfa Wasserman Alera/Axcel	7	118.71	4.23	83.00-155.00	3.53	7	194.86	7.43	137.00-254.00	3.79	7	302.71	10.23	212.00-394.00	3.37	7	25.71	1.89	18.00-34.00	7.31	7	110.29	4.99	77.00-143.00	4.55
PG: Siemens Dimension/AR/ES/RxL/Xpand/EXL / IFCC	5	116.00	0.00	81.00-151.00	0.00	5	204.75	3.50	144.00-267.00	1.71	5	322.75	7.14	226.00-420.00	2.20	5	18.25	1.50	13.00-23.00	8.33	5	110.25	2.50	77.00-143.00	2.27
PG: Beckman AU	7	106.86	4.45	75.00-139.00	4.21	7	187.14	7.49	131.00-243.00	4.01	7	296.71	10.03	208.00-386.00	3.37	6	16.17	0.98	11.00-21.00	6.25	6	101.83	3.43	71.00-133.00	3.33
SG: Select Beckman Coulter Instruments (Chem)	9	109.89	7.36	77.00-143.00	6.73	9	190.78	10.59	134.00-248.00	5.55	9	301.78	14.28	211.00-393.00	4.74	8	17.00	1.77	12.00-22.00	10.59	8	104.75	6.59	74.00-137.00	6.29
SG: Abaxis Piccolo Instruments	13	144.00	6.11	101.00-187.00	4.24	14	223.57	6.43	157.00-291.00	2.86															
SG: International Abaxis Piccolo Instruments	13	144.00	6.11	101.00-187.00	4.24	14	223.57	6.43	157.00-291.00	2.86															
PG: Intl Abaxis Piccolo	5	136.67	2.08	96.00-178.00	1.53	5	221.33	0.58	155.00-287.00	0.27	5	344.67	4.51	242.00-449.00	1.30	5	37.67	3.51	27.00-49.00	9.21	5	132.33	3.21	92.00-172.00	2.42
SG: Select Gen Chem Instruments	22	115.77	8.88	81.00-151.00	7.67	22	195.55	14.04	137.00-255.00	7.14	22	307.55	20.80	216.00-400.00	6.75	21	20.33	4.64	14.00-26.00	23.00	21	108.52	7.99	76.00-142.00	7.34
MG: Method Group	43	125.00	15.00	88.00-163.00	12.00	44	205.00	17.50	144.00-267.00	8.54	32	307.00	25.40	215.00-399.00	8.27	29	20.00	4.10	14.00-26.00	20.50	31	111.00	10.20	78.00-144.00	9.19
MG: Method Group	43	125.00	15.00	88.00-163.00	12.00	44	205.00	17.50	144.00-267.00	8.54	32	307.00	25.40	215.00-399.00	8.27	29	20.00	4.10	14.00-26.00	20.50	31	111.00	10.20	78.00-144.00	9.19
GGT	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%
PG: Abaxis Piccolo (Waived)	22	33.64	1.22	30.00-38.00	3.53	22	116.23	3.46	106.00-127.00	3.02															
PG: Beckman AU	7	29.71	1.25	26.00-34.00	4.33	7	102.71	4.03	91.00-115.00	3.88	7	162.57	6.24	144.00-182.00	3.80	6	13.00	0.89	10.00-16.00	6.92	6	58.17	1.72	53.00-63.00	2.93
SG: Select Beckman Coulter Instruments (Chem)	9	29.33	1.32	25.00-33.00	4.48	9	109.00	12.96	70.00-148.00	11.93	9	171.56	18.63	116.00-228.00	10.81	8	14.25	2.49	7.00-22.00	17.86	8	62.50	8.16	38.00-88.00	13.02
SG: Abaxis Piccolo Instruments	24	33.58	1.21	30.00-38.00	3.53	24	116.33	3.41	106.00-126.00	2.93															
PG: Intl OCD Vitros 250,350,400,500,750,950	5	38.00	1.00	35.00-41.00	2.63	5	191.67	4.04	180.00-204.00	2.08	5	286.67	5.13	272.00-302.00	1.78	5	23.50	0.71	22.00-26.00	2.92	5	111.50	3.54	102.00-123.00	3.13
SG: International Abaxis Piccolo Instruments	25	33.56	1.19	30.00-38.00	3.53	25	116.40	3.35	106.00-126.00	2.93															
SG: Select Gen Chem Instruments	15	31.20	3.43	21.00-41.00	10.97	15	108.87	14.90	64.00-154.00	13.67	15	169.53	22.91	101.00-239.00	13.47	14	15.21	2.42	8.00-22.00	16.00	14	63.00	8.46	38.00-89.00	13.49
AG: All Group	44	32.70	6.13	15.00-51.00	18.48	44	116.20	13.89	74.00-158.00	11.98	22	181.36	30.27	90.00-272.00	16.74	21	16.57	4.07	5.00-29.00	24.12	21	67.38	11.64	32.00-102.00	17.31
LDH	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%
PG: Siemens Dimension/AR/ES/RxL/Xpand/EXL / IFCC	5	163.00	0.00	130.00-196.00	0.00	5	460.00	0.00	368.00-552.00	0.00	5	664.00	0.00	531.00-797.00	0.00	5	140.00	0.00	112.00-168.00	0.00	5	302.00	0.00	242.00-362.00	0.00
PG: Intl OCD Vitros 250,350,400,500,750,950	5	428.67	22.37	343.00-515.00	5.22	5	1351.33	49.44	1081.00-1621.00	3.66	5	1905.67	81.59	1525.00-2287.00	4.28	5	459.00	8.48	367.00-551.00	1.85	5	931.50	40.31	746.00-1118.00	4.32
SG: Select Gen Chem Instruments	10	135.80	15.57	109.00-163.00	11.47	10	425.60	48.72	341.00-511.00	11.43	10	609.20	71.33	487.00-731.00	11.71	9	130.00	11.81	104.00-156.00	9.08	9	279.56	27.99	224.00-336.00	10.00
MG: Method Group	11	138.00	16.50	110.00-166.00	11.96	11	431.00	49.60	345.00-5																



MG: Method Group	64	0.20	0.12	0.00-0.55	58.00	65	1.25	0.29	0.37-2.13	23.52	65	1.95	1.04	0.00-5.08	53.49	65	0.25	0.17	0.00-0.75	66.40	65	0.73	0.21	0.10-1.36	28.90
Lactic Acid	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%
PG: Abaxis Piccolo (Non-Waived)	5	0.80	0.00	0.80-0.80	0.00	5	4.00	0.10	3.70-4.30	2.50	5	6.80	0.10	6.50-7.10	1.47	5	0.30	0.00	0.30-0.30	0.00	5	2.23	0.06	2.00-2.40	2.73
AG: All Group	6	0.82	0.04	0.70-0.90	5.00	6	4.07	0.12	3.70-4.50	2.93	6	6.78	0.20	6.20-7.40	3.09	6	0.30	0.00	0.30-0.30	0.00	6	2.27	0.10	2.00-2.60	4.35
Magnesium	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%
PG: Alfa Wasserman Alera/Axcel	8	1.39	0.22	1.00-1.80	15.71	8	3.03	0.21	2.30-3.80	7.00	8	4.31	0.27	3.20-5.40	6.28	8	1.18	0.18	0.90-1.50	15.00	8	2.19	0.19	1.70-2.80	8.64
PG: Siemens Dimension/AR/ES/RxL/Xpand/EXL	6	1.30	0.06	1.00-1.60	4.62	6	3.20	0.17	2.40-4.00	5.31	6	4.62	0.10	3.40-5.80	2.17	6	1.10	0.06	0.80-1.40	5.45	6	2.22	0.07	1.70-2.80	3.64
PG: Beckman AU	8	1.39	0.04	1.00-1.80	2.86	8	3.20	0.08	2.40-4.00	2.50	8	4.51	0.08	3.40-5.60	1.78	7	1.19	0.04	0.90-1.50	3.33	7	2.23	0.05	1.70-2.80	2.27
SG: Select Beckman Coulter Instruments (Chem)	11	1.37	0.05	1.00-1.80	3.57	11	3.13	0.14	2.30-3.90	4.52	11	4.47	0.10	3.40-5.60	2.22	10	1.17	0.05	0.90-1.50	4.17	10	2.16	0.11	1.70-2.80	5.00
PG: Intl Abaxis Piccolo	5	1.39	0.02	1.00-1.80	1.43	5	3.14	0.11	2.30-3.90	3.55	5	4.59	0.18	3.40-5.80	3.91	5	1.24	0.05	0.90-1.50	4.17	5	2.11	0.04	1.60-2.60	1.90
SG: Select Gen Chem Instruments	25	1.46	0.22	1.10-1.90	14.67	26	3.12	0.18	2.30-3.90	5.81	26	4.45	0.23	3.30-5.50	5.23	25	1.15	0.14	0.90-1.50	11.67	25	2.18	0.14	1.70-2.80	6.36
AG: All Group	40	1.40	0.20	1.00-1.80	14.29	41	3.12	0.20	2.30-3.90	6.45	41	4.46	0.25	3.40-5.60	5.56	41	1.16	0.12	0.90-1.50	10.00	40	2.17	0.14	1.70-2.80	6.82
Uric Acid	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%
PG: Abaxis Piccolo (Waived)	6	1.23	0.19	1.00-1.40	15.83	7	5.83	0.17	4.80-6.80	2.93															
PG: Alfa Wasserman Alera/Axcel	12	3.20	0.42	2.70-3.70	13.13	12	7.24	0.37	6.00-8.40	5.14	12	10.38	0.43	8.60-12.20	4.13	11	2.83	0.22	2.30-3.30	7.86	12	5.03	0.23	4.20-5.80	4.60
PG: Siemens Dimension/AR/ES/RxL/Xpand/EXL	11	1.82	0.10	1.50-2.10	5.56	11	6.65	0.19	5.50-7.70	2.88	11	9.86	0.32	8.20-11.60	3.23	11	1.65	0.07	1.40-2.00	4.12	11	4.16	0.12	3.50-4.90	2.86
PG: Beckman AU	11	1.87	0.16	1.60-2.20	8.42	11	6.95	0.12	5.80-8.20	1.71	11	10.35	0.28	8.50-12.10	2.72	10	1.77	0.07	1.50-2.10	3.89	10	4.43	0.08	3.70-5.10	1.82
SG: Select Beckman Coulter Instruments (Chem)	15	1.82	0.18	1.50-2.10	10.00	15	6.89	0.19	5.70-8.10	2.75	15	10.23	0.36	8.50-11.90	3.63	13	1.75	0.07	1.50-2.10	3.89	14	4.39	0.12	3.70-5.10	2.73
SG: Select Roche Instruments	5	1.64	0.06	1.30-1.90	3.13	5	6.82	0.16	5.60-8.00	2.35	5	10.20	0.21	8.50-11.90	2.06	5	1.62	0.04	1.30-1.90	2.50	5	4.30	0.07	3.60-5.00	1.63
SG: Select OCD Vitros Instruments	7	1.84	0.08	1.50-2.10	4.44	7	6.96	0.11	5.80-8.20	1.57	7	10.60	0.08	8.80-12.40	0.75	7	1.56	0.05	1.30-1.90	3.13	7	4.29	0.07	3.60-5.00	1.63
SG: Abaxis Piccolo Instruments	6	1.23	0.19	1.00-1.40	15.83	7	5.83	0.17	4.80-6.80	2.93															
SG: Select Abbott Architect Instruments	5	1.76	0.06	1.50-2.10	2.78	5	7.40	0.12	6.10-8.70	1.62	5	11.18	0.19	9.30-13.10	1.70	5	1.74	0.09	1.40-2.00	5.29	5	4.56	0.06	3.80-5.40	1.09
SG: International OCD Vitros Instruments	10	1.83	0.07	1.50-2.10	3.89	10	6.96	0.11	5.80-8.20	1.57	10	10.60	0.12	8.80-12.40	1.13	9	1.57	0.05	1.30-1.90	3.13	9	4.30	0.07	3.60-5.00	1.63
AG: All Group	62	1.80	0.34	1.50-2.10	18.89	71	6.86	0.46	5.70-8.10	6.67	64	10.35	0.49	8.60-12.20	4.71	51	1.71	0.20	1.40-2.00	11.76	62	4.46	0.33	3.70-5.30	7.33
Iron	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%
PG: Siemens Dimension/AR/ES/RxL/Xpand/EXL	5	171.80	1.30	138.00-206.00	0.76	5	133.80	1.92	107.00-161.00	1.42	5	203.40	2.07	162.00-244.00	1.03	5	29.80	1.30	24.00-36.00	4.33	5	82.40	1.52	66.00-98.00	1.83
PG: Beckman AU	7	186.43	11.49	149.00-223.00	6.18	7	147.29	12.49	118.00-176.00	8.50	7	225.86	19.69	181.00-271.00	8.72	7	31.86	2.27	26.00-38.00	7.19	7	90.29	5.59	72.00-108.00	6.22
SG: Select Beckman Coulter Instruments (Chem)	10	185.90	9.68	149.00-223.00	5.22	10	146.60	10.92	118.00-176.00	7.41	10	223.80	16.84	179.00-269.00	7.50	10	30.80	2.70	25.00-37.00	8.71	10	89.00	5.08	71.00-107.00	5.73
SG: Select Abbott Architect Instruments	6	187.33	3.98	150.00-224.00	2.14	6	140.67	2.81	113.00-169.00	1.99	6	212.67	3.93	170.00-256.00	1.83	6	32.00	1.67	26.00-38.00	5.31	6	86.17	3.25	69.00-103.00	3.84
SG: International OCD Vitros Instruments	5	216.00	10.16	173.00-259.00	4.72	5	152.75	5.44	122.00-184.00	3.53	5	258.25	10.11	206.00-310.00	3.91	5	19.00	2.65	15.00-23.00	13.68	5	86.33	1.53	69.00-103.00	1.74
SG: Select Gen Chem Instruments	24	182.77	9.92	146.00-220.00	5.41	24	142.07	8.71	114.00-170.00	6.13	24	214.43	14.49	171.00-257.00	6.78	24	31.90	3.43	26.00-38.00	10.63	24	86.57	5.19	70.00-104.00	5.98
AG: All Group	35	184.44	15.04	147.00-221.00	8.15	35	141.88	8.88	114.00-170.00	6.27	35	217.64	19.64	174.00-262.00	8.99	31	30.86	3.34	25.00-37.00	10.65	34	85.81	4.68	69.00-103.00	5.47
TIBC	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%
SG: Select Beckman Coulter Instruments (Chem)	5	600.00	18.40	544.80-655.20	3.07	5	271.80	50.74	119.60-424.00	18.67	5	387.40	75.77	160.10-614.70	19.56	5	102.60	6.11	84.30-120.90	5.96	5	186.20	30.36	95.10-277.30	16.31
SG: Select Siemens Instruments	5	556.25	11.90	520.60-592.00	2.14	5	215.25	2.99	206.30-224.30	1.39	5	309.00	12.19	272.40-345.60	3.94	5	72.50	6.61	52.70-92.30	9.12	5	140.75	3.30	130.90-150.70	2.34
PG: Intl OCD Vitros 250,350,400,500,750,950	5	622.33	7.57	599.60-645.00	1.22	5	339.00	28.28	254.20-423.80	8.34	5	564.33	34.56	460.60-668.00	6.12	5	111.33	16.29	62.40-160.20	14.64	5	194.67	9.61	165.90-223.50	4.94
SG: Select Gen Chem Instruments	9	579.61	33.64	478.70-680.50	5.80	9	260.67	50.45	109.30-412.10	19.35	9	349.33	77.12	117.90-580.70	22.08	9	115.17	27.55	32.50-197.90	23.91	9	191.03	37.22	79.30-302.70	19.49
AG: All Group	15	578.30	39.57	459.60-697.00	6.84	15	256.93	55.09	91.60-422.20	21.44	16	368.94	99.14	71.50-666.30	26.87	15	105.03	30.85	12.40-197.60	29.38	15	179.82	38.63	63.90-295.70	21.48
UIBC	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%
AG: All Group	5	394.80	34.05	292.70-496.90	8.62	5	140.00	13.34	100.00-180.00	9.53	5	183.80	17.01	132.80-234.80	9.25	5	73.40	9.89	43.70-103.10	13.47	5	107.40	12.88	68.80-146.00	11.99
746 - Auto Diff 6 (Abbott 5-Part Diff)																									
White Blood Cell Count	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%
PG: Abbott Cell-Dyn Ruby	8	3.10	0.11	2.60-3.60	3.55	8	19.92	0.39	16.90-22.90	1.96	8	7.64	0.22	6.50-8.70	2.89	8	20.09	0.52	17.10-23.10	2.59	8	7.64	0.25	6.50-8.70	3.29
Lymphocyte	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%
PG: Abbott Cell-Dyn Ruby	8	36.61	3.68	25.60-47.60	10.05	8	15.30	0.89	12.60-18.00	5.82	8	27.06	0.96	24.20-30.00	3.54	8	15.36	0.47	14.00-16.80	3.05	8	27.28	0.76	25.00-29.60	2.78
Monocyte	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%
PG: Abbott Cell-Dyn Ruby	8	9.24	3.16	0.00-18.72	34.20	8	3.93	0.42	2.68-5.18	10.64	8	5.57	0.70	3.47-7.67	12.55	8	3.70	0.25	2.95-4.45	6.73	8	5.50	0.44	4.18-6.82	8.00
Granulocyte	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%
PG: Abbott Cell-Dyn Ruby	8	49.26	2.64	41.40-57.20	5.35	8	76.91	0.45	75.60-78.30	0.59	8	62.39	0.66	60.40-64.40	1.06	8	77.26	0.20	76.70-77.90	0.26	8	62.40	0.76	60.10-64.70	1.22
Eosinophil	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%
PG: Abbott Cell-Dyn Ruby	8	4.40	0.61	2.60-6.20	13.86	8	3.34	0.16	2.80-3.80	4.85	8	4.96	0.39	3.90-6.10	7.60	8	3.32	0.23	2.60-4.00	6.97	8	4.76	0.29	3.90-5.70	6.04



AG: All Group	10	7.04	0.57	5.33-8.75	8.10	10	3.02	0.47	1.61-4.43	15.56															
600 - Cholestech Chemistry	CT-1																								
CholesterolTotal (Waived)	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%
PG: Cholestech LDX (Waived)	71	286.33	38.08	172.00-400.00	13.32																				
Cholesterol HDL (Waived)	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%
PG: Cholestech LDX (Waived)	62	70.60	10.42	40.00-102.00	14.65																				
Cholesterol LDL (Waived)	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%
PG: Cholestech LDX (Waived)	42	189.04	36.02	81.00-297.00	19.05																				
Triglycerides (Waived)	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%
PG: Cholestech LDX (Waived)	58	135.83	6.77	116.00-156.00	5.00																				
Glucose (Waived)	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%
PG: Cholestech LDX (Waived)	16	127.69	6.56	108.00-148.00	5.16																				
763 - CoaguChek XS Plus/XS Pro (2 specimens);764 - CoaguChek XS Plus/XS Pro (5 specimens)	CXP-1					CXP-2					CXP-3					CXP-4					CXP-5				
Protime W Blood - Waived	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%
PG: CoaguChek XS Plus	20	37.04	1.87	31.40-42.60	5.05	20	14.28	1.07	11.10-17.50	7.48															
INR	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%
PG: CoaguChek XS	7	3.06	0.11	2.80-3.40	3.55	7	1.16	0.05	1.00-1.40	4.17															
PG: CoaguChek XS Plus	113	3.09	0.19	2.50-3.70	6.13	113	1.18	0.09	0.90-1.50	7.50	47	1.79	0.13	1.40-2.20	7.22	47	1.17	0.12	0.80-1.60	10.00	47	1.79	0.14	1.40-2.20	7.78
PG: CoaguChek XS Pro	13	3.11	0.13	2.70-3.50	4.19	13	1.18	0.04	1.10-1.30	3.33	10	1.84	0.05	1.60-2.00	2.78	10	1.18	0.06	1.00-1.40	5.00	10	1.86	0.07	1.70-2.10	3.68
Protime - Whole Blood	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%
PG: CoaguChek XS Plus	19	37.79	1.40	32.10-43.50	3.70	19	13.97	0.30	11.90-16.10	2.21	7	21.91	0.43	18.60-25.20	1.96	7	13.97	0.24	11.90-16.10	1.71	7	21.93	0.56	18.60-25.20	2.56
PG: CoaguChek XS Pro	10	37.13	1.65	31.50-42.70	4.45	10	14.24	0.31	12.10-16.30	2.18	10	22.01	0.59	18.70-25.30	2.68	10	13.99	0.56	11.90-16.10	4.00	10	22.41	0.70	19.00-25.80	3.13
745 - Auto Diff 5 (Beckman Coulter 5-Part Diff)	DIF-1					DIF-2					DIF-3					DIF-4					DIF-5				
White Blood Cell Count	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%
SG: Select Beckman Coulter Instruments (5-part Diff)	5	4.18	0.13	3.60-4.80	3.10	5	20.80	0.14	17.70-23.90	0.67	5	9.20	0.25	7.80-10.60	2.61	5	20.53	0.17	17.40-23.60	0.83	5	9.28	0.25	7.90-10.70	2.69
Lymphocyte	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%
SG: Select Beckman Coulter Instruments (5-part Diff)	5	38.78	1.21	35.20-42.40	3.12	5	20.60	0.85	18.10-23.20	4.13	5	30.70	1.33	26.70-34.70	4.33	5	21.00	0.49	19.50-22.50	2.33	5	29.95	1.32	26.00-34.00	4.40
Monocyte	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%
SG: Select Beckman Coulter Instruments (5-part Diff)	5	1.13	0.53	0.00-2.71	46.46	5	7.45	0.57	5.73-9.17	7.70	5	4.23	0.90	1.53-6.93	21.28	5	6.93	0.58	5.19-8.67	8.37	5	4.28	0.92	1.51-7.05	21.54
Neutrophil	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%
SG: Select Beckman Coulter Instruments (5-part Diff)	5	51.98	0.81	49.60-54.40	1.56	5	66.22	0.67	64.20-68.20	1.01	5	58.45	0.93	55.70-61.30	1.59	5	66.70	0.73	64.50-68.90	1.11	5	58.97	0.65	57.10-60.90	1.08
Eosinophil	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%
SG: Select Beckman Coulter Instruments (5-part Diff)	5	8.10	0.36	7.00-9.20	4.44	5	5.63	0.86	3.00-8.20	15.36	5	6.53	0.35	5.50-7.50	5.38	5	5.30	0.88	2.70-7.90	16.60	5	6.68	0.48	5.30-8.10	7.16
Basophil	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%
SG: Select Beckman Coulter Instruments (5-part Diff)	5	0.10	0.00	0.10-0.10	0.00	5	0.13	0.06	0.00-0.30	44.62	5	0.13	0.06	0.00-0.30	44.62	5	0.10	0.00	0.10-0.10	0.00	5	0.13	0.05	0.00-0.28	38.46
Red Blood Cell Count	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%
SG: Select Beckman Coulter Instruments (5-part Diff)	5	2.63	0.07	2.47-2.79	2.81	5	5.41	0.17	5.09-5.73	3.14	5	4.19	0.06	3.94-4.44	1.38	5	5.38	0.14	5.06-5.70	2.58	5	4.17	0.07	3.92-4.42	1.68
Hemoglobin	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%
SG: Select Beckman Coulter Instruments (5-part Diff)	5	6.35	0.10	6.00-6.80	1.56	5	17.08	0.22	15.90-18.30	1.29	5	12.20	0.22	11.30-13.10	1.80	5	17.00	0.12	15.80-18.20	0.71	5	12.10	0.22	11.30-12.90	1.82
Hematocrit	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%
SG: Select Beckman Coulter Instruments (5-part Diff)	5	19.97	0.22	18.80-21.20	1.10	5	54.10	0.84	50.90-57.30	1.55	5	38.08	0.54	35.80-40.40	1.42	5	54.08	0.83	50.90-57.30	1.53	5	37.95	0.62	35.70-40.30	1.63
Platelet Count	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%
SG: Select Beckman Coulter Instruments (5-part Diff)	5	76.00	2.94	57.00-95.00	3.82	5	472.50	20.21	355.00-591.00	4.27	5	238.50	7.51	179.00-299.00	3.14	5	464.25	18.95	348.00-580.00	4.07	5	231.75	15.86	174.00-290.00	6.85
MCV	Labs																								



PG: Siemens DCA Vantage	106	9.48	0.80	7.10-11.90	8.42																				
SG: Select Siemens (Bayer) DCA Instruments	118	9.62	0.40	8.40-10.80	4.17																				
PG: Bayer A1CNow+	33	6.79	0.92	4.00-9.60	13.53																				
SG: Select A1C Now Instruments	37	6.81	0.86	4.20-9.40	12.65																				
PG: TOSOH G8	6	10.13	0.20	9.50-10.70	1.98																				
PG: PTS A1CNow	19	6.99	1.05	3.80-10.20	15.00																				
SG: Select Glycohemoglobin Instruments	37	11.32	0.74	9.10-13.50	6.55																				
761 - High Sensitivity CRP	HCR-1					HCR-2																			
hsCRP	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%
AG: All Group	5	0.27	0.02	0.22-0.32	6.67	5	0.09	0.01	0.06-0.12	10.00															
602 - Hemoglobin and Hematocrit - Waived;614 - Auto Diff II (All Other 3-part);610 - Regular Hematology;601 - Waived/PPM Tests II;609 - Waived/PPM Tests I	HD-1					HD-2					HD-3					HD-4					HD-5				
	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%
Granulocyte	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%
PG: Abbott Cell-Dyn 1700	10	85.01	1.21	81.40-88.60	1.42	10	39.22	2.74	31.00-47.40	6.99	10	35.18	2.00	29.20-41.20	5.68	10	69.55	1.67	64.60-74.60	2.40	10	82.73	0.83	80.20-85.20	1.00
PG: Abbott Cell-Dyn 1800	11	85.11	1.72	79.90-90.30	2.02	11	39.12	2.36	32.00-46.20	6.04	11	35.25	1.81	29.90-40.70	5.13	11	70.17	1.61	65.40-75.00	2.29	11	82.42	0.75	80.20-84.70	0.91
PG: CDS Medonic 620	9	81.69	0.82	79.20-84.20	1.00	8	29.85	2.00	23.90-35.90	6.69	9	26.81	0.99	23.80-29.80	3.69	9	61.41	2.48	54.00-68.80	4.04	9	79.20	0.43	77.90-80.50	0.54
PG: Beckman Coulter AcT Diff/Diff 2	167	82.78	4.99	67.80-97.80	6.03	166	33.72	3.07	24.50-42.90	9.11	166	31.91	2.74	23.70-40.10	8.59	166	65.60	5.49	49.10-82.10	8.37	166	80.75	7.19	59.20-102.40	8.90
PG: Horiba Micros/45/60	52	89.58	0.77	87.30-91.90	0.86	52	50.99	2.17	44.50-57.50	4.25	52	43.14	2.07	36.90-49.30	4.80	52	77.33	1.00	74.30-80.30	1.29	52	84.86	0.64	83.00-86.80	0.75
SG: Select Abbott Cell Dyn Instruments (3-part Diff)	112	86.81	8.95	60.00-113.60	10.31	112	45.56	5.58	28.90-62.30	12.24	112	38.56	4.23	25.90-51.30	10.96	112	73.59	6.77	53.30-93.90	9.20	112	82.22	6.52	62.60-101.80	7.93
PG: Drew Scientific D3	13	87.95	0.60	86.20-89.80	0.68	13	43.02	1.31	39.00-47.00	3.07	13	37.47	0.91	34.80-40.20	2.43	13	72.93	0.88	70.30-75.50	1.21	13	83.72	0.79	81.30-86.10	0.94
PG: Abbott Cell-Dyn Emerald	91	87.22	9.87	57.60-116.80	11.32	91	47.04	5.03	31.90-62.10	10.70	91	39.33	4.25	26.50-52.00	10.81	91	74.45	7.21	52.80-96.00	9.69	91	82.14	7.23	60.40-103.80	8.81
PG: CDS M Series	141	81.43	0.95	78.60-84.30	1.17	141	30.10	2.38	23.00-37.20	7.91	142	27.90	1.89	22.20-33.60	6.77	142	61.30	1.66	56.30-66.30	2.71	142	78.90	5.30	63.00-94.80	6.73
SG: Select CDS Instruments	150	81.45	0.94	78.60-84.20	1.15	149	30.08	2.36	23.00-37.20	7.84	151	27.84	1.86	22.20-33.40	6.69	151	61.30	1.71	56.20-66.40	2.79	151	78.92	5.14	63.50-94.30	6.51
AG: All Group	495	84.13	5.89	66.40-101.80	7.00	481	37.18	7.70	14.10-60.30	20.70	481	33.32	5.26	17.50-49.10	15.80	495	67.52	7.33	45.50-89.50	10.86	495	81.03	6.20	62.40-99.60	7.65
Hematocrit	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%
PG: Abbott Cell-Dyn 1700	11	23.16	0.59	21.80-24.60	2.59	11	29.47	0.81	27.70-31.30	2.75	11	35.95	0.67	33.80-38.20	1.86	11	21.67	0.61	20.40-23.00	2.81	11	53.56	1.20	50.40-56.80	2.24
PG: Abbott Cell-Dyn 1800	8	24.25	0.83	22.80-25.80	3.42	11	30.14	0.61	28.30-31.90	2.03	11	36.75	1.08	34.50-38.90	2.94	11	22.52	0.82	21.10-23.90	3.64	11	53.85	2.32	50.60-57.00	4.31
PG: CDS Medonic 620	10	22.78	1.06	21.40-24.20	4.69	10	27.78	0.62	26.10-29.50	2.23	10	34.27	0.67	32.20-36.40	1.95	10	20.63	0.62	19.40-21.80	3.01	10	54.47	1.45	51.20-57.80	2.66
PG: Beckman Coulter AcT Diff/Diff 2	175	23.37	0.82	22.00-24.80	3.50	176	29.12	1.03	27.40-30.80	3.54	176	35.47	1.81	33.40-37.60	5.10	175	21.38	0.64	20.10-22.70	2.99	176	53.19	3.67	50.00-56.40	6.90
PG: Horiba Micros/45/60	52	21.18	0.61	19.90-22.50	2.88	52	27.37	0.75	25.80-29.00	2.74	52	33.33	1.18	31.30-35.30	3.54	52	19.27	0.77	18.10-20.50	3.99	52	51.24	2.22	48.10-54.30	4.34
SG: Select Abbott Cell Dyn Instruments (3-part Diff)	110	23.92	0.83	22.50-25.30	3.47	112	29.95	0.98	28.10-31.70	3.28	114	36.34	1.19	34.10-38.50	3.28	114	22.10	0.82	20.80-23.40	3.71	114	54.74	1.81	51.40-58.00	3.31
PG: Drew Scientific D3	12	23.28	0.89	21.90-24.70	3.82	12	29.31	1.04	27.50-31.10	3.55	11	35.75	1.11	33.60-37.80	3.11	13	21.47	1.01	20.20-22.80	4.70	13	53.08	2.29	49.90-56.30	4.31
PG: Abbott Cell-Dyn Emerald	91	23.99	0.81	22.60-25.40	3.38	90	29.98	1.03	28.20-31.80	3.43	92	36.34	1.25	34.10-38.50	3.44	92	22.10	0.82	20.80-23.40	3.71	92	54.99	1.74	51.70-58.30	3.16
PG: CDS M Series	144	22.34	0.55	21.00-23.60	2.47	142	27.82	0.59	26.10-29.50	2.12	144	34.13	0.73	32.10-36.10	2.14	144	20.50	0.78	19.30-21.70	3.80	144	54.68	1.28	51.40-58.00	2.36
SG: Select CDS Instruments	154	22.37	0.60	21.10-23.70	2.68	152	27.82	0.59	26.10-29.50	2.12	154	34.13	0.72	32.10-36.10	2.11										



PG: CDS M Series	141	73.24	1.29	69.30-77.10	1.76	141	71.83	2.35	64.80-78.80	3.27	141	79.81	1.23	76.10-83.50	1.54	141	80.62	1.36	76.50-84.70	1.69	139	92.75	1.42	88.50-97.10	1.53
SG: Select CDS Instruments	150	73.21	1.32	69.20-77.20	1.80	150	71.80	2.32	64.80-78.80	3.23	150	79.78	1.29	75.90-83.70	1.62	150	80.60	1.41	76.40-84.80	1.75	148	92.70	1.52	88.10-97.30	1.64
AG: All Group	499	75.22	3.12	65.80-84.60	4.15	498	74.13	2.96	65.20-83.00	3.99	498	82.41	3.10	73.10-91.70	3.76	498	82.74	3.42	72.40-93.00	4.14	495	93.59	3.15	84.10-103.00	3.37
Monocyte/Mixed	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%
PG: Abbott Cell-Dyn 1700	10	7.18	0.78	4.90-9.50	10.83	10	21.54	1.61	16.70-26.30	7.49	10	19.37	2.65	11.40-27.30	13.66	10	12.12	0.58	10.40-13.80	4.79	10	5.81	0.47	4.40-7.20	8.10
PG: Abbott Cell-Dyn 1800	10	7.61	0.57	5.90-9.30	7.50	10	25.31	1.05	22.10-28.50	4.15	10	24.50	1.47	20.10-28.90	6.00	10	14.17	0.57	12.50-15.90	4.01	10	7.34	0.22	6.60-8.00	3.01
PG: CDS Medonic 620	9	6.91	0.95	4.10-9.80	13.77	9	11.58	1.29	7.70-15.50	11.12	9	10.24	1.33	6.20-14.20	13.04	9	9.69	1.51	5.10-14.30	15.67	9	5.04	0.59	3.20-6.80	11.80
PG: Beckman Coulter AcT Diff/Diff 2	168	7.30	1.24	3.60-11.00	16.99	167	19.33	2.61	11.50-27.10	13.52	166	14.21	1.79	8.80-19.60	12.61	167	11.54	1.41	7.30-15.70	12.26	167	4.37	0.87	1.80-7.00	19.77
PG: Horiba Micros/45/60	52	3.37	0.29	2.50-4.30	8.53	52	19.07	1.47	14.70-23.50	7.70	52	21.93	1.64	17.00-26.80	7.49	52	8.69	0.59	6.90-10.50	6.90	52	4.90	0.41	3.70-6.10	8.37
SG: Select Abbott Cell Dyn Instruments (3-part Diff)	89	2.56	0.30	1.70-3.50	11.54	88	12.29	1.69	7.20-17.40	13.74	101	13.77	2.68	5.80-21.80	19.42	91	6.61	1.18	3.10-10.10	17.88	86	3.29	0.49	1.80-4.80	14.85
PG: Drew Scientific D3	13	3.29	0.34	2.30-4.30	10.30	13	13.74	0.59	11.90-15.50	4.31	13	12.78	0.93	10.00-15.60	7.27	13	7.45	0.57	5.70-9.10	7.70	13	3.25	0.29	2.30-4.10	9.06
PG: Abbott Cell-Dyn Emerald	91	2.51	0.42	1.20-3.80	16.80	91	12.21	2.00	6.20-18.20	16.39	91	12.95	2.21	6.40-19.60	17.00	91	6.44	0.99	3.40-9.40	15.47	91	3.30	0.56	1.60-5.00	16.97
PG: CDS M Series	140	7.01	0.56	5.30-8.70	8.00	141	10.66	1.79	5.30-16.10	16.73	141	10.03	1.45	5.70-14.30	14.50	141	9.76	0.99	6.80-12.80	10.10	141	4.99	0.54	3.40-6.60	10.80
SG: Select CDS Instruments	149	7.01	0.59	5.20-8.80	8.43	150	10.72	1.77	5.40-16.00	16.54	150	10.05	1.44	5.70-14.30	14.40	150	9.75	1.02	6.70-12.90	10.41	150	4.99	0.54	3.40-6.60	10.80
AG: All Group	469	5.77	2.17	0.00-12.30	37.41	467	15.20	4.36	2.20-28.30	28.62	480	13.65	3.79	2.20-25.00	27.87	471	9.64	2.11	3.30-15.90	21.98	458	4.38	0.79	2.00-6.80	17.95
Platelet Count	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%
PG: Abbott Cell-Dyn 1700	11	398.64	27.88	299.00-499.00	6.99	11	113.27	10.40	85.00-141.00	9.20	11	258.18	19.02	194.00-323.00	7.36	11	73.55	8.03	56.00-93.00	10.81	11	589.45	38.85	442.00-736.00	6.60
PG: Abbott Cell-Dyn 1800	11	388.09	49.75	291.00-485.00	12.81	11	115.45	6.47	86.00-144.00	5.65	11	267.45	14.86	200.00-334.00	5.58	11	79.73	7.27	60.00-100.00	9.13	11	602.82	46.32	452.00-754.00	7.68
PG: CDS Medonic 620	9	378.67	16.75	284.00-474.00	4.41	10	101.40	7.14	76.00-126.00	7.03	10	245.90	7.81	185.00-308.00	3.17	10	75.80	4.02	57.00-95.00	5.26	10	532.30	18.37	399.00-665.00	3.46
PG: Beckman Coulter AcT Diff/Diff 2	177	410.59	28.07	308.00-514.00	6.84	176	111.80	7.94	84.00-140.00	7.05	176	253.13	12.98	190.00-316.00	5.14	175	80.14	7.97	60.00-100.00	10.00	176	569.84	46.05	428.00-713.00	8.07
PG: Horiba Micros/45/60	52	411.65	20.49	309.00-515.00	4.98	52	118.79	10.03	89.00-149.00	8.40	52	254.21	22.79	191.00-318.00	8.98	52	84.35	7.78	63.00-105.00	9.29	52	534.42	26.80	401.00-668.00	5.02
SG: Select Abbott Cell Dyn Instruments (3-part Diff)	114	400.04	31.24	300.00-500.00	7.80	114	120.27	13.44	90.00-150.00	11.17	114	254.59	17.89	191.00-319.00	7.02	112	82.21	14.57	62.00-103.00	17.80	114	545.67	43.46	410.00-683.00	7.97
PG: Drew Scientific D3	13	395.62	22.91	297.00-495.00	5.78	13	114.69	10.79	86.00-144.00	9.39	13	245.54	16.74	185.00-308.00	6.79	13	78.85	9.92	59.00-99.00	12.53	13	523.62	32.67	393.00-655.00	6.24
PG: Abbott Cell-Dyn Emerald	92	401.63	28.88	302.00-503.00	7.19	92	121.68	14.04	92.00-153.00	11.48	92	252.62	17.55	190.00-316.00	6.96	90	83.58	15.48	63.00-105.00	18.45	92	533.60	34.20	401.00-668.00	6.40
PG: CDS M Series	144	369.63	16.02	278.00-463.00	4.32	142	101.63	6.63	77.00-128.00	6.47	144	242.31	11.00	182.00-303.00	4.55	144	72.94	4.79	55.00-91.00	6.58	144	519.69	20.13	390.00-650.00	3.87
SG: Select CDS Instruments	153	370.16	16.14	278.00-463.00	4.35	152	101.62	6.64	77.00-128.00	6.47	154	242.54	10.84	182.00-304.00	4.44	154	73.12	4.78	55.00-91.00	6.58	154	520.51	20.20	391.00-651.00	3.88
AG: All Group	510	395.75	30.43	297.00-495.00	7.68	508	111.44	11.85	83.00-139.00	10.63	510	250.17	15.82	188.00-313.00	6.32	507	78.84	10.01	59.00-99.00	12.66	510	544.61	42.31	409.00-681.00	7.76
RDW	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%
PG: Abbott Cell-Dyn 1700	10	18.24	0.55	16.60-19.80	3.02	10	18.09	0.54	16.50-19.70	2.98	10	16.40	0.77	14.10-18.70	4.63	10	15.77	0.64	13.90-17.70	4.05	10	15.39	0.47	14.00-16.80	3.05
PG: Abbott Cell-Dyn 1800	10	17.47	0.32	16.50-18.50	1.83	10	17.13	0.50	15.60-18.60	2.87	10	15.67	0.69	13.60-17.80	4.39	10	15.34	0.54	13.70-16.90	3.53	10	14.48	0.32	13.50-15.50	2.21
PG: CDS Medonic 620	9	14.61	0.68	12.60-16.60	4.66	9	14.60	0.95	11.80-17.40	6.51	9	13.67	0.73	11.50-15.90	5.33	9	12.60	0.74	10.40-14.80	5.87	9	13.74	0.76	11.40-16.00	5.55
PG: Beckman Coulter AcT Diff/Diff 2	168	15.85	0.37	14.70-16.90	2.34	168	15.66	0.32	14.70-16.70	2.04	168	14.11	0.50	12.60-15.60	3.48	168	13.98	0.30	13.10-14.90	2.14	168	13.47	0.48	12.10-14.90	3.56
PG: Horiba Micros/45/60	45	13.72	0.38	12.60-14.80	2.77	45	13.79	0.45	12.50-15.20	3.26	45	12.95	0.43	11.70-14.30	3.31	45	12.70	0.38	11.60-13.80	2.99	45	12.28	0.32	11.30-13.30	2.60
SG: Select Abbott Cell Dyn Instruments (3-part Diff)	111	15.34	1.48	10.90-19.70	9.67	111	15.48	1.30	11.60-19.40	8.39	111	14.25	1.16	10.80-17.60	8.10	111	13.86	1.03	10.80-17.00	7.41	111	13.37	1.05	10.30-16.60	7.84
PG: Drew Scientific D3	13	15.03	0.91	12.30-17.70	6.00	13	15.39	0.89	12.70-18.10	5.78	13	14.08	0.51	12.60-15.60	3.62	13	13.48	0.55	11.80-15.20	4.07	13	13.15	0.69	11.10-15.30	5.23
PG: Abbott Cell-Dyn Emerald	91	14.79	0.94	12.00-17.60	6.35	91	15.01	0.86	12.40-17.60	5.73	91	13.85	0.79	11.50-16.30	5.68	91	13.48	0.67	11.50-15.50	4.96	91	13.03	0.78	10.70-15.30	6.00
PG: CDS M Series	141	14.50	1.09	11.20-17.80	7.52	141	14.41	0.92	11.60-17.20	6.39	141	13.54	0.87	10.90-16.10	6.44	141	12.70	0.90	10.00-15.40	7.09	141	13.55	1.12	10.10-16.90	8.30
SG: Select CDS Instruments	150	14.51	1.07	11.30-17.70	7.38	150	14.42	0.92	11.60-17.20	6.39	150	13.55	0.86	11.00-16.20	6.32	150	12.70	0.89	10.00-15.40	7.01	150	13.56	1.10	10.30-16.90	8.09
AG: All Group	488	15.11	1.20	11.50-18.70	7.95	488	15.07	1.09	11.80-18.40	7.22	488	13.87	0.90	11.20-16.60	6.47	488	13.43	0.95	10.60-16.30	7.09	488	13.36	0.93	10.60-16.20	6.94
Red Blood Cell Count	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%
PG: Abbott Cell-Dyn 1700	11	3.12	0.08	2.93-3.31	2.47	11	3.97	0.11	3.73-4.21	2.72	11	4.33	0.11	4.07-4.59	2.45	11	2.61	0.06	2.45-2.77	2.45	11	5.70	0.16	5.36-6.04	2.82
PG: Abbott Cell-Dyn 1800	7	3.10	0.12	2.91-3.29	3.97	11	3.89	0.10	3.66-4.12	2.52	10	4.25	0.13	4.00-4.50	3.06	11	2.61	0.11	2.45-2.77	4.02	11	5.48	0.23	5.15-5.81	4.27
PG: CDS Medonic 620	10	3.13	0.11	2.94-3.32	3.61	10	3.90	0.06	3.67-4.13	1.62	10	4.32	0.05	4.06-4.58	1.23	10	2.57	0.05	2.42-2.72	1.87	10	5.93	0.05	5.57-6.29	0.79
PG: Beckman Coulter AcT Diff/Diff 2	174	3.09	0.09	2.90-3.28	2.88	176	3.92	0.14	3.68-4.16	3.47	176	4.30	0.22	4.04-4.56	5.07	174	2.59	0.07	2.43-2.75	2.51	176	5.74	0.35	5.40-6.08	6.15
PG: Horiba Micros/45/60	52	3.03	0.07	2.85-3.21	2.15	52	3.85	0.09	3.62-4.08	2.23	52	4.22	0.10	3.97-4.47	2.30	52	2.50	0.09	2.35-2.65	3.76	52	5.72	0.23	5.38-6.06	3.93
SG: Select Abbott Cell Dyn Instruments (3-part Diff)	109	3.03	0.10	2.85-3.21	3.33	112	3.84	0.12	3.61-4.07	3.15	113	4.19	0.14	3.94-4.44	3.22										



PG: Horiba Micros/45/60	44	1.90	0.13	1.50-2.30	6.84	44	1.12	0.14	0.70-1.50	12.73	44	1.74	0.21	1.10-2.30	12.35	45	1.34	0.13	0.90-1.70	10.00	45	2.05	0.16	1.60-2.60	7.62
Monocyte/Mixed Abs #	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%
PG: Horiba Micros/45/60	45	0.88	0.07	0.70-1.10	7.78	45	0.70	0.06	0.50-0.90	8.57	45	1.09	0.11	0.80-1.40	10.00	45	0.82	0.07	0.60-1.00	8.75	45	0.96	0.10	0.70-1.30	10.00
MCH	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%
PG: Horiba Micros/45/60	45	27.87	0.83	25.40-30.40	2.97	45	25.71	0.48	24.30-27.10	1.87	45	28.81	0.59	27.00-30.60	2.05	45	30.06	0.63	28.20-32.00	2.06	45	31.93	0.72	29.70-34.10	2.26
MCHC	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%
PG: Horiba Micros/45/60	45	39.90	1.27	36.10-43.70	3.18	45	36.23	0.86	33.60-38.80	2.38	45	36.50	0.83	34.00-39.00	2.27	45	39.05	1.00	36.00-42.00	2.56	45	35.68	0.98	32.80-38.60	2.75
Hematocrit (Waived)	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%
PG: Microhematocrit, Spun	18	20.83	2.01	14.80-26.80	9.66	11	25.09	1.58	20.40-29.80	6.29															
SG: Select Microhematocrit Instruments	27	21.66	2.50	14.20-29.20	11.52																				
AG: All Group	27	21.66	2.50	14.20-29.20	11.52																				
Hemoglobin (Waived)	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%
PG: HemoCue	45	7.81	0.65	5.80-9.80	8.33																				
PG: HemoCue 201	97	7.80	0.69	5.70-9.90	8.85																				
PG: Alere HemoPoint H2	39	7.80	0.98	4.90-10.70	12.56	5	10.03	0.49	8.50-11.50	4.90															
SG: Select Waived Hemoglobin Instruments	183	7.80	0.74	5.60-10.00	9.49																				
AG: All Group	186	7.79	0.75	5.50-10.10	9.62																				
608 - Hemocue Glucose/Hemoglobin	HQ-1																								
Glucose Whole Blood (Hemocue)	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%
PG: HemoCue	40	109.43	15.35	63.00-155.00	14.13																				
PG: HemoCue 201	137	111.03	16.17	62.00-160.00	14.59																				
Hemoglobin (Hemocue)	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%
PG: HemoCue	23	12.99	0.78	10.70-15.30	6.00																				
PG: HemoCue 201	66	12.84	0.68	10.80-14.80	5.31																				
606 - Hemocue WBC	HQW-1					HQW-2					HQW-3					HQW-4					HQW-5				
White Blood Cell Count	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%
PG: Hemocue WBC System	11	8.22	0.36	7.00-9.40	4.39	11	21.94	1.01	18.60-25.20	4.61	11	2.82	0.29	2.40-3.20	10.36	11	21.71	0.64	18.40-25.00	2.90	11	2.86	0.16	2.50-3.30	5.52
700 - i-Stat Chemistry;721 - i-Stat Chemistry & Blood Gases;722 - i-Stat Chemistry - Waived;723 - i-Stat Blood Gases	IST-1					IST-2					IST-3					IST-4					IST-5				
Blood Gas - pCO2	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%
PG: Abbott i-Stat	8	42.84	1.80	38.00-48.00	4.19	8	45.16	1.95	40.00-50.00	4.44	8	49.17	1.76	44.00-54.00	3.67	8	69.64	3.92	64.00-76.00	5.57	8	43.11	1.08	38.00-48.00	2.56
SG: International Abbott i-Stat Instruments	12	42.81	1.52	38.00-48.00	3.49	12	45.28	1.59	40.00-50.00	3.56	12	49.20	1.66	44.00-54.00	3.47	12	70.43	3.37	64.00-76.00	4.86	12	43.33	0.96	38.00-48.00	2.33
Blood Gas - pH	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%
PG: Abbott i-Stat	8	7.29	0.01	7.25-7.33	0.07	8	7.31	0.01	7.27-7.35	0.11	8	7.32	0.01	7.28-7.36	0.08	8	6.92	0.01	6.88-6.96	0.22	8	7.29	0.01	7.25-7.33	0.10
SG: International Abbott i-Stat Instruments	12	7.29	0.01	7.25-7.33	0.08	12	7.31	0.01	7.27-7.35	0.10	12	7.32	0.01	7.28-7.36	0.07	12	6.92	0.01	6.88-6.96	0.22	12	7.29	0.01	7.25-7.33	0.08
Blood Gas - pO2	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%
PG: Abbott i-Stat	8	77.88	4.26	65.00-91.00	5.51	8	168.00	8.52	143.00-194.00	5.06	8	193.63	13.51	154.00-235.00	6.96	8	62.00	6.33	43.00-81.00	10.16	7	75.86	3.89	64.00-88.00	5.13
SG: International Abbott i-Stat Instruments	12	77.83	3.66	67.00-89.00	4.74	12	170.25	7.76	147.00-193.00	4.59	12	195.25	11.25	161.00-229.00	5.79	12	59.08	6.87	38.00-80.00	11.69	11	75.91	3.56	65.00-87.00	4.74
Chloride (CL-)	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%
PG: Abbott i-Stat / Thermal	5	80.00	0.71	76.00-84.00	0.89	5	97.40	0.55	92.00-102.00	0.56	5	92.00	0.71	87.00-97.00	0.77	5	112.00	0.71	106.40-117.60	0.63	5	79.60	0.55	76.00-84.00	0.69
PG: Intl Abbott i-Stat / Thermal	5	80.00	0.00	76.00-84.00	0.00	5	97.75	0.50	93.00-103.00	0.51	5	91.25	0.50	86.00-96.00	0.55	5	112.25	0.50	106.70-117.90	0.45	5	80.00	0.00	76.00-84.00	0.00
PG: Abbott i-Stat (Waived) / Thermal	18	79.44	1.04	75.00-83.00	1.31																				
CO2 (Total CO2 Bicarb)	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%
PG: Abbott i-Stat / Thermal	6	23.50	1.05	20.40-26.60	4.47	6	25.17	0.75	22.90-27.40	2.98	6	28.00	1.55	23.40-32.60	5.54	6	18.17	0.98	15.30-21.10	5.38	6	23.00	1.26	19.20-26.80	5.48
PG: Intl Abbott i-Stat / Thermal	5	23.00	1.41	18.80-27.20	6.13	5	24.50	1.29	20.60-28.40	5.27	5	26.50	1.29	22.60-30.40	4.87	5	17.00	1.16	13.60-20.40	6.76	5	23.25	0.96	20.40-26.20	4.12
PG: Abbott i-Stat (Waived) / Thermal	20	23.50	0.95	20.60-26.40	4.04																				
Glucose	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%
PG: Abbott i-Stat / Thermal	8	50.63	0.74	45.00-57.00	1.46	8	197.25	1.17	177.00-217.00	0.61	8	87.13	1.13	78.40-95.80	1.30	8	36.50	0.54	31.00-43.00	1.45	8	50.75	1.03	45.00-57.00	2.05
PG: Intl Abbott i-Stat / Thermal	5	51.36	1.04	45.00-57.00	2.02	5	195.07	1.73	176.00-215.00	0.89	5	86.95	0.90	78.20-95.60	1.04	5	37.39	0.90	31.00-43.00	2.43	5	50.91	0.90	45.00-57.00	1.77
PG: Abbott i-Stat (Waived) / Thermal	20	50.90	1.17	45.00-57.00	2.35																				
Hematocrit	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%
PG: Abbott i-Stat / Thermal	5	15.40	0.55	14.50-16.30	3.57	5	32.80	0.84	30.80-34.80	2.56	5	36.80	0.84	34.60-39.00	2.28	5	19.00	0.71	17.90-20.10	3.74	5	15.20	0.45	14.30-16.10	2.96
PG: Intl Abbott i-Stat / Thermal	5	15.00	0.00	14.10-15.90	0.00	5	32.00	0.00	30.10-33.90	0.00	5	36.00	0.00	33.80-38.20	0.00	5	19.00	0.00	17.90-20.10	0.00	5	15.00	0.00	14.10-15.90	0.00
PG: Abbott i-Stat (Waived) / Thermal	5	15.20	0.45	14.30-16.10	2.96																				
Hemoglobin	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%
PG: Abbott i-Stat / Thermal	5	5.14	0.15	4.70-5.50	2.94	5	11.24	0.45	10.40-12.00	4.02	5	12.52	0.34	11.60-13.40	2.64	5	6.52	0.32	6.00-7.00	4.92	5	5.08	0.04	4.70-5.50	0.78
PG: Intl Abbott i-Stat / Thermal	5	5.10	0.00	4.70-5.50	0.00	5	10.90	0.00	10.10-11.70	0.00	5	12.20	0.00	11.30-13.10	0.00	5	6.50	0.00	6.00-7.00	0.00	5	5.10	0.00	4.70-5.50	0.00
PG: Abbott i-Stat (Waived) / Thermal	5	5.16	0.13	4.80-5.60	2.50																				
Ionized Calcium	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%
PG: Abbott i-Stat / Thermal	6	0.60	0.01	0.58-0.62	1.00	6	0.97	0.02	0.92-1.02	1.86	6	1.30	0.01	1.28-1.32	0.46	6	2.48	0.05	2.34-2.62	1.94	6	0.60	0.01	0.58-0.62	1.00



PG: Abbott i-Stat / Thermal	8	126.13	0.35	122.10-130.10	0.32	8	155.88	0.64	152.00-160.00	0.41	8	137.88	0.64	134.00-142.00	0.46	8	131.75	0.46	128.00-136.00	0.38	8	126.38	0.52	122.40-130.40	0.41
PG: Intl Abbott i-Stat / Thermal	5	126.00	0.00	122.00-130.00	0.00	5	155.00	0.00	151.00-159.00	0.00	5	137.50	0.58	134.00-142.00	0.43	5	131.75	0.50	128.00-136.00	0.38	5	125.75	0.50	121.80-129.80	0.40
PG: Abbott i-Stat (Waived) / Thermal	20	126.50	0.69	122.50-130.50	0.55																				
Urea Nitrogen (BUN)	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%
PG: Abbott i-Stat / Thermal	5	9.00	0.00	7.00-11.00	0.00	5	53.00	0.00	48.20-57.80	0.00	5	34.00	0.00	30.90-37.10	0.00	5	5.00	0.00	3.00-7.00	0.00	5	9.00	0.00	7.00-11.00	0.00
PG: Intl Abbott i-Stat / Thermal	5	8.89	0.14	6.90-10.90	1.57	5	53.15	0.27	48.40-58.00	0.51	5	34.17	0.23	31.10-37.30	0.67	5	5.39	0.14	3.40-7.40	2.59	5	9.03	0.14	7.00-11.00	1.56
PG: Abbott i-Stat (Waived) / Thermal	20	9.05	0.22	7.10-11.10	2.42																				
HCO3	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%
PG: Abbott i-Stat	5	20.33	0.64	18.40-22.20	3.15	5	22.60	1.28	18.80-26.40	5.66	5	24.93	0.90	22.20-27.60	3.61	5	14.17	0.74	12.00-16.40	5.21	5	20.30	0.61	18.50-22.10	3.00
PG: Intl Abbott i-Stat	5	20.65	0.60	18.90-22.50	2.90	5	22.65	0.38	21.60-23.80	1.67	5	24.90	0.85	22.30-27.40	3.41	5	14.23	0.26	13.40-15.00	1.83	5	20.90	0.27	20.10-21.70	1.29
Lactate	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%
PG: Abbott i-Stat	5	6.38	0.00	6.38-6.38	0.00	5	2.22	0.00	2.22-2.22	0.00	5	0.72	0.00	0.72-0.72	0.00	5	4.95	0.00	4.95-4.95	0.00	5	6.38	0.00	6.38-6.38	0.00
PG: Intl Abbott i-Stat	5	6.51	0.10	6.20-6.82	1.60	5	2.29	0.02	2.24-2.34	0.79	5	0.77	0.01	0.73-0.81	1.69	5	5.09	0.02	5.02-5.16	0.47	5	6.48	0.05	6.33-6.63	0.76
Anion Gap	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%
PG: Abbott i-Stat (Waived) / Thermal	7	27.00	1.29	23.00-31.00	4.81																				
Creatinine	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%
PG: Abbott i-Stat / Thermal	10	0.37	0.05	0.10-0.70	12.50	11	4.64	0.13	3.90-5.30	2.83	11	2.48	0.09	2.10-2.90	3.60	11	6.97	0.26	6.00-8.10	3.71	10	0.36	0.05	0.10-0.70	12.50
PG: Intl Abbott i-Stat / Thermal	5	0.37	0.03	0.10-0.70	7.50	5	4.53	0.06	3.80-5.20	1.33	5	2.43	0.07	2.00-2.80	2.92	5	6.75	0.33	5.80-7.80	4.85	5	0.37	0.03	0.10-0.70	7.50
PG: Abbott i-Stat (Waived) / Thermal	23	0.36	0.05	0.10-0.70	12.50																				
716 - Blood Lead - Waived Methods;709 - Blood Lead	LE-1					LE-2					LE-3					LE-4					LE-5				
Blood Lead (Waived)	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%
PG: ESA LeadCare II (Waived)	49	20.71	2.41	13.50-27.90	11.64																				
662 - Neonatal Bilirubin (5 specimens);663 - Neonatal Bilirubin Add-On	NB-1					NB-2					NB-3					NB-4					NB-5				
Neonatal Bilirubin (Total)	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%
PG: Reichert (Leica) Unistat Bilirubinometer	25	12.07	0.56	9.70-14.50	4.63	5	0.10	0.00	0.00-0.50	0.00	25	19.04	1.69	15.20-22.80	8.89	25	17.62	0.49	14.10-21.10	2.78	25	7.63	0.19	6.10-9.10	2.50
AG: All Group	33	11.90	0.67	9.50-14.30	5.63	7	0.17	0.10	0.00-0.60	50.00	31	18.90	1.59	15.10-22.70	8.41	31	17.19	1.28	13.80-20.60	7.44	31	7.62	0.19	6.10-9.10	2.50
748 - i-Stat Prottime	PTI-1					PTI-2					PTI-3					PTI-4					PTI-5				
Prottime - Whole Blood	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%
PG: Abbott i-Stat Prottime	6	25.23	1.29	21.40-29.00	5.12	6	14.63	0.77	12.40-16.80	5.27	6	16.07	4.14	13.70-18.50	25.71	6	26.45	2.56	22.50-30.50	9.66	6	25.88	0.79	22.00-29.80	3.05
PG: Abbott Coaousense	5	34.90	2.39	29.70-40.10	6.85	5	15.25	1.77	13.00-17.60	11.57	5	14.50	0.36	12.30-16.70	2.48	5	36.83	5.10	31.30-42.30	13.86	5	35.57	2.82	30.30-40.90	7.92
SG: International Abbott i-Stat Instruments	7	25.67	1.65	21.80-29.60	6.42	7	14.73	0.74	12.50-16.90	5.03	6	14.53	0.46	12.30-16.70	3.17	7	26.67	2.41	22.70-30.70	9.03	7	26.16	1.02	22.30-30.10	3.89
SG: International Abbott Coaousense Instruments	5	35.90	2.79	30.50-41.30	7.77	5	15.03	1.30	12.80-17.30	8.73	5	14.72	0.54	12.50-16.90	3.67	5	37.38	4.59	31.80-43.00	12.25	5	36.13	2.56	30.70-41.50	7.09
INR	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%
PG: Abbott i-Stat Prottime	7	2.16	0.11	1.90-2.50	5.00	7	1.24	0.08	1.00-1.40	6.67	6	1.22	0.04	1.10-1.30	3.33	7	2.29	0.23	1.60-3.00	10.00	7	2.23	0.08	2.00-2.40	3.64
PG: Abbott Coaousense	5	3.03	0.48	1.60-4.40	16.00	5	1.30	0.14	0.90-1.70	10.77	5	1.20	0.00	1.20-1.20	0.00	5	2.95	0.37	1.90-4.10	12.33	5	3.05	0.47	1.70-4.50	15.16
SG: International Abbott i-Stat Instruments	12	2.20	0.14	1.80-2.60	5.91	12	1.23	0.08	1.00-1.40	6.67	11	1.19	0.07	1.00-1.40	5.83	12	2.28	0.19	1.70-2.90	8.26	12	2.24	0.08	2.00-2.40	3.64
PG: Intl Abbott i-Stat Prottime	5	2.26	0.15	1.80-2.80	6.52	5	1.22	0.08	1.00-1.40	6.67	5	1.16	0.09	0.90-1.50	7.50	5	2.26	0.13	1.90-2.70	5.65	5	2.26	0.09	2.00-2.60	3.91
SG: International Abbott Coaousense Instruments	5	3.04	0.42	1.70-4.30	14.00	5	1.27	0.12	0.90-1.70	9.23	5	1.23	0.05	1.00-1.40	4.17	5	3.00	0.34	2.00-4.00	11.33	5	3.04	0.40	1.80-4.20	13.33
686 - Special Chemistry;667 - PSA Add-on	SC-1					SC-2																			
CEA	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range											



PG: Beckman ACCESS	10	1171.00	58.66	995.00-1347.00	5.01	10	141.50	11.32	108.00-176.00	7.96															
SG: Select Beckman Coulter Instruments (Chem)	15	1199.13	103.82	888.00-1510.00	8.66	15	141.20	11.82	106.00-176.00	8.37															
SG: Select Abbott Architect Instruments	11	1669.27	112.40	1332.00-2006.00	6.73	11	183.00	11.35	149.00-217.00	6.17															
AG: All Group	38	1430.82	288.19	566.00-2296.00	20.14	38	173.71	33.88	72.00-276.00	19.48															
Progesterone	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%
AG: All Group	12	39.17	7.84	15.70-62.70	20.00	13	2.41	0.48	1.00-3.80	20.00															
Prolactin	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%
AG: All Group	10	55.94	4.89	41.20-70.60	8.75	10	4.11	0.84	1.60-6.60	20.49															
PSA	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%
PG: Beckman ACCESS	20	30.71	1.46	26.34-35.08	4.74	19	0.57	0.03	0.49-0.65	4.91															
PG: Siemens Dimension/AR/ES/RxL/Xpand/EXL	11	27.88	1.78	22.55-33.21	6.37	11	0.45	0.06	0.27-0.63	13.33															
PG: Qualigen FastPack PSA	7	42.90	6.20	24.31-61.49	14.45	7	0.43	0.05	0.29-0.57	10.70															
PG: TOSOH AIA	13	22.47	1.04	19.34-25.60	4.64	12	0.49	0.03	0.40-0.58	6.12															
SG: Select Beckman Coulter Instruments (Chem)	25	30.41	1.66	25.44-35.38	5.45	24	0.56	0.03	0.47-0.65	5.54															
SG: Select Siemens Instruments	13	27.19	2.65	19.23-35.15	9.76	13	0.45	0.07	0.25-0.65	14.67															
SG: Select OCD Vitros Instruments	5	15.18	1.00	12.19-18.17	6.57	5	0.51	0.03	0.43-0.59	5.49															
PG: Roche Cobas e Series	5	29.82	2.14	23.40-36.24	7.18	5	0.59	0.03	0.50-0.68	4.92															
SG: Select Abbott Architect Instruments	11	24.13	1.76	18.86-29.40	7.28	11	0.38	0.03	0.29-0.47	7.63															
PG: NanoEnTek Frend System	9	25.00	0.00	25.00-25.00	0.00	9	0.44	0.06	0.26-0.62	13.86															
Testosterone	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%
PG: Beckman ACCESS	13	13.43	0.61	11.60-15.20	4.55	13	1.10	0.07	0.90-1.30	6.36															
PG: TOSOH AIA	7	17.64	1.62	12.70-22.50	9.20	7	1.22	0.18	0.70-1.70	15.00															
SG: Select Beckman Coulter Instruments (Chem)	18	13.09	0.79	10.70-15.50	6.03	18	1.09	0.06	0.90-1.30	5.45															
SG: Select Abbott Architect Instruments	6	11.73	2.53	4.10-19.30	21.62	6	1.02	0.05	0.90-1.10	5.00															
AG: All Group	43	14.04	2.73	5.80-22.20	19.50	39	1.18	0.19	0.60-1.80	15.83															
685 - Special Chemistry II	SPC-1					SPC-2																			
Insulin	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%
SG: Select Siemens Immulite Instruments	5	35.60	4.81	21.20-50.00	13.51	5	3.35	0.92	0.60-6.20	27.06															
AG: All Group	5	100.60	22.11	34.30-166.90	21.98	5	7.33	0.73	5.10-9.50	10.00															
PTH	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%
SG: Select Beckman Coulter Instruments (Chem)	5	115.50	18.37	60.40-170.60	15.90	5	15.40	2.62	7.50-23.30	17.01															
AG: All Group	8	124.95	16.11	76.70-173.30	12.89	9	17.10	4.17	4.60-29.60	24.39															
Vitamin D	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%
PG: Beckman ACCESS	17	93.89	9.83	64.39-123.39	10.47	16	14.82	2.79	6.44-23.20	18.86															
PG: Horiba Pentra 400	5	13.70	0.00	13.70-13.70	0.00	5	1.40	0.00	1.40-1.40	0.00															
PG: TOSOH AIA	6	114.13	7.46	91.73-136.53	6.54	6	30.80	1.42	26.55-35.05	4.60															
SG: Select Beckman Coulter Instruments (Chem)	21	97.75	12.14	61.32-134.18	12.42	19	14.96	2.75	6.71-23.21	18.39															
PG: Roche Cobas e Series	5	60.00	0.00	60.00-60.00	0.00	5	11.42	1.44	7.11-15.73	12.58															
SG: Select Abbott Architect Instruments	14	92.31	4.50	78.82-105.80	4.87	14	13.50	0.70	11.38-15.62	5.22															
PG: NanoEnTek Frend System	5	96.00	0.00	96.00-96.00	0.00	5	19.46	3.43	9.18-29.74	17.61															
PG: Intl OCD Vitros ECI/ECIQ	5	60.90	68.02	0.00-264.97	111.70	5	19.00	0.28	18.15-19.85	1.49															
AG: All Group	56	91.79	17.52	39.24-144.34	19.08	51	15.64	4.17	3.12-28.16	26.68															
752 - Auto Diff 8 (Sysmex 5-part Diff)	SXT-1					SXT-2					SXT-3					SXT-4					SXT-5				
White Blood Cell Count	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%
PG: Sysmex XS-1000i	14	3.82	0.08	3.20-4.40	2.11	14	21.21	0.30	18.00-24.40	1.42	14	7.65	0.17	6.50-8.90	2.21	14	21.17	0.							



PG: Sysmex XS-1000i	14	2.41	0.02	2.27-2.55	0.75	14	5.59	0.08	5.25-5.93	1.48	14	4.69	0.04	4.41-4.97	0.85	14	5.60	0.05	5.26-5.94	0.84	14	4.69	0.11	4.41-4.97	2.24
SG: Select Sysmex Instruments (5-part Diff)	21	2.41	0.03	2.27-2.55	1.33	21	5.61	0.08	5.27-5.95	1.50	21	4.69	0.05	4.41-4.97	1.13	21	5.62	0.06	5.28-5.96	1.05	21	4.69	0.09	4.41-4.97	1.94
PG: Sysmex XN Series	6	2.39	0.03	2.25-2.53	1.21	6	5.65	0.07	5.31-5.99	1.15	6	4.69	0.06	4.41-4.97	1.30	6	5.64	0.06	5.30-5.98	1.13	6	4.67	0.05	4.39-4.95	1.07
Hemoglobin	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%
PG: Sysmex XS-1000i	14	5.81	0.08	5.40-6.20	1.38	14	17.52	0.11	16.30-18.70	0.63	14	13.26	0.13	12.40-14.20	0.98	14	17.51	0.15	16.30-18.70	0.86	14	13.29	0.28	12.40-14.20	2.11
SG: Select Sysmex Instruments (5-part Diff)	21	5.86	0.10	5.50-6.30	1.69	21	17.49	0.13	16.30-18.70	0.74	21	13.27	0.12	12.40-14.20	0.90	21	17.49	0.14	16.30-18.70	0.80	21	13.28	0.24	12.40-14.20	1.80
PG: Sysmex XN Series	6	5.97	0.05	5.60-6.40	0.83	6	17.42	0.15	16.20-18.60	0.86	6	13.32	0.10	12.40-14.20	0.75	6	17.48	0.13	16.30-18.70	0.74	6	13.30	0.09	12.40-14.20	0.68
Hematocrit	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%
PG: Sysmex XS-1000i	14	18.04	0.17	16.90-19.10	0.94	14	52.84	0.66	49.60-56.00	1.25	14	40.13	0.32	37.70-42.50	0.80	14	52.75	0.55	49.60-56.00	1.04	14	40.10	0.75	37.70-42.50	1.87
SG: Select Sysmex Instruments (5-part Diff)	21	17.97	0.35	16.90-19.10	1.94	21	52.71	0.87	49.50-55.90	1.65	21	40.07	0.58	37.70-42.50	1.45	21	52.63	0.81	49.40-55.80	1.54	21	40.00	0.75	37.60-42.40	1.88
PG: Sysmex XN Series	6	17.67	0.40	16.60-18.80	2.26	6	52.27	1.19	49.20-55.40	2.28	6	39.75	0.83	37.40-42.20	2.09	6	52.18	1.16	49.10-55.30	2.20	6	39.65	0.71	37.20-42.00	1.79
MCV	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%
PG: Sysmex XS-1000i	14	75.01	0.53	73.40-76.60	0.71	14	94.42	0.52	92.80-96.00	0.55	14	85.66	0.50	84.20-87.20	0.58	14	94.29	0.59	92.50-96.10	0.63	14	85.49	0.59	83.70-87.30	0.69
SG: Select Sysmex Instruments (5-part Diff)	21	74.66	0.95	71.90-77.50	1.27	21	93.84	1.17	90.30-97.30	1.25	21	85.39	0.90	82.70-88.10	1.05	21	93.75	1.17	90.30-97.30	1.25	21	85.30	0.90	82.60-88.00	1.06
PG: Sysmex XN Series	6	73.80	1.28	70.00-77.60	1.73	6	92.52	1.30	88.60-96.40	1.41	6	84.75	1.37	80.70-88.90	1.62	6	92.47	1.32	88.50-96.50	1.43	6	84.85	1.40	80.60-89.00	1.65
Platelet Count	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%
PG: Sysmex XS-1000i	14	65.93	3.36	50.00-83.00	5.15	14	355.07	8.42	266.00-444.00	2.37	14	210.14	4.52	158.00-263.00	2.14	14	354.43	8.60	266.00-443.00	2.43	14	208.79	7.80	157.00-261.00	3.73
SG: Select Sysmex Instruments (5-part Diff)	21	63.86	4.55	48.00-80.00	7.19	21	363.67	14.67	273.00-455.00	4.04	21	212.19	5.80	159.00-265.00	2.74	21	361.05	12.23	271.00-451.00	3.38	21	210.52	7.00	158.00-264.00	3.32
PG: Sysmex XN Series	6	58.67	2.88	44.00-74.00	4.92	6	380.83	7.47	286.00-476.00	1.97	6	215.00	5.69	161.00-269.00	2.65	6	373.67	5.92	281.00-468.00	1.58	6	213.17	2.64	160.00-266.00	1.22
RDW	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%
PG: Sysmex XS-1000i	14	16.81	0.19	16.20-17.40	1.13	14	14.51	0.20	13.90-15.10	1.38	14	15.22	0.10	14.90-15.50	0.66	14	14.51	0.24	13.80-15.20	1.66	14	15.22	0.10	14.90-15.50	0.66
SG: Select Sysmex Instruments (5-part Diff)	20	16.91	0.28	16.10-17.70	1.66	20	14.40	0.27	13.60-15.20	1.81	20	15.18	0.13	14.80-15.60	0.86	20	14.38	0.30	13.50-15.30	2.08	20	15.20	0.12	14.90-15.50	0.72
PG: Sysmex XN Series	5	17.26	0.11	17.00-17.60	0.64	5	14.06	0.11	13.80-14.40	0.78	5	15.06	0.17	14.60-15.60	1.13	5	14.02	0.13	13.60-14.40	0.93	5	15.16	0.13	14.80-15.60	0.86
613 - Auto Diff I (Sysmex 3-part)	SYX-1					SYX-2					SYX-3					SYX-4					SYX-5				
White Blood Cell Count	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%
PG: Sysmex K-800,1000,4500,KX-21	15	2.62	0.09	2.20-3.00	3.46	15	19.91	0.55	16.90-22.90	2.76	15	7.91	0.23	6.70-9.10	2.91	15	19.88	0.53	16.90-22.90	2.66	15	7.93	0.26	6.70-9.10	3.29
PG: Sysmex pocH-100i	24	2.63	0.10	2.20-3.00	3.85	24	19.93	0.48	16.90-22.90	2.41	24	7.93	0.22	6.70-9.10	2.78	24	20.06	0.48	17.10-23.10	2.39	24	7.93	0.20	6.70-9.10	2.53
SG: Select Sysmex Instruments (3-part Diff)	120	2.69	0.19	2.30-3.10	7.04	120	20.58	0.79	17.50-23.70	3.83	120	8.20	0.31	7.00-9.40	3.78	120	20.54	0.65	17.40-23.60	3.17	119	8.20	0.32	7.00-9.40	3.90
PG: Sysmex XP-300	81	2.72	0.22	2.30-3.10	8.15	81	20.90	0.71	17.80-24.00	3.40	81	8.33	0.26	7.10-9.50	3.13	81	20.80	0.55	17.70-23.90	2.64	80	8.34	0.27	7.10-9.50	3.25
SG: International Sysmex Instruments	28	2.63	0.09	2.20-3.00	3.46	28	19.90	0.47	16.90-22.90	2.36	28	7.93	0.22	6.70-9.10	2.78	28	20.05	0.46	17.00-23.00	2.30	28	7.93	0.20	6.70-9.10	2.53
Neutrophil	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%
PG: Sysmex K-800,1000,4500,KX-21	15	70.61	1.32	66.60-74.60	1.87																				

Print Date: 4/18/2019



SG: Select Medisense Instruments	32	100.16	9.29	72.00-128.00	9.30	5	45.25	2.06	39.00-51.00	4.67															
SG: Select Roche WBG Meters	289	113.21	6.24	94.00-132.00	5.49	49	58.45	2.68	50.00-66.00	4.66															
SG: Select Bayer WBG Meters	166	91.45	17.22	39.00-143.00	18.90	22	42.45	6.75	22.00-62.00	15.95															
PG: Nova Biomedical StatStrip Xpress	8	100.63	4.90	86.00-116.00	4.85																				
SG: Select True Track WBG Instruments	7	269.00	20.94	206.00-332.00	7.77	5	138.00	0.00	138.00-138.00	0.00															
PG: Nova Biomedical StatStrip Glucose meter	5	101.00	4.69	87.00-115.00	4.65																				
SG: Select Nova Biomedical StatStrip Instruments	12	100.75	4.62	87.00-115.00	4.55																				
PG: Bayer Contour Next	18	90.06	10.18	59.00-121.00	11.33	6	44.17	9.85	15.00-73.00	22.27															
PG: Roche AccuChek Performa	153	114.29	3.19	104.00-124.00	2.81	8	57.50	2.00	52.00-64.00	3.45															
PG: LifeScan OneTouch Verio	5	104.50	5.51	89.00-122.00	5.24																				
SG: Select TrueResult Meters	6	96.33	3.44	86.00-106.00	3.54																				
PG: True Metrix Blood Glucose Meter	8	92.75	11.97	57.00-129.00	12.90																				
PG: McKesson True Metrix Pro Blood Glucose Meter	31	97.90	9.40	70.00-126.00	9.59																				
SG: Select Arkray Whole Blood Glucose	5	120.60	10.97	88.00-154.00	9.09																				
PG: EvenCare G3 Blood Glucose System	5	105.00	6.21	86.00-124.00	5.90																				
SG: Select EvenCare Whole Blood Glucose	7	107.86	7.60	85.00-131.00	7.04																				
AG: All Group	631	105.77	12.54	69.00-144.00	11.79	95	51.60	8.52	27.00-78.00	16.35	20	84.15	11.20	50.00-118.00	13.33	20	189.75	21.68	125.00-255.00	11.42	20	46.95	7.57	24.00-70.00	16.17



789 - Group B Screen	Labs	Response	Number	Percent
		Strep B Ag Detect		
AD-1		* Denotes Correct Response		
AG: All Methods	1			
		*Positive	1	100.00%
Cepheid GeneXpert	1	*Positive	1	100.00%
AD-2		* Denotes Correct Response		
AG: All Methods	1			
		*Positive	1	100.00%
Cepheid GeneXpert	1	*Positive	1	100.00%
AD-3		* Denotes Correct Response		
AG: All Methods	1			
		*Negative	1	100.00%
Cepheid GeneXpert	1	*Negative	1	100.00%
AD-4		* Denotes Correct Response		
AG: All Methods	1			
		*Negative	1	100.00%
Cepheid GeneXpert	1	*Negative	1	100.00%
AD-5		* Denotes Correct Response		
AG: All Methods	1			
		*Positive	1	100.00%
Cepheid GeneXpert	1	*Positive	1	100.00%
760 - ANA (ELISA)	Labs	Response	Number	Percent
		ANA (ELISA) (International)		
AE-1		* Denotes Correct Response		
AG: All Methods	1			
		*Positive	1	100.00%
Dynex DS2	1	*Positive	1	100.00%
AE-2		* Denotes Correct Response		
AG: All Methods	1			
		*Positive	1	100.00%
Dynex DS2	1	*Positive	1	100.00%
AE-3		* Denotes Correct Response		
AG: All Methods	1			
		*Negative	1	100.00%
Dynex DS2	1	*Negative	1	100.00%
AE-4		* Denotes Correct Response		
AG: All Methods	1			
		*Negative	1	100.00%
Dynex DS2	1	*Negative	1	100.00%
AE-5		* Denotes Correct Response		
AG: All Methods	1			
		*Negative	1	100.00%
Dynex DS2	1	*Negative	1	100.00%
765 - Allergen Testing	Labs	Response	Number	Percent
		Ash - White Alrgn		
AL-1		* Denotes Correct Response		
PG: Hitachi CLA-1	3			
		*0	1	33.33%
		*1	1	33.33%
		*1/0	1	33.33%
AL-2		* Denotes Correct Response		
PG: Hitachi CLA-1	3			
		*0	3	100.00%
		Birch/Alder Mix Alrgn		
AL-1		* Denotes Correct Response		
PG: Hitachi CLA-1	4			
		*1	3	75.00%
		*2	1	25.00%
AL-2		* Denotes Correct Response		
PG: Hitachi CLA-1	4			
		0	2	50.00%
		1	1	25.00%
		2	1	25.00%
		Cedar - Mountain Alrgn		
AL-1		* Denotes Correct Response		
PG: Hitachi CLA-1	4			



		0	2	50.00%
		1/0	1	25.00%
		2	1	25.00%
AL-2		* Denotes Correct Response		
PG: Hitachi CLA-1	4			
		*0	3	75.00%
		*1/0	1	25.00%
		Cottonwood - Eastern Alrgn		
AL-1		* Denotes Correct Response		
PG: Hitachi CLA-1	4			
		*0	2	50.00%
		*1	1	25.00%
		*1/0	1	25.00%
AL-2		* Denotes Correct Response		
PG: Hitachi CLA-1	4			
		*0	4	100.00%
		Elm - White Alrgn		
AL-1		* Denotes Correct Response		
PG: Hitachi CLA-1	4			
		*0	2	50.00%
		*1	1	25.00%
		*1/0	1	25.00%
AL-2		* Denotes Correct Response		
PG: Hitachi CLA-1	4			
		*0	3	75.00%
		*1/0	1	25.00%
		Mulberry Mix Alrgn		
AL-1		* Denotes Correct Response		
PG: Hitachi CLA-1	3			
		*0	3	100.00%
AL-2		* Denotes Correct Response		
PG: Hitachi CLA-1	3			
		*0	3	100.00%
		Pine Mix Alrgn		
AL-1		* Denotes Correct Response		
PG: Hitachi CLA-1	2			
		*0	2	100.00%
AL-2		* Denotes Correct Response		
PG: Hitachi CLA-1	2			
		*0	2	100.00%
		Privet Alrgn		
AL-1		* Denotes Correct Response		
PG: Hitachi CLA-1	2			
		*0	1	50.00%
		*1/0	1	50.00%
AL-2		* Denotes Correct Response		
PG: Hitachi CLA-1	2			
		*0	2	100.00%
		Sycamore - American Alrgn		
AL-1		* Denotes Correct Response		
PG: Hitachi CLA-1	2			
		*0	1	50.00%
		*1	1	50.00%
AL-2		* Denotes Correct Response		
PG: Hitachi CLA-1	2			
		*0	1	50.00%
		*1/0	1	50.00%
		Cocklebur Alrgn		
AL-1		* Denotes Correct Response		
PG: Hitachi CLA-1	2			
		*0	2	100.00%
AL-2		* Denotes Correct Response		
PG: Hitachi CLA-1	2			
		*0	2	100.00%
		Lambs Quarters Alrgn		
AL-1		* Denotes Correct Response		
PG: Hitachi CLA-1	2			



		*1	1	50.00%
		*1/0	1	50.00%
AL-2		* Denotes Correct Response		
PG: Hitachi CLA-1	2			
		*0	2	100.00%
		Marshelder - Rough Alrgn		
AL-1		* Denotes Correct Response		
PG: Hitachi CLA-1	2			
		*0	1	50.00%
		*1/0	1	50.00%
AL-2		* Denotes Correct Response		
PG: Hitachi CLA-1	2			
		*0	2	100.00%
		Pigweed Alrgn		
AL-1		* Denotes Correct Response		
PG: Hitachi CLA-1	4			
		0	2	50.00%
		1	1	25.00%
		2	1	25.00%
AL-2		* Denotes Correct Response		
PG: Hitachi CLA-1	4			
		*0	3	75.00%
		*1	1	25.00%
		Ragweed - Short Alrgn		
AL-1		* Denotes Correct Response		
PG: Hitachi CLA-1	4			
		*1	1	25.00%
		*2	3	75.00%
AL-2		* Denotes Correct Response		
PG: Hitachi CLA-1	4			
		*1	1	25.00%
		*2	2	50.00%
		0	1	25.00%
		Sheep Sorrel Alrgn		
AL-1		* Denotes Correct Response		
PG: Hitachi CLA-1	4			
		*0	1	25.00%
		*1	3	75.00%
AL-2		* Denotes Correct Response		
PG: Hitachi CLA-1	4			
		0	2	50.00%
		1	1	25.00%
		2	1	25.00%
		Waterhemp Alrgn		
AL-1		* Denotes Correct Response		
PG: Hitachi CLA-1	2			
		*2	1	50.00%
		*3	1	50.00%
AL-2		* Denotes Correct Response		
PG: Hitachi CLA-1	2			
		0	1	50.00%
		3	1	50.00%
		Mite D Farinae Alrgn		
AL-1		* Denotes Correct Response		
PG: Hitachi CLA-1	5			
		*4	5	100.00%
AL-2		* Denotes Correct Response		
PG: Hitachi CLA-1	5			
		*4	4	80.00%
		0	1	20.00%
		Cladosporium Alrgn		
AL-1		* Denotes Correct Response		
PG: Hitachi CLA-1	4			
		*0	2	50.00%
		*1	1	25.00%
		*1/0	1	25.00%
AL-2		* Denotes Correct Response		



PG: Hitachi CLA-1	4			
		*0	3	75.00%
		*1	1	25.00%
		Tomato Alrgn		
AL-1		* Denotes Correct Response		
PG: Hitachi CLA-1	3			
		*0	2	66.67%
		*1	1	33.33%
AL-2		* Denotes Correct Response		
PG: Hitachi CLA-1	3			
		0	1	33.33%
		1	1	33.33%
		2	1	33.33%
		Wheat Alrgn		
AL-1		* Denotes Correct Response		
PG: Hitachi CLA-1	4			
		0	1	25.00%
		1	1	25.00%
		2	1	25.00%
		3	1	25.00%
AL-2		* Denotes Correct Response		
PG: Hitachi CLA-1	4			
		*1	1	25.00%
		*2	1	25.00%
		*3	2	50.00%
		Vegetable Mix Alrgn		
AL-1		* Denotes Correct Response		
PG: Hitachi CLA-1	2			
		*0	1	50.00%
		*1	1	50.00%
AL-2		* Denotes Correct Response		
PG: Hitachi CLA-1	2			
		*1	2	100.00%
		Tuna Alrgn		
AL-1		* Denotes Correct Response		
PG: Hitachi CLA-1	2			
		*0	1	50.00%
		*1/0	1	50.00%
AL-2		* Denotes Correct Response		
PG: Hitachi CLA-1	2			
		*1	1	50.00%
		*1/0	1	50.00%
		Shellfish Mix Alrgn		
AL-1		* Denotes Correct Response		
PG: Hitachi CLA-1	3			
		0	1	33.33%
		4	2	66.67%
AL-2		* Denotes Correct Response		
PG: Hitachi CLA-1	3			
		*4	3	100.00%
		Rice Alrgn		
AL-1		* Denotes Correct Response		
PG: Hitachi CLA-1	4			
		0	1	25.00%
		1	1	25.00%
		2	2	50.00%
AL-2		* Denotes Correct Response		
PG: Hitachi CLA-1	4			
		*1	2	50.00%
		*2	2	50.00%
		Potato Alrgn		
AL-1		* Denotes Correct Response		
PG: Hitachi CLA-1	2			
		0	1	50.00%
		2	1	50.00%
AL-2		* Denotes Correct Response		
PG: Hitachi CLA-1	2			



		*2	2	100.00%
		Orange Alrgn		
AL-1		* Denotes Correct Response		
PG: Hitachi CLA-1	3			
		0	1	33.33%
		1	1	33.33%
		2	1	33.33%
AL-2		* Denotes Correct Response		
PG: Hitachi CLA-1	3			
		*2	2	66.67%
		*3	1	33.33%
		Onion Mix Alrgn		
AL-1		* Denotes Correct Response		
PG: Hitachi CLA-1	3			
		*0	1	33.33%
		*1	2	66.67%
AL-2		* Denotes Correct Response		
PG: Hitachi CLA-1	3			
		0	1	33.33%
		1	1	33.33%
		2	1	33.33%
		Oat Alrgn		
AL-1		* Denotes Correct Response		
PG: Hitachi CLA-1	2			
		0	1	50.00%
		2	1	50.00%
AL-2		* Denotes Correct Response		
PG: Hitachi CLA-1	2			
		*2	2	100.00%
		Milk Alrgn		
AL-1		* Denotes Correct Response		
PG: Hitachi CLA-1	4			
		0	1	25.00%
		4	3	75.00%
AL-2		* Denotes Correct Response		
PG: Hitachi CLA-1	4			
		*4	4	100.00%
		Garlic Alrgn		
AL-1		* Denotes Correct Response		
PG: Hitachi CLA-1	2			
		0	1	50.00%
		2	1	50.00%
AL-2		* Denotes Correct Response		
PG: Hitachi CLA-1	2			
		*2	2	100.00%
		Egg - whole Alrgn		
AL-1		* Denotes Correct Response		
PG: Hitachi CLA-1	4			
		0	1	25.00%
		4	3	75.00%
AL-2		* Denotes Correct Response		
PG: Hitachi CLA-1	4			
		*3	1	25.00%
		*4	3	75.00%
		Corn Alrgn		
AL-1		* Denotes Correct Response		
PG: Hitachi CLA-1	3			
		*0	1	33.33%
		*1	1	33.33%
		*1/0	1	33.33%
AL-2		* Denotes Correct Response		
PG: Hitachi CLA-1	3			
		*1	2	66.67%
		*1/0	1	33.33%
		Chocolate Alrgn		
AL-1		* Denotes Correct Response		
PG: Hitachi CLA-1	2			



		*0	2	100.00%
AL-2		* Denotes Correct Response		
PG: Hitachi CLA-1	2			
		*0	2	100.00%
		Chicken Alrgn		
AL-1		* Denotes Correct Response		
PG: Hitachi CLA-1	3			
		*0	3	100.00%
AL-2		* Denotes Correct Response		
PG: Hitachi CLA-1	3			
		*0	2	66.67%
		*1/0	1	33.33%
		Beef Alrgn		
AL-1		* Denotes Correct Response		
PG: Hitachi CLA-1	2			
		*0	2	100.00%
AL-2		* Denotes Correct Response		
PG: Hitachi CLA-1	2			
		1/0	1	50.00%
		2	1	50.00%
		Bean - White Alrgn		
AL-1		* Denotes Correct Response		
PG: Hitachi CLA-1	2			
		*0	2	100.00%
AL-2		* Denotes Correct Response		
PG: Hitachi CLA-1	2			
		*1	2	100.00%
		Barley Alrgn		
AL-1		* Denotes Correct Response		
PG: Hitachi CLA-1	2			
		*0	1	50.00%
		*1	1	50.00%
AL-2		* Denotes Correct Response		
PG: Hitachi CLA-1	2			
		*0	1	50.00%
		*1/0	1	50.00%
		Apple Alrgn		
AL-1		* Denotes Correct Response		
PG: Hitachi CLA-1	2			
		*0	1	50.00%
		*1/0	1	50.00%
AL-2		* Denotes Correct Response		
PG: Hitachi CLA-1	2			
		*1	1	50.00%
		*1/0	1	50.00%
		Almond Alrgn		
AL-1		* Denotes Correct Response		
PG: Hitachi CLA-1	2			
		*0	1	50.00%
		*1/0	1	50.00%
AL-2		* Denotes Correct Response		
PG: Hitachi CLA-1	2			
		*1	1	50.00%
		*2	1	50.00%
		Alternaria alternata Alrgn		
AL-1		* Denotes Correct Response		
PG: Hitachi CLA-1	5			
		*1	2	40.00%
		*1/0	3	60.00%
AL-2		* Denotes Correct Response		
PG: Hitachi CLA-1	5			
		*0	2	40.00%
		*1	1	20.00%
		*1/0	2	40.00%
		Aspergillus fumigatus Alrgn		
AL-1		* Denotes Correct Response		
PG: Hitachi CLA-1	4			



		*0	3	75.00%
		*1/0	1	25.00%
AL-2		* Denotes Correct Response		
PG: Hitachi CLA-1	4			
		*0	3	75.00%
		*1/0	1	25.00%
		Bermuda Grass Alrgn		
AL-1		* Denotes Correct Response		
PG: Hitachi CLA-1	4			
		*1	1	25.00%
		*1/0	1	25.00%
		*2	2	50.00%
AL-2		* Denotes Correct Response		
PG: Hitachi CLA-1	5			
		*1	3	60.00%
		*2	1	20.00%
		0	1	20.00%
		Cat Epithelium Alrgn		
AL-1		* Denotes Correct Response		
PG: Hitachi CLA-1	2			
		*4	2	100.00%
AL-2		* Denotes Correct Response		
PG: Hitachi CLA-1	2			
		*4	2	100.00%
		Dog Epithelium Alrgn		
AL-1		* Denotes Correct Response		
PG: Hitachi CLA-1	2			
		*3	2	100.00%
AL-2		* Denotes Correct Response		
PG: Hitachi CLA-1	2			
		*3	2	100.00%
		Dog Dander Alrgn		
AL-1		* Denotes Correct Response		
PG: Hitachi CLA-1	3			
		*2	1	33.33%
		*3	1	33.33%
		*4	1	33.33%
AL-2		* Denotes Correct Response		
PG: Hitachi CLA-1	3			
		0	1	33.33%
		2	1	33.33%
		3	1	33.33%
		House Dust Mite Alrgn		
AL-1		* Denotes Correct Response		
PG: Hitachi CLA-1	4			
		*1	3	75.00%
		*1/0	1	25.00%
AL-2		* Denotes Correct Response		
PG: Hitachi CLA-1	5			
		*1	3	60.00%
		*2	1	20.00%
		0	1	20.00%
		Timothy Grass Alrgn		
AL-1		* Denotes Correct Response		
PG: Hitachi CLA-1	4			
		*4	4	100.00%
AL-2		* Denotes Correct Response		
PG: Hitachi CLA-1	4			
		0	1	25.00%
		4	3	75.00%
		Beech, American Alrgn		
AL-1		* Denotes Correct Response		
PG: Hitachi CLA-1	2			
		1/0	1	50.00%
		2	1	50.00%
AL-2		* Denotes Correct Response		
PG: Hitachi CLA-1	2			



		*0	1	50.00%
		*1	1	50.00%
		Box Elder, Maple Alrgn		
AL-1		* Denotes Correct Response		
PG: Hitachi CLA-1	3			
		*0	1	33.33%
		*1/0	2	66.67%
AL-2		* Denotes Correct Response		
PG: Hitachi CLA-1	3			
		*0	3	100.00%
		Peanut (Arachis hyogaea) Alrgn		
AL-1		* Denotes Correct Response		
PG: Hitachi CLA-1	4			
		0	1	25.00%
		4	3	75.00%
AL-2		* Denotes Correct Response		
PG: Hitachi CLA-1	4			
		*4	4	100.00%
		Cat Dander Alrgn		
AL-1		* Denotes Correct Response		
PG: Hitachi CLA-1	3			
		*4	3	100.00%
AL-2		* Denotes Correct Response		
PG: Hitachi CLA-1	3			
		0	1	33.33%
		4	2	66.67%
		Soybean Alrgn		
AL-1		* Denotes Correct Response		
PG: Hitachi CLA-1	4			
		0	1	25.00%
		3	3	75.00%
AL-2		* Denotes Correct Response		
PG: Hitachi CLA-1	4			
		*3	4	100.00%
		Oak, White Alrgn		
AL-1		* Denotes Correct Response		
PG: Hitachi CLA-1	2			
		*2	2	100.00%
AL-2		* Denotes Correct Response		
PG: Hitachi CLA-1	2			
		0	1	50.00%
		2	1	50.00%
		Penicillium Alrgn		
AL-1		* Denotes Correct Response		
PG: Hitachi CLA-1	3			
		*0	2	66.67%
		*1	1	33.33%
AL-2		* Denotes Correct Response		
PG: Hitachi CLA-1	3			
		*0	2	66.67%
		*1/0	1	33.33%
		Walnut/Hickory/Pecan Mix Alrgn		
AL-1		* Denotes Correct Response		
PG: Hitachi CLA-1	4			
		*0	2	50.00%
		*1	1	25.00%
		*1/0	1	25.00%
AL-2		* Denotes Correct Response		
PG: Hitachi CLA-1	4			
		*0	3	75.00%
		*1/0	1	25.00%
		Yeast - Bakers Alrgn		
AL-1		* Denotes Correct Response		
PG: Hitachi CLA-1	4			
		*0	3	75.00%
		*1	1	25.00%
AL-2		* Denotes Correct Response		



PG: Hitachi CLA-1	4			
		*0	3	75.00%
		*1	1	25.00%
683 - Antistreptolysin O (ASO)	Labs	Response	Number	Percent
		ASO (Qual) (International)		
AS-1		* Denotes Correct Response		
AG: All Methods	1			
		*Positive	1	100.00%
Wampole ColorCard ASO	1	*Positive	1	100.00%
AS-2		* Denotes Correct Response		
AG: All Methods	1			
		*Negative	1	100.00%
Wampole ColorCard ASO	1	*Negative	1	100.00%
AS-3		* Denotes Correct Response		
AG: All Methods	1			
		*Negative	1	100.00%
Wampole ColorCard ASO	1	*Negative	1	100.00%
AS-4		* Denotes Correct Response		
AG: All Methods	1			
		*Positive	1	100.00%
Wampole ColorCard ASO	1	*Positive	1	100.00%
AS-5		* Denotes Correct Response		
AG: All Methods	1			
		*Positive	1	100.00%
Wampole ColorCard ASO	1	*Positive	1	100.00%
782 - Complete Bacteriology	Labs	Response	Number	Percent
		Wound Culture		
BA-1		* Denotes Correct Response		
AG: All Methods	2			
		*Escherichia coli	2	100.00%
		*Staphylococcus epidermidis	1	50.00%
		No Secondary Organism Present	1	50.00%
Cultures	2	*Escherichia coli	2	100.00%
Cultures		*Staphylococcus epidermidis	1	50.00%
Cultures		No Secondary Organism Present	1	50.00%
975 - Bacterial Vaginosis	Labs	Response	Number	Percent
		Bacterial Vaginosa (Sialidase)		
BV-1		* Denotes Correct Response		
AG: All Methods	97			
		*Positive	83	85.57%
		Negative	14	14.43%
Sekisui Diagnostics OSOM BVBlue	97	*Positive	83	85.57%
Sekisui Diagnostics OSOM BVBlue		Negative	14	14.43%
BV-2		* Denotes Correct Response		
AG: All Methods	97			
		*Positive	84	86.60%
		Negative	13	13.40%
Sekisui Diagnostics OSOM BVBlue	97	*Positive	84	86.60%
Sekisui Diagnostics OSOM BVBlue		Negative	13	13.40%
661 - Urine Colony Count;784 - Presump ID/Colony Count - Urine;785 - Presump ID/Colony Count-Urine & Strep Combo;775 - Urine Culture;778 - Urine & Throat Cultures;779 - Urine & Genital Cultures;781 - Urine, Throat & Genital Cultures;782 - Complete Bacteri	Labs	Response	Number	Percent
		Colony Count		
CC-1		* Denotes Correct Response		
AG: All Methods	31			
		*10,000 - 100,000	13	41.94%
		*Greater than 100,000	16	51.61%
		Less than 10,000	1	3.23%
		No growth	1	3.23%
Calibrated Loop - Manual	10	*10,000 - 100,000	4	40.00%
Uri-Chek (Troy Bio)	6	*10,000 - 100,000	4	66.67%
Uricult (Orion)	13	*10,000 - 100,000	5	38.46%
Calibrated Loop - Manual	10	*Greater than 100,000	6	60.00%
Uri-Chek (Troy Bio)	6	*Greater than 100,000	2	33.33%
Uricult (Orion)	13	*Greater than 100,000	6	46.15%
Uri-Three (Culture Kits)	1	*Greater than 100,000	1	100.00%



Uricult (Lifesign)	1	*Greater than 100,000	1	100.00%
Uricult (Orion)	13	Less than 10,000	1	7.69%
Uricult (Orion)		No growth	1	7.69%
784 - Presump ID/Colony Count - Urine;785 - Presump ID/Colony Count-Urine & Strep Combo	Labs	Response	Number	Percent
		Colony Count		
CC-1		* Denotes Correct Response		
AG: All Methods	31			
		*10,000 - 100,000	13	41.94%
		*Greater than 100,000	16	51.61%
		Less than 10,000	1	3.23%
		No growth	1	3.23%
Calibrated Loop - Manual	10	*10,000 - 100,000	4	40.00%
Uri-Chek (Troy Bio)	6	*10,000 - 100,000	4	66.67%
Uricult (Orion)	13	*10,000 - 100,000	5	38.46%
Calibrated Loop - Manual	10	*Greater than 100,000	6	60.00%
Uri-Chek (Troy Bio)	6	*Greater than 100,000	2	33.33%
Uricult (Orion)	13	*Greater than 100,000	6	46.15%
Uri-Three (Culture Kits	1	*Greater than 100,000	1	100.00%
Uricult (Lifesign)	1	*Greater than 100,000	1	100.00%
Uricult (Orion)	13	Less than 10,000	1	7.69%
Uricult (Orion)		No growth	1	7.69%
CC-2		* Denotes Correct Response		
AG: All Methods	5			
		*Greater than 100,000	5	100.00%
Calibrated Loop - Manual	4	*Greater than 100,000	4	100.00%
Uri-Chek (Troy Bio)	1	*Greater than 100,000	1	100.00%
CC-3		* Denotes Correct Response		
AG: All Methods	5			
		*Greater than 100,000	5	100.00%
Calibrated Loop - Manual	4	*Greater than 100,000	4	100.00%
Uri-Chek (Troy Bio)	1	*Greater than 100,000	1	100.00%
CC-4		* Denotes Correct Response		
AG: All Methods	5			
		*Greater than 100,000	5	100.00%
Calibrated Loop - Manual	4	*Greater than 100,000	4	100.00%
Uri-Chek (Troy Bio)	1	*Greater than 100,000	1	100.00%
CC-5		* Denotes Correct Response		
AG: All Methods	5			
		*Greater than 100,000	5	100.00%
Calibrated Loop - Manual	4	*Greater than 100,000	4	100.00%
Uri-Chek (Troy Bio)	1	*Greater than 100,000	1	100.00%
		Urine Presump Id		
CC-1		* Denotes Correct Response		
AG: All Methods	2			
		*Gram Negative Organsim	1	50.00%
		*Presumptive Pseudomonas sp.	1	50.00%
Calibrated Loop - Manual	2	*Gram Negative Organsim	1	50.00%
Calibrated Loop - Manual		*Presumptive Pseudomonas sp.	1	50.00%
CC-2		* Denotes Correct Response		
AG: All Methods	2			
		*Gram Positive Organism	1	50.00%
		*Yeast	1	50.00%
Calibrated Loop - Manual	2	*Gram Positive Organism	1	50.00%
Calibrated Loop - Manual		*Yeast	1	50.00%
CC-3		* Denotes Correct Response		
AG: All Methods	2			
		*Gram Negative Organsim	1	50.00%
		*Presumptive E. Coli	1	50.00%
Calibrated Loop - Manual	2	*Gram Negative Organsim	1	50.00%
Calibrated Loop - Manual		*Presumptive E. Coli	1	50.00%
CC-4		* Denotes Correct Response		
AG: All Methods	2			
		*Gram Negative Organsim	1	50.00%
		*Presumptive Enterobacter sp.	1	50.00%
Calibrated Loop - Manual	2	*Gram Negative Organsim	1	50.00%
Calibrated Loop - Manual		*Presumptive Enterobacter sp.	1	50.00%
CC-5		* Denotes Correct Response		



AG: All Methods	2			
		*Presumptive Proteus sp.	1	50.00%
		Gram Positive Organism	1	50.00%
Calibrated Loop - Manual	2	*Presumptive Proteus sp.	1	50.00%
Calibrated Loop - Manual		Gram Positive Organism	1	50.00%
		Susceptibility Cefepime		
CC-1		* Denotes Correct Response		
AG: All Methods	1			
		*Sensitive	1	100.00%
Susceptibility	1	*Sensitive	1	100.00%
		Susceptibility Levofloxacin		
CC-1		* Denotes Correct Response		
AG: All Methods	3			
		*Intermediate	2	66.67%
		*Resistant	1	33.33%
Susceptibility	3	*Intermediate	2	66.67%
Susceptibility		*Resistant	1	33.33%
		Susceptibility Piperacillin/Tazobactam		
CC-1		* Denotes Correct Response		
AG: All Methods	1			
		*Sensitive	1	100.00%
Susceptibility	1	*Sensitive	1	100.00%
		Susceptibility Tobramycin		
CC-1		* Denotes Correct Response		
AG: All Methods	2			
		*Sensitive	2	100.00%
Susceptibility	2	*Sensitive	2	100.00%
964 - C. difficile Toxin/Ag (2 specimens);963 - C. difficile Toxin/Ag (5 specimens)	Labs	Response	Number	Percent
		C. difficile Antigen Detection		
CDF-1		* Denotes Correct Response		
AG: All Methods	1			
		*Positive	1	100.00%
Cepheid GeneXpert - Xpert C. diff	1	*Positive	1	100.00%
CDF-2		* Denotes Correct Response		
AG: All Methods	1			
		*Negative	1	100.00%
Cepheid GeneXpert - Xpert C. diff	1	*Negative	1	100.00%
795 - Vaginal Pathogens II (2 specimens);794 - Vaginal Pathogens II (5 specimens);977 - Trichomonas Ag Det (Waved)	Labs	Response	Number	Percent
		T. vaginalis - Non Culture Tech		
CGT-1		* Denotes Correct Response		
AG: All Methods	2			
		*Positive	2	100.00%
Cepheid GeneXpert	1	*Positive	1	100.00%
BD Max	1	*Positive	1	100.00%
CGT-2		* Denotes Correct Response		
AG: All Methods	2			
		*Negative	2	100.00%
Cepheid GeneXpert	1	*Negative	1	100.00%
BD Max	1	*Negative	1	100.00%
		Chlamydia - Non Culture Tech		
CGT-1		* Denotes Correct Response		
AG: All Methods	18			
		*Negative	18	100.00%
Quidel QuickVue Chlamydia	2	*Negative	2	100.00%
Cepheid GeneXpert	14	*Negative	14	100.00%
BD Viper XTR	1	*Negative	1	100.00%
BD Max	1	*Negative	1	100.00%
CGT-2		* Denotes Correct Response		
AG: All Methods	18			
		*Negative	18	100.00%
Quidel QuickVue Chlamydia	2	*Negative	2	100.00%
Cepheid GeneXpert	14	*Negative	14	100.00%
BD Viper XTR	1	*Negative	1	100.00%
BD Max	1	*Negative	1	100.00%
		GC - Non Culture Tech		



CGT-1		* Denotes Correct Response		
AG: All Methods	16			
		*Positive	16	100.00%
Cepheid GeneXpert	14	*Positive	14	100.00%
BD Viper XTR	1	*Positive	1	100.00%
BD Max	1	*Positive	1	100.00%
CGT-2		* Denotes Correct Response		
AG: All Methods	16			
		*Negative	16	100.00%
Cepheid GeneXpert	14	*Negative	14	100.00%
BD Viper XTR	1	*Negative	1	100.00%
BD Max	1	*Negative	1	100.00%
794 - Vaginal Pathogens II (5 specimens)	Labs	Response	Number	Percent
		Chlamydia - Non Culture Tech		
CGT-1		* Denotes Correct Response		
AG: All Methods	18			
		*Negative	18	100.00%
Quidel QuickVue Chlamydia	2	*Negative	2	100.00%
Cepheid GeneXpert	14	*Negative	14	100.00%
BD Viper XTR	1	*Negative	1	100.00%
BD Max	1	*Negative	1	100.00%
CGT-2		* Denotes Correct Response		
AG: All Methods	18			
		*Negative	18	100.00%
Quidel QuickVue Chlamydia	2	*Negative	2	100.00%
Cepheid GeneXpert	14	*Negative	14	100.00%
BD Viper XTR	1	*Negative	1	100.00%
BD Max	1	*Negative	1	100.00%
CGT-3		* Denotes Correct Response		
AG: All Methods	18			
		*Positive	18	100.00%
Quidel QuickVue Chlamydia	2	*Positive	2	100.00%
Cepheid GeneXpert	14	*Positive	14	100.00%
BD Viper XTR	1	*Positive	1	100.00%
BD Max	1	*Positive	1	100.00%
CGT-4		* Denotes Correct Response		
AG: All Methods	18			
		*Positive	18	100.00%
Quidel QuickVue Chlamydia	2	*Positive	2	100.00%
Cepheid GeneXpert	14	*Positive	14	100.00%
BD Viper XTR	1	*Positive	1	100.00%
BD Max	1	*Positive	1	100.00%
CGT-5		* Denotes Correct Response		
AG: All Methods	18			
		*Positive	18	100.00%
Quidel QuickVue Chlamydia	2	*Positive	2	100.00%
Cepheid GeneXpert	14	*Positive	14	100.00%
BD Viper XTR	1	*Positive	1	100.00%
BD Max	1	*Positive	1	100.00%
		GC - Non Culture Tech		
CGT-1		* Denotes Correct Response		
AG: All Methods	16			
		*Positive	16	100.00%
Cepheid GeneXpert	14	*Positive	14	100.00%
BD Viper XTR	1	*Positive	1	100.00%
BD Max	1	*Positive	1	100.00%
CGT-2		* Denotes Correct Response		
AG: All Methods	16			
		*Negative	16	100.00%
Cepheid GeneXpert	14	*Negative	14	100.00%
BD Viper XTR	1	*Negative	1	100.00%
BD Max	1	*Negative	1	100.00%
CGT-3		* Denotes Correct Response		
AG: All Methods	15			
		*Negative	9	60.00%
		Positive	6	40.00%
Cepheid GeneXpert	13	*Negative	7	53.85%
BD Viper XTR	1	*Negative	1	100.00%



BD Max	1	*Negative	1	100.00%
Cepheid GeneXpert	13	Positive	6	46.15%
CGT-4		* Denotes Correct Response		
AG: All Methods	16			
		*Positive	16	100.00%
Cepheid GeneXpert	14	*Positive	14	100.00%
BD Viper XTR	1	*Positive	1	100.00%
BD Max	1	*Positive	1	100.00%
CGT-5		* Denotes Correct Response		
AG: All Methods	16			
		*Negative	16	100.00%
Cepheid GeneXpert	14	*Negative	14	100.00%
BD Viper XTR	1	*Negative	1	100.00%
BD Max	1	*Negative	1	100.00%
		T. vaginalis - Non Culture Tech		
CGT-1		* Denotes Correct Response		
AG: All Methods	2			
		*Positive	2	100.00%
Cepheid GeneXpert	1	*Positive	1	100.00%
BD Max	1	*Positive	1	100.00%
CGT-2		* Denotes Correct Response		
AG: All Methods	2			
		*Negative	2	100.00%
Cepheid GeneXpert	1	*Negative	1	100.00%
BD Max	1	*Negative	1	100.00%
CGT-3		* Denotes Correct Response		
AG: All Methods	2			
		*Negative	2	100.00%
Cepheid GeneXpert	1	*Negative	1	100.00%
BD Max	1	*Negative	1	100.00%
CGT-4		* Denotes Correct Response		
AG: All Methods	2			
		*Negative	2	100.00%
Cepheid GeneXpert	1	*Negative	1	100.00%
BD Max	1	*Negative	1	100.00%
CGT-5		* Denotes Correct Response		
AG: All Methods	2			
		*Positive	2	100.00%
Cepheid GeneXpert	1	*Positive	1	100.00%
BD Max	1	*Positive	1	100.00%
		Trichomonas Aq Det (Waived)		
CGT-1		* Denotes Correct Response		
AG: All Methods	96			
		*Positive	87	90.63%
		Negative	9	9.38%
Sekisui Diagnostics OSOM Trichomonas	96	*Positive	87	90.63%
Sekisui Diagnostics OSOM Trichomonas		Negative	9	9.38%
CGT-2		* Denotes Correct Response		
AG: All Methods	95			
		*Negative	85	89.47%
		Positive	10	10.53%
Sekisui Diagnostics OSOM Trichomonas	95	*Negative	85	89.47%
Sekisui Diagnostics OSOM Trichomonas		Positive	10	10.53%
639 - Complete Urinalysis Package;638 - Urine Sediment;641 - Clinical Microscopy;601 - Waived/PPM Tests II;609 - Waived/PPM Tests I	Labs	Response	Number	Percent
		Urine Sediment Id		
CM-1		* Denotes Correct Response		
AG: All Methods	634			
		*Waxy cast	467	73.66%
		Broad cast	20	3.15%
		Cellular (RTE) cast	6	0.95%
		Cellular (WBC) cast	3	0.47%
		Fatty (fat globules) cast	2	0.32%
		Fiber (exogenous)/fecal contamination	3	0.47%
		Granular cast	8	1.26%
		Hyaline cast	122	19.24%
		RBC/Blood/Hgb cast	1	0.16%



		Renal tubular epithelial (RTE)	1	0.16%
		RTE with fat globules (oval fat bodies)	1	0.16%
Microscopy	634	*Waxy cast	467	73.66%
Microscopy		Broad cast	20	3.15%
Microscopy		Cellular (RTE) cast	6	0.95%
Microscopy		Cellular (WBC) cast	3	0.47%
Microscopy		Fatty (fat globules) cast	2	0.32%
Microscopy		Fiber (exogenous)/fecal contamination	3	0.47%
Microscopy		Granular cast	8	1.26%
Microscopy		Hyaline cast	122	19.24%
Microscopy		RBC/Blood/Hgb cast	1	0.16%
Microscopy		Renal tubular epithelial (RTE)	1	0.16%
Microscopy		RTE with fat globules (oval fat bodies)	1	0.16%
CM-2		* Denotes Correct Response		
AG: All Methods	633			
		*Renal tubular epithelial (RTE)	486	76.78%
		Cellular (RTE) cast	2	0.32%
		RBC/Blood/Hgb cast	1	0.16%
		Red blood cell (RBC) (Erythrocyte)	6	0.95%
		RTE with fat globules (oval fat bodies)	7	1.11%
		Transitional epithelial cell (urothelial)	88	13.90%
		White blood cell (WBC) (Leukocyte)	43	6.79%
Microscopy	633	*Renal tubular epithelial (RTE)	486	76.78%
Microscopy		Cellular (RTE) cast	2	0.32%
Microscopy		RBC/Blood/Hgb cast	1	0.16%
Microscopy		Red blood cell (RBC) (Erythrocyte)	6	0.95%
Microscopy		RTE with fat globules (oval fat bodies)	7	1.11%
Microscopy		Transitional epithelial cell (urothelial)	88	13.90%
Microscopy		White blood cell (WBC) (Leukocyte)	43	6.79%
CM-3		* Denotes Correct Response		
AG: All Methods	633			
		*Ident Unknown - Referred	1	0.16%
		*Squamous epithelial cell	625	98.74%
		Calcium oxalate crystal	1	0.16%
		Renal tubular epithelial (RTE)	1	0.16%
		RTE with fat globules (oval fat bodies)	1	0.16%
		Transitional epithelial cell (urothelial)	4	0.63%
Microscopy	633	*Ident Unknown - Referred	1	0.16%
Microscopy		*Squamous epithelial cell	625	98.74%
Microscopy		Calcium oxalate crystal	1	0.16%
Microscopy		Renal tubular epithelial (RTE)	1	0.16%
Microscopy		RTE with fat globules (oval fat bodies)	1	0.16%
Microscopy		Transitional epithelial cell (urothelial)	4	0.63%
CM-4		* Denotes Correct Response		
AG: All Methods	633			
		*Ident Unknown - Referred	1	0.16%
		*Starch granules	423	66.82%
		Ammonium biurate crystal	2	0.32%
		Ammonium magnesium phosphate (Triple phosphate)	2	0.32%
		Amorphous urates/phosphates	1	0.16%
		Calcium carbonate	6	0.95%
		Calcium oxalate crystal	153	24.17%
		Calcium phosphate crystal	2	0.32%
		Cholesterol crystal	2	0.32%
		Cystine crystal	4	0.63%
		Fat droplets or globules	11	1.74%
		Hippuric acid	2	0.32%
		Leucine crystal	4	0.63%
		Red blood cell (RBC) (Erythrocyte)	2	0.32%
		RTE with fat globules (oval fat bodies)	2	0.32%
		Squamous epithelial cell	1	0.16%
		Sulfonamide crystal	1	0.16%
		Uric acid crystal	13	2.05%
		Yeast/fungi	1	0.16%
Microscopy	633	*Ident Unknown - Referred	1	0.16%
Microscopy		*Starch granules	423	66.82%
Microscopy		Ammonium biurate crystal	2	0.32%
Microscopy		Ammonium magnesium phosphate (Triple phosphate)	2	0.32%
Microscopy		Amorphous urates/phosphates	1	0.16%



Microscopy		Calcium carbonate	6	0.95%
Microscopy		Calcium oxalate crystal	153	24.17%
Microscopy		Calcium phosphate crystal	2	0.32%
Microscopy		Cholesterol crystal	2	0.32%
Microscopy		Cystine crystal	4	0.63%
Microscopy		Fat droplets or globules	11	1.74%
Microscopy		Hippuric acid	2	0.32%
Microscopy		Leucine crystal	4	0.63%
Microscopy		Red blood cell (RBC) (Erythrocyte)	2	0.32%
Microscopy		RTE with fat globules (oval fat bodies)	2	0.32%
Microscopy		Squamous epithelial cell	1	0.16%
Microscopy		Sulfonamide crystal	1	0.16%
Microscopy		Uric acid crystal	13	2.05%
Microscopy		Yeast/fungi	1	0.16%
		Vaginal Wet Prep		
CM-5		* Denotes Correct Response		
AG: All Methods	875			
		*Yeast/fungi present	860	98.29%
		Bacteria present	1	0.11%
		Clue cell	2	0.23%
		Fiber (exogenous)/fecal contamination	1	0.11%
		Hyaline cast	1	0.11%
		Mucous strand	1	0.11%
		Spermatozoa present	1	0.11%
		Squamous epithelial cell	7	0.80%
		Starch granule	1	0.11%
Microscopy	875	*Yeast/fungi present	860	98.29%
Microscopy		Bacteria present	1	0.11%
Microscopy		Clue cell	2	0.23%
Microscopy		Fiber (exogenous)/fecal contamination	1	0.11%
Microscopy		Hyaline cast	1	0.11%
Microscopy		Mucous strand	1	0.11%
Microscopy		Spermatozoa present	1	0.11%
Microscopy		Squamous epithelial cell	7	0.80%
Microscopy		Starch granule	1	0.11%
CM-6		* Denotes Correct Response		
AG: All Methods	875			
		*Squamous epithelial cell	703	80.34%
		Bacteria absent	1	0.11%
		Bacteria present	5	0.57%
		Clue cell	150	17.14%
		Red Blood Cell (RBC) (Erythrocyte)	1	0.11%
		White Blood Cell (WBC) (Leukocyte)	7	0.80%
		Yeast/fungi present	8	0.91%
Microscopy	875	*Squamous epithelial cell	703	80.34%
Microscopy		Bacteria absent	1	0.11%
Microscopy		Bacteria present	5	0.57%
Microscopy		Clue cell	150	17.14%
Microscopy		Red Blood Cell (RBC) (Erythrocyte)	1	0.11%
Microscopy		White Blood Cell (WBC) (Leukocyte)	7	0.80%
Microscopy		Yeast/fungi present	8	0.91%
		Skin Scraping		
CM-7		* Denotes Correct Response		
AG: All Methods	589			
		*Scabies present	578	98.13%
		*Yeast/fungi absent	10	1.70%
		No scabies present	1	0.17%
Microscopy	589	*Scabies present	578	98.13%
Microscopy		*Yeast/fungi absent	10	1.70%
Microscopy		No scabies present	1	0.17%
		Pinworm Prep		
CM-8		* Denotes Correct Response		
AG: All Methods	276			
		*Pinworm/pinworm eggs present	274	99.28%
		No pinworm/pinworm eggs present	2	0.72%
Microscopy	276	*Pinworm/pinworm eggs present	274	99.28%
Microscopy		No pinworm/pinworm eggs present	2	0.72%
683 - Antistreptolysin O (ASO);684 - C-Reactive Protein	Labs	Response	Number	Percent



		CRP Qualitative		
CR-1		* Denotes Correct Response		
AG: All Methods	1			
		*Positive	1	100.00%
Inverness Immunex	1	*Positive	1	100.00%
684 - C-Reactive Protein	Labs	Response	Number	Percent
CR-1		CRP Qualitative		
		* Denotes Correct Response		
AG: All Methods	1			
		*Positive	1	100.00%
Inverness Immunex	1	*Positive	1	100.00%
CR-2		* Denotes Correct Response		
AG: All Methods	1			
		*Positive	1	100.00%
Inverness Immunex	1	*Positive	1	100.00%
772 - Dermatophyte Screen	Labs	Response	Number	Percent
		Dermatophyte Screen		
DM-1		* Denotes Correct Response		
AG: All Methods	13			
		*Positive, Dermatophytes Present	7	53.85%
		Negative, No Dermatophytes Present	6	46.15%
Dermatophytes	13	*Positive, Dermatophytes Present	7	53.85%
Dermatophytes		Negative, No Dermatophytes Present	6	46.15%
DM-2		* Denotes Correct Response		
AG: All Methods	13			
		*Negative, No Dermatophytes Present	13	100.00%
Dermatophytes	13	*Negative, No Dermatophytes Present	13	100.00%
DM-3		* Denotes Correct Response		
AG: All Methods	13			
		*Negative, No Dermatophytes Present	13	100.00%
Dermatophytes	13	*Negative, No Dermatophytes Present	13	100.00%
DM-4		* Denotes Correct Response		
AG: All Methods	13			
		*Positive, Dermatophytes Present	13	100.00%
Dermatophytes	13	*Positive, Dermatophytes Present	13	100.00%
DM-5		* Denotes Correct Response		
AG: All Methods	13			
		*Positive, Dermatophytes Present	13	100.00%
Dermatophytes	13	*Positive, Dermatophytes Present	13	100.00%
777 - Genital Culture:796 - Genital Culture (2 specimens)	Labs	Response	Number	Percent
		Genital Culture		
GC-1		* Denotes Correct Response		
AG: All Methods	2			
		*Presumptive positive for GC, GR(-) diplococci, oxidase positive organism	1	50.00%
		*Presumptive positive for GC, growth referred for ID/confirmation	1	50.00%
Cultures	2	*Presumptive positive for GC, GR(-) diplococci, oxidase positive organism	1	50.00%
Cultures		*Presumptive positive for GC, growth referred for ID/confirmation	1	50.00%
GC-2		* Denotes Correct Response		
AG: All Methods	2			
		*No growth, GC selective media	2	100.00%
Cultures	2	*No growth, GC selective media	2	100.00%
GC-3		* Denotes Correct Response		
AG: All Methods	2			
		*No growth, GC selective media	2	100.00%
Cultures	2	*No growth, GC selective media	2	100.00%
GC-4		* Denotes Correct Response		
AG: All Methods	4			
		*Group B Streptococcus not present	1	25.00%
		*Negative for GC	1	25.00%
		*No growth, GC selective media	2	50.00%
Cultures	4	*Group B Streptococcus not present	1	25.00%
Cultures		*Negative for GC	1	25.00%
Cultures		*No growth, GC selective media	2	50.00%
GC-5		* Denotes Correct Response		
AG: All Methods	3			



		*Positive for GC	1	33.33%
		*Presumptive positive for GC, GR(-) diplococci, oxidase positive organism	1	33.33%
		*Presumptive positive for GC, growth referred for ID/confirmation	1	33.33%
Cultures	3	*Positive for GC	1	33.33%
Cultures		*Presumptive positive for GC, GR(-) diplococci, oxidase positive organism	1	33.33%
Cultures		*Presumptive positive for GC, growth referred for ID/confirmation	1	33.33%
668 - Gram Stain	Labs	Response	Number	Percent
		Gram Stain		
GS-1		* Denotes Correct Response		
AG: All Methods	7			
		*Gram Negative	6	85.71%
		Gram Positive	1	14.29%
Gram Stain	7	*Gram Negative	6	85.71%
Gram Stain		Gram Positive	1	14.29%
GS-2		* Denotes Correct Response		
AG: All Methods	7			
		*Gram Positive	7	100.00%
Gram Stain	7	*Gram Positive	7	100.00%
GS-3		* Denotes Correct Response		
AG: All Methods	7			
		*Gram Negative	6	85.71%
		Gram Positive	1	14.29%
Gram Stain	7	*Gram Negative	6	85.71%
Gram Stain		Gram Positive	1	14.29%
GS-4		* Denotes Correct Response		
AG: All Methods	7			
		*Gram Negative	6	85.71%
		Gram Positive	1	14.29%
Gram Stain	7	*Gram Negative	6	85.71%
Gram Stain		Gram Positive	1	14.29%
GS-5		* Denotes Correct Response		
AG: All Methods	7			
		*Gram Positive	7	100.00%
Gram Stain	7	*Gram Positive	7	100.00%
		Gram Stain Morph		
GS-1		* Denotes Correct Response		
AG: All Methods	6			
		*Rods/Bacilli	6	100.00%
Gram Stain	6	*Rods/Bacilli	6	100.00%
GS-2		* Denotes Correct Response		
AG: All Methods	6			
		*Cocci	4	66.67%
		*Cocci in Chains	1	16.67%
		*Cocci in Pairs	1	16.67%
Gram Stain	6	*Cocci	4	66.67%
Gram Stain		*Cocci in Chains	1	16.67%
Gram Stain		*Cocci in Pairs	1	16.67%
GS-3		* Denotes Correct Response		
AG: All Methods	6			
		*Cocci	4	66.67%
		*Diplococci	2	33.33%
Gram Stain	6	*Cocci	4	66.67%
Gram Stain		*Diplococci	2	33.33%
GS-4		* Denotes Correct Response		
AG: All Methods	6			
		*Rods/Bacilli	6	100.00%
Gram Stain	6	*Rods/Bacilli	6	100.00%
GS-5		* Denotes Correct Response		
AG: All Methods	6			
		*Cocci	4	66.67%
		*Cocci in Pairs	2	33.33%
Gram Stain	6	*Cocci	4	66.67%
Gram Stain		*Cocci in Pairs	2	33.33%
978 - Gastrointestinal Viral Detection	Labs	Response	Number	Percent
		Rotavirus - Aq Detect		



GVD-1		* Denotes Correct Response		
AG: All Methods	8			
		*Negative	8	100.00%
Meridian Immunocard Stat! Rotavirus	4	*Negative	4	100.00%
SAS Rota-Test	4	*Negative	4	100.00%
GVD-2		* Denotes Correct Response		
AG: All Methods	8			
		*Positive	7	87.50%
		Negative	1	12.50%
Meridian Immunocard Stat! Rotavirus	4	*Positive	4	100.00%
SAS Rota-Test	4	*Positive	3	75.00%
SAS Rota-Test		Negative	1	25.00%
GVD-3		* Denotes Correct Response		
AG: All Methods	8			
		*Negative	8	100.00%
Meridian Immunocard Stat! Rotavirus	4	*Negative	4	100.00%
SAS Rota-Test	4	*Negative	4	100.00%
GVD-4		* Denotes Correct Response		
AG: All Methods	8			
		*Positive	8	100.00%
Meridian Immunocard Stat! Rotavirus	4	*Positive	4	100.00%
SAS Rota-Test	4	*Positive	4	100.00%
GVD-5		* Denotes Correct Response		
AG: All Methods	8			
		*Negative	8	100.00%
Meridian Immunocard Stat! Rotavirus	4	*Negative	4	100.00%
SAS Rota-Test	4	*Negative	4	100.00%
678 - Serum hCG	Labs	Response	Number	Percent
HCG-1		* Denotes Correct Response		
AG: All Methods	78			
		*Positive	78	100.00%
Beck Coulter ICON 25	7	*Positive	7	100.00%
Fisher HC Sure-Vue hCG STAT	2	*Positive	2	100.00%
Sekisui Diagnostics OSOM	2	*Positive	2	100.00%
Henry Schein One Step + (hCG)	5	*Positive	5	100.00%
Immunostics Immuno/hCG	1	*Positive	1	100.00%
Inverness Signify HCG	2	*Positive	2	100.00%
LifeSign Status hCG	1	*Positive	1	100.00%
McKesson Medi-Lab hCG - Urine	1	*Positive	1	100.00%
Quidel QuickVue + hCG	2	*Positive	2	100.00%
Quidel QuickVue One-Step hCG Combo (U/S)	27	*Positive	27	100.00%
Select Diag hCG Combo	3	*Positive	3	100.00%
SP Brand (Cardinal) HCG Combo	14	*Positive	14	100.00%
Wampole Clearview	1	*Positive	1	100.00%
Polymedco Poly stat hCG Serum	1	*Positive	1	100.00%
Stanbio True 20 Plus	1	*Positive	1	100.00%
Consult Diag hCG Combo Cassette	7	*Positive	7	100.00%
Alere hCG Combo Cassette (25 mIU/mL)	1	*Positive	1	100.00%
HCG-2		* Denotes Correct Response		
AG: All Methods	78			
		*Negative	78	100.00%
Beck Coulter ICON 25	7	*Negative	7	100.00%
Fisher HC Sure-Vue hCG STAT	2	*Negative	2	100.00%
Sekisui Diagnostics OSOM	2	*Negative	2	100.00%
Henry Schein One Step + (hCG)	5	*Negative	5	100.00%
Immunostics Immuno/hCG	1	*Negative	1	100.00%
Inverness Signify HCG	2	*Negative	2	100.00%
LifeSign Status hCG	1	*Negative	1	100.00%
McKesson Medi-Lab hCG - Urine	1	*Negative	1	100.00%
Quidel QuickVue + hCG	2	*Negative	2	100.00%
Quidel QuickVue One-Step hCG Combo (U/S)	27	*Negative	27	100.00%
Select Diag hCG Combo	3	*Negative	3	100.00%
SP Brand (Cardinal) HCG Combo	14	*Negative	14	100.00%
Wampole Clearview	1	*Negative	1	100.00%
Polymedco Poly stat hCG Serum	1	*Negative	1	100.00%
Stanbio True 20 Plus	1	*Negative	1	100.00%
Consult Diag hCG Combo Cassette	7	*Negative	7	100.00%



Alere hCG Combo Cassette (25 mIU/mL)	1	*Negative	1	100.00%
HCG-3		* Denotes Correct Response		
AG: All Methods	78			
		*Negative	78	100.00%
Beck Coulter ICON 25	7	*Negative	7	100.00%
Fisher HC Sure-Vue hCG STAT	2	*Negative	2	100.00%
Sekisui Diagnostics OSOM	2	*Negative	2	100.00%
Henry Schein One Step + (hCG)	5	*Negative	5	100.00%
Immunostics Immuno/hCG	1	*Negative	1	100.00%
Inverness Signify HCG	2	*Negative	2	100.00%
LifeSign Status hCG	1	*Negative	1	100.00%
McKesson Medi-Lab hCG - Urine	1	*Negative	1	100.00%
Quidel QuickVue + hCG	2	*Negative	2	100.00%
Quidel QuickVue One-Step hCG Combo (U/S)	27	*Negative	27	100.00%
Select Diag hCG Combo	3	*Negative	3	100.00%
SP Brand (Cardinal) HCG Combo	14	*Negative	14	100.00%
Wampole Clearview	1	*Negative	1	100.00%
Polymedco Poly stat hCG Serum	1	*Negative	1	100.00%
Stanbio True 20 Plus	1	*Negative	1	100.00%
Consult Diag hCG Combo Cassette	7	*Negative	7	100.00%
Alere hCG Combo Cassette (25 mIU/mL)	1	*Negative	1	100.00%
HCG-4		* Denotes Correct Response		
AG: All Methods	78			
		*Positive	78	100.00%
Beck Coulter ICON 25	7	*Positive	7	100.00%
Fisher HC Sure-Vue hCG STAT	2	*Positive	2	100.00%
Sekisui Diagnostics OSOM	2	*Positive	2	100.00%
Henry Schein One Step + (hCG)	5	*Positive	5	100.00%
Immunostics Immuno/hCG	1	*Positive	1	100.00%
Inverness Signify HCG	2	*Positive	2	100.00%
LifeSign Status hCG	1	*Positive	1	100.00%
McKesson Medi-Lab hCG - Urine	1	*Positive	1	100.00%
Quidel QuickVue + hCG	2	*Positive	2	100.00%
Quidel QuickVue One-Step hCG Combo (U/S)	27	*Positive	27	100.00%
Select Diag hCG Combo	3	*Positive	3	100.00%
SP Brand (Cardinal) HCG Combo	14	*Positive	14	100.00%
Wampole Clearview	1	*Positive	1	100.00%
Polymedco Poly stat hCG Serum	1	*Positive	1	100.00%
Stanbio True 20 Plus	1	*Positive	1	100.00%
Consult Diag hCG Combo Cassette	7	*Positive	7	100.00%
Alere hCG Combo Cassette (25 mIU/mL)	1	*Positive	1	100.00%
HCG-5		* Denotes Correct Response		
AG: All Methods	78			
		*Positive	78	100.00%
Beck Coulter ICON 25	7	*Positive	7	100.00%
Fisher HC Sure-Vue hCG STAT	2	*Positive	2	100.00%
Sekisui Diagnostics OSOM	2	*Positive	2	100.00%
Henry Schein One Step + (hCG)	5	*Positive	5	100.00%
Immunostics Immuno/hCG	1	*Positive	1	100.00%
Inverness Signify HCG	2	*Positive	2	100.00%
LifeSign Status hCG	1	*Positive	1	100.00%
McKesson Medi-Lab hCG - Urine	1	*Positive	1	100.00%
Quidel QuickVue + hCG	2	*Positive	2	100.00%
Quidel QuickVue One-Step hCG Combo (U/S)	27	*Positive	27	100.00%
Select Diag hCG Combo	3	*Positive	3	100.00%
SP Brand (Cardinal) HCG Combo	14	*Positive	14	100.00%
Wampole Clearview	1	*Positive	1	100.00%
Polymedco Poly stat hCG Serum	1	*Positive	1	100.00%
Stanbio True 20 Plus	1	*Positive	1	100.00%
Consult Diag hCG Combo Cassette	7	*Positive	7	100.00%
Alere hCG Combo Cassette (25 mIU/mL)	1	*Positive	1	100.00%
		HCG - Qual Serum (International)		
HCG-1		* Denotes Correct Response		
AG: All Methods	3			
		*Positive	3	100.00%
Henry Schein One Step + (hCG)	1	*Positive	1	100.00%
SAS Ultra HCG	2	*Positive	2	100.00%
HCG-2		* Denotes Correct Response		
AG: All Methods	3			



		*Negative	3	100.00%
Henry Schein One Step + (hCG)	1	*Negative	1	100.00%
SAS Ultra HCG	2	*Negative	2	100.00%
HCG-3		* Denotes Correct Response		
AG: All Methods	3			
		*Negative	3	100.00%
Henry Schein One Step + (hCG)	1	*Negative	1	100.00%
SAS Ultra HCG	2	*Negative	2	100.00%
HCG-4		* Denotes Correct Response		
AG: All Methods	3			
		*Positive	3	100.00%
Henry Schein One Step + (hCG)	1	*Positive	1	100.00%
SAS Ultra HCG	2	*Positive	2	100.00%
HCG-5		* Denotes Correct Response		
AG: All Methods	3			
		*Positive	3	100.00%
Henry Schein One Step + (hCG)	1	*Positive	1	100.00%
SAS Ultra HCG	2	*Positive	2	100.00%
933 - HIV Antibody/Antigen Combo (2 specimens):934 - HIV Antibody/Antigen Combo (5 specimens)	Labs	Response	Number	Percent
		HIV Antibody (Waived)		
HIV-1		* Denotes Correct Response		
AG: All Methods	98			
		*Positive	97	98.98%
		Negative	1	1.02%
Alere Clearview HIV 1/2 STAT PAK (WB)	45	*Positive	45	100.00%
Orasure OraQuick Anti-HIV	10	*Positive	9	90.00%
Trinity Biotech Uni-Gold Recombigen HIV	2	*Positive	2	100.00%
Alere Determine HIV Ag/Ab	30	*Positive	30	100.00%
Chembio HIV 1/2 STAT-PAK Dipstick	3	*Positive	3	100.00%
Chembio HIV 1/2 STAT-PAK Cassette	3	*Positive	3	100.00%
bioLytical INSTI HIV Antibody	5	*Positive	5	100.00%
Orasure OraQuick Anti-HIV	10	Negative	1	10.00%
HIV-2		* Denotes Correct Response		
AG: All Methods	98			
		*Negative	93	94.90%
		Positive	5	5.10%
Alere Clearview HIV 1/2 STAT PAK (WB)	45	*Negative	45	100.00%
Orasure OraQuick Anti-HIV	10	*Negative	9	90.00%
Trinity Biotech Uni-Gold Recombigen HIV	2	*Negative	2	100.00%
Alere Determine HIV Ag/Ab	30	*Negative	26	86.67%
Chembio HIV 1/2 STAT-PAK Dipstick	3	*Negative	3	100.00%
Chembio HIV 1/2 STAT-PAK Cassette	3	*Negative	3	100.00%
bioLytical INSTI HIV Antibody	5	*Negative	5	100.00%
Orasure OraQuick Anti-HIV	10	Positive	1	10.00%
Alere Determine HIV Ag/Ab	30	Positive	4	13.33%
		HIV Antibody (International)		
HIV-1		* Denotes Correct Response		
AG: All Methods	5			
		*Positive	5	100.00%
Alere Clearview HIV 1/2 STAT PAK (WB)	1	*Positive	1	100.00%
Intl OCD Vitros ECI/ECiQ	3	*Positive	3	100.00%
Intl Chembio HIV 1/2 STAT-PAK Dipstick	1	*Positive	1	100.00%
HIV-2		* Denotes Correct Response		
AG: All Methods	5			
		*Negative	5	100.00%
Alere Clearview HIV 1/2 STAT PAK (WB)	1	*Negative	1	100.00%
Intl OCD Vitros ECI/ECiQ	3	*Negative	3	100.00%
Intl Chembio HIV 1/2 STAT-PAK Dipstick	1	*Negative	1	100.00%
		HIV Antigen (p24) (Waived)		
HIV-1		* Denotes Correct Response		
AG: All Methods	30			
		*Negative	29	96.67%
		Positive	1	3.33%
Alere Determine HIV Ag/Ab	30	*Negative	29	96.67%
Alere Determine HIV Ag/Ab		Positive	1	3.33%
HIV-2		* Denotes Correct Response		
AG: All Methods	30			



		*Positive	30	100.00%
Alere Determine HIV Ag/Ab	30	*Positive	30	100.00%
934 - HIV Antibody/Antigen Combo (5 specimens)	Labs	Response	Number	Percent
		HIV Antibody (Waived)		
HIV-1		* Denotes Correct Response		
AG: All Methods	98			
		*Positive	97	98.98%
		Negative	1	1.02%
Alere Clearview HIV 1/2 STAT PAK (WB)	45	*Positive	45	100.00%
Orasure OraQuick Anti-HIV	10	*Positive	9	90.00%
Trinity Biotech Uni-Gold Recombigen HIV	2	*Positive	2	100.00%
Alere Determine HIV Ag/Ab	30	*Positive	30	100.00%
Chembio HIV 1/2 STAT-PAK Dipstick	3	*Positive	3	100.00%
Chembio HIV 1/2 STAT-PAK Cassette	3	*Positive	3	100.00%
bioLytical INSTI HIV Antibody	5	*Positive	5	100.00%
Orasure OraQuick Anti-HIV	10	Negative	1	10.00%
HIV-2		* Denotes Correct Response		
AG: All Methods	98			
		*Negative	93	94.90%
		Positive	5	5.10%
Alere Clearview HIV 1/2 STAT PAK (WB)	45	*Negative	45	100.00%
Orasure OraQuick Anti-HIV	10	*Negative	9	90.00%
Trinity Biotech Uni-Gold Recombigen HIV	2	*Negative	2	100.00%
Alere Determine HIV Ag/Ab	30	*Negative	26	86.67%
Chembio HIV 1/2 STAT-PAK Dipstick	3	*Negative	3	100.00%
Chembio HIV 1/2 STAT-PAK Cassette	3	*Negative	3	100.00%
bioLytical INSTI HIV Antibody	5	*Negative	5	100.00%
Orasure OraQuick Anti-HIV	10	Positive	1	10.00%
Alere Determine HIV Ag/Ab	30	Positive	4	13.33%
HIV-3		* Denotes Correct Response		
AG: All Methods	8			
		*Positive	8	100.00%
Alere Clearview HIV 1/2 STAT PAK (WB)	4	*Positive	4	100.00%
Orasure OraQuick Anti-HIV	1	*Positive	1	100.00%
Trinity Biotech Uni-Gold Recombigen HIV	1	*Positive	1	100.00%
Alere Determine HIV Ag/Ab	1	*Positive	1	100.00%
bioLytical INSTI HIV Antibody	1	*Positive	1	100.00%
HIV-4		* Denotes Correct Response		
AG: All Methods	8			
		*Negative	8	100.00%
Alere Clearview HIV 1/2 STAT PAK (WB)	4	*Negative	4	100.00%
Orasure OraQuick Anti-HIV	1	*Negative	1	100.00%
Trinity Biotech Uni-Gold Recombigen HIV	1	*Negative	1	100.00%
Alere Determine HIV Ag/Ab	1	*Negative	1	100.00%
bioLytical INSTI HIV Antibody	1	*Negative	1	100.00%
HIV-5		* Denotes Correct Response		
AG: All Methods	8			
		*Negative	8	100.00%
Alere Clearview HIV 1/2 STAT PAK (WB)	4	*Negative	4	100.00%
Orasure OraQuick Anti-HIV	1	*Negative	1	100.00%
Trinity Biotech Uni-Gold Recombigen HIV	1	*Negative	1	100.00%
Alere Determine HIV Ag/Ab	1	*Negative	1	100.00%
bioLytical INSTI HIV Antibody	1	*Negative	1	100.00%
		HIV Antibody (International)		
HIV-1		* Denotes Correct Response		
AG: All Methods	5			
		*Positive	5	100.00%
Alere Clearview HIV 1/2 STAT PAK (WB)	1	*Positive	1	100.00%
Intl OCD Vitros ECI/ECiQ	3	*Positive	3	100.00%
Intl Chembio HIV 1/2 STAT-PAK Dipstick	1	*Positive	1	100.00%
HIV-2		* Denotes Correct Response		
AG: All Methods	5			
		*Negative	5	100.00%
Alere Clearview HIV 1/2 STAT PAK (WB)	1	*Negative	1	100.00%
Intl OCD Vitros ECI/ECiQ	3	*Negative	3	100.00%
Intl Chembio HIV 1/2 STAT-PAK Dipstick	1	*Negative	1	100.00%
HIV-3		* Denotes Correct Response		
AG: All Methods	3			



		*Positive	3	100.00%
Intl OCD Vitros ECI/ECiQ	3	*Positive	3	100.00%
HIV-4		* Denotes Correct Response		
AG: All Methods	3			
		*Negative	3	100.00%
Intl OCD Vitros ECI/ECiQ	3	*Negative	3	100.00%
HIV-5		* Denotes Correct Response		
AG: All Methods	2			
		*Negative	2	100.00%
Intl OCD Vitros ECI/ECiQ	2	*Negative	2	100.00%
		HIV Antigen (p24) (Waived)		
HIV-1		* Denotes Correct Response		
AG: All Methods	30			
		*Negative	29	96.67%
		Positive	1	3.33%
Alere Determine HIV Ag/Ab	30	*Negative	29	96.67%
Alere Determine HIV Ag/Ab		Positive	1	3.33%
HIV-2		* Denotes Correct Response		
AG: All Methods	30			
		*Positive	30	100.00%
Alere Determine HIV Ag/Ab	30	*Positive	30	100.00%
HIV-3		* Denotes Correct Response		
AG: All Methods	1			
		*Negative	1	100.00%
Alere Determine HIV Ag/Ab	1	*Negative	1	100.00%
HIV-4		* Denotes Correct Response		
AG: All Methods	1			
		*Negative	1	100.00%
Alere Determine HIV Ag/Ab	1	*Negative	1	100.00%
HIV-5		* Denotes Correct Response		
AG: All Methods	1			
		*Negative	1	100.00%
Alere Determine HIV Ag/Ab	1	*Negative	1	100.00%
644 - H. pylori	Labs	Response	Number	Percent
		H. pylori Ab Detect		
HP-1		* Denotes Correct Response		
AG: All Methods	64			
		*Negative	63	98.44%
		Positive	1	1.56%
Beck Coulter FlexSure	3	*Negative	3	100.00%
Beck Coulter ICON HP	2	*Negative	2	100.00%
McKesson Medi-Lab Perf H. pylori	7	*Negative	7	100.00%
Polymedco Poly stat (WB - Waived)	1	*Negative	1	100.00%
Henry Schein One Step+ (H. pylori)	10	*Negative	10	100.00%
Select Diag H. pylori Cassette	2	*Negative	2	100.00%
Quidel QuickVue H. pylori gII (S/P)	3	*Negative	3	100.00%
Quidel QuickVue H. pylori gII (WB)	7	*Negative	6	85.71%
SP Brand (Cardinal) HP Rapid Test	9	*Negative	9	100.00%
Alere Clearview II S/P (mod)	2	*Negative	2	100.00%
Alere Clearview II WB (waived)	3	*Negative	3	100.00%
LifeSign Status H. pylori	2	*Negative	2	100.00%
Consult Diag H pylori	12	*Negative	12	100.00%
Jant Accutest H pylori	1	*Negative	1	100.00%
Quidel QuickVue H. pylori gII (WB)	7	Positive	1	14.29%
929 - Genital Viral Detection I	Labs	Response	Number	Percent
		HSV 1		
HVD-1		* Denotes Correct Response		
AG: All Methods	1			
		*Positive	1	100.00%
Hologic Aptima HSV 1 & 2	1	*Positive	1	100.00%
HVD-2		* Denotes Correct Response		
AG: All Methods	1			
		*Negative	1	100.00%
Hologic Aptima HSV 1 & 2	1	*Negative	1	100.00%
HVD-3		* Denotes Correct Response		
AG: All Methods	1			
		*Positive	1	100.00%
Hologic Aptima HSV 1 & 2	1	*Positive	1	100.00%



HVD-4		* Denotes Correct Response		
AG: All Methods	1			
		*Negative	1	100.00%
Hologic Aptima HSV 1 & 2	1	*Negative	1	100.00%
HVD-5		* Denotes Correct Response		
AG: All Methods	1			
		*Negative	1	100.00%
Hologic Aptima HSV 1 & 2	1	*Negative	1	100.00%
HSV 2				
HVD-1		* Denotes Correct Response		
AG: All Methods	1			
		*Negative	1	100.00%
Hologic Aptima HSV 1 & 2	1	*Negative	1	100.00%
HVD-2		* Denotes Correct Response		
AG: All Methods	1			
		*Positive	1	100.00%
Hologic Aptima HSV 1 & 2	1	*Positive	1	100.00%
HVD-3		* Denotes Correct Response		
AG: All Methods	1			
		*Negative	1	100.00%
Hologic Aptima HSV 1 & 2	1	*Negative	1	100.00%
HVD-4		* Denotes Correct Response		
AG: All Methods	1			
		*Positive	1	100.00%
Hologic Aptima HSV 1 & 2	1	*Positive	1	100.00%
HVD-5		* Denotes Correct Response		
AG: All Methods	1			
		*Negative	1	100.00%
Hologic Aptima HSV 1 & 2	1	*Negative	1	100.00%
603 - Blood Cell Identification:610 - Regular Hematology	Labs	Response	Number	Percent
		Blood Cell Identification		
ID-1		* Denotes Correct Response		
AG: All Methods	61			
		*Eosinophil, any stage	61	100.00%
Blood Cell Microscopy	61	*Eosinophil, any stage	61	100.00%
ID-2		* Denotes Correct Response		
AG: All Methods	61			
		*Seg/band neutrophil with toxic granulation	2	3.28%
		*Segmented or band neutrophil	57	93.44%
		Seg/band neutrophil with Dohle body	2	3.28%
Blood Cell Microscopy	61	*Seg/band neutrophil with toxic granulation	2	3.28%
Blood Cell Microscopy		*Segmented or band neutrophil	57	93.44%
Blood Cell Microscopy		Seg/band neutrophil with Dohle body	2	3.28%
ID-3		* Denotes Correct Response		
AG: All Methods	60			
		*Immature/abnormal cell - Referred	1	1.67%
		*Target cell (codocyte)	59	98.33%
Blood Cell Microscopy	60	*Immature/abnormal cell - Referred	1	1.67%
Blood Cell Microscopy		*Target cell (codocyte)	59	98.33%
ID-4		* Denotes Correct Response		
AG: All Methods	60			
		*Immature/abnormal cell - Referred	1	1.67%
		*Metarubricyte	1	1.67%
		*Nucleated red cell	57	95.00%
		Red Blood Cell, normal (RBC) (Erythrocyte)	1	1.67%
Blood Cell Microscopy	60	*Immature/abnormal cell - Referred	1	1.67%
Blood Cell Microscopy		*Metarubricyte	1	1.67%
Blood Cell Microscopy		*Nucleated red cell	57	95.00%
Blood Cell Microscopy		Red Blood Cell, normal (RBC) (Erythrocyte)	1	1.67%
ID-5		* Denotes Correct Response		
AG: All Methods	61			
		*Platelet, giant (macrothrombocyte)	2	3.28%
		*Platelet, normal	59	96.72%
Blood Cell Microscopy	61	*Platelet, giant (macrothrombocyte)	2	3.28%
Blood Cell Microscopy		*Platelet, normal	59	96.72%
ID-6		* Denotes Correct Response		
AG: All Methods	60			
		*Fungii, intracellular	31	51.67%



		Bacteria, extracellular	1	1.67%
		Identification Unknown - Referred	10	16.67%
		Immature/abnormal cell - Referred	3	5.00%
		Leukocyte with phagocytized bacteria	5	8.33%
		Other, not listed	6	10.00%
		Platelet satellitism	1	1.67%
		Platelet, normal	1	1.67%
		Seg/band neutrophil with Dohle body	1	1.67%
		Seg/band neutrophil with vacuoles	1	1.67%
Blood Cell Microscopy	60	*Fungii, intracellular	31	51.67%
Blood Cell Microscopy		Bacteria, extracellular	1	1.67%
Blood Cell Microscopy		Identification Unknown - Referred	10	16.67%
Blood Cell Microscopy		Immature/abnormal cell - Referred	3	5.00%
Blood Cell Microscopy		Leukocyte with phagocytized bacteria	5	8.33%
Blood Cell Microscopy		Other, not listed	6	10.00%
Blood Cell Microscopy		Platelet satellitism	1	1.67%
Blood Cell Microscopy		Platelet, normal	1	1.67%
Blood Cell Microscopy		Seg/band neutrophil with Dohle body	1	1.67%
Blood Cell Microscopy		Seg/band neutrophil with vacuoles	1	1.67%
		Blood Cell Identification (International)		
ID-1		* Denotes Correct Response		
AG: All Methods	2			
		*Eosinophil, any stage	2	100.00%
Intl Blood Cell Microscopy	2	*Eosinophil, any stage	2	100.00%
ID-2		* Denotes Correct Response		
AG: All Methods	2			
		*Segmented or band neutrophil	2	100.00%
Intl Blood Cell Microscopy	2	*Segmented or band neutrophil	2	100.00%
ID-3		* Denotes Correct Response		
AG: All Methods	2			
		*Target cell (codocyte)	2	100.00%
Intl Blood Cell Microscopy	2	*Target cell (codocyte)	2	100.00%
ID-4		* Denotes Correct Response		
AG: All Methods	2			
		*Nucleated red cell	2	100.00%
Intl Blood Cell Microscopy	2	*Nucleated red cell	2	100.00%
ID-5		* Denotes Correct Response		
AG: All Methods	2			
		*Platelet, giant (macrothrombocyte)	2	100.00%
Intl Blood Cell Microscopy	2	*Platelet, giant (macrothrombocyte)	2	100.00%
ID-6		* Denotes Correct Response		
AG: All Methods	2			
		*Fungii, intracellular	1	50.00%
		Identification Unknown - Referred	1	50.00%
Intl Blood Cell Microscopy	2	*Fungii, intracellular	1	50.00%
Intl Blood Cell Microscopy		Identification Unknown - Referred	1	50.00%
609 - Waived/PPM Tests I;675 - Infectious Mononucleosis;753 - Infectious Mononucleosis (Waived) (1 specimen);754 - Infectious Mono/RF Combo	Labs	Response	Number	Percent
		Infectious Mono (Waived)		
IM-1		* Denotes Correct Response		
AG: All Methods	473			
		*Negative	469	99.15%
		Positive	4	0.85%
Alere Acceava (S)	1	*Negative	1	100.00%
Alere Acceava (WB)	56	*Negative	54	96.43%
Sekisui Diagnostics OSOM Mono (WB)	187	*Negative	186	99.47%
Henry Schein One-Step+ Mono	45	*Negative	45	100.00%
Instant Technologies i-Screen (WB)	1	*Negative	1	100.00%
Inverness Signify Mono	4	*Negative	4	100.00%
LifeSign Status Inf Mono	29	*Negative	29	100.00%
McKesson Medi-Lab Perf - Inf Mono	10	*Negative	10	100.00%
Polymedco Poly stat Mono - (WB)	6	*Negative	6	100.00%
Quidel QuickVue+ (WB)	10	*Negative	10	100.00%
SP Brand (Cardinal) Rapid Test	30	*Negative	30	100.00%
Alere Clearview Mono Cassette	18	*Negative	17	94.44%
Pulse Scientific	2	*Negative	2	100.00%
Diag Test Group Clarity Mono	2	*Negative	2	100.00%
Fisher HC Sure-View Mono	6	*Negative	6	100.00%



Immunostics Inc. ImmunoDetector Mono (WB)	2	*Negative	2	100.00%
Meridian Immunocard STAT! Mono (WB)	6	*Negative	6	100.00%
NDC ProAdvantage Mono	1	*Negative	1	100.00%
Consult Diag Mono Cassette	42	*Negative	42	100.00%
Alere Acceava (WB) Mono Cassette	15	*Negative	15	100.00%
Alere Acceava (WB)	56	Positive	2	3.57%
Sekisui Diagnostics OSOM Mono (WB)	187	Positive	1	0.53%
Alere Clearview Mono Cassette	18	Positive	1	5.56%
		Infectious Mononucleosis		
IM-1		* Denotes Correct Response		
AG: All Methods	45			
		*Negative	45	100.00%
Beck Coulter ICON (WB)	1	*Negative	1	100.00%
Alere Acceava (WB)	1	*Negative	1	100.00%
Sekisui Diagnostics OSOM Mono (WB)	9	*Negative	9	100.00%
LifeSign Status Inf Mono	2	*Negative	2	100.00%
Lifesign Unistep	1	*Negative	1	100.00%
McKesson Medi-Lab Perf - Inf Mono	2	*Negative	2	100.00%
Polymedco Poly stat Mono (S/P)	2	*Negative	2	100.00%
Quidel QuickVue+ (S/P)	2	*Negative	2	100.00%
Quidel QuickVue+ (WB)	8	*Negative	8	100.00%
Seradyn Color Slide II	3	*Negative	3	100.00%
SP Brand (Cardinal) Rapid Test	1	*Negative	1	100.00%
Alere Clearview Mono Cassette	8	*Negative	8	100.00%
Sekisui Diagnostics OSOM Mono Non-Waived	1	*Negative	1	100.00%
Consult Diag Mono Cassette	4	*Negative	4	100.00%
675 - Infectious Mononucleosis;754 - Infectious Mono/RF Combo	Labs	Response	Number	Percent
		Infectious Mononucleosis		
IM-1		* Denotes Correct Response		
AG: All Methods	45			
		*Negative	45	100.00%
Beck Coulter ICON (WB)	1	*Negative	1	100.00%
Alere Acceava (WB)	1	*Negative	1	100.00%
Sekisui Diagnostics OSOM Mono (WB)	9	*Negative	9	100.00%
LifeSign Status Inf Mono	2	*Negative	2	100.00%
Lifesign Unistep	1	*Negative	1	100.00%
McKesson Medi-Lab Perf - Inf Mono	2	*Negative	2	100.00%
Polymedco Poly stat Mono (S/P)	2	*Negative	2	100.00%
Quidel QuickVue+ (S/P)	2	*Negative	2	100.00%
Quidel QuickVue+ (WB)	8	*Negative	8	100.00%
Seradyn Color Slide II	3	*Negative	3	100.00%
SP Brand (Cardinal) Rapid Test	1	*Negative	1	100.00%
Alere Clearview Mono Cassette	8	*Negative	8	100.00%
Sekisui Diagnostics OSOM Mono Non-Waived	1	*Negative	1	100.00%
Consult Diag Mono Cassette	4	*Negative	4	100.00%
IM-2		* Denotes Correct Response		
AG: All Methods	45			
		*Negative	45	100.00%
Beck Coulter ICON (WB)	1	*Negative	1	100.00%
Alere Acceava (WB)	1	*Negative	1	100.00%
Sekisui Diagnostics OSOM Mono (WB)	9	*Negative	9	100.00%
LifeSign Status Inf Mono	2	*Negative	2	100.00%
Lifesign Unistep	1	*Negative	1	100.00%
McKesson Medi-Lab Perf - Inf Mono	2	*Negative	2	100.00%
Polymedco Poly stat Mono (S/P)	2	*Negative	2	100.00%
Quidel QuickVue+ (S/P)	2	*Negative	2	100.00%
Quidel QuickVue+ (WB)	8	*Negative	8	100.00%
Seradyn Color Slide II	3	*Negative	3	100.00%
SP Brand (Cardinal) Rapid Test	1	*Negative	1	100.00%
Alere Clearview Mono Cassette	8	*Negative	8	100.00%
Sekisui Diagnostics OSOM Mono Non-Waived	1	*Negative	1	100.00%
Consult Diag Mono Cassette	4	*Negative	4	100.00%
IM-3		* Denotes Correct Response		
AG: All Methods	45			
		*Positive	45	100.00%
Beck Coulter ICON (WB)	1	*Positive	1	100.00%
Alere Acceava (WB)	1	*Positive	1	100.00%
Sekisui Diagnostics OSOM Mono (WB)	9	*Positive	9	100.00%



LifeSign Status Inf Mono	2	*Positive	2	100.00%
Lifesign Unistep	1	*Positive	1	100.00%
McKesson Medi-Lab Perf - Inf Mono	2	*Positive	2	100.00%
Polymedco Poly stat Mono (S/P)	2	*Positive	2	100.00%
Quidel QuickVue+ (S/P)	2	*Positive	2	100.00%
Quidel QuickVue+ (WB)	8	*Positive	8	100.00%
Seradyn Color Slide II	3	*Positive	3	100.00%
SP Brand (Cardinal) Rapid Test	1	*Positive	1	100.00%
Alere Clearview Mono Cassette	8	*Positive	8	100.00%
Sekisui Diagnostics OSOM Mono Non-Waived	1	*Positive	1	100.00%
Consult Diag Mono Cassette	4	*Positive	4	100.00%
IM-4		* Denotes Correct Response		
AG: All Methods	45			
		*Positive	45	100.00%
Beck Coulter ICON (WB)	1	*Positive	1	100.00%
Alere Acceava (WB)	1	*Positive	1	100.00%
Sekisui Diagnostics OSOM Mono (WB)	9	*Positive	9	100.00%
LifeSign Status Inf Mono	2	*Positive	2	100.00%
Lifesign Unistep	1	*Positive	1	100.00%
McKesson Medi-Lab Perf - Inf Mono	2	*Positive	2	100.00%
Polymedco Poly stat Mono (S/P)	2	*Positive	2	100.00%
Quidel QuickVue+ (S/P)	2	*Positive	2	100.00%
Quidel QuickVue+ (WB)	8	*Positive	8	100.00%
Seradyn Color Slide II	3	*Positive	3	100.00%
SP Brand (Cardinal) Rapid Test	1	*Positive	1	100.00%
Alere Clearview Mono Cassette	8	*Positive	8	100.00%
Sekisui Diagnostics OSOM Mono Non-Waived	1	*Positive	1	100.00%
Consult Diag Mono Cassette	4	*Positive	4	100.00%
IM-5		* Denotes Correct Response		
AG: All Methods	45			
		*Negative	45	100.00%
Beck Coulter ICON (WB)	1	*Negative	1	100.00%
Alere Acceava (WB)	1	*Negative	1	100.00%
Sekisui Diagnostics OSOM Mono (WB)	9	*Negative	9	100.00%
LifeSign Status Inf Mono	2	*Negative	2	100.00%
Lifesign Unistep	1	*Negative	1	100.00%
McKesson Medi-Lab Perf - Inf Mono	2	*Negative	2	100.00%
Polymedco Poly stat Mono (S/P)	2	*Negative	2	100.00%
Quidel QuickVue+ (S/P)	2	*Negative	2	100.00%
Quidel QuickVue+ (WB)	8	*Negative	8	100.00%
Seradyn Color Slide II	3	*Negative	3	100.00%
SP Brand (Cardinal) Rapid Test	1	*Negative	1	100.00%
Alere Clearview Mono Cassette	8	*Negative	8	100.00%
Sekisui Diagnostics OSOM Mono Non-Waived	1	*Negative	1	100.00%
Consult Diag Mono Cassette	4	*Negative	4	100.00%
		Infectious Mono (Waived)		
IM-1		* Denotes Correct Response		
AG: All Methods	473			
		*Negative	469	99.15%
		Positive	4	0.85%
Alere Acceava (S)	1	*Negative	1	100.00%
Alere Acceava (WB)	56	*Negative	54	96.43%
Sekisui Diagnostics OSOM Mono (WB)	187	*Negative	186	99.47%
Henry Schein One-Step+ Mono	45	*Negative	45	100.00%
Instant Technologies i-Screen (WB)	1	*Negative	1	100.00%
Inverness Signify Mono	4	*Negative	4	100.00%
LifeSign Status Inf Mono	29	*Negative	29	100.00%
McKesson Medi-Lab Perf - Inf Mono	10	*Negative	10	100.00%
Polymedco Poly stat Mono - (WB)	6	*Negative	6	100.00%
Quidel QuickVue+ (WB)	10	*Negative	10	100.00%
SP Brand (Cardinal) Rapid Test	30	*Negative	30	100.00%
Alere Clearview Mono Cassette	18	*Negative	17	94.44%
Pulse Scientific	2	*Negative	2	100.00%
Diag Test Group Clarity Mono	2	*Negative	2	100.00%
Fisher HC Sure-View Mono	6	*Negative	6	100.00%
Immunostics Inc. ImmunoDetector Mono (WB)	2	*Negative	2	100.00%
Meridian Immunocard STAT! Mono (WB)	6	*Negative	6	100.00%
NDC ProAdvantage Mono	1	*Negative	1	100.00%
Consult Diag Mono Cassette	42	*Negative	42	100.00%



Alere Acceava (WB) Mono Cassette	15	*Negative	15	100.00%
Alere Acceava (WB)	56	Positive	2	3.57%
Sekisui Diagnostics OSOM Mono (WB)	187	Positive	1	0.53%
Alere Clearview Mono Cassette	18	Positive	1	5.56%
		Infectious Mononucleosis (International)		
IM-1		* Denotes Correct Response		
AG: All Methods	2			
		*Negative	2	100.00%
Alere Clearview Mono Cassette	1	*Negative	1	100.00%
Fisher HC Sure-Vue Mono	1	*Negative	1	100.00%
IM-2		* Denotes Correct Response		
AG: All Methods	2			
		*Negative	2	100.00%
Alere Clearview Mono Cassette	1	*Negative	1	100.00%
Fisher HC Sure-Vue Mono	1	*Negative	1	100.00%
IM-3		* Denotes Correct Response		
AG: All Methods	2			
		*Positive	2	100.00%
Alere Clearview Mono Cassette	1	*Positive	1	100.00%
Fisher HC Sure-Vue Mono	1	*Positive	1	100.00%
IM-4		* Denotes Correct Response		
AG: All Methods	2			
		*Positive	2	100.00%
Alere Clearview Mono Cassette	1	*Positive	1	100.00%
Fisher HC Sure-Vue Mono	1	*Positive	1	100.00%
IM-5		* Denotes Correct Response		
AG: All Methods	2			
		*Negative	2	100.00%
Alere Clearview Mono Cassette	1	*Negative	1	100.00%
Fisher HC Sure-Vue Mono	1	*Negative	1	100.00%
637 - KOH - Glass Slides (Skin Hair Nail)	Labs	Response	Number	Percent
		KOH Preparation		
KOH-1		* Denotes Correct Response		
AG: All Methods	19			
		*No fungal elements/yeast observed	18	94.74%
		Fungal elements/yeast observed	1	5.26%
KOH Glass Slide	19	*No fungal elements/yeast observed	18	94.74%
KOH Glass Slide		Fungal elements/yeast observed	1	5.26%
KOH-2		* Denotes Correct Response		
AG: All Methods	19			
		*Fungal elements/yeast observed	19	100.00%
KOH Glass Slide	19	*Fungal elements/yeast observed	19	100.00%
768 - MRSA (Molecular) (2 specimens);690 - MRSA (Molecular) (5 specimens)	Labs	Response	Number	Percent
		MRSA (Molecular)		
MRSA-1		* Denotes Correct Response		
AG: All Methods	1			
		*Positive	1	100.00%
Cepheid GeneXpert - Xpert MRSA	1	*Positive	1	100.00%
MRSA-2		* Denotes Correct Response		
AG: All Methods	1			
		*Positive	1	100.00%
Cepheid GeneXpert - Xpert MRSA	1	*Positive	1	100.00%
758 - Mycoplasma Antibody	Labs	Response	Number	Percent
		Mycoplasma Antibody		
MY-1		* Denotes Correct Response		
AG: All Methods	2			
		*Negative	2	100.00%
Meridian ImmunoCard	2	*Negative	2	100.00%
791 - Fecal Occult Blood Add-On;601 - Waived/PPM Tests II;609 - Waived/PPM Tests I;647 - Mini-Waived Package	Labs	Response	Number	Percent
		Fecal Occult Blood		
OC-1		* Denotes Correct Response		
MG: Guaiac methods	689			
		*Positive	670	97.24%
		Negative	19	2.76%
OC-1		* Denotes Correct Response		



MG: Immunochemical methods	100			
		*Positive	95	95.00%
		Negative	5	5.00%
950 - Post Vasectomy Sperm Analysis	Labs	Response	Number	Percent
		Sperm Count (Post Vasectomy)		
PVAS-1		* Denotes Correct Response		
AG: All Methods	14			
		*Positive for Sperm	14	100.00%
Hemocytometer	4	*Positive for Sperm	4	100.00%
Post Vasectomy Microscopy	10	*Positive for Sperm	10	100.00%
PVAS-2		* Denotes Correct Response		
AG: All Methods	14			
		*Negative for Sperm	14	100.00%
Hemocytometer	4	*Negative for Sperm	4	100.00%
Post Vasectomy Microscopy	10	*Negative for Sperm	10	100.00%
754 - Infectious Mono/RF Combo;653 - Rheumatoid Factor	Labs	Response	Number	Percent
		Rheumatoid Factor		
RF-1		* Denotes Correct Response		
AG: All Methods	10			
		*Positive	9	90.00%
		N/A (Negative)	1	10.00%
ASI	1	*Positive	1	100.00%
Siemens RapiTex	1	*Positive	1	100.00%
Seradyn Color Slide	2	*Positive	1	50.00%
Wampole ColorCard RF	2	*Positive	2	100.00%
Wampole Rheumatex	2	*Positive	2	100.00%
Fisher HC Sure-Vue RF	1	*Positive	1	100.00%
Biokit Rhemajet RF	1	*Positive	1	100.00%
Seradyn Color Slide	2	N/A (Negative)	1	50.00%
RF-2		* Denotes Correct Response		
AG: All Methods	10			
		*Positive	10	100.00%
ASI	1	*Positive	1	100.00%
Siemens RapiTex	1	*Positive	1	100.00%
Seradyn Color Slide	2	*Positive	2	100.00%
Wampole ColorCard RF	2	*Positive	2	100.00%
Wampole Rheumatex	2	*Positive	2	100.00%
Fisher HC Sure-Vue RF	1	*Positive	1	100.00%
Biokit Rhemajet RF	1	*Positive	1	100.00%
RF-3		* Denotes Correct Response		
AG: All Methods	10			
		*N/A (Negative)	9	90.00%
		Positive	1	10.00%
ASI	1	*N/A (Negative)	1	100.00%
Siemens RapiTex	1	*N/A (Negative)	1	100.00%
Seradyn Color Slide	2	*N/A (Negative)	1	50.00%
Wampole ColorCard RF	2	*N/A (Negative)	2	100.00%
Wampole Rheumatex	2	*N/A (Negative)	2	100.00%
Fisher HC Sure-Vue RF	1	*N/A (Negative)	1	100.00%
Biokit Rhemajet RF	1	*N/A (Negative)	1	100.00%
Seradyn Color Slide	2	Positive	1	50.00%
RF-4		* Denotes Correct Response		
AG: All Methods	10			
		*N/A (Negative)	9	90.00%
		Positive	1	10.00%
ASI	1	*N/A (Negative)	1	100.00%
Siemens RapiTex	1	*N/A (Negative)	1	100.00%
Seradyn Color Slide	2	*N/A (Negative)	1	50.00%
Wampole ColorCard RF	2	*N/A (Negative)	2	100.00%
Wampole Rheumatex	2	*N/A (Negative)	2	100.00%
Fisher HC Sure-Vue RF	1	*N/A (Negative)	1	100.00%
Biokit Rhemajet RF	1	*N/A (Negative)	1	100.00%
Seradyn Color Slide	2	Positive	1	50.00%
RF-5		* Denotes Correct Response		
AG: All Methods	10			
		*Positive	9	90.00%
		N/A (Negative)	1	10.00%



ASI	1	*Positive	1	100.00%
Siemens RapiTex	1	*Positive	1	100.00%
Seradyn Color Slide	2	*Positive	1	50.00%
Wampole ColorCard RF	2	*Positive	2	100.00%
Wampole Rheumatex	2	*Positive	2	100.00%
Fisher HC Sure-Vue RF	1	*Positive	1	100.00%
Biokit Rhemajet RF	1	*Positive	1	100.00%
Seradyn Color Slide	2	N/A (Negative)	1	50.00%
		Rheumatoid Factor (International)		
RF-1		* Denotes Correct Response		
AG: All Methods	3			
		*Positive	3	100.00%
Fisher HC Sure-Vue RF	2	*Positive	2	100.00%
Immunostics Rhuemacol/LEX	1	*Positive	1	100.00%
RF-2		* Denotes Correct Response		
AG: All Methods	3			
		*Positive	3	100.00%
Fisher HC Sure-Vue RF	2	*Positive	2	100.00%
Immunostics Rhuemacol/LEX	1	*Positive	1	100.00%
RF-3		* Denotes Correct Response		
AG: All Methods	3			
		*Negative	3	100.00%
Fisher HC Sure-Vue RF	2	*Negative	2	100.00%
Immunostics Rhuemacol/LEX	1	*Negative	1	100.00%
RF-4		* Denotes Correct Response		
AG: All Methods	3			
		*Negative	3	100.00%
Fisher HC Sure-Vue RF	2	*Negative	2	100.00%
Immunostics Rhuemacol/LEX	1	*Negative	1	100.00%
RF-5		* Denotes Correct Response		
AG: All Methods	3			
		*Positive	3	100.00%
Fisher HC Sure-Vue RF	2	*Positive	2	100.00%
Immunostics Rhuemacol/LEX	1	*Positive	1	100.00%
645 - Fetal Membrane Rupture	Labs	Response	Number	Percent
		Fetal Membranes Rupture PAMG-1		
ROM-1		* Denotes Correct Response		
AG: All Methods	17			
		*Negative	17	100.00%
AmniSure ROM Test	17	*Negative	17	100.00%
ROM-2		* Denotes Correct Response		
AG: All Methods	17			
		*Positive	17	100.00%
AmniSure ROM Test	17	*Positive	17	100.00%
601 - Waived/PPM Tests II;609 - Waived/PPM Tests I;785 - Presump ID/Colony Count-Urine & Strep Combo;691 - Group A Strep Antigen (5 specimens);783 - Group A Strep Antigen/Add-on (2 specimens);665 - Group A Strep Antigen (Waived) (2 specimens)	Labs	Response	Number	Percent
		Strep A Ag Det. (Waived)		
RS-1		* Denotes Correct Response		
AG: All Methods	852			
		*Negative	845	99.18%
		Positive	7	0.82%
Abbott Signify Strep A - Waived	6	*Negative	6	100.00%
Beck Coulter ICON DS Strep	1	*Negative	1	100.00%
Alere Acceava Strep A Dipstick	114	*Negative	114	100.00%
Fisher HC Sure-Vue Strep A (Waived)	9	*Negative	9	100.00%
Sekisui Diagnostics OSOM	33	*Negative	33	100.00%
Sekisui Diagnostics OSOM Ultra Strep	64	*Negative	64	100.00%
Henry Schein OneStep + Rapid Strep Cass	6	*Negative	6	100.00%
Henry Schein OneStep+ Rapid Strep Dip	70	*Negative	68	97.14%
Innovacon One Step Strep Ag	1	*Negative	1	100.00%
Instant Technologies I Strep	1	*Negative	1	100.00%
Inverness Signify Strep A Dipstick	9	*Negative	9	100.00%
McKesson Strep A Cas	2	*Negative	2	100.00%
McKesson Strep A Dip	34	*Negative	34	100.00%
McKesson Strep A Test - Twist	2	*Negative	2	100.00%
Permaxim Rediscreen Strep A Rapid	1	*Negative	1	100.00%



Polymedco Poly Stat Strep A (W)	1	*Negative	1	100.00%
PSS Rapid Test Strip Strep A	2	*Negative	1	50.00%
Quidel QuickVue Dipstick Strep	82	*Negative	82	100.00%
Quidel QuickVue In-Line - 25 Test	58	*Negative	58	100.00%
Quidel QuickVue+ Strep A	32	*Negative	32	100.00%
Select Diag Strep A Dip	3	*Negative	3	100.00%
SP (Cardinal) Rapid Strep Cassette	4	*Negative	4	100.00%
SP Brand (Cardinal) Rapid Strep A Dipstick	68	*Negative	68	100.00%
Alere Clearview Exact II Dip	39	*Negative	39	100.00%
Sekisui Diagnostics OSOM Ultra	7	*Negative	6	85.71%
Consult Diag Strep A Dipstick	101	*Negative	98	97.03%
Consult Diag Strep A Twist	1	*Negative	1	100.00%
MooreBrand Strep A Dipstick	4	*Negative	4	100.00%
NDC ProAdvantage Strep Strips	1	*Negative	1	100.00%
Meridian ImmunoCard STAT Strep A	1	*Negative	1	100.00%
LifeSign Bio Strep A Test STRip	1	*Negative	1	100.00%
Quidel Sofia/Sofia 2 Strep A FIA	34	*Negative	34	100.00%
BD Veritor Strep A	30	*Negative	30	100.00%
Alere ID NOW Strep A	15	*Negative	15	100.00%
Medline Strep A Dipstick	4	*Negative	4	100.00%
Roche Cobas Liat Strep A	11	*Negative	11	100.00%
Henry Schein OneStep+ Rapid Strep Dip	70	Positive	2	2.86%
PSS Rapid Test Strip Strep A	2	Positive	1	50.00%
Sekisui Diagnostics OSOM Ultra	7	Positive	1	14.29%
Consult Diag Strep A Dipstick	101	Positive	3	2.97%
RS-2		* Denotes Correct Response		
AG: All Methods	145			
		*Positive	139	95.86%
		Negative	6	4.14%
Abbott Signify Strep A - Waived	1	*Positive	1	100.00%
Alere Acceava Strep A Dipstick	7	*Positive	7	100.00%
Fisher HC Sure-Vue Strep A (Waived)	1	*Positive	1	100.00%
Sekisui Diagnostics OSOM	8	*Positive	8	100.00%
Sekisui Diagnostics OSOM Ultra Strep	24	*Positive	24	100.00%
Henry Schein OneStep + Rapid Strep Cass	1	*Positive	1	100.00%
Henry Schein OneStep+ Rapid Strep Dip	12	*Positive	10	83.33%
Innovacon One Step Strep Ag	1	*Positive	1	100.00%
McKesson Strep A Dip	6	*Positive	6	100.00%
Quidel QuickVue Dipstick Strep	14	*Positive	13	92.86%
Quidel QuickVue In-Line - 25 Test	14	*Positive	14	100.00%
Quidel QuickVue+ Strep A	4	*Positive	4	100.00%
SP Brand (Cardinal) Rapid Strep A Dipstick	13	*Positive	13	100.00%
Sekisui Diagnostics OSOM Ultra	1	*Positive	1	100.00%
Consult Diag Strep A Dipstick	13	*Positive	12	92.31%
MooreBrand Strep A Dipstick	2	*Positive	1	50.00%
BD Veritor Strep A	7	*Positive	7	100.00%
Alere ID NOW Strep A	4	*Positive	3	75.00%
Medline Strep A Dipstick	2	*Positive	2	100.00%
Roche Cobas Liat Strep A	10	*Positive	10	100.00%
Henry Schein OneStep+ Rapid Strep Dip	12	Negative	2	16.67%
Quidel QuickVue Dipstick Strep	14	Negative	1	7.14%
Consult Diag Strep A Dipstick	13	Negative	1	7.69%
MooreBrand Strep A Dipstick	2	Negative	1	50.00%
Alere ID NOW Strep A	4	Negative	1	25.00%
		Strep A Ag Det.		
RS-1		* Denotes Correct Response		
AG: All Methods	162			
		*Negative	161	99.38%
		Positive	1	0.62%
Abbott Signify Strep A - Waived	1	*Negative	1	100.00%
BD Check Strep A (Dip)	2	*Negative	2	100.00%
Beck Coulter ICON DS Strep	1	*Negative	1	100.00%
Alere Acceava Strep A Dipstick	3	*Negative	3	100.00%
Fisher HC Sure-Vue Select Strep Ag	1	*Negative	1	100.00%
Sekisui Diagnostics OSOM	14	*Negative	14	100.00%
Sekisui Diagnostics OSOM Ultra Strep	5	*Negative	5	100.00%
Henry Schein OneStep+ Rapid Strep Dip	7	*Negative	7	100.00%
Immunostics Immuno/Strep A Det	2	*Negative	2	100.00%
Inverness Signify Strep A Dipstick	1	*Negative	1	100.00%



McKesson Strep A Cas	2	*Negative	2	100.00%
McKesson Strep A Dip	9	*Negative	9	100.00%
Permaxim Rediscreen Strep A Rapid	1	*Negative	1	100.00%
PSS Rapid Test Strip Strep A	1	*Negative	1	100.00%
Quidel QuickVue Dipstick Strep	8	*Negative	8	100.00%
Quidel QuickVue In-Line - 25 Test	9	*Negative	9	100.00%
Quidel QuickVue+ Strep A	17	*Negative	17	100.00%
Select Diag Strep A Dip	1	*Negative	1	100.00%
SP Brand (Cardinal) Rapid Strep A Dipstick	1	*Negative	1	100.00%
Stanbio QuStick Strep A	2	*Negative	2	100.00%
Consult Diag Strep A Dipstick	18	*Negative	18	100.00%
Meridian ImmunoCard STAT Strep A	1	*Negative	1	100.00%
Quidel Sofia/Sofia 2 Strep A FIA	10	*Negative	10	100.00%
BD Veritor Strep A	6	*Negative	6	100.00%
Quidel Solana	6	*Negative	6	100.00%
Cepheid GeneXpert - Xpert Strep A	33	*Negative	32	96.97%
Cepheid GeneXpert - Xpert Strep A		Positive	1	3.03%
RS-2		* Denotes Correct Response		
AG: All Methods	162			
		*Positive	161	99.38%
		Negative	1	0.62%
Abbott Signify Strep A - Waived	1	*Positive	1	100.00%
BD Check Strep A (Dip)	2	*Positive	2	100.00%
Beck Coulter ICON DS Strep	1	*Positive	1	100.00%
Alere Acceava Strep A Dipstick	3	*Positive	3	100.00%
Fisher HC Sure-Vue Select Strep Ag	1	*Positive	1	100.00%
Sekisui Diagnostics OSOM	14	*Positive	14	100.00%
Sekisui Diagnostics OSOM Ultra Strep	5	*Positive	5	100.00%
Henry Schein OneStep+ Rapid Strep Dip	7	*Positive	7	100.00%
Immunostics Immuno/Strep A Det	2	*Positive	2	100.00%
Inverness Signify Strep A Dipstick	1	*Positive	1	100.00%
McKesson Strep A Cas	2	*Positive	2	100.00%
McKesson Strep A Dip	9	*Positive	9	100.00%
Permaxim Rediscreen Strep A Rapid	1	*Positive	1	100.00%
PSS Rapid Test Strip Strep A	1	*Positive	1	100.00%
Quidel QuickVue Dipstick Strep	8	*Positive	8	100.00%
Quidel QuickVue In-Line - 25 Test	9	*Positive	9	100.00%
Quidel QuickVue+ Strep A	17	*Positive	17	100.00%
Select Diag Strep A Dip	1	*Positive	1	100.00%
SP Brand (Cardinal) Rapid Strep A Dipstick	1	*Positive	1	100.00%
Stanbio QuStick Strep A	2	*Positive	2	100.00%
Consult Diag Strep A Dipstick	18	*Positive	18	100.00%
Meridian ImmunoCard STAT Strep A	1	*Positive	1	100.00%
Quidel Sofia/Sofia 2 Strep A FIA	10	*Positive	10	100.00%
BD Veritor Strep A	6	*Positive	6	100.00%
Quidel Solana	6	*Positive	6	100.00%
Cepheid GeneXpert - Xpert Strep A	33	*Positive	32	96.97%
Cepheid GeneXpert - Xpert Strep A		Negative	1	3.03%
		Strep A Ag Det. (International)		
RS-1		* Denotes Correct Response		
AG: All Methods	1			
		*Negative	1	100.00%
Alere Clearview Exact II Cas	1	*Negative	1	100.00%
RS-2		* Denotes Correct Response		
AG: All Methods	1			
		*Positive	1	100.00%
Alere Clearview Exact II Cas	1	*Positive	1	100.00%
691 - Group A Strep Antigen (5 specimens)	Labs	Response	Number	Percent
		Strep A Ag Det.		
RS-1		* Denotes Correct Response		
AG: All Methods	162			
		*Negative	161	99.38%
		Positive	1	0.62%
Abbott Signify Strep A - Waived	1	*Negative	1	100.00%
BD Check Strep A (Dip)	2	*Negative	2	100.00%
Beck Coulter ICON DS Strep	1	*Negative	1	100.00%
Alere Acceava Strep A Dipstick	3	*Negative	3	100.00%
Fisher HC Sure-Vue Select Strep Ag	1	*Negative	1	100.00%
Sekisui Diagnostics OSOM	14	*Negative	14	100.00%



Sekisui Diagnostics OSOM Ultra Strep	5	*Negative	5	100.00%
Henry Schein OneStep+ Rapid Strep Dip	7	*Negative	7	100.00%
Immunostics Immuno/Strep A Det	2	*Negative	2	100.00%
Inverness Signify Strep A Dipstick	1	*Negative	1	100.00%
McKesson Strep A Cas	2	*Negative	2	100.00%
McKesson Strep A Dip	9	*Negative	9	100.00%
Permaxim Rediscreen Strep A Rapid	1	*Negative	1	100.00%
PSS Rapid Test Strip Strep A	1	*Negative	1	100.00%
Quidel QuickVue Dipstick Strep	8	*Negative	8	100.00%
Quidel QuickVue In-Line - 25 Test	9	*Negative	9	100.00%
Quidel QuickVue+ Strep A	17	*Negative	17	100.00%
Select Diag Strep A Dip	1	*Negative	1	100.00%
SP Brand (Cardinal) Rapid Strep A Dipstick	1	*Negative	1	100.00%
Stanbio QuStick Strep A	2	*Negative	2	100.00%
Consult Diag Strep A Dipstick	18	*Negative	18	100.00%
Meridian ImmunoCard STAT Strep A	1	*Negative	1	100.00%
Quidel Sofia/Sofia 2 Strep A FIA	10	*Negative	10	100.00%
BD Veritor Strep A	6	*Negative	6	100.00%
Quidel Solana	6	*Negative	6	100.00%
Cepheid GeneXpert - Xpert Strep A	33	*Negative	32	96.97%
Cepheid GeneXpert - Xpert Strep A		Positive	1	3.03%
RS-2		* Denotes Correct Response		
AG: All Methods	162			
		*Positive	161	99.38%
		Negative	1	0.62%
Abbott Signify Strep A - Waived	1	*Positive	1	100.00%
BD Check Strep A (Dip)	2	*Positive	2	100.00%
Beck Coulter ICON DS Strep	1	*Positive	1	100.00%
Alere Acceava Strep A Dipstick	3	*Positive	3	100.00%
Fisher HC Sure-Vue Select Strep Ag	1	*Positive	1	100.00%
Sekisui Diagnostics OSOM	14	*Positive	14	100.00%
Sekisui Diagnostics OSOM Ultra Strep	5	*Positive	5	100.00%
Henry Schein OneStep+ Rapid Strep Dip	7	*Positive	7	100.00%
Immunostics Immuno/Strep A Det	2	*Positive	2	100.00%
Inverness Signify Strep A Dipstick	1	*Positive	1	100.00%
McKesson Strep A Cas	2	*Positive	2	100.00%
McKesson Strep A Dip	9	*Positive	9	100.00%
Permaxim Rediscreen Strep A Rapid	1	*Positive	1	100.00%
PSS Rapid Test Strip Strep A	1	*Positive	1	100.00%
Quidel QuickVue Dipstick Strep	8	*Positive	8	100.00%
Quidel QuickVue In-Line - 25 Test	9	*Positive	9	100.00%
Quidel QuickVue+ Strep A	17	*Positive	17	100.00%
Select Diag Strep A Dip	1	*Positive	1	100.00%
SP Brand (Cardinal) Rapid Strep A Dipstick	1	*Positive	1	100.00%
Stanbio QuStick Strep A	2	*Positive	2	100.00%
Consult Diag Strep A Dipstick	18	*Positive	18	100.00%
Meridian ImmunoCard STAT Strep A	1	*Positive	1	100.00%
Quidel Sofia/Sofia 2 Strep A FIA	10	*Positive	10	100.00%
BD Veritor Strep A	6	*Positive	6	100.00%
Quidel Solana	6	*Positive	6	100.00%
Cepheid GeneXpert - Xpert Strep A	33	*Positive	32	96.97%
Cepheid GeneXpert - Xpert Strep A		Negative	1	3.03%
RS-3		* Denotes Correct Response		
AG: All Methods	117			
		*Positive	117	100.00%
Abbott Signify Strep A - Waived	1	*Positive	1	100.00%
BD Check Strep A (Dip)	2	*Positive	2	100.00%
Beck Coulter ICON DS Strep	1	*Positive	1	100.00%
Alere Acceava Strep A Dipstick	3	*Positive	3	100.00%
Fisher HC Sure-Vue Select Strep Ag	1	*Positive	1	100.00%
Sekisui Diagnostics OSOM	14	*Positive	14	100.00%
Sekisui Diagnostics OSOM Ultra Strep	4	*Positive	4	100.00%
Henry Schein OneStep+ Rapid Strep Dip	7	*Positive	7	100.00%
Immunostics Immuno/Strep A Det	2	*Positive	2	100.00%
Inverness Signify Strep A Dipstick	1	*Positive	1	100.00%
McKesson Strep A Cas	1	*Positive	1	100.00%
McKesson Strep A Dip	4	*Positive	4	100.00%
Permaxim Rediscreen Strep A Rapid	1	*Positive	1	100.00%
PSS Rapid Test Strip Strep A	1	*Positive	1	100.00%



Quidel QuickVue Dipstick Strep	7	*Positive	7	100.00%
Quidel QuickVue In-Line - 25 Test	9	*Positive	9	100.00%
Quidel QuickVue+ Strep A	16	*Positive	16	100.00%
SP Brand (Cardinal) Rapid Strep A Dipstick	1	*Positive	1	100.00%
Consult Diag Strep A Dipstick	14	*Positive	14	100.00%
Meridian ImmunoCard STAT Strep A	1	*Positive	1	100.00%
Quidel Sofia/Sofia 2 Strep A FIA	5	*Positive	5	100.00%
BD Veritor Strep A	1	*Positive	1	100.00%
Quidel Solana	4	*Positive	4	100.00%
Cepheid GeneXpert - Xpert Strep A	16	*Positive	16	100.00%
RS-4		* Denotes Correct Response		
AG: All Methods	117			
		*Positive	116	99.15%
		Negative	1	0.85%
Abbott Signify Strep A - Waived	1	*Positive	1	100.00%
BD Check Strep A (Dip)	2	*Positive	2	100.00%
Beck Coulter ICON DS Strep	1	*Positive	1	100.00%
Alere Acceava Strep A Dipstick	3	*Positive	3	100.00%
Fisher HC Sure-Vue Select Strep Ag	1	*Positive	1	100.00%
Sekisui Diagnostics OSOM	14	*Positive	14	100.00%
Sekisui Diagnostics OSOM Ultra Strep	4	*Positive	4	100.00%
Henry Schein OneStep+ Rapid Strep Dip	7	*Positive	7	100.00%
Immunostics Immuno/Strep A Det	2	*Positive	2	100.00%
Inverness Signify Strep A Dipstick	1	*Positive	1	100.00%
McKesson Strep A Cas	1	*Positive	1	100.00%
McKesson Strep A Dip	4	*Positive	4	100.00%
Permaxim Rediscreen Strep A Rapid	1	*Positive	1	100.00%
PSS Rapid Test Strip Strep A	1	*Positive	1	100.00%
Quidel QuickVue Dipstick Strep	7	*Positive	7	100.00%
Quidel QuickVue In-Line - 25 Test	9	*Positive	9	100.00%
Quidel QuickVue+ Strep A	16	*Positive	15	93.75%
SP Brand (Cardinal) Rapid Strep A Dipstick	1	*Positive	1	100.00%
Consult Diag Strep A Dipstick	14	*Positive	14	100.00%
Meridian ImmunoCard STAT Strep A	1	*Positive	1	100.00%
Quidel Sofia/Sofia 2 Strep A FIA	5	*Positive	5	100.00%
BD Veritor Strep A	1	*Positive	1	100.00%
Quidel Solana	4	*Positive	4	100.00%
Cepheid GeneXpert - Xpert Strep A	16	*Positive	16	100.00%
Quidel QuickVue+ Strep A	16	Negative	1	6.25%
RS-5		* Denotes Correct Response		
AG: All Methods	117			
		*Negative	114	97.44%
		Positive	3	2.56%
Abbott Signify Strep A - Waived	1	*Negative	1	100.00%
BD Check Strep A (Dip)	2	*Negative	2	100.00%
Beck Coulter ICON DS Strep	1	*Negative	1	100.00%
Alere Acceava Strep A Dipstick	3	*Negative	3	100.00%
Fisher HC Sure-Vue Select Strep Ag	1	*Negative	1	100.00%
Sekisui Diagnostics OSOM	14	*Negative	13	92.86%
Sekisui Diagnostics OSOM Ultra Strep	4	*Negative	4	100.00%
Henry Schein OneStep+ Rapid Strep Dip	7	*Negative	7	100.00%
Immunostics Immuno/Strep A Det	2	*Negative	2	100.00%
Inverness Signify Strep A Dipstick	1	*Negative	1	100.00%
McKesson Strep A Cas	1	*Negative	1	100.00%
McKesson Strep A Dip	4	*Negative	4	100.00%
Permaxim Rediscreen Strep A Rapid	1	*Negative	1	100.00%
PSS Rapid Test Strip Strep A	1	*Negative	1	100.00%
Quidel QuickVue Dipstick Strep	7	*Negative	7	100.00%
Quidel QuickVue In-Line - 25 Test	9	*Negative	8	88.89%
Quidel QuickVue+ Strep A	16	*Negative	15	93.75%
SP Brand (Cardinal) Rapid Strep A Dipstick	1	*Negative	1	100.00%
Consult Diag Strep A Dipstick	14	*Negative	14	100.00%
Meridian ImmunoCard STAT Strep A	1	*Negative	1	100.00%
Quidel Sofia/Sofia 2 Strep A FIA	5	*Negative	5	100.00%
BD Veritor Strep A	1	*Negative	1	100.00%
Quidel Solana	4	*Negative	4	100.00%
Cepheid GeneXpert - Xpert Strep A	16	*Negative	16	100.00%
Sekisui Diagnostics OSOM	14	Positive	1	7.14%
Quidel QuickVue In-Line - 25 Test	9	Positive	1	11.11%



Quidel QuickVue+ Strep A	16	Positive	1	6.25%
		Strep A Ag Det. (Waived)		
RS-1		* Denotes Correct Response		
AG: All Methods	852			
		*Negative	845	99.18%
		Positive	7	0.82%
Abbott Signify Strep A - Waived	6	*Negative	6	100.00%
Beck Coulter ICON DS Strep	1	*Negative	1	100.00%
Alere Acceava Strep A Dipstick	114	*Negative	114	100.00%
Fisher HC Sure-Vue Strep A (Waived)	9	*Negative	9	100.00%
Sekisui Diagnostics OSOM	33	*Negative	33	100.00%
Sekisui Diagnostics OSOM Ultra Strep	64	*Negative	64	100.00%
Henry Schein OneStep + Rapid Strep Cass	6	*Negative	6	100.00%
Henry Schein OneStep+ Rapid Strep Dip	70	*Negative	68	97.14%
Innovacon One Step Strep Ag	1	*Negative	1	100.00%
Instant Technologies i Strep	1	*Negative	1	100.00%
Inverness Signify Strep A Dipstick	9	*Negative	9	100.00%
McKesson Strep A Cas	2	*Negative	2	100.00%
McKesson Strep A Dip	34	*Negative	34	100.00%
McKesson Strep A Test - Twist	2	*Negative	2	100.00%
Permaxim Rediscreen Strep A Rapid	1	*Negative	1	100.00%
Polymedco Poly Stat Strep A (W)	1	*Negative	1	100.00%
PSS Rapid Test Strip Strep A	2	*Negative	1	50.00%
Quidel QuickVue Dipstick Strep	82	*Negative	82	100.00%
Quidel QuickVue In-Line - 25 Test	58	*Negative	58	100.00%
Quidel QuickVue+ Strep A	32	*Negative	32	100.00%
Select Diag Strep A Dip	3	*Negative	3	100.00%
SP (Cardinal) Rapid Strep Cassette	4	*Negative	4	100.00%
SP Brand (Cardinal) Rapid Strep A Dipstick	68	*Negative	68	100.00%
Alere Clearview Exact II Dip	39	*Negative	39	100.00%
Sekisui Diagnostics OSOM Ultra	7	*Negative	6	85.71%
Consult Diag Strep A Dipstick	101	*Negative	98	97.03%
Consult Diag Strep A Twist	1	*Negative	1	100.00%
MooreBrand Strep A Dipstick	4	*Negative	4	100.00%
NDC ProAdvantage Strep Strips	1	*Negative	1	100.00%
Meridian ImmunoCard STAT Strep A	1	*Negative	1	100.00%
LifeSign Bio Strep A Test STRip	1	*Negative	1	100.00%
Quidel Sofia/Sofia 2 Strep A FIA	34	*Negative	34	100.00%
BD Veritor Strep A	30	*Negative	30	100.00%
Alere ID NOW Strep A	15	*Negative	15	100.00%
Medline Strep A Dipstick	4	*Negative	4	100.00%
Roche Cobas Liat Strep A	11	*Negative	11	100.00%
Henry Schein OneStep+ Rapid Strep Dip	70	Positive	2	2.86%
PSS Rapid Test Strip Strep A	2	Positive	1	50.00%
Sekisui Diagnostics OSOM Ultra	7	Positive	1	14.29%
Consult Diag Strep A Dipstick	101	Positive	3	2.97%
RS-2		* Denotes Correct Response		
AG: All Methods	145			
		*Positive	139	95.86%
		Negative	6	4.14%
Abbott Signify Strep A - Waived	1	*Positive	1	100.00%
Alere Acceava Strep A Dipstick	7	*Positive	7	100.00%
Fisher HC Sure-Vue Strep A (Waived)	1	*Positive	1	100.00%
Sekisui Diagnostics OSOM	8	*Positive	8	100.00%
Sekisui Diagnostics OSOM Ultra Strep	24	*Positive	24	100.00%
Henry Schein OneStep + Rapid Strep Cass	1	*Positive	1	100.00%
Henry Schein OneStep+ Rapid Strep Dip	12	*Positive	10	83.33%
Innovacon One Step Strep Ag	1	*Positive	1	100.00%
McKesson Strep A Dip	6	*Positive	6	100.00%
Quidel QuickVue Dipstick Strep	14	*Positive	13	92.86%
Quidel QuickVue In-Line - 25 Test	14	*Positive	14	100.00%
Quidel QuickVue+ Strep A	4	*Positive	4	100.00%
SP Brand (Cardinal) Rapid Strep A Dipstick	13	*Positive	13	100.00%
Sekisui Diagnostics OSOM Ultra	1	*Positive	1	100.00%
Consult Diag Strep A Dipstick	13	*Positive	12	92.31%
MooreBrand Strep A Dipstick	2	*Positive	1	50.00%
BD Veritor Strep A	7	*Positive	7	100.00%
Alere ID NOW Strep A	4	*Positive	3	75.00%
Medline Strep A Dipstick	2	*Positive	2	100.00%



Roche Cobas Liat Strep A	10	*Positive	10	100.00%
Henry Schein OneStep+ Rapid Strep Dip	12	Negative	2	16.67%
Quidel QuickVue Dipstick Strep	14	Negative	1	7.14%
Consult Diag Strep A Dipstick	13	Negative	1	7.69%
MooreBrand Strep A Dipstick	2	Negative	1	50.00%
Alere ID NOW Strep A	4	Negative	1	25.00%
RS-3		* Denotes Correct Response		
AG: All Methods	3			
		*Positive	3	100.00%
Quidel QuickVue In-Line - 25 Test	1	*Positive	1	100.00%
Medline Strep A Dipstick	2	*Positive	2	100.00%
RS-4		* Denotes Correct Response		
AG: All Methods	3			
		*Positive	3	100.00%
Quidel QuickVue In-Line - 25 Test	1	*Positive	1	100.00%
Medline Strep A Dipstick	2	*Positive	2	100.00%
RS-5		* Denotes Correct Response		
AG: All Methods	3			
		*Negative	3	100.00%
Quidel QuickVue In-Line - 25 Test	1	*Negative	1	100.00%
Medline Strep A Dipstick	2	*Negative	2	100.00%
		Strep A Aq Det. (International)		
RS-1		* Denotes Correct Response		
AG: All Methods	1			
		*Negative	1	100.00%
Alere Clearview Exact II Cas	1	*Negative	1	100.00%
RS-2		* Denotes Correct Response		
AG: All Methods	1			
		*Positive	1	100.00%
Alere Clearview Exact II Cas	1	*Positive	1	100.00%
RS-3		* Denotes Correct Response		
AG: All Methods	1			
		*Positive	1	100.00%
Alere Clearview Exact II Cas	1	*Positive	1	100.00%
RS-4		* Denotes Correct Response		
AG: All Methods	1			
		*Positive	1	100.00%
Alere Clearview Exact II Cas	1	*Positive	1	100.00%
RS-5		* Denotes Correct Response		
AG: All Methods	1			
		*Negative	1	100.00%
Alere Clearview Exact II Cas	1	*Negative	1	100.00%
739 - Sickle Cell Screen	Labs	Response	Number	Percent
		Sickle Cell Screen		
SS-1		* Denotes Correct Response		
AG: All Methods	12			
		*Negative	12	100.00%
Pacific Hemostasis Sickle Screen	3	*Negative	3	100.00%
Streck SickieDex	9	*Negative	9	100.00%
679 - Syphilis Serology	Labs	Response	Number	Percent
		RPR		
SY-1		* Denotes Correct Response		
AG: All Methods	10			
		*Reactive	10	100.00%
ASI RPR	2	*Reactive	2	100.00%
Becton Dickinson BBL	4	*Reactive	4	100.00%
New Horizons	1	*Reactive	1	100.00%
Stanbio Quicktest	1	*Reactive	1	100.00%
Wampole Impact	1	*Reactive	1	100.00%
Intl ASI RPR	1	*Reactive	1	100.00%
SY-2		* Denotes Correct Response		
AG: All Methods	10			
		*Nonreactive	10	100.00%
ASI RPR	2	*Nonreactive	2	100.00%
Becton Dickinson BBL	4	*Nonreactive	4	100.00%
New Horizons	1	*Nonreactive	1	100.00%
Stanbio Quicktest	1	*Nonreactive	1	100.00%
Wampole Impact	1	*Nonreactive	1	100.00%



Intl ASI RPR	1	*Nonreactive	1	100.00%
SY-3		* Denotes Correct Response		
AG: All Methods	10			
		*Reactive	10	100.00%
ASI RPR	2	*Reactive	2	100.00%
Becton Dickinson BBL	4	*Reactive	4	100.00%
New Horizons	1	*Reactive	1	100.00%
Stanbio Quicktest	1	*Reactive	1	100.00%
Wampole Impact	1	*Reactive	1	100.00%
Intl ASI RPR	1	*Reactive	1	100.00%
SY-4		* Denotes Correct Response		
AG: All Methods	10			
		*Reactive	9	90.00%
		*Weakly Reactive	1	10.00%
ASI RPR	2	*Reactive	2	100.00%
Becton Dickinson BBL	4	*Reactive	4	100.00%
New Horizons	1	*Reactive	1	100.00%
Stanbio Quicktest	1	*Reactive	1	100.00%
Intl ASI RPR	1	*Reactive	1	100.00%
Wampole Impact	1	*Weakly Reactive	1	100.00%
SY-5		* Denotes Correct Response		
AG: All Methods	10			
		*Nonreactive	10	100.00%
ASI RPR	2	*Nonreactive	2	100.00%
Becton Dickinson BBL	4	*Nonreactive	4	100.00%
New Horizons	1	*Nonreactive	1	100.00%
Stanbio Quicktest	1	*Nonreactive	1	100.00%
Wampole Impact	1	*Nonreactive	1	100.00%
Intl ASI RPR	1	*Nonreactive	1	100.00%
		RPR/VDRL (titer)		
SY-1		* Denotes Correct Response		
AG: All Methods	3			
		*1:1	1	33.33%
		*1:2	1	33.33%
		*1:4	1	33.33%
New Horizons	1	*1:1	1	100.00%
Becton Dickinson BBL	1	*1:2	1	100.00%
ASI RPR	1	*1:4	1	100.00%
SY-2		* Denotes Correct Response		
AG: All Methods	3			
		*Nonreactive	3	100.00%
ASI RPR	1	*Nonreactive	1	100.00%
Becton Dickinson BBL	1	*Nonreactive	1	100.00%
New Horizons	1	*Nonreactive	1	100.00%
SY-3		* Denotes Correct Response		
AG: All Methods	3			
		*1:2	1	33.33%
		*1:4	2	66.67%
New Horizons	1	*1:2	1	100.00%
ASI RPR	1	*1:4	1	100.00%
Becton Dickinson BBL	1	*1:4	1	100.00%
SY-4		* Denotes Correct Response		
AG: All Methods	3			
		*1:1	1	33.33%
		*1:2	2	66.67%
New Horizons	1	*1:1	1	100.00%
ASI RPR	1	*1:2	1	100.00%
Becton Dickinson BBL	1	*1:2	1	100.00%
SY-5		* Denotes Correct Response		
AG: All Methods	3			
		*Nonreactive	3	100.00%
ASI RPR	1	*Nonreactive	1	100.00%
Becton Dickinson BBL	1	*Nonreactive	1	100.00%
New Horizons	1	*Nonreactive	1	100.00%
		RPR (International)		
SY-1		* Denotes Correct Response		
AG: All Methods	3			
		*Reactive	3	100.00%



BD Macro-Vue RPR	1	*Reactive	1	100.00%
Fisher HC Sure-Vue RPR	1	*Reactive	1	100.00%
Wampole Impact	1	*Reactive	1	100.00%
SY-2		* Denotes Correct Response		
AG: All Methods	3			
		*Nonreactive	3	100.00%
BD Macro-Vue RPR	1	*Nonreactive	1	100.00%
Fisher HC Sure-Vue RPR	1	*Nonreactive	1	100.00%
Wampole Impact	1	*Nonreactive	1	100.00%
SY-3		* Denotes Correct Response		
AG: All Methods	3			
		*Reactive	3	100.00%
BD Macro-Vue RPR	1	*Reactive	1	100.00%
Fisher HC Sure-Vue RPR	1	*Reactive	1	100.00%
Wampole Impact	1	*Reactive	1	100.00%
SY-4		* Denotes Correct Response		
AG: All Methods	3			
		*Reactive	3	100.00%
BD Macro-Vue RPR	1	*Reactive	1	100.00%
Fisher HC Sure-Vue RPR	1	*Reactive	1	100.00%
Wampole Impact	1	*Reactive	1	100.00%
SY-5		* Denotes Correct Response		
AG: All Methods	3			
		*Nonreactive	3	100.00%
BD Macro-Vue RPR	1	*Nonreactive	1	100.00%
Fisher HC Sure-Vue RPR	1	*Nonreactive	1	100.00%
Wampole Impact	1	*Nonreactive	1	100.00%
776 - Throat Culture;780 - Throat & Genital Cultures	Labs	Response	Number	Percent
		Throat Culture		
TC-1		* Denotes Correct Response		
AG: All Methods	61			
		*Beta hemolysis, referred for ID/confirmation	2	3.28%
		*Growth, referred for identification	1	1.64%
		*Positive for Group A Streptococcus	31	50.82%
		*Presumptive Positive for Group A Streptococcus	22	36.07%
		*Presumptive Streptococcus species	2	3.28%
		*Streptococcus pyogenes	1	1.64%
		Negative for Group A Streptococcus	1	1.64%
		Streptococcus, beta-hemolytic, not Group A	1	1.64%
Cultures	61	*Beta hemolysis, referred for ID/confirmation	2	3.28%
Cultures		*Growth, referred for identification	1	1.64%
Cultures		*Positive for Group A Streptococcus	31	50.82%
Cultures		*Presumptive Positive for Group A Streptococcus	22	36.07%
Cultures		*Presumptive Streptococcus species	2	3.28%
Cultures		*Streptococcus pyogenes	1	1.64%
Cultures		Negative for Group A Streptococcus	1	1.64%
Cultures		Streptococcus, beta-hemolytic, not Group A	1	1.64%
TC-2		* Denotes Correct Response		
AG: All Methods	60			
		*Beta hemolysis, referred for ID/confirmation	2	3.33%
		*Negative for Group A Streptococcus	50	83.33%
		*Presumptive Streptococcus species	1	1.67%
		*Streptococcus, beta-hemolytic, not Group A	3	5.00%
		Positive for Group A Streptococcus	3	5.00%
		Presumptive Positive for Group A Streptococcus	1	1.67%
Cultures	60	*Beta hemolysis, referred for ID/confirmation	2	3.33%
Cultures		*Negative for Group A Streptococcus	50	83.33%
Cultures		*Presumptive Streptococcus species	1	1.67%
Cultures		*Streptococcus, beta-hemolytic, not Group A	3	5.00%
Cultures		Positive for Group A Streptococcus	3	5.00%
Cultures		Presumptive Positive for Group A Streptococcus	1	1.67%
TC-3		* Denotes Correct Response		
AG: All Methods	62			
		*Beta hemolysis, referred for ID/confirmation	3	4.84%
		*Positive for Group A Streptococcus	30	48.39%
		*Presumptive Positive for Group A Streptococcus	24	38.71%
		*Presumptive Streptococcus species	2	3.23%
		*Streptococcus pyogenes	1	1.61%



		Negative for Group A Streptococcus	2	3.23%
Cultures	62	*Beta hemolysis, referred for ID/confirmation	3	4.84%
Cultures		*Positive for Group A Streptococcus	30	48.39%
Cultures		*Presumptive Positive for Group A Streptococcus	24	38.71%
Cultures		*Presumptive Streptococcus species	2	3.23%
Cultures		*Streptococcus pyogenes	1	1.61%
Cultures		Negative for Group A Streptococcus	2	3.23%
TC-4		* Denotes Correct Response		
AG: All Methods	76			
		*Beta hemolysis, referred for ID/confirmation	6	7.89%
		*Growth, referred for identification	1	1.32%
		*Positive for Group A Streptococcus	40	52.63%
		*Presumptive Positive for Group A Streptococcus	24	31.58%
		*Presumptive Streptococcus species	2	2.63%
		*Streptococcus pyogenes	1	1.32%
		*Streptococcus, beta-hemolytic, Group A	1	1.32%
		Negative for Group A Streptococcus	1	1.32%
Cultures	76	*Beta hemolysis, referred for ID/confirmation	6	7.89%
Cultures		*Growth, referred for identification	1	1.32%
Cultures		*Positive for Group A Streptococcus	40	52.63%
Cultures		*Presumptive Positive for Group A Streptococcus	24	31.58%
Cultures		*Presumptive Streptococcus species	2	2.63%
Cultures		*Streptococcus pyogenes	1	1.32%
Cultures		*Streptococcus, beta-hemolytic, Group A	1	1.32%
Cultures		Negative for Group A Streptococcus	1	1.32%
TC-5		* Denotes Correct Response		
AG: All Methods	76			
		*Negative for Group A Streptococcus	59	77.63%
		*No beta hemolysis	3	3.95%
		*No Growth - Strep Selective medium	11	14.47%
		Positive for Group A Streptococcus	3	3.95%
Cultures	76	*Negative for Group A Streptococcus	59	77.63%
Cultures		*No beta hemolysis	3	3.95%
Cultures		*No Growth - Strep Selective medium	11	14.47%
Cultures		Positive for Group A Streptococcus	3	3.95%
640 - Urinalysis;642 - Urinalysis/hCG;639 - Complete Urinalysis Package;601 - Waived/PPM Tests II;609 - Waived/PPM Tests I;647 - Mini-Waived Package	Labs	Response	Number	Percent
		Urine Sp Gravity		
UA-1		* Denotes Correct Response		
PG: Siemens Clinitek 10 / 100	12			
		*1.015	4	33.33%
		*1.020	4	33.33%
		*1.025	3	25.00%
		1.010	1	8.33%
UA-1		* Denotes Correct Response		
PG: Siemens Clinitek 50	69			
		*1.015	10	14.49%
		*1.020	33	47.83%
		*1.025	26	37.68%
UA-1		* Denotes Correct Response		
PG: Siemens Clinitek 500	5			
		*1.020	3	60.00%
		*1.025	2	40.00%
UA-1		* Denotes Correct Response		
PG: Siemens Clinitek Status	660			
		*1.020	190	28.79%
		*1.025	448	67.88%
		*1.030	16	2.42%
		1.010	1	0.15%
		1.015	5	0.76%
UA-1		* Denotes Correct Response		
PG: Siemens Clinitek Status / Siemens - Plus Sys / Siemens Urine Dipstick	13			
		*1.020	2	15.38%
		*1.025	11	84.62%
UA-1		* Denotes Correct Response		
PG: Siemens Multistix - Visual	207			
		*1.010	16	7.73%



		*1.015	97	46.86%
		*1.020	48	23.19%
		*1.025	41	19.81%
		1.000	3	1.45%
		1.005	1	0.48%
		1.030	1	0.48%
UA-1		* Denotes Correct Response		
PG: Henry Schein Urispec	5			
		*1.015	5	100.00%
UA-1		* Denotes Correct Response		
PG: Consult Diagnostics Urine Reagent Strips	5			
		*1.015	2	40.00%
		*1.020	3	60.00%
UA-1		* Denotes Correct Response		
PG: Diag Test Group Clarity Urochek	25			
		*1.015	4	16.00%
		*1.020	15	60.00%
		*1.025	6	24.00%
UA-1		* Denotes Correct Response		
PG: Roche Chemstrips - Visual	84			
		*1.010	24	28.57%
		*1.015	53	63.10%
		1.005	3	3.57%
		1.020	1	1.19%
		1.025	3	3.57%
UA-1		* Denotes Correct Response		
PG: Roche Urisys / Roche - Conventional / Roche Urine Dipstick	2			
		*1.010	1	50.00%
		*1.015	1	50.00%
UA-1		* Denotes Correct Response		
PG: UriScan Reagent Strips	6			
		*1.010	1	16.67%
		*1.015	3	50.00%
		*1.020	1	16.67%
		1.025	1	16.67%
UA-1		* Denotes Correct Response		
PG: Medi-Lab Performance Urine Strip - Visual	8			
		*1.015	2	25.00%
		*1.020	5	62.50%
		*1.025	1	12.50%
UA-1		* Denotes Correct Response		
PG: Siemens Clinitek Advantus	23			
		*1.015	1	4.35%
		*1.020	20	86.96%
		*1.025	2	8.70%
UA-1		* Denotes Correct Response		
PG: Consult Diag Urinalysis Test Strips	7			
		*1.015	2	28.57%
		*1.020	4	57.14%
		*1.025	1	14.29%
UA-1		* Denotes Correct Response		
PG: Consult Diag Urine Analyzer	8			
		*1.015	1	12.50%
		*1.020	6	75.00%
		*1.025	1	12.50%
UA-1		* Denotes Correct Response		
PG: 120 Urine Analyzer	43			
		*1.015	8	18.60%
		*1.020	31	72.09%
		*1.025	4	9.30%
UA-1		* Denotes Correct Response		
PG: McKesson 10 SG Urine Dipstick	5			
		*1.010	1	20.00%
		*1.015	2	40.00%
		*1.020	2	40.00%
UA-1		* Denotes Correct Response		
SG: Select Siemens Clinitek/Visual	977			
		*1.010	18	1.84%



		*1.015	117	11.98%
		*1.020	298	30.50%
		*1.025	523	53.53%
		*1.030	17	1.74%
		1.000	3	0.31%
		1.005	1	0.10%
		1.010	18	1.84%
		1.015	117	11.98%
		1.030	17	1.74%
UA-1		* Denotes Correct Response		
SG: Select Siemens Clinitek/Visual / Siemens - Plus Sys / Siemens Urine Dipstick	14			
		*1.020	3	21.43%
		*1.025	11	78.57%
UA-1		* Denotes Correct Response		
SG: Select Roche Urinalysis Inst/Visual	93			
		*1.010	25	26.88%
		*1.015	59	63.44%
		*1.020	2	2.15%
		*1.025	4	4.30%
		1.005	3	3.23%
		1.020	2	2.15%
		1.025	4	4.30%
UA-1		* Denotes Correct Response		
SG: Select Consult Diag Urinalysis Inst/Visual	20			
		*1.015	5	25.00%
		*1.020	13	65.00%
		*1.025	2	10.00%
UA-1		* Denotes Correct Response		
SG: Select McKesson Urinalysis Inst/Visual	48			
		*1.010	1	2.08%
		*1.015	10	20.83%
		*1.020	33	68.75%
		*1.025	4	8.33%
UA-1		* Denotes Correct Response		
SG: Select Henry Schein Urinalysis Inst/Visual	8			
		*1.015	6	75.00%
		*1.025	2	25.00%
UA-1		* Denotes Correct Response		
SG: Select Medi-Lab Performance UA Inst/Visual	9			
		*1.015	3	33.33%
		*1.020	5	55.56%
		*1.025	1	11.11%
UA-1		* Denotes Correct Response		
MG: Method Group	1206			
		*1.010	48	3.98%
		*1.015	213	17.66%
		*1.020	374	31.01%
		*1.025	547	45.36%
		*1.030	17	1.41%
		1.000	3	0.25%
		1.005	4	0.33%
		1.010	48	3.98%
		1.015	213	17.66%
		1.020	374	31.01%
		1.025	547	45.36%
		1.030	17	1.41%
UA-1		* Denotes Correct Response		
MG: Method Group	2			
		*1.010	1	50.00%
		*1.015	1	50.00%
UA-1		* Denotes Correct Response		
MG: Method Group	14			
		*1.020	3	21.43%
		*1.025	11	78.57%
UA-1		* Denotes Correct Response		
MG: Method Group	0			
		Urine Bilirubin		



UA-1		* Denotes Correct Response		
PG: Siemens Clinitek 10 / 100	12			
		*Negative	12	100.00%
UA-1		* Denotes Correct Response		
PG: Siemens Clinitek 50	69			
		*Negative	69	100.00%
UA-1		* Denotes Correct Response		
PG: Siemens Clinitek 500	5			
		*Negative	5	100.00%
UA-1		* Denotes Correct Response		
PG: Siemens Clinitek Status	654			
		*Negative	654	100.00%
UA-1		* Denotes Correct Response		
PG: Siemens Clinitek Status / Siemens - Plus Sys / Siemens Urine Dipstick	13			
		*Negative	13	100.00%
UA-1		* Denotes Correct Response		
PG: Siemens Multistix - Visual	204			
		*Negative	199	97.55%
		Large (3+)	1	0.49%
		Moderate (2+)	2	0.98%
		Trace/Small (1+)	2	0.98%
UA-1		* Denotes Correct Response		
PG: Germaine Airstick	4			
		*Negative	4	100.00%
UA-1		* Denotes Correct Response		
PG: Henry Schein Urispec	5			
		*Negative	5	100.00%
UA-1		* Denotes Correct Response		
PG: Consult Diagnostics Urine Reagent Strips	5			
		*Negative	5	100.00%
UA-1		* Denotes Correct Response		
PG: Diag Test Group Clarity Urochek	24			
		*Negative	24	100.00%
UA-1		* Denotes Correct Response		
PG: Roche Chemstrips - Visual	85			
		*Negative	85	100.00%
UA-1		* Denotes Correct Response		
PG: Roche Urisys / Roche - Conventional / Roche Urine Dipstick	2			
		*Negative	2	100.00%
UA-1		* Denotes Correct Response		
PG: UriScan Reagent Strips	7			
		*Negative	7	100.00%
UA-1		* Denotes Correct Response		
PG: Medi-Lab Performance Urine Strip - Visual	8			
		*Negative	8	100.00%
UA-1		* Denotes Correct Response		
PG: Siemens Clinitek Advantus	23			
		*Negative	23	100.00%
UA-1		* Denotes Correct Response		
PG: Consult Diag Urinalysis Test Strips	7			
		*Negative	7	100.00%
UA-1		* Denotes Correct Response		
PG: Consult Diag Urine Analyzer	8			
		*Negative	8	100.00%
UA-1		* Denotes Correct Response		
PG: 120 Urine Analyzer	43			
		*Negative	43	100.00%
UA-1		* Denotes Correct Response		
PG: McKesson 10 SG Urine Dipstick	5			
		*Negative	5	100.00%
UA-1		* Denotes Correct Response		
SG: Select Siemens Clinitek/Visual	968			
		*Negative	963	99.48%
		Large (3+)	1	0.10%
		Moderate (2+)	2	0.21%
		Trace/Small (1+)	2	0.21%
UA-1		* Denotes Correct Response		



SG: Select Siemens Clinitek/Visual / Siemens - Plus Sys / Siemens Urine Dipstick	14			
		*Negative	14	100.00%
UA-1		* Denotes Correct Response		
SG: Select Roche Urinalysis Inst/Visual	94			
		*Negative	94	100.00%
UA-1		* Denotes Correct Response		
SG: Select Consult Diag Urinalysis Inst/Visual	20			
		*Negative	20	100.00%
UA-1		* Denotes Correct Response		
SG: Select McKesson Urinalysis Inst/Visual	48			
		*Negative	48	100.00%
UA-1		* Denotes Correct Response		
SG: Select Henry Schein Urinalysis Inst/Visual	8			
		*Negative	8	100.00%
UA-1		* Denotes Correct Response		
SG: Select Medi-Lab Performance UA Inst/Visual	9			
		*Negative	9	100.00%
UA-1		* Denotes Correct Response		
MG: Method Group	1197			
		*Negative	1191	99.50%
		Large (3+)	1	0.08%
		Moderate (2+)	2	0.17%
		Trace/Small (1+)	3	0.25%
UA-1		* Denotes Correct Response		
MG: Method Group	2			
		*Negative	2	100.00%
UA-1		* Denotes Correct Response		
MG: Method Group	14			
		*Negative	14	100.00%
UA-1		* Denotes Correct Response		
MG: Method Group	0			
		Urine Ketones		
UA-1		* Denotes Correct Response		
PG: Siemens Clinitek 10 / 100	12			
		*Negative	12	100.00%
UA-1		* Denotes Correct Response		
PG: Siemens Clinitek 50	69			
		*Negative	69	100.00%
UA-1		* Denotes Correct Response		
PG: Siemens Clinitek 500	5			
		*Negative	5	100.00%
UA-1		* Denotes Correct Response		
PG: Siemens Clinitek Status	659			
		*Negative	659	100.00%
UA-1		* Denotes Correct Response		
PG: Siemens Clinitek Status / Siemens - Plus Sys / Siemens Urine Dipstick	13			
		*Negative	13	100.00%
UA-1		* Denotes Correct Response		
PG: Siemens Multistix - Visual	209			
		*Negative	202	96.65%
		Large (3+) / 80 mg/dL	1	0.48%
		Small (1+) / 15mg/dL	3	1.44%
		Trace/5 mg/dL	3	1.44%
UA-1		* Denotes Correct Response		
PG: Henry Schein Urispec	5			
		*Negative	5	100.00%
UA-1		* Denotes Correct Response		
PG: Consult Diagnostics Urine Reagent Strips	5			
		*Negative	5	100.00%
UA-1		* Denotes Correct Response		
PG: Diag Test Group Clarity Urochek	25			
		*Negative	25	100.00%
UA-1		* Denotes Correct Response		
PG: Roche Chemstrips - Visual	86			
		*Negative	85	98.84%
		Small (1+) / 15mg/dL	1	1.16%



UA-1		* Denotes Correct Response		
PG: Roche Urisys / Roche - Conventional / Roche Urine Dipstick	2			
		*Negative	2	100.00%
UA-1		* Denotes Correct Response		
PG: Medi-Lab Performance Urine Strip - Visual	8			
		*Negative	8	100.00%
UA-1		* Denotes Correct Response		
PG: Siemens Clinitek Advantus	23			
		*Negative	23	100.00%
UA-1		* Denotes Correct Response		
PG: Consult Diag Urinalysis Test Strips	7			
		*Negative	7	100.00%
UA-1		* Denotes Correct Response		
PG: Consult Diag Urine Analyzer	8			
		*Negative	8	100.00%
UA-1		* Denotes Correct Response		
PG: 120 Urine Analyzer	43			
		*Negative	43	100.00%
UA-1		* Denotes Correct Response		
PG: McKesson 10 SG Urine Dipstick	5			
		*Negative	5	100.00%
UA-1		* Denotes Correct Response		
SG: Select Siemens Clinitek/Visual	978			
		*Negative	971	99.28%
		Large (3+) / 80 mg/dL	1	0.10%
		Small (1+) / 15mg/dL	3	0.31%
		Trace/5 mg/dL	3	0.31%
UA-1		* Denotes Correct Response		
SG: Select Siemens Clinitek/Visual / Siemens - Plus Sys / Siemens Urine Dipstick	14			
		*Negative	14	100.00%
UA-1		* Denotes Correct Response		
SG: Select Roche Urinalysis Inst/Visual	95			
		*Negative	94	98.95%
		Small (1+) / 15mg/dL	1	1.05%
UA-1		* Denotes Correct Response		
SG: Select Consult Diag Urinalysis Inst/Visual	20			
		*Negative	20	100.00%
UA-1		* Denotes Correct Response		
SG: Select McKesson Urinalysis Inst/Visual	48			
		*Negative	48	100.00%
UA-1		* Denotes Correct Response		
SG: Select Henry Schein Urinalysis Inst/Visual	8			
		*Negative	8	100.00%
UA-1		* Denotes Correct Response		
SG: Select Medi-Lab Performance UA Inst/Visual	9			
		*Negative	9	100.00%
UA-1		* Denotes Correct Response		
MG: Method Group	1210			
		*Negative	1200	99.17%
		Large (3+) / 80 mg/dL	1	0.08%
		Moderate (2+) / 40 mg/dL	1	0.08%
		Small (1+) / 15mg/dL	4	0.33%
		Trace/5 mg/dL	4	0.33%
UA-1		* Denotes Correct Response		
MG: Method Group	2			
		*Negative	2	100.00%
UA-1		* Denotes Correct Response		
MG: Method Group	14			
		*Negative	14	100.00%
UA-1		* Denotes Correct Response		
MG: Method Group	0			
		Urine pH		
UA-1		* Denotes Correct Response		
PG: Siemens Clinitek 10 / 100	12			
		*7	9	75.00%
		*7.5	3	25.00%
UA-1		* Denotes Correct Response		



PG: Siemens Clinitek 50	69			
		*7	69	100.00%
UA-1		* Denotes Correct Response		
PG: Siemens Clinitek 500	5			
		*7	5	100.00%
UA-1		* Denotes Correct Response		
PG: Siemens Clinitek Status	660			
		*7	653	98.94%
		*7.5	5	0.76%
		6	1	0.15%
		6.5	1	0.15%
UA-1		* Denotes Correct Response		
PG: Siemens Clinitek Status / Siemens - Plus Sys / Siemens Urine Dipstick	13			
		*7	13	100.00%
UA-1		* Denotes Correct Response		
PG: Siemens Multistix - Visual	208			
		*7	106	50.96%
		*7.5	96	46.15%
		6.5	5	2.40%
		8 or greater	1	0.48%
UA-1		* Denotes Correct Response		
PG: Henry Schein Urispec	5			
		*7	5	100.00%
UA-1		* Denotes Correct Response		
PG: Consult Diagnostics Urine Reagent Strips	5			
		*6.5	1	20.00%
		*7	3	60.00%
		6	1	20.00%
UA-1		* Denotes Correct Response		
PG: Diag Test Group Clarity Urochek	25			
		*6.5	3	12.00%
		*7	22	88.00%
UA-1		* Denotes Correct Response		
PG: Roche Chemstrips - Visual	85			
		*7	81	95.29%
		6	1	1.18%
		6.5	1	1.18%
		7.5	1	1.18%
		8 or greater	1	1.18%
UA-1		* Denotes Correct Response		
PG: Roche Urisys / Roche - Conventional / Roche Urine Dipstick	2			
		*7	2	100.00%
UA-1		* Denotes Correct Response		
PG: UriScan Reagent Strips	8			
		*6.5	1	12.50%
		*7	6	75.00%
		*7.5	1	12.50%
UA-1		* Denotes Correct Response		
PG: Siemens Clinitek Advantus	23			
		*7	15	65.22%
		*7.5	8	34.78%
UA-1		* Denotes Correct Response		
PG: Consult Diag Urinalysis Test Strips	7			
		*6.5	2	28.57%
		*7	3	42.86%
		*7.5	1	14.29%
		6	1	14.29%
UA-1		* Denotes Correct Response		
SG: Select Siemens Clinitek/Visual	979			
		*7	859	87.74%
		*7.5	112	11.44%
		6	1	0.10%
		6.5	6	0.61%
		8 or greater	1	0.10%
UA-1		* Denotes Correct Response		
SG: Select Siemens Clinitek/Visual / Siemens - Plus Sys / Siemens Urine Dipstick	14			
		*7	14	100.00%



UA-1		* Denotes Correct Response		
SG: Select Roche Urinalysis Inst/Visual	93			
		*7	89	95.70%
		6	1	1.08%
		6.5	1	1.08%
		7.5	1	1.08%
		8 or greater	1	1.08%
UA-1		* Denotes Correct Response		
SG: Select McKesson Urinalysis Inst/Visual	48			
		*6.5	2	4.17%
		*7	44	91.67%
		*7.5	2	4.17%
UA-1		* Denotes Correct Response		
SG: Select Henry Schein Urinalysis Inst/Visual	8			
		*7	8	100.00%
UA-1		* Denotes Correct Response		
SG: Select Medi-Lab Performance UA Inst/Visual	9			
		*6.5	1	11.11%
		*7	8	88.89%
UA-1		* Denotes Correct Response		
MG: Method Group	1210			
		*6.5	19	1.57%
		*7	1068	88.26%
		*7.5	117	9.67%
		6	4	0.33%
		6.5	19	1.57%
		7.5	117	9.67%
		8 or greater	2	0.17%
UA-1		* Denotes Correct Response		
MG: Method Group	2			
		*7	2	100.00%
UA-1		* Denotes Correct Response		
MG: Method Group	14			
		*7	14	100.00%
UA-1		* Denotes Correct Response		
MG: Method Group	0			
UA-1		Urine Protein - Qual		
UA-1		* Denotes Correct Response		
PG: Siemens Clinitek 10 / 100	12			
		*100 (2+)	8	66.67%
		*300-500 (3+)	3	25.00%
		Negative	1	8.33%
UA-1		* Denotes Correct Response		
PG: Siemens Clinitek 50	69			
		*100 (2+)	27	39.13%
		*300-500 (3+)	42	60.87%
UA-1		* Denotes Correct Response		
PG: Siemens Clinitek 500	5			
		*100 (2+)	2	40.00%
		*300-500 (3+)	3	60.00%
UA-1		* Denotes Correct Response		
PG: Siemens Clinitek Status	660			
		*100 (2+)	253	38.33%
		*300-500 (3+)	405	61.36%
		30 (1+)	1	0.15%
		Negative	1	0.15%
UA-1		* Denotes Correct Response		
PG: Siemens Clinitek Status / Siemens - Plus Sys / Siemens Urine Dipstick	13			
		*2+	2	15.38%
		*3+	11	84.62%
UA-1		* Denotes Correct Response		
PG: Siemens Multistix - Visual	210			
		*100 (2+)	112	53.33%
		*300-500 (3+)	89	42.38%
		>/= 1000 (4+)	4	1.90%
		30 (1+)	3	1.43%
		Negative	2	0.95%



UA-1		* Denotes Correct Response		
PG: Henry Schein Urispec	5			
		*100 (2+)	4	80.00%
		*30 (1+)	1	20.00%
UA-1		* Denotes Correct Response		
PG: Consult Diagnostics Urine Reagent Strips	5			
		*100 (2+)	4	80.00%
		*300-500 (3+)	1	20.00%
UA-1		* Denotes Correct Response		
PG: Diag Test Group Clarity Urochek	25			
		*100 (2+)	14	56.00%
		*30 (1+)	3	12.00%
		*300-500 (3+)	8	32.00%
UA-1		* Denotes Correct Response		
PG: Roche Chemstrips - Visual	90			
		*100 (2+)	71	78.89%
		*300-500 (3+)	17	18.89%
		30 (1+)	2	2.22%
UA-1		* Denotes Correct Response		
PG: Roche Urisys / Roche - Conventional / Roche Urine Dipstick	2			
		*100 mg/dL	1	50.00%
		*30 mg/dL	1	50.00%
UA-1		* Denotes Correct Response		
PG: UriScan Reagent Strips	8			
		*100 (2+)	3	37.50%
		*30 (1+)	3	37.50%
		*300-500 (3+)	1	12.50%
		>/= 1000 (4+)	1	12.50%
UA-1		* Denotes Correct Response		
PG: Medi-Lab Performance Urine Strip - Visual	8			
		*100 (2+)	5	62.50%
		*300-500 (3+)	3	37.50%
UA-1		* Denotes Correct Response		
PG: Siemens Clinitek Advantus	23			
		*100 (2+)	22	95.65%
		Trace	1	4.35%
UA-1		* Denotes Correct Response		
PG: Consult Diag Urinalysis Test Strips	7			
		*100 (2+)	3	42.86%
		*300-500 (3+)	3	42.86%
		30 (1+)	1	14.29%
UA-1		* Denotes Correct Response		
PG: Consult Diag Urine Analyzer	8			
		*100 (2+)	6	75.00%
		*300-500 (3+)	2	25.00%
UA-1		* Denotes Correct Response		
PG: 120 Urine Analyzer	43			
		*100 (2+)	38	88.37%
		*30 (1+)	4	9.30%
		300-500 (3+)	1	2.33%
UA-1		* Denotes Correct Response		
PG: McKesson 10 SG Urine Dipstick	5			
		*100 (2+)	3	60.00%
		*300-500 (3+)	2	40.00%
UA-1		* Denotes Correct Response		
SG: Select Siemens Clinitek/Visual	984			
		*100 (2+)	429	43.60%
		*300-500 (3+)	542	55.08%
		>/= 1000 (4+)	4	0.41%
		30 (1+)	4	0.41%
		Negative	4	0.41%
		Trace	1	0.10%
UA-1		* Denotes Correct Response		
SG: Select Siemens Clinitek/Visual / Siemens - Plus Sys / Siemens Urine Dipstick	14			
		*2+	3	21.43%
		*3+	11	78.57%
UA-1		* Denotes Correct Response		



SG: Select Roche Urinalysis Inst/Visual	99			
		*100 (2+)	77	77.78%
		*30 (1+)	3	3.03%
		*300-500 (3+)	19	19.19%
		30 (1+)	3	3.03%
UA-1		* Denotes Correct Response		
SG: Select Consult Diag Urinalysis Inst/Visual	20			
		*100 (2+)	13	65.00%
		*300-500 (3+)	6	30.00%
		30 (1+)	1	5.00%
UA-1		* Denotes Correct Response		
SG: Select McKesson Urinalysis Inst/Visual	48			
		*100 (2+)	41	85.42%
		*30 (1+)	4	8.33%
		*300-500 (3+)	3	6.25%
		300-500 (3+)	3	6.25%
UA-1		* Denotes Correct Response		
SG: Select Henry Schein Urinalysis Inst/Visual	8			
		*100 (2+)	6	75.00%
		*30 (1+)	1	12.50%
		*300-500 (3+)	1	12.50%
UA-1		* Denotes Correct Response		
SG: Select Medi-Lab Performance UA Inst/Visual	9			
		*100 (2+)	6	66.67%
		*300-500 (3+)	3	33.33%
UA-1		* Denotes Correct Response		
MG: Method Group	1221			
		*100 (2+)	605	49.55%
		*30 (1+)	20	1.64%
		*300-500 (3+)	585	47.91%
		>/= 1000 (4+)	6	0.49%
		30 (1+)	20	1.64%
		300-500 (3+)	585	47.91%
		Negative	4	0.33%
		Trace	1	0.08%
UA-1		* Denotes Correct Response		
MG: Method Group	2			
		*100 mg/dL	1	50.00%
		*30 mg/dL	1	50.00%
UA-1		* Denotes Correct Response		
MG: Method Group	14			
		*2+	3	21.43%
		*3+	11	78.57%
UA-1		* Denotes Correct Response		
MG: Method Group	0			
		Urine Blood/Hemoglobin		
UA-1		* Denotes Correct Response		
PG: Siemens Clinitek 10 / 100	12			
		*Large (3+)	9	75.00%
		*Moderate (2+)	2	16.67%
		Small (1+)	1	8.33%
UA-1		* Denotes Correct Response		
PG: Siemens Clinitek 50	69			
		*Large (3+)	65	94.20%
		*Moderate (2+)	3	4.35%
		Negative	1	1.45%
UA-1		* Denotes Correct Response		
PG: Siemens Clinitek 500	5			
		*Large (3+)	5	100.00%
UA-1		* Denotes Correct Response		
PG: Siemens Clinitek Status	659			
		*Large (3+)	623	94.54%
		*Moderate (2+)	32	4.86%
		Negative	4	0.61%
UA-1		* Denotes Correct Response		
PG: Siemens Clinitek Status / Siemens - Plus Sys / Siemens Urine Dipstick	13			
		*3+	13	100.00%



UA-1		* Denotes Correct Response		
PG: Siemens Multistix - Visual	208			
		*Large (3+)	202	97.12%
		*Moderate (2+)	6	2.88%
UA-1		* Denotes Correct Response		
PG: Henry Schein Urispec	5			
		*Large (3+)	4	80.00%
		*Moderate (2+)	1	20.00%
UA-1		* Denotes Correct Response		
PG: Consult Diagnostics Urine Reagent Strips	5			
		*Large (3+)	5	100.00%
UA-1		* Denotes Correct Response		
PG: Diag Test Group Clarity Urochek	25			
		*Large (3+)	25	100.00%
UA-1		* Denotes Correct Response		
PG: Roche Chemstrips - Visual	91			
		*Large (3+)	85	93.41%
		*Moderate (2+)	5	5.49%
		Small (1+)	1	1.10%
UA-1		* Denotes Correct Response		
PG: Roche Urisys / Roche - Conventional / Roche Urine Dipstick	2			
		*250 Ery/uL	2	100.00%
UA-1		* Denotes Correct Response		
PG: Medi-Lab Performance Urine Strip - Visual	8			
		*Large (3+)	8	100.00%
UA-1		* Denotes Correct Response		
PG: Siemens Clinitek Advantus	23			
		*Large (3+)	23	100.00%
UA-1		* Denotes Correct Response		
PG: Consult Diag Urinalysis Test Strips	7			
		*Large (3+)	7	100.00%
UA-1		* Denotes Correct Response		
PG: Consult Diag Urine Analyzer	8			
		*Large (3+)	8	100.00%
UA-1		* Denotes Correct Response		
PG: 120 Urine Analyzer	43			
		*Large (3+)	42	97.67%
		*Moderate (2+)	1	2.33%
UA-1		* Denotes Correct Response		
PG: McKesson 10 SG Urine Dipstick	5			
		*Large (3+)	5	100.00%
UA-1		* Denotes Correct Response		
SG: Select Siemens Clinitek/Visual	978			
		*Large (3+)	929	94.99%
		*Moderate (2+)	43	4.40%
		Negative	5	0.51%
		Small (1+)	1	0.10%
UA-1		* Denotes Correct Response		
SG: Select Siemens Clinitek/Visual / Siemens - Plus Sys / Siemens Urine Dipstick	14			
		*3+	14	100.00%
UA-1		* Denotes Correct Response		
SG: Select Roche Urinalysis Inst/Visual	100			
		*Large (3+)	94	94.00%
		*Moderate (2+)	5	5.00%
		Small (1+)	1	1.00%
UA-1		* Denotes Correct Response		
SG: Select Consult Diag Urinalysis Inst/Visual	20			
		*Large (3+)	20	100.00%
UA-1		* Denotes Correct Response		
SG: Select McKesson Urinalysis Inst/Visual	48			
		*Large (3+)	47	97.92%
		*Moderate (2+)	1	2.08%
UA-1		* Denotes Correct Response		
SG: Select Henry Schein Urinalysis Inst/Visual	8			
		*Large (3+)	7	87.50%
		*Moderate (2+)	1	12.50%
UA-1		* Denotes Correct Response		



SG: Select Medi-Lab Performance UA Inst/Visual	9			
		*Large (3+)	9	100.00%
UA-1		* Denotes Correct Response		
MG: Method Group	1215			
		*Large (3+)	1154	94.98%
		*Moderate (2+)	52	4.28%
		Negative	6	0.49%
		Small (1+)	3	0.25%
UA-1		* Denotes Correct Response		
MG: Method Group	2			
		*250 Ery/uL	2	100.00%
UA-1		* Denotes Correct Response		
MG: Method Group	14			
		*3+	14	100.00%
UA-1		* Denotes Correct Response		
MG: Method Group	0			
		Urine Glucose (Reducing Sub)		
UA-1		* Denotes Correct Response		
PG: Siemens Clinitek 10 / 100	12			
		*Negative	12	100.00%
UA-1		* Denotes Correct Response		
PG: Siemens Clinitek 50	69			
		*Negative	69	100.00%
UA-1		* Denotes Correct Response		
PG: Siemens Clinitek 500	5			
		*Negative	5	100.00%
UA-1		* Denotes Correct Response		
PG: Siemens Clinitek Status	657			
		*Negative	656	99.85%
		<100 mg/dL	1	0.15%
UA-1		* Denotes Correct Response		
PG: Siemens Clinitek Status / Siemens - Plus Sys / Siemens Urine Dipstick	13			
		*Negative	13	100.00%
UA-1		* Denotes Correct Response		
PG: Siemens Multistix - Visual	210			
		*Negative	210	100.00%
UA-1		* Denotes Correct Response		
PG: Henry Schein Urispec	5			
		*Negative	5	100.00%
UA-1		* Denotes Correct Response		
PG: Consult Diagnostics Urine Reagent Strips	5			
		*Negative	5	100.00%
UA-1		* Denotes Correct Response		
PG: Diag Test Group Clarity Urochek	25			
		*Negative	25	100.00%
UA-1		* Denotes Correct Response		
PG: Roche Chemstrips - Visual	91			
		*Negative	91	100.00%
UA-1		* Denotes Correct Response		
PG: Roche Urisys / Roche - Conventional / Roche Urine Dipstick	2			
		*Normal	2	100.00%
UA-1		* Denotes Correct Response		
PG: UriScan Reagent Strips	8			
		*Negative	8	100.00%
UA-1		* Denotes Correct Response		
PG: Medi-Lab Performance Urine Strip - Visual	8			
		*Negative	8	100.00%
UA-1		* Denotes Correct Response		
PG: Siemens Clinitek Advantus	23			
		*Negative	23	100.00%
UA-1		* Denotes Correct Response		
PG: Consult Diag Urinalysis Test Strips	7			
		*Negative	7	100.00%
UA-1		* Denotes Correct Response		
PG: Consult Diag Urine Analyzer	8			
		*Negative	8	100.00%



UA-1		* Denotes Correct Response		
PG: 120 Urine Analyzer	42			
		*Negative	41	97.62%
		100 mg/dL	1	2.38%
UA-1		* Denotes Correct Response		
PG: McKesson 10 SG Urine Dipstick	5			
		*Negative	5	100.00%
UA-1		* Denotes Correct Response		
SG: Select Siemens Clinitek/Visual	981			
		*Negative	980	99.90%
		<100 mg/dL	1	0.10%
UA-1		* Denotes Correct Response		
SG: Select Siemens Clinitek/Visual / Siemens - Plus Sys / Siemens Urine Dipstick	14			
		*Negative	14	100.00%
UA-1		* Denotes Correct Response		
SG: Select Roche Urinalysis Inst/Visual	100			
		*Negative	100	100.00%
UA-1		* Denotes Correct Response		
SG: Select Consult Diag Urinalysis Inst/Visual	20			
		*Negative	20	100.00%
UA-1		* Denotes Correct Response		
SG: Select McKesson Urinalysis Inst/Visual	47			
		*Negative	46	97.87%
		100 mg/dL	1	2.13%
UA-1		* Denotes Correct Response		
SG: Select Henry Schein Urinalysis Inst/Visual	8			
		*Negative	8	100.00%
UA-1		* Denotes Correct Response		
SG: Select Medi-Lab Performance UA Inst/Visual	9			
		*Negative	9	100.00%
UA-1		* Denotes Correct Response		
MG: Method Group	1218			
		*Negative	1216	99.84%
		<100 mg/dL	1	0.08%
		100 mg/dL	1	0.08%
UA-1		* Denotes Correct Response		
MG: Method Group	2			
		*Normal	2	100.00%
UA-1		* Denotes Correct Response		
MG: Method Group	14			
		*Negative	14	100.00%
UA-1		* Denotes Correct Response		
MG: Method Group	0			
		Urine Leuk Esterase		
UA-1		* Denotes Correct Response		
PG: Siemens Clinitek 10 / 100	12			
		*Negative	12	100.00%
UA-1		* Denotes Correct Response		
PG: Siemens Clinitek 50	69			
		*Negative	69	100.00%
UA-1		* Denotes Correct Response		
PG: Siemens Clinitek 500	5			
		*Negative	5	100.00%
UA-1		* Denotes Correct Response		
PG: Siemens Clinitek Status	656			
		*Negative	655	99.85%
		Trace	1	0.15%
UA-1		* Denotes Correct Response		
PG: Siemens Clinitek Status / Siemens - Plus Sys / Siemens Urine Dipstick	13			
		*Negative	13	100.00%
UA-1		* Denotes Correct Response		
PG: Siemens Multistix - Visual	209			
		*Negative	209	100.00%
UA-1		* Denotes Correct Response		
PG: Henry Schein Urispec	5			
		*Negative	5	100.00%



UA-1		* Denotes Correct Response		
PG: Consult Diagnostics Urine Reagent Strips	5			
		*Negative	5	100.00%
UA-1		* Denotes Correct Response		
PG: Diag Test Group Clarity Urochek	25			
		*Negative	25	100.00%
UA-1		* Denotes Correct Response		
PG: Roche Chemstrips - Visual	91			
		*Negative	91	100.00%
UA-1		* Denotes Correct Response		
PG: Roche Urisys / Roche - Conventional / Roche Urine Dipstick	2			
		*Negative	2	100.00%
UA-1		* Denotes Correct Response		
PG: UriScan Reagent Strips	8			
		*Negative	7	87.50%
		Moderate (2+)	1	12.50%
UA-1		* Denotes Correct Response		
PG: Medi-Lab Performance Urine Strip - Visual	8			
		*Negative	8	100.00%
UA-1		* Denotes Correct Response		
PG: Siemens Clinitek Advantus	23			
		*Negative	23	100.00%
UA-1		* Denotes Correct Response		
PG: Consult Diag Urinalysis Test Strips	7			
		*Negative	7	100.00%
UA-1		* Denotes Correct Response		
PG: Consult Diag Urine Analyzer	8			
		*Negative	8	100.00%
UA-1		* Denotes Correct Response		
PG: 120 Urine Analyzer	43			
		*Negative	43	100.00%
UA-1		* Denotes Correct Response		
PG: McKesson 10 SG Urine Dipstick	5			
		*Negative	5	100.00%
UA-1		* Denotes Correct Response		
SG: Select Siemens Clinitek/Visual	976			
		*Negative	975	99.90%
		Trace	1	0.10%
UA-1		* Denotes Correct Response		
SG: Select Siemens Clinitek/Visual / Siemens - Plus Sys / Siemens Urine Dipstick	14			
		*Negative	14	100.00%
UA-1		* Denotes Correct Response		
SG: Select Roche Urinalysis Inst/Visual	100			
		*Negative	100	100.00%
UA-1		* Denotes Correct Response		
SG: Select Consult Diag Urinalysis Inst/Visual	20			
		*Negative	20	100.00%
UA-1		* Denotes Correct Response		
SG: Select McKesson Urinalysis Inst/Visual	48			
		*Negative	48	100.00%
UA-1		* Denotes Correct Response		
SG: Select Henry Schein Urinalysis Inst/Visual	8			
		*Negative	8	100.00%
UA-1		* Denotes Correct Response		
SG: Select Medi-Lab Performance UA Inst/Visual	9			
		*Negative	9	100.00%
UA-1		* Denotes Correct Response		
MG: Method Group	1213			
		*Negative	1211	99.84%
		Moderate (2+)	1	0.08%
		Trace	1	0.08%
UA-1		* Denotes Correct Response		
MG: Method Group	2			
		*Negative	2	100.00%
UA-1		* Denotes Correct Response		
MG: Method Group	14			
		*Negative	14	100.00%



UA-1		* Denotes Correct Response		
MG: Method Group	0			
		Urine Nitrite		
UA-1		* Denotes Correct Response		
PG: Siemens Clinitek 10 / 100	12			
		*Negative	12	100.00%
UA-1		* Denotes Correct Response		
PG: Siemens Clinitek 50	69			
		*Negative	68	98.55%
		Positive	1	1.45%
UA-1		* Denotes Correct Response		
PG: Siemens Clinitek 500	5			
		*Negative	5	100.00%
UA-1		* Denotes Correct Response		
PG: Siemens Clinitek Status	660			
		*Negative	657	99.55%
		Positive	3	0.45%
UA-1		* Denotes Correct Response		
PG: Siemens Clinitek Status / Siemens - Plus Sys / Siemens Urine Dipstick	13			
		*Negative	13	100.00%
UA-1		* Denotes Correct Response		
PG: Siemens Multistix - Visual	207			
		*Negative	207	100.00%
UA-1		* Denotes Correct Response		
PG: Henry Schein Urispec	5			
		*Negative	5	100.00%
UA-1		* Denotes Correct Response		
PG: Consult Diagnostics Urine Reagent Strips	5			
		*Negative	5	100.00%
UA-1		* Denotes Correct Response		
PG: Diag Test Group Clarity Urochek	24			
		*Negative	24	100.00%
UA-1		* Denotes Correct Response		
PG: Roche Chemstrips - Visual	91			
		*Negative	91	100.00%
UA-1		* Denotes Correct Response		
PG: Roche Urisys / Roche - Conventional / Roche Urine Dipstick	2			
		*Negative	2	100.00%
UA-1		* Denotes Correct Response		
PG: UriScan Reagent Strips	8			
		*Negative	7	87.50%
		Positive	1	12.50%
UA-1		* Denotes Correct Response		
PG: Medi-Lab Performance Urine Strip - Visual	8			
		*Negative	8	100.00%
UA-1		* Denotes Correct Response		
PG: Siemens Clinitek Advantus	23			
		*Negative	23	100.00%
UA-1		* Denotes Correct Response		
PG: Consult Diag Urinalysis Test Strips	7			
		*Negative	7	100.00%
UA-1		* Denotes Correct Response		
PG: Consult Diag Urine Analyzer	8			
		*Negative	8	100.00%
UA-1		* Denotes Correct Response		
PG: 120 Urine Analyzer	43			
		*Negative	43	100.00%
UA-1		* Denotes Correct Response		
PG: McKesson 10 SG Urine Dipstick	5			
		*Negative	5	100.00%
UA-1		* Denotes Correct Response		
SG: Select Siemens Clinitek/Visual	978			
		*Negative	974	99.59%
		Positive	4	0.41%
UA-1		* Denotes Correct Response		
SG: Select Siemens Clinitek/Visual / Siemens - Plus Sys / Siemens Urine Dipstick	14			



		*Negative	14	100.00%
UA-1		* Denotes Correct Response		
SG: Select Roche Urinalysis Inst/Visual	100			
		*Negative	99	99.00%
		Positive	1	1.00%
UA-1		* Denotes Correct Response		
SG: Select Consult Diag Urinalysis Inst/Visual	20			
		*Negative	20	100.00%
UA-1		* Denotes Correct Response		
SG: Select McKesson Urinalysis Inst/Visual	48			
		*Negative	48	100.00%
UA-1		* Denotes Correct Response		
SG: Select Henry Schein Urinalysis Inst/Visual	8			
		*Negative	8	100.00%
UA-1		* Denotes Correct Response		
SG: Select Medi-Lab Performance UA Inst/Visual	9			
		*Negative	9	100.00%
UA-1		* Denotes Correct Response		
MG: Method Group	1214			
		*Negative	1208	99.51%
		Positive	6	0.49%
UA-1		* Denotes Correct Response		
MG: Method Group	2			
		*Negative	2	100.00%
UA-1		* Denotes Correct Response		
MG: Method Group	14			
		*Negative	14	100.00%
UA-1		* Denotes Correct Response		
MG: Method Group	0			
		Urine Urobilinogen		
UA-1		* Denotes Correct Response		
PG: Siemens Clinitek 10 / 100	12			
		*0.2/Normal	12	100.00%
UA-1		* Denotes Correct Response		
PG: Siemens Clinitek 50	69			
		*0.2/Normal	69	100.00%
UA-1		* Denotes Correct Response		
PG: Siemens Clinitek 500	5			
		*0.2/Normal	5	100.00%
UA-1		* Denotes Correct Response		
PG: Siemens Clinitek Status	652			
		*0.2/Normal	647	99.23%
		2.0 mg/dL	5	0.77%
UA-1		* Denotes Correct Response		
PG: Siemens Clinitek Status / Siemens - Plus Sys / Siemens Urine Dipstick	6			
		*3.2 umol/L	6	100.00%
UA-1		* Denotes Correct Response		
PG: Siemens Clinitek Status / Siemens - Plus Sys / Siemens Urine Urobilinogen	6			
		*0.2 EU/dL	6	100.00%
UA-1		* Denotes Correct Response		
PG: Siemens Multistix - Visual	200			
		*0.2/Normal	196	98.00%
		1 mg/dL	2	1.00%
		2.0 mg/dL	2	1.00%
UA-1		* Denotes Correct Response		
PG: Consult Diagnostics Urine Reagent Strips	5			
		*0.2/Normal	5	100.00%
UA-1		* Denotes Correct Response		
PG: Diag Test Group Clarity Urochek	25			
		*0.2/Normal	24	96.00%
		2.0 mg/dL	1	4.00%
UA-1		* Denotes Correct Response		
PG: Roche Chemstrips - Visual	85			
		*0.2/Normal	85	100.00%
UA-1		* Denotes Correct Response		
PG: Roche Urisys / Roche - Conventional / Roche Urine Dipstick	2			



		*Normal	2	100.00%
UA-1		* Denotes Correct Response		
PG: UriScan Reagent Strips	7			
		*0.2/Normal	6	85.71%
		>= 8 mg/dL	1	14.29%
UA-1		* Denotes Correct Response		
PG: Medi-Lab Performance Urine Strip - Visual	8			
		*0.2/Normal	8	100.00%
UA-1		* Denotes Correct Response		
PG: Siemens Clinitek Advantus	23			
		*0.2/Normal	23	100.00%
UA-1		* Denotes Correct Response		
PG: Consult Diag Urinalysis Test Strips	7			
		*0.2/Normal	7	100.00%
UA-1		* Denotes Correct Response		
PG: Consult Diag Urine Analyzer	8			
		*0.2/Normal	8	100.00%
UA-1		* Denotes Correct Response		
PG: 120 Urine Analyzer	42			
		*0.2/Normal	42	100.00%
UA-1		* Denotes Correct Response		
PG: McKesson 10 SG Urine Dipstick	5			
		*0.2/Normal	5	100.00%
UA-1		* Denotes Correct Response		
SG: Select Siemens Clinitek/Visual	962			
		*0.2/Normal	953	99.06%
		1 mg/dL	2	0.21%
		2.0 mg/dL	7	0.73%
UA-1		* Denotes Correct Response		
SG: Select Siemens Clinitek/Visual / Siemens - Plus Sys / Siemens Urine Dipstick	7			
		*3.2 umol/L	7	100.00%
UA-1		* Denotes Correct Response		
SG: Select Roche Urinalysis Inst/Visual	94			
		*0.2/Normal	93	98.94%
		1 mg/dL	1	1.06%
UA-1		* Denotes Correct Response		
SG: Select Consult Diag Urinalysis Inst/Visual	20			
		*0.2/Normal	20	100.00%
UA-1		* Denotes Correct Response		
SG: Select McKesson Urinalysis Inst/Visual	47			
		*0.2/Normal	47	100.00%
UA-1		* Denotes Correct Response		
SG: Select Henry Schein Urinalysis Inst/Visual	7			
		*0.2/Normal	7	100.00%
UA-1		* Denotes Correct Response		
SG: Select Medi-Lab Performance UA Inst/Visual	9			
		*0.2/Normal	9	100.00%
UA-1		* Denotes Correct Response		
MG: Method Group	1190			
		*0.2/Normal	1177	98.91%
		>= 8 mg/dL	2	0.17%
		1 mg/dL	3	0.25%
		2.0 mg/dL	8	0.67%
UA-1		* Denotes Correct Response		
MG: Method Group	2			
		*Normal	2	100.00%
UA-1		* Denotes Correct Response		
MG: Method Group	13			
		*0.2 EU/dL	6	46.15%
		*3.2 umol/L	7	53.85%
UA-1		* Denotes Correct Response		
MG: Method Group	0			
		Urine Microalbumin		
UA-1		* Denotes Correct Response		
PG: Siemens Clinitek 50	8			
		*150 mg/L	8	100.00%
UA-1		* Denotes Correct Response		



PG: Siemens Clinitek Status	50			
		*150 mg/L	48	96.00%
		100 mg/L (3+)	1	2.00%
		Negative	1	2.00%
UA-1		* Denotes Correct Response		
PG: Roche Chemstrip Micral	27			
		*100 mg/L (3+)	16	59.26%
		*150 mg/L	3	11.11%
		*50 mg/L (2+)	7	25.93%
		20/30 mg/L (Small	1	3.70%
UA-1		* Denotes Correct Response		
SG: Semi-quant Methods - Ur Microalbumin-Creat	58			
		*150 mg/L	56	96.55%
		100 mg/L (3+)	1	1.72%
		Negative	1	1.72%
UA-1		* Denotes Correct Response		
SG: Select Siemens Clinitek/Visual	63			
		*100 mg/L (3+)	2	3.17%
		*150 mg/L	58	92.06%
		*50 mg/L (2+)	1	1.59%
		100 mg/L (3+)	2	3.17%
		20/30 mg/L (Small	1	1.59%
		Negative	1	1.59%
UA-1		* Denotes Correct Response		
MG: Method Group	92			
		*100 mg/L (3+)	19	20.65%
		*150 mg/L	62	67.39%
		*50 mg/L (2+)	8	8.70%
		100 mg/L (3+)	19	20.65%
		20/30 mg/L (Small	2	2.17%
		Negative	1	1.09%
		Fluid pH		
UA-1		* Denotes Correct Response		
PG: Nitrazine Paper	12			
		*7	2	16.67%
		*7.5	10	83.33%
		Urine hCG		
UA-1		* Denotes Correct Response		
AG: All Methods	1062			
		*Positive	1053	99.15%
		Negative	9	0.85%
ASI ProPhase Plus	1	*Positive	1	100.00%
Siemens Clinitek Status	13	*Positive	13	100.00%
Beck Coulter ICON 20 HCG	7	*Positive	7	100.00%
Beck Coulter ICON 25 hCG	65	*Positive	64	98.46%
Alere Acceava HCG Basic II	24	*Positive	23	95.83%
Alere Acceava hCG Combo	5	*Positive	5	100.00%
Alere Acceava hCG-Urine	21	*Positive	20	95.24%
Conception Tech One Step hCG	1	*Positive	1	100.00%
Diag Test Group Clarity hCG Plus Urine Cassettes	6	*Positive	6	100.00%
Fisher HC Sure-Vue hCG STAT	4	*Positive	4	100.00%
Fisher HC Sure-Vue Urine hCG	16	*Positive	16	100.00%
Fisher HC Sure-Vure U/S hCG	5	*Positive	5	100.00%
Sekisui Diagnostics OSOM Urine	12	*Positive	12	100.00%
Sekisui Diagnostics OSOM Card Pregnancy	31	*Positive	31	100.00%
Sekisui Diagnostics OSOM hCG Combo	39	*Positive	39	100.00%
Henry Schein One Step + (hCG)	54	*Positive	54	100.00%
Immunostics Detector Combo	2	*Positive	2	100.00%
Immunostics hCG Detector-urine	3	*Positive	3	100.00%
Instant Tech. i Pregnancy	3	*Positive	3	100.00%
Inverness Clearview Combo II	2	*Positive	2	100.00%
Inv Clearview Urine - Cassette	4	*Positive	4	100.00%
Inverness Signify hCG combo	2	*Positive	2	100.00%
Inverness Signify hCG urine	4	*Positive	4	100.00%
LifeSign Status hCG	1	*Positive	1	100.00%
McKesson Medi-Lab hCG - Urine	32	*Positive	31	96.88%
Meridian Immunocard Stat!	1	*Positive	1	100.00%
Permaxim RediScreen	3	*Positive	3	100.00%



Quidel QuickVue + hCG	22	*Positive	22	100.00%
Quidel QuickVue One-Step hCG Combo U/S	58	*Positive	58	100.00%
Quidel QuickVue One-Step Urine	138	*Positive	137	99.28%
Quidel RapidVue hCG - Urine	8	*Positive	8	100.00%
SAS Ultra HCG	2	*Positive	2	100.00%
Select Diag hCG Cassette (urine)	19	*Positive	19	100.00%
Select Diag hCG Combo	10	*Positive	10	100.00%
Select Diag hCG Dipstick	6	*Positive	5	83.33%
SP (Cardinal) Brand combo	26	*Positive	26	100.00%
SP (Cardinal) Brand-cassette	28	*Positive	28	100.00%
SP Brand (Cardinal) HCG Combo	58	*Positive	58	100.00%
Stanbio Qupid	5	*Positive	5	100.00%
Stanbio QuStick	1	*Positive	1	100.00%
Stanbio TRUE hCG	5	*Positive	5	100.00%
MooreBrand hCG Cassette	1	*Positive	1	100.00%
Nova Scientific Century S/U hCG	2	*Positive	2	100.00%
Diag Test Group Clarity hCG Urine Strips	1	*Positive	1	100.00%
Jant AccuTest	4	*Positive	4	100.00%
Instant-View Urine Pregnancy	1	*Positive	1	100.00%
Becton Dickinson Precise One-Step	1	*Positive	1	100.00%
Consult Diag hCG Urine Cassette	137	*Positive	135	98.54%
CLIAwaived Inc Pregnancy Urine Dip-Strip	4	*Positive	4	100.00%
Germaine AIM Step Combo Pregnancy	1	*Positive	1	100.00%
NDC ProAdvantage Urine hCG Cassette	5	*Positive	4	80.00%
Consult Diag hCG Urine Dipstick	47	*Positive	47	100.00%
SP (Cardinal) Brand Dipstick	11	*Positive	11	100.00%
Consult Diag hCG Combo Cassette (Urine)	10	*Positive	10	100.00%
McKesson Medi-Lab Performance hCG Urine Cassette	7	*Positive	7	100.00%
Alere hCG Combo Cassette (25 mIU/mL)	2	*Positive	2	100.00%
Alere hCG Cassette (25 mIU/mL)	12	*Positive	12	100.00%
Alere hCG Dipstick	7	*Positive	7	100.00%
Medline Pregnancy Urine Dip-Strip	2	*Positive	2	100.00%
Medline Pregnancy Urine Cassette Test	6	*Positive	6	100.00%
McKesson hCG Urine Test Dipstick	7	*Positive	7	100.00%
Quidel QuickVue hCG Urine Test	18	*Positive	18	100.00%
MooreBrand hCG Dipstick	1	*Positive	1	100.00%
Siemens Clinitest hCG	1	*Positive	1	100.00%
Henry Schein One Step + HCG Combo	2	*Positive	2	100.00%
Henry Schein One Step hCG Urine Cassette	22	*Positive	22	100.00%
Medline hCG Test Strip 25 mIU/mL	1	*Positive	1	100.00%
NCS hCG One Step Pregnancy Test	1	*Positive	1	100.00%
MediChoice Urine hCG Cassette	1	*Positive	1	100.00%
Beck Coulter ICON 25 hCG	65	Negative	1	1.54%
Alere Acceava HCG Basic II	24	Negative	1	4.17%
Alere Acceava hCG-Urine	21	Negative	1	4.76%
McKesson Medi-Lab hCG - Urine	32	Negative	1	3.13%
Quidel QuickVue One-Step Urine	138	Negative	1	0.72%
Select Diag hCG Dipstick	6	Negative	1	16.67%
Consult Diag hCG Urine Cassette	137	Negative	2	1.46%
NDC ProAdvantage Urine hCG Cassette	5	Negative	1	20.00%
		Urine hCG (non-waived)		
UA-1		* Denotes Correct Response		
PG: Quidel Sofia/Sofia 2 hCG (urine)	1			
		*Positive	1	100.00%
661 - Urine Colony Count;784 - Presump ID/Colony Count - Urine;785 - Presump ID/Colony Count-Urine & Strep Combo;775 - Urine Culture;778 - Urine & Throat Cultures;779 - Urine & Genital Cultures;781 - Urine, Throat & Genital Cultures;782 - Complete Bacteri	Labs	Response	Number	Percent
		Susceptibility Cefepime		
		* Denotes Correct Response		
UC-1				
AG: All Methods	1			
		*Sensitive	1	100.00%
Susceptibility	1	*Sensitive	1	100.00%
		Susceptibility Levofloxacin		
		* Denotes Correct Response		
UC-1				
AG: All Methods	3			
		*Intermediate	2	66.67%
		*Resistant	1	33.33%



Susceptibility	3	*Intermediate	2	66.67%
Susceptibility		*Resistant	1	33.33%
		Susceptibility Piperacillin/Tazobactam		
UC-1		* Denotes Correct Response		
AG: All Methods	1			
		*Sensitive	1	100.00%
Susceptibility	1	*Sensitive	1	100.00%
		Susceptibility Tobramycin		
UC-1		* Denotes Correct Response		
AG: All Methods	2			
		*Sensitive	2	100.00%
Susceptibility	2	*Sensitive	2	100.00%
		Urine Culture		
UC-1		* Denotes Correct Response		
AG: All Methods	31			
		*Gram Negative organism	3	9.68%
		*Growth, referred for identification	13	41.94%
		*Presumptive Pseudomonas species	2	6.45%
		*Pseudomonas aeruginosa	2	6.45%
		*Pseudomonas species, NOS	2	6.45%
		*Uricult - Growth (Positive)	8	25.81%
		Candida species, NOS	1	3.23%
Cultures	31	*Gram Negative organism	3	9.68%
Cultures		*Growth, referred for identification	13	41.94%
Cultures		*Presumptive Pseudomonas species	2	6.45%
Cultures		*Pseudomonas aeruginosa	2	6.45%
Cultures		*Pseudomonas species, NOS	2	6.45%
Cultures		*Uricult - Growth (Positive)	8	25.81%
Cultures		Candida species, NOS	1	3.23%
UC-2		* Denotes Correct Response		
AG: All Methods	31			
		*Candida species, NOS	1	3.23%
		*Gram Positive organism	2	6.45%
		*Growth, referred for identification	15	48.39%
		*Uricult - Growth (Positive)	7	22.58%
		*Yeast, NOS	5	16.13%
		Uricult - No Growth (Negative)	1	3.23%
Cultures	31	*Candida species, NOS	1	3.23%
Cultures		*Gram Positive organism	2	6.45%
Cultures		*Growth, referred for identification	15	48.39%
Cultures		*Uricult - Growth (Positive)	7	22.58%
Cultures		*Yeast, NOS	5	16.13%
Cultures		Uricult - No Growth (Negative)	1	3.23%
UC-3		* Denotes Correct Response		
AG: All Methods	30			
		*Escherichia coli	2	6.67%
		*Gram Negative organism	4	13.33%
		*Growth, referred for identification	13	43.33%
		*Presumptive E. coli	2	6.67%
		*Uricult - Growth (Positive)	8	26.67%
		Yeast, NOS	1	3.33%
Cultures	30	*Escherichia coli	2	6.67%
Cultures		*Gram Negative organism	4	13.33%
Cultures		*Growth, referred for identification	13	43.33%
Cultures		*Presumptive E. coli	2	6.67%
Cultures		*Uricult - Growth (Positive)	8	26.67%
Cultures		Yeast, NOS	1	3.33%
UC-4		* Denotes Correct Response		
AG: All Methods	14			
		*Gram Negative organism	1	7.14%
		*Growth, referred for identification	4	28.57%
		*Serratia marcescens	2	14.29%
		*Uricult - Growth (Positive)	5	35.71%
		Candida species, NOS	1	7.14%
		Presumptive Klebsiella species	1	7.14%
Cultures	14	*Gram Negative organism	1	7.14%
Cultures		*Growth, referred for identification	4	28.57%
Cultures		*Serratia marcescens	2	14.29%



Cultures		*Uricult - Growth (Positive)	5	35.71%
Cultures		Candida species, NOS	1	7.14%
Cultures		Presumptive Klebsiella species	1	7.14%
UC-5		* Denotes Correct Response		
AG: All Methods	13			
		*Gram Positive organism	2	15.38%
		*Growth, referred for identification	3	23.08%
		*Presumptive Staphylococcus species	1	7.69%
		*Staphylococcus aureus	2	15.38%
		*Uricult - Growth (Positive)	5	38.46%
Cultures	13	*Gram Positive organism	2	15.38%
Cultures		*Growth, referred for identification	3	23.08%
Cultures		*Presumptive Staphylococcus species	1	7.69%
Cultures		*Staphylococcus aureus	2	15.38%
Cultures		*Uricult - Growth (Positive)	5	38.46%
703 - Urine Drug Screen;713 - Urine Drug Screen (50 ml sample volume)	Labs	Response	Number	Percent
		Barbiturates		
UDS-1		* Denotes Correct Response		
AG: All Methods	137			
		*Negative	136	99.27%
		Positive	1	0.73%
AlereTriage	3	*Negative	3	100.00%
E-Screen/Alere	26	*Negative	26	100.00%
Express Diag. Express Check	1	*Negative	1	100.00%
Exp Results Integ Multi-Drug Screen Cup	1	*Negative	1	100.00%
Alere iCup Drugs of Abuse	13	*Negative	13	100.00%
Instant Technologies i Screen	1	*Negative	1	100.00%
MedTox Profile II	2	*Negative	2	100.00%
Beckman AU	1	*Negative	1	100.00%
Innovacon Multi-Clin	9	*Negative	9	100.00%
MedTox Profile III ER	1	*Negative	1	100.00%
CLIAwaived Inc Rapid Drug Cup	4	*Negative	4	100.00%
MedTox Profile V	2	*Negative	2	100.00%
Branan Med Corp QuickTox	1	*Negative	1	100.00%
Jant Accutest Multi Drug Screen	2	*Negative	2	100.00%
Amedica Biotech AmediCheck Instant Test Cup	4	*Negative	4	100.00%
Nova Century Scientific Nova Screen DOA Cassette	1	*Negative	1	100.00%
American Screening Corp Reveal Multi-Drug Test Cup	4	*Negative	4	100.00%
Mindray BS-200	2	*Negative	2	100.00%
Medica Easy RA	2	*Negative	2	100.00%
Alere iCassette	2	*Negative	2	100.00%
Roche Hitachi 717	1	*Negative	1	100.00%
Alere USD UScreen Cup	9	*Negative	9	100.00%
Alere iScreen Drugs of Abuse Card	1	*Negative	1	100.00%
Noble Medical Insta-Screen Test Cup	3	*Negative	3	100.00%
Hemosure First Sign UDS Dipstick	15	*Negative	14	93.33%
Siemens Viva/ProE System	2	*Negative	2	100.00%
Premier Biotech Bio-Cup	12	*Negative	12	100.00%
UCP Biosciences Home Drug Screening Test Cup	1	*Negative	1	100.00%
Abbott Architect c and ci Series	2	*Negative	2	100.00%
Consult Diag Drugs of Abuse Cups	1	*Negative	1	100.00%
Branan Med Corp ToxCup	6	*Negative	6	100.00%
American Screening Discover Plus Cup	1	*Negative	1	100.00%
McKesson Drugs of Abuse Test Cup	1	*Negative	1	100.00%
Hemosure First Sign UDS Dipstick	15	Positive	1	6.67%
UDS-2		* Denotes Correct Response		
AG: All Methods	138			
		*Negative	137	99.28%
		Positive	1	0.72%
AlereTriage	3	*Negative	3	100.00%
E-Screen/Alere	26	*Negative	26	100.00%
Express Diag. Express Check	1	*Negative	1	100.00%
Exp Results Integ Multi-Drug Screen Cup	1	*Negative	1	100.00%
Alere iCup Drugs of Abuse	13	*Negative	13	100.00%
Instant Technologies i Screen	1	*Negative	1	100.00%
MedTox Profile II	2	*Negative	2	100.00%
Beckman AU	1	*Negative	1	100.00%
Innovacon Multi-Clin	9	*Negative	9	100.00%



MedTox Profile III ER	1	*Negative	1	100.00%
CLIAwaived Inc Rapid Drug Cup	5	*Negative	5	100.00%
MedTox Profile V	2	*Negative	2	100.00%
Branan Med Corp QuickTox	1	*Negative	1	100.00%
Jant Accutest Multi Drug Screen	2	*Negative	2	100.00%
Amedica Biotech AmediCheck Instant Test Cup	4	*Negative	4	100.00%
Nova Century Scientific Nova Screen DOA Cassette	1	*Negative	1	100.00%
American Screening Corp Reveal Multi-Drug Test Cup	4	*Negative	4	100.00%
Mindray BS-200	2	*Negative	2	100.00%
Medica Easy RA	2	*Negative	2	100.00%
Alere iCassette	2	*Negative	2	100.00%
Roche Hitachi 717	1	*Negative	1	100.00%
Alere USD UScreen Cup	9	*Negative	9	100.00%
Alere iScreen Drugs of Abuse Card	1	*Negative	1	100.00%
Noble Medical Insta-Screen Test Cup	3	*Negative	3	100.00%
Hemosure First Sign UDS Dipstick	15	*Negative	14	93.33%
Siemens Viva/ProE System	2	*Negative	2	100.00%
Premier Biotech Bio-Cup	12	*Negative	12	100.00%
UCP Biosciences Home Drug Screening Test Cup	1	*Negative	1	100.00%
Abbott Architect c and ci Series	2	*Negative	2	100.00%
Consult Diag Drugs of Abuse Cups	1	*Negative	1	100.00%
Branan Med Corp ToxCup	6	*Negative	6	100.00%
American Screening Discover Plus Cup	1	*Negative	1	100.00%
McKesson Drugs of Abuse Test Cup	1	*Negative	1	100.00%
Hemosure First Sign UDS Dipstick	15	Positive	1	6.67%
		Benzodiazepines		
UDS-1		* Denotes Correct Response		
AG: All Methods	148			
		*Positive	141	95.27%
		Negative	7	4.73%
AlereTriage	3	*Positive	3	100.00%
Siemens Dim/AR/ES/RxL/Xpand/EXL	1	*Positive	1	100.00%
E-Screen/Alere	26	*Positive	25	96.15%
Express Diag. Express Check	1	*Positive	1	100.00%
Exp Results Integ Multi-Drug Screen Cup	1	*Positive	1	100.00%
Alere iCup Drugs of Abuse	16	*Positive	15	93.75%
Instant Technologies i Screen	1	*Positive	1	100.00%
MedTox Profile II	2	*Positive	2	100.00%
Noble Med Split Specimen Cup	1	*Positive	1	100.00%
Beckman AU	1	*Positive	1	100.00%
Innovacon Multi-Clin	9	*Positive	9	100.00%
MedTox Profile III ER	1	*Positive	1	100.00%
CLIAwaived Inc Rapid Drug Cup	6	*Positive	6	100.00%
MedTox Profile V	2	*Positive	2	100.00%
Branan Med Corp QuickTox	1	*Positive	1	100.00%
Jant Accutest Multi Drug Screen	2	*Positive	2	100.00%
Amedica Biotech AmediCheck Instant Test Cup	4	*Positive	4	100.00%
Nova Century Scientific Nova Screen DOA Cassette	1	*Positive	1	100.00%
American Screening Corp Reveal Multi-Drug Test Cup	4	*Positive	4	100.00%
Mindray BS-200	2	*Positive	2	100.00%
Medica Easy RA	2	*Positive	2	100.00%
Alere iCassette	2	*Positive	1	50.00%
Roche Hitachi 717	1	*Positive	1	100.00%
Alere USD UScreen Cup	9	*Positive	8	88.89%
ThermoScientific Indiko Plus	1	*Positive	1	100.00%
Alere iScreen Drugs of Abuse Card	1	*Positive	1	100.00%
Noble Medical Insta-Screen Test Cup	3	*Positive	2	66.67%
Hemosure First Sign UDS Dipstick	15	*Positive	14	93.33%
Siemens Viva/ProE System	4	*Positive	4	100.00%
TFS DRI Assay	1	*Positive	1	100.00%
Premier Biotech Bio-Cup	12	*Positive	12	100.00%
UCP Biosciences Home Drug Screening Test Cup	1	*Positive	1	100.00%
Abbott Architect c and ci Series	2	*Positive	2	100.00%
Consult Diag Drugs of Abuse Cups	1	*Positive	1	100.00%
Branan Med Corp ToxCup	6	*Positive	5	83.33%
American Screening Discover Plus Cup	1	*Positive	1	100.00%
McKesson Drugs of Abuse Test Cup	1	*Positive	1	100.00%
E-Screen/Alere	26	Negative	1	3.85%
Alere iCup Drugs of Abuse	16	Negative	1	6.25%



Alere iCassette	2	Negative	1	50.00%
Alere USD UScreen Cup	9	Negative	1	11.11%
Noble Medical Insta-Screen Test Cup	3	Negative	1	33.33%
Hemosure First Sign UDS Dipstick	15	Negative	1	6.67%
Branan Med Corp ToxCup	6	Negative	1	16.67%
UDS-2		* Denotes Correct Response		
AG: All Methods	148			
		*Negative	145	97.97%
		Positive	3	2.03%
AlereTriage	3	*Negative	3	100.00%
Siemens Dim/AR/ES/RxL/Xpand/EXL	1	*Negative	1	100.00%
E-Screen/Alere	26	*Negative	26	100.00%
Express Diag. Express Check	1	*Negative	1	100.00%
Exp Results Integ Multi-Drug Screen Cup	1	*Negative	1	100.00%
Alere iCup Drugs of Abuse	16	*Negative	16	100.00%
Instant Technologies i Screen	1	*Negative	1	100.00%
MedTox Profile II	2	*Negative	2	100.00%
Noble Med Split Specimen Cup	1	*Negative	1	100.00%
Beckman AU	1	*Negative	1	100.00%
Innovacon Multi-Clin	9	*Negative	9	100.00%
MedTox Profile III ER	1	*Negative	1	100.00%
CLIAwaived Inc Rapid Drug Cup	6	*Negative	6	100.00%
MedTox Profile V	2	*Negative	2	100.00%
Branan Med Corp QuickTox	1	*Negative	1	100.00%
Jant Accutest Multi Drug Screen	2	*Negative	1	50.00%
Amedica Biotech AmediCheck Instant Test Cup	4	*Negative	4	100.00%
Nova Century Scientific Nova Screen DOA Cassette	1	*Negative	1	100.00%
American Screening Corp Reveal Multi-Drug Test Cup	4	*Negative	4	100.00%
Mindray BS-200	2	*Negative	2	100.00%
Medica Easy RA	2	*Negative	2	100.00%
Alere iCassette	2	*Negative	2	100.00%
Roche Hitachi 717	1	*Negative	1	100.00%
Alere USD UScreen Cup	9	*Negative	8	88.89%
ThermoScientific Indiko Plus	1	*Negative	1	100.00%
Alere iScreen Drugs of Abuse Card	1	*Negative	1	100.00%
Noble Medical Insta-Screen Test Cup	3	*Negative	3	100.00%
Hemosure First Sign UDS Dipstick	15	*Negative	14	93.33%
Siemens Viva/ProE System	4	*Negative	4	100.00%
TFS DRI Assay	1	*Negative	1	100.00%
Premier Biotech Bio-Cup	12	*Negative	12	100.00%
UCP Biosciences Home Drug Screening Test Cup	1	*Negative	1	100.00%
Abbott Architect c and ci Series	2	*Negative	2	100.00%
Consult Diag Drugs of Abuse Cups	1	*Negative	1	100.00%
Branan Med Corp ToxCup	6	*Negative	6	100.00%
American Screening Discover Plus Cup	1	*Negative	1	100.00%
McKesson Drugs of Abuse Test Cup	1	*Negative	1	100.00%
Jant Accutest Multi Drug Screen	2	Positive	1	50.00%
Alere USD UScreen Cup	9	Positive	1	11.11%
Hemosure First Sign UDS Dipstick	15	Positive	1	6.67%
		Amphetamines		
UDS-1		* Denotes Correct Response		
AG: All Methods	152			
		*Negative	148	97.37%
		Positive	4	2.63%
AlereTriage	3	*Negative	3	100.00%
Siemens Dim/AR/ES/RxL/Xpand/EXL	2	*Negative	2	100.00%
E-Screen/Alere	23	*Negative	23	100.00%
Express Diag. Express Check	1	*Negative	1	100.00%
Exp Results Integ Multi-Drug Screen Cup	1	*Negative	1	100.00%
Alere iCup Drugs of Abuse	19	*Negative	17	89.47%
Instant Technologies i Screen	1	*Negative	1	100.00%
MedTox Profile II	3	*Negative	3	100.00%
Noble Med Split Specimen Cup	1	*Negative	1	100.00%
Beckman AU	1	*Negative	1	100.00%
Innovacon Multi-Clin	8	*Negative	8	100.00%
MedTox Profile III ER	1	*Negative	1	100.00%
CLIAwaived Inc Rapid Drug Cup	6	*Negative	6	100.00%
MedTox Profile V	2	*Negative	2	100.00%
Branan Med Corp QuickTox	1	*Negative	1	100.00%



Jant Accutest Multi Drug Screen	2	*Negative	2	100.00%
Amedica Biotech AmediCheck Instant Test Cup	4	*Negative	4	100.00%
Nova Century Scientific Nova Screen DOA Cassette	1	*Negative	1	100.00%
American Screening Corp Reveal Multi-Drug Test Cup	4	*Negative	4	100.00%
Mindray BS-200	2	*Negative	2	100.00%
Medica Easy RA	2	*Negative	2	100.00%
Alere iCassette	4	*Negative	3	75.00%
Roche Hitachi 717	1	*Negative	1	100.00%
Alere USD UScreen Cup	9	*Negative	9	100.00%
ThermoScientific Indiko Plus	1	*Negative	1	100.00%
Alere iScreen Drugs of Abuse Card	1	*Negative	1	100.00%
Noble Medical Insta-Screen Test Cup	3	*Negative	3	100.00%
Hemosure First Sign UDS Dipstick	15	*Negative	14	93.33%
Siemens Viva/ProE System	5	*Negative	5	100.00%
TFS DRI Assay	1	*Negative	1	100.00%
Premier Biotech Bio-Cup	12	*Negative	12	100.00%
UCP Biosciences Home Drug Screening Test Cup	1	*Negative	1	100.00%
Abbott Architect c and ci Series	2	*Negative	2	100.00%
Consult Diag Drugs of Abuse Cups	1	*Negative	1	100.00%
Branan Med Corp ToxCup	6	*Negative	6	100.00%
American Screening Discover Plus Cup	1	*Negative	1	100.00%
McKesson Drugs of Abuse Test Cup	1	*Negative	1	100.00%
Alere iCup Drugs of Abuse	19	Positive	2	10.53%
Alere iCassette	4	Positive	1	25.00%
Hemosure First Sign UDS Dipstick	15	Positive	1	6.67%
UDS-2		* Denotes Correct Response		
AG: All Methods	152			
		*Negative	151	99.34%
		Positive	1	0.66%
AlereTriage	3	*Negative	3	100.00%
Siemens Dim/AR/ES/RxL/Xpand/EXL	2	*Negative	2	100.00%
E-Screen/Alere	23	*Negative	23	100.00%
Express Diag. Express Check	1	*Negative	1	100.00%
Exp Results Integ Multi-Drug Screen Cup	1	*Negative	1	100.00%
Alere iCup Drugs of Abuse	19	*Negative	19	100.00%
Instant Technologies i Screen	1	*Negative	1	100.00%
MedTox Profile II	3	*Negative	3	100.00%
Noble Med Split Specimen Cup	1	*Negative	1	100.00%
Beckman AU	1	*Negative	1	100.00%
Innovacon Multi-Clin	8	*Negative	8	100.00%
MedTox Profile III ER	1	*Negative	1	100.00%
CLIAwaived Inc Rapid Drug Cup	6	*Negative	6	100.00%
MedTox Profile V	2	*Negative	2	100.00%
Branan Med Corp QuickTox	1	*Negative	1	100.00%
Jant Accutest Multi Drug Screen	2	*Negative	2	100.00%
Amedica Biotech AmediCheck Instant Test Cup	4	*Negative	4	100.00%
Nova Century Scientific Nova Screen DOA Cassette	1	*Negative	1	100.00%
American Screening Corp Reveal Multi-Drug Test Cup	4	*Negative	4	100.00%
Mindray BS-200	2	*Negative	2	100.00%
Medica Easy RA	2	*Negative	2	100.00%
Alere iCassette	4	*Negative	4	100.00%
Roche Hitachi 717	1	*Negative	1	100.00%
Alere USD UScreen Cup	9	*Negative	9	100.00%
ThermoScientific Indiko Plus	1	*Negative	1	100.00%
Alere iScreen Drugs of Abuse Card	1	*Negative	1	100.00%
Noble Medical Insta-Screen Test Cup	3	*Negative	3	100.00%
Hemosure First Sign UDS Dipstick	15	*Negative	14	93.33%
Siemens Viva/ProE System	5	*Negative	5	100.00%
TFS DRI Assay	1	*Negative	1	100.00%
Premier Biotech Bio-Cup	12	*Negative	12	100.00%
UCP Biosciences Home Drug Screening Test Cup	1	*Negative	1	100.00%
Abbott Architect c and ci Series	2	*Negative	2	100.00%
Consult Diag Drugs of Abuse Cups	1	*Negative	1	100.00%
Branan Med Corp ToxCup	6	*Negative	6	100.00%
American Screening Discover Plus Cup	1	*Negative	1	100.00%
McKesson Drugs of Abuse Test Cup	1	*Negative	1	100.00%
Hemosure First Sign UDS Dipstick	15	Positive	1	6.67%
		Cannabinoids (THC)		
UDS-1		* Denotes Correct Response		



AG: All Methods	183			
		*Negative	180	98.36%
		Positive	3	1.64%
AlereTriage	3	*Negative	3	100.00%
Siemens Dim/AR/ES/RxL/Xpand/EXL	2	*Negative	2	100.00%
E-Screen/Alere	51	*Negative	51	100.00%
Express Diag. Express Check	1	*Negative	1	100.00%
Exp Results Integ Multi-Drug Screen Cup	1	*Negative	1	100.00%
Alere iCup Drugs of Abuse	19	*Negative	18	94.74%
Instant Technologies i Screen	2	*Negative	2	100.00%
Integrated E-Z Split Key Cup	1	*Negative	1	100.00%
MedTox Profile II	3	*Negative	3	100.00%
Noble Med Split Specimen Cup	1	*Negative	1	100.00%
Beckman AU	1	*Negative	1	100.00%
Innovacon Multi-Clin	9	*Negative	9	100.00%
MedTox Profile III ER	1	*Negative	1	100.00%
CLIAwaived Inc Rapid Drug Cup	6	*Negative	5	83.33%
MedTox Profile V	2	*Negative	2	100.00%
Branan Med Corp QuickTox	1	*Negative	1	100.00%
Jant Accutest Multi Drug Screen	2	*Negative	2	100.00%
Amedica Biotech AmediCheck Instant Test Cup	4	*Negative	4	100.00%
Nova Century Scientific Nova Screen DOA Cassette	1	*Negative	1	100.00%
American Screening Corp Reveal Multi-Drug Test Cup	4	*Negative	4	100.00%
Mindray BS-200	2	*Negative	2	100.00%
Medica Easy RA	2	*Negative	2	100.00%
Alere iCassette	4	*Negative	4	100.00%
Roche Hitachi 717	1	*Negative	1	100.00%
Alere USD UScreen Cup	9	*Negative	9	100.00%
ThermoScientific Indiko Plus	1	*Negative	1	100.00%
Alere iScreen Drugs of Abuse Card	1	*Negative	1	100.00%
Noble Medical Insta-Screen Test Cup	3	*Negative	3	100.00%
CRL Stat Cup	1	*Negative	1	100.00%
Hemosure First Sign UDS Dipstick	15	*Negative	14	93.33%
Siemens Viva/ProE System	4	*Negative	4	100.00%
TFS DRI Assay	1	*Negative	1	100.00%
Premier Biotech Bio-Cup	12	*Negative	12	100.00%
UCP Biosciences Home Drug Screening Test Cup	1	*Negative	1	100.00%
Abbott Architect c and ci Series	2	*Negative	2	100.00%
Consult Diag Drugs of Abuse Cups	1	*Negative	1	100.00%
Branan Med Corp ToxCup	6	*Negative	6	100.00%
American Screening Discover Plus Cup	1	*Negative	1	100.00%
McKesson Drugs of Abuse Test Cup	1	*Negative	1	100.00%
Alere iCup Drugs of Abuse	19	Positive	1	5.26%
CLIAwaived Inc Rapid Drug Cup	6	Positive	1	16.67%
Hemosure First Sign UDS Dipstick	15	Positive	1	6.67%
UDS-2		* Denotes Correct Response		
AG: All Methods	182			
		*Negative	181	99.45%
		Positive	1	0.55%
AlereTriage	3	*Negative	3	100.00%
Siemens Dim/AR/ES/RxL/Xpand/EXL	2	*Negative	2	100.00%
E-Screen/Alere	51	*Negative	51	100.00%
Express Diag. Express Check	1	*Negative	1	100.00%
Exp Results Integ Multi-Drug Screen Cup	1	*Negative	1	100.00%
Alere iCup Drugs of Abuse	19	*Negative	19	100.00%
Instant Technologies i Screen	2	*Negative	2	100.00%
Integrated E-Z Split Key Cup	1	*Negative	1	100.00%
MedTox Profile II	3	*Negative	3	100.00%
Noble Med Split Specimen Cup	1	*Negative	1	100.00%
Beckman AU	1	*Negative	1	100.00%
Innovacon Multi-Clin	9	*Negative	9	100.00%
MedTox Profile III ER	1	*Negative	1	100.00%
CLIAwaived Inc Rapid Drug Cup	6	*Negative	6	100.00%
MedTox Profile V	2	*Negative	2	100.00%
Branan Med Corp QuickTox	1	*Negative	1	100.00%
Jant Accutest Multi Drug Screen	2	*Negative	2	100.00%
Amedica Biotech AmediCheck Instant Test Cup	4	*Negative	4	100.00%
Nova Century Scientific Nova Screen DOA Cassette	1	*Negative	1	100.00%
American Screening Corp Reveal Multi-Drug Test Cup	4	*Negative	4	100.00%



Mindray BS-200	2	*Negative	2	100.00%
Medica Easy RA	2	*Negative	2	100.00%
Alere iCassette	4	*Negative	4	100.00%
Roche Hitachi 717	1	*Negative	1	100.00%
Alere USD UScreen Cup	9	*Negative	9	100.00%
ThermoScientific Indiko Plus	1	*Negative	1	100.00%
Alere iScreen Drugs of Abuse Card	1	*Negative	1	100.00%
Noble Medical Insta-Screen Test Cup	3	*Negative	3	100.00%
CRL Stat Cup	1	*Negative	1	100.00%
Hemosure First Sign UDS Dipstick	15	*Negative	14	93.33%
Siemens Viva/ProE System	3	*Negative	3	100.00%
TFS DRI Assay	1	*Negative	1	100.00%
Premier Biotech Bio-Cup	12	*Negative	12	100.00%
UCP Biosciences Home Drug Screening Test Cup	1	*Negative	1	100.00%
Abbott Architect c and ci Series	2	*Negative	2	100.00%
Consult Diag Drugs of Abuse Cups	1	*Negative	1	100.00%
Branan Med Corp ToxCup	6	*Negative	6	100.00%
American Screening Discover Plus Cup	1	*Negative	1	100.00%
McKesson Drugs of Abuse Test Cup	1	*Negative	1	100.00%
Hemosure First Sign UDS Dipstick	15	Positive	1	6.67%
		Cocaine Metabolites		
UDS-1		* Denotes Correct Response		
AG: All Methods	183			
		*Negative	181	98.91%
		Positive	2	1.09%
AlereTriage	3	*Negative	3	100.00%
Siemens Dim/AR/ES/RxL/Xpand/EXL	2	*Negative	2	100.00%
E-Screen/Alere	51	*Negative	51	100.00%
Express Diag. Express Check	1	*Negative	1	100.00%
Exp Results Integ Multi-Drug Screen Cup	1	*Negative	1	100.00%
Alere iCup Drugs of Abuse	19	*Negative	18	94.74%
Instant Technologies i Screen	2	*Negative	2	100.00%
Integrated E-Z Split Key Cup	1	*Negative	1	100.00%
MedTox Profile II	3	*Negative	3	100.00%
Noble Med Split Specimen Cup	1	*Negative	1	100.00%
Beckman AU	1	*Negative	1	100.00%
Innovacon Multi-Clin	9	*Negative	9	100.00%
MedTox Profile III ER	1	*Negative	1	100.00%
CLIAwaived Inc Rapid Drug Cup	6	*Negative	6	100.00%
MedTox Profile V	2	*Negative	2	100.00%
Branan Med Corp QuickTox	1	*Negative	1	100.00%
Jant Accutest Multi Drug Screen	2	*Negative	2	100.00%
Amedica Biotech AmediCheck Instant Test Cup	4	*Negative	4	100.00%
Nova Century Scientific Nova Screen DOA Cassette	1	*Negative	1	100.00%
American Screening Corp Reveal Multi-Drug Test Cup	4	*Negative	4	100.00%
Mindray BS-200	2	*Negative	2	100.00%
Medica Easy RA	2	*Negative	2	100.00%
Alere iCassette	4	*Negative	4	100.00%
Roche Hitachi 717	1	*Negative	1	100.00%
Alere USD UScreen Cup	9	*Negative	9	100.00%
Alere iScreen Drugs of Abuse Card	1	*Negative	1	100.00%
Noble Medical Insta-Screen Test Cup	3	*Negative	3	100.00%
CRL Stat Cup	1	*Negative	1	100.00%
Hemosure First Sign UDS Dipstick	15	*Negative	14	93.33%
Siemens Viva/ProE System	5	*Negative	5	100.00%
TFS DRI Assay	1	*Negative	1	100.00%
Premier Biotech Bio-Cup	12	*Negative	12	100.00%
UCP Biosciences Home Drug Screening Test Cup	1	*Negative	1	100.00%
Abbott Architect c and ci Series	2	*Negative	2	100.00%
Consult Diag Drugs of Abuse Cups	1	*Negative	1	100.00%
Branan Med Corp ToxCup	6	*Negative	6	100.00%
American Screening Discover Plus Cup	1	*Negative	1	100.00%
McKesson Drugs of Abuse Test Cup	1	*Negative	1	100.00%
Alere iCup Drugs of Abuse	19	Positive	1	5.26%
Hemosure First Sign UDS Dipstick	15	Positive	1	6.67%
UDS-2		* Denotes Correct Response		
AG: All Methods	183			
		*Negative	182	99.45%
		Positive	1	0.55%



AlereTriage	3	*Negative	3	100.00%
Siemens Dim/AR/ES/RxL/Xpand/EXL	2	*Negative	2	100.00%
E-Screen/Alere	51	*Negative	51	100.00%
Express Diag. Express Check	1	*Negative	1	100.00%
Exp Results Integ Multi-Drug Screen Cup	1	*Negative	1	100.00%
Alere iCup Drugs of Abuse	19	*Negative	19	100.00%
Instant Technologies i Screen	2	*Negative	2	100.00%
Integrated E-Z Split Key Cup	1	*Negative	1	100.00%
MedTox Profile II	3	*Negative	3	100.00%
Noble Med Split Specimen Cup	1	*Negative	1	100.00%
Beckman AU	1	*Negative	1	100.00%
Innovacon Multi-Clin	9	*Negative	9	100.00%
MedTox Profile III ER	1	*Negative	1	100.00%
CLIAwaived Inc Rapid Drug Cup	6	*Negative	6	100.00%
MedTox Profile V	2	*Negative	2	100.00%
Branan Med Corp QuickTox	1	*Negative	1	100.00%
Jant Accutest Multi Drug Screen	2	*Negative	2	100.00%
Amedica Biotech AmediCheck Instant Test Cup	4	*Negative	4	100.00%
Nova Century Scientific Nova Screen DOA Cassette	1	*Negative	1	100.00%
American Screening Corp Reveal Multi-Drug Test Cup	4	*Negative	4	100.00%
Mindray BS-200	2	*Negative	2	100.00%
Medica Easy RA	2	*Negative	2	100.00%
Alere iCassette	4	*Negative	4	100.00%
Roche Hitachi 717	1	*Negative	1	100.00%
Alere USD UScreen Cup	9	*Negative	9	100.00%
Alere iScreen Drugs of Abuse Card	1	*Negative	1	100.00%
Noble Medical Insta-Screen Test Cup	3	*Negative	3	100.00%
CRL Stat Cup	1	*Negative	1	100.00%
Hemosure First Sign UDS Dipstick	15	*Negative	14	93.33%
Siemens Viva/ProE System	5	*Negative	5	100.00%
TFS DRI Assay	1	*Negative	1	100.00%
Premier Biotech Bio-Cup	12	*Negative	12	100.00%
UCP Biosciences Home Drug Screening Test Cup	1	*Negative	1	100.00%
Abbott Architect c and ci Series	2	*Negative	2	100.00%
Consult Diag Drugs of Abuse Cups	1	*Negative	1	100.00%
Branan Med Corp ToxCup	6	*Negative	6	100.00%
American Screening Discover Plus Cup	1	*Negative	1	100.00%
McKesson Drugs of Abuse Test Cup	1	*Negative	1	100.00%
Hemosure First Sign UDS Dipstick	15	Positive	1	6.67%
		Ethanol		
UDS-1		* Denotes Correct Response		
AG: All Methods	8			
		*Negative	8	100.00%
Mindray BS-200	2	*Negative	2	100.00%
Medica Easy RA	1	*Negative	1	100.00%
Roche Hitachi 717	1	*Negative	1	100.00%
Siemens Viva/ProE System	2	*Negative	2	100.00%
TFS DRI Assay	1	*Negative	1	100.00%
Abbott Architect c and ci Series	1	*Negative	1	100.00%
UDS-2		* Denotes Correct Response		
AG: All Methods	8			
		*Positive	5	62.50%
		Negative	3	37.50%
Mindray BS-200	2	*Positive	1	50.00%
Medica Easy RA	1	*Positive	1	100.00%
Siemens Viva/ProE System	2	*Positive	1	50.00%
TFS DRI Assay	1	*Positive	1	100.00%
Abbott Architect c and ci Series	1	*Positive	1	100.00%
Mindray BS-200	2	Negative	1	50.00%
Roche Hitachi 717	1	Negative	1	100.00%
Siemens Viva/ProE System	2	Negative	1	50.00%
		Acetaminophen		
UDS-1		* Denotes Correct Response		
AG: All Methods	5			
		*Negative	5	100.00%
AlereTriage	3	*Negative	3	100.00%
E-Screen/Alere	1	*Negative	1	100.00%
Abbott Architect c and ci Series	1	*Negative	1	100.00%
UDS-2		* Denotes Correct Response		



AG: All Methods	5			
		*Positive	4	80.00%
		Negative	1	20.00%
AlereTriage	3	*Positive	3	100.00%
Abbott Architect c and ci Series	1	*Positive	1	100.00%
E-Screen/Alere	1	Negative	1	100.00%
		MDMA (Ecstasy)		
UDS-1		* Denotes Correct Response		
AG: All Methods	95			
		*Negative	93	97.89%
		Positive	2	2.11%
E-Screen/Alere	4	*Negative	4	100.00%
Alere iCup Drugs of Abuse	13	*Negative	12	92.31%
Innovacon Multi-Clin	9	*Negative	9	100.00%
CLIAwaived Inc Rapid Drug Cup	5	*Negative	5	100.00%
Branan Med Corp QuickTox	1	*Negative	1	100.00%
Jant Accutest Multi Drug Screen	2	*Negative	2	100.00%
Amedica Biotech AmediCheck Instant Test Cup	4	*Negative	4	100.00%
Nova Century Scientific Nova Screen DOA Cassette	1	*Negative	1	100.00%
American Screening Corp Reveal Multi-Drug Test Cup	2	*Negative	2	100.00%
Medica Easy RA	1	*Negative	1	100.00%
Alere iCassette	2	*Negative	2	100.00%
Alere USD UScreen Cup	9	*Negative	9	100.00%
Alere iScreen Drugs of Abuse Card	1	*Negative	1	100.00%
Noble Medical Insta-Screen Test Cup	3	*Negative	3	100.00%
Hemosure First Sign UDS Dipstick	15	*Negative	14	93.33%
Siemens Viva/ProE System	1	*Negative	1	100.00%
Premier Biotech Bio-Cup	12	*Negative	12	100.00%
UCP Biosciences Home Drug Screening Test Cup	1	*Negative	1	100.00%
Consult Diag Drugs of Abuse Cups	1	*Negative	1	100.00%
Branan Med Corp ToxCup	6	*Negative	6	100.00%
American Screening Discover Plus Cup	1	*Negative	1	100.00%
McKesson Drugs of Abuse Test Cup	1	*Negative	1	100.00%
Alere iCup Drugs of Abuse	13	Positive	1	7.69%
Hemosure First Sign UDS Dipstick	15	Positive	1	6.67%
UDS-2		* Denotes Correct Response		
AG: All Methods	95			
		*Negative	93	97.89%
		Positive	2	2.11%
E-Screen/Alere	4	*Negative	4	100.00%
Alere iCup Drugs of Abuse	13	*Negative	12	92.31%
Innovacon Multi-Clin	9	*Negative	9	100.00%
CLIAwaived Inc Rapid Drug Cup	5	*Negative	5	100.00%
Branan Med Corp QuickTox	1	*Negative	1	100.00%
Jant Accutest Multi Drug Screen	2	*Negative	2	100.00%
Amedica Biotech AmediCheck Instant Test Cup	4	*Negative	4	100.00%
Nova Century Scientific Nova Screen DOA Cassette	1	*Negative	1	100.00%
American Screening Corp Reveal Multi-Drug Test Cup	2	*Negative	2	100.00%
Medica Easy RA	1	*Negative	1	100.00%
Alere iCassette	2	*Negative	2	100.00%
Alere USD UScreen Cup	9	*Negative	9	100.00%
Alere iScreen Drugs of Abuse Card	1	*Negative	1	100.00%
Noble Medical Insta-Screen Test Cup	3	*Negative	3	100.00%
Hemosure First Sign UDS Dipstick	15	*Negative	14	93.33%
Siemens Viva/ProE System	1	*Negative	1	100.00%
Premier Biotech Bio-Cup	12	*Negative	12	100.00%
UCP Biosciences Home Drug Screening Test Cup	1	*Negative	1	100.00%
Consult Diag Drugs of Abuse Cups	1	*Negative	1	100.00%
Branan Med Corp ToxCup	6	*Negative	6	100.00%
American Screening Discover Plus Cup	1	*Negative	1	100.00%
McKesson Drugs of Abuse Test Cup	1	*Negative	1	100.00%
Alere iCup Drugs of Abuse	13	Positive	1	7.69%
Hemosure First Sign UDS Dipstick	15	Positive	1	6.67%
		Methadone		
UDS-1		* Denotes Correct Response		
AG: All Methods	119			
		*Negative	118	99.16%
		Positive	1	0.84%
AlereTriage	3	*Negative	3	100.00%



Siemens Dim/AR/ES/RxL/Xpand/EXL	1	*Negative	1	100.00%
E-Screen/Alere	26	*Negative	26	100.00%
Express Diag. Express Check	1	*Negative	1	100.00%
Exp Results Integ Multi-Drug Screen Cup	1	*Negative	1	100.00%
Alere iCup Drugs of Abuse	15	*Negative	15	100.00%
MedTox Profile II	1	*Negative	1	100.00%
Beckman AU	1	*Negative	1	100.00%
MedTox Profile III ER	1	*Negative	1	100.00%
CLIAwaived Inc Rapid Drug Cup	6	*Negative	6	100.00%
MedTox Profile V	2	*Negative	2	100.00%
Branan Med Corp QuickTox	1	*Negative	1	100.00%
Jant Accutest Multi Drug Screen	2	*Negative	2	100.00%
Amedica Biotech AmediCheck Instant Test Cup	4	*Negative	4	100.00%
Nova Century Scientific Nova Screen DOA Cassette	1	*Negative	1	100.00%
American Screening Corp Reveal Multi-Drug Test Cup	4	*Negative	4	100.00%
Mindray BS-200	2	*Negative	2	100.00%
Medica Easy RA	2	*Negative	2	100.00%
Roche Hitachi 717	1	*Negative	1	100.00%
Alere USD UScreen Cup	9	*Negative	9	100.00%
Alere iScreen Drugs of Abuse Card	1	*Negative	1	100.00%
Noble Medical Insta-Screen Test Cup	3	*Negative	3	100.00%
Hemosure First Sign UDS Dipstick	15	*Negative	14	93.33%
Siemens Viva/ProE System	3	*Negative	3	100.00%
TFS DRI Assay	1	*Negative	1	100.00%
UCP Biosciences Home Drug Screening Test Cup	1	*Negative	1	100.00%
Abbott Architect c and ci Series	2	*Negative	2	100.00%
Consult Diag Drugs of Abuse Cups	1	*Negative	1	100.00%
Branan Med Corp ToxCup	6	*Negative	6	100.00%
American Screening Discover Plus Cup	1	*Negative	1	100.00%
McKesson Drugs of Abuse Test Cup	1	*Negative	1	100.00%
Hemosure First Sign UDS Dipstick	15	Positive	1	6.67%
UDS-2		* Denotes Correct Response		
AG: All Methods	119			
		*Negative	118	99.16%
		Positive	1	0.84%
AlereTriage	3	*Negative	3	100.00%
Siemens Dim/AR/ES/RxL/Xpand/EXL	1	*Negative	1	100.00%
E-Screen/Alere	26	*Negative	26	100.00%
Express Diag. Express Check	1	*Negative	1	100.00%
Exp Results Integ Multi-Drug Screen Cup	1	*Negative	1	100.00%
Alere iCup Drugs of Abuse	15	*Negative	15	100.00%
MedTox Profile II	1	*Negative	1	100.00%
Beckman AU	1	*Negative	1	100.00%
MedTox Profile III ER	1	*Negative	1	100.00%
CLIAwaived Inc Rapid Drug Cup	6	*Negative	6	100.00%
MedTox Profile V	2	*Negative	2	100.00%
Branan Med Corp QuickTox	1	*Negative	1	100.00%
Jant Accutest Multi Drug Screen	2	*Negative	2	100.00%
Amedica Biotech AmediCheck Instant Test Cup	4	*Negative	4	100.00%
Nova Century Scientific Nova Screen DOA Cassette	1	*Negative	1	100.00%
American Screening Corp Reveal Multi-Drug Test Cup	4	*Negative	4	100.00%
Mindray BS-200	2	*Negative	2	100.00%
Medica Easy RA	2	*Negative	2	100.00%
Roche Hitachi 717	1	*Negative	1	100.00%
Alere USD UScreen Cup	9	*Negative	9	100.00%
Alere iScreen Drugs of Abuse Card	1	*Negative	1	100.00%
Noble Medical Insta-Screen Test Cup	3	*Negative	3	100.00%
Hemosure First Sign UDS Dipstick	15	*Negative	14	93.33%
Siemens Viva/ProE System	3	*Negative	3	100.00%
TFS DRI Assay	1	*Negative	1	100.00%
UCP Biosciences Home Drug Screening Test Cup	1	*Negative	1	100.00%
Abbott Architect c and ci Series	2	*Negative	2	100.00%
Consult Diag Drugs of Abuse Cups	1	*Negative	1	100.00%
Branan Med Corp ToxCup	6	*Negative	6	100.00%
American Screening Discover Plus Cup	1	*Negative	1	100.00%
McKesson Drugs of Abuse Test Cup	1	*Negative	1	100.00%
Hemosure First Sign UDS Dipstick	15	Positive	1	6.67%
		Methamphetamine		
UDS-1		* Denotes Correct Response		



AG: All Methods	150			
		*Negative	148	98.67%
		Positive	2	1.33%
AlereTriage	3	*Negative	3	100.00%
E-Screen/Alere	50	*Negative	50	100.00%
Express Diag. Express Check	1	*Negative	1	100.00%
Exp Results Integ Multi-Drug Screen Cup	1	*Negative	1	100.00%
Alere iCup Drugs of Abuse	18	*Negative	17	94.44%
Instant Technologies i Screen	1	*Negative	1	100.00%
Integrated E-Z Split Key Cup	1	*Negative	1	100.00%
MedTox Profile II	2	*Negative	2	100.00%
Noble Med Split Specimen Cup	1	*Negative	1	100.00%
MedTox Profile III ER	1	*Negative	1	100.00%
CLIAwaived Inc Rapid Drug Cup	6	*Negative	6	100.00%
MedTox Profile V	2	*Negative	2	100.00%
Branan Med Corp QuickTox	1	*Negative	1	100.00%
Jant Accutest Multi Drug Screen	2	*Negative	2	100.00%
Amedica Biotech AmediCheck Instant Test Cup	4	*Negative	4	100.00%
Nova Century Scientific Nova Screen DOA Cassette	1	*Negative	1	100.00%
American Screening Corp Reveal Multi-Drug Test Cup	4	*Negative	4	100.00%
Alere USD UScreen Cup	9	*Negative	9	100.00%
Alere iScreen Drugs of Abuse Card	1	*Negative	1	100.00%
Noble Medical Insta-Screen Test Cup	3	*Negative	3	100.00%
CRL Stat Cup	1	*Negative	1	100.00%
Hemosure First Sign UDS Dipstick	15	*Negative	14	93.33%
Premier Biotech Bio-Cup	12	*Negative	12	100.00%
UCP Biosciences Home Drug Screening Test Cup	1	*Negative	1	100.00%
Consult Diag Drugs of Abuse Cups	1	*Negative	1	100.00%
Branan Med Corp ToxCup	6	*Negative	6	100.00%
American Screening Discover Plus Cup	1	*Negative	1	100.00%
McKesson Drugs of Abuse Test Cup	1	*Negative	1	100.00%
Alere iCup Drugs of Abuse	18	Positive	1	5.56%
Hemosure First Sign UDS Dipstick	15	Positive	1	6.67%
UDS-2		* Denotes Correct Response		
AG: All Methods	150			
		*Negative	149	99.33%
		Positive	1	0.67%
AlereTriage	3	*Negative	3	100.00%
E-Screen/Alere	50	*Negative	50	100.00%
Express Diag. Express Check	1	*Negative	1	100.00%
Exp Results Integ Multi-Drug Screen Cup	1	*Negative	1	100.00%
Alere iCup Drugs of Abuse	18	*Negative	18	100.00%
Instant Technologies i Screen	1	*Negative	1	100.00%
Integrated E-Z Split Key Cup	1	*Negative	1	100.00%
MedTox Profile II	2	*Negative	2	100.00%
Noble Med Split Specimen Cup	1	*Negative	1	100.00%
MedTox Profile III ER	1	*Negative	1	100.00%
CLIAwaived Inc Rapid Drug Cup	6	*Negative	6	100.00%
MedTox Profile V	2	*Negative	2	100.00%
Branan Med Corp QuickTox	1	*Negative	1	100.00%
Jant Accutest Multi Drug Screen	2	*Negative	2	100.00%
Amedica Biotech AmediCheck Instant Test Cup	4	*Negative	4	100.00%
Nova Century Scientific Nova Screen DOA Cassette	1	*Negative	1	100.00%
American Screening Corp Reveal Multi-Drug Test Cup	4	*Negative	4	100.00%
Alere USD UScreen Cup	9	*Negative	9	100.00%
Alere iScreen Drugs of Abuse Card	1	*Negative	1	100.00%
Noble Medical Insta-Screen Test Cup	3	*Negative	3	100.00%
CRL Stat Cup	1	*Negative	1	100.00%
Hemosure First Sign UDS Dipstick	15	*Negative	14	93.33%
Premier Biotech Bio-Cup	12	*Negative	12	100.00%
UCP Biosciences Home Drug Screening Test Cup	1	*Negative	1	100.00%
Consult Diag Drugs of Abuse Cups	1	*Negative	1	100.00%
Branan Med Corp ToxCup	6	*Negative	6	100.00%
American Screening Discover Plus Cup	1	*Negative	1	100.00%
McKesson Drugs of Abuse Test Cup	1	*Negative	1	100.00%
Hemosure First Sign UDS Dipstick	15	Positive	1	6.67%
		Methaqualone		
UDS-1		* Denotes Correct Response		
AG: All Methods	1			



		*Negative	1	100.00%
Beckman AU	1	*Negative	1	100.00%
UDS-2		* Denotes Correct Response		
AG: All Methods	1			
		*Positive	1	100.00%
Beckman AU	1	*Positive	1	100.00%
		Opiates		
UDS-1		* Denotes Correct Response		
AG: All Methods	181			
		*Negative	164	90.61%
		Positive	17	9.39%
AlereTriage	3	*Negative	3	100.00%
Siemens Dim/AR/ES/RxL/Xpand/EXL	2	*Negative	2	100.00%
E-Screen/Alere	51	*Negative	51	100.00%
Exp Results Integ Multi-Drug Screen Cup	1	*Negative	1	100.00%
Alere iCup Drugs of Abuse	19	*Negative	15	78.95%
Instant Technologies i Screen	2	*Negative	2	100.00%
Integrated E-Z Split Key Cup	1	*Negative	1	100.00%
MedTox Profile II	2	*Negative	2	100.00%
Noble Med Split Specimen Cup	1	*Negative	1	100.00%
Beckman AU	1	*Negative	1	100.00%
Innovacon Multi-Clin	9	*Negative	9	100.00%
CLIAwaived Inc Rapid Drug Cup	6	*Negative	5	83.33%
Branan Med Corp QuickTox	1	*Negative	1	100.00%
Jant Accutest Multi Drug Screen	2	*Negative	2	100.00%
Amedica Biotech AmediCheck Instant Test Cup	4	*Negative	4	100.00%
Nova Century Scientific Nova Screen DOA Cassette	1	*Negative	1	100.00%
American Screening Corp Reveal Multi-Drug Test Cup	4	*Negative	4	100.00%
Mindray BS-200	2	*Negative	2	100.00%
Medica Easy RA	2	*Negative	2	100.00%
Alere iCassette	4	*Negative	4	100.00%
Roche Hitachi 717	1	*Negative	1	100.00%
Alere USD UScreen Cup	8	*Negative	8	100.00%
ThermoScientific Indiko Plus	1	*Negative	1	100.00%
Alere iScreen Drugs of Abuse Card	1	*Negative	1	100.00%
Noble Medical Insta-Screen Test Cup	3	*Negative	3	100.00%
CRL Stat Cup	1	*Negative	1	100.00%
Hemosure First Sign UDS Dipstick	15	*Negative	14	93.33%
Siemens Viva/ProE System	4	*Negative	4	100.00%
TFS DRI Assay	1	*Negative	1	100.00%
Premier Biotech Bio-Cup	12	*Negative	7	58.33%
UCP Biosciences Home Drug Screening Test Cup	1	*Negative	1	100.00%
Abbott Architect c and ci Series	2	*Negative	2	100.00%
Consult Diag Drugs of Abuse Cups	1	*Negative	1	100.00%
Branan Med Corp ToxCup	6	*Negative	4	66.67%
American Screening Discover Plus Cup	1	*Negative	1	100.00%
McKesson Drugs of Abuse Test Cup	1	*Negative	1	100.00%
Express Diag. Express Check	1	Positive	1	100.00%
Alere iCup Drugs of Abuse	19	Positive	4	21.05%
MedTox Profile III ER	1	Positive	1	100.00%
CLIAwaived Inc Rapid Drug Cup	6	Positive	1	16.67%
MedTox Profile V	2	Positive	2	100.00%
Hemosure First Sign UDS Dipstick	15	Positive	1	6.67%
Premier Biotech Bio-Cup	12	Positive	5	41.67%
Branan Med Corp ToxCup	6	Positive	2	33.33%
UDS-2		* Denotes Correct Response		
AG: All Methods	182			
		*Negative	181	99.45%
		Positive	1	0.55%
AlereTriage	3	*Negative	3	100.00%
Siemens Dim/AR/ES/RxL/Xpand/EXL	2	*Negative	2	100.00%
E-Screen/Alere	51	*Negative	51	100.00%
Express Diag. Express Check	1	*Negative	1	100.00%
Exp Results Integ Multi-Drug Screen Cup	1	*Negative	1	100.00%
Alere iCup Drugs of Abuse	19	*Negative	19	100.00%
Instant Technologies i Screen	2	*Negative	2	100.00%
Integrated E-Z Split Key Cup	1	*Negative	1	100.00%
MedTox Profile II	2	*Negative	2	100.00%
Noble Med Split Specimen Cup	1	*Negative	1	100.00%



Beckman AU	1	*Negative	1	100.00%
Innovacon Multi-Clin	9	*Negative	9	100.00%
MedTox Profile III ER	1	*Negative	1	100.00%
CLIAwaived Inc Rapid Drug Cup	6	*Negative	6	100.00%
MedTox Profile V	2	*Negative	2	100.00%
Branan Med Corp QuickTox	1	*Negative	1	100.00%
Jant Accutest Multi Drug Screen	2	*Negative	2	100.00%
Amedica Biotech AmediCheck Instant Test Cup	4	*Negative	4	100.00%
Nova Century Scientific Nova Screen DOA Cassette	1	*Negative	1	100.00%
American Screening Corp Reveal Multi-Drug Test Cup	4	*Negative	4	100.00%
Mindray BS-200	2	*Negative	2	100.00%
Medica Easy RA	2	*Negative	2	100.00%
Alere iCassette	4	*Negative	4	100.00%
Roche Hitachi 717	1	*Negative	1	100.00%
Alere USD UScreen Cup	8	*Negative	8	100.00%
ThermoScientific Indiko Plus	1	*Negative	1	100.00%
Alere iScreen Drugs of Abuse Card	1	*Negative	1	100.00%
Noble Medical Insta-Screen Test Cup	3	*Negative	3	100.00%
CRL Stat Cup	1	*Negative	1	100.00%
Hemosure First Sign UDS Dipstick	15	*Negative	14	93.33%
Siemens Viva/ProE System	5	*Negative	5	100.00%
TFS DRI Assay	1	*Negative	1	100.00%
Premier Biotech Bio-Cup	12	*Negative	12	100.00%
UCP Biosciences Home Drug Screening Test Cup	1	*Negative	1	100.00%
Abbott Architect c and ci Series	2	*Negative	2	100.00%
Consult Diag Drugs of Abuse Cups	1	*Negative	1	100.00%
Branan Med Corp ToxCup	6	*Negative	6	100.00%
American Screening Discover Plus Cup	1	*Negative	1	100.00%
McKesson Drugs of Abuse Test Cup	1	*Negative	1	100.00%
Hemosure First Sign UDS Dipstick	15	Positive	1	6.67%
		Oxycodone		
UDS-1		* Denotes Correct Response		
AG: All Methods	112			
		*Negative	107	95.54%
		Positive	5	4.46%
E-Screen/Alere	18	*Negative	17	94.44%
Exp Results Integ Multi-Drug Screen Cup	1	*Negative	1	100.00%
Alere iCup Drugs of Abuse	15	*Negative	13	86.67%
MedTox Profile II	1	*Negative	1	100.00%
Beckman AU	1	*Negative	1	100.00%
Innovacon Multi-Clin	9	*Negative	9	100.00%
MedTox Profile III ER	1	*Negative	1	100.00%
CLIAwaived Inc Rapid Drug Cup	6	*Negative	5	83.33%
MedTox Profile V	2	*Negative	2	100.00%
Branan Med Corp QuickTox	1	*Negative	1	100.00%
Jant Accutest Multi Drug Screen	2	*Negative	2	100.00%
Amedica Biotech AmediCheck Instant Test Cup	4	*Negative	4	100.00%
Nova Century Scientific Nova Screen DOA Cassette	1	*Negative	1	100.00%
American Screening Corp Reveal Multi-Drug Test Cup	4	*Negative	4	100.00%
Mindray BS-200	2	*Negative	2	100.00%
Medica Easy RA	2	*Negative	2	100.00%
Alere iCassette	2	*Negative	2	100.00%
Roche Hitachi 717	1	*Negative	1	100.00%
Alere USD UScreen Cup	9	*Negative	8	88.89%
Alere iScreen Drugs of Abuse Card	1	*Negative	1	100.00%
Noble Medical Insta-Screen Test Cup	3	*Negative	3	100.00%
Hemosure First Sign UDS Cup	1	*Negative	1	100.00%
Siemens Viva/ProE System	2	*Negative	2	100.00%
TFS DRI Assay	1	*Negative	1	100.00%
Premier Biotech Bio-Cup	12	*Negative	12	100.00%
UCP Biosciences Home Drug Screening Test Cup	1	*Negative	1	100.00%
Consult Diag Drugs of Abuse Cups	1	*Negative	1	100.00%
Branan Med Corp ToxCup	6	*Negative	6	100.00%
American Screening Discover Plus Cup	1	*Negative	1	100.00%
McKesson Drugs of Abuse Test Cup	1	*Negative	1	100.00%
E-Screen/Alere	18	Positive	1	5.56%
Alere iCup Drugs of Abuse	15	Positive	2	13.33%
CLIAwaived Inc Rapid Drug Cup	6	Positive	1	16.67%
Alere USD UScreen Cup	9	Positive	1	11.11%



UDS-2		* Denotes Correct Response		
AG: All Methods	112			
		*Positive	111	99.11%
		Negative	1	0.89%
E-Screen/Alere	18	*Positive	18	100.00%
Exp Results Integ Multi-Drug Screen Cup	1	*Positive	1	100.00%
Alere iCup Drugs of Abuse	15	*Positive	15	100.00%
MedTox Profile II	1	*Positive	1	100.00%
Beckman AU	1	*Positive	1	100.00%
Innovacon Multi-Clin	9	*Positive	9	100.00%
MedTox Profile III ER	1	*Positive	1	100.00%
CLIAwaived Inc Rapid Drug Cup	6	*Positive	6	100.00%
MedTox Profile V	2	*Positive	2	100.00%
Branan Med Corp QuickTox	1	*Positive	1	100.00%
Jant Accutest Multi Drug Screen	2	*Positive	2	100.00%
Amedica Biotech AmediCheck Instant Test Cup	4	*Positive	4	100.00%
Nova Century Scientific Nova Screen DOA Cassette	1	*Positive	1	100.00%
American Screening Corp Reveal Multi-Drug Test Cup	4	*Positive	4	100.00%
Mindray BS-200	2	*Positive	2	100.00%
Medica Easy RA	2	*Positive	2	100.00%
Alere iCassette	2	*Positive	2	100.00%
Roche Hitachi 717	1	*Positive	1	100.00%
Alere USD UScreen Cup	9	*Positive	8	88.89%
Alere iScreen Drugs of Abuse Card	1	*Positive	1	100.00%
Noble Medical Insta-Screen Test Cup	3	*Positive	3	100.00%
Hemosure First Sign UDS Cup	1	*Positive	1	100.00%
Siemens Viva/ProE System	2	*Positive	2	100.00%
TFS DRI Assay	1	*Positive	1	100.00%
Premier Biotech Bio-Cup	12	*Positive	12	100.00%
UCP Biosciences Home Drug Screening Test Cup	1	*Positive	1	100.00%
Consult Diag Drugs of Abuse Cups	1	*Positive	1	100.00%
Branan Med Corp ToxCup	6	*Positive	6	100.00%
American Screening Discover Plus Cup	1	*Positive	1	100.00%
McKesson Drugs of Abuse Test Cup	1	*Positive	1	100.00%
Alere USD UScreen Cup	9	Negative	1	11.11%
		PCP (Phencyclidine)		
UDS-1		* Denotes Correct Response		
AG: All Methods	163			
		*Positive	156	95.71%
		Negative	7	4.29%
AlereTriage	3	*Positive	3	100.00%
Siemens Dim/AR/ES/RxL/Xpand/EXL	1	*Positive	1	100.00%
E-Screen/Alere	51	*Positive	50	98.04%
Express Diag. Express Check	1	*Positive	1	100.00%
Exp Results Integ Multi-Drug Screen Cup	1	*Positive	1	100.00%
Alere iCup Drugs of Abuse	13	*Positive	12	92.31%
Instant Technologies i Screen	2	*Positive	2	100.00%
Integrated E-Z Split Key Cup	1	*Positive	1	100.00%
MedTox Profile II	3	*Positive	3	100.00%
Beckman AU	1	*Positive	1	100.00%
Innovacon Multi-Clin	9	*Positive	9	100.00%
MedTox Profile III ER	1	*Positive	1	100.00%
CLIAwaived Inc Rapid Drug Cup	6	*Positive	6	100.00%
MedTox Profile V	2	*Positive	2	100.00%
Branan Med Corp QuickTox	1	*Positive	1	100.00%
Jant Accutest Multi Drug Screen	2	*Positive	2	100.00%
Amedica Biotech AmediCheck Instant Test Cup	1	*Positive	1	100.00%
Nova Century Scientific Nova Screen DOA Cassette	1	*Positive	1	100.00%
American Screening Corp Reveal Multi-Drug Test Cup	4	*Positive	4	100.00%
Mindray BS-200	1	*Positive	1	100.00%
Medica Easy RA	1	*Positive	1	100.00%
Alere iCassette	4	*Positive	4	100.00%
Roche Hitachi 717	1	*Positive	1	100.00%
Alere USD UScreen Cup	6	*Positive	6	100.00%
Alere iScreen Drugs of Abuse Card	1	*Positive	1	100.00%
Noble Medical Insta-Screen Test Cup	3	*Positive	2	66.67%
Hemosure First Sign UDS Dipstick	15	*Positive	14	93.33%
TFS DRI Assay	1	*Positive	1	100.00%
Premier Biotech Bio-Cup	12	*Positive	12	100.00%



UCP Biosciences Home Drug Screening Test Cup	1	*Positive	1	100.00%
Abbott Architect c and ci Series	2	*Positive	2	100.00%
Consult Diag Drugs of Abuse Cups	1	*Positive	1	100.00%
Branan Med Corp ToxCup	6	*Positive	6	100.00%
McKesson Drugs of Abuse Test Cup	1	*Positive	1	100.00%
E-Screen/Alere	51	Negative	1	1.96%
Alere iCup Drugs of Abuse	13	Negative	1	7.69%
Noble Medical Insta-Screen Test Cup	3	Negative	1	33.33%
CRL Stat Cup	1	Negative	1	100.00%
Hemosure First Sign UDS Dipstick	15	Negative	1	6.67%
Siemens Viva/ProE System	1	Negative	1	100.00%
American Screening Discover Plus Cup	1	Negative	1	100.00%
UDS-2		* Denotes Correct Response		
AG: All Methods	163			
		*Negative	162	99.39%
		Positive	1	0.61%
AlereTriage	3	*Negative	3	100.00%
Siemens Dim/AR/ES/RxL/Xpand/EXL	1	*Negative	1	100.00%
E-Screen/Alere	51	*Negative	51	100.00%
Express Diag. Express Check	1	*Negative	1	100.00%
Exp Results Integ Multi-Drug Screen Cup	1	*Negative	1	100.00%
Alere iCup Drugs of Abuse	13	*Negative	13	100.00%
Instant Technologies i Screen	2	*Negative	2	100.00%
Integrated E-Z Split Key Cup	1	*Negative	1	100.00%
MedTox Profile II	3	*Negative	3	100.00%
Beckman AU	1	*Negative	1	100.00%
Innovacon Multi-Clin	9	*Negative	9	100.00%
MedTox Profile III ER	1	*Negative	1	100.00%
CLIAwaived Inc Rapid Drug Cup	6	*Negative	6	100.00%
MedTox Profile V	2	*Negative	2	100.00%
Branan Med Corp QuickTox	1	*Negative	1	100.00%
Jant Accutest Multi Drug Screen	2	*Negative	2	100.00%
Amedica Biotech AmediCheck Instant Test Cup	1	*Negative	1	100.00%
Nova Century Scientific Nova Screen DOA Cassette	1	*Negative	1	100.00%
American Screening Corp Reveal Multi-Drug Test Cup	4	*Negative	4	100.00%
Mindray BS-200	1	*Negative	1	100.00%
Medica Easy RA	1	*Negative	1	100.00%
Alere iCassette	4	*Negative	4	100.00%
Roche Hitachi 717	1	*Negative	1	100.00%
Alere USD UScreen Cup	6	*Negative	6	100.00%
Alere iScreen Drugs of Abuse Card	1	*Negative	1	100.00%
Noble Medical Insta-Screen Test Cup	3	*Negative	3	100.00%
CRL Stat Cup	1	*Negative	1	100.00%
Hemosure First Sign UDS Dipstick	15	*Negative	14	93.33%
Siemens Viva/ProE System	1	*Negative	1	100.00%
TFS DRI Assay	1	*Negative	1	100.00%
Premier Biotech Bio-Cup	12	*Negative	12	100.00%
UCP Biosciences Home Drug Screening Test Cup	1	*Negative	1	100.00%
Abbott Architect c and ci Series	2	*Negative	2	100.00%
Consult Diag Drugs of Abuse Cups	1	*Negative	1	100.00%
Branan Med Corp ToxCup	6	*Negative	6	100.00%
American Screening Discover Plus Cup	1	*Negative	1	100.00%
McKesson Drugs of Abuse Test Cup	1	*Negative	1	100.00%
Hemosure First Sign UDS Dipstick	15	Positive	1	6.67%
		Propoxyphene		
UDS-1		* Denotes Correct Response		
AG: All Methods	39			
		*Negative	37	94.87%
		Positive	2	5.13%
E-Screen/Alere	4	*Negative	4	100.00%
Exp Results Integ Multi-Drug Screen Cup	1	*Negative	1	100.00%
Alere iCup Drugs of Abuse	11	*Negative	9	81.82%
MedTox Profile II	1	*Negative	1	100.00%
Beckman AU	1	*Negative	1	100.00%
Innovacon Multi-Clin	7	*Negative	7	100.00%
CLIAwaived Inc Rapid Drug Cup	1	*Negative	1	100.00%
MedTox Profile V	2	*Negative	2	100.00%
Alere iCassette	2	*Negative	2	100.00%
Abbott Architect c and ci Series	1	*Negative	1	100.00%



Consult Diag Drugs of Abuse Cups	1	*Negative	1	100.00%
Branan Med Corp ToxCup	6	*Negative	6	100.00%
McKesson Drugs of Abuse Test Cup	1	*Negative	1	100.00%
Alere iCup Drugs of Abuse	11	Positive	2	18.18%
UDS-2		* Denotes Correct Response		
AG: All Methods	39			
		*Negative	39	100.00%
E-Screen/Alere	4	*Negative	4	100.00%
Exp Results Integ Multi-Drug Screen Cup	1	*Negative	1	100.00%
Alere iCup Drugs of Abuse	11	*Negative	11	100.00%
MedTox Profile II	1	*Negative	1	100.00%
Beckman AU	1	*Negative	1	100.00%
Innovacon Multi-Clin	7	*Negative	7	100.00%
CLIAwaived Inc Rapid Drug Cup	1	*Negative	1	100.00%
MedTox Profile V	2	*Negative	2	100.00%
Alere iCassette	2	*Negative	2	100.00%
Abbott Architect c and ci Series	1	*Negative	1	100.00%
Consult Diag Drugs of Abuse Cups	1	*Negative	1	100.00%
Branan Med Corp ToxCup	6	*Negative	6	100.00%
McKesson Drugs of Abuse Test Cup	1	*Negative	1	100.00%
		Tricyclic Antidepressants		
UDS-1		* Denotes Correct Response		
AG: All Methods	53			
		*Negative	53	100.00%
AlereTriage	3	*Negative	3	100.00%
Alere iCup Drugs of Abuse	12	*Negative	12	100.00%
Innovacon Multi-Clin	8	*Negative	8	100.00%
CLIAwaived Inc Rapid Drug Cup	1	*Negative	1	100.00%
MedTox Profile V	2	*Negative	2	100.00%
Branan Med Corp QuickTox	1	*Negative	1	100.00%
Jant Accutest Multi Drug Screen	2	*Negative	2	100.00%
Amedica Biotech AmediCheck Instant Test Cup	1	*Negative	1	100.00%
Nova Century Scientific Nova Screen DOA Cassette	1	*Negative	1	100.00%
Alere iCassette	2	*Negative	2	100.00%
Alere USD UScreen Cup	6	*Negative	6	100.00%
Alere iScreen Drugs of Abuse Card	1	*Negative	1	100.00%
Noble Medical Insta-Screen Test Cup	3	*Negative	3	100.00%
UCP Biosciences Home Drug Screening Test Cup	1	*Negative	1	100.00%
Consult Diag Drugs of Abuse Cups	1	*Negative	1	100.00%
Branan Med Corp ToxCup	6	*Negative	6	100.00%
American Screening Discover Plus Cup	1	*Negative	1	100.00%
McKesson Drugs of Abuse Test Cup	1	*Negative	1	100.00%
UDS-2		* Denotes Correct Response		
AG: All Methods	53			
		*Negative	53	100.00%
AlereTriage	3	*Negative	3	100.00%
Alere iCup Drugs of Abuse	12	*Negative	12	100.00%
Innovacon Multi-Clin	8	*Negative	8	100.00%
CLIAwaived Inc Rapid Drug Cup	1	*Negative	1	100.00%
MedTox Profile V	2	*Negative	2	100.00%
Branan Med Corp QuickTox	1	*Negative	1	100.00%
Jant Accutest Multi Drug Screen	2	*Negative	2	100.00%
Amedica Biotech AmediCheck Instant Test Cup	1	*Negative	1	100.00%
Nova Century Scientific Nova Screen DOA Cassette	1	*Negative	1	100.00%
Alere iCassette	2	*Negative	2	100.00%
Alere USD UScreen Cup	6	*Negative	6	100.00%
Alere iScreen Drugs of Abuse Card	1	*Negative	1	100.00%
Noble Medical Insta-Screen Test Cup	3	*Negative	3	100.00%
UCP Biosciences Home Drug Screening Test Cup	1	*Negative	1	100.00%
Consult Diag Drugs of Abuse Cups	1	*Negative	1	100.00%
Branan Med Corp ToxCup	6	*Negative	6	100.00%
American Screening Discover Plus Cup	1	*Negative	1	100.00%
McKesson Drugs of Abuse Test Cup	1	*Negative	1	100.00%
		Buprenorphine		
UDS-1		* Denotes Correct Response		
AG: All Methods	58			
		*Negative	57	98.28%
		Positive	1	1.72%
Exp Results Integ Multi-Drug Screen Cup	1	*Negative	1	100.00%



Alere iCup Drugs of Abuse	12	*Negative	12	100.00%
CLIAwaived Inc Rapid Drug Cup	5	*Negative	5	100.00%
MedTox Profile V	1	*Negative	1	100.00%
Jant Accutest Multi Drug Screen	1	*Negative	1	100.00%
Mindray BS-200	2	*Negative	2	100.00%
Medica Easy RA	2	*Negative	2	100.00%
Alere USD UScreen Cup	2	*Negative	2	100.00%
Alere iScreen Drugs of Abuse Card	1	*Negative	1	100.00%
Noble Medical Insta-Screen Test Cup	3	*Negative	3	100.00%
Hemosure First Sign UDS Cup	1	*Negative	1	100.00%
Siemens Viva/ProE System	5	*Negative	4	80.00%
TFS DRI Assay	1	*Negative	1	100.00%
Premier Biotech Bio-Cup	12	*Negative	12	100.00%
Consult Diag Drugs of Abuse Cups	1	*Negative	1	100.00%
Branan Med Corp ToxCup	6	*Negative	6	100.00%
American Screening Discover Plus Cup	1	*Negative	1	100.00%
McKesson Drugs of Abuse Test Cup	1	*Negative	1	100.00%
Siemens Viva/ProE System	5	Positive	1	20.00%
UDS-2		* Denotes Correct Response		
AG: All Methods	58			
		*Negative	58	100.00%
Exp Results Integ Multi-Drug Screen Cup	1	*Negative	1	100.00%
Alere iCup Drugs of Abuse	12	*Negative	12	100.00%
CLIAwaived Inc Rapid Drug Cup	5	*Negative	5	100.00%
MedTox Profile V	1	*Negative	1	100.00%
Jant Accutest Multi Drug Screen	1	*Negative	1	100.00%
Mindray BS-200	2	*Negative	2	100.00%
Medica Easy RA	2	*Negative	2	100.00%
Alere USD UScreen Cup	2	*Negative	2	100.00%
Alere iScreen Drugs of Abuse Card	1	*Negative	1	100.00%
Noble Medical Insta-Screen Test Cup	3	*Negative	3	100.00%
Hemosure First Sign UDS Cup	1	*Negative	1	100.00%
Siemens Viva/ProE System	5	*Negative	5	100.00%
TFS DRI Assay	1	*Negative	1	100.00%
Premier Biotech Bio-Cup	12	*Negative	12	100.00%
Consult Diag Drugs of Abuse Cups	1	*Negative	1	100.00%
Branan Med Corp ToxCup	6	*Negative	6	100.00%
American Screening Discover Plus Cup	1	*Negative	1	100.00%
McKesson Drugs of Abuse Test Cup	1	*Negative	1	100.00%
		6-Monoacetylmorphine (6 AM)		
UDS-1		* Denotes Correct Response		
AG: All Methods	4			
		*Negative	4	100.00%
Siemens Viva/ProE System	4	*Negative	4	100.00%
UDS-2		* Denotes Correct Response		
AG: All Methods	4			
		*Negative	4	100.00%
Siemens Viva/ProE System	4	*Negative	4	100.00%
		EDDP		
UDS-1		* Denotes Correct Response		
AG: All Methods	15			
		*Negative	15	100.00%
Siemens Viva/ProE System	3	*Negative	3	100.00%
Premier Biotech Bio-Cup	12	*Negative	12	100.00%
UDS-2		* Denotes Correct Response		
AG: All Methods	15			
		*Negative	15	100.00%
Siemens Viva/ProE System	3	*Negative	3	100.00%
Premier Biotech Bio-Cup	12	*Negative	12	100.00%
756 - HSV Antibodies	Labs	Response	Number	Percent
		HSV-2 Antibody		
VA-1		* Denotes Correct Response		
AG: All Methods	1			
		*Negative	1	100.00%
Biokit HSV-II	1	*Negative	1	100.00%
VA-2		* Denotes Correct Response		
AG: All Methods	1			
		*Positive	1	100.00%



Biokit HSV-II	1	*Positive	1	100.00%
770 - Viral Antigen Detection;792 - Viral Antigen Detection - Waived Methods	Labs	Response	Number	Percent
		Influenza A - Aq Detect		
VAD-1		* Denotes Correct Response		
AG: All Methods	554			
		*Negative	536	96.75%
		Positive	18	3.25%
Binax Now Flu A	3	*Negative	3	100.00%
Alere BinaxNOW Influenza A+B Card	44	*Negative	43	97.73%
Sekisui Diagnostics OSOM Influenza A & B	16	*Negative	16	100.00%
Quidel QuickVue A+B	42	*Negative	40	95.24%
Quidel Sofia/Sofia 2 Influenza A+B	187	*Negative	182	97.33%
Alere Influenza A & B	7	*Negative	6	85.71%
BD Veritor Flu A+B	68	*Negative	66	97.06%
Henry Schein OneStep + Influenza A&B (distinguish)	8	*Negative	8	100.00%
Medline Influenza A and B	1	*Negative	1	100.00%
Cepheid GeneXpert - Xpert Flu	21	*Negative	21	100.00%
Alere ID NOW Influenza A & B	22	*Negative	20	90.91%
Roche Cobas Liat Influenza A/B	7	*Negative	7	100.00%
Sekisui Diagnostics OSOM Ultra Influenza A & B	84	*Negative	81	96.43%
Consult Diag Influenza A & B	21	*Negative	21	100.00%
Lifesign Status Influenza A & B	23	*Negative	21	91.30%
Alere BinaxNOW Influenza A+B Card	44	Positive	1	2.27%
Quidel QuickVue A+B	42	Positive	2	4.76%
Quidel Sofia/Sofia 2 Influenza A+B	187	Positive	5	2.67%
Alere Influenza A & B	7	Positive	1	14.29%
BD Veritor Flu A+B	68	Positive	2	2.94%
Alere ID NOW Influenza A & B	22	Positive	2	9.09%
Sekisui Diagnostics OSOM Ultra Influenza A & B	84	Positive	3	3.57%
Lifesign Status Influenza A & B	23	Positive	2	8.70%
VAD-2		* Denotes Correct Response		
AG: All Methods	556			
		*Positive	543	97.66%
		Negative	13	2.34%
Binax Now Flu A	3	*Positive	3	100.00%
Alere BinaxNOW Influenza A+B Card	44	*Positive	43	97.73%
Sekisui Diagnostics OSOM Influenza A & B	16	*Positive	15	93.75%
Quidel QuickVue A+B	42	*Positive	41	97.62%
Quidel Sofia/Sofia 2 Influenza A+B	188	*Positive	184	97.87%
Alere Influenza A & B	7	*Positive	6	85.71%
BD Veritor Flu A+B	68	*Positive	68	100.00%
Henry Schein OneStep + Influenza A&B (distinguish)	8	*Positive	8	100.00%
Medline Influenza A and B	1	*Positive	1	100.00%
Cepheid GeneXpert - Xpert Flu	21	*Positive	21	100.00%
Alere ID NOW Influenza A & B	22	*Positive	21	95.45%
Roche Cobas Liat Influenza A/B	7	*Positive	7	100.00%
Sekisui Diagnostics OSOM Ultra Influenza A & B	84	*Positive	82	97.62%
Consult Diag Influenza A & B	21	*Positive	21	100.00%
Lifesign Status Influenza A & B	24	*Positive	22	91.67%
Alere BinaxNOW Influenza A+B Card	44	Negative	1	2.27%
Sekisui Diagnostics OSOM Influenza A & B	16	Negative	1	6.25%
Quidel QuickVue A+B	42	Negative	1	2.38%
Quidel Sofia/Sofia 2 Influenza A+B	188	Negative	4	2.13%
Alere Influenza A & B	7	Negative	1	14.29%
Alere ID NOW Influenza A & B	22	Negative	1	4.55%
Sekisui Diagnostics OSOM Ultra Influenza A & B	84	Negative	2	2.38%
Lifesign Status Influenza A & B	24	Negative	2	8.33%
VAD-3		* Denotes Correct Response		
AG: All Methods	48			
		*Negative	48	100.00%
Alere BinaxNOW Influenza A+B Card	5	*Negative	5	100.00%
Sekisui Diagnostics OSOM Influenza A & B	8	*Negative	8	100.00%
Quidel QuickVue A+B	12	*Negative	12	100.00%
Quidel Sofia/Sofia 2 Influenza A+B	7	*Negative	7	100.00%
BD Veritor Flu A+B	5	*Negative	5	100.00%
Henry Schein OneStep + Influenza A&B (distinguish)	1	*Negative	1	100.00%
Cepheid GeneXpert - Xpert Flu	6	*Negative	6	100.00%
Sekisui Diagnostics OSOM Ultra Influenza A & B	1	*Negative	1	100.00%



Consult Diag Influenza A & B	3	*Negative	3	100.00%
VAD-4		* Denotes Correct Response		
AG: All Methods	48			
		*Positive	48	100.00%
Alere BinaxNOW Influenza A+B Card	5	*Positive	5	100.00%
Sekisui Diagnostics OSOM Influenza A & B	8	*Positive	8	100.00%
Quidel QuickVue A+B	12	*Positive	12	100.00%
Quidel Sofia/Sofia 2 Influenza A+B	7	*Positive	7	100.00%
BD Veritor Flu A+B	5	*Positive	5	100.00%
Henry Schein OneStep + Influenza A&B (distinguish)	1	*Positive	1	100.00%
Cepheid GeneXpert - Xpert Flu	6	*Positive	6	100.00%
Sekisui Diagnostics OSOM Ultra Influenza A & B	1	*Positive	1	100.00%
Consult Diag Influenza A & B	3	*Positive	3	100.00%
VAD-5		* Denotes Correct Response		
AG: All Methods	48			
		*Negative	47	97.92%
		Positive	1	2.08%
Alere BinaxNOW Influenza A+B Card	5	*Negative	5	100.00%
Sekisui Diagnostics OSOM Influenza A & B	8	*Negative	8	100.00%
Quidel QuickVue A+B	12	*Negative	12	100.00%
Quidel Sofia/Sofia 2 Influenza A+B	7	*Negative	7	100.00%
BD Veritor Flu A+B	5	*Negative	5	100.00%
Henry Schein OneStep + Influenza A&B (distinguish)	1	*Negative	1	100.00%
Cepheid GeneXpert - Xpert Flu	6	*Negative	5	83.33%
Sekisui Diagnostics OSOM Ultra Influenza A & B	1	*Negative	1	100.00%
Consult Diag Influenza A & B	3	*Negative	3	100.00%
Cepheid GeneXpert - Xpert Flu	6	Positive	1	16.67%
		Influenza B - Ag Detect		
VAD-1		* Denotes Correct Response		
AG: All Methods	556			
		*Positive	546	98.20%
		Negative	10	1.80%
Binax Now Flu B	3	*Positive	3	100.00%
Alere BinaxNOW Influenza A+B Card	44	*Positive	43	97.73%
Sekisui Diagnostics OSOM Influenza A & B	16	*Positive	16	100.00%
Quidel QuickVue A+B	42	*Positive	41	97.62%
Quidel Sofia/Sofia 2 Influenza A+B	188	*Positive	184	97.87%
Alere Influenza A & B	7	*Positive	6	85.71%
BD Veritor Flu A+B	68	*Positive	68	100.00%
Henry Schein OneStep + Influenza A&B (distinguish)	8	*Positive	8	100.00%
Medline Influenza A and B	1	*Positive	1	100.00%
Cepheid GeneXpert - Xpert Flu	21	*Positive	21	100.00%
Alere ID NOW Influenza A & B	22	*Positive	22	100.00%
Roche Cobas Liat Influenza A/B	7	*Positive	7	100.00%
Sekisui Diagnostics OSOM Ultra Influenza A & B	84	*Positive	83	98.81%
Consult Diag Influenza A & B	21	*Positive	21	100.00%
Lifesign Status Influenza A & B	24	*Positive	22	91.67%
Alere BinaxNOW Influenza A+B Card	44	Negative	1	2.27%
Quidel QuickVue A+B	42	Negative	1	2.38%
Quidel Sofia/Sofia 2 Influenza A+B	188	Negative	4	2.13%
Alere Influenza A & B	7	Negative	1	14.29%
Sekisui Diagnostics OSOM Ultra Influenza A & B	84	Negative	1	1.19%
Lifesign Status Influenza A & B	24	Negative	2	8.33%
VAD-2		* Denotes Correct Response		
AG: All Methods	554			
		*Negative	541	97.65%
		Positive	13	2.35%
Binax Now Flu B	3	*Negative	3	100.00%
Alere BinaxNOW Influenza A+B Card	44	*Negative	43	97.73%
Sekisui Diagnostics OSOM Influenza A & B	16	*Negative	15	93.75%
Quidel QuickVue A+B	42	*Negative	41	97.62%
Quidel Sofia/Sofia 2 Influenza A+B	187	*Negative	182	97.33%
Alere Influenza A & B	7	*Negative	6	85.71%
BD Veritor Flu A+B	68	*Negative	68	100.00%
Henry Schein OneStep + Influenza A&B (distinguish)	8	*Negative	8	100.00%
Medline Influenza A and B	1	*Negative	1	100.00%
Cepheid GeneXpert - Xpert Flu	21	*Negative	21	100.00%
Alere ID NOW Influenza A & B	22	*Negative	21	95.45%
Roche Cobas Liat Influenza A/B	7	*Negative	7	100.00%



Sekisui Diagnostics OSOM Ultra Influenza A & B	84	*Negative	82	97.62%
Consult Diag Influenza A & B	21	*Negative	21	100.00%
Lifesign Status Influenza A & B	23	*Negative	22	95.65%
Alere BinaxNOW Influenza A+B Card	44	Positive	1	2.27%
Sekisui Diagnostics OSOM Influenza A & B	16	Positive	1	6.25%
Quidel QuickVue A+B	42	Positive	1	2.38%
Quidel Sofia/Sofia 2 Influenza A+B	187	Positive	5	2.67%
Alere Influenza A & B	7	Positive	1	14.29%
Alere ID NOW Influenza A & B	22	Positive	1	4.55%
Sekisui Diagnostics OSOM Ultra Influenza A & B	84	Positive	2	2.38%
Lifesign Status Influenza A & B	23	Positive	1	4.35%
VAD-3		* Denotes Correct Response		
AG: All Methods	48			
		*Negative	47	97.92%
		Positive	1	2.08%
Alere BinaxNOW Influenza A+B Card	5	*Negative	5	100.00%
Sekisui Diagnostics OSOM Influenza A & B	8	*Negative	8	100.00%
Quidel QuickVue A+B	12	*Negative	11	91.67%
Quidel Sofia/Sofia 2 Influenza A+B	7	*Negative	7	100.00%
BD Veritor Flu A+B	5	*Negative	5	100.00%
Henry Schein OneStep + Influenza A&B (distinguish)	1	*Negative	1	100.00%
Cepheid GeneXpert - Xpert Flu	6	*Negative	6	100.00%
Sekisui Diagnostics OSOM Ultra Influenza A & B	1	*Negative	1	100.00%
Consult Diag Influenza A & B	3	*Negative	3	100.00%
Quidel QuickVue A+B	12	Positive	1	8.33%
VAD-4		* Denotes Correct Response		
AG: All Methods	48			
		*Negative	48	100.00%
Alere BinaxNOW Influenza A+B Card	5	*Negative	5	100.00%
Sekisui Diagnostics OSOM Influenza A & B	8	*Negative	8	100.00%
Quidel QuickVue A+B	12	*Negative	12	100.00%
Quidel Sofia/Sofia 2 Influenza A+B	7	*Negative	7	100.00%
BD Veritor Flu A+B	5	*Negative	5	100.00%
Henry Schein OneStep + Influenza A&B (distinguish)	1	*Negative	1	100.00%
Cepheid GeneXpert - Xpert Flu	6	*Negative	6	100.00%
Sekisui Diagnostics OSOM Ultra Influenza A & B	1	*Negative	1	100.00%
Consult Diag Influenza A & B	3	*Negative	3	100.00%
VAD-5		* Denotes Correct Response		
AG: All Methods	48			
		*Positive	48	100.00%
Alere BinaxNOW Influenza A+B Card	5	*Positive	5	100.00%
Sekisui Diagnostics OSOM Influenza A & B	8	*Positive	8	100.00%
Quidel QuickVue A+B	12	*Positive	12	100.00%
Quidel Sofia/Sofia 2 Influenza A+B	7	*Positive	7	100.00%
BD Veritor Flu A+B	5	*Positive	5	100.00%
Henry Schein OneStep + Influenza A&B (distinguish)	1	*Positive	1	100.00%
Cepheid GeneXpert - Xpert Flu	6	*Positive	6	100.00%
Sekisui Diagnostics OSOM Ultra Influenza A & B	1	*Positive	1	100.00%
Consult Diag Influenza A & B	3	*Positive	3	100.00%
		Influenza A/B - Aq Detect		
VAD-1		* Denotes Correct Response		
AG: All Methods	17			
		*Positive	17	100.00%
Alere BinaxNOW Influenza A+B Card	1	*Positive	1	100.00%
Quidel QuickVue	8	*Positive	8	100.00%
Henry Schein OneStep + Influenza A/B	8	*Positive	8	100.00%
VAD-2		* Denotes Correct Response		
AG: All Methods	17			
		*Positive	17	100.00%
Alere BinaxNOW Influenza A+B Card	1	*Positive	1	100.00%
Quidel QuickVue	8	*Positive	8	100.00%
Henry Schein OneStep + Influenza A/B	8	*Positive	8	100.00%
		RSV - Aq Detect		
VAD-1		* Denotes Correct Response		
AG: All Methods	278			
		*Positive	274	98.56%
		Negative	4	1.44%
BD Directigen RSV	1	*Positive	1	100.00%



BD EZ RSV	1	*Positive	1	100.00%
Alere BinaxNOW RSV	154	*Positive	154	100.00%
Quidel QuickVue RSV	27	*Positive	27	100.00%
Wampole Clearview RSV	3	*Positive	2	66.67%
SAS RSV Alert	1	*Positive	1	100.00%
Quidel QuickVue RSV 10	1	*Positive	1	100.00%
Quidel Sofia/Sofia 2 RSV	36	*Positive	34	94.44%
BD Veritor RSV	36	*Positive	35	97.22%
Cepheid GeneXpert - Xpert RSV	14	*Positive	14	100.00%
Alere ID NOW RSV	4	*Positive	4	100.00%
Wampole Clearview RSV	3	Negative	1	33.33%
Quidel Sofia/Sofia 2 RSV	36	Negative	2	5.56%
BD Veritor RSV	36	Negative	1	2.78%
VAD-2		* Denotes Correct Response		
AG: All Methods	279			
		*Negative	265	94.98%
		Positive	14	5.02%
BD Directigen RSV	1	*Negative	1	100.00%
BD EZ RSV	1	*Negative	1	100.00%
Alere BinaxNOW RSV	155	*Negative	147	94.84%
Quidel QuickVue RSV	27	*Negative	27	100.00%
Wampole Clearview RSV	3	*Negative	2	66.67%
SAS RSV Alert	1	*Negative	1	100.00%
Quidel QuickVue RSV 10	1	*Negative	1	100.00%
Quidel Sofia/Sofia 2 RSV	36	*Negative	33	91.67%
BD Veritor RSV	36	*Negative	34	94.44%
Cepheid GeneXpert - Xpert RSV	14	*Negative	14	100.00%
Alere ID NOW RSV	4	*Negative	4	100.00%
Alere BinaxNOW RSV	155	Positive	8	5.16%
Wampole Clearview RSV	3	Positive	1	33.33%
Quidel Sofia/Sofia 2 RSV	36	Positive	3	8.33%
BD Veritor RSV	36	Positive	2	5.56%
VAD-3		* Denotes Correct Response		
AG: All Methods	23			
		*Positive	22	95.65%
		Negative	1	4.35%
BD EZ RSV	1	*Positive	1	100.00%
Alere BinaxNOW RSV	9	*Positive	9	100.00%
Quidel QuickVue RSV	4	*Positive	4	100.00%
Wampole Clearview RSV	2	*Positive	1	50.00%
SAS RSV Alert	1	*Positive	1	100.00%
Quidel QuickVue RSV 10	1	*Positive	1	100.00%
Quidel Sofia/Sofia 2 RSV	3	*Positive	3	100.00%
Cepheid GeneXpert - Xpert RSV	2	*Positive	2	100.00%
Wampole Clearview RSV	2	Negative	1	50.00%
VAD-4		* Denotes Correct Response		
AG: All Methods	23			
		*Negative	22	95.65%
		Positive	1	4.35%
BD EZ RSV	1	*Negative	1	100.00%
Alere BinaxNOW RSV	9	*Negative	9	100.00%
Quidel QuickVue RSV	4	*Negative	4	100.00%
Wampole Clearview RSV	2	*Negative	1	50.00%
SAS RSV Alert	1	*Negative	1	100.00%
Quidel QuickVue RSV 10	1	*Negative	1	100.00%
Quidel Sofia/Sofia 2 RSV	3	*Negative	3	100.00%
Cepheid GeneXpert - Xpert RSV	2	*Negative	2	100.00%
Wampole Clearview RSV	2	Positive	1	50.00%
VAD-5		* Denotes Correct Response		
AG: All Methods	23			
		*Negative	22	95.65%
		Positive	1	4.35%
BD EZ RSV	1	*Negative	1	100.00%
Alere BinaxNOW RSV	9	*Negative	8	88.89%
Quidel QuickVue RSV	4	*Negative	4	100.00%
Wampole Clearview RSV	2	*Negative	2	100.00%
SAS RSV Alert	1	*Negative	1	100.00%
Quidel QuickVue RSV 10	1	*Negative	1	100.00%
Quidel Sofia/Sofia 2 RSV	3	*Negative	3	100.00%



Cepheid GeneXpert - Xpert RSV	2	*Negative	2	100.00%
Alere BinaxNOW RSV	9	Positive	1	11.11%
		Adenovirus		
VAD-1		* Denotes Correct Response		
AG: All Methods	9			
		*Negative	8	88.89%
		Positive	1	11.11%
SAS Adeno Test	9	*Negative	8	88.89%
SAS Adeno Test		Positive	1	11.11%
VAD-2		* Denotes Correct Response		
AG: All Methods	9			
		*Positive	9	100.00%
SAS Adeno Test	9	*Positive	9	100.00%
VAD-3		* Denotes Correct Response		
AG: All Methods	9			
		*Negative	8	88.89%
		Positive	1	11.11%
SAS Adeno Test	9	*Negative	8	88.89%
SAS Adeno Test		Positive	1	11.11%
VAD-4		* Denotes Correct Response		
AG: All Methods	9			
		*Positive	9	100.00%
SAS Adeno Test	9	*Positive	9	100.00%
VAD-5		* Denotes Correct Response		
AG: All Methods	9			
		*Negative	9	100.00%
SAS Adeno Test	9	*Negative	9	100.00%
932 - Hepatitis C Antibody (Waived);767 - Hepatitis Markers	Labs	Response	Number	Percent
		HCV (Waived)		
VM-1		* Denotes Correct Response		
AG: All Methods	4			
		*Negative	4	100.00%
OraSure QraQuick HCV Rapid Antibody Test	4	*Negative	4	100.00%
VM-2		* Denotes Correct Response		
AG: All Methods	4			
		*Positive	4	100.00%
OraSure QraQuick HCV Rapid Antibody Test	4	*Positive	4	100.00%
		HBsAg		
VM-1		* Denotes Correct Response		
AG: All Methods	1			
		*Positive	1	100.00%
Roche Cobas e Series	1	*Positive	1	100.00%
VM-2		* Denotes Correct Response		
AG: All Methods	1			
		*Positive	1	100.00%
Roche Cobas e Series	1	*Positive	1	100.00%
767 - Hepatitis Markers	Labs	Response	Number	Percent
		HBsAg		
VM-1		* Denotes Correct Response		
AG: All Methods	1			
		*Positive	1	100.00%
Roche Cobas e Series	1	*Positive	1	100.00%
VM-2		* Denotes Correct Response		
AG: All Methods	1			
		*Positive	1	100.00%
Roche Cobas e Series	1	*Positive	1	100.00%
VM-3		* Denotes Correct Response		
AG: All Methods	2			
		*Negative	2	100.00%
Roche Cobas e Series	1	*Negative	1	100.00%
Abbott Architect ci4100	1	*Negative	1	100.00%
VM-4		* Denotes Correct Response		
AG: All Methods	2			
		*Negative	2	100.00%
Roche Cobas e Series	1	*Negative	1	100.00%
Abbott Architect ci4100	1	*Negative	1	100.00%
VM-5		* Denotes Correct Response		



AG: All Methods	2			
		*Negative	2	100.00%
Roche Cobas e Series	1	*Negative	1	100.00%
Abbott Architect ci4100	1	*Negative	1	100.00%
		HCV		
VM-1		* Denotes Correct Response		
AG: All Methods	2			
		*Negative	2	100.00%
Abbott Architect ci4100	2	*Negative	2	100.00%
VM-2		* Denotes Correct Response		
AG: All Methods	2			
		*Positive	2	100.00%
Abbott Architect ci4100	2	*Positive	2	100.00%
VM-3		* Denotes Correct Response		
AG: All Methods	2			
		*Negative	2	100.00%
Abbott Architect ci4100	2	*Negative	2	100.00%
VM-4		* Denotes Correct Response		
AG: All Methods	2			
		*Negative	2	100.00%
Abbott Architect ci4100	2	*Negative	2	100.00%
VM-5		* Denotes Correct Response		
AG: All Methods	2			
		*Positive	2	100.00%
Abbott Architect ci4100	2	*Positive	2	100.00%
		Anti-HBc IgM		
VM-1		* Denotes Correct Response		
AG: All Methods	1			
		*Negative	1	100.00%
Abbott Architect ci4100	1	*Negative	1	100.00%
VM-2		* Denotes Correct Response		
AG: All Methods	1			
		*Negative	1	100.00%
Abbott Architect ci4100	1	*Negative	1	100.00%
VM-3		* Denotes Correct Response		
AG: All Methods	1			
		*Negative	1	100.00%
Abbott Architect ci4100	1	*Negative	1	100.00%
VM-4		* Denotes Correct Response		
AG: All Methods	1			
		*Negative	1	100.00%
Abbott Architect ci4100	1	*Negative	1	100.00%
VM-5		* Denotes Correct Response		
AG: All Methods	1			
		*Negative	1	100.00%
Abbott Architect ci4100	1	*Negative	1	100.00%
		HAV IgM		
VM-1		* Denotes Correct Response		
AG: All Methods	1			
		*Negative	1	100.00%
Abbott Architect ci4100	1	*Negative	1	100.00%
VM-2		* Denotes Correct Response		
AG: All Methods	1			
		*Negative	1	100.00%
Abbott Architect ci4100	1	*Negative	1	100.00%
VM-3		* Denotes Correct Response		
AG: All Methods	1			
		*Negative	1	100.00%
Abbott Architect ci4100	1	*Negative	1	100.00%
VM-4		* Denotes Correct Response		
AG: All Methods	1			
		*Negative	1	100.00%
Abbott Architect ci4100	1	*Negative	1	100.00%
VM-5		* Denotes Correct Response		
AG: All Methods	1			
		*Negative	1	100.00%
Abbott Architect ci4100	1	*Negative	1	100.00%



		Anti-HBc Total (International)		
VM-1		* Denotes Correct Response		
AG: All Methods	3			
		*Positive	3	100.00%
Intl OCD Vitros ECi/ECiQ	3	*Positive	3	100.00%
VM-2		* Denotes Correct Response		
AG: All Methods	3			
		*Positive	3	100.00%
Intl OCD Vitros ECi/ECiQ	3	*Positive	3	100.00%
VM-3		* Denotes Correct Response		
AG: All Methods	3			
		*Negative	3	100.00%
Intl OCD Vitros ECi/ECiQ	3	*Negative	3	100.00%
VM-4		* Denotes Correct Response		
AG: All Methods	3			
		*Positive	3	100.00%
Intl OCD Vitros ECi/ECiQ	3	*Positive	3	100.00%
VM-5		* Denotes Correct Response		
AG: All Methods	3			
		*Positive	3	100.00%
Intl OCD Vitros ECi/ECiQ	3	*Positive	3	100.00%
		HBsAb (International)		
VM-1		* Denotes Correct Response		
AG: All Methods	2			
		*Negative	2	100.00%
Intl OCD Vitros ECi/ECiQ	2	*Negative	2	100.00%
VM-2		* Denotes Correct Response		
AG: All Methods	2			
		*Negative	2	100.00%
Intl OCD Vitros ECi/ECiQ	2	*Negative	2	100.00%
VM-3		* Denotes Correct Response		
AG: All Methods	2			
		*Positive	2	100.00%
Intl OCD Vitros ECi/ECiQ	2	*Positive	2	100.00%
VM-4		* Denotes Correct Response		
AG: All Methods	2			
		*Positive	2	100.00%
Intl OCD Vitros ECi/ECiQ	2	*Positive	2	100.00%
VM-5		* Denotes Correct Response		
AG: All Methods	2			
		*Negative	2	100.00%
Intl OCD Vitros ECi/ECiQ	2	*Negative	2	100.00%
		HBsAg (International)		
VM-1		* Denotes Correct Response		
AG: All Methods	3			
		*Positive	3	100.00%
Intl OCD Vitros ECi/ECiQ	3	*Positive	3	100.00%
VM-2		* Denotes Correct Response		
AG: All Methods	3			
		*Positive	3	100.00%
Intl OCD Vitros ECi/ECiQ	3	*Positive	3	100.00%
VM-3		* Denotes Correct Response		
AG: All Methods	3			
		*Negative	3	100.00%
Intl OCD Vitros ECi/ECiQ	3	*Negative	3	100.00%
VM-4		* Denotes Correct Response		
AG: All Methods	3			
		*Negative	3	100.00%
Intl OCD Vitros ECi/ECiQ	3	*Negative	3	100.00%
VM-5		* Denotes Correct Response		
AG: All Methods	3			
		*Negative	3	100.00%
Intl OCD Vitros ECi/ECiQ	3	*Negative	3	100.00%
		HCV (International)		
VM-1		* Denotes Correct Response		
AG: All Methods	4			
		*Negative	4	100.00%



Intl OCD Vitros ECi/ECiQ	4	*Negative	4	100.00%
VM-2		* Denotes Correct Response		
AG: All Methods	4			
		*Positive	4	100.00%
Intl OCD Vitros ECi/ECiQ	4	*Positive	4	100.00%
VM-3		* Denotes Correct Response		
AG: All Methods	4			
		*Negative	4	100.00%
Intl OCD Vitros ECi/ECiQ	4	*Negative	4	100.00%
VM-4		* Denotes Correct Response		
AG: All Methods	5			
		*Negative	4	80.00%
		Positive	1	20.00%
Intl OCD Vitros ECi/ECiQ	5	*Negative	4	80.00%
Intl OCD Vitros ECi/ECiQ		Positive	1	20.00%
VM-5		* Denotes Correct Response		
AG: All Methods	4			
		*Positive	4	100.00%
Intl OCD Vitros ECi/ECiQ	4	*Positive	4	100.00%
771 - Vaginal Pathogens I	Labs	Response	Number	Percent
		Candida sp. - Non Culture		
VP-1		* Denotes Correct Response		
AG: All Methods	48			
		*Negative	47	97.92%
		Positive	1	2.08%
Affirm VP	47	*Negative	47	100.00%
BD Max	1	Positive	1	100.00%
VP-2		* Denotes Correct Response		
AG: All Methods	48			
		*Positive	47	97.92%
		Negative	1	2.08%
Affirm VP	47	*Positive	46	97.87%
BD Max	1	*Positive	1	100.00%
Affirm VP	47	Negative	1	2.13%
VP-3		* Denotes Correct Response		
AG: All Methods	48			
		*Negative	48	100.00%
Affirm VP	47	*Negative	47	100.00%
BD Max	1	*Negative	1	100.00%
VP-4		* Denotes Correct Response		
AG: All Methods	48			
		*Negative	48	100.00%
Affirm VP	47	*Negative	47	100.00%
BD Max	1	*Negative	1	100.00%
VP-5		* Denotes Correct Response		
AG: All Methods	48			
		*Positive	46	95.83%
		Negative	2	4.17%
Affirm VP	47	*Positive	45	95.74%
BD Max	1	*Positive	1	100.00%
Affirm VP	47	Negative	2	4.26%
		G. vaginalis - Non Culture		
VP-1		* Denotes Correct Response		
AG: All Methods	48			
		*Negative	48	100.00%
Affirm VP	47	*Negative	47	100.00%
BD Max	1	*Negative	1	100.00%
VP-2		* Denotes Correct Response		
AG: All Methods	48			
		*Negative	47	97.92%
		Positive	1	2.08%
Affirm VP	47	*Negative	46	97.87%
BD Max	1	*Negative	1	100.00%
Affirm VP	47	Positive	1	2.13%
VP-3		* Denotes Correct Response		
AG: All Methods	48			
		*Positive	48	100.00%
Affirm VP	47	*Positive	47	100.00%



BD Max	1	*Positive	1	100.00%
VP-4		* Denotes Correct Response		
AG: All Methods	48			
		*Positive	48	100.00%
Affirm VP	47	*Positive	47	100.00%
BD Max	1	*Positive	1	100.00%
VP-5		* Denotes Correct Response		
AG: All Methods	48			
		*Negative	46	95.83%
		Positive	2	4.17%
Affirm VP	47	*Negative	45	95.74%
BD Max	1	*Negative	1	100.00%
Affirm VP	47	Positive	2	4.26%
		T. vaginalis - Non Culture Tech		
VP-1		* Denotes Correct Response		
AG: All Methods	48			
		*Positive	48	100.00%
Affirm VP	47	*Positive	47	100.00%
BD Max	1	*Positive	1	100.00%
VP-2		* Denotes Correct Response		
AG: All Methods	48			
		*Negative	48	100.00%
Affirm VP	47	*Negative	47	100.00%
BD Max	1	*Negative	1	100.00%
VP-3		* Denotes Correct Response		
AG: All Methods	48			
		*Negative	48	100.00%
Affirm VP	47	*Negative	47	100.00%
BD Max	1	*Negative	1	100.00%
VP-4		* Denotes Correct Response		
AG: All Methods	48			
		*Negative	48	100.00%
Affirm VP	47	*Negative	47	100.00%
BD Max	1	*Negative	1	100.00%
VP-5		* Denotes Correct Response		
AG: All Methods	48			
		*Positive	47	97.92%
		Negative	1	2.08%
Affirm VP	47	*Positive	46	97.87%
BD Max	1	*Positive	1	100.00%
Affirm VP	47	Negative	1	2.13%