



SG: Select Roche Instruments	6	286.17	16.29	200.00-372.00	5.70	6	69.32	4.24	48.00-90.00	6.09	6	194.95	10.63	137.00-254.00	5.44	6	224.45	15.11	157.00-291.00	6.74	6	102.63	6.85	72.00-134.00	6.60
SG: Select OCD Vitros Instruments	9	204.89	8.42	144.00-267.00	4.10	9	67.00	2.60	47.00-87.00	3.88	9	161.00	6.86	113.00-209.00	4.29	9	178.78	6.65	125.00-233.00	3.69	9	95.44	3.75	67.00-124.00	3.89
SG: Abaxis Piccolo Instruments	39	217.92	27.68	153.00-283.00	12.71	39	62.69	4.03	44.00-82.00	6.35															
SG: Select Abbott Architect Instruments	9	293.11	12.42	205.00-381.00	4.23	9	70.78	2.95	50.00-92.00	4.08	9	198.22	8.07	139.00-257.00	4.09	9	229.11	10.41	160.00-298.00	4.54	9	105.00	5.45	74.00-137.00	5.24
SG: International OCD Vitros Instruments	12	204.83	8.58	144.00-267.00	4.20	12	66.42	3.17	46.00-86.00	4.85	12	160.00	7.35	112.00-208.00	4.56	12	178.08	6.72	125.00-231.00	3.76	12	94.92	3.40	67.00-124.00	3.58
SG: International Abaxis Piccolo Instruments	40	217.63	27.39	153.00-283.00	12.57	40	62.65	3.99	44.00-82.00	6.35															
SG: Select Gen Chem Instruments	58	272.50	23.23	191.00-355.00	8.50	58	65.87	5.98	46.00-86.00	9.09	58	185.71	16.57	130.00-242.00	8.92	58	212.54	19.65	149.00-277.00	9.20	58	97.56	8.99	69.00-127.00	9.18
MG: Method Group	118	252.00	38.70	176.00-328.00	15.36	118	65.00	5.60	46.00-85.00	8.62	83	184.00	18.40	129.00-239.00	10.00	83	210.00	23.10	147.00-273.00	11.00	83	98.00	8.90	69.00-127.00	9.08
MG: IFCC	9	301.00	23.00	211.00-391.00	7.64	9	72.00	6.30	50.00-94.00	8.75	9	205.00	15.80	144.00-267.00	7.71	9	238.00	18.60	167.00-309.00	7.82	9	107.00	7.70	75.00-139.00	7.20
ALT/SGPT	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%
PG: Abaxis Piccolo (Non-Waived)	6	195.00	8.48	156.00-234.00	4.36	6	44.33	4.80	35.00-53.00	10.91	6	129.33	5.13	103.00-155.00	3.95	6	150.00	6.00	120.00-180.00	4.00	6	67.67	6.50	54.00-82.00	9.56
PG: Abaxis Piccolo (Waived)	35	188.89	25.19	151.00-227.00	13.33	35	44.89	2.49	36.00-54.00	5.56															
PG: Abbott Architect	6	222.00	1.10	178.00-266.00	0.50	6	44.00	1.79	35.00-53.00	4.09	6	147.50	1.52	118.00-178.00	1.01	6	173.17	3.37	138.00-208.00	1.97	6	71.00	2.45	57.00-85.00	3.38
PG: Alfa Wasserman Alera/Axcel	23	179.91	9.81	144.00-216.00	5.44	23	34.17	3.61	27.00-41.00	10.59	23	120.39	3.47	96.00-144.00	2.92	23	140.78	4.21	113.00-169.00	2.98	23	56.83	2.13	46.00-68.00	3.68
PG: Siemens Dimension/AR/ES/RxL/Xpand/EXL / IFCC	9	225.89	4.34	181.00-271.00	1.90	9	47.56	1.51	38.00-58.00	3.13	9	152.22	2.95	122.00-182.00	1.91	9	177.11	6.35	142.00-212.00	3.62	9	75.00	2.24	60.00-90.00	2.93
PG: Beckman AU	19	194.63	6.05	156.00-234.00	3.08	19	38.42	1.22	30.00-46.00	3.16	19	129.26	3.44	103.00-155.00	2.64	19	151.53	4.56	122.00-182.00	3.03	19	61.89	1.70	50.00-74.00	2.74
PG: Roche Cobas Integra	5	202.20	16.36	162.00-242.00	8.12	5	40.40	3.13	32.00-48.00	7.75	5	134.80	10.40	108.00-162.00	7.70	5	157.00	13.34	126.00-188.00	8.47	5	64.80	5.26	52.00-78.00	8.15
PG: OCD Vitros 250,350,400,500,750,950	7	219.43	5.53	175.00-263.00	2.51	7	62.00	4.24	50.00-74.00	6.77	7	153.71	3.35	123.00-185.00	2.21	7	174.43	3.73	139.00-209.00	2.13	7	84.43	2.37	67.00-101.00	2.86
SG: Select Beckman Coulter Instruments (Chem)	21	196.48	8.18	157.00-235.00	4.18	21	39.52	3.67	32.00-48.00	9.25	21	130.81	5.88	105.00-157.00	4.50	21	153.19	6.82	122.00-184.00	4.44	21	63.05	3.98	50.00-76.00	6.35
SG: Select Siemens Instruments	5	233.25	10.87	186.00-280.00	4.68	5	48.00	2.45	38.00-58.00	5.00	5	157.50	6.25	126.00-190.00	3.92	5	182.75	7.80	146.00-220.00	4.26	5	76.25	2.99	61.00-91.00	3.95
SG: Select Roche Instruments	6	203.33	14.90	162.00-244.00	7.34	6	40.50	2.81	33.00-49.00	6.83	6	135.67	9.54	109.00-163.00	6.99	6	158.17	12.27	126.00-190.00	7.78	6	65.17	4.79	52.00-78.00	7.38
SG: Select OCD Vitros Instruments	8	219.38	5.13	175.00-263.00	2.33	8	61.75	3.99	50.00-74.00	6.45	8	153.88	3.14	123.00-185.00	2.01	8	174.38	3.46	139.00-209.00	2.01	8	84.63	2.26	68.00-102.00	2.71
SG: Abaxis Piccolo Instruments	41	189.78	23.52	152.00-228.00	12.37	41	44.80	2.87	36.00-54.00	6.44	6	129.33	5.13	103.00-155.00	3.95	6	150.00	6.00	120.00-180.00	4.00	6	67.67	6.50	54.00-82.00	9.56
SG: Select Abbott Architect Instruments	9	221.67	1.66	178.00-266.00	0.77	9	43.78	1.56	35.00-53.00	3.64	9	147.67	1.58	118.00-178.00	1.08	9	173.00	2.78	138.00-208.00	1.62	9	71.00	2.12	57.00-85.00	2.96
SG: International OCD Vitros Instruments	10	218.90	4.72	175.00-263.00	2.15	10	61.20	4.16	49.00-73.00	6.89	10	152.80	3.65	122.00-184.00	2.35	10	174.30	3.20	139.00-209.00	1.84	10	83.90	3.14	67.00-101.00	3.69
SG: International Abaxis Piccolo Instruments	42	189.88	23.24	152.00-228.00	12.21	42	44.79	2.83	36.00-54.00	6.22	7	129.43	4.68	103.00-155.00	3.64	7	150.29	5.53	120.00-180.00	3.67	7	67.29	6.02	54.00-80.00	8.96
SG: Select Gen Chem Instruments	57	194.94	17.22	156.00-234.00	8.82	57	38.40	5.07	30.00-46.00	13.42	57	130.31	11.13	104.00-156.00	8.54	57	152.23	12.83	122.00-182.00	8.42	57	62.38	6.25	50.00-74.00	10.00
MG: IFCC	12	225.00	10.40	180.00-270.00	4.62	12	47.00	3.20	38.00-56.00	6.81	12	151.00	7.50	121.00-181.00	4.97	12	176.00	9.80	141.00-211.00	5.57	12	74.00	3.90	59.00-89.00	5.27
MG: IFCC	12	225.00	10.40	180.00-270.00	4.62	12	47.00	3.20	38.00-56.00	6.81	12	151.00	7.50	121.00-181.00	4.97	12	176.00	9.80	141.00-211.00	5.57	12	74.00	3.90	59.00-89.00	5.27
AG: All Group	130	199.03	21.96	159.00-239.00	11.06	122	42.32	5.30	34.00-50.00	12.62	95	136.17	13.76	109.00-163.00	10.15	95	158.52	15.48	127.00-191.00	9.75	95	66.80	8.99	54.00-80.00	13.43
AST/SGOT	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%
PG: Abaxis Piccolo (Non-Waived)	6	243.67	8.73	195.00-293.00	3.57	6	124.00	6.54	99.00-149.00	5.24	6	190.83	7.41	153.00-229.00	3.87	6	208.83	8.21	167.00-251.00	3.92	6	141.17	6.85	113.00-169.00	4.89
PG: Abaxis Piccolo (Waived)	35	239.46	21.05	191.00-287.00	8.79	35	123.74	3.94	99.00-149.00	3.15															
PG: Abbott Architect	6	248.83	2.99	199.00-299.00	1.20	6	124.17	0.75	99.00-149.00	0.65	6	194.67	2.88	156.00-234.00	1.49	6	214.17	2.14	171.00-257.00	0.98	6	141.17	1.47	113.00-169.00	1.06
PG: Alfa Wasserman Alera/Axcel	23	225.26	6.29	180.00-270.00	2.80	23	113.43	6.42	90.00-136.00	5.66	23	178.57	4.63	143.00-215.00	2.57	23	194.52	4.93	156.00-234.00	2.51	23	129.43	3.16	103.00-155.00	2.48
PG: Siemens Dimension/AR/ES/RxL/Xpand/EXL	12	265.50	8.85	213.00-319.00	3.31	12	126.75	3.67	102.00-152.00	2.91	12	207.08	7.57	166.00-248.00	3.67	12	227.00	7.91	182.00-272.00	3.48	12	147.83	4.93	118.00-178.00	3.31
PG: Beckman AU	19	222.63	6.19	178.00-268.00	2.78	19	111.42	3.82	89.00-133.00	3.42	19	174.79	5.03	140.00-210.00	2.86	19	190.37	5.17	152.00-228.00	2.74	19	127.21	5.54	102.00-152.00	4.33
PG: Roche Cobas Integra	5	246.00	24.62	197.00-295.00	10.00	5	123.80	12.40	99.00-149.00	10.00	5	193.00	19.74	154.00-232.00	10.21	5	210.80	21.34	169.00-253.00	10.09	5	140.40	14.84	112.00-168.00	10.57
PG: OCD Vitros 250,350,400,500,750,950	7	298.29	8.86	238.00-358.00	2.99	7	127.43	2.44	102.00-152.00	1.89	7	218.71	5.85	175.00-263.00	2.69	7	250.86	9.51	201.00-301.00	3.78	7	149.14	3.02	119.00-179.00	2.01
SG: Select Beckman Coulter Instruments (Chem)	21	224.67	8.73	180.00-270.00	3.87	21	112.71	5.46	90.00-136.00	4.87	21	176.48	7.17	141.00-211.00	4.09	21	192.24	7.70	154.00-230.00	4.01	21	128.57	6.79	103.00-155.00	5.27
SG: Select Siemens Instruments	13	266.54	9.26	214.00-320.00	3.48	13	127.85	5.29	102.00-154.00	4.14	13	208.08	8.09	166.00-250.00	3.89	13	227.92	8.27	182.00-274.00	3.64	13	149.23	6.91	119.00-179.00	4.63
SG: Select Roche Instruments	6	244.50	22.32	196.00-294.00	9.10	6	123.17	11.20	98.00-148.00	9.11	6	192.00	17.82	154.00-230.00	9.27	6	210.00	19.18	168.00-252.00	9.14	6	139.67	13.40	112.00-168.00	9.57
SG: Select OCD Vitros Instruments	9	296.11	11.10	237.00-355.00	3.75	9	127.11	2.67	102.00-152.00	2.13	9	217.78	6.30	174.00-262.00	2.89	9	249.22	9.82	199.00-299.00	3.94	9	148.33	3.81	118.00-178.00	2.57
SG: Abaxis Piccolo Instruments	41	240.07	19.70	192.00-288.00	8.21	41	123.78	4.31	99.00-149.00	3.47	6	190.83	7.41	153.00-229.00	3.87	6	208.83	8.21	167.00-251.00	3.92	6	141.17	6.85	113.00-169.00	4.89
SG: Select Abbott Architect Instruments	9	248.00	2.96	198.00-298.00	1.21	9	123.56	1.67	99.00-149.00	1.37	9	194.33	2.74	155.00-233.00	1.39	9	213.44	2.65	170.00-256.00	1.27	9	140.56	1.94	113.00-169.00	1.35
SG: International OCD Vitros Instruments	12	297.92	10.31	238.00-358.00	3.46	12	127.67																		



Granulocyte Abs #	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%
PG: Horiba Micros/45/60	46	0.90	0.06	0.70-1.10	6.67	46	12.91	0.34	11.90-13.90	2.64	46	5.22	0.17	4.70-5.70	3.27	46	20.69	0.61	18.90-22.50	2.95	45	5.80	0.17	5.30-6.30	2.93
Lymphocyte Abs #	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%
PG: Horiba Micros/45/60	46	0.71	0.09	0.50-0.90	11.43	46	2.73	0.23	2.00-3.40	8.52	46	1.78	0.14	1.40-2.20	7.78	46	2.81	0.23	2.10-3.50	8.21	45	1.62	0.16	1.10-2.10	10.00
Monocyte/Mixed Abs #	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%
PG: Horiba Micros/45/60	45	0.42	0.06	0.20-0.60	15.00	46	1.26	0.12	0.90-1.70	9.23	46	0.66	0.08	0.50-0.90	11.43	46	1.11	0.09	0.80-1.40	8.18	46	0.77	0.10	0.50-1.10	12.50
MCH	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%
PG: Horiba Micros/45/60	46	25.30	0.48	23.90-26.70	1.90	46	29.50	0.39	28.30-30.70	1.32	46	29.78	0.40	28.60-31.00	1.34	46	29.67	0.46	28.30-31.10	1.55	46	29.70	0.42	28.40-31.00	1.41
MCHC	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%
PG: Horiba Micros/45/60	46	36.87	0.88	34.30-39.50	2.38	46	36.20	0.79	33.80-38.60	2.18	46	36.77	0.61	35.00-38.60	1.66	46	36.14	0.91	33.40-38.80	2.52	46	37.03	0.68	35.00-39.00	1.84
Hematocrit (Waived)	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%
PG: Microhematocrit, Spun	14	16.64	1.74	11.40-21.80	10.48	8	38.25	2.12	31.90-44.70	5.54															
AG: All Group	20	16.95	1.78	11.50-22.30	10.59	8	38.25	2.12	31.90-44.70	5.54															
Hemoglobin (Waived)	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%
PG: HemoCue	40	6.04	0.26	5.20-6.80	4.33																				
PG: HemoCue 201	99	6.03	0.37	4.90-7.10	6.17																				
PG: Alere HemoPoint H2	42	6.17	0.17	5.70-6.70	2.74																				
SG: Select Waived Hemoglobin Instruments	184	6.05	0.39	4.80-7.20	6.50																				
AG: All Group	189	6.13	0.91	3.40-8.80	14.92	7	15.53	1.30	11.60-19.40	8.39															
608 - Hemocue Glucose/Hemoglobin	HQ-3																								
Glucose Whole Blood (Hemocue)	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%
PG: HemoCue	35	48.27	9.18	20.00-76.00	19.17																				
PG: HemoCue 201	140	49.01	9.59	20.00-78.00	19.59																				
Hemoglobin (Hemocue)	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%
PG: HemoCue	21	4.90	0.43	3.60-6.20	8.78																				
PG: HemoCue 201	80	5.03	0.26	4.20-5.80	5.20																				
606 - Hemocue WBC	HQW-11					HQW-12					HQW-13					HQW-14					HQW-15				
White Blood Cell Count	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%
PG: Hemocue WBC System	11	21.39	0.86	18.20-24.60	4.07	11	2.83	0.15	2.40-3.20	5.36	11	7.77	0.38	6.60-9.00	4.87	11	21.75	1.12	18.40-25.00	5.16	11	22.16	2.71	18.90-25.50	12.21
700 - i-Stat Chemistry;721 - i-Stat Chemistry & Blood Gases;722 - i-STAT Chemistry (1 specimen);723 - i-Stat Blood Gases	IST-11					IST-12					IST-13					IST-14					IST-15				
Blood Gas - pCO2	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%
PG: Abbott i-Stat	9	19.14	1.30	14.00-24.00	6.84	9	27.92	0.94	23.00-33.00	3.21	9	24.07	1.39	19.00-29.00	5.83	9	56.40	1.55	51.00-61.00	2.68	9	19.63	0.50	15.00-25.00	2.50
SG: International Abbott i-Stat Instruments	13	19.41	1.14	14.00-24.00	5.79	13	27.87	1.22	23.00-33.00	4.29	13	24.12	1.17	19.00-29.00	5.00	13	56.51	1.42	52.00-62.00	2.46	13	19.52	0.51	15.00-25.00	2.50
Blood Gas - pH	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%
PG: Abbott i-Stat	9	7.48	0.00	7.44-7.52	0.05	9	7.50	0.00	7.46-7.54	0.05	9	7.49	0.01	7.45-7.53	0.09	9	7.20	0.01	7.16-7.24	0.17	9	7.48	0.01	7.44-7.52	0.07
SG: International Abbott i-Stat Instruments	13	7.48	0.00	7.44-7.52	0.05	13	7.50	0.00	7.46-7.54	0.05	13	7.49	0.01	7.45-7.53	0.08	13	7.19	0.01	7.15-7.23	0.17	13	7.48	0.01	7.44-7.52	0.07
Blood Gas - pO2	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%
PG: Abbott i-Stat	8	161.50	6.80	142.00-182.00	4.20	9	120.56	7.78	98.00-144.00	6.45	9	138.44	4.77	124.00-152.00	3.48	8	84.75	6.78	65.00-105.00	8.00	9	162.56	7.05	142.00-184.00	4.36
SG: International Abbott i-Stat Instruments	12	165.50	8.26	141.00-191.00	5.00	13	119.23	6.78	99.00-139.00	5.71	13	140.15	5.11	125.00-155.00	3.64	12	83.00	6.27	64.00-102.00	7.59	13	164.00	6.42	145.00-183.00	3.90
Chloride (CL-)	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%	Labs	Mean	SD	Range	CV.%
PG: Abbott i-Stat / Thermal	6	103.50	0.55	98.30-108.70	0.53	6	116.00	0.63	110.00-122.00	0.54	6	84.67	0.82	81.00-89.00	0.97	6	75.17	0.75	71.00-79.00	1.07	6	103.17	0.41	98.00-108.00	0.40
SG: International Abbott i-Stat Instruments / Thermal	35	102.97	0.75	97.80-108.20	0.73	10	115.80	0.63	110.00-122.00	0.54	10														



PG: Sysmex poch-100i	21	6.25	0.11	5.80-6.60	1.77	21	17.75	1.31	16.50-18.90	7.40	20	12.12	0.21	11.30-12.90	1.74	21	18.09	0.27	16.80-19.40	1.49	21	12.12	0.19	11.30-12.90	1.57
SG: Select Sysmex Instruments (3-part Diff)	121	6.03	0.17	5.60-6.40	2.83	124	17.88	0.71	16.60-19.20	3.97	122	12.03	0.19	11.20-12.80	1.58	124	17.97	0.63	16.70-19.30	3.50	123	12.10	0.43	11.30-12.90	3.55
PG: Sysmex XP-300	85	5.98	0.14	5.60-6.40	2.33	88	17.91	0.56	16.60-19.20	3.13	87	12.01	0.19	11.20-12.80	1.58	88	17.93	0.73	16.60-19.20	4.08	87	12.09	0.50	11.30-12.90	4.13
SG: International Sysmex Instruments	25	6.23	0.12	5.80-6.60	1.94	25	17.80	1.20	16.60-19.00	6.74	24	12.10	0.20	11.30-12.90	1.65	25	18.08	0.26	16.80-19.40	1.44	25	12.10	0.19	11.30-12.90	1.57
Lymphocyte	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%
PG: Sysmex K-800,1000,4500,KX-21	14	11.31	1.16	7.80-14.80	10.27	15	63.61	0.52	62.00-65.20	0.82	15	31.65	0.79	29.20-34.00	2.50	15	63.63	0.52	62.00-65.20	0.82	15	31.79	0.61	30.00-33.60	1.92
PG: Sysmex poch-100i	19	10.55	1.03	7.40-13.60	9.81	18	63.49	0.63	61.60-65.40	0.99	19	30.03	1.37	25.90-34.10	4.57	20	63.51	0.73	61.30-65.70	1.15	20	30.33	1.10	27.00-33.60	3.63
SG: Select Sysmex Instruments (3-part Diff)	115	11.58	1.32	7.60-15.60	11.38	120	63.56	0.68	61.60-65.60	1.07	119	31.41	1.52	26.80-36.00	4.84	121	63.62	0.88	61.00-66.20	1.38	120	31.46	1.44	27.20-35.80	4.57
PG: Sysmex XP-300	82	11.87	1.29	8.00-15.80	10.84	87	63.57	0.72	61.40-65.80	1.13	85	31.68	1.50	27.20-36.20	4.73	86	63.64	0.96	60.70-66.50	1.51	85	31.67	1.49	27.20-36.20	4.70
SG: International Sysmex Instruments	23	10.55	1.01	7.50-13.50	9.62	22	63.38	0.70	61.30-65.50	1.10	23	30.08	1.26	26.30-33.90	4.19	24	63.38	0.76	61.10-65.70	1.20	24	30.17	1.22	26.50-33.90	4.04
Hematocrit	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%
PG: Sysmex K-800,1000,4500,KX-21	15	17.10	0.34	16.10-18.10	1.99	15	47.40	0.62	44.60-50.20	1.31	15	32.02	0.47	30.10-33.90	1.47	15	47.65	0.88	44.70-50.50	1.85	15	31.95	0.54	30.00-33.80	1.69
PG: Sysmex poch-100i	21	18.30	0.43	17.20-19.40	2.35	21	49.70	3.81	46.70-52.70	7.67	20	34.15	0.75	32.10-36.10	2.20	21	50.42	1.08	47.40-53.40	2.14	21	34.19	0.70	32.10-36.30	2.05
SG: Select Sysmex Instruments (3-part Diff)	121	17.48	0.51	16.40-18.60	2.91	124	48.04	2.06	45.10-50.90	4.29	122	32.47	0.92	30.60-34.50	2.83	124	48.45	2.75	45.60-51.40	5.65	122	32.61	1.28	30.60-34.60	3.93
PG: Sysmex XP-300	85	17.35	0.32	16.30-18.30	1.85	88	47.76	1.33	44.90-50.70	2.78	87	32.16	0.48	30.30-34.10	1.52	88	48.12	3.02	45.20-51.00	6.28	86	32.35	1.19	30.40-34.20	3.68
SG: International Sysmex Instruments	25	18.29	0.40	17.20-19.40	2.19	25	49.82	3.50	46.80-52.80	7.03	24	34.14	0.73	32.10-36.10	2.14	25	50.42	1.07	47.40-53.40	2.12	25	34.18	0.68	32.10-36.30	1.99
Mono/Mixed	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%
PG: Sysmex K-800,1000,4500,KX-21	15	18.49	1.47	14.10-22.90	7.95	15	12.80	0.78	10.50-15.10	6.09	15	17.55	1.02	14.40-20.60	5.83	15	13.02	0.69	11.00-15.00	5.23	15	17.42	1.14	14.00-20.80	6.55
PG: Sysmex poch-100i	20	16.91	1.43	12.60-21.20	8.46	20	12.13	1.11	8.80-15.40	9.17	20	15.99	1.30	12.10-19.90	8.19	19	12.04	0.82	9.50-14.50	6.83	19	15.99	0.97	13.10-18.90	6.06
SG: Select Sysmex Instruments (3-part Diff)	116	18.16	1.66	13.30-23.10	9.07	122	12.52	1.19	8.90-16.10	9.52	121	17.17	1.68	12.20-22.20	9.77	119	12.68	0.86	10.10-15.30	6.77	116	17.07	1.15	13.70-20.60	6.73
PG: Sysmex XP-300	81	18.42	1.61	13.60-23.20	8.75	87	12.57	1.25	8.80-16.40	9.92	86	17.38	1.74	12.20-22.60	10.00	85	12.76	0.84	10.30-15.30	6.56	82	17.26	1.06	14.10-20.50	6.13
SG: International Sysmex Instruments	24	16.82	1.56	12.10-21.50	9.29	24	12.11	1.07	8.90-15.30	8.84	24	15.87	1.29	12.00-19.80	8.11	23	12.12	0.80	9.70-14.50	6.61	23	15.93	0.91	13.20-18.60	5.72
MCV	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%
PG: Sysmex K-800,1000,4500,KX-21	15	70.31	1.16	66.80-73.80	1.65	15	84.30	1.16	80.80-87.80	1.38	15	79.99	1.06	76.80-83.20	1.33	15	84.48	1.10	81.20-87.80	1.30	15	79.85	1.13	76.50-83.30	1.41
PG: Sysmex poch-100i	16	73.70	1.29	69.80-77.60	1.75	16	88.22	1.61	83.40-93.00	1.83	16	84.07	1.98	78.20-90.00	2.35	16	88.21	1.13	84.80-91.60	1.28	16	83.76	1.26	80.00-87.60	1.50
SG: Select Sysmex Instruments (3-part Diff)	117	71.45	1.44	67.10-75.70	2.02	118	85.38	1.60	80.60-90.20	1.87	118	80.93	1.72	75.70-86.10	2.13	117	85.36	1.55	80.80-90.10	1.81	118	80.91	1.54	76.30-85.50	1.90
PG: Sysmex XP-300	86	71.23	1.03	68.10-74.30	1.45	87	85.04	0.99	82.00-88.00	1.16	87	80.52	0.99	77.50-83.50	1.23	86	84.99	1.03	81.90-88.10	1.21	87	80.57	0.95	77.70-83.50	1.19
SG: International Sysmex Instruments	20	73.79	1.23	70.10-77.50	1.67	20	88.36	1.51	83.90-92.90	1.71	20	84.00	1.79	78.60-89.40	2.13	20	88.24	1.06	85.00-91.40	1.20	20	83.76	1.16	80.30-87.30	1.38
Platelet Count	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%
PG: Sysmex K-800,1000,4500,KX-21	15	62.87</																							



SG: Select Siemens (Bayer) DCA Instruments	15	68.85	4.10	57.00-81.00	5.94	15	71.57	2.90	63.00-81.00	4.03															
PG: Afinion ACR	10	65.85	3.36	56.00-76.00	5.15	10	64.78	2.39	58.00-72.00	3.69															
SG: Select Gen Chem Instruments	13	60.18	5.59	43.00-77.00	9.33	13	61.15	7.03	40.00-82.00	11.48															
AG: All Group	44	64.80	5.65	48.00-82.00	8.77	44	65.66	6.37	47.00-85.00	9.70															
728 - Urine Chemistry	URC-5					URC-6																			
Creatinine - Urine	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%
AG: All Group	5	185.87	1.40	181.70-190.10	0.75	5	62.73	0.83	60.20-65.20	1.32															
Total Protein - Urine	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%
AG: All Group	5	98.95	7.21	77.30-120.50	7.29	5	11.28	4.60	0.00-25.10	40.71															
601 - Waived/PPM Tests II;609 - Waived/PPM Tests I;664 - Whole Blood Glucose;669 - Whole Blood Glucose	WBG-11					WBG-12					WBG-13					WBG-14					WBG-15				
Glucose Wh Blood (Waived)	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%	Labs	Mean	SD	Range	CV,%
PG: Abbott Freestyle Lite	11	158.91	9.14	132.00-186.00	5.72																				
PG: Bayer Ascensia Contour	136	134.65	15.72	88.00-182.00	11.63	14	46.71	5.92	29.00-65.00	12.55															
PG: Bayer Ascensia Elite XL	8	137.50	18.57	82.00-194.00	13.48																				
PG: Lifescan One Touch Ultra	31	176.19	25.78	99.00-253.00	14.66																				
PG: McKesson/SunMark True Track	5	408.00	7.21	386.00-430.00	1.76																				
PG: Roche Accu-Chek Advantage/Comp/GTS	18	165.78	13.28	126.00-206.00	8.01																				
PG: Roche AccuChek Aviva	21	168.62	13.20	129.00-209.00	7.81																				
PG: Roche Accu-Chek II	6	162.50	9.83	134.00-192.00	6.01																				
PG: Roche Accu-Chek, others	16	169.25	20.42	108.00-230.00	12.07																				
SG: Select Lifescan One Touch Instruments	40	174.40	23.91	102.00-246.00	13.74																				
SG: Select Abbott Freestyle Instruments	14	156.43	9.86	126.00-186.00	6.35																				
PG: Roche AccuChek Inform	73	170.71	6.76	151.00-191.00	3.98	35	69.89	2.07	64.00-76.00	3.00															
PG: Consult Diag Quintet Blood Glucose	18	168.94	15.77	122.00-216.00	9.35																				
PG: Medisense Precision Xceed Pro	26	157.96	9.98	128.00-188.00	6.33																				
SG: Select Medisense Instruments	31	158.23	9.71	129.00-187.00	6.14																				
SG: Select Roche WBG Meters	283	168.87	17.22	117.00-221.00	10.18	48	69.73	2.15	63.00-77.00	3.14	6	275.50	8.94	249.00-303.00	3.22	6	163.33	10.97	130.00-196.00	6.75	6	70.00	2.19	63.00-77.00	3.14
SG: Select Bayer WBG Meters	168	134.88	15.71	88.00-182.00	11.63	21	47.48	6.98	26.00-68.00	14.89															
PG: Nova Biomedical StatStrip Xpress	10	147.80	8.54	123.00-174.00	5.74																				
PG: Nova Biomedical StatStrip Glucose meter	5	144.60	0.55	144.00-147.00	0.34																				
SG: Select Nova Biomedical StatStrip Instruments	15	146.73	7.03	126.00-168.00	4.76																				
PG: Bayer Contour Next	19	135.00	12.76	97.00-173.00	9.48	6	51.33	7.29	29.00-73.00	14.31															
PG: Roche AccuChek Performa	149	168.58	21.19	105.00-233.00	12.54	5	69.60	2.07	64.00-76.00	3.00															
PG: True Metrix Blood Glucose Meter	13	149.15	7.23	127.00-171.00	4.83																				
PG: McKesson True Metrix Pro Blood Glucose Meter	39	147.69	10.14	118.00-178.00	6.82	6	55.50	2.59	48.00-64.00	4.64															
SG: Select Arkray Whole Blood Glucose	6	184.17	29.68	95.00-273.00	16.14																				
PG: EvenCare G3 Blood Glucose System	6	171.50	9.91	142.00-202.00	5.76																				
SG: Select EvenCare Whole Blood Glucose	9	169.67	8.72	144.00-196.00	5.12																				
AG: All Group	659	159.34	29.67	70.00-248.00	18.68						20	258.60	27.91	175.00-343.00	10.77	20	153.70	16.46	105.00-204.00	10.71	20	56.15	10.70	24.00-88.00	19.11



789 - Group B Screen	Labs	Response	Number	Percent
		Strep B Ag Detect		
AD-11		* Denotes Correct Response		
AG: All Methods	1			
		*Positive	1	100.00%
Cepheid GeneXpert	1	*Positive	1	100.00%
AD-12		* Denotes Correct Response		
AG: All Methods	1			
		*Negative	1	100.00%
Cepheid GeneXpert	1	*Negative	1	100.00%
AD-13		* Denotes Correct Response		
AG: All Methods	1			
		*Positive	1	100.00%
Cepheid GeneXpert	1	*Positive	1	100.00%
AD-14		* Denotes Correct Response		
AG: All Methods	1			
		*Negative	1	100.00%
Cepheid GeneXpert	1	*Negative	1	100.00%
AD-15		* Denotes Correct Response		
AG: All Methods	1			
		*Positive	1	100.00%
Cepheid GeneXpert	1	*Positive	1	100.00%
760 - Anti-nuclear Antibodies (ANA)	Labs	Response	Number	Percent
		ANA (ELISA) (International)		
AE-11		* Denotes Correct Response		
AG: All Methods	1			
		*Negative	1	100.00%
Dynex DS2	1	*Negative	1	100.00%
AE-12		* Denotes Correct Response		
AG: All Methods	1			
		*Negative	1	100.00%
Dynex DS2	1	*Negative	1	100.00%
AE-13		* Denotes Correct Response		
AG: All Methods	1			
		*Positive	1	100.00%
Dynex DS2	1	*Positive	1	100.00%
AE-14		* Denotes Correct Response		
AG: All Methods	1			
		*Negative	1	100.00%
Dynex DS2	1	*Negative	1	100.00%
AE-15		* Denotes Correct Response		
AG: All Methods	1			
		*Positive	1	100.00%
Dynex DS2	1	*Positive	1	100.00%
765 - Allergen Testing	Labs	Response	Number	Percent
		Ash - White Alrgn		
AL-5		* Denotes Correct Response		
PG: Hitachi CLA-1	1			
		0	1	100.00%
AL-6		* Denotes Correct Response		
PG: Hitachi CLA-1	2			
		*0	1	50.00%
		*1	1	50.00%
		Birch/Alder Mix Alrgn		
AL-5		* Denotes Correct Response		
PG: Hitachi CLA-1	2			
		*1	2	100.00%
AL-6		* Denotes Correct Response		
PG: Hitachi CLA-1	3			
		*1	1	33.33%
		*2	2	66.67%
		Cedar - Mountain Alrgn		
AL-5		* Denotes Correct Response		
PG: Hitachi CLA-1	2			
		*1/0	2	100.00%
AL-6		* Denotes Correct Response		
PG: Hitachi CLA-1	3			



		*1	2	66.67%
		*2	1	33.33%
		Cottonwood - Eastern Alrgn		
AL-5		* Denotes Correct Response		
PG: Hitachi CLA-1	2			
		*0	1	50.00%
		*1/0	1	50.00%
AL-6		* Denotes Correct Response		
PG: Hitachi CLA-1	3			
		*1	1	33.33%
		*1/0	2	66.67%
		Elm - White Alrgn		
AL-5		* Denotes Correct Response		
PG: Hitachi CLA-1	2			
		*0	1	50.00%
		*1/0	1	50.00%
AL-6		* Denotes Correct Response		
PG: Hitachi CLA-1	3			
		*1	2	66.67%
		*2	1	33.33%
		Mulberry Mix Alrgn		
AL-5		* Denotes Correct Response		
PG: Hitachi CLA-1	2			
		*0	1	50.00%
		*1/0	1	50.00%
AL-6		* Denotes Correct Response		
PG: Hitachi CLA-1	3			
		*1	2	66.67%
		*2	1	33.33%
		Pine Mix Alrgn		
AL-5		* Denotes Correct Response		
PG: Hitachi CLA-1	1			
		0	1	100.00%
AL-6		* Denotes Correct Response		
PG: Hitachi CLA-1	2			
		*1	1	50.00%
		*1/0	1	50.00%
		Privet Alrgn		
AL-5		* Denotes Correct Response		
PG: Hitachi CLA-1	1			
		1/0	1	100.00%
AL-6		* Denotes Correct Response		
PG: Hitachi CLA-1	2			
		*0	1	50.00%
		*1	1	50.00%
		Sycamore - American Alrgn		
AL-5		* Denotes Correct Response		
PG: Hitachi CLA-1	1			
		1	1	100.00%
AL-6		* Denotes Correct Response		
PG: Hitachi CLA-1	2			
		*1	1	50.00%
		*2	1	50.00%
		Cocklebur Alrgn		
AL-5		* Denotes Correct Response		
PG: Hitachi CLA-1	1			
		0	1	100.00%
AL-6		* Denotes Correct Response		
PG: Hitachi CLA-1	2			
		*1	1	50.00%
		*1/0	1	50.00%
		Lambs Quarters Alrgn		
AL-5		* Denotes Correct Response		
PG: Hitachi CLA-1	1			
		1	1	100.00%
AL-6		* Denotes Correct Response		
PG: Hitachi CLA-1	2			



		*1	1	50.00%
		*2	1	50.00%
		Marshelder - Rough Alrgn		
AL-5		* Denotes Correct Response		
PG: Hitachi CLA-1	1			
		1/0	1	100.00%
AL-6		* Denotes Correct Response		
PG: Hitachi CLA-1	2			
		*1	1	50.00%
		*2	1	50.00%
		Pigweed Alrgn		
AL-5		* Denotes Correct Response		
PG: Hitachi CLA-1	2			
		*1/0	2	100.00%
AL-6		* Denotes Correct Response		
PG: Hitachi CLA-1	3			
		*1	3	100.00%
		Ragweed - Short Alrgn		
AL-5		* Denotes Correct Response		
PG: Hitachi CLA-1	3			
		*1	1	33.33%
		*2	1	33.33%
		*3	1	33.33%
AL-6		* Denotes Correct Response		
PG: Hitachi CLA-1	4			
		*1	1	25.00%
		*2	1	25.00%
		*3	1	25.00%
		*4	1	25.00%
		Sheep Sorrel Alrgn		
AL-5		* Denotes Correct Response		
PG: Hitachi CLA-1	2			
		*2	1	50.00%
		*3	1	50.00%
AL-6		* Denotes Correct Response		
PG: Hitachi CLA-1	3			
		*1	1	33.33%
		*2	1	33.33%
		*3	1	33.33%
		Waterhemp Alrgn		
AL-5		* Denotes Correct Response		
PG: Hitachi CLA-1	1			
		2	1	100.00%
AL-6		* Denotes Correct Response		
PG: Hitachi CLA-1	2			
		*2	1	50.00%
		*3	1	50.00%
		Mite D Farinae Alrgn		
AL-5		* Denotes Correct Response		
PG: Hitachi CLA-1	2			
		*4	2	100.00%
AL-6		* Denotes Correct Response		
PG: Hitachi CLA-1	4			
		*4	3	75.00%
		3	1	25.00%
		Cladosporium Alrgn		
AL-5		* Denotes Correct Response		
PG: Hitachi CLA-1	2			
		*1	2	100.00%
AL-6		* Denotes Correct Response		
PG: Hitachi CLA-1	4			
		*3	3	75.00%
		*4	1	25.00%
		Tomato Alrgn		
AL-5		* Denotes Correct Response		
PG: Hitachi CLA-1	2			
		*1/0	2	100.00%



AL-6		* Denotes Correct Response		
PG: Hitachi CLA-1	1			
		1/0	1	100.00%
		Wheat Alrgn		
AL-5		* Denotes Correct Response		
PG: Hitachi CLA-1	2			
		*1	2	100.00%
AL-6		* Denotes Correct Response		
PG: Hitachi CLA-1	1			
		1	1	100.00%
		Vegetable Mix Alrgn		
AL-5		* Denotes Correct Response		
PG: Hitachi CLA-1	2			
		*1/0	2	100.00%
AL-6		* Denotes Correct Response		
PG: Hitachi CLA-1	1			
		1/0	1	100.00%
		Tuna Alrgn		
AL-5		* Denotes Correct Response		
PG: Hitachi CLA-1	2			
		*0	2	100.00%
AL-6		* Denotes Correct Response		
PG: Hitachi CLA-1	1			
		0	1	100.00%
		Shellfish Mix Alrgn		
AL-5		* Denotes Correct Response		
PG: Hitachi CLA-1	2			
		*4	2	100.00%
AL-6		* Denotes Correct Response		
PG: Hitachi CLA-1	1			
		4	1	100.00%
		Rice Alrgn		
AL-5		* Denotes Correct Response		
PG: Hitachi CLA-1	2			
		*1	1	50.00%
		*1/0	1	50.00%
AL-6		* Denotes Correct Response		
PG: Hitachi CLA-1	1			
		1	1	100.00%
		Potato Alrgn		
AL-5		* Denotes Correct Response		
PG: Hitachi CLA-1	2			
		*1	1	50.00%
		*2	1	50.00%
AL-6		* Denotes Correct Response		
PG: Hitachi CLA-1	1			
		2	1	100.00%
		Pork Alrgn		
AL-5		* Denotes Correct Response		
PG: Hitachi CLA-1	2			
		*2	1	50.00%
		*3	1	50.00%
AL-6		* Denotes Correct Response		
PG: Hitachi CLA-1	1			
		2	1	100.00%
		Orange Alrgn		
AL-5		* Denotes Correct Response		
PG: Hitachi CLA-1	2			
		*2	2	100.00%
AL-6		* Denotes Correct Response		
PG: Hitachi CLA-1	1			
		2	1	100.00%
		Onion Mix Alrgn		
AL-5		* Denotes Correct Response		
PG: Hitachi CLA-1	2			
		*1	2	100.00%
AL-6		* Denotes Correct Response		



PG: Hitachi CLA-1	1			
		1	1	100.00%
		Oat Alrgn		
AL-5		* Denotes Correct Response		
PG: Hitachi CLA-1	2			
		*1	1	50.00%
		*2	1	50.00%
AL-6		* Denotes Correct Response		
PG: Hitachi CLA-1	1			
		1/0	1	100.00%
		Garlic Alrgn		
AL-5		* Denotes Correct Response		
PG: Hitachi CLA-1	2			
		*2	2	100.00%
AL-6		* Denotes Correct Response		
PG: Hitachi CLA-1	1			
		1	1	100.00%
		Egg - whole Alrgn		
AL-5		* Denotes Correct Response		
PG: Hitachi CLA-1	3			
		*3	2	66.67%
		*4	1	33.33%
AL-6		* Denotes Correct Response		
PG: Hitachi CLA-1	1			
		1	1	100.00%
		Corn Alrgn		
AL-5		* Denotes Correct Response		
PG: Hitachi CLA-1	2			
		*1	1	50.00%
		*1/0	1	50.00%
AL-6		* Denotes Correct Response		
PG: Hitachi CLA-1	1			
		1/0	1	100.00%
		Chocolate Alrgn		
AL-5		* Denotes Correct Response		
PG: Hitachi CLA-1	2			
		*0	2	100.00%
AL-6		* Denotes Correct Response		
PG: Hitachi CLA-1	1			
		0	1	100.00%
		Chicken Alrgn		
AL-5		* Denotes Correct Response		
PG: Hitachi CLA-1	2			
		*0	1	50.00%
		*1/0	1	50.00%
AL-6		* Denotes Correct Response		
PG: Hitachi CLA-1	1			
		0	1	100.00%
		Beef Alrgn		
AL-5		* Denotes Correct Response		
PG: Hitachi CLA-1	2			
		*1/0	2	100.00%
AL-6		* Denotes Correct Response		
PG: Hitachi CLA-1	1			
		2	1	100.00%
		Bean - White Alrgn		
AL-5		* Denotes Correct Response		
PG: Hitachi CLA-1	2			
		*1	1	50.00%
		*1/0	1	50.00%
AL-6		* Denotes Correct Response		
PG: Hitachi CLA-1	1			
		1	1	100.00%
		Barley Alrgn		
AL-5		* Denotes Correct Response		
PG: Hitachi CLA-1	2			
		*0	1	50.00%



		*1/0	1	50.00%
AL-6		* Denotes Correct Response		
PG: Hitachi CLA-1	1			
		0	1	100.00%
		Apple Alrgn		
AL-5		* Denotes Correct Response		
PG: Hitachi CLA-1	3			
		*0	1	33.33%
		*1	1	33.33%
		1/0	1	33.33%
AL-6		* Denotes Correct Response		
PG: Hitachi CLA-1	1			
		2	1	100.00%
		Almond Alrgn		
AL-5		* Denotes Correct Response		
PG: Hitachi CLA-1	2			
		*1	2	100.00%
AL-6		* Denotes Correct Response		
PG: Hitachi CLA-1	1			
		1	1	100.00%
		Alternaria alternata Alrgn		
AL-5		* Denotes Correct Response		
PG: Hitachi CLA-1	2			
		*2	2	100.00%
AL-6		* Denotes Correct Response		
PG: Hitachi CLA-1	3			
		*2	2	66.67%
		*3	1	33.33%
		Aspergillus fumigatus Alrgn		
AL-5		* Denotes Correct Response		
PG: Hitachi CLA-1	2			
		*1	2	100.00%
AL-6		* Denotes Correct Response		
PG: Hitachi CLA-1	3			
		*2	3	100.00%
		Bermuda Grass Alrgn		
AL-5		* Denotes Correct Response		
PG: Hitachi CLA-1	2			
		*2	2	100.00%
AL-6		* Denotes Correct Response		
PG: Hitachi CLA-1	3			
		*2	2	66.67%
		*3	1	33.33%
		Cat Epithelium Alrgn		
AL-5		* Denotes Correct Response		
PG: Hitachi CLA-1	2			
		*4	1	50.00%
		3	1	50.00%
AL-6		* Denotes Correct Response		
PG: Hitachi CLA-1	1			
		4	1	100.00%
		Common Ragweed Alrgn		
AL-5		* Denotes Correct Response		
PG: Hitachi CLA-1	1			
		3	1	100.00%
		Dog Epithelium Alrgn		
AL-5		* Denotes Correct Response		
PG: Hitachi CLA-1	1			
		3	1	100.00%
AL-6		* Denotes Correct Response		
PG: Hitachi CLA-1	1			
		*4	1	100.00%
		Dog Dander Alrgn		
AL-5		* Denotes Correct Response		
PG: Hitachi CLA-1	1			
		3	1	100.00%
AL-6		* Denotes Correct Response		



PG: Hitachi CLA-1	3			
		*2	1	33.33%
		*3	2	66.67%
		House Dust Mite Alrgn		
AL-5		* Denotes Correct Response		
PG: Hitachi CLA-1	3			
		*1	1	33.33%
		*3	1	33.33%
		2	1	33.33%
AL-6		* Denotes Correct Response		
PG: Hitachi CLA-1	4			
		*1	2	50.00%
		*3	1	25.00%
		2	1	25.00%
		Timothy Grass Alrgn		
AL-5		* Denotes Correct Response		
PG: Hitachi CLA-1	2			
		*4	2	100.00%
AL-6		* Denotes Correct Response		
PG: Hitachi CLA-1	4			
		*3	1	25.00%
		*4	3	75.00%
		Box Elder, Maple Alrgn		
AL-5		* Denotes Correct Response		
PG: Hitachi CLA-1	1			
		1	1	100.00%
AL-6		* Denotes Correct Response		
PG: Hitachi CLA-1	2			
		*1	1	50.00%
		*2	1	50.00%
		Peanut (Arachis hyogaea) Alrgn		
AL-5		* Denotes Correct Response		
PG: Hitachi CLA-1	2			
		*4	2	100.00%
AL-6		* Denotes Correct Response		
PG: Hitachi CLA-1	2			
		*4	1	50.00%
		3	1	50.00%
		Cat Dander Alrgn		
AL-5		* Denotes Correct Response		
PG: Hitachi CLA-1	2			
		*4	1	50.00%
		3	1	50.00%
AL-6		* Denotes Correct Response		
PG: Hitachi CLA-1	2			
		*4	2	100.00%
		Soybean Alrgn		
AL-5		* Denotes Correct Response		
PG: Hitachi CLA-1	2			
		*3	2	100.00%
AL-6		* Denotes Correct Response		
PG: Hitachi CLA-1	1			
		2	1	100.00%
		Cockroach Mix Alrgn		
AL-5		* Denotes Correct Response		
PG: Hitachi CLA-1	1			
		1	1	100.00%
AL-6		* Denotes Correct Response		
PG: Hitachi CLA-1	2			
		*2	2	100.00%
		Oak, White Alrgn		
AL-5		* Denotes Correct Response		
PG: Hitachi CLA-1	1			
		2	1	100.00%
AL-6		* Denotes Correct Response		
PG: Hitachi CLA-1	2			
		*3	2	100.00%



		Penicillium Alrgn		
AL-5		* Denotes Correct Response		
PG: Hitachi CLA-1	2			
		0	1	50.00%
		1/0	1	50.00%
AL-6		* Denotes Correct Response		
PG: Hitachi CLA-1	2			
		*1	2	100.00%
		Walnut/Hickory/Pecan Mix Alrgn		
AL-5		* Denotes Correct Response		
PG: Hitachi CLA-1	2			
		*1/0	2	100.00%
AL-6		* Denotes Correct Response		
PG: Hitachi CLA-1	2			
		*1	1	50.00%
		*2	1	50.00%
		Yeast - Bakers Alrgn		
AL-5		* Denotes Correct Response		
PG: Hitachi CLA-1	2			
		*0	2	100.00%
AL-6		* Denotes Correct Response		
PG: Hitachi CLA-1	1			
		0	1	100.00%
683 - Antistreptolysin O (ASO)	Labs	Response	Number	Percent
		ASO (Qual) (International)		
AS-11		* Denotes Correct Response		
PG: Wampole ColorCard ASO	1			
		*Negative	1	100.00%
AS-12		* Denotes Correct Response		
PG: Wampole ColorCard ASO	1			
		*Positive	1	100.00%
AS-13		* Denotes Correct Response		
PG: Wampole ColorCard ASO	1			
		*Positive	1	100.00%
AS-14		* Denotes Correct Response		
PG: Wampole ColorCard ASO	1			
		*Positive	1	100.00%
AS-15		* Denotes Correct Response		
PG: Wampole ColorCard ASO	1			
		*Negative	1	100.00%
782 - Complete Bacteriology	Labs	Response	Number	Percent
		Stool Culture		
BA-7		* Denotes Correct Response		
AG: All Methods	1			
		*Campylobacter jejuni	1	100.00%
		*Escherichia coli	1	100.00%
Cultures	1	*Campylobacter jejuni	1	100.00%
Cultures		*Escherichia coli	1	100.00%
		Wound Culture		
BA-8		* Denotes Correct Response		
AG: All Methods	2			
		*Enterobacter aerogenes	1	50.00%
		*Gram Negative organism	1	50.00%
Cultures	2	*Enterobacter aerogenes	1	50.00%
Cultures		*Gram Negative organism	1	50.00%
975 - Bacterial Vaginosis	Labs	Response	Number	Percent
		Bacterial Vaginosa (Sialidase)		
BV-5		* Denotes Correct Response		
AG: All Methods	98			
		*Positive	84	85.71%
		Negative	14	14.29%
Sekisui Diagnostics OSOM BVBlue	98	*Positive	84	85.71%
Sekisui Diagnostics OSOM BVBlue		Negative	14	14.29%
BV-6		* Denotes Correct Response		
AG: All Methods	98			
		*Negative	88	89.80%
		Positive	10	10.20%



Sekisui Diagnostics OSOM BVBlue	98	*Negative	88	89.80%
Sekisui Diagnostics OSOM BVBlue		Positive	10	10.20%
661 - Urine Colony Count;784 - Presump ID/Colony Count - Urine;775 - Urine Culture;778 - Urine & Throat Cultures;779 - Urine & Genital Cultures;781 - Urine, Throat & Genital Cultures;782 - Complete Bacteriology;790 - Urine Susceptibility Add-On	Labs	Response	Number	Percent
		Colony Count		
CC-11		* Denotes Correct Response		
AG: All Methods	31			
		*10,000 - 100,000	2	6.45%
		*Greater than 100,000	29	93.55%
Calibrated Loop - Manual	10	*10,000 - 100,000	1	10.00%
Uricult (Orion)	12	*10,000 - 100,000	1	8.33%
Calibrated Loop - Manual	10	*Greater than 100,000	9	90.00%
Healthlink URI-ID	1	*Greater than 100,000	1	100.00%
Uri-Chek (Troy Bio)	6	*Greater than 100,000	6	100.00%
Uricult (Orion)	12	*Greater than 100,000	11	91.67%
Uri-Three (Culture Kits)	1	*Greater than 100,000	1	100.00%
Uricult (Lifesign)	1	*Greater than 100,000	1	100.00%
784 - Presump ID/Colony Count - Urine	Labs	Response	Number	Percent
		Colony Count		
CC-11		* Denotes Correct Response		
AG: All Methods	31			
		*10,000 - 100,000	2	6.45%
		*Greater than 100,000	29	93.55%
Calibrated Loop - Manual	10	*10,000 - 100,000	1	10.00%
Uricult (Orion)	12	*10,000 - 100,000	1	8.33%
Calibrated Loop - Manual	10	*Greater than 100,000	9	90.00%
Healthlink URI-ID	1	*Greater than 100,000	1	100.00%
Uri-Chek (Troy Bio)	6	*Greater than 100,000	6	100.00%
Uricult (Orion)	12	*Greater than 100,000	11	91.67%
Uri-Three (Culture Kits)	1	*Greater than 100,000	1	100.00%
Uricult (Lifesign)	1	*Greater than 100,000	1	100.00%
CC-12		* Denotes Correct Response		
AG: All Methods	5			
		*Greater than 100,000	5	100.00%
Calibrated Loop - Manual	4	*Greater than 100,000	4	100.00%
Uri-Chek (Troy Bio)	1	*Greater than 100,000	1	100.00%
CC-14		* Denotes Correct Response		
AG: All Methods	5			
		*10,000 - 100,000	1	20.00%
		*Greater than 100,000	4	80.00%
Calibrated Loop - Manual	4	*10,000 - 100,000	1	25.00%
Calibrated Loop - Manual		*Greater than 100,000	3	75.00%
Uri-Chek (Troy Bio)	1	*Greater than 100,000	1	100.00%
CC-15		* Denotes Correct Response		
AG: All Methods	5			
		*Greater than 100,000	5	100.00%
Calibrated Loop - Manual	4	*Greater than 100,000	4	100.00%
Uri-Chek (Troy Bio)	1	*Greater than 100,000	1	100.00%
		Urine Presump Id		
CC-11		* Denotes Correct Response		
AG: All Methods	3			
		*Gram Negative Organsim	1	33.33%
		*Presumptive Enterobacter sp.	2	66.67%
Calibrated Loop - Manual	3	*Gram Negative Organsim	1	33.33%
Calibrated Loop - Manual		*Presumptive Enterobacter sp.	2	66.67%
CC-12		* Denotes Correct Response		
AG: All Methods	3			
		*Gram Negative Organsim	1	33.33%
		*Presumptive Klebsiella sp.	2	66.67%
Calibrated Loop - Manual	3	*Gram Negative Organsim	1	33.33%
Calibrated Loop - Manual		*Presumptive Klebsiella sp.	2	66.67%
CC-14		* Denotes Correct Response		
AG: All Methods	3			
		*Gram Negative Organsim	1	33.33%
		*Presumptive E. Coli	2	66.67%
Calibrated Loop - Manual	3	*Gram Negative Organsim	1	33.33%



Calibrated Loop - Manual		*Presumptive E. Coli	2	66.67%
CC-15		* Denotes Correct Response		
AG: All Methods	3			
		*Growth, Referred for ID	1	33.33%
		*Presumptive Morganela sp	1	33.33%
		Gram Positive Organism	1	33.33%
Calibrated Loop - Manual	3	*Growth, Referred for ID	1	33.33%
Calibrated Loop - Manual		*Presumptive Morganela sp	1	33.33%
Calibrated Loop - Manual		Gram Positive Organism	1	33.33%
		Susceptibility Amox/Clavulanate		
CC-11		* Denotes Correct Response		
AG: All Methods	3			
		*Resistant	3	100.00%
Susceptibility	3	*Resistant	3	100.00%
		Susceptibility Ampicillin		
CC-11		* Denotes Correct Response		
AG: All Methods	5			
		*Intermediate	1	20.00%
		*Resistant	4	80.00%
Susceptibility	5	*Intermediate	1	20.00%
Susceptibility		*Resistant	4	80.00%
		Susceptibility Ampicillin/Sulbactam		
CC-11		* Denotes Correct Response		
AG: All Methods	1			
		*Intermediate	1	100.00%
Susceptibility	1	*Intermediate	1	100.00%
		Susceptibility Cefazolin		
CC-11		* Denotes Correct Response		
AG: All Methods	1			
		*Resistant	1	100.00%
Susceptibility	1	*Resistant	1	100.00%
		Susceptibility Cefepime		
CC-11		* Denotes Correct Response		
AG: All Methods	1			
		*Sensitive	1	100.00%
Susceptibility	1	*Sensitive	1	100.00%
		Susceptibility Cefpodoxime		
CC-11		* Denotes Correct Response		
AG: All Methods	1			
		*Sensitive	1	100.00%
Susceptibility	1	*Sensitive	1	100.00%
		Susceptibility Ceftazidime		
CC-11		* Denotes Correct Response		
AG: All Methods	1			
		*Sensitive	1	100.00%
Susceptibility	1	*Sensitive	1	100.00%
		Susceptibility Ceftriaxone		
CC-11		* Denotes Correct Response		
AG: All Methods	3			
		*Sensitive	3	100.00%
Susceptibility	3	*Sensitive	3	100.00%
		Susceptibility Cephalothin		
CC-11		* Denotes Correct Response		
AG: All Methods	4			
		*Resistant	4	100.00%
Susceptibility	4	*Resistant	4	100.00%
		Susceptibility Ciprofloxacin		
CC-11		* Denotes Correct Response		
AG: All Methods	6			
		*Sensitive	6	100.00%
Susceptibility	6	*Sensitive	6	100.00%
		Susceptibility Gentamicin		
CC-11		* Denotes Correct Response		
AG: All Methods	4			
		*Sensitive	4	100.00%
Susceptibility	4	*Sensitive	4	100.00%
		Susceptibility Imipenem		



CC-11		* Denotes Correct Response		
AG: All Methods	2			
		*Sensitive	2	100.00%
Susceptibility	2	*Sensitive	2	100.00%
		Susceptibility Levofloxacin		
CC-11		* Denotes Correct Response		
AG: All Methods	3			
		*Sensitive	3	100.00%
Susceptibility	3	*Sensitive	3	100.00%
		Susceptibility Minocycline		
CC-11		* Denotes Correct Response		
AG: All Methods	1			
		*Sensitive	1	100.00%
Susceptibility	1	*Sensitive	1	100.00%
		Susceptibility Piperacillin		
CC-11		* Denotes Correct Response		
AG: All Methods	1			
		*Sensitive	1	100.00%
Susceptibility	1	*Sensitive	1	100.00%
		Susceptibility Piperacillin/Tazobactam		
CC-11		* Denotes Correct Response		
AG: All Methods	1			
		*Sensitive	1	100.00%
Susceptibility	1	*Sensitive	1	100.00%
		Susceptibility Sulfamethoxazole		
CC-11		* Denotes Correct Response		
AG: All Methods	1			
		*Sensitive	1	100.00%
Susceptibility	1	*Sensitive	1	100.00%
		Susceptibility Tetracycline		
CC-11		* Denotes Correct Response		
AG: All Methods	4			
		*Intermediate	1	25.00%
		*Sensitive	3	75.00%
Susceptibility	4	*Intermediate	1	25.00%
Susceptibility		*Sensitive	3	75.00%
		Susceptibility Tobramycin		
CC-11		* Denotes Correct Response		
AG: All Methods	2			
		*Sensitive	2	100.00%
Susceptibility	2	*Sensitive	2	100.00%
		Susceptibility Trimethoprim		
CC-11		* Denotes Correct Response		
AG: All Methods	1			
		*Sensitive	1	100.00%
Susceptibility	1	*Sensitive	1	100.00%
		Susceptibility Trimeth./Sulfameth.		
CC-11		* Denotes Correct Response		
AG: All Methods	4			
		*Sensitive	4	100.00%
Susceptibility	4	*Sensitive	4	100.00%
964 - C. difficile Toxin/Ag (2 specimens);963 - C. difficile Toxin/Ag (5 specimens)	Labs	Response	Number	Percent
		C. difficile Toxin Detect		
CDF-11		* Denotes Correct Response		
AG: All Methods	1			
		*Negative	1	100.00%
Cepheid GeneXpert - Xpert C. diff	1	*Negative	1	100.00%
CDF-12		* Denotes Correct Response		
AG: All Methods	1			
		*Negative	1	100.00%
Cepheid GeneXpert - Xpert C. diff	1	*Negative	1	100.00%
795 - Vaginal Pathogens II (2 specimens);794 - Vaginal Pathogens II (5 specimens);977 - Trichomonas Ag Det (Wavied)	Labs	Response	Number	Percent
		T. vaginalis - Non Culture Tech		
CGT-11		* Denotes Correct Response		



AG: All Methods	3			
		*Positive	3	100.00%
Cepheid GeneXpert	2	*Positive	2	100.00%
BD Max	1	*Positive	1	100.00%
CGT-12		* Denotes Correct Response		
AG: All Methods	3			
		*Positive	3	100.00%
Cepheid GeneXpert	2	*Positive	2	100.00%
BD Max	1	*Positive	1	100.00%
		Chlamydia - Non Culture Tech		
CGT-11		* Denotes Correct Response		
AG: All Methods	24			
		*Negative	24	100.00%
Quidel QuickVue Chlamydia	2	*Negative	2	100.00%
Cepheid GeneXpert	20	*Negative	20	100.00%
BD Viper XTR	1	*Negative	1	100.00%
BD Max	1	*Negative	1	100.00%
CGT-12		* Denotes Correct Response		
AG: All Methods	24			
		*Positive	24	100.00%
Quidel QuickVue Chlamydia	2	*Positive	2	100.00%
Cepheid GeneXpert	20	*Positive	20	100.00%
BD Viper XTR	1	*Positive	1	100.00%
BD Max	1	*Positive	1	100.00%
		GC - Non Culture Tech		
CGT-11		* Denotes Correct Response		
AG: All Methods	22			
		*Positive	22	100.00%
Cepheid GeneXpert	20	*Positive	20	100.00%
BD Viper XTR	1	*Positive	1	100.00%
BD Max	1	*Positive	1	100.00%
CGT-12		* Denotes Correct Response		
AG: All Methods	22			
		*Negative	22	100.00%
Cepheid GeneXpert	20	*Negative	20	100.00%
BD Viper XTR	1	*Negative	1	100.00%
BD Max	1	*Negative	1	100.00%
794 - Vaginal Pathogens II (5 specimens)	Labs	Response	Number	Percent
		Chlamydia - Non Culture Tech		
CGT-11		* Denotes Correct Response		
AG: All Methods	24			
		*Negative	24	100.00%
Quidel QuickVue Chlamydia	2	*Negative	2	100.00%
Cepheid GeneXpert	20	*Negative	20	100.00%
BD Viper XTR	1	*Negative	1	100.00%
BD Max	1	*Negative	1	100.00%
CGT-12		* Denotes Correct Response		
AG: All Methods	24			
		*Positive	24	100.00%
Quidel QuickVue Chlamydia	2	*Positive	2	100.00%
Cepheid GeneXpert	20	*Positive	20	100.00%
BD Viper XTR	1	*Positive	1	100.00%
BD Max	1	*Positive	1	100.00%
CGT-13		* Denotes Correct Response		
AG: All Methods	23			
		*Positive	23	100.00%
Quidel QuickVue Chlamydia	2	*Positive	2	100.00%
Cepheid GeneXpert	19	*Positive	19	100.00%
BD Viper XTR	1	*Positive	1	100.00%
BD Max	1	*Positive	1	100.00%
CGT-14		* Denotes Correct Response		
AG: All Methods	23			
		*Positive	23	100.00%
Quidel QuickVue Chlamydia	2	*Positive	2	100.00%
Cepheid GeneXpert	19	*Positive	19	100.00%
BD Viper XTR	1	*Positive	1	100.00%
BD Max	1	*Positive	1	100.00%
CGT-15		* Denotes Correct Response		



AG: All Methods	23			
		*Negative	22	95.65%
		Positive	1	4.35%
Quidel QuickVue Chlamydia	2	*Negative	2	100.00%
Cepheid GeneXpert	19	*Negative	18	94.74%
BD Viper XTR	1	*Negative	1	100.00%
BD Max	1	*Negative	1	100.00%
Cepheid GeneXpert	19	Positive	1	5.26%
		GC - Non Culture Tech		
CGT-11		* Denotes Correct Response		
AG: All Methods	22			
		*Positive	22	100.00%
Cepheid GeneXpert	20	*Positive	20	100.00%
BD Viper XTR	1	*Positive	1	100.00%
BD Max	1	*Positive	1	100.00%
CGT-12		* Denotes Correct Response		
AG: All Methods	22			
		*Negative	22	100.00%
Cepheid GeneXpert	20	*Negative	20	100.00%
BD Viper XTR	1	*Negative	1	100.00%
BD Max	1	*Negative	1	100.00%
CGT-13		* Denotes Correct Response		
AG: All Methods	21			
		*Positive	21	100.00%
Cepheid GeneXpert	19	*Positive	19	100.00%
BD Viper XTR	1	*Positive	1	100.00%
BD Max	1	*Positive	1	100.00%
CGT-14		* Denotes Correct Response		
AG: All Methods	21			
		*Negative	21	100.00%
Cepheid GeneXpert	19	*Negative	19	100.00%
BD Viper XTR	1	*Negative	1	100.00%
BD Max	1	*Negative	1	100.00%
CGT-15		* Denotes Correct Response		
AG: All Methods	21			
		*Positive	20	95.24%
		Negative	1	4.76%
Cepheid GeneXpert	19	*Positive	18	94.74%
BD Viper XTR	1	*Positive	1	100.00%
BD Max	1	*Positive	1	100.00%
Cepheid GeneXpert	19	Negative	1	5.26%
		T. vaginalis - Non Culture Tech		
CGT-11		* Denotes Correct Response		
AG: All Methods	3			
		*Positive	3	100.00%
Cepheid GeneXpert	2	*Positive	2	100.00%
BD Max	1	*Positive	1	100.00%
CGT-12		* Denotes Correct Response		
AG: All Methods	3			
		*Positive	3	100.00%
Cepheid GeneXpert	2	*Positive	2	100.00%
BD Max	1	*Positive	1	100.00%
CGT-13		* Denotes Correct Response		
AG: All Methods	3			
		*Negative	3	100.00%
Cepheid GeneXpert	2	*Negative	2	100.00%
BD Max	1	*Negative	1	100.00%
CGT-14		* Denotes Correct Response		
AG: All Methods	3			
		*Negative	3	100.00%
Cepheid GeneXpert	2	*Negative	2	100.00%
BD Max	1	*Negative	1	100.00%
CGT-15		* Denotes Correct Response		
AG: All Methods	1			
		*Negative	1	100.00%
BD Max	1	*Negative	1	100.00%
		Trichomonas Ag Det (Waived)		
CGT-11		* Denotes Correct Response		



AG: All Methods	95			
		*Positive	86	90.53%
		Negative	9	9.47%
Sekisui Diagnostics OSOM Trichomonas	95	*Positive	86	90.53%
Sekisui Diagnostics OSOM Trichomonas		Negative	9	9.47%
CGT-12		* Denotes Correct Response		
AG: All Methods	95			
		*Positive	83	87.37%
		Negative	12	12.63%
Sekisui Diagnostics OSOM Trichomonas	95	*Positive	83	87.37%
Sekisui Diagnostics OSOM Trichomonas		Negative	12	12.63%
612 - MultiChem 5;616 - MultiChem 10;617 - MultiChem 15;618 - MultiChem 20;619 - MultiChem 21+;676 - Endocrinology	Labs	Response	Number	Percent
		TSH (Semi-Quant)		
CH-11		* Denotes Correct Response		
PG: Biomedical Thyro Test	1			
		<= 5	1	100.00%
CH-12		* Denotes Correct Response		
PG: Biomedical Thyro Test	1			
		<= 5	1	100.00%
CH-13		* Denotes Correct Response		
PG: Biomedical Thyro Test	1			
		<= 5	1	100.00%
CH-14		* Denotes Correct Response		
PG: Biomedical Thyro Test	1			
		<= 5	1	100.00%
CH-15		* Denotes Correct Response		
PG: Biomedical Thyro Test	1			
		*<= 5	1	100.00%
639 - Complete Urinalysis Package;638 - Urine Sediment;641 - Clinical Microscopy;601 - Waived/PPM Tests II;609 - Waived/PPM Tests I	Labs	Response	Number	Percent
		Urine Sediment Id		
CM-19		* Denotes Correct Response		
AG: All Methods	601			
		*Ammonium magnesium phosphate (Triple phosphate)	567	94.34%
		*Ident Unknown - Referred	1	0.17%
		Ammonium biurate crystal	3	0.50%
		Calcium oxalate crystal	9	1.50%
		Calcium phosphate crystal	4	0.67%
		Cholesterol crystal	2	0.33%
		Hippuric acid	3	0.50%
		Hyaline cast	1	0.17%
		Red blood cell (RBC) (Erythrocyte)	5	0.83%
		Uric acid crystal	5	0.83%
		Waxy cast	1	0.17%
Microscopy	601	*Ammonium magnesium phosphate (Triple phosphate)	567	94.34%
Microscopy		*Ident Unknown - Referred	1	0.17%
Microscopy		Ammonium biurate crystal	3	0.50%
Microscopy		Calcium oxalate crystal	9	1.50%
Microscopy		Calcium phosphate crystal	4	0.67%
Microscopy		Cholesterol crystal	2	0.33%
Microscopy		Hippuric acid	3	0.50%
Microscopy		Hyaline cast	1	0.17%
Microscopy		Red blood cell (RBC) (Erythrocyte)	5	0.83%
Microscopy		Uric acid crystal	5	0.83%
Microscopy		Waxy cast	1	0.17%
CM-20		* Denotes Correct Response		
AG: All Methods	605			
		*Red blood cell (RBC) (Erythrocyte)	572	94.55%
		Ammonium magnesium phosphate (Triple phosphate)	4	0.66%
		Calcium oxalate crystal	1	0.17%
		Fat droplets or globules	2	0.33%
		Leucine crystal	3	0.50%
		RBC/Blood/Hgb cast	2	0.33%
		Renal tubular epithelial (RTE)	3	0.50%
		Starch granules	1	0.17%
		Transitional epithelial cell (urothelial)	1	0.17%



		White blood cell (WBC) (Leukocyte)	16	2.64%
Microscopy	605	*Red blood cell (RBC) (Erythrocyte)	572	94.55%
Microscopy		Ammonium magnesium phosphate (Triple phosphate)	4	0.66%
Microscopy		Calcium oxalate crystal	1	0.17%
Microscopy		Fat droplets or globules	2	0.33%
Microscopy		Leucine crystal	3	0.50%
Microscopy		RBC/Blood/Hgb cast	2	0.33%
Microscopy		Renal tubular epithelial (RTE)	3	0.50%
Microscopy		Starch granules	1	0.17%
Microscopy		Transitional epithelial cell (urothelial)	1	0.17%
Microscopy		White blood cell (WBC) (Leukocyte)	16	2.64%
CM-21		* Denotes Correct Response		
AG: All Methods	600			
		*Granular cast	508	84.67%
		Bacteria	1	0.17%
		Cellular (RTE) cast	3	0.50%
		Cellular (WBC) cast	60	10.00%
		Fatty (fat globules) cast	1	0.17%
		Hyaline cast	9	1.50%
		RBC/Blood/Hgb cast	9	1.50%
		RTE with fat globules (oval fat bodies)	1	0.17%
		Squamous epithelial cell	1	0.17%
		Waxy cast	1	0.17%
		White blood cell (WBC) (Leukocyte)	6	1.00%
Microscopy	600	*Granular cast	508	84.67%
Microscopy		Bacteria	1	0.17%
Microscopy		Cellular (RTE) cast	3	0.50%
Microscopy		Cellular (WBC) cast	60	10.00%
Microscopy		Fatty (fat globules) cast	1	0.17%
Microscopy		Hyaline cast	9	1.50%
Microscopy		RBC/Blood/Hgb cast	9	1.50%
Microscopy		RTE with fat globules (oval fat bodies)	1	0.17%
Microscopy		Squamous epithelial cell	1	0.17%
Microscopy		Waxy cast	1	0.17%
Microscopy		White blood cell (WBC) (Leukocyte)	6	1.00%
CM-22		* Denotes Correct Response		
AG: All Methods	603			
		*Ident Unknown - Referred	1	0.17%
		*White blood cell (WBC) (Leukocyte)	497	82.42%
		Bacteria	1	0.17%
		Cellular (RTE) cast	1	0.17%
		Cellular (WBC) cast	4	0.66%
		Fat droplets or globules	9	1.49%
		Granular cast	5	0.83%
		RBC/Blood/Hgb cast	2	0.33%
		Renal tubular epithelial (RTE)	47	7.79%
		RTE with fat globules (oval fat bodies)	19	3.15%
		Starch granules	1	0.17%
		Transitional epithelial cell (urothelial)	16	2.65%
Microscopy	603	*Ident Unknown - Referred	1	0.17%
Microscopy		*White blood cell (WBC) (Leukocyte)	497	82.42%
Microscopy		Bacteria	1	0.17%
Microscopy		Cellular (RTE) cast	1	0.17%
Microscopy		Cellular (WBC) cast	4	0.66%
Microscopy		Fat droplets or globules	9	1.49%
Microscopy		Granular cast	5	0.83%
Microscopy		RBC/Blood/Hgb cast	2	0.33%
Microscopy		Renal tubular epithelial (RTE)	47	7.79%
Microscopy		RTE with fat globules (oval fat bodies)	19	3.15%
Microscopy		Starch granules	1	0.17%
Microscopy		Transitional epithelial cell (urothelial)	16	2.65%
		Vaginal Wet Prep		
CM-23		* Denotes Correct Response		
AG: All Methods	870			
		*White Blood Cell (WBC) (Leukocyte)	823	94.60%
		Bacteria absent	1	0.11%
		Bacteria present	3	0.34%
		Clue cell	21	2.41%
		Red Blood Cell (RBC) (Erythrocyte)	13	1.49%



		Squamous epithelial cell	4	0.46%
		Trichomonas vaginalis	4	0.46%
		Yeast/fungi present	1	0.11%
Microscopy	870	*White Blood Cell (WBC) (Leukocyte)	823	94.60%
Microscopy		Bacteria absent	1	0.11%
Microscopy		Bacteria present	3	0.34%
Microscopy		Clue cell	21	2.41%
Microscopy		Red Blood Cell (RBC) (Erythrocyte)	13	1.49%
Microscopy		Squamous epithelial cell	4	0.46%
Microscopy		Trichomonas vaginalis	4	0.46%
Microscopy		Yeast/fungi present	1	0.11%
CM-24		* Denotes Correct Response		
AG: All Methods	870			
		*Clue cell	827	95.06%
		Bacteria present	4	0.46%
		Squamous epithelial cell	20	2.30%
		White Blood Cell (WBC) (Leukocyte)	14	1.61%
		Yeast/fungi present	5	0.57%
Microscopy	870	*Clue cell	827	95.06%
Microscopy		Bacteria present	4	0.46%
Microscopy		Squamous epithelial cell	20	2.30%
Microscopy		White Blood Cell (WBC) (Leukocyte)	14	1.61%
Microscopy		Yeast/fungi present	5	0.57%
		Skin Scraping		
CM-25		* Denotes Correct Response		
AG: All Methods	595			
		*No scabies present	1	0.17%
		*Yeast/fungi present	593	99.66%
		Scabies present	1	0.17%
Microscopy	595	*No scabies present	1	0.17%
Microscopy		*Yeast/fungi present	593	99.66%
Microscopy		Scabies present	1	0.17%
		Fern Test		
CM-26		* Denotes Correct Response		
AG: All Methods	358			
		*Positive Fern Test	358	100.00%
Microscopy	358	*Positive Fern Test	358	100.00%
772 - Dermatophyte Screen	Labs	Response	Number	Percent
		Dermatophyte Screen		
DM-11		* Denotes Correct Response		
AG: All Methods	12			
		*Negative, No Dermatophytes Present	11	91.67%
		Positive, Dermatophytes Present	1	8.33%
Dermatophytes	12	*Negative, No Dermatophytes Present	11	91.67%
Dermatophytes		Positive, Dermatophytes Present	1	8.33%
DM-12		* Denotes Correct Response		
AG: All Methods	12			
		*Positive, Dermatophytes Present	12	100.00%
Dermatophytes	12	*Positive, Dermatophytes Present	12	100.00%
DM-13		* Denotes Correct Response		
AG: All Methods	12			
		*Positive, Dermatophytes Present	12	100.00%
Dermatophytes	12	*Positive, Dermatophytes Present	12	100.00%
DM-14		* Denotes Correct Response		
AG: All Methods	12			
		*Positive, Dermatophytes Present	3	25.00%
		Negative, No Dermatophytes Present	9	75.00%
Dermatophytes	12	*Positive, Dermatophytes Present	3	25.00%
Dermatophytes		Negative, No Dermatophytes Present	9	75.00%
DM-15		* Denotes Correct Response		
AG: All Methods	12			
		*Positive, Dermatophytes Present	12	100.00%
Dermatophytes	12	*Positive, Dermatophytes Present	12	100.00%
777 - Genital Culture	Labs	Response	Number	Percent
		Genital Culture		
GC-11		* Denotes Correct Response		
AG: All Methods	2			
		*No growth, GC selective media	2	100.00%



Cultures	2	*No growth, GC selective media	2	100.00%
GC-12		* Denotes Correct Response		
AG: All Methods	2			
		*No growth, GC selective media	2	100.00%
Cultures	2	*No growth, GC selective media	2	100.00%
GC-13		* Denotes Correct Response		
AG: All Methods	2			
		*Presumptive positive for GC, GR(-) diplococci, oxidase positive organism	1	50.00%
		*Presumptive positive for GC, growth referred for ID/confirmation	1	50.00%
Cultures	2	*Presumptive positive for GC, GR(-) diplococci, oxidase positive organism	1	50.00%
Cultures		*Presumptive positive for GC, growth referred for ID/confirmation	1	50.00%
GC-14		* Denotes Correct Response		
AG: All Methods	4			
		*Group B Streptococcus present	1	25.00%
		*Negative for GC	1	25.00%
		*No growth, GC selective media	2	50.00%
Cultures	4	*Group B Streptococcus present	1	25.00%
Cultures		*Negative for GC	1	25.00%
Cultures		*No growth, GC selective media	2	50.00%
GC-15		* Denotes Correct Response		
AG: All Methods	3			
		*Negative for GC	2	66.67%
		*No Gram Negative diplococci present	1	33.33%
Cultures	3	*Negative for GC	2	66.67%
Cultures		*No Gram Negative diplococci present	1	33.33%
630 - Gastric Occult Blood & pH	Labs	Response	Number	Percent
		Gastric Occult Blood		
GOB-5		* Denotes Correct Response		
AG: All Methods	1			
		*Positive	1	100.00%
Beckman Coulter Gastrocult	1	*Positive	1	100.00%
GOB-6		* Denotes Correct Response		
AG: All Methods	1			
		*Negative	1	100.00%
Beckman Coulter Gastrocult	1	*Negative	1	100.00%
668 - Gram Stain	Labs	Response	Number	Percent
		Gram Stain		
GS-11		* Denotes Correct Response		
AG: All Methods	7			
		*Gram Negative	7	100.00%
Gram Stain	7	*Gram Negative	7	100.00%
GS-12		* Denotes Correct Response		
AG: All Methods	7			
		*Gram Positive	7	100.00%
Gram Stain	7	*Gram Positive	7	100.00%
GS-13		* Denotes Correct Response		
AG: All Methods	7			
		*Gram Positive	7	100.00%
Gram Stain	7	*Gram Positive	7	100.00%
GS-14		* Denotes Correct Response		
AG: All Methods	7			
		*Gram Positive	6	85.71%
		Gram Negative	1	14.29%
Gram Stain	7	*Gram Positive	6	85.71%
Gram Stain		Gram Negative	1	14.29%
GS-15		* Denotes Correct Response		
AG: All Methods	7			
		*Gram Negative	6	85.71%
		Gram Positive	1	14.29%
Gram Stain	7	*Gram Negative	6	85.71%
Gram Stain		Gram Positive	1	14.29%
		Gram Stain Morph		
GS-11		* Denotes Correct Response		
AG: All Methods	6			
		*Rods/Bacilli	6	100.00%



Gram Stain	6	*Rods/Bacilli	6	100.00%
GS-12		* Denotes Correct Response		
AG: All Methods	6			
		*Cocci	5	83.33%
		*Cocci in Chains	1	16.67%
Gram Stain	6	*Cocci	5	83.33%
Gram Stain		*Cocci in Chains	1	16.67%
GS-13		* Denotes Correct Response		
AG: All Methods	6			
		*Cocci	4	66.67%
		*Cocci in Chains	1	16.67%
		Cocco-Bacilli	1	16.67%
Gram Stain	6	*Cocci	4	66.67%
Gram Stain		*Cocci in Chains	1	16.67%
Gram Stain		Cocco-Bacilli	1	16.67%
GS-14		* Denotes Correct Response		
AG: All Methods	6			
		*Cocci	4	66.67%
		*Cocci in Chains	2	33.33%
Gram Stain	6	*Cocci	4	66.67%
Gram Stain		*Cocci in Chains	2	33.33%
GS-15		* Denotes Correct Response		
AG: All Methods	6			
		*Cocco-Bacilli	2	33.33%
		*Rods/Bacilli	4	66.67%
Gram Stain	6	*Cocco-Bacilli	2	33.33%
Gram Stain		*Rods/Bacilli	4	66.67%
978 - Gastrointestinal Viral Detection	Labs	Response	Number	Percent
		Rotavirus - Ag Detect		
GVD-11		* Denotes Correct Response		
AG: All Methods	10			
		*Positive	9	90.00%
		Negative	1	10.00%
Meridian Immunocard Stat! Rotavirus	5	*Positive	4	80.00%
SAS Rota-Test	5	*Positive	5	100.00%
Meridian Immunocard Stat! Rotavirus	5	Negative	1	20.00%
GVD-12		* Denotes Correct Response		
AG: All Methods	10			
		*Negative	10	100.00%
Meridian Immunocard Stat! Rotavirus	5	*Negative	5	100.00%
SAS Rota-Test	5	*Negative	5	100.00%
GVD-13		* Denotes Correct Response		
AG: All Methods	10			
		*Negative	10	100.00%
Meridian Immunocard Stat! Rotavirus	5	*Negative	5	100.00%
SAS Rota-Test	5	*Negative	5	100.00%
GVD-14		* Denotes Correct Response		
AG: All Methods	10			
		*Positive	9	90.00%
		Negative	1	10.00%
Meridian Immunocard Stat! Rotavirus	5	*Positive	4	80.00%
SAS Rota-Test	5	*Positive	5	100.00%
Meridian Immunocard Stat! Rotavirus	5	Negative	1	20.00%
GVD-15		* Denotes Correct Response		
AG: All Methods	10			
		*Positive	8	80.00%
		Negative	2	20.00%
Meridian Immunocard Stat! Rotavirus	5	*Positive	4	80.00%
SAS Rota-Test	5	*Positive	4	80.00%
Meridian Immunocard Stat! Rotavirus	5	Negative	1	20.00%
SAS Rota-Test	5	Negative	1	20.00%
678 - Serum hCG	Labs	Response	Number	Percent
		HCG - Qual Serum		
HCG-11		* Denotes Correct Response		
AG: All Methods	74			
		*Positive	74	100.00%
Beck Coulter ICON 25	4	*Positive	4	100.00%
Fisher HC Sure-View hCG STAT	2	*Positive	2	100.00%



Sekisui Diagnostics OSOM	2	*Positive	2	100.00%
Henry Schein One Step + (hCG)	6	*Positive	6	100.00%
Immunostics Immuno/hCG	1	*Positive	1	100.00%
Inverness Signify HCG	1	*Positive	1	100.00%
LifeSign Status hCG	1	*Positive	1	100.00%
McKesson Medi-Lab hCG - Urine	1	*Positive	1	100.00%
Quidel QuickVue + hCG	2	*Positive	2	100.00%
Quidel QuickVue One-Step hCG Combo (U/S)	22	*Positive	22	100.00%
Select Diag hCG Combo	3	*Positive	3	100.00%
SP Brand (Cardinal) HCG Combo	17	*Positive	17	100.00%
Wampole Clearview	1	*Positive	1	100.00%
Polymedco Poly stat hCG Serum	1	*Positive	1	100.00%
Stanbio True 20 Plus	1	*Positive	1	100.00%
Consult Diag hCG Combo Cassette	8	*Positive	8	100.00%
Alere hCG Combo Cassette (25 mIU/mL)	1	*Positive	1	100.00%
HCG-12		* Denotes Correct Response		
AG: All Methods	74			
		*Negative	74	100.00%
Beck Coulter ICON 25	4	*Negative	4	100.00%
Fisher HC Sure-Vue hCG STAT	2	*Negative	2	100.00%
Sekisui Diagnostics OSOM	2	*Negative	2	100.00%
Henry Schein One Step + (hCG)	6	*Negative	6	100.00%
Immunostics Immuno/hCG	1	*Negative	1	100.00%
Inverness Signify HCG	1	*Negative	1	100.00%
LifeSign Status hCG	1	*Negative	1	100.00%
McKesson Medi-Lab hCG - Urine	1	*Negative	1	100.00%
Quidel QuickVue + hCG	2	*Negative	2	100.00%
Quidel QuickVue One-Step hCG Combo (U/S)	22	*Negative	22	100.00%
Select Diag hCG Combo	3	*Negative	3	100.00%
SP Brand (Cardinal) HCG Combo	17	*Negative	17	100.00%
Wampole Clearview	1	*Negative	1	100.00%
Polymedco Poly stat hCG Serum	1	*Negative	1	100.00%
Stanbio True 20 Plus	1	*Negative	1	100.00%
Consult Diag hCG Combo Cassette	8	*Negative	8	100.00%
Alere hCG Combo Cassette (25 mIU/mL)	1	*Negative	1	100.00%
HCG-13		* Denotes Correct Response		
AG: All Methods	74			
		*Negative	74	100.00%
Beck Coulter ICON 25	4	*Negative	4	100.00%
Fisher HC Sure-Vue hCG STAT	2	*Negative	2	100.00%
Sekisui Diagnostics OSOM	2	*Negative	2	100.00%
Henry Schein One Step + (hCG)	6	*Negative	6	100.00%
Immunostics Immuno/hCG	1	*Negative	1	100.00%
Inverness Signify HCG	1	*Negative	1	100.00%
LifeSign Status hCG	1	*Negative	1	100.00%
McKesson Medi-Lab hCG - Urine	1	*Negative	1	100.00%
Quidel QuickVue + hCG	2	*Negative	2	100.00%
Quidel QuickVue One-Step hCG Combo (U/S)	22	*Negative	22	100.00%
Select Diag hCG Combo	3	*Negative	3	100.00%
SP Brand (Cardinal) HCG Combo	17	*Negative	17	100.00%
Wampole Clearview	1	*Negative	1	100.00%
Polymedco Poly stat hCG Serum	1	*Negative	1	100.00%
Stanbio True 20 Plus	1	*Negative	1	100.00%
Consult Diag hCG Combo Cassette	8	*Negative	8	100.00%
Alere hCG Combo Cassette (25 mIU/mL)	1	*Negative	1	100.00%
HCG-14		* Denotes Correct Response		
AG: All Methods	74			
		*Positive	74	100.00%
Beck Coulter ICON 25	4	*Positive	4	100.00%
Fisher HC Sure-Vue hCG STAT	2	*Positive	2	100.00%
Sekisui Diagnostics OSOM	2	*Positive	2	100.00%
Henry Schein One Step + (hCG)	6	*Positive	6	100.00%
Immunostics Immuno/hCG	1	*Positive	1	100.00%
Inverness Signify HCG	1	*Positive	1	100.00%
LifeSign Status hCG	1	*Positive	1	100.00%
McKesson Medi-Lab hCG - Urine	1	*Positive	1	100.00%
Quidel QuickVue + hCG	2	*Positive	2	100.00%
Quidel QuickVue One-Step hCG Combo (U/S)	22	*Positive	22	100.00%
Select Diag hCG Combo	3	*Positive	3	100.00%



SP Brand (Cardinal) HCG Combo	17	*Positive	17	100.00%
Wampole Clearview	1	*Positive	1	100.00%
Polymedco Poly stat hCG Serum	1	*Positive	1	100.00%
Stanbio True 20 Plus	1	*Positive	1	100.00%
Consult Diag hCG Combo Cassette	8	*Positive	8	100.00%
Alere hCG Combo Cassette (25 mIU/mL)	1	*Positive	1	100.00%
HCG-15		* Denotes Correct Response		
AG: All Methods	74			
		*Positive	73	98.65%
		Negative	1	1.35%
Beck Coulter ICON 25	4	*Positive	4	100.00%
Fisher HC Sure-Vue hCG STAT	2	*Positive	2	100.00%
Sekisui Diagnostics OSOM	2	*Positive	2	100.00%
Henry Schein One Step + (hCG)	6	*Positive	5	83.33%
Immunostics Immuno/hCG	1	*Positive	1	100.00%
Inverness Signify HCG	1	*Positive	1	100.00%
LifeSign Status hCG	1	*Positive	1	100.00%
McKesson Medi-Lab hCG - Urine	1	*Positive	1	100.00%
Quidel QuickVue + hCG	2	*Positive	2	100.00%
Quidel QuickVue One-Step hCG Combo (U/S)	22	*Positive	22	100.00%
Select Diag hCG Combo	3	*Positive	3	100.00%
SP Brand (Cardinal) HCG Combo	17	*Positive	17	100.00%
Wampole Clearview	1	*Positive	1	100.00%
Polymedco Poly stat hCG Serum	1	*Positive	1	100.00%
Stanbio True 20 Plus	1	*Positive	1	100.00%
Consult Diag hCG Combo Cassette	8	*Positive	8	100.00%
Alere hCG Combo Cassette (25 mIU/mL)	1	*Positive	1	100.00%
Henry Schein One Step + (hCG)	6	Negative	1	16.67%
		HCG - Qual Serum (International)		
HCG-11		* Denotes Correct Response		
AG: All Methods	3			
		*Positive	3	100.00%
Henry Schein One Step + (hCG)	1	*Positive	1	100.00%
SAS Ultra HCG	2	*Positive	2	100.00%
HCG-12		* Denotes Correct Response		
AG: All Methods	3			
		*Negative	3	100.00%
Henry Schein One Step + (hCG)	1	*Negative	1	100.00%
SAS Ultra HCG	2	*Negative	2	100.00%
HCG-13		* Denotes Correct Response		
AG: All Methods	3			
		*Negative	3	100.00%
Henry Schein One Step + (hCG)	1	*Negative	1	100.00%
SAS Ultra HCG	2	*Negative	2	100.00%
HCG-14		* Denotes Correct Response		
AG: All Methods	3			
		*Positive	3	100.00%
Henry Schein One Step + (hCG)	1	*Positive	1	100.00%
SAS Ultra HCG	2	*Positive	2	100.00%
HCG-15		* Denotes Correct Response		
AG: All Methods	3			
		*Positive	3	100.00%
Henry Schein One Step + (hCG)	1	*Positive	1	100.00%
SAS Ultra HCG	2	*Positive	2	100.00%
933 - HIV Antibody/Antigen Combo (2 specimens):934 - HIV Antibody/Antigen Combo (5 specimens)	Labs	Response	Number	Percent
		HIV Antibody (Waived)		
HIV-11		* Denotes Correct Response		
AG: All Methods	83			
		*Positive	80	96.39%
		Negative	3	3.61%
Alere Clearview Complete HIV 1/2 (mod)	1	*Positive	1	100.00%
Alere Clearview HIV 1/2 STAT-PAK (WB)	42	*Positive	40	95.24%
Orasure OraQuick Anti-HIV	11	*Positive	11	100.00%
Trinity Biotech Uni-Gold Recombigen HIV	1	*Positive	1	100.00%
Alere Determine HIV Ag/Ab	12	*Positive	11	91.67%
Chembio HIV 1/2 STAT-PAK Dipstick	7	*Positive	7	100.00%
Chembio HIV 1/2 STAT-PAK Cassette	4	*Positive	4	100.00%



bioLytical INSTI HIV Antibody	5	*Positive	5	100.00%
Alere Clearview HIV 1/2 STAT PAK (WB)	42	Negative	2	4.76%
Alere Determine HIV Ag/Ab	12	Negative	1	8.33%
HIV-12		* Denotes Correct Response		
AG: All Methods	83			
		*Negative	80	96.39%
		Positive	3	3.61%
Alere Clearview HIV 1/2 STAT PAK (WB)	42	*Negative	41	97.62%
Orasure OraQuick Anti-HIV	11	*Negative	11	100.00%
Trinity Biotech Uni-Gold Recombigen HIV	1	*Negative	1	100.00%
Alere Determine HIV Ag/Ab	12	*Negative	11	91.67%
Chembio HIV 1/2 STAT-PAK Dipstick	7	*Negative	7	100.00%
Chembio HIV 1/2 STAT-PAK Cassette	4	*Negative	4	100.00%
bioLytical INSTI HIV Antibody	5	*Negative	5	100.00%
Alere Clearview Complete HIV 1/2 (mod)	1	Positive	1	100.00%
Alere Clearview HIV 1/2 STAT PAK (WB)	42	Positive	1	2.38%
Alere Determine HIV Ag/Ab	12	Positive	1	8.33%
		HIV Antibody (International)		
HIV-11		* Denotes Correct Response		
AG: All Methods	3			
		*Positive	3	100.00%
Alere Clearview HIV 1/2 STAT PAK (WB)	1	*Positive	1	100.00%
Intl OCD Vitros ECi/ECiQ	1	*Positive	1	100.00%
Intl Chembio HIV 1/2 STAT-PAK Dipstick	1	*Positive	1	100.00%
HIV-12		* Denotes Correct Response		
AG: All Methods	3			
		*Negative	3	100.00%
Alere Clearview HIV 1/2 STAT PAK (WB)	1	*Negative	1	100.00%
Intl OCD Vitros ECi/ECiQ	1	*Negative	1	100.00%
Intl Chembio HIV 1/2 STAT-PAK Dipstick	1	*Negative	1	100.00%
		HIV Antigen (p24) (Waived)		
HIV-11		* Denotes Correct Response		
AG: All Methods	12			
		*Negative	11	91.67%
		Positive	1	8.33%
Alere Determine HIV Ag/Ab	12	*Negative	11	91.67%
Alere Determine HIV Ag/Ab		Positive	1	8.33%
HIV-12		* Denotes Correct Response		
AG: All Methods	12			
		*Positive	11	91.67%
		Negative	1	8.33%
Alere Determine HIV Ag/Ab	12	*Positive	11	91.67%
Alere Determine HIV Ag/Ab		Negative	1	8.33%
934 - HIV Antibody/Antigen Combo (5 specimens)	Labs	Response	Number	Percent
		HIV Antibody (Waived)		
HIV-11		* Denotes Correct Response		
AG: All Methods	83			
		*Positive	80	96.39%
		Negative	3	3.61%
Alere Clearview Complete HIV 1/2 (mod)	1	*Positive	1	100.00%
Alere Clearview HIV 1/2 STAT PAK (WB)	42	*Positive	40	95.24%
Orasure OraQuick Anti-HIV	11	*Positive	11	100.00%
Trinity Biotech Uni-Gold Recombigen HIV	1	*Positive	1	100.00%
Alere Determine HIV Ag/Ab	12	*Positive	11	91.67%
Chembio HIV 1/2 STAT-PAK Dipstick	7	*Positive	7	100.00%
Chembio HIV 1/2 STAT-PAK Cassette	4	*Positive	4	100.00%
bioLytical INSTI HIV Antibody	5	*Positive	5	100.00%
Alere Clearview HIV 1/2 STAT PAK (WB)	42	Negative	2	4.76%
Alere Determine HIV Ag/Ab	12	Negative	1	8.33%
HIV-12		* Denotes Correct Response		
AG: All Methods	83			
		*Negative	80	96.39%
		Positive	3	3.61%
Alere Clearview HIV 1/2 STAT PAK (WB)	42	*Negative	41	97.62%
Orasure OraQuick Anti-HIV	11	*Negative	11	100.00%
Trinity Biotech Uni-Gold Recombigen HIV	1	*Negative	1	100.00%
Alere Determine HIV Ag/Ab	12	*Negative	11	91.67%
Chembio HIV 1/2 STAT-PAK Dipstick	7	*Negative	7	100.00%



Chembio HIV 1/2 STAT-PAK Cassette	4	*Negative	4	100.00%
bioLytical INSTI HIV Antibody	5	*Negative	5	100.00%
Alere Clearview Complete HIV 1/2 (mod)	1	Positive	1	100.00%
Alere Clearview HIV 1/2 STAT PAK (WB)	42	Positive	1	2.38%
Alere Determine HIV Ag/Ab	12	Positive	1	8.33%
HIV-13		* Denotes Correct Response		
AG: All Methods	9			
		*Negative	9	100.00%
Alere Clearview Complete HIV 1/2 (mod)	1	*Negative	1	100.00%
Alere Clearview HIV 1/2 STAT PAK (WB)	4	*Negative	4	100.00%
Orasure OraQuick Anti-HIV	1	*Negative	1	100.00%
Trinity Biotech Uni-Gold Recombigen HIV	1	*Negative	1	100.00%
Chembio HIV 1/2 STAT-PAK Cassette	1	*Negative	1	100.00%
bioLytical INSTI HIV Antibody	1	*Negative	1	100.00%
HIV-14		* Denotes Correct Response		
AG: All Methods	9			
		*Negative	9	100.00%
Alere Clearview Complete HIV 1/2 (mod)	1	*Negative	1	100.00%
Alere Clearview HIV 1/2 STAT PAK (WB)	4	*Negative	4	100.00%
Orasure OraQuick Anti-HIV	1	*Negative	1	100.00%
Trinity Biotech Uni-Gold Recombigen HIV	1	*Negative	1	100.00%
Chembio HIV 1/2 STAT-PAK Cassette	1	*Negative	1	100.00%
bioLytical INSTI HIV Antibody	1	*Negative	1	100.00%
HIV-15		* Denotes Correct Response		
AG: All Methods	9			
		*Positive	9	100.00%
Alere Clearview Complete HIV 1/2 (mod)	1	*Positive	1	100.00%
Alere Clearview HIV 1/2 STAT PAK (WB)	4	*Positive	4	100.00%
Orasure OraQuick Anti-HIV	1	*Positive	1	100.00%
Trinity Biotech Uni-Gold Recombigen HIV	1	*Positive	1	100.00%
Chembio HIV 1/2 STAT-PAK Cassette	1	*Positive	1	100.00%
bioLytical INSTI HIV Antibody	1	*Positive	1	100.00%
		HIV Antibody (International)		
HIV-11		* Denotes Correct Response		
AG: All Methods	3			
		*Positive	3	100.00%
Alere Clearview HIV 1/2 STAT PAK (WB)	1	*Positive	1	100.00%
Intl OCD Vitros ECI/ECiQ	1	*Positive	1	100.00%
Intl Chembio HIV 1/2 STAT-PAK Dipstick	1	*Positive	1	100.00%
HIV-12		* Denotes Correct Response		
AG: All Methods	3			
		*Negative	3	100.00%
Alere Clearview HIV 1/2 STAT PAK (WB)	1	*Negative	1	100.00%
Intl OCD Vitros ECI/ECiQ	1	*Negative	1	100.00%
Intl Chembio HIV 1/2 STAT-PAK Dipstick	1	*Negative	1	100.00%
HIV-13		* Denotes Correct Response		
AG: All Methods	1			
		*Negative	1	100.00%
Intl OCD Vitros ECI/ECiQ	1	*Negative	1	100.00%
HIV-14		* Denotes Correct Response		
AG: All Methods	1			
		*Negative	1	100.00%
Intl OCD Vitros ECI/ECiQ	1	*Negative	1	100.00%
HIV-15		* Denotes Correct Response		
AG: All Methods	1			
		*Positive	1	100.00%
Intl OCD Vitros ECI/ECiQ	1	*Positive	1	100.00%
		HIV Antibody		
HIV-11		* Denotes Correct Response		
AG: All Methods	2			
		*Positive	2	100.00%
OCD Vitros 5600	2	*Positive	2	100.00%
HIV-12		* Denotes Correct Response		
AG: All Methods	2			
		*Negative	1	50.00%
		Positive	1	50.00%
OCD Vitros 5600	2	*Negative	1	50.00%
OCD Vitros 5600		Positive	1	50.00%



HIV-13		* Denotes Correct Response		
AG: All Methods	2			
		*Negative	2	100.00%
OCD Vitros 5600	2	*Negative	2	100.00%
HIV-14		* Denotes Correct Response		
AG: All Methods	2			
		*Negative	2	100.00%
OCD Vitros 5600	2	*Negative	2	100.00%
HIV-15		* Denotes Correct Response		
AG: All Methods	2			
		*Positive	2	100.00%
OCD Vitros 5600	2	*Positive	2	100.00%
HIV Ab/P24 Ag Combo (International)				
HIV-11		* Denotes Correct Response		
AG: All Methods	1			
		*Positive	1	100.00%
Intl OCD Vitros ECI/ECiQ	1	*Positive	1	100.00%
HIV-12		* Denotes Correct Response		
AG: All Methods	2			
		*Positive	2	100.00%
Intl OCD Vitros ECI/ECiQ	2	*Positive	2	100.00%
HIV-13		* Denotes Correct Response		
AG: All Methods	2			
		*Negative	2	100.00%
Intl OCD Vitros ECI/ECiQ	2	*Negative	2	100.00%
HIV-14		* Denotes Correct Response		
AG: All Methods	2			
		*Negative	2	100.00%
Intl OCD Vitros ECI/ECiQ	2	*Negative	2	100.00%
HIV-15		* Denotes Correct Response		
AG: All Methods	2			
		*Positive	2	100.00%
Intl OCD Vitros ECI/ECiQ	2	*Positive	2	100.00%
HIV Antigen (p24) (Waived)				
HIV-11		* Denotes Correct Response		
AG: All Methods	12			
		*Negative	11	91.67%
		Positive	1	8.33%
Alere Determine HIV Ag/Ab	12	*Negative	11	91.67%
Alere Determine HIV Ag/Ab		Positive	1	8.33%
HIV-12		* Denotes Correct Response		
AG: All Methods	12			
		*Positive	11	91.67%
		Negative	1	8.33%
Alere Determine HIV Ag/Ab	12	*Positive	11	91.67%
Alere Determine HIV Ag/Ab		Negative	1	8.33%
644 - H. pylori	Labs	Response	Number	Percent
H. pylori Ab Detect				
HP-3		* Denotes Correct Response		
AG: All Methods	63			
		*Positive	61	96.83%
		Negative	2	3.17%
Beck Coulter FlexSure	3	*Positive	3	100.00%
Beck Coulter ICON HP	2	*Positive	2	100.00%
McKesson Medi-Lab Perf H. pylori	6	*Positive	5	83.33%
Henry Schein One Step+ (H. pylori)	11	*Positive	11	100.00%
Select Diag H. pylori Cassette	2	*Positive	2	100.00%
Quidel QuickVue H. pylori gII (S/P)	3	*Positive	3	100.00%
Quidel QuickVue H. pylori gII (WB)	7	*Positive	7	100.00%
SP Brand (Cardinal) HP Rapid Test	8	*Positive	8	100.00%
Alere Clearview II S/P (mod)	2	*Positive	2	100.00%
Alere Clearview II WB (waived)	3	*Positive	3	100.00%
LifeSign Status H. pylori	2	*Positive	2	100.00%
Consult Diag H pylori	12	*Positive	11	91.67%
Jant Accutest H pylori	1	*Positive	1	100.00%
Medline H. pylori Cassette	1	*Positive	1	100.00%
McKesson Medi-Lab Perf H. pylori	6	Negative	1	16.67%
Consult Diag H pylori	12	Negative	1	8.33%



929 - Genital Viral Detection I	Labs	Response	Number	Percent
		HSV 1		
HVD-11		* Denotes Correct Response		
AG: All Methods	2			
		*Positive	2	100.00%
Hologic Aptima HSV 1 & 2	2	*Positive	2	100.00%
HVD-12		* Denotes Correct Response		
AG: All Methods	2			
		*Positive	1	50.00%
		Negative	1	50.00%
Hologic Aptima HSV 1 & 2	2	*Positive	1	50.00%
Hologic Aptima HSV 1 & 2		Negative	1	50.00%
HVD-13		* Denotes Correct Response		
AG: All Methods	2			
		*Negative	2	100.00%
Hologic Aptima HSV 1 & 2	2	*Negative	2	100.00%
HVD-14		* Denotes Correct Response		
AG: All Methods	2			
		*Negative	2	100.00%
Hologic Aptima HSV 1 & 2	2	*Negative	2	100.00%
HVD-15		* Denotes Correct Response		
AG: All Methods	2			
		*Negative	2	100.00%
Hologic Aptima HSV 1 & 2	2	*Negative	2	100.00%
		HSV 2		
HVD-11		* Denotes Correct Response		
AG: All Methods	1			
		*Negative	1	100.00%
Hologic Aptima HSV 1 & 2	1	*Negative	1	100.00%
HVD-12		* Denotes Correct Response		
AG: All Methods	1			
		*Negative	1	100.00%
Hologic Aptima HSV 1 & 2	1	*Negative	1	100.00%
HVD-13		* Denotes Correct Response		
AG: All Methods	1			
		*Positive	1	100.00%
Hologic Aptima HSV 1 & 2	1	*Positive	1	100.00%
HVD-14		* Denotes Correct Response		
AG: All Methods	1			
		*Negative	1	100.00%
Hologic Aptima HSV 1 & 2	1	*Negative	1	100.00%
HVD-15		* Denotes Correct Response		
AG: All Methods	1			
		*Positive	1	100.00%
Hologic Aptima HSV 1 & 2	1	*Positive	1	100.00%
603 - Blood Cell Identification:610 - Regular Hematology	Labs	Response	Number	Percent
		Blood Cell Identification		
ID-13		* Denotes Correct Response		
AG: All Methods	62			
		*Segmented or band neutrophil	62	100.00%
Blood Cell Microscopy	62	*Segmented or band neutrophil	62	100.00%
ID-14		* Denotes Correct Response		
AG: All Methods	62			
		*Red Blood Cell, normal (RBC) (Erythrocyte)	59	95.16%
		Spherocyte	3	4.84%
Blood Cell Microscopy	62	*Red Blood Cell, normal (RBC) (Erythrocyte)	59	95.16%
Blood Cell Microscopy		Spherocyte	3	4.84%
ID-15		* Denotes Correct Response		
AG: All Methods	62			
		*Identification Unknown - Referred	3	4.84%
		*Monocyte	57	91.94%
		Lymphocyte, reactive (atypical, variant, Downey)	2	3.23%
Blood Cell Microscopy	62	*Identification Unknown - Referred	3	4.84%
Blood Cell Microscopy		*Monocyte	57	91.94%
Blood Cell Microscopy		Lymphocyte, reactive (atypical, variant, Downey)	2	3.23%
ID-16		* Denotes Correct Response		
AG: All Methods	62			



		*Identification Unknown - Referred	1	1.61%
		*Immature/abnormal cell - Referred	3	4.84%
		*Lymphocyte, reactive (atypical, variant, Downey)	57	91.94%
		Monoblast	1	1.61%
Blood Cell Microscopy	62	*Identification Unknown - Referred	1	1.61%
Blood Cell Microscopy		*Immature/abnormal cell - Referred	3	4.84%
Blood Cell Microscopy		*Lymphocyte, reactive (atypical, variant, Downey)	57	91.94%
Blood Cell Microscopy		Monoblast	1	1.61%
ID-17		* Denotes Correct Response		
AG: All Methods	62			
		*Lymphocyte	62	100.00%
Blood Cell Microscopy	62	*Lymphocyte	62	100.00%
ID-18		* Denotes Correct Response		
AG: All Methods	62			
		*Platelet satellitism	59	95.16%
		Immature/abnormal cell - Referred	1	1.61%
		Other, not listed	1	1.61%
		Platelet clumping	1	1.61%
Blood Cell Microscopy	62	*Platelet satellitism	59	95.16%
Blood Cell Microscopy		Immature/abnormal cell - Referred	1	1.61%
Blood Cell Microscopy		Other, not listed	1	1.61%
Blood Cell Microscopy		Platelet clumping	1	1.61%
		Blood Cell Identification (International)		
ID-13		* Denotes Correct Response		
AG: All Methods	2			
		*Segmented or band neutrophil	1	50.00%
		Degenerated neutrophil	1	50.00%
Intl Blood Cell Microscopy	2	*Segmented or band neutrophil	1	50.00%
Intl Blood Cell Microscopy		Degenerated neutrophil	1	50.00%
ID-14		* Denotes Correct Response		
AG: All Methods	2			
		*Red Blood Cell, normal (RBC) (Erythrocyte)	2	100.00%
Intl Blood Cell Microscopy	2	*Red Blood Cell, normal (RBC) (Erythrocyte)	2	100.00%
ID-15		* Denotes Correct Response		
AG: All Methods	2			
		*Monocyte	2	100.00%
Intl Blood Cell Microscopy	2	*Monocyte	2	100.00%
ID-16		* Denotes Correct Response		
AG: All Methods	2			
		*Lymphocyte, reactive (atypical, variant, Downey)	2	100.00%
Intl Blood Cell Microscopy	2	*Lymphocyte, reactive (atypical, variant, Downey)	2	100.00%
ID-17		* Denotes Correct Response		
AG: All Methods	2			
		*Lymphocyte	2	100.00%
Intl Blood Cell Microscopy	2	*Lymphocyte	2	100.00%
ID-18		* Denotes Correct Response		
AG: All Methods	2			
		*Platelet satellitism	2	100.00%
Intl Blood Cell Microscopy	2	*Platelet satellitism	2	100.00%
609 - Waived/PPM Tests I;675 - Infectious Mononucleosis;753 - Infectious Mononucleosis (Waived) (1 specimen);754 - Infectious Mono/RF Combo	Labs	Response	Number	Percent
		Infectious Mono (Waived)		
IM-11		* Denotes Correct Response		
AG: All Methods	469			
		*Positive	464	98.93%
		Negative	5	1.07%
Alere Acceava (S)	1	*Positive	1	100.00%
Alere Acceava (WB)	53	*Positive	51	96.23%
Sekisui Diagnostics OSOM Mono (WB)	184	*Positive	182	98.91%
Henry Schein One-Step+ Mono	50	*Positive	50	100.00%
Instant Technologies i-Screen (WB)	1	*Positive	1	100.00%
Inverness Signify Mono	3	*Positive	3	100.00%
LifeSign Status Inf Mono (WB)	31	*Positive	31	100.00%
McKesson Medi-Lab Perf - Inf Mono	8	*Positive	8	100.00%
Polymedco Poly stat Mono - (WB)	6	*Positive	6	100.00%
Quidel QuickVue+ (WB)	9	*Positive	9	100.00%
SP Brand (Cardinal) Rapid Test	30	*Positive	30	100.00%



Alere Clearview Mono Cassette	18	*Positive	18	100.00%
Pulse Scientific	2	*Positive	2	100.00%
Diag Test Group Clarity Mono	2	*Positive	2	100.00%
Fisher HC Sure-View Mono	6	*Positive	6	100.00%
Immunostics Inc. ImmunoDetector Mono (WB)	2	*Positive	2	100.00%
Meridian Immunocard STAT! Mono (WB)	5	*Positive	4	80.00%
NDC ProAdvantage Mono	2	*Positive	2	100.00%
Consult Diag Mono Cassette (WB)	41	*Positive	41	100.00%
Alere Acceava (WB) Mono Cassette	15	*Positive	15	100.00%
Alere Acceava (WB)	53	Negative	2	3.77%
Sekisui Diagnostics OSOM Mono (WB)	184	Negative	2	1.09%
Meridian Immunocard STAT! Mono (WB)	5	Negative	1	20.00%
		Infectious Mononucleosis		
IM-11		* Denotes Correct Response		
AG: All Methods	45			
		*Positive	45	100.00%
Beck Coulter ICON (WB)	1	*Positive	1	100.00%
Alere Acceava (WB)	1	*Positive	1	100.00%
Sekisui Diagnostics OSOM Mono (WB)	9	*Positive	9	100.00%
LifeSign Status Inf Mono (WB)	1	*Positive	1	100.00%
Lifesign Unistep	1	*Positive	1	100.00%
McKesson Medi-Lab Perf - Inf Mono	2	*Positive	2	100.00%
Polymedco Poly stat Mono (S/P)	2	*Positive	2	100.00%
Quidel QuickVue+ (S/P)	2	*Positive	2	100.00%
Quidel QuickVue+ (WB)	8	*Positive	8	100.00%
Seradyn Color Slide II	3	*Positive	3	100.00%
SP Brand (Cardinal) Rapid Test	1	*Positive	1	100.00%
Alere Clearview Mono Cassette	8	*Positive	8	100.00%
Sekisui Diagnostics OSOM Mono Non-Waived	1	*Positive	1	100.00%
Consult Diag Mono Cassette (WB)	3	*Positive	3	100.00%
Consult Diag Mono Cassette (Serum)	1	*Positive	1	100.00%
LifeSign Status Inf Mono (S/P)	1	*Positive	1	100.00%
675 - Infectious Mononucleosis;754 - Infectious Mono/RF Combo	Labs	Response	Number	Percent
		Infectious Mononucleosis		
IM-11		* Denotes Correct Response		
AG: All Methods	45			
		*Positive	45	100.00%
Beck Coulter ICON (WB)	1	*Positive	1	100.00%
Alere Acceava (WB)	1	*Positive	1	100.00%
Sekisui Diagnostics OSOM Mono (WB)	9	*Positive	9	100.00%
LifeSign Status Inf Mono (WB)	1	*Positive	1	100.00%
Lifesign Unistep	1	*Positive	1	100.00%
McKesson Medi-Lab Perf - Inf Mono	2	*Positive	2	100.00%
Polymedco Poly stat Mono (S/P)	2	*Positive	2	100.00%
Quidel QuickVue+ (S/P)	2	*Positive	2	100.00%
Quidel QuickVue+ (WB)	8	*Positive	8	100.00%
Seradyn Color Slide II	3	*Positive	3	100.00%
SP Brand (Cardinal) Rapid Test	1	*Positive	1	100.00%
Alere Clearview Mono Cassette	8	*Positive	8	100.00%
Sekisui Diagnostics OSOM Mono Non-Waived	1	*Positive	1	100.00%
Consult Diag Mono Cassette (WB)	3	*Positive	3	100.00%
Consult Diag Mono Cassette (Serum)	1	*Positive	1	100.00%
LifeSign Status Inf Mono (S/P)	1	*Positive	1	100.00%
IM-12		* Denotes Correct Response		
AG: All Methods	45			
		*Negative	45	100.00%
Beck Coulter ICON (WB)	1	*Negative	1	100.00%
Alere Acceava (WB)	1	*Negative	1	100.00%
Sekisui Diagnostics OSOM Mono (WB)	9	*Negative	9	100.00%
LifeSign Status Inf Mono (WB)	1	*Negative	1	100.00%
Lifesign Unistep	1	*Negative	1	100.00%
McKesson Medi-Lab Perf - Inf Mono	2	*Negative	2	100.00%
Polymedco Poly stat Mono (S/P)	2	*Negative	2	100.00%
Quidel QuickVue+ (S/P)	2	*Negative	2	100.00%
Quidel QuickVue+ (WB)	8	*Negative	8	100.00%
Seradyn Color Slide II	3	*Negative	3	100.00%
SP Brand (Cardinal) Rapid Test	1	*Negative	1	100.00%
Alere Clearview Mono Cassette	8	*Negative	8	100.00%



Sekisui Diagnostics OSOM Mono Non-Waived	1	*Negative	1	100.00%
Consult Diag Mono Cassette (WB)	3	*Negative	3	100.00%
Consult Diag Mono Cassette (Serum)	1	*Negative	1	100.00%
LifeSign Status Inf Mono (S/P)	1	*Negative	1	100.00%
IM-13		* Denotes Correct Response		
AG: All Methods	45			
		*Positive	45	100.00%
Beck Coulter ICON (WB)	1	*Positive	1	100.00%
Alere Acceava (WB)	1	*Positive	1	100.00%
Sekisui Diagnostics OSOM Mono (WB)	9	*Positive	9	100.00%
LifeSign Status Inf Mono (WB)	1	*Positive	1	100.00%
Lifesign Unistep	1	*Positive	1	100.00%
McKesson Medi-Lab Perf - Inf Mono	2	*Positive	2	100.00%
Polymedco Poly stat Mono (S/P)	2	*Positive	2	100.00%
Quidel QuickVue+ (S/P)	2	*Positive	2	100.00%
Quidel QuickVue+ (WB)	8	*Positive	8	100.00%
Seradyn Color Slide II	3	*Positive	3	100.00%
SP Brand (Cardinal) Rapid Test	1	*Positive	1	100.00%
Alere Clearview Mono Cassette	8	*Positive	8	100.00%
Sekisui Diagnostics OSOM Mono Non-Waived	1	*Positive	1	100.00%
Consult Diag Mono Cassette (WB)	3	*Positive	3	100.00%
Consult Diag Mono Cassette (Serum)	1	*Positive	1	100.00%
LifeSign Status Inf Mono (S/P)	1	*Positive	1	100.00%
IM-14		* Denotes Correct Response		
AG: All Methods	45			
		*Negative	45	100.00%
Beck Coulter ICON (WB)	1	*Negative	1	100.00%
Alere Acceava (WB)	1	*Negative	1	100.00%
Sekisui Diagnostics OSOM Mono (WB)	9	*Negative	9	100.00%
LifeSign Status Inf Mono (WB)	1	*Negative	1	100.00%
Lifesign Unistep	1	*Negative	1	100.00%
McKesson Medi-Lab Perf - Inf Mono	2	*Negative	2	100.00%
Polymedco Poly stat Mono (S/P)	2	*Negative	2	100.00%
Quidel QuickVue+ (S/P)	2	*Negative	2	100.00%
Quidel QuickVue+ (WB)	8	*Negative	8	100.00%
Seradyn Color Slide II	3	*Negative	3	100.00%
SP Brand (Cardinal) Rapid Test	1	*Negative	1	100.00%
Alere Clearview Mono Cassette	8	*Negative	8	100.00%
Sekisui Diagnostics OSOM Mono Non-Waived	1	*Negative	1	100.00%
Consult Diag Mono Cassette (WB)	3	*Negative	3	100.00%
Consult Diag Mono Cassette (Serum)	1	*Negative	1	100.00%
LifeSign Status Inf Mono (S/P)	1	*Negative	1	100.00%
IM-15		* Denotes Correct Response		
AG: All Methods	45			
		*Negative	45	100.00%
Beck Coulter ICON (WB)	1	*Negative	1	100.00%
Alere Acceava (WB)	1	*Negative	1	100.00%
Sekisui Diagnostics OSOM Mono (WB)	9	*Negative	9	100.00%
LifeSign Status Inf Mono (WB)	1	*Negative	1	100.00%
Lifesign Unistep	1	*Negative	1	100.00%
McKesson Medi-Lab Perf - Inf Mono	2	*Negative	2	100.00%
Polymedco Poly stat Mono (S/P)	2	*Negative	2	100.00%
Quidel QuickVue+ (S/P)	2	*Negative	2	100.00%
Quidel QuickVue+ (WB)	8	*Negative	8	100.00%
Seradyn Color Slide II	3	*Negative	3	100.00%
SP Brand (Cardinal) Rapid Test	1	*Negative	1	100.00%
Alere Clearview Mono Cassette	8	*Negative	8	100.00%
Sekisui Diagnostics OSOM Mono Non-Waived	1	*Negative	1	100.00%
Consult Diag Mono Cassette (WB)	3	*Negative	3	100.00%
Consult Diag Mono Cassette (Serum)	1	*Negative	1	100.00%
LifeSign Status Inf Mono (S/P)	1	*Negative	1	100.00%
		Infectious Mono (Waived)		
IM-11		* Denotes Correct Response		
AG: All Methods	469			
		*Positive	464	98.93%
		Negative	5	1.07%
Alere Acceava (S)	1	*Positive	1	100.00%
Alere Acceava (WB)	53	*Positive	51	96.23%
Sekisui Diagnostics OSOM Mono (WB)	184	*Positive	182	98.91%



Henry Schein One-Step+ Mono	50	*Positive	50	100.00%
Instant Technologies i-Screen (WB)	1	*Positive	1	100.00%
Inverness Signify Mono	3	*Positive	3	100.00%
LifeSign Status Inf Mono (WB)	31	*Positive	31	100.00%
McKesson Medi-Lab Perf - Inf Mono	8	*Positive	8	100.00%
Polymedco Poly stat Mono - (WB)	6	*Positive	6	100.00%
Quidel QuickVue+ (WB)	9	*Positive	9	100.00%
SP Brand (Cardinal) Rapid Test	30	*Positive	30	100.00%
Alere Clearview Mono Cassette	18	*Positive	18	100.00%
Pulse Scientific	2	*Positive	2	100.00%
Diag Test Group Clarity Mono	2	*Positive	2	100.00%
Fisher HC Sure-View Mono	6	*Positive	6	100.00%
Immunostics Inc. ImmunoDetector Mono (WB)	2	*Positive	2	100.00%
Meridian Immunocard STAT! Mono (WB)	5	*Positive	4	80.00%
NDC ProAdvantage Mono	2	*Positive	2	100.00%
Consult Diag Mono Cassette (WB)	41	*Positive	41	100.00%
Alere Acceava (WB) Mono Cassette	15	*Positive	15	100.00%
Alere Acceava (WB)	53	Negative	2	3.77%
Sekisui Diagnostics OSOM Mono (WB)	184	Negative	2	1.09%
Meridian Immunocard STAT! Mono (WB)	5	Negative	1	20.00%
		Infectious Mononucleosis (International)		
IM-11		* Denotes Correct Response		
AG: All Methods	2			
		*Positive	2	100.00%
Alere Clearview Mono Cassette	1	*Positive	1	100.00%
Fisher HC Sure-View Mono	1	*Positive	1	100.00%
IM-12		* Denotes Correct Response		
AG: All Methods	2			
		*Negative	2	100.00%
Alere Clearview Mono Cassette	1	*Negative	1	100.00%
Fisher HC Sure-View Mono	1	*Negative	1	100.00%
IM-13		* Denotes Correct Response		
AG: All Methods	2			
		*Positive	2	100.00%
Alere Clearview Mono Cassette	1	*Positive	1	100.00%
Fisher HC Sure-View Mono	1	*Positive	1	100.00%
IM-14		* Denotes Correct Response		
AG: All Methods	2			
		*Negative	2	100.00%
Alere Clearview Mono Cassette	1	*Negative	1	100.00%
Fisher HC Sure-View Mono	1	*Negative	1	100.00%
IM-15		* Denotes Correct Response		
AG: All Methods	2			
		*Negative	2	100.00%
Alere Clearview Mono Cassette	1	*Negative	1	100.00%
Fisher HC Sure-View Mono	1	*Negative	1	100.00%
637 - KOH - Glass Slides (Skin Hair Nail)	Labs	Response	Number	Percent
		KOH Preparation		
KOH-5		* Denotes Correct Response		
AG: All Methods	20			
		*Fungal elements/yeast observed	19	95.00%
		No fungal elements/yeast observed	1	5.00%
KOH Glass Slide	20	*Fungal elements/yeast observed	19	95.00%
KOH Glass Slide		No fungal elements/yeast observed	1	5.00%
KOH-6		* Denotes Correct Response		
AG: All Methods	20			
		*No fungal elements/yeast observed	16	80.00%
		Fungal elements/yeast observed	4	20.00%
KOH Glass Slide	20	*No fungal elements/yeast observed	16	80.00%
KOH Glass Slide		Fungal elements/yeast observed	4	20.00%
654 - Lyme Serology	Labs	Response	Number	Percent
		Lyme Antibody IgG		
LY-3		* Denotes Correct Response		
PG: Quidel Sofia/Sofia 2 Lyme FIA	4			
		*Positive	4	100.00%
		Lyme Antibody IgM		
LY-3		* Denotes Correct Response		
AG: All Methods	4			



		*Positive	4	100.00%
Quidel Sofia/Sofia 2 Lyme FIA	4	*Positive	4	100.00%
768 - MRSA Molecular (2 specimens);690 - MRSA (Molecular) (5 specimens)	Labs	Response	Number	Percent
		MRSA (Molecular)		
MRSA-11		* Denotes Correct Response		
AG: All Methods	1			
		*Positive	1	100.00%
Cepheid GeneXpert - Xpert MRSA	1	*Positive	1	100.00%
MRSA-12		* Denotes Correct Response		
AG: All Methods	1			
		*Negative	1	100.00%
Cepheid GeneXpert - Xpert MRSA	1	*Negative	1	100.00%
758 - Mycoplasma Antibody	Labs	Response	Number	Percent
		Mycoplasma Antibody		
MY-3		* Denotes Correct Response		
AG: All Methods	2			
		*Positive	2	100.00%
Meridian ImmunoCard	2	*Positive	2	100.00%
791 - Fecal Occult Blood Add-On;601 - Waived/PPM Tests II;609 - Waived/PPM Tests I;647 - Urinalysis/hCG/Fecal Occult Blood	Labs	Response	Number	Percent
		Fecal Occult Blood		
OC-3		* Denotes Correct Response		
PG: Occult Blood / Guaiac methods	666			
		*Positive	641	96.25%
		Negative	25	3.75%
OC-3		* Denotes Correct Response		
PG: Occult Blood / Immunochemical methods	113			
		*Positive	110	97.35%
		Negative	3	2.65%
950 - Post Vasectomy Sperm Analysis	Labs	Response	Number	Percent
		Sperm Count (Post Vasectomy)		
PVAS-5		* Denotes Correct Response		
AG: All Methods	19			
		*Positive for Sperm	18	94.74%
		Negative for Sperm	1	5.26%
Hemocytometer	8	*Positive for Sperm	8	100.00%
Post Vasectomy Microscopy	11	*Positive for Sperm	10	90.91%
Post Vasectomy Microscopy		Negative for Sperm	1	9.09%
PVAS-6		* Denotes Correct Response		
AG: All Methods	19			
		*Negative for Sperm	18	94.74%
		Positive for Sperm	1	5.26%
Hemocytometer	8	*Negative for Sperm	8	100.00%
Post Vasectomy Microscopy	11	*Negative for Sperm	10	90.91%
Post Vasectomy Microscopy		Positive for Sperm	1	9.09%
754 - Infectious Mono/RF Combo;653 - Rheumatoid Factor	Labs	Response	Number	Percent
		Rheumatoid Factor		
RF-11		* Denotes Correct Response		
AG: All Methods	11			
		*Positive	11	100.00%
ASI	2	*Positive	2	100.00%
Siemens RapiTex	1	*Positive	1	100.00%
Seradyn Color Slide	1	*Positive	1	100.00%
Stanbio	1	*Positive	1	100.00%
Wampole ColorCard RF	2	*Positive	2	100.00%
Wampole Rheumatex	2	*Positive	2	100.00%
Fisher HC Sure-Vue RF	1	*Positive	1	100.00%
Biokit Rhemajet RF	1	*Positive	1	100.00%
RF-12		* Denotes Correct Response		
AG: All Methods	11			
		*N/A (Negative)	9	81.82%
		Positive	2	18.18%
ASI	2	*N/A (Negative)	1	50.00%
Siemens RapiTex	1	*N/A (Negative)	1	100.00%
Seradyn Color Slide	1	*N/A (Negative)	1	100.00%



Stanbio	1	*N/A (Negative)	1	100.00%
Wampole ColorCard RF	2	*N/A (Negative)	2	100.00%
Wampole Rheumatex	2	*N/A (Negative)	1	50.00%
Fisher HC Sure-Vue RF	1	*N/A (Negative)	1	100.00%
Biokit Rhemajet RF	1	*N/A (Negative)	1	100.00%
ASI	2	Positive	1	50.00%
Wampole Rheumatex	2	Positive	1	50.00%
RF-13		* Denotes Correct Response		
AG: All Methods	11			
		*N/A (Negative)	9	81.82%
		Positive	2	18.18%
ASI	2	*N/A (Negative)	1	50.00%
Siemens RapiTex	1	*N/A (Negative)	1	100.00%
Seradyn Color Slide	1	*N/A (Negative)	1	100.00%
Stanbio	1	*N/A (Negative)	1	100.00%
Wampole ColorCard RF	2	*N/A (Negative)	2	100.00%
Wampole Rheumatex	2	*N/A (Negative)	1	50.00%
Fisher HC Sure-Vue RF	1	*N/A (Negative)	1	100.00%
Biokit Rhemajet RF	1	*N/A (Negative)	1	100.00%
ASI	2	Positive	1	50.00%
Wampole Rheumatex	2	Positive	1	50.00%
RF-14		* Denotes Correct Response		
AG: All Methods	11			
		*Positive	11	100.00%
ASI	2	*Positive	2	100.00%
Siemens RapiTex	1	*Positive	1	100.00%
Seradyn Color Slide	1	*Positive	1	100.00%
Stanbio	1	*Positive	1	100.00%
Wampole ColorCard RF	2	*Positive	2	100.00%
Wampole Rheumatex	2	*Positive	2	100.00%
Fisher HC Sure-Vue RF	1	*Positive	1	100.00%
Biokit Rhemajet RF	1	*Positive	1	100.00%
RF-15		* Denotes Correct Response		
AG: All Methods	11			
		*N/A (Negative)	9	81.82%
		Positive	2	18.18%
ASI	2	*N/A (Negative)	1	50.00%
Siemens RapiTex	1	*N/A (Negative)	1	100.00%
Seradyn Color Slide	1	*N/A (Negative)	1	100.00%
Stanbio	1	*N/A (Negative)	1	100.00%
Wampole ColorCard RF	2	*N/A (Negative)	2	100.00%
Wampole Rheumatex	2	*N/A (Negative)	1	50.00%
Fisher HC Sure-Vue RF	1	*N/A (Negative)	1	100.00%
Biokit Rhemajet RF	1	*N/A (Negative)	1	100.00%
ASI	2	Positive	1	50.00%
Wampole Rheumatex	2	Positive	1	50.00%
		Rheumatoid Factor (International)		
RF-11		* Denotes Correct Response		
AG: All Methods	3			
		*Positive	3	100.00%
Fisher HC Sure-Vue RF	2	*Positive	2	100.00%
Immunostics Rheumacol/LEX	1	*Positive	1	100.00%
RF-12		* Denotes Correct Response		
AG: All Methods	3			
		*Negative	3	100.00%
Fisher HC Sure-Vue RF	2	*Negative	2	100.00%
Immunostics Rheumacol/LEX	1	*Negative	1	100.00%
RF-13		* Denotes Correct Response		
AG: All Methods	3			
		*Negative	3	100.00%
Fisher HC Sure-Vue RF	2	*Negative	2	100.00%
Immunostics Rheumacol/LEX	1	*Negative	1	100.00%
RF-14		* Denotes Correct Response		
AG: All Methods	3			
		*Positive	3	100.00%
Fisher HC Sure-Vue RF	2	*Positive	2	100.00%
Immunostics Rheumacol/LEX	1	*Positive	1	100.00%
RF-15		* Denotes Correct Response		



AG: All Methods	3			
		*Negative	3	100.00%
Fisher HC Sure-Vue RF	2	*Negative	2	100.00%
Immunostics Rheumacol/LEX	1	*Negative	1	100.00%
		Rheumatoid Factor (Semi Quant)		
RF-11		* Denotes Correct Response		
AG: All Methods	1			
		*80 IU/mL	1	100.00%
Fisher HC Sure-Vue RF (Semi Quant)	1	*80 IU/mL	1	100.00%
RF-12		* Denotes Correct Response		
AG: All Methods	1			
		*Negative	1	100.00%
Fisher HC Sure-Vue RF (Semi Quant)	1	*Negative	1	100.00%
RF-13		* Denotes Correct Response		
AG: All Methods	1			
		*Negative	1	100.00%
Fisher HC Sure-Vue RF (Semi Quant)	1	*Negative	1	100.00%
RF-14		* Denotes Correct Response		
AG: All Methods	1			
		*320 IU/mL	1	100.00%
Fisher HC Sure-Vue RF (Semi Quant)	1	*320 IU/mL	1	100.00%
RF-15		* Denotes Correct Response		
AG: All Methods	1			
		*Negative	1	100.00%
Fisher HC Sure-Vue RF (Semi Quant)	1	*Negative	1	100.00%
645 - Fetal Membrane Rupture	Labs	Response	Number	Percent
		Fetal Membranes Rupture PAMG-1		
ROM-5		* Denotes Correct Response		
AG: All Methods	18			
		*Positive	18	100.00%
AmniSure ROM Test	17	*Positive	17	100.00%
Clinical Innovations ROM Plus	1	*Positive	1	100.00%
ROM-6		* Denotes Correct Response		
AG: All Methods	18			
		*Negative	18	100.00%
AmniSure ROM Test	17	*Negative	17	100.00%
Clinical Innovations ROM Plus	1	*Negative	1	100.00%
601 - Waived/PPM Tests II;609 - Waived/PPM Tests I;691 - Group A Strep Antigen (5 specimens);783 - Group A Strep Antigen/Add-on (2 specimens);665 - Group A Strep Antigen (Waived) (2 specimens)	Labs	Response	Number	Percent
		Strep A Ag Det. (Waived)		
RS-11		* Denotes Correct Response		
AG: All Methods	852			
		*Negative	842	98.83%
		Positive	10	1.17%
Abbott Signify Strep A - Cassette	1	*Negative	1	100.00%
Abbott Signify Strep A - Waived	6	*Negative	6	100.00%
Beck Coulter ICON DS Strep	1	*Negative	1	100.00%
Alere Acceava Strep A Dipstick	110	*Negative	109	99.09%
Fisher HC Sure-Vue Strep A (Waived)	8	*Negative	8	100.00%
Sekisui Diagnostics OSOM	27	*Negative	26	96.30%
Sekisui Diagnostics OSOM Ultra Strep	62	*Negative	61	98.39%
Henry Schein OneStep + Rapid Strep Cass	6	*Negative	6	100.00%
Henry Schein OneStep+ Rapid Strep Dip	82	*Negative	80	97.56%
Innovacon One Step Strep Ag	1	*Negative	1	100.00%
Instant Technologies i Strep	1	*Negative	1	100.00%
Inverness Signify Strep A Dipstick	8	*Negative	8	100.00%
McKesson Strep A Cas	3	*Negative	3	100.00%
McKesson Strep A Dip	35	*Negative	35	100.00%
McKesson Strep A Test - Twist	1	*Negative	1	100.00%
Permaxim Rediscreen Strep A Rapid	1	*Negative	1	100.00%
Polymedco Poly Stat Strep A (W)	1	*Negative	1	100.00%
PSS Rapid Test Strip Strep A	2	*Negative	2	100.00%
Quidel QuickVue Dipstick Strep	85	*Negative	84	98.82%
Quidel QuickVue In-Line - 25 Test	67	*Negative	67	100.00%
Quidel QuickVue+ Strep A	24	*Negative	24	100.00%
Select Diag Strep A Dip	3	*Negative	3	100.00%
SP (Cardinal) Rapid Strep Cassette	4	*Negative	4	100.00%



SP Brand (Cardinal) Rapid Strep A Dipstick	55	*Negative	55	100.00%
Alere Clearview Exact II Dip	37	*Negative	35	94.59%
Sekisui Diagnostics OSOM Ultra	7	*Negative	7	100.00%
Consult Diag Strep A Dipstick	103	*Negative	101	98.06%
Consult Diag Strep A Twist	1	*Negative	1	100.00%
MooreBrand Strep A Dipstick	2	*Negative	2	100.00%
NDC ProAdvantage Strep Strips	1	*Negative	1	100.00%
LifeSign Bio Strep A Test STRip	1	*Negative	1	100.00%
Quidel Sofia/Sofia 2 Strep A FIA	46	*Negative	46	100.00%
BD Veritor Strep A	32	*Negative	32	100.00%
Alere ID NOW Strep A	18	*Negative	18	100.00%
Medline Strep A Dipstick	4	*Negative	4	100.00%
Roche Cobas Liat Strep A	1	*Negative	1	100.00%
Cepheid GeneXpert - Xpert Strep A	3	*Negative	3	100.00%
Lifesign Status Strep A Plus	1	*Negative	1	100.00%
Consult Diag Strep A Cassette	1	*Negative	1	100.00%
Alere Acceava Strep A Dipstick	110	Positive	1	0.91%
Sekisui Diagnostics OSOM	27	Positive	1	3.70%
Sekisui Diagnostics OSOM Ultra Strep	62	Positive	1	1.61%
Henry Schein OneStep+ Rapid Strep Dip	82	Positive	2	2.44%
Quidel QuickVue Dipstick Strep	85	Positive	1	1.18%
Alere Clearview Exact II Dip	37	Positive	2	5.41%
Consult Diag Strep A Dipstick	103	Positive	2	1.94%
RS-12		* Denotes Correct Response		
AG: All Methods	146			
		*Positive	141	96.58%
		Negative	5	3.42%
Abbott Signify Strep A - Cassette	1	*Positive	1	100.00%
Abbott Signify Strep A - Waived	1	*Positive	1	100.00%
Alere Acceava Strep A Dipstick	7	*Positive	6	85.71%
Fisher HC Sure-Vue Strep A (Waived)	1	*Positive	1	100.00%
Sekisui Diagnostics OSOM	5	*Positive	5	100.00%
Sekisui Diagnostics OSOM Ultra Strep	24	*Positive	24	100.00%
Henry Schein OneStep+ Rapid Strep Dip	16	*Positive	15	93.75%
Innovacon One Step Strep Ag	1	*Positive	1	100.00%
McKesson Strep A Dip	6	*Positive	6	100.00%
Quidel QuickVue Dipstick Strep	14	*Positive	14	100.00%
Quidel QuickVue In-Line - 25 Test	23	*Positive	22	95.65%
Quidel QuickVue+ Strep A	4	*Positive	4	100.00%
Select Diag Strep A Dip	1	*Positive	1	100.00%
SP Brand (Cardinal) Rapid Strep A Dipstick	12	*Positive	12	100.00%
Sekisui Diagnostics OSOM Ultra	1	*Positive	1	100.00%
Consult Diag Strep A Dipstick	11	*Positive	9	81.82%
MooreBrand Strep A Dipstick	1	*Positive	1	100.00%
BD Veritor Strep A	10	*Positive	10	100.00%
Alere ID NOW Strep A	6	*Positive	6	100.00%
Roche Cobas Liat Strep A	1	*Positive	1	100.00%
Alere Acceava Strep A Dipstick	7	Negative	1	14.29%
Henry Schein OneStep+ Rapid Strep Dip	16	Negative	1	6.25%
Quidel QuickVue In-Line - 25 Test	23	Negative	1	4.35%
Consult Diag Strep A Dipstick	11	Negative	2	18.18%
		Strep A Ag Det.		
RS-11		* Denotes Correct Response		
AG: All Methods	173			
		*Negative	171	98.84%
		Positive	2	1.16%
Abbott Signify Strep A - Waived	2	*Negative	1	50.00%
BD Check Strep A (Dip)	2	*Negative	2	100.00%
Beck Coulter ICON DS Strep	1	*Negative	1	100.00%
Alere Acceava Strep A Dipstick	2	*Negative	2	100.00%
Fisher HC Sure-Vue Select Strep Ag	1	*Negative	1	100.00%
Sekisui Diagnostics OSOM	15	*Negative	15	100.00%
Sekisui Diagnostics OSOM Ultra Strep	5	*Negative	5	100.00%
Henry Schein OneStep+ Rapid Strep Dip	9	*Negative	9	100.00%
Inverness Signify Strep A Dipstick	1	*Negative	1	100.00%
McKesson Strep A Cas	2	*Negative	2	100.00%
McKesson Strep A Dip	4	*Negative	4	100.00%
Permaxim Rediscreen Strep A Rapid	1	*Negative	1	100.00%
PSS Rapid Test Strip Strep A	1	*Negative	1	100.00%



Quidel QuickVue Dipstick Strep	8	*Negative	8	100.00%
Quidel QuickVue In-Line - 25 Test	11	*Negative	11	100.00%
Quidel QuickVue+ Strep A	17	*Negative	17	100.00%
Select Diag Strep A Dip	1	*Negative	1	100.00%
SP Brand (Cardinal) Rapid Strep A Dipstick	1	*Negative	1	100.00%
Stanbio QuStick Strep A	2	*Negative	1	50.00%
Consult Diag Strep A Dipstick	16	*Negative	16	100.00%
Meridian ImmunoCard STAT Strep A	1	*Negative	1	100.00%
Quidel Sofia/Sofia 2 Strep A FIA	9	*Negative	9	100.00%
BD Veritor Strep A	6	*Negative	6	100.00%
Quidel Solana	6	*Negative	6	100.00%
Cepheid GeneXpert - Xpert Strep A	49	*Negative	49	100.00%
Abbott Signify Strep A - Waived	2	Positive	1	50.00%
Stanbio QuStick Strep A	2	Positive	1	50.00%
RS-12		* Denotes Correct Response		
AG: All Methods	173			
		*Positive	170	98.27%
		Negative	3	1.73%
Abbott Signify Strep A - Waived	2	*Positive	2	100.00%
BD Check Strep A (Dip)	2	*Positive	2	100.00%
Beck Coulter ICON DS Strep	1	*Positive	1	100.00%
Alere Acceava Strep A Dipstick	2	*Positive	2	100.00%
Fisher HC Sure-Vue Select Strep Ag	1	*Positive	1	100.00%
Sekisui Diagnostics OSOM	15	*Positive	15	100.00%
Sekisui Diagnostics OSOM Ultra Strep	5	*Positive	5	100.00%
Henry Schein OneStep+ Rapid Strep Dip	9	*Positive	9	100.00%
Inverness Signify Strep A Dipstick	1	*Positive	1	100.00%
McKesson Strep A Cas	2	*Positive	2	100.00%
McKesson Strep A Dip	4	*Positive	4	100.00%
Permaxim Rediscreen Strep A Rapid	1	*Positive	1	100.00%
Quidel QuickVue Dipstick Strep	8	*Positive	8	100.00%
Quidel QuickVue In-Line - 25 Test	11	*Positive	11	100.00%
Quidel QuickVue+ Strep A	17	*Positive	17	100.00%
Select Diag Strep A Dip	1	*Positive	1	100.00%
SP Brand (Cardinal) Rapid Strep A Dipstick	1	*Positive	1	100.00%
Stanbio QuStick Strep A	2	*Positive	1	50.00%
Consult Diag Strep A Dipstick	16	*Positive	16	100.00%
Meridian ImmunoCard STAT Strep A	1	*Positive	1	100.00%
Quidel Sofia/Sofia 2 Strep A FIA	9	*Positive	8	88.89%
BD Veritor Strep A	6	*Positive	6	100.00%
Quidel Solana	6	*Positive	6	100.00%
Cepheid GeneXpert - Xpert Strep A	49	*Positive	49	100.00%
PSS Rapid Test Strip Strep A	1	Negative	1	100.00%
Stanbio QuStick Strep A	2	Negative	1	50.00%
Quidel Sofia/Sofia 2 Strep A FIA	9	Negative	1	11.11%
691 - Group A Strep Antigen (5 specimens)	Labs	Response	Number	Percent
		Strep A Ag Det.		
RS-11		* Denotes Correct Response		
AG: All Methods	173			
		*Negative	171	98.84%
		Positive	2	1.16%
Abbott Signify Strep A - Waived	2	*Negative	1	50.00%
BD Check Strep A (Dip)	2	*Negative	2	100.00%
Beck Coulter ICON DS Strep	1	*Negative	1	100.00%
Alere Acceava Strep A Dipstick	2	*Negative	2	100.00%
Fisher HC Sure-Vue Select Strep Ag	1	*Negative	1	100.00%
Sekisui Diagnostics OSOM	15	*Negative	15	100.00%
Sekisui Diagnostics OSOM Ultra Strep	5	*Negative	5	100.00%
Henry Schein OneStep+ Rapid Strep Dip	9	*Negative	9	100.00%
Inverness Signify Strep A Dipstick	1	*Negative	1	100.00%
McKesson Strep A Cas	2	*Negative	2	100.00%
McKesson Strep A Dip	4	*Negative	4	100.00%
Permaxim Rediscreen Strep A Rapid	1	*Negative	1	100.00%
PSS Rapid Test Strip Strep A	1	*Negative	1	100.00%
Quidel QuickVue Dipstick Strep	8	*Negative	8	100.00%
Quidel QuickVue In-Line - 25 Test	11	*Negative	11	100.00%
Quidel QuickVue+ Strep A	17	*Negative	17	100.00%
Select Diag Strep A Dip	1	*Negative	1	100.00%
SP Brand (Cardinal) Rapid Strep A Dipstick	1	*Negative	1	100.00%



Stanbio QuStick Strep A	2	*Negative	1	50.00%
Consult Diag Strep A Dipstick	16	*Negative	16	100.00%
Meridian ImmunoCard STAT Strep A	1	*Negative	1	100.00%
Quidel Sofia/Sofia 2 Strep A FIA	9	*Negative	9	100.00%
BD Veritor Strep A	6	*Negative	6	100.00%
Quidel Solana	6	*Negative	6	100.00%
Cepheid GeneXpert - Xpert Strep A	49	*Negative	49	100.00%
Abbott Signify Strep A - Waived	2	Positive	1	50.00%
Stanbio QuStick Strep A	2	Positive	1	50.00%
RS-12		* Denotes Correct Response		
AG: All Methods	173			
		*Positive	170	98.27%
		Negative	3	1.73%
Abbott Signify Strep A - Waived	2	*Positive	2	100.00%
BD Check Strep A (Dip)	2	*Positive	2	100.00%
Beck Coulter ICON DS Strep	1	*Positive	1	100.00%
Alere Acceava Strep A Dipstick	2	*Positive	2	100.00%
Fisher HC Sure-Vue Select Strep Ag	1	*Positive	1	100.00%
Sekisui Diagnostics OSOM	15	*Positive	15	100.00%
Sekisui Diagnostics OSOM Ultra Strep	5	*Positive	5	100.00%
Henry Schein OneStep+ Rapid Strep Dip	9	*Positive	9	100.00%
Inverness Signify Strep A Dipstick	1	*Positive	1	100.00%
McKesson Strep A Cas	2	*Positive	2	100.00%
McKesson Strep A Dip	4	*Positive	4	100.00%
Permaxim Rediscreen Strep A Rapid	1	*Positive	1	100.00%
Quidel QuickVue Dipstick Strep	8	*Positive	8	100.00%
Quidel QuickVue In-Line - 25 Test	11	*Positive	11	100.00%
Quidel QuickVue+ Strep A	17	*Positive	17	100.00%
Select Diag Strep A Dip	1	*Positive	1	100.00%
SP Brand (Cardinal) Rapid Strep A Dipstick	1	*Positive	1	100.00%
Stanbio QuStick Strep A	2	*Positive	1	50.00%
Consult Diag Strep A Dipstick	16	*Positive	16	100.00%
Meridian ImmunoCard STAT Strep A	1	*Positive	1	100.00%
Quidel Sofia/Sofia 2 Strep A FIA	9	*Positive	8	88.89%
BD Veritor Strep A	6	*Positive	6	100.00%
Quidel Solana	6	*Positive	6	100.00%
Cepheid GeneXpert - Xpert Strep A	49	*Positive	49	100.00%
PSS Rapid Test Strip Strep A	1	Negative	1	100.00%
Stanbio QuStick Strep A	2	Negative	1	50.00%
Quidel Sofia/Sofia 2 Strep A FIA	9	Negative	1	11.11%
RS-13		* Denotes Correct Response		
AG: All Methods	119			
		*Positive	118	99.16%
		Negative	1	0.84%
Abbott Signify Strep A - Waived	2	*Positive	2	100.00%
BD Check Strep A (Dip)	2	*Positive	2	100.00%
Beck Coulter ICON DS Strep	1	*Positive	1	100.00%
Alere Acceava Strep A Dipstick	2	*Positive	2	100.00%
Fisher HC Sure-Vue Select Strep Ag	1	*Positive	1	100.00%
Sekisui Diagnostics OSOM	15	*Positive	15	100.00%
Sekisui Diagnostics OSOM Ultra Strep	4	*Positive	4	100.00%
Henry Schein OneStep+ Rapid Strep Dip	9	*Positive	9	100.00%
Inverness Signify Strep A Dipstick	1	*Positive	1	100.00%
McKesson Strep A Cas	1	*Positive	1	100.00%
McKesson Strep A Dip	4	*Positive	4	100.00%
Permaxim Rediscreen Strep A Rapid	1	*Positive	1	100.00%
PSS Rapid Test Strip Strep A	1	*Positive	1	100.00%
Quidel QuickVue Dipstick Strep	7	*Positive	7	100.00%
Quidel QuickVue In-Line - 25 Test	11	*Positive	11	100.00%
Quidel QuickVue+ Strep A	16	*Positive	15	93.75%
SP Brand (Cardinal) Rapid Strep A Dipstick	1	*Positive	1	100.00%
Consult Diag Strep A Dipstick	13	*Positive	13	100.00%
Meridian ImmunoCard STAT Strep A	1	*Positive	1	100.00%
Quidel Sofia/Sofia 2 Strep A FIA	4	*Positive	4	100.00%
BD Veritor Strep A	1	*Positive	1	100.00%
Quidel Solana	4	*Positive	4	100.00%
Cepheid GeneXpert - Xpert Strep A	17	*Positive	17	100.00%
Quidel QuickVue+ Strep A	16	Negative	1	6.25%
RS-14		* Denotes Correct Response		



AG: All Methods	119			
		*Negative	118	99.16%
		Positive	1	0.84%
Abbott Signify Strep A - Waived	2	*Negative	2	100.00%
BD Check Strep A (Dip)	2	*Negative	2	100.00%
Beck Coulter ICON DS Strep	1	*Negative	1	100.00%
Alere Acceava Strep A Dipstick	2	*Negative	2	100.00%
Fisher HC Sure-Vue Select Strep Ag	1	*Negative	1	100.00%
Sekisui Diagnostics OSOM	15	*Negative	15	100.00%
Sekisui Diagnostics OSOM Ultra Strep	4	*Negative	3	75.00%
Henry Schein OneStep+ Rapid Strep Dip	9	*Negative	9	100.00%
Inverness Signify Strep A Dipstick	1	*Negative	1	100.00%
McKesson Strep A Cas	1	*Negative	1	100.00%
McKesson Strep A Dip	4	*Negative	4	100.00%
Permaxim Rediscreeen Strep A Rapid	1	*Negative	1	100.00%
PSS Rapid Test Strip Strep A	1	*Negative	1	100.00%
Quidel QuickVue Dipstick Strep	7	*Negative	7	100.00%
Quidel QuickVue In-Line - 25 Test	11	*Negative	11	100.00%
Quidel QuickVue+ Strep A	16	*Negative	16	100.00%
SP Brand (Cardinal) Rapid Strep A Dipstick	1	*Negative	1	100.00%
Consult Diag Strep A Dipstick	13	*Negative	13	100.00%
Meridian ImmunoCard STAT Strep A	1	*Negative	1	100.00%
Quidel Sofia/Sofia 2 Strep A FIA	4	*Negative	4	100.00%
BD Veritor Strep A	1	*Negative	1	100.00%
Quidel Solana	4	*Negative	4	100.00%
Cepheid GeneXpert - Xpert Strep A	17	*Negative	17	100.00%
Sekisui Diagnostics OSOM Ultra Strep	4	Positive	1	25.00%
RS-15		* Denotes Correct Response		
AG: All Methods	119			
		*Positive	117	98.32%
		Negative	2	1.68%
Abbott Signify Strep A - Waived	2	*Positive	2	100.00%
BD Check Strep A (Dip)	2	*Positive	2	100.00%
Beck Coulter ICON DS Strep	1	*Positive	1	100.00%
Alere Acceava Strep A Dipstick	2	*Positive	2	100.00%
Fisher HC Sure-Vue Select Strep Ag	1	*Positive	1	100.00%
Sekisui Diagnostics OSOM	15	*Positive	15	100.00%
Sekisui Diagnostics OSOM Ultra Strep	4	*Positive	3	75.00%
Henry Schein OneStep+ Rapid Strep Dip	9	*Positive	9	100.00%
Inverness Signify Strep A Dipstick	1	*Positive	1	100.00%
McKesson Strep A Cas	1	*Positive	1	100.00%
McKesson Strep A Dip	4	*Positive	4	100.00%
Permaxim Rediscreeen Strep A Rapid	1	*Positive	1	100.00%
Quidel QuickVue Dipstick Strep	7	*Positive	7	100.00%
Quidel QuickVue In-Line - 25 Test	11	*Positive	11	100.00%
Quidel QuickVue+ Strep A	16	*Positive	16	100.00%
SP Brand (Cardinal) Rapid Strep A Dipstick	1	*Positive	1	100.00%
Consult Diag Strep A Dipstick	13	*Positive	13	100.00%
Meridian ImmunoCard STAT Strep A	1	*Positive	1	100.00%
Quidel Sofia/Sofia 2 Strep A FIA	4	*Positive	4	100.00%
BD Veritor Strep A	1	*Positive	1	100.00%
Quidel Solana	4	*Positive	4	100.00%
Cepheid GeneXpert - Xpert Strep A	17	*Positive	17	100.00%
Sekisui Diagnostics OSOM Ultra Strep	4	Negative	1	25.00%
PSS Rapid Test Strip Strep A	1	Negative	1	100.00%
		Strep A Ag Det. (Waived)		
RS-11		* Denotes Correct Response		
AG: All Methods	852			
		*Negative	842	98.83%
		Positive	10	1.17%
Abbott Signify Strep A - Cassette	1	*Negative	1	100.00%
Abbott Signify Strep A - Waived	6	*Negative	6	100.00%
Beck Coulter ICON DS Strep	1	*Negative	1	100.00%
Alere Acceava Strep A Dipstick	110	*Negative	109	99.09%
Fisher HC Sure-Vue Strep A (Waived)	8	*Negative	8	100.00%
Sekisui Diagnostics OSOM	27	*Negative	26	96.30%
Sekisui Diagnostics OSOM Ultra Strep	62	*Negative	61	98.39%
Henry Schein OneStep + Rapid Strep Cass	6	*Negative	6	100.00%
Henry Schein OneStep+ Rapid Strep Dip	82	*Negative	80	97.56%



Innovacon One Step Strep Ag	1	*Negative	1	100.00%
Instant Technologies i Strep	1	*Negative	1	100.00%
Inverness Signify Strep A Dipstick	8	*Negative	8	100.00%
McKesson Strep A Cas	3	*Negative	3	100.00%
McKesson Strep A Dip	35	*Negative	35	100.00%
McKesson Strep A Test - Twist	1	*Negative	1	100.00%
Permaxim Rediscreen Strep A Rapid	1	*Negative	1	100.00%
Polymedco Poly Stat Strep A (W)	1	*Negative	1	100.00%
PSS Rapid Test Strip Strep A	2	*Negative	2	100.00%
Quidel QuickVue Dipstick Strep	85	*Negative	84	98.82%
Quidel QuickVue In-Line - 25 Test	67	*Negative	67	100.00%
Quidel QuickVue+ Strep A	24	*Negative	24	100.00%
Select Diag Strep A Dip	3	*Negative	3	100.00%
SP (Cardinal) Rapid Strep Cassette	4	*Negative	4	100.00%
SP Brand (Cardinal) Rapid Strep A Dipstick	55	*Negative	55	100.00%
Alere Clearview Exact II Dip	37	*Negative	35	94.59%
Sekisui Diagnostics OSOM Ultra	7	*Negative	7	100.00%
Consult Diag Strep A Dipstick	103	*Negative	101	98.06%
Consult Diag Strep A Twist	1	*Negative	1	100.00%
MooreBrand Strep A Dipstick	2	*Negative	2	100.00%
NDC ProAdvantage Strep Strips	1	*Negative	1	100.00%
LifeSign Bio Strep A Test STRip	1	*Negative	1	100.00%
Quidel Sofia/Sofia 2 Strep A FIA	46	*Negative	46	100.00%
BD Veritor Strep A	32	*Negative	32	100.00%
Alere ID NOW Strep A	18	*Negative	18	100.00%
Medline Strep A Dipstick	4	*Negative	4	100.00%
Roche Cobas Liat Strep A	1	*Negative	1	100.00%
Cepheid GeneXpert - Xpert Strep A	3	*Negative	3	100.00%
Lifesign Status Strep A Plus	1	*Negative	1	100.00%
Consult Diag Strep A Cassette	1	*Negative	1	100.00%
Alere Acceava Strep A Dipstick	110	Positive	1	0.91%
Sekisui Diagnostics OSOM	27	Positive	1	3.70%
Sekisui Diagnostics OSOM Ultra Strep	62	Positive	1	1.61%
Henry Schein OneStep+ Rapid Strep Dip	82	Positive	2	2.44%
Quidel QuickVue Dipstick Strep	85	Positive	1	1.18%
Alere Clearview Exact II Dip	37	Positive	2	5.41%
Consult Diag Strep A Dipstick	103	Positive	2	1.94%
RS-12		* Denotes Correct Response		
AG: All Methods	146			
		*Positive	141	96.58%
		Negative	5	3.42%
Abbott Signify Strep A - Cassette	1	*Positive	1	100.00%
Abbott Signify Strep A - Waived	1	*Positive	1	100.00%
Alere Acceava Strep A Dipstick	7	*Positive	6	85.71%
Fisher HC Sure-Vue Strep A (Waived)	1	*Positive	1	100.00%
Sekisui Diagnostics OSOM	5	*Positive	5	100.00%
Sekisui Diagnostics OSOM Ultra Strep	24	*Positive	24	100.00%
Henry Schein OneStep+ Rapid Strep Dip	16	*Positive	15	93.75%
Innovacon One Step Strep Ag	1	*Positive	1	100.00%
McKesson Strep A Dip	6	*Positive	6	100.00%
Quidel QuickVue Dipstick Strep	14	*Positive	14	100.00%
Quidel QuickVue In-Line - 25 Test	23	*Positive	22	95.65%
Quidel QuickVue+ Strep A	4	*Positive	4	100.00%
Select Diag Strep A Dip	1	*Positive	1	100.00%
SP Brand (Cardinal) Rapid Strep A Dipstick	12	*Positive	12	100.00%
Sekisui Diagnostics OSOM Ultra	1	*Positive	1	100.00%
Consult Diag Strep A Dipstick	11	*Positive	9	81.82%
MooreBrand Strep A Dipstick	1	*Positive	1	100.00%
BD Veritor Strep A	10	*Positive	10	100.00%
Alere ID NOW Strep A	6	*Positive	6	100.00%
Roche Cobas Liat Strep A	1	*Positive	1	100.00%
Alere Acceava Strep A Dipstick	7	Negative	1	14.29%
Henry Schein OneStep+ Rapid Strep Dip	16	Negative	1	6.25%
Quidel QuickVue In-Line - 25 Test	23	Negative	1	4.35%
Consult Diag Strep A Dipstick	11	Negative	2	18.18%
RS-13		* Denotes Correct Response		
AG: All Methods	1			
		*Positive	1	100.00%
Alere ID NOW Strep A	1	*Positive	1	100.00%



RS-14		* Denotes Correct Response		
AG: All Methods	2			
		*Negative	2	100.00%
Quidel QuickVue In-Line - 25 Test	1	*Negative	1	100.00%
Alere ID NOW Strep A	1	*Negative	1	100.00%
RS-15		* Denotes Correct Response		
AG: All Methods	3			
		*Positive	2	66.67%
		Negative	1	33.33%
Quidel QuickVue In-Line - 25 Test	2	*Positive	1	50.00%
Alere ID NOW Strep A	1	*Positive	1	100.00%
Quidel QuickVue In-Line - 25 Test	2	Negative	1	50.00%
		Strep A Ag Det. (International)		
RS-11		* Denotes Correct Response		
AG: All Methods	1			
		*Negative	1	100.00%
Alere Clearview Exact II Cas	1	*Negative	1	100.00%
RS-12		* Denotes Correct Response		
AG: All Methods	1			
		*Positive	1	100.00%
Alere Clearview Exact II Cas	1	*Positive	1	100.00%
RS-13		* Denotes Correct Response		
AG: All Methods	1			
		*Positive	1	100.00%
Alere Clearview Exact II Cas	1	*Positive	1	100.00%
RS-14		* Denotes Correct Response		
AG: All Methods	1			
		*Negative	1	100.00%
Alere Clearview Exact II Cas	1	*Negative	1	100.00%
RS-15		* Denotes Correct Response		
AG: All Methods	1			
		*Positive	1	100.00%
Alere Clearview Exact II Cas	1	*Positive	1	100.00%
739 - Sickie Cell Screen	Labs	Response	Number	Percent
		Sickle Cell Screen		
SS-3		* Denotes Correct Response		
AG: All Methods	2			
		*Negative	2	100.00%
Pacific Hemostasis Sickie Screen	1	*Negative	1	100.00%
Streck SickieDex	1	*Negative	1	100.00%
679 - Syphilis Serology	Labs	Response	Number	Percent
		RPR		
SY-11		* Denotes Correct Response		
AG: All Methods	10			
		*Reactive	10	100.00%
ASI RPR	3	*Reactive	3	100.00%
Becton Dickinson BBL	4	*Reactive	4	100.00%
BD Macro-Vue RPR	1	*Reactive	1	100.00%
Wampole Impact	1	*Reactive	1	100.00%
Intl ASI RPR	1	*Reactive	1	100.00%
SY-12		* Denotes Correct Response		
AG: All Methods	10			
		*Nonreactive	10	100.00%
ASI RPR	3	*Nonreactive	3	100.00%
Becton Dickinson BBL	4	*Nonreactive	4	100.00%
BD Macro-Vue RPR	1	*Nonreactive	1	100.00%
Wampole Impact	1	*Nonreactive	1	100.00%
Intl ASI RPR	1	*Nonreactive	1	100.00%
SY-13		* Denotes Correct Response		
AG: All Methods	10			
		*Reactive	10	100.00%
ASI RPR	3	*Reactive	3	100.00%
Becton Dickinson BBL	4	*Reactive	4	100.00%
BD Macro-Vue RPR	1	*Reactive	1	100.00%
Wampole Impact	1	*Reactive	1	100.00%
Intl ASI RPR	1	*Reactive	1	100.00%
SY-14		* Denotes Correct Response		
AG: All Methods	10			



		*Nonreactive	10	100.00%
ASI RPR	3	*Nonreactive	3	100.00%
Becton Dickinson BBL	4	*Nonreactive	4	100.00%
BD Macro-Vue RPR	1	*Nonreactive	1	100.00%
Wampole Impact	1	*Nonreactive	1	100.00%
Intl ASI RPR	1	*Nonreactive	1	100.00%
SY-15		* Denotes Correct Response		
AG: All Methods	10			
		*Reactive	10	100.00%
ASI RPR	3	*Reactive	3	100.00%
Becton Dickinson BBL	4	*Reactive	4	100.00%
BD Macro-Vue RPR	1	*Reactive	1	100.00%
Wampole Impact	1	*Reactive	1	100.00%
Intl ASI RPR	1	*Reactive	1	100.00%
		RPR/VDRL (titer)		
SY-11		* Denotes Correct Response		
AG: All Methods	3			
		*1:1	1	33.33%
		*1:2	2	66.67%
BD Macro-Vue RPR	1	*1:1	1	100.00%
ASI RPR	1	*1:2	1	100.00%
Becton Dickinson BBL	1	*1:2	1	100.00%
SY-12		* Denotes Correct Response		
AG: All Methods	3			
		*N/A (Negative)	1	33.33%
		*Nonreactive	2	66.67%
ASI RPR	1	*N/A (Negative)	1	100.00%
Becton Dickinson BBL	1	*Nonreactive	1	100.00%
BD Macro-Vue RPR	1	*Nonreactive	1	100.00%
SY-13		* Denotes Correct Response		
AG: All Methods	3			
		*1:2	1	33.33%
		*1:4	2	66.67%
BD Macro-Vue RPR	1	*1:2	1	100.00%
ASI RPR	1	*1:4	1	100.00%
Becton Dickinson BBL	1	*1:4	1	100.00%
SY-14		* Denotes Correct Response		
AG: All Methods	3			
		*N/A (Negative)	1	33.33%
		*Nonreactive	2	66.67%
ASI RPR	1	*N/A (Negative)	1	100.00%
Becton Dickinson BBL	1	*Nonreactive	1	100.00%
BD Macro-Vue RPR	1	*Nonreactive	1	100.00%
SY-15		* Denotes Correct Response		
AG: All Methods	3			
		*1:2	3	100.00%
ASI RPR	1	*1:2	1	100.00%
Becton Dickinson BBL	1	*1:2	1	100.00%
BD Macro-Vue RPR	1	*1:2	1	100.00%
		RPR (International)		
SY-11		* Denotes Correct Response		
AG: All Methods	3			
		*Reactive	3	100.00%
BD Macro-Vue RPR	1	*Reactive	1	100.00%
Fisher HC Sure-Vue RPR	1	*Reactive	1	100.00%
Wampole Impact	1	*Reactive	1	100.00%
SY-12		* Denotes Correct Response		
AG: All Methods	3			
		*Nonreactive	3	100.00%
BD Macro-Vue RPR	1	*Nonreactive	1	100.00%
Fisher HC Sure-Vue RPR	1	*Nonreactive	1	100.00%
Wampole Impact	1	*Nonreactive	1	100.00%
SY-13		* Denotes Correct Response		
AG: All Methods	3			
		*Reactive	3	100.00%
BD Macro-Vue RPR	1	*Reactive	1	100.00%
Fisher HC Sure-Vue RPR	1	*Reactive	1	100.00%
Wampole Impact	1	*Reactive	1	100.00%



SY-14		* Denotes Correct Response		
AG: All Methods	3			
		*Nonreactive	3	100.00%
BD Macro-Vue RPR	1	*Nonreactive	1	100.00%
Fisher HC Sure-Vue RPR	1	*Nonreactive	1	100.00%
Wampole Impact	1	*Nonreactive	1	100.00%
SY-15		* Denotes Correct Response		
AG: All Methods	3			
		*Reactive	3	100.00%
BD Macro-Vue RPR	1	*Reactive	1	100.00%
Fisher HC Sure-Vue RPR	1	*Reactive	1	100.00%
Wampole Impact	1	*Reactive	1	100.00%
776 - Throat Culture;780 - Throat & Genital Cultures	Labs	Response	Number	Percent
Throat Culture				
TC-11		* Denotes Correct Response		
AG: All Methods	57			
		*Beta hemolysis, referred for ID/confirmation	2	3.51%
		*Positive for Group A Streptococcus	25	43.86%
		*Presumptive Positive for Group A Streptococcus	24	42.11%
		*Presumptive Streptococcus species	1	1.75%
		*Streptococcus, beta-hemolytic, Group A	1	1.75%
		Negative for Group A Streptococcus	4	7.02%
Cultures	57	*Beta hemolysis, referred for ID/confirmation	2	3.51%
Cultures		*Positive for Group A Streptococcus	25	43.86%
Cultures		*Presumptive Positive for Group A Streptococcus	24	42.11%
Cultures		*Presumptive Streptococcus species	1	1.75%
Cultures		*Streptococcus, beta-hemolytic, Group A	1	1.75%
Cultures		Negative for Group A Streptococcus	4	7.02%
TC-12		* Denotes Correct Response		
AG: All Methods	57			
		*Growth, referred for identification	1	1.75%
		*Negative for Group A Streptococcus	46	80.70%
		*No beta hemolysis	3	5.26%
		Positive for Group A Streptococcus	4	7.02%
		Presumptive Positive for Group A Streptococcus	3	5.26%
Cultures	57	*Growth, referred for identification	1	1.75%
Cultures		*Negative for Group A Streptococcus	46	80.70%
Cultures		*No beta hemolysis	3	5.26%
Cultures		Positive for Group A Streptococcus	4	7.02%
Cultures		Presumptive Positive for Group A Streptococcus	3	5.26%
TC-13		* Denotes Correct Response		
AG: All Methods	58			
		*Beta hemolysis, referred for ID/confirmation	2	3.45%
		*Positive for Group A Streptococcus	25	43.10%
		*Presumptive Positive for Group A Streptococcus	25	43.10%
		*Presumptive Streptococcus species	1	1.72%
		*Streptococcus, beta-hemolytic, Group A	1	1.72%
		Negative for Group A Streptococcus	4	6.90%
Cultures	58	*Beta hemolysis, referred for ID/confirmation	2	3.45%
Cultures		*Positive for Group A Streptococcus	25	43.10%
Cultures		*Presumptive Positive for Group A Streptococcus	25	43.10%
Cultures		*Presumptive Streptococcus species	1	1.72%
Cultures		*Streptococcus, beta-hemolytic, Group A	1	1.72%
Cultures		Negative for Group A Streptococcus	4	6.90%
TC-14		* Denotes Correct Response		
AG: All Methods	70			
		*Beta hemolysis, referred for ID/confirmation	4	5.71%
		*Growth, referred for identification	1	1.43%
		*Negative for Group A Streptococcus	56	80.00%
		*No Growth - Strep Selective medium	3	4.29%
		*Presumptive Streptococcus species	1	1.43%
		*Streptococcus, beta-hemolytic, not Group A	2	2.86%
		Positive for Group A Streptococcus	3	4.29%
Cultures	70	*Beta hemolysis, referred for ID/confirmation	4	5.71%
Cultures		*Growth, referred for identification	1	1.43%
Cultures		*Negative for Group A Streptococcus	56	80.00%
Cultures		*No Growth - Strep Selective medium	3	4.29%
Cultures		*Presumptive Streptococcus species	1	1.43%



Cultures		*Streptococcus, beta-hemolytic, not Group A	2	2.86%
Cultures		Positive for Group A Streptococcus	3	4.29%
TC-15		* Denotes Correct Response		
AG: All Methods	70			
		*Growth, referred for identification	1	1.43%
		*Negative for Group A Streptococcus	57	81.43%
		*No beta hemolysis	4	5.71%
		*No Growth - Strep Selective medium	8	11.43%
Cultures	70	*Growth, referred for identification	1	1.43%
Cultures		*Negative for Group A Streptococcus	57	81.43%
Cultures		*No beta hemolysis	4	5.71%
Cultures		*No Growth - Strep Selective medium	8	11.43%
640 - Urinalysis;642 - Urinalysis/hCG;639 - Complete Urinalysis Package;601 - Waived/PPM Tests II;609 - Waived/PPM Tests I;647 - Urinalysis/hCG/Fecal Occult Blood	Labs	Response	Number	Percent
		Urine Sp Gravity		
UA-3		* Denotes Correct Response		
PG: Siemens Clinitek 10 / 100	11			
		*1.010	3	27.27%
		*1.015	4	36.36%
		*1.020	2	18.18%
		1.005	1	9.09%
		1.025	1	9.09%
UA-3		* Denotes Correct Response		
PG: Siemens Clinitek 50	63			
		*1.010	11	17.46%
		*1.015	28	44.44%
		*1.020	24	38.10%
UA-3		* Denotes Correct Response		
PG: Siemens Clinitek Status	686			
		*1.010	2	0.29%
		*1.015	296	43.15%
		*1.020	382	55.69%
		1.025	4	0.58%
		1.030	2	0.29%
UA-3		* Denotes Correct Response		
PG: Siemens Clinitek Status / Siemens - Plus Sys / Siemens Urine Dipstick	16			
		*1.015	8	50.00%
		*1.020	7	43.75%
		1.025	1	6.25%
UA-3		* Denotes Correct Response		
PG: Siemens Multistix - Visual	193			
		*1.005	29	15.03%
		*1.010	92	47.67%
		*1.015	37	19.17%
		*1.020	31	16.06%
		1.000	3	1.55%
		1.025	1	0.52%
UA-3		* Denotes Correct Response		
PG: Diaq Test Group Clarity Urochek	24			
		*1.010	16	66.67%
		*1.015	3	12.50%
		*1.020	4	16.67%
		1.000	1	4.17%
UA-3		* Denotes Correct Response		
PG: Roche Chemstrips - Visual	56			
		*1.005	10	17.86%
		*1.010	28	50.00%
		*1.015	10	17.86%
		1.000	6	10.71%
		1.020	1	1.79%
		1.025	1	1.79%
UA-3		* Denotes Correct Response		
PG: Roche Urisys	30			
		*1.005	1	3.33%
		*1.010	17	56.67%
		*1.015	10	33.33%



		*1.020	1	3.33%
		1.025	1	3.33%
UA-3		* Denotes Correct Response		
PG: Roche Urisys / Roche - Conventional / Roche Urine Dipstick	2			
		*1.015	1	50.00%
		*1.020	1	50.00%
UA-3		* Denotes Correct Response		
PG: Roche Urisys / Roche - Arbitrary / Roche Urine Arbitrary Result	1			
		*1.010	1	100.00%
UA-3		* Denotes Correct Response		
PG: Siemens Clinitek Advantus	25			
		*1.010	4	16.00%
		*1.015	21	84.00%
UA-3		* Denotes Correct Response		
PG: 120 Urine Analyzer	44			
		*1.010	22	50.00%
		*1.015	22	50.00%
UA-3		* Denotes Correct Response		
SG: Select Siemens Clinitek/Visual	984			
		*1.005	30	3.05%
		*1.010	114	11.59%
		*1.015	387	39.33%
		*1.020	442	44.92%
		1.000	3	0.30%
		1.005	30	3.05%
		1.025	6	0.61%
		1.030	2	0.20%
UA-3		* Denotes Correct Response		
SG: Select Roche Urinalysis Inst/Visual	92			
		*1.005	12	13.04%
		*1.010	49	53.26%
		*1.015	21	22.83%
		*1.020	2	2.17%
		1.000	6	6.52%
		1.020	2	2.17%
		1.025	2	2.17%
UA-3		* Denotes Correct Response		
SG: Select Consult Diag Urinalysis Inst/Visual	21			
		*1.010	13	61.90%
		*1.015	6	28.57%
		*1.020	2	9.52%
UA-3		* Denotes Correct Response		
SG: Select McKesson Urinalysis Inst/Visual	51			
		*1.010	26	50.98%
		*1.015	25	49.02%
UA-3		* Denotes Correct Response		
SG: Select Henry Schein Urinalysis Inst/Visual	9			
		*1.010	1	11.11%
		*1.015	7	77.78%
		*1.020	1	11.11%
UA-3		* Denotes Correct Response		
SG: Select Medi-Lab Performance UA Inst/Visual	8			
		*1.005	2	25.00%
		*1.010	5	62.50%
		*1.015	1	12.50%
UA-3		* Denotes Correct Response		
MG: Method Group	1211			
		*1.005	49	4.05%
		*1.010	231	19.08%
		*1.015	457	37.74%
		*1.020	452	37.32%
		1.000	10	0.83%
		1.005	49	4.05%
		1.020	452	37.32%
		1.025	9	0.74%
		1.030	3	0.25%
UA-3		* Denotes Correct Response		



MG: Method Group	2			
		*1.015	1	50.00%
		*1.020	1	50.00%
UA-3		* Denotes Correct Response		
MG: Siemens - Plus Sys	17			
		*1.015	9	52.94%
		*1.020	7	41.18%
		1.025	1	5.88%
UA-3		* Denotes Correct Response		
MG: Siemens - Plus Sys	1			
		*1.010	1	100.00%
		Urine Bilirubin		
UA-3		* Denotes Correct Response		
PG: Siemens Clinitek 10 / 100	11			
		*Large (3+)	4	36.36%
		*Moderate (2+)	7	63.64%
UA-3		* Denotes Correct Response		
PG: Siemens Clinitek 50	63			
		*Moderate (2+)	53	84.13%
		*Trace/Small (1+)	8	12.70%
		Large (3+)	1	1.59%
		Negative	1	1.59%
UA-3		* Denotes Correct Response		
PG: Siemens Clinitek Status	679			
		*Large (3+)	18	2.65%
		*Moderate (2+)	578	85.13%
		*Trace/Small (1+)	78	11.49%
		Negative	5	0.74%
UA-3		* Denotes Correct Response		
PG: Siemens Clinitek Status / Siemens - Plus Sys / Siemens Urine Dipstick	16			
		*1+	6	37.50%
		*2+	10	62.50%
UA-3		* Denotes Correct Response		
PG: Siemens Multistix - Visual	187			
		*Large (3+)	78	41.71%
		*Moderate (2+)	68	36.36%
		*Trace/Small (1+)	33	17.65%
		Negative	8	4.28%
UA-3		* Denotes Correct Response		
PG: Diag Test Group Clarity Urochek	23			
		*Large (3+)	6	26.09%
		*Moderate (2+)	10	43.48%
		*Trace/Small (1+)	6	26.09%
		Negative	1	4.35%
UA-3		* Denotes Correct Response		
PG: Roche Chemstrips - Visual	57			
		*Large (3+)	10	17.54%
		*Moderate (2+)	28	49.12%
		*Trace/Small (1+)	13	22.81%
		Negative	6	10.53%
UA-3		* Denotes Correct Response		
PG: Roche Urisys	30			
		*Large (3+)	6	20.00%
		*Moderate (2+)	3	10.00%
		*Trace/Small (1+)	21	70.00%
UA-3		* Denotes Correct Response		
PG: Roche Urisys / Roche - Conventional / Roche Urine Dipstick	2			
		*1 mg/dL	1	50.00%
		*3 mg/dL	1	50.00%
UA-3		* Denotes Correct Response		
PG: Roche Urisys / Roche - Arbitrary / Roche Urine Arbitrary Result	1			
		*++	1	100.00%
UA-3		* Denotes Correct Response		
PG: Siemens Clinitek Advantus	25			
		*Moderate (2+)	23	92.00%
		*Trace/Small (1+)	2	8.00%



UA-3		* Denotes Correct Response		
PG: 120 Urine Analyzer	44			
		*Moderate (2+)	14	31.82%
		*Trace/Small (1+)	29	65.91%
		Negative	1	2.27%
UA-3		* Denotes Correct Response		
SG: Select Siemens Clinitek/Visual	971			
		*Large (3+)	101	10.40%
		*Moderate (2+)	733	75.49%
		*Trace/Small (1+)	123	12.67%
		Large (3+)	101	10.40%
		Negative	14	1.44%
UA-3		* Denotes Correct Response		
SG: Select Roche Urinalysis Inst/Visual	93			
		*Large (3+)	18	19.35%
		*Moderate (2+)	33	35.48%
		*Trace/Small (1+)	35	37.63%
		Negative	7	7.53%
UA-3		* Denotes Correct Response		
SG: Select Roche Urinalysis Inst/Visual / Roche - Arbitrary / Roche Urine Arbitrary Result	1			
		*++	1	100.00%
UA-3		* Denotes Correct Response		
SG: Select Consult Diag Urinalysis Inst/Visual	21			
		*Large (3+)	1	4.76%
		*Moderate (2+)	11	52.38%
		*Trace/Small (1+)	7	33.33%
		Negative	2	9.52%
UA-3		* Denotes Correct Response		
SG: Select McKesson Urinalysis Inst/Visual	51			
		*Large (3+)	2	3.92%
		*Moderate (2+)	15	29.41%
		*Trace/Small (1+)	32	62.75%
		Negative	2	3.92%
UA-3		* Denotes Correct Response		
SG: Select Henry Schein Urinalysis Inst/Visual	9			
		*Large (3+)	4	44.44%
		*Moderate (2+)	4	44.44%
		*Trace/Small (1+)	1	11.11%
UA-3		* Denotes Correct Response		
SG: Select Medi-Lab Performance UA Inst/Visual	8			
		*Trace/Small (1+)	7	87.50%
		Negative	1	12.50%
UA-3		* Denotes Correct Response		
MG: Method Group	1198			
		*Large (3+)	138	11.52%
		*Moderate (2+)	811	67.70%
		*Trace/Small (1+)	220	18.36%
		Large (3+)	138	11.52%
		Negative	29	2.42%
UA-3		* Denotes Correct Response		
MG: Method Group	2			
		*1 mg/dL	1	50.00%
		*3 mg/dL	1	50.00%
UA-3		* Denotes Correct Response		
MG: Siemens - Plus Sys	17			
		*1+	7	41.18%
		*2+	10	58.82%
UA-3		* Denotes Correct Response		
MG: Siemens - Plus Sys	1			
		*++	1	100.00%
		Urine Ketones		
UA-3		* Denotes Correct Response		
PG: Siemens Clinitek 10 / 100	11			
		*Large (3+) / 80 mg/dL	11	100.00%
UA-3		* Denotes Correct Response		
PG: Siemens Clinitek 50	63			
		*>=160 mg/dL	38	60.32%



		*Large (3+) / 80 mg/dL	25	39.68%
UA-3		* Denotes Correct Response		
PG: Siemens Clinitek Status	686			
		*>=160 mg/dL	609	88.78%
		*Large (3+) / 80 mg/dL	73	10.64%
		Moderate (2+) / 40 mg/dL	3	0.44%
		Small (1+) / 15mg/dL	1	0.15%
UA-3		* Denotes Correct Response		
PG: Siemens Clinitek Status / Siemens - Plus Sys / Siemens Urine Dipstick	16			
		*3+	4	25.00%
		*4+	12	75.00%
UA-3		* Denotes Correct Response		
PG: Siemens Multistix - Visual	195			
		*Large (3+) / 80 mg/dL	192	98.46%
		Moderate (2+) / 40 mg/dL	3	1.54%
UA-3		* Denotes Correct Response		
PG: Diag Test Group Clarity Urochek	24			
		*Large (3+) / 80 mg/dL	24	100.00%
UA-3		* Denotes Correct Response		
PG: Roche Chemstrips - Visual	58			
		*Large (3+) / 80 mg/dL	57	98.28%
		Moderate (2+) / 40 mg/dL	1	1.72%
UA-3		* Denotes Correct Response		
PG: Roche Urisys	30			
		*Large (3+) / 80 mg/dL	27	90.00%
		Moderate (2+) / 40 mg/dL	2	6.67%
		Small (1+) / 15mg/dL	1	3.33%
UA-3		* Denotes Correct Response		
PG: Roche Urisys / Roche - Conventional / Roche Urine Dipstick	2			
		*150 mg/dL	2	100.00%
UA-3		* Denotes Correct Response		
PG: Roche Urisys / Roche - Arbitrary / Roche Urine Arbitrary Result	1			
		*+++	1	100.00%
UA-3		* Denotes Correct Response		
PG: Siemens Clinitek Advantus	24			
		*Large (3+) / 80 mg/dL	24	100.00%
UA-3		* Denotes Correct Response		
PG: 120 Urine Analyzer	44			
		*Large (3+) / 80 mg/dL	44	100.00%
UA-3		* Denotes Correct Response		
SG: Select Siemens Clinitek/Visual	985			
		*>=160 mg/dL	647	65.69%
		*Large (3+) / 80 mg/dL	331	33.60%
		Moderate (2+) / 40 mg/dL	6	0.61%
		Small (1+) / 15mg/dL	1	0.10%
UA-3		* Denotes Correct Response		
SG: Select Roche Urinalysis Inst/Visual	94			
		*Large (3+) / 80 mg/dL	88	93.62%
		Moderate (2+) / 40 mg/dL	3	3.19%
		Negative	2	2.13%
		Small (1+) / 15mg/dL	1	1.06%
UA-3		* Denotes Correct Response		
SG: Select Consult Diag Urinalysis Inst/Visual	21			
		*Large (3+) / 80 mg/dL	21	100.00%
UA-3		* Denotes Correct Response		
SG: Select McKesson Urinalysis Inst/Visual	51			
		*Large (3+) / 80 mg/dL	51	100.00%
UA-3		* Denotes Correct Response		
SG: Select Henry Schein Urinalysis Inst/Visual	9			
		*Large (3+) / 80 mg/dL	8	88.89%
		Moderate (2+) / 40 mg/dL	1	11.11%
UA-3		* Denotes Correct Response		
SG: Select Medi-Lab Performance UA Inst/Visual	8			
		*Large (3+) / 80 mg/dL	8	100.00%
UA-3		* Denotes Correct Response		
MG: Method Group	1215			



		*≥160 mg/dL	647	53.25%
		*Large (3+) / 80 mg/dL	554	45.60%
		Moderate (2+) / 40 mg/dL	10	0.82%
		Negative	2	0.16%
		Small (1+) / 15mg/dL	2	0.16%
UA-3		* Denotes Correct Response		
MG: Method Group	2			
		*150 mg/dL	2	100.00%
UA-3		* Denotes Correct Response		
MG: Siemens - Plus Sys	17			
		*3+	5	29.41%
		*4+	12	70.59%
UA-3		* Denotes Correct Response		
MG: Siemens - Plus Sys	1			
		*+++	1	100.00%
		Urine pH		
UA-3		* Denotes Correct Response		
PG: Siemens Clinitek 10 / 100	11			
		*8 or greater	10	90.91%
		7	1	9.09%
UA-3		* Denotes Correct Response		
PG: Siemens Clinitek 50	63			
		*7.5	5	7.94%
		*8 or greater	58	92.06%
UA-3		* Denotes Correct Response		
PG: Siemens Clinitek Status	685			
		*7.5	41	5.99%
		*8 or greater	641	93.58%
		5 or less	2	0.29%
		7	1	0.15%
UA-3		* Denotes Correct Response		
PG: Siemens Clinitek Status / Siemens - Plus Sys / Siemens Urine Dipstick	16			
		*8.0	4	25.00%
		*8.5	10	62.50%
		7.5	2	12.50%
UA-3		* Denotes Correct Response		
PG: Siemens Multistix - Visual	194			
		*7.5	17	8.76%
		*8 or greater	173	89.18%
		5 or less	1	0.52%
		6	1	0.52%
		6.5	1	0.52%
		7	1	0.52%
UA-3		* Denotes Correct Response		
PG: Diag Test Group Clarity Urochek	24			
		*7.5	2	8.33%
		*8 or greater	22	91.67%
UA-3		* Denotes Correct Response		
PG: Roche Chemstrips - Visual	56			
		*7.5	1	1.79%
		*8 or greater	51	91.07%
		7	4	7.14%
UA-3		* Denotes Correct Response		
PG: Roche Urisys	30			
		*8 or greater	27	90.00%
		7	3	10.00%
UA-3		* Denotes Correct Response		
PG: Roche Urisys / Roche - Conventional / Roche Urine Dipstick	2			
		*8	2	100.00%
UA-3		* Denotes Correct Response		
PG: Roche Urisys / Roche - Arbitrary / Roche Urine Arbitrary Result	1			
		*8.0	1	100.00%
UA-3		* Denotes Correct Response		
PG: Siemens Clinitek Advantus	25			
		*7.5	2	8.00%
		*8 or greater	23	92.00%



UA-3		* Denotes Correct Response		
PG: 120 Urine Analyzer	44			
		*7.5	9	20.45%
		*8 or greater	35	79.55%
UA-3		* Denotes Correct Response		
SG: Select Siemens Clinitek/Visual	985			
		*7.5	65	6.60%
		*8 or greater	912	92.59%
		5 or less	3	0.30%
		6	1	0.10%
		6.5	1	0.10%
		7	3	0.30%
UA-3		* Denotes Correct Response		
SG: Select Roche Urinalysis Inst/Visual	92			
		*7.5	1	1.09%
		*8 or greater	83	90.22%
		7	8	8.70%
UA-3		* Denotes Correct Response		
SG: Select Consult Diag Urinalysis Inst/Visual	21			
		*7.5	4	19.05%
		*8 or greater	17	80.95%
UA-3		* Denotes Correct Response		
SG: Select McKesson Urinalysis Inst/Visual	51			
		*7.5	12	23.53%
		*8 or greater	39	76.47%
UA-3		* Denotes Correct Response		
SG: Select Henry Schein Urinalysis Inst/Visual	9			
		*8 or greater	9	100.00%
UA-3		* Denotes Correct Response		
SG: Select Medi-Lab Performance UA Inst/Visual	8			
		*7.5	2	25.00%
		*8 or greater	6	75.00%
UA-3		* Denotes Correct Response		
MG: Method Group	1213			
		*7.5	92	7.58%
		*8 or greater	1103	90.93%
		5 or less	3	0.25%
		6	1	0.08%
		6.5	1	0.08%
		7	13	1.07%
UA-3		* Denotes Correct Response		
MG: Method Group	2			
		*8	2	100.00%
UA-3		* Denotes Correct Response		
MG: Siemens - Plus Sys	17			
		*8.0	5	29.41%
		*8.5	10	58.82%
		7.5	2	11.76%
UA-3		* Denotes Correct Response		
MG: Siemens - Plus Sys	1			
		*8.0	1	100.00%
UA-3		Urine Protein - Qual		
UA-3		* Denotes Correct Response		
PG: Siemens Clinitek 10 / 100	11			
		*100 (2+)	7	63.64%
		*300-500 (3+)	4	36.36%
UA-3		* Denotes Correct Response		
PG: Siemens Clinitek 50	63			
		*100 (2+)	18	28.57%
		*300-500 (3+)	45	71.43%
UA-3		* Denotes Correct Response		
PG: Siemens Clinitek Status	687			
		*100 (2+)	214	31.15%
		*300-500 (3+)	466	67.83%
		>/= 1000 (4+)	2	0.29%
		30 (1+)	2	0.29%
		Negative	3	0.44%
UA-3		* Denotes Correct Response		



PG: Siemens Clinitek Status / Siemens - Plus Sys / Siemens Urine Dipstick	16			
		*2+	3	18.75%
		*3+	13	81.25%
UA-3		* Denotes Correct Response		
PG: Siemens Multistix - Visual	196			
		*100 (2+)	98	50.00%
		*300-500 (3+)	85	43.37%
		>/= 1000 (4+)	7	3.57%
		30 (1+)	4	2.04%
		Negative	2	1.02%
UA-3		* Denotes Correct Response		
PG: Diag Test Group Clarity Urochek	24			
		*100 (2+)	15	62.50%
		*300-500 (3+)	8	33.33%
		>/= 1000 (4+)	1	4.17%
UA-3		* Denotes Correct Response		
PG: Roche Chemstrips - Visual	63			
		*100 (2+)	46	73.02%
		*300-500 (3+)	14	22.22%
		30 (1+)	3	4.76%
UA-3		* Denotes Correct Response		
PG: Roche Urisys	30			
		*100 (2+)	20	66.67%
		*300-500 (3+)	6	20.00%
		30 (1+)	4	13.33%
UA-3		* Denotes Correct Response		
PG: Roche Urisys / Roche - Conventional / Roche Urine Dipstick	2			
		*500 mg/dL	1	50.00%
		30 mg/dL	1	50.00%
UA-3		* Denotes Correct Response		
PG: Roche Urisys / Roche - Arbitrary / Roche Urine Arbitrary Result	1			
		*++	1	100.00%
UA-3		* Denotes Correct Response		
PG: Siemens Clinitek Advantus	25			
		*100 (2+)	25	100.00%
UA-3		* Denotes Correct Response		
PG: 120 Urine Analyzer	44			
		*100 (2+)	38	86.36%
		*300-500 (3+)	6	13.64%
UA-3		* Denotes Correct Response		
SG: Select Siemens Clinitek/Visual	992			
		*100 (2+)	367	37.00%
		*300-500 (3+)	604	60.89%
		>/= 1000 (4+)	9	0.91%
		30 (1+)	7	0.71%
		Negative	5	0.50%
UA-3		* Denotes Correct Response		
SG: Select Roche Urinalysis Inst/Visual	99			
		*100 (2+)	70	70.71%
		*300-500 (3+)	21	21.21%
		30 (1+)	7	7.07%
		Negative	1	1.01%
UA-3		* Denotes Correct Response		
SG: Select Consult Diag Urinalysis Inst/Visual	21			
		*100 (2+)	12	57.14%
		*300-500 (3+)	9	42.86%
UA-3		* Denotes Correct Response		
SG: Select McKesson Urinalysis Inst/Visual	51			
		*100 (2+)	43	84.31%
		*300-500 (3+)	8	15.69%
UA-3		* Denotes Correct Response		
SG: Select Henry Schein Urinalysis Inst/Visual	9			
		*100 (2+)	6	66.67%
		*300-500 (3+)	2	22.22%
		30 (1+)	1	11.11%
UA-3		* Denotes Correct Response		
SG: Select Medi-Lab Performance UA Inst/Visual	8			



		*100 (2+)	6	75.00%
		*300-500 (3+)	1	12.50%
		30 (1+)	1	12.50%
UA-3		* Denotes Correct Response		
MG: Method Group	1228			
		*100 (2+)	536	43.65%
		*300-500 (3+)	657	53.50%
		>/= 1000 (4+)	11	0.90%
		30 (1+)	18	1.47%
		Negative	6	0.49%
UA-3		* Denotes Correct Response		
MG: Method Group	2			
		*500 mg/dL	1	50.00%
		30 mg/dL	1	50.00%
UA-3		* Denotes Correct Response		
MG: Siemens - Plus Sys	17			
		*2+	4	23.53%
		*3+	13	76.47%
UA-3		* Denotes Correct Response		
MG: Siemens - Plus Sys	1			
		*++	1	100.00%
		Urine Blood/Hemoglobin		
UA-3		* Denotes Correct Response		
PG: Siemens Clinitek 10 / 100	11			
		*Large (3+)	9	81.82%
		*Moderate (2+)	2	18.18%
UA-3		* Denotes Correct Response		
PG: Siemens Clinitek 50	63			
		*Large (3+)	43	68.25%
		*Moderate (2+)	18	28.57%
		Negative	1	1.59%
		Trace	1	1.59%
UA-3		* Denotes Correct Response		
PG: Siemens Clinitek Status	684			
		*Large (3+)	550	80.41%
		*Moderate (2+)	130	19.01%
		Negative	2	0.29%
		Small (1+)	1	0.15%
		Trace	1	0.15%
UA-3		* Denotes Correct Response		
PG: Siemens Clinitek Status / Siemens - Plus Sys / Siemens Urine Dipstick	16			
		*2+	1	6.25%
		*3+	15	93.75%
UA-3		* Denotes Correct Response		
PG: Siemens Multistix - Visual	194			
		*Large (3+)	182	93.81%
		*Moderate (2+)	9	4.64%
		Negative	2	1.03%
		Small (1+)	1	0.52%
UA-3		* Denotes Correct Response		
PG: Diag Test Group Clarity Urochek	24			
		*Large (3+)	23	95.83%
		Negative	1	4.17%
UA-3		* Denotes Correct Response		
PG: Roche Chemstrips - Visual	63			
		*Large (3+)	53	84.13%
		*Moderate (2+)	9	14.29%
		Negative	1	1.59%
UA-3		* Denotes Correct Response		
PG: Roche Urisys	30			
		*Large (3+)	25	83.33%
		*Moderate (2+)	5	16.67%
UA-3		* Denotes Correct Response		
PG: Roche Urisys / Roche - Conventional / Roche Urine Dipstick	2			
		*250 Ery/uL	2	100.00%
UA-3		* Denotes Correct Response		



PG: Roche Urisys / Roche - Arbitrary / Roche Urine Arbitrary Result	1			
		*++	1	100.00%
UA-3		* Denotes Correct Response		
PG: Siemens Clinitek Advantus	25			
		*Large (3+)	23	92.00%
		*Moderate (2+)	2	8.00%
UA-3		* Denotes Correct Response		
PG: 120 Urine Analyzer	44			
		*Large (3+)	41	93.18%
		*Moderate (2+)	2	4.55%
		Negative	1	2.27%
UA-3		* Denotes Correct Response		
SG: Select Siemens Clinitek/Visual	984			
		*Large (3+)	814	82.72%
		*Moderate (2+)	161	16.36%
		Negative	5	0.51%
		Small (1+)	2	0.20%
		Trace	2	0.20%
UA-3		* Denotes Correct Response		
SG: Select Roche Urinalysis Inst/Visual	99			
		*Large (3+)	84	84.85%
		*Moderate (2+)	14	14.14%
		Negative	1	1.01%
UA-3		* Denotes Correct Response		
SG: Select Consult Diag Urinalysis Inst/Visual	21			
		*Large (3+)	18	85.71%
		*Moderate (2+)	3	14.29%
UA-3		* Denotes Correct Response		
SG: Select McKesson Urinalysis Inst/Visual	51			
		*Large (3+)	47	92.16%
		*Moderate (2+)	3	5.88%
		Negative	1	1.96%
UA-3		* Denotes Correct Response		
SG: Select Henry Schein Urinalysis Inst/Visual	9			
		*Large (3+)	6	66.67%
		*Moderate (2+)	3	33.33%
UA-3		* Denotes Correct Response		
SG: Select Medi-Lab Performance UA Inst/Visual	8			
		*Large (3+)	8	100.00%
UA-3		* Denotes Correct Response		
MG: Method Group	1218			
		*Large (3+)	1017	83.50%
		*Moderate (2+)	189	15.52%
		Negative	8	0.66%
		Small (1+)	2	0.16%
		Trace	2	0.16%
UA-3		* Denotes Correct Response		
MG: Method Group	2			
		*250 Ery/uL	2	100.00%
UA-3		* Denotes Correct Response		
MG: Siemens - Plus Sys	17			
		*2+	1	5.88%
		*3+	16	94.12%
UA-3		* Denotes Correct Response		
MG: Siemens - Plus Sys	1			
		*++	1	100.00%
		Urine Glucose (Reducing Sub)		
UA-3		* Denotes Correct Response		
PG: Siemens Clinitek 10 / 100	11			
		*Negative	11	100.00%
UA-3		* Denotes Correct Response		
PG: Siemens Clinitek 50	63			
		*Negative	63	100.00%
UA-3		* Denotes Correct Response		
PG: Siemens Clinitek Status	684			
		*Negative	682	99.71%
		100 mg/dL	1	0.15%



		250 mg/dL	1	0.15%
UA-3		* Denotes Correct Response		
PG: Siemens Clinitek Status / Siemens - Plus Sys / Siemens Urine Dipstick	16			
		*Negative	16	100.00%
UA-3		* Denotes Correct Response		
PG: Siemens Multistix - Visual	196			
		*Negative	192	97.96%
		100 mg/dL	4	2.04%
UA-3		* Denotes Correct Response		
PG: Diag Test Group Clarity Urochek	24			
		*Negative	24	100.00%
UA-3		* Denotes Correct Response		
PG: Roche Chemstrips - Visual	63			
		*Negative	63	100.00%
UA-3		* Denotes Correct Response		
PG: Roche Urisys	30			
		*Negative	29	96.67%
		250 mg/dL	1	3.33%
UA-3		* Denotes Correct Response		
PG: Roche Urisys / Roche - Conventional / Roche Urine Dipstick	2			
		*Normal	2	100.00%
UA-3		* Denotes Correct Response		
PG: Roche Urisys / Roche - Arbitrary / Roche Urine Arbitrary Result	1			
		*Normal	1	100.00%
UA-3		* Denotes Correct Response		
PG: Siemens Clinitek Advantus	25			
		*Negative	25	100.00%
UA-3		* Denotes Correct Response		
PG: 120 Urine Analyzer	44			
		*Negative	44	100.00%
UA-3		* Denotes Correct Response		
SG: Select Siemens Clinitek/Visual	989			
		*Negative	983	99.39%
		100 mg/dL	5	0.51%
		250 mg/dL	1	0.10%
UA-3		* Denotes Correct Response		
SG: Select Roche Urinalysis Inst/Visual	99			
		*Negative	98	98.99%
		250 mg/dL	1	1.01%
UA-3		* Denotes Correct Response		
SG: Select Consult Diag Urinalysis Inst/Visual	21			
		*Negative	21	100.00%
UA-3		* Denotes Correct Response		
SG: Select McKesson Urinalysis Inst/Visual	51			
		*Negative	51	100.00%
UA-3		* Denotes Correct Response		
SG: Select Henry Schein Urinalysis Inst/Visual	9			
		*Negative	9	100.00%
UA-3		* Denotes Correct Response		
SG: Select Medi-Lab Performance UA Inst/Visual	8			
		*Negative	8	100.00%
UA-3		* Denotes Correct Response		
MG: Method Group	1223			
		*Negative	1216	99.43%
		100 mg/dL	5	0.41%
		250 mg/dL	2	0.16%
UA-3		* Denotes Correct Response		
MG: Method Group	2			
		*Normal	2	100.00%
UA-3		* Denotes Correct Response		
MG: Siemens - Plus Sys	17			
		*Negative	17	100.00%
UA-3		* Denotes Correct Response		
MG: Siemens - Plus Sys	1			
		*Normal	1	100.00%
		Urine Leuk Esterase		



UA-3		* Denotes Correct Response		
PG: Siemens Clinitek 10 / 100	11			
		*Negative	11	100.00%
UA-3		* Denotes Correct Response		
PG: Siemens Clinitek 50	63			
		*Negative	62	98.41%
		Trace	1	1.59%
UA-3		* Denotes Correct Response		
PG: Siemens Clinitek Status	685			
		*Negative	682	99.56%
		Small (1+)	3	0.44%
UA-3		* Denotes Correct Response		
PG: Siemens Clinitek Status / Siemens - Plus Sys / Siemens Urine Dipstick	16			
		*Negative	16	100.00%
UA-3		* Denotes Correct Response		
PG: Siemens Multistix - Visual	195			
		*Negative	191	97.95%
		Small (1+)	1	0.51%
		Trace	3	1.54%
UA-3		* Denotes Correct Response		
PG: Diag Test Group Clarity Urochek	24			
		*Negative	24	100.00%
UA-3		* Denotes Correct Response		
PG: Roche Chemstrips - Visual	63			
		*Negative	57	90.48%
		Small (1+)	2	3.17%
		Trace	4	6.35%
UA-3		* Denotes Correct Response		
PG: Roche Urisys	30			
		*Negative	28	93.33%
		Small (1+)	2	6.67%
UA-3		* Denotes Correct Response		
PG: Roche Urisys / Roche - Conventional / Roche Urine Dipstick	2			
		*Negative	2	100.00%
UA-3		* Denotes Correct Response		
PG: Roche Urisys / Roche - Arbitrary / Roche Urine Arbitrary Result	1			
		*Negative	1	100.00%
UA-3		* Denotes Correct Response		
PG: Siemens Clinitek Advantus	25			
		*Negative	25	100.00%
UA-3		* Denotes Correct Response		
PG: 120 Urine Analyzer	44			
		*Negative	44	100.00%
UA-3		* Denotes Correct Response		
SG: Select Siemens Clinitek/Visual	986			
		*Negative	978	99.19%
		Small (1+)	4	0.41%
		Trace	4	0.41%
UA-3		* Denotes Correct Response		
SG: Select Roche Urinalysis Inst/Visual	99			
		*Negative	90	90.91%
		Small (1+)	5	5.05%
		Trace	4	4.04%
UA-3		* Denotes Correct Response		
SG: Select Consult Diag Urinalysis Inst/Visual	21			
		*Negative	20	95.24%
		Small (1+)	1	4.76%
UA-3		* Denotes Correct Response		
SG: Select McKesson Urinalysis Inst/Visual	51			
		*Negative	49	96.08%
		Small (1+)	1	1.96%
		Trace	1	1.96%
UA-3		* Denotes Correct Response		
SG: Select Henry Schein Urinalysis Inst/Visual	9			
		*Negative	9	100.00%
UA-3		* Denotes Correct Response		



SG: Select Medi-Lab Performance UA Inst/Visual	8			
		*Negative	8	100.00%
UA-3		* Denotes Correct Response		
MG: Method Group	1220			
		*Negative	1199	98.28%
		Large (3+)	1	0.08%
		Small (1+)	11	0.90%
		Trace	9	0.74%
UA-3		* Denotes Correct Response		
MG: Method Group	2			
		*Negative	2	100.00%
UA-3		* Denotes Correct Response		
MG: Siemens - Plus Sys	17			
		*Negative	17	100.00%
UA-3		* Denotes Correct Response		
MG: Siemens - Plus Sys	1			
		*Negative	1	100.00%
		Urine Nitrite		
UA-3		* Denotes Correct Response		
PG: Siemens Clinitek 10 / 100	11			
		*Negative	11	100.00%
UA-3		* Denotes Correct Response		
PG: Siemens Clinitek 50	63			
		*Negative	63	100.00%
UA-3		* Denotes Correct Response		
PG: Siemens Clinitek Status	686			
		*Negative	685	99.85%
		Positive	1	0.15%
UA-3		* Denotes Correct Response		
PG: Siemens Clinitek Status / Siemens - Plus Sys / Siemens Urine Dipstick	16			
		*Negative	16	100.00%
UA-3		* Denotes Correct Response		
PG: Siemens Multistix - Visual	194			
		*Negative	191	98.45%
		Positive	3	1.55%
UA-3		* Denotes Correct Response		
PG: Diag Test Group Clarity Urochek	23			
		*Negative	23	100.00%
UA-3		* Denotes Correct Response		
PG: Roche Chemstrips - Visual	63			
		*Negative	62	98.41%
		Positive	1	1.59%
UA-3		* Denotes Correct Response		
PG: Roche Urisys	30			
		*Negative	29	96.67%
		Positive	1	3.33%
UA-3		* Denotes Correct Response		
PG: Roche Urisys / Roche - Conventional / Roche Urine Dipstick	2			
		*Negative	2	100.00%
UA-3		* Denotes Correct Response		
PG: Roche Urisys / Roche - Arbitrary / Roche Urine Arbitrary Result	1			
		*Negative	1	100.00%
UA-3		* Denotes Correct Response		
PG: Siemens Clinitek Advantus	25			
		*Negative	25	100.00%
UA-3		* Denotes Correct Response		
PG: 120 Urine Analyzer	44			
		*Negative	44	100.00%
UA-3		* Denotes Correct Response		
SG: Select Siemens Clinitek/Visual	986			
		*Negative	982	99.59%
		Positive	4	0.41%
UA-3		* Denotes Correct Response		
SG: Select Roche Urinalysis Inst/Visual	99			
		*Negative	97	97.98%
		Positive	2	2.02%



UA-3		* Denotes Correct Response		
SG: Select Consult Diag Urinalysis Inst/Visual	21			
		*Negative	20	95.24%
		Positive	1	4.76%
UA-3		* Denotes Correct Response		
SG: Select McKesson Urinalysis Inst/Visual	51			
		*Negative	51	100.00%
UA-3		* Denotes Correct Response		
SG: Select Henry Schein Urinalysis Inst/Visual	9			
		*Negative	9	100.00%
UA-3		* Denotes Correct Response		
SG: Select Medi-Lab Performance UA Inst/Visual	8			
		*Negative	8	100.00%
UA-3		* Denotes Correct Response		
MG: Method Group	1219			
		*Negative	1212	99.43%
		Positive	7	0.57%
UA-3		* Denotes Correct Response		
MG: Method Group	2			
		*Negative	2	100.00%
UA-3		* Denotes Correct Response		
MG: Siemens - Plus Sys	17			
		*Negative	17	100.00%
UA-3		* Denotes Correct Response		
MG: Siemens - Plus Sys	1			
		*Negative	1	100.00%
UA-3		Urine Urobilinogen		
		* Denotes Correct Response		
PG: Siemens Clinitek 10 / 100	11			
		*≥/ = 8 mg/dL	1	9.09%
		*1 mg/dL	2	18.18%
		*2.0 mg/dL	3	27.27%
		*4 mg/dL	3	27.27%
		0.2/Normal	2	18.18%
UA-3		* Denotes Correct Response		
PG: Siemens Clinitek Status / Siemens - Plus Sys / Siemens Urine Dipstick	7			
		*16 umol/L	1	14.29%
		*66 umol/L	5	71.43%
		3.2 umol/L	1	14.29%
UA-3		* Denotes Correct Response		
PG: Siemens Clinitek Status / Siemens - Plus Sys / Siemens Urine Urobilinogen	8			
		*2.0 EU/dL	1	12.50%
		*4.0 EU/dL	7	87.50%
UA-3		* Denotes Correct Response		
PG: Cole Taylor Strips	1			
		*0.2/Normal	1	100.00%
UA-3		* Denotes Correct Response		
PG: Germaine Aimstick	4			
		*0.2/Normal	3	75.00%
		*2.0 mg/dL	1	25.00%
UA-3		* Denotes Correct Response		
PG: Hypoguard DiaScreen - Visual	1			
		*≥/ = 8 mg/dL	1	100.00%
UA-3		* Denotes Correct Response		
PG: Iris Diagnostics vChem Strips	1			
		*2.0 mg/dL	1	100.00%
UA-3		* Denotes Correct Response		
PG: Biosys Uriscan Optima II Urine Analyzer	1			
		*0.2/Normal	1	100.00%
UA-3		* Denotes Correct Response		
PG: Quidel QuickVue UrinChek - Visual	1			
		*≥/ = 8 mg/dL	1	100.00%
UA-3		* Denotes Correct Response		
PG: Diag Test Group Clarity Urochek	24			
		*0.2/Normal	4	16.67%
		*1 mg/dL	7	29.17%



		*2.0 mg/dL	4	16.67%
		*4 mg/dL	7	29.17%
		>/= 8 mg/dL	2	8.33%
UA-3		* Denotes Correct Response		
PG: Roche Urisys / Roche - Conventional / Roche Urine Dipstick	2			
		*4 mg/dL	1	50.00%
		*Normal	1	50.00%
UA-3		* Denotes Correct Response		
PG: Roche Urisys / Roche - Arbitrary / Roche Urine Arbitrary Result	1			
		*+	1	100.00%
UA-3		* Denotes Correct Response		
PG: UriScan Reagent Strips	6			
		*0.2/Normal	3	50.00%
		*1 mg/dL	2	33.33%
		*4 mg/dL	1	16.67%
UA-3		* Denotes Correct Response		
PG: 120 Urine Analyzer	43			
		*0.2/Normal	28	65.12%
		*1 mg/dL	14	32.56%
		*2.0 mg/dL	1	2.33%
UA-3		* Denotes Correct Response		
PG: Medline Urinalysis Reagent Strips	2			
		*2.0 mg/dL	1	50.00%
		*4 mg/dL	1	50.00%
UA-3		* Denotes Correct Response		
PG: Moore Medical Urine Reagent Strips	1			
		*0.2/Normal	1	100.00%
UA-3		* Denotes Correct Response		
PG: LW Scientific Urine Reagent Strip	1			
		*4 mg/dL	1	100.00%
UA-3		* Denotes Correct Response		
SG: Select Siemens Clinitek/Visual	969			
		*>/= 8 mg/dL	53	5.47%
		*1 mg/dL	18	1.86%
		*2.0 mg/dL	120	12.38%
		*4 mg/dL	752	77.61%
		0.2/Normal	26	2.68%
		1 mg/dL	18	1.86%
UA-3		* Denotes Correct Response		
SG: Select Siemens Clinitek/Visual / Siemens - Plus Sys / Siemens Urine Dipstick	8			
		*16 umol/L	1	12.50%
		*33 umol/L	1	12.50%
		*66 umol/L	5	62.50%
		3.2 umol/L	1	12.50%
UA-3		* Denotes Correct Response		
SG: Select Roche Urinalysis Inst/Visual	93			
		*0.2/Normal	22	23.66%
		*1 mg/dL	51	54.84%
		*2.0 mg/dL	3	3.23%
		*4 mg/dL	13	13.98%
		>/= 8 mg/dL	4	4.30%
UA-3		* Denotes Correct Response		
SG: Select Consult Diag Urinalysis Inst/Visual	21			
		*0.2/Normal	8	38.10%
		*1 mg/dL	4	19.05%
		*2.0 mg/dL	4	19.05%
		*4 mg/dL	5	23.81%
UA-3		* Denotes Correct Response		
SG: Select McKesson Urinalysis Inst/Visual	50			
		*0.2/Normal	30	60.00%
		*1 mg/dL	16	32.00%
		*2.0 mg/dL	3	6.00%
		>/= 8 mg/dL	1	2.00%
UA-3		* Denotes Correct Response		
SG: Select Henry Schein Urinalysis Inst/Visual	8			
		*0.2/Normal	6	75.00%
		*4 mg/dL	2	25.00%



UA-3		* Denotes Correct Response		
SG: Select Medi-Lab Performance UA Inst/Visual	8			
		*0.2/Normal	3	37.50%
		*2.0 mg/dL	5	62.50%
UA-3		* Denotes Correct Response		
PG: URS - Urine Reagent Strips	3			
		*0.2/Normal	1	33.33%
		*1 mg/dL	1	33.33%
		*2.0 mg/dL	1	33.33%
		Urine Microalbumin		
UA-3		* Denotes Correct Response		
PG: Siemens Clinitek 50	9			
		*150 mg/L	9	100.00%
UA-3		* Denotes Correct Response		
PG: Siemens Clinitek Status	59			
		*150 mg/L	58	98.31%
		Negative	1	1.69%
UA-3		* Denotes Correct Response		
PG: Roche Chemstrip Micral	22			
		*100 mg/L (3+)	17	77.27%
		*150 mg/L	3	13.64%
		*50 mg/L (2+)	2	9.09%
UA-3		* Denotes Correct Response		
SG: Select Siemens Clinitek/Visual	72			
		*100 mg/L (3+)	1	1.39%
		*150 mg/L	69	95.83%
		20/30 mg/L (Small	1	1.39%
		Negative	1	1.39%
UA-3		* Denotes Correct Response		
MG: Method Group	96			
		*100 mg/L (3+)	19	19.79%
		*150 mg/L	73	76.04%
		*50 mg/L (2+)	2	2.08%
		20/30 mg/L (Small	1	1.04%
		Negative	1	1.04%
		Fluid pH		
UA-3		* Denotes Correct Response		
PG: Nitrazine Paper	14			
		*7	4	28.57%
		*7.5	7	50.00%
		*8 or greater	3	21.43%
		Urine hCG		
UA-3		* Denotes Correct Response		
AG: All Methods	1089			
		*Positive	1079	99.08%
		Negative	10	0.92%
ASI ProPhase Plus	1	*Positive	1	100.00%
Siemens Clinitek 50	2	*Positive	2	100.00%
Siemens Clinitek Status	13	*Positive	13	100.00%
Beck Coulter ICON 20 HCG	6	*Positive	6	100.00%
Beck Coulter ICON 25 hCG	70	*Positive	70	100.00%
Alere Acceava HCG Basic II	24	*Positive	24	100.00%
Alere Acceava hCG Combo	5	*Positive	5	100.00%
Alere Acceava hCG-Urine	20	*Positive	20	100.00%
Conception Tech One Step hCG	1	*Positive	1	100.00%
Diag Test Group Clarity hCG Plus Urine Cassettes	6	*Positive	6	100.00%
Fisher HC Sure-Vue hCG STAT	4	*Positive	4	100.00%
Fisher HC Sure-Vue Urine hCG	23	*Positive	23	100.00%
Fisher HC Sure-Vure U/S hCG	5	*Positive	5	100.00%
Sekisui Diagnostics OSOM Urine	11	*Positive	11	100.00%
Sekisui Diagnostics OSOM Card Pregnancy	30	*Positive	30	100.00%
Sekisui Diagnostics OSOM hCG Combo	40	*Positive	40	100.00%
Henry Schein One Step + (hCG)	56	*Positive	55	98.21%
Immunostics Detector Combo	3	*Positive	3	100.00%
Immunostics hCG Detector-urine	2	*Positive	2	100.00%
Instant Tech. i Pregnancy	3	*Positive	3	100.00%
Inverness Clearview Combo II	1	*Positive	1	100.00%
Inv Clearview Urine - Cassette	3	*Positive	2	66.67%



Inverness Signify hCG combo	1	*Positive	1	100.00%
Inverness Signify hCG urine	3	*Positive	3	100.00%
LifeSign Status hCG	1	*Positive	1	100.00%
McKesson Medi-Lab hCG - Urine	32	*Positive	32	100.00%
Meridian Immunocard Stat!	1	*Positive	1	100.00%
Permaxim RediScreen	4	*Positive	4	100.00%
Quidel QuickVue + hCG	20	*Positive	20	100.00%
Quidel QuickVue One-Step hCG Combo U/S	55	*Positive	54	98.18%
Quidel QuickVue One-Step Urine	141	*Positive	140	99.29%
Quidel RapidVue hCG - Urine	6	*Positive	6	100.00%
SAS Ultra HCG	2	*Positive	2	100.00%
Select Diag hCG Cassette (urine)	19	*Positive	19	100.00%
Select Diag hCG Combo	10	*Positive	9	90.00%
Select Diag hCG Dipstick	7	*Positive	7	100.00%
SP (Cardinal) Brand combo	24	*Positive	24	100.00%
SP (Cardinal) Brand-cassette	29	*Positive	27	93.10%
SP Brand (Cardinal) HCG Combo	64	*Positive	64	100.00%
Stanbio Qupid	5	*Positive	5	100.00%
Stanbio QuStick	1	*Positive	1	100.00%
Stanbio TRUE hCG	5	*Positive	4	80.00%
MooreBrand hCG Cassette	2	*Positive	2	100.00%
Nova Scientific Century S/U hCG	2	*Positive	2	100.00%
Diag Test Group Clarity hCG Urine Strips	1	*Positive	1	100.00%
Jant AccuTest	4	*Positive	4	100.00%
Becton Dickinson Precise One-Step	1	*Positive	1	100.00%
Consult Diag hCG Urine Cassette	141	*Positive	141	100.00%
CLIAwaived Inc Pregnancy Urine Dip-Strip	4	*Positive	4	100.00%
Germaine AIM Step Combo Pregnancy	1	*Positive	1	100.00%
NDC ProAdvantage Urine hCG Cassette	5	*Positive	5	100.00%
Consult Diag hCG Urine Dipstick	52	*Positive	52	100.00%
SP (Cardinal) Brand Dipstick	12	*Positive	11	91.67%
Consult Diag hCG Combo Cassette (Urine)	11	*Positive	11	100.00%
McKesson Medi-Lab Performance hCG Urine Cassette	7	*Positive	7	100.00%
Alere hCG Combo Cassette (25 mIU/mL)	2	*Positive	2	100.00%
Alere hCG Cassette (25 mIU/mL)	12	*Positive	12	100.00%
Alere hCG Combo Cassette (20/10 mIU/mL)	1	*Positive	1	100.00%
Alere hCG Dipstick	7	*Positive	7	100.00%
Medline Pregnancy Urine Dip-Strip	2	*Positive	2	100.00%
Medline Pregnancy Urine Cassette Test	4	*Positive	4	100.00%
McKesson hCG Urine Test Dipstick	3	*Positive	3	100.00%
Medline Pregnancy Combo Cassette (urine)	1	*Positive	1	100.00%
Quidel QuickVue hCG Urine Test	20	*Positive	20	100.00%
MooreBrand hCG Dipstick	1	*Positive	1	100.00%
Siemens Clinitest hCG Pregnancy Test	2	*Positive	2	100.00%
Henry Schein One Step + HCG Combo	3	*Positive	2	66.67%
Henry Schein One Step hCG Urine Cassette	26	*Positive	26	100.00%
Medline hCG Test Strip 25 mIU/mL	1	*Positive	1	100.00%
NCS hCG One Step Pregnancy Test	1	*Positive	1	100.00%
MediChoice Urine hCG Cassette	1	*Positive	1	100.00%
Henry Schein One Step + (hCG)	56	Negative	1	1.79%
Inv Clearview Urine - Cassette	3	Negative	1	33.33%
Quidel QuickVue One-Step hCG Combo U/S	55	Negative	1	1.82%
Quidel QuickVue One-Step Urine	141	Negative	1	0.71%
Select Diag hCG Combo	10	Negative	1	10.00%
SP (Cardinal) Brand-cassette	29	Negative	2	6.90%
Stanbio TRUE hCG	5	Negative	1	20.00%
SP (Cardinal) Brand Dipstick	12	Negative	1	8.33%
Henry Schein One Step + HCG Combo	3	Negative	1	33.33%
		Urine hCG (non-waived)		
UA-3		* Denotes Correct Response		
PG: Quidel Sofia/Sofia 2 hCG (urine)	1			
		*Positive	1	100.00%
940 - Urine Adulteration	Labs	Response	Number	Percent
		Adulteraton - Creatinine (Semi-quant)		
UAD-5		* Denotes Correct Response		
PG: Confirm Biosciences T Cup	2			
		*50 mg/dL	1	50.00%
		100 mg/dL	1	50.00%
UAD-6		* Denotes Correct Response		



PG: Confirm Biosciences T Cup	2			
		*Negative	1	50.00%
		100 mg/dL	1	50.00%
		Adulteration - pH (Semi-quant)		
UAD-5		* Denotes Correct Response		
PG: Confirm Biosciences T Cup	2			
		*3	1	50.00%
		7	1	50.00%
UAD-6		* Denotes Correct Response		
PG: Confirm Biosciences T Cup	2			
		*7	2	100.00%
		Adulteration - SG (Semi-quant)		
UAD-5		* Denotes Correct Response		
PG: Confirm Biosciences T Cup	2			
		*1.005	2	100.00%
UAD-6		* Denotes Correct Response		
PG: Confirm Biosciences T Cup	2			
		*1.000	1	50.00%
		1.015	1	50.00%
661 - Urine Colony Count;784 - Presump ID/Colony Count - Urine;775 - Urine Culture;778 - Urine & Throat Cultures;779 - Urine & Genital Cultures;781 - Urine, Throat & Genital Cultures;782 - Complete Bacteriology;790 - Urine Susceptibility Add-On	Labs	Response	Number	Percent
		Susceptibility Cefazolin		
UC-11		* Denotes Correct Response		
AG: All Methods	1			
		*Resistant	1	100.00%
Susceptibility	1	*Resistant	1	100.00%
		Susceptibility Cefepime		
UC-11		* Denotes Correct Response		
AG: All Methods	1			
		*Sensitive	1	100.00%
Susceptibility	1	*Sensitive	1	100.00%
		Susceptibility Cefpodoxime		
UC-11		* Denotes Correct Response		
AG: All Methods	1			
		*Sensitive	1	100.00%
Susceptibility	1	*Sensitive	1	100.00%
		Susceptibility Ceftriaxone		
UC-11		* Denotes Correct Response		
AG: All Methods	1			
		*Sensitive	1	100.00%
Susceptibility	1	*Sensitive	1	100.00%
		Susceptibility Ceftriaxone		
UC-11		* Denotes Correct Response		
AG: All Methods	3			
		*Sensitive	3	100.00%
Susceptibility	3	*Sensitive	3	100.00%
		Susceptibility Cephalothin		
UC-11		* Denotes Correct Response		
AG: All Methods	4			
		*Resistant	4	100.00%
Susceptibility	4	*Resistant	4	100.00%
		Susceptibility Ciprofloxacin		
UC-11		* Denotes Correct Response		
AG: All Methods	6			
		*Sensitive	6	100.00%
Susceptibility	6	*Sensitive	6	100.00%
		Susceptibility Gentamicin		
UC-11		* Denotes Correct Response		
AG: All Methods	4			
		*Sensitive	4	100.00%
Susceptibility	4	*Sensitive	4	100.00%
		Susceptibility Imipenem		
UC-11		* Denotes Correct Response		
AG: All Methods	2			
		*Sensitive	2	100.00%



Susceptibility	2	*Sensitive	2	100.00%
		Susceptibility Levofloxacin		
UC-11		* Denotes Correct Response		
AG: All Methods	3			
		*Sensitive	3	100.00%
Susceptibility	3	*Sensitive	3	100.00%
		Susceptibility Minocycline		
UC-11		* Denotes Correct Response		
AG: All Methods	1			
		*Sensitive	1	100.00%
Susceptibility	1	*Sensitive	1	100.00%
		Susceptibility Piperacillin		
UC-11		* Denotes Correct Response		
AG: All Methods	1			
		*Sensitive	1	100.00%
Susceptibility	1	*Sensitive	1	100.00%
		Susceptibility Piperacillin/Tazobactam		
UC-11		* Denotes Correct Response		
AG: All Methods	1			
		*Sensitive	1	100.00%
Susceptibility	1	*Sensitive	1	100.00%
		Susceptibility Sulfamethoxazole		
UC-11		* Denotes Correct Response		
AG: All Methods	1			
		*Sensitive	1	100.00%
Susceptibility	1	*Sensitive	1	100.00%
		Susceptibility Tetracycline		
UC-11		* Denotes Correct Response		
AG: All Methods	4			
		*Intermediate	1	25.00%
		*Sensitive	3	75.00%
Susceptibility	4	*Intermediate	1	25.00%
Susceptibility		*Sensitive	3	75.00%
		Susceptibility Tobramycin		
UC-11		* Denotes Correct Response		
AG: All Methods	2			
		*Sensitive	2	100.00%
Susceptibility	2	*Sensitive	2	100.00%
		Susceptibility Trimethoprim		
UC-11		* Denotes Correct Response		
AG: All Methods	1			
		*Sensitive	1	100.00%
Susceptibility	1	*Sensitive	1	100.00%
		Susceptibility Trimeth./Sulfameth.		
UC-11		* Denotes Correct Response		
AG: All Methods	4			
		*Sensitive	4	100.00%
Susceptibility	4	*Sensitive	4	100.00%
		Urine Culture		
UC-11		* Denotes Correct Response		
AG: All Methods	27			
		*Enterobacter cloacae	1	3.70%
		*Gram Negative organism	5	18.52%
		*Growth, referred for identification	13	48.15%
		*Presumptive Enterobacter species	2	7.41%
		*Uricult - Growth (Positive)	6	22.22%
Cultures	27	*Enterobacter cloacae	1	3.70%
Cultures		*Gram Negative organism	5	18.52%
Cultures		*Growth, referred for identification	13	48.15%
Cultures		*Presumptive Enterobacter species	2	7.41%
Cultures		*Uricult - Growth (Positive)	6	22.22%
UC-12		* Denotes Correct Response		
AG: All Methods	27			
		*Gram Negative organism	4	14.81%
		*Growth, referred for identification	13	48.15%
		*Klebsiella oxytoca	1	3.70%
		*Presumptive Klebsiella species	1	3.70%



		*Uricult - Growth (Positive)	5	18.52%
		Escherichia coli	1	3.70%
		Presumptive E. coli	1	3.70%
		Uricult - No Growth (Negative)	1	3.70%
Cultures	27	*Gram Negative organism	4	14.81%
Cultures		*Growth, referred for identification	13	48.15%
Cultures		*Klebsiella oxytoca	1	3.70%
Cultures		*Presumptive Klebsiella species	1	3.70%
Cultures		*Uricult - Growth (Positive)	5	18.52%
Cultures		Escherichia coli	1	3.70%
Cultures		Presumptive E. coli	1	3.70%
Cultures		Uricult - No Growth (Negative)	1	3.70%
UC-13		* Denotes Correct Response		
AG: All Methods	26			
		*Enterococcus faecalis	1	3.85%
		*Gram Positive organism	3	11.54%
		*Growth, referred for identification	12	46.15%
		*Presumptive Enterococcus species	1	3.85%
		*Uricult - Growth (Positive)	4	15.38%
		No Growth	3	11.54%
		Uricult - No Growth (Negative)	2	7.69%
Cultures	26	*Enterococcus faecalis	1	3.85%
Cultures		*Gram Positive organism	3	11.54%
Cultures		*Growth, referred for identification	12	46.15%
Cultures		*Presumptive Enterococcus species	1	3.85%
Cultures		*Uricult - Growth (Positive)	4	15.38%
Cultures		No Growth	3	11.54%
Cultures		Uricult - No Growth (Negative)	2	7.69%
UC-14		* Denotes Correct Response		
AG: All Methods	11			
		*Escherichia coli	1	9.09%
		*Gram Negative organism	3	27.27%
		*Growth, referred for identification	2	18.18%
		*Presumptive E. coli	1	9.09%
		*Uricult - Growth (Positive)	4	36.36%
Cultures	11	*Escherichia coli	1	9.09%
Cultures		*Gram Negative organism	3	27.27%
Cultures		*Growth, referred for identification	2	18.18%
Cultures		*Presumptive E. coli	1	9.09%
Cultures		*Uricult - Growth (Positive)	4	36.36%
UC-15		* Denotes Correct Response		
AG: All Methods	11			
		*Gram Negative organism	3	27.27%
		*Growth, referred for identification	2	18.18%
		*Morganella morganii	1	9.09%
		*Uricult - Growth (Positive)	4	36.36%
		Presumptive Proteus species	1	9.09%
Cultures	11	*Gram Negative organism	3	27.27%
Cultures		*Growth, referred for identification	2	18.18%
Cultures		*Morganella morganii	1	9.09%
Cultures		*Uricult - Growth (Positive)	4	36.36%
Cultures		Presumptive Proteus species	1	9.09%
		Susceptibility Amox/Clavulanate		
UC-11		* Denotes Correct Response		
AG: All Methods	3			
		*Resistant	3	100.00%
Susceptibility	3	*Resistant	3	100.00%
		Susceptibility Ampicillin		
UC-11		* Denotes Correct Response		
AG: All Methods	5			
		*Intermediate	1	20.00%
		*Resistant	4	80.00%
Susceptibility	5	*Intermediate	1	20.00%
Susceptibility		*Resistant	4	80.00%
		Susceptibility Ampicillin/Sulbactam		
UC-11		* Denotes Correct Response		
AG: All Methods	1			
		*Intermediate	1	100.00%



Susceptibility	1	*Intermediate	1	100.00%
703 - Urine Drug Screen 10mL;713 - Urine Drug Screen (50 ml sample volume)	Labs	Response	Number	Percent
		Barbiturates		
UDS-5		* Denotes Correct Response		
AG: All Methods	132			
		*Negative	130	98.48%
		Positive	2	1.52%
AlereTriage	4	*Negative	4	100.00%
E-Screen/Alere	25	*Negative	25	100.00%
Express Diag. Express Check	1	*Negative	1	100.00%
Alere iCup Drugs of Abuse	12	*Negative	12	100.00%
Instant Technologies i Screen	1	*Negative	1	100.00%
MedTox Profile II	2	*Negative	2	100.00%
Beckman AU	1	*Negative	1	100.00%
Innovacon Multi-Clin	5	*Negative	5	100.00%
CLIAwaived Inc Rapid Drug Cup	2	*Negative	2	100.00%
MedTox Profile V	2	*Negative	2	100.00%
Jant Accutest Multi Drug Screen	1	*Negative	1	100.00%
Amedica Biotech AmediCheck Instant Test Cup	4	*Negative	4	100.00%
Nova Century Scientific Nova Screen DOA Cassette	1	*Negative	1	100.00%
American Screening Corp Reveal Multi-Drug Test Cup	4	*Negative	4	100.00%
Mindray BS-200	2	*Negative	2	100.00%
Medica Easy RA	1	*Negative	1	100.00%
Alere iCassette	3	*Negative	3	100.00%
Roche Hitachi 717	1	*Negative	1	100.00%
Alere USD UScreen Cup	8	*Negative	8	100.00%
Alere iScreen Drugs of Abuse Card	1	*Negative	1	100.00%
Noble Medical Insta-Screen Test Cup	3	*Negative	3	100.00%
Hemosure First Sign UDS Dipstick	14	*Negative	14	100.00%
Siemens Viva/ProE System	1	*Negative	1	100.00%
Premier Biotech Bio-Cup	12	*Negative	12	100.00%
UCP Biosciences Home Drug Screening Test Cup	1	*Negative	1	100.00%
Abbott Architect c and ci Series	1	*Negative	1	100.00%
Consult Diag Drugs of Abuse Cups	2	*Negative	2	100.00%
Branan Med Corp ToxCup	5	*Negative	5	100.00%
American Screening Discover Plus Cup	1	*Negative	1	100.00%
McKesson Drugs of Abuse Test Cup	2	*Negative	1	50.00%
Henry Schein OneStep+ DOA	5	*Negative	5	100.00%
Siemens Atellica	1	*Negative	1	100.00%
Confirm Biosciences T Cup	1	*Negative	1	100.00%
Healgen Sci Drug Test Cup	1	*Negative	1	100.00%
ASC One Screen Drug Test Card	1	Positive	1	100.00%
McKesson Drugs of Abuse Test Cup	2	Positive	1	50.00%
UDS-6		* Denotes Correct Response		
AG: All Methods	132			
		*Positive	129	97.73%
		Negative	3	2.27%
AlereTriage	4	*Positive	4	100.00%
E-Screen/Alere	25	*Positive	25	100.00%
Express Diag. Express Check	1	*Positive	1	100.00%
Alere iCup Drugs of Abuse	12	*Positive	12	100.00%
Instant Technologies i Screen	1	*Positive	1	100.00%
MedTox Profile II	2	*Positive	2	100.00%
Beckman AU	1	*Positive	1	100.00%
Innovacon Multi-Clin	5	*Positive	5	100.00%
CLIAwaived Inc Rapid Drug Cup	2	*Positive	2	100.00%
MedTox Profile V	2	*Positive	2	100.00%
Jant Accutest Multi Drug Screen	1	*Positive	1	100.00%
Amedica Biotech AmediCheck Instant Test Cup	4	*Positive	4	100.00%
Nova Century Scientific Nova Screen DOA Cassette	1	*Positive	1	100.00%
American Screening Corp Reveal Multi-Drug Test Cup	4	*Positive	4	100.00%
Mindray BS-200	2	*Positive	2	100.00%
Medica Easy RA	1	*Positive	1	100.00%
Alere iCassette	3	*Positive	3	100.00%
Roche Hitachi 717	1	*Positive	1	100.00%
Alere USD UScreen Cup	8	*Positive	8	100.00%
Alere iScreen Drugs of Abuse Card	1	*Positive	1	100.00%
Noble Medical Insta-Screen Test Cup	3	*Positive	2	66.67%



Hemosure First Sign UDS Dipstick	14	*Positive	14	100.00%
Siemens Viva/ProE System	1	*Positive	1	100.00%
Premier Biotech Bio-Cup	12	*Positive	12	100.00%
UCP Biosciences Home Drug Screening Test Cup	1	*Positive	1	100.00%
Abbott Architect c and ci Series	1	*Positive	1	100.00%
Consult Diag Drugs of Abuse Cups	2	*Positive	2	100.00%
Branan Med Corp ToxCup	5	*Positive	5	100.00%
American Screening Discover Plus Cup	1	*Positive	1	100.00%
McKesson Drugs of Abuse Test Cup	2	*Positive	2	100.00%
Henry Schein OneStep+ DOA	5	*Positive	4	80.00%
Siemens Atellica	1	*Positive	1	100.00%
Confirm Biosciences T Cup	1	*Positive	1	100.00%
Healgen Sci Drug Test Cup	1	*Positive	1	100.00%
ASC One Screen Drug Test Card	1	Negative	1	100.00%
Noble Medical Insta-Screen Test Cup	3	Negative	1	33.33%
Henry Schein OneStep+ DOA	5	Negative	1	20.00%
		Benzodiazepines		
UDS-5		* Denotes Correct Response		
AG: All Methods	143			
		*Negative	143	100.00%
AlereTriage	4	*Negative	4	100.00%
Siemens Dim/AR/ES/RxL/Xpand/EXL	1	*Negative	1	100.00%
E-Screen/Alere	25	*Negative	25	100.00%
Express Diag. Express Check	1	*Negative	1	100.00%
Alere iCup Drugs of Abuse	15	*Negative	15	100.00%
Instant Technologies i Screen	1	*Negative	1	100.00%
MedTox Profile II	2	*Negative	2	100.00%
Noble Med Split Specimen Cup	1	*Negative	1	100.00%
Beckman AU	1	*Negative	1	100.00%
Innovacon Multi-Clin	5	*Negative	5	100.00%
CLIAwaived Inc Rapid Drug Cup	3	*Negative	3	100.00%
ASC One Screen Drug Test Card	1	*Negative	1	100.00%
MedTox Profile V	2	*Negative	2	100.00%
Jant Accutest Multi Drug Screen	1	*Negative	1	100.00%
Amedica Biotech AmediCheck Instant Test Cup	4	*Negative	4	100.00%
Nova Century Scientific Nova Screen DOA Cassette	1	*Negative	1	100.00%
American Screening Corp Reveal Multi-Drug Test Cup	4	*Negative	4	100.00%
Mindray BS-200	2	*Negative	2	100.00%
Medica Easy RA	1	*Negative	1	100.00%
Alere iCassette	3	*Negative	3	100.00%
Roche Hitachi 717	1	*Negative	1	100.00%
Alere USD UScreen Cup	8	*Negative	8	100.00%
TFS Indiko/Indiko Plus	3	*Negative	3	100.00%
Alere iScreen Drugs of Abuse Card	1	*Negative	1	100.00%
Noble Medical Insta-Screen Test Cup	3	*Negative	3	100.00%
Hemosure First Sign UDS Dipstick	14	*Negative	14	100.00%
Siemens Viva/ProE System	3	*Negative	3	100.00%
TFS DRI Assay	1	*Negative	1	100.00%
Premier Biotech Bio-Cup	12	*Negative	12	100.00%
UCP Biosciences Home Drug Screening Test Cup	1	*Negative	1	100.00%
Abbott Architect c and ci Series	1	*Negative	1	100.00%
Consult Diag Drugs of Abuse Cups	2	*Negative	2	100.00%
Branan Med Corp ToxCup	4	*Negative	4	100.00%
American Screening Discover Plus Cup	1	*Negative	1	100.00%
McKesson Drugs of Abuse Test Cup	2	*Negative	2	100.00%
Henry Schein OneStep+ DOA	5	*Negative	5	100.00%
Siemens Atellica	1	*Negative	1	100.00%
Confirm Biosciences T Cup	1	*Negative	1	100.00%
Healgen Sci Drug Test Cup	1	*Negative	1	100.00%
UDS-6		* Denotes Correct Response		
AG: All Methods	143			
		*Negative	143	100.00%
AlereTriage	4	*Negative	4	100.00%
Siemens Dim/AR/ES/RxL/Xpand/EXL	1	*Negative	1	100.00%
E-Screen/Alere	25	*Negative	25	100.00%
Express Diag. Express Check	1	*Negative	1	100.00%
Alere iCup Drugs of Abuse	15	*Negative	15	100.00%
Instant Technologies i Screen	1	*Negative	1	100.00%
MedTox Profile II	2	*Negative	2	100.00%



Noble Med Split Specimen Cup	1	*Negative	1	100.00%
Beckman AU	1	*Negative	1	100.00%
Innovacon Multi-Clin	5	*Negative	5	100.00%
CLIAwaived Inc Rapid Drug Cup	3	*Negative	3	100.00%
ASC One Screen Drug Test Card	1	*Negative	1	100.00%
MedTox Profile V	2	*Negative	2	100.00%
Jant Accutest Multi Drug Screen	1	*Negative	1	100.00%
Amedica Biotech AmediCheck Instant Test Cup	4	*Negative	4	100.00%
Nova Century Scientific Nova Screen DOA Cassette	1	*Negative	1	100.00%
American Screening Corp Reveal Multi-Drug Test Cup	4	*Negative	4	100.00%
Mindray BS-200	2	*Negative	2	100.00%
Medica Easy RA	1	*Negative	1	100.00%
Alere iCassette	3	*Negative	3	100.00%
Roche Hitachi 717	1	*Negative	1	100.00%
Alere USD UScreen Cup	8	*Negative	8	100.00%
TFS Indiko/Indiko Plus	3	*Negative	3	100.00%
Alere iScreen Drugs of Abuse Card	1	*Negative	1	100.00%
Noble Medical Insta-Screen Test Cup	3	*Negative	3	100.00%
Hemosure First Sign UDS Dipstick	14	*Negative	14	100.00%
Siemens Viva/ProE System	3	*Negative	3	100.00%
TFS DRI Assay	1	*Negative	1	100.00%
Premier Biotech Bio-Cup	12	*Negative	12	100.00%
UCP Biosciences Home Drug Screening Test Cup	1	*Negative	1	100.00%
Abbott Architect c and ci Series	1	*Negative	1	100.00%
Consult Diag Drugs of Abuse Cups	2	*Negative	2	100.00%
Branan Med Corp ToxCup	4	*Negative	4	100.00%
American Screening Discover Plus Cup	1	*Negative	1	100.00%
McKesson Drugs of Abuse Test Cup	2	*Negative	2	100.00%
Henry Schein OneStep+ DOA	5	*Negative	5	100.00%
Siemens Atellica	1	*Negative	1	100.00%
Confirm Biosciences T Cup	1	*Negative	1	100.00%
Healgen Sci Drug Test Cup	1	*Negative	1	100.00%
		Amphetamines		
UDS-5		* Denotes Correct Response		
AG: All Methods	149			
		*Negative	149	100.00%
AlereTriage	4	*Negative	4	100.00%
Siemens Dim/AR/ES/RxL/Xpand/EXL	1	*Negative	1	100.00%
E-Screen/Alere	22	*Negative	22	100.00%
Express Diag. Express Check	1	*Negative	1	100.00%
Alere iCup Drugs of Abuse	19	*Negative	19	100.00%
Instant Technologies i Screen	1	*Negative	1	100.00%
MedTox Profile II	3	*Negative	3	100.00%
Noble Med Split Specimen Cup	1	*Negative	1	100.00%
Beckman AU	1	*Negative	1	100.00%
Innovacon Multi-Clin	5	*Negative	5	100.00%
CLIAwaived Inc Rapid Drug Cup	3	*Negative	3	100.00%
ASC One Screen Drug Test Card	1	*Negative	1	100.00%
MedTox Profile V	2	*Negative	2	100.00%
Jant Accutest Multi Drug Screen	1	*Negative	1	100.00%
Amedica Biotech AmediCheck Instant Test Cup	4	*Negative	4	100.00%
Nova Century Scientific Nova Screen DOA Cassette	1	*Negative	1	100.00%
American Screening Corp Reveal Multi-Drug Test Cup	4	*Negative	4	100.00%
Mindray BS-200	2	*Negative	2	100.00%
Medica Easy RA	1	*Negative	1	100.00%
Alere iCassette	6	*Negative	6	100.00%
Roche Hitachi 717	1	*Negative	1	100.00%
Alere USD UScreen Cup	8	*Negative	8	100.00%
TFS Indiko/Indiko Plus	3	*Negative	3	100.00%
Alere iScreen Drugs of Abuse Card	1	*Negative	1	100.00%
Noble Medical Insta-Screen Test Cup	3	*Negative	3	100.00%
Hemosure First Sign UDS Dipstick	14	*Negative	14	100.00%
Siemens Viva/ProE System	4	*Negative	4	100.00%
TFS DRI Assay	1	*Negative	1	100.00%
Premier Biotech Bio-Cup	12	*Negative	12	100.00%
UCP Biosciences Home Drug Screening Test Cup	1	*Negative	1	100.00%
Abbott Architect c and ci Series	1	*Negative	1	100.00%
Consult Diag Drugs of Abuse Cups	2	*Negative	2	100.00%
Branan Med Corp ToxCup	5	*Negative	5	100.00%



American Screening Discover Plus Cup	1	*Negative	1	100.00%
McKesson Drugs of Abuse Test Cup	2	*Negative	2	100.00%
Henry Schein OneStep+ DOA	5	*Negative	5	100.00%
Siemens Atellica	1	*Negative	1	100.00%
Healgen Sci Drug Test Cup	1	*Negative	1	100.00%
UDS-6		* Denotes Correct Response		
AG: All Methods	147			
		*Negative	147	100.00%
AlereTriage	4	*Negative	4	100.00%
Siemens Dim/AR/ES/RxL/Xpand/EXL	1	*Negative	1	100.00%
E-Screen/Alere	22	*Negative	22	100.00%
Express Diag. Express Check	1	*Negative	1	100.00%
Alere iCup Drugs of Abuse	18	*Negative	18	100.00%
Instant Technologies i Screen	1	*Negative	1	100.00%
MedTox Profile II	3	*Negative	3	100.00%
Noble Med Split Specimen Cup	1	*Negative	1	100.00%
Beckman AU	1	*Negative	1	100.00%
Innovacon Multi-Clin	5	*Negative	5	100.00%
CLIAwaived Inc Rapid Drug Cup	3	*Negative	3	100.00%
ASC One Screen Drug Test Card	1	*Negative	1	100.00%
MedTox Profile V	2	*Negative	2	100.00%
Jant Accutest Multi Drug Screen	1	*Negative	1	100.00%
Amedica Biotech AmediCheck Instant Test Cup	4	*Negative	4	100.00%
Nova Century Scientific Nova Screen DOA Cassette	1	*Negative	1	100.00%
American Screening Corp Reveal Multi-Drug Test Cup	4	*Negative	4	100.00%
Mindray BS-200	2	*Negative	2	100.00%
Medica Easy RA	1	*Negative	1	100.00%
Alere iCassette	6	*Negative	6	100.00%
Roche Hitachi 717	1	*Negative	1	100.00%
Alere USD UScreen Cup	7	*Negative	7	100.00%
TFS Indiko/Indiko Plus	3	*Negative	3	100.00%
Alere iScreen Drugs of Abuse Card	1	*Negative	1	100.00%
Noble Medical Insta-Screen Test Cup	3	*Negative	3	100.00%
Hemosure First Sign UDS Dipstick	14	*Negative	14	100.00%
Siemens Viva/ProE System	4	*Negative	4	100.00%
TFS DRI Assay	1	*Negative	1	100.00%
Premier Biotech Bio-Cup	12	*Negative	12	100.00%
UCP Biosciences Home Drug Screening Test Cup	1	*Negative	1	100.00%
Abbott Architect c and ci Series	1	*Negative	1	100.00%
Consult Diag Drugs of Abuse Cups	2	*Negative	2	100.00%
Branan Med Corp ToxCup	5	*Negative	5	100.00%
American Screening Discover Plus Cup	1	*Negative	1	100.00%
McKesson Drugs of Abuse Test Cup	2	*Negative	2	100.00%
Henry Schein OneStep+ DOA	5	*Negative	5	100.00%
Siemens Atellica	1	*Negative	1	100.00%
Healgen Sci Drug Test Cup	1	*Negative	1	100.00%
		Cannabinoids (THC)		
UDS-5		* Denotes Correct Response		
AG: All Methods	177			
		*Negative	174	98.31%
		Positive	3	1.69%
AlereTriage	4	*Negative	4	100.00%
Siemens Dim/AR/ES/RxL/Xpand/EXL	1	*Negative	1	100.00%
E-Screen/Alere	49	*Negative	49	100.00%
Express Diag. Express Check	1	*Negative	1	100.00%
Alere iCup Drugs of Abuse	19	*Negative	19	100.00%
Instant Technologies i Screen	2	*Negative	2	100.00%
Integrated E-Z Split Key Cup	1	*Negative	1	100.00%
MedTox Profile II	3	*Negative	3	100.00%
Beckman AU	1	*Negative	1	100.00%
Innovacon Multi-Clin	5	*Negative	5	100.00%
CLIAwaived Inc Rapid Drug Cup	3	*Negative	3	100.00%
MedTox Profile V	2	*Negative	2	100.00%
Jant Accutest Multi Drug Screen	1	*Negative	1	100.00%
Amedica Biotech AmediCheck Instant Test Cup	4	*Negative	4	100.00%
Nova Century Scientific Nova Screen DOA Cassette	1	*Negative	1	100.00%
American Screening Corp Reveal Multi-Drug Test Cup	4	*Negative	4	100.00%
Mindray BS-200	2	*Negative	2	100.00%
Medica Easy RA	1	*Negative	1	100.00%



Alere iCassette	6	*Negative	6	100.00%
Roche Hitachi 717	1	*Negative	1	100.00%
Alere USD UScreen Cup	8	*Negative	8	100.00%
TFS Indiko/Indiko Plus	2	*Negative	2	100.00%
Alere iScreen Drugs of Abuse Card	1	*Negative	1	100.00%
Noble Medical Insta-Screen Test Cup	3	*Negative	3	100.00%
Hemosure First Sign UDS Dipstick	14	*Negative	14	100.00%
Siemens Viva/ProE System	3	*Negative	3	100.00%
TFS DRI Assay	1	*Negative	1	100.00%
Premier Biotech Bio-Cup	12	*Negative	12	100.00%
UCP Biosciences Home Drug Screening Test Cup	1	*Negative	1	100.00%
Abbott Architect c and ci Series	1	*Negative	1	100.00%
Consult Diag Drugs of Abuse Cups	2	*Negative	2	100.00%
Branan Med Corp ToxCup	5	*Negative	5	100.00%
American Screening Discover Plus Cup	1	*Negative	1	100.00%
McKesson Drugs of Abuse Test Cup	2	*Negative	1	50.00%
Henry Schein OneStep+ DOA	5	*Negative	5	100.00%
Siemens Atellica	1	*Negative	1	100.00%
Confirm Biosciences T Cup	1	*Negative	1	100.00%
Healgen Sci Drug Test Cup	1	*Negative	1	100.00%
Noble Med Split Specimen Cup	1	Positive	1	100.00%
ASC One Screen Drug Test Card	1	Positive	1	100.00%
McKesson Drugs of Abuse Test Cup	2	Positive	1	50.00%
UDS-6		* Denotes Correct Response		
AG: All Methods	177			
		*Positive	173	97.74%
		Negative	4	2.26%
AlereTriage	4	*Positive	4	100.00%
Siemens Dim/AR/ES/RxL/Xpand/EXL	1	*Positive	1	100.00%
E-Screen/Alere	49	*Positive	48	97.96%
Express Diag. Express Check	1	*Positive	1	100.00%
Alere iCup Drugs of Abuse	19	*Positive	19	100.00%
Instant Technologies i Screen	2	*Positive	1	50.00%
Integrated E-Z Split Key Cup	1	*Positive	1	100.00%
MedTox Profile II	3	*Positive	3	100.00%
Beckman AU	1	*Positive	1	100.00%
Innovacon Multi-Clin	5	*Positive	5	100.00%
CLIAwaived Inc Rapid Drug Cup	3	*Positive	3	100.00%
MedTox Profile V	2	*Positive	2	100.00%
Jant Accutest Multi Drug Screen	1	*Positive	1	100.00%
Amedica Biotech AmediCheck Instant Test Cup	4	*Positive	4	100.00%
Nova Century Scientific Nova Screen DOA Cassette	1	*Positive	1	100.00%
American Screening Corp Reveal Multi-Drug Test Cup	4	*Positive	4	100.00%
Mindray BS-200	2	*Positive	2	100.00%
Medica Easy RA	1	*Positive	1	100.00%
Alere iCassette	6	*Positive	6	100.00%
Roche Hitachi 717	1	*Positive	1	100.00%
Alere USD UScreen Cup	8	*Positive	8	100.00%
TFS Indiko/Indiko Plus	2	*Positive	2	100.00%
Alere iScreen Drugs of Abuse Card	1	*Positive	1	100.00%
Noble Medical Insta-Screen Test Cup	3	*Positive	3	100.00%
Hemosure First Sign UDS Dipstick	14	*Positive	14	100.00%
Siemens Viva/ProE System	3	*Positive	3	100.00%
TFS DRI Assay	1	*Positive	1	100.00%
Premier Biotech Bio-Cup	12	*Positive	12	100.00%
UCP Biosciences Home Drug Screening Test Cup	1	*Positive	1	100.00%
Abbott Architect c and ci Series	1	*Positive	1	100.00%
Consult Diag Drugs of Abuse Cups	2	*Positive	2	100.00%
Branan Med Corp ToxCup	5	*Positive	5	100.00%
American Screening Discover Plus Cup	1	*Positive	1	100.00%
McKesson Drugs of Abuse Test Cup	2	*Positive	2	100.00%
Henry Schein OneStep+ DOA	5	*Positive	5	100.00%
Siemens Atellica	1	*Positive	1	100.00%
Confirm Biosciences T Cup	1	*Positive	1	100.00%
Healgen Sci Drug Test Cup	1	*Positive	1	100.00%
E-Screen/Alere	49	Negative	1	2.04%
Instant Technologies i Screen	2	Negative	1	50.00%
Noble Med Split Specimen Cup	1	Negative	1	100.00%
ASC One Screen Drug Test Card	1	Negative	1	100.00%



		Cocaine Metabolites		
UDS-5		* Denotes Correct Response		
AG: All Methods	178			
		*Positive	173	97.19%
		Negative	5	2.81%
AlereTriage	4	*Positive	4	100.00%
Siemens Dim/AR/ES/RxL/Xpand/EXL	1	*Positive	1	100.00%
E-Screen/Alere	49	*Positive	47	95.92%
Express Diag. Express Check	1	*Positive	1	100.00%
Alere iCup Drugs of Abuse	19	*Positive	19	100.00%
Instant Technologies i Screen	2	*Positive	2	100.00%
Integrated E-Z Split Key Cup	1	*Positive	1	100.00%
MedTox Profile II	3	*Positive	3	100.00%
Beckman AU	1	*Positive	1	100.00%
Innovacon Multi-Clin	5	*Positive	5	100.00%
CLIAwaived Inc Rapid Drug Cup	3	*Positive	3	100.00%
MedTox Profile V	2	*Positive	2	100.00%
Jant Accutest Multi Drug Screen	1	*Positive	1	100.00%
Amedica Biotech AmediCheck Instant Test Cup	4	*Positive	4	100.00%
Nova Century Scientific Nova Screen DOA Cassette	1	*Positive	1	100.00%
American Screening Corp Reveal Multi-Drug Test Cup	4	*Positive	4	100.00%
Mindray BS-200	2	*Positive	2	100.00%
Medica Easy RA	1	*Positive	1	100.00%
Alere iCassette	6	*Positive	6	100.00%
Roche Hitachi 717	1	*Positive	1	100.00%
Alere USD UScreen Cup	8	*Positive	8	100.00%
TFS Indiko/Indiko Plus	2	*Positive	2	100.00%
Alere iScreen Drugs of Abuse Card	1	*Positive	1	100.00%
Noble Medical Insta-Screen Test Cup	3	*Positive	3	100.00%
Hemosure First Sign UDS Dipstick	14	*Positive	14	100.00%
Siemens Viva/ProE System	4	*Positive	4	100.00%
TFS DRI Assay	1	*Positive	1	100.00%
Premier Biotech Bio-Cup	12	*Positive	12	100.00%
UCP Biosciences Home Drug Screening Test Cup	1	*Positive	1	100.00%
Abbott Architect c and ci Series	1	*Positive	1	100.00%
Consult Diag Drugs of Abuse Cups	2	*Positive	2	100.00%
Branan Med Corp ToxCup	5	*Positive	5	100.00%
American Screening Discover Plus Cup	1	*Positive	1	100.00%
McKesson Drugs of Abuse Test Cup	2	*Positive	1	50.00%
Henry Schein OneStep+ DOA	5	*Positive	5	100.00%
Siemens Atellica	1	*Positive	1	100.00%
Confirm Biosciences T Cup	1	*Positive	1	100.00%
Healgen Sci Drug Test Cup	1	*Positive	1	100.00%
E-Screen/Alere	49	Negative	2	4.08%
Noble Med Split Specimen Cup	1	Negative	1	100.00%
ASC One Screen Drug Test Card	1	Negative	1	100.00%
McKesson Drugs of Abuse Test Cup	2	Negative	1	50.00%
UDS-6		* Denotes Correct Response		
AG: All Methods	178			
		*Negative	175	98.31%
		Positive	3	1.69%
AlereTriage	4	*Negative	4	100.00%
Siemens Dim/AR/ES/RxL/Xpand/EXL	1	*Negative	1	100.00%
E-Screen/Alere	49	*Negative	48	97.96%
Express Diag. Express Check	1	*Negative	1	100.00%
Alere iCup Drugs of Abuse	19	*Negative	19	100.00%
Instant Technologies i Screen	2	*Negative	2	100.00%
Integrated E-Z Split Key Cup	1	*Negative	1	100.00%
MedTox Profile II	3	*Negative	3	100.00%
Beckman AU	1	*Negative	1	100.00%
Innovacon Multi-Clin	5	*Negative	5	100.00%
CLIAwaived Inc Rapid Drug Cup	3	*Negative	3	100.00%
MedTox Profile V	2	*Negative	2	100.00%
Jant Accutest Multi Drug Screen	1	*Negative	1	100.00%
Amedica Biotech AmediCheck Instant Test Cup	4	*Negative	4	100.00%
Nova Century Scientific Nova Screen DOA Cassette	1	*Negative	1	100.00%
American Screening Corp Reveal Multi-Drug Test Cup	4	*Negative	4	100.00%
Mindray BS-200	2	*Negative	2	100.00%
Medica Easy RA	1	*Negative	1	100.00%



Alere iCassette	6	*Negative	6	100.00%
Roche Hitachi 717	1	*Negative	1	100.00%
Alere USD UScreen Cup	8	*Negative	8	100.00%
TFS Indiko/Indiko Plus	2	*Negative	2	100.00%
Alere iScreen Drugs of Abuse Card	1	*Negative	1	100.00%
Noble Medical Insta-Screen Test Cup	3	*Negative	3	100.00%
Hemosure First Sign UDS Dipstick	14	*Negative	14	100.00%
Siemens Viva/ProE System	4	*Negative	4	100.00%
TFS DRI Assay	1	*Negative	1	100.00%
Premier Biotech Bio-Cup	12	*Negative	12	100.00%
UCP Biosciences Home Drug Screening Test Cup	1	*Negative	1	100.00%
Abbott Architect c and ci Series	1	*Negative	1	100.00%
Consult Diag Drugs of Abuse Cups	2	*Negative	2	100.00%
Branan Med Corp ToxCup	5	*Negative	5	100.00%
American Screening Discover Plus Cup	1	*Negative	1	100.00%
McKesson Drugs of Abuse Test Cup	2	*Negative	2	100.00%
Henry Schein OneStep+ DOA	5	*Negative	5	100.00%
Siemens Atellica	1	*Negative	1	100.00%
Confirm Biosciences T Cup	1	*Negative	1	100.00%
Healgen Sci Drug Test Cup	1	*Negative	1	100.00%
E-Screen/Alere	49	Positive	1	2.04%
Noble Med Split Specimen Cup	1	Positive	1	100.00%
ASC One Screen Drug Test Card	1	Positive	1	100.00%
		Ethanol		
UDS-5		* Denotes Correct Response		
AG: All Methods	7			
		*Negative	6	85.71%
		Positive	1	14.29%
Mindray BS-200	2	*Negative	2	100.00%
Roche Hitachi 717	1	*Negative	1	100.00%
TFS Indiko/Indiko Plus	1	*Negative	1	100.00%
Siemens Viva/ProE System	1	*Negative	1	100.00%
TFS DRI Assay	1	*Negative	1	100.00%
Siemens Atellica	1	Positive	1	100.00%
UDS-6		* Denotes Correct Response		
AG: All Methods	7			
		*Negative	6	85.71%
		Positive	1	14.29%
Mindray BS-200	2	*Negative	2	100.00%
Roche Hitachi 717	1	*Negative	1	100.00%
TFS Indiko/Indiko Plus	1	*Negative	1	100.00%
Siemens Viva/ProE System	1	*Negative	1	100.00%
TFS DRI Assay	1	*Negative	1	100.00%
Siemens Atellica	1	Positive	1	100.00%
		Acetaminophen		
UDS-5		* Denotes Correct Response		
AG: All Methods	6			
		*Negative	6	100.00%
AlereTriage	3	*Negative	3	100.00%
E-Screen/Alere	1	*Negative	1	100.00%
Noble Med Split Specimen Cup	1	*Negative	1	100.00%
Jant Accutest Multi Drug Screen	1	*Negative	1	100.00%
UDS-6		* Denotes Correct Response		
AG: All Methods	6			
		*Negative	6	100.00%
AlereTriage	3	*Negative	3	100.00%
E-Screen/Alere	1	*Negative	1	100.00%
Noble Med Split Specimen Cup	1	*Negative	1	100.00%
Jant Accutest Multi Drug Screen	1	*Negative	1	100.00%
		MDMA (Ecstasy)		
UDS-5		* Denotes Correct Response		
AG: All Methods	88			
		*Negative	88	100.00%
E-Screen/Alere	4	*Negative	4	100.00%
Alere iCup Drugs of Abuse	11	*Negative	11	100.00%
Innovacon Multi-Clin	4	*Negative	4	100.00%
CLIAwaived Inc Rapid Drug Cup	3	*Negative	3	100.00%
Jant Accutest Multi Drug Screen	1	*Negative	1	100.00%



Amedica Biotech AmediCheck Instant Test Cup	4	*Negative	4	100.00%
Nova Century Scientific Nova Screen DOA Cassette	1	*Negative	1	100.00%
American Screening Corp Reveal Multi-Drug Test Cup	2	*Negative	2	100.00%
Alere iCassette	3	*Negative	3	100.00%
Alere USD UScreen Cup	8	*Negative	8	100.00%
Alere iScreen Drugs of Abuse Card	1	*Negative	1	100.00%
Noble Medical Insta-Screen Test Cup	3	*Negative	3	100.00%
Hemosure First Sign UDS Dipstick	14	*Negative	14	100.00%
Premier Biotech Bio-Cup	12	*Negative	12	100.00%
UCP Biosciences Home Drug Screening Test Cup	1	*Negative	1	100.00%
Consult Diag Drugs of Abuse Cups	2	*Negative	2	100.00%
Branan Med Corp ToxCup	5	*Negative	5	100.00%
American Screening Discover Plus Cup	1	*Negative	1	100.00%
McKesson Drugs of Abuse Test Cup	2	*Negative	2	100.00%
Henry Schein OneStep+ DOA	5	*Negative	5	100.00%
Healgen Sci Drug Test Cup	1	*Negative	1	100.00%
UDS-6		* Denotes Correct Response		
AG: All Methods	89			
		*Negative	89	100.00%
E-Screen/Alere	4	*Negative	4	100.00%
Alere iCup Drugs of Abuse	11	*Negative	11	100.00%
Innovacon Multi-Clin	5	*Negative	5	100.00%
CLIAwaived Inc Rapid Drug Cup	3	*Negative	3	100.00%
Jant Accutest Multi Drug Screen	1	*Negative	1	100.00%
Amedica Biotech AmediCheck Instant Test Cup	4	*Negative	4	100.00%
Nova Century Scientific Nova Screen DOA Cassette	1	*Negative	1	100.00%
American Screening Corp Reveal Multi-Drug Test Cup	2	*Negative	2	100.00%
Alere iCassette	3	*Negative	3	100.00%
Alere USD UScreen Cup	8	*Negative	8	100.00%
Alere iScreen Drugs of Abuse Card	1	*Negative	1	100.00%
Noble Medical Insta-Screen Test Cup	3	*Negative	3	100.00%
Hemosure First Sign UDS Dipstick	14	*Negative	14	100.00%
Premier Biotech Bio-Cup	12	*Negative	12	100.00%
UCP Biosciences Home Drug Screening Test Cup	1	*Negative	1	100.00%
Consult Diag Drugs of Abuse Cups	2	*Negative	2	100.00%
Branan Med Corp ToxCup	5	*Negative	5	100.00%
American Screening Discover Plus Cup	1	*Negative	1	100.00%
McKesson Drugs of Abuse Test Cup	2	*Negative	2	100.00%
Henry Schein OneStep+ DOA	5	*Negative	5	100.00%
Healgen Sci Drug Test Cup	1	*Negative	1	100.00%
		Methadone		
UDS-5		* Denotes Correct Response		
AG: All Methods	115			
		*Positive	99	86.09%
		Negative	16	13.91%
AlereTriage	3	*Positive	3	100.00%
Siemens Dim/AR/ES/RxL/Xpand/EXL	1	*Positive	1	100.00%
E-Screen/Alere	25	*Positive	25	100.00%
Express Diag. Express Check	1	*Positive	1	100.00%
Alere iCup Drugs of Abuse	14	*Positive	10	71.43%
MedTox Profile II	1	*Positive	1	100.00%
Beckman AU	1	*Positive	1	100.00%
CLIAwaived Inc Rapid Drug Cup	3	*Positive	3	100.00%
MedTox Profile V	2	*Positive	2	100.00%
Jant Accutest Multi Drug Screen	1	*Positive	1	100.00%
Amedica Biotech AmediCheck Instant Test Cup	4	*Positive	3	75.00%
Nova Century Scientific Nova Screen DOA Cassette	1	*Positive	1	100.00%
American Screening Corp Reveal Multi-Drug Test Cup	4	*Positive	4	100.00%
Mindray BS-200	2	*Positive	2	100.00%
Medica Easy RA	1	*Positive	1	100.00%
Roche Hitachi 717	1	*Positive	1	100.00%
Alere USD UScreen Cup	8	*Positive	6	75.00%
Hemosure First Sign UDS Dipstick	14	*Positive	14	100.00%
Siemens Viva/ProE System	2	*Positive	2	100.00%
TFS DRI Assay	1	*Positive	1	100.00%
UCP Biosciences Home Drug Screening Test Cup	1	*Positive	1	100.00%
Abbott Architect c and ci Series	1	*Positive	1	100.00%
Consult Diag Drugs of Abuse Cups	2	*Positive	2	100.00%
Branan Med Corp ToxCup	5	*Positive	4	80.00%



American Screening Discover Plus Cup	1	*Positive	1	100.00%
McKesson Drugs of Abuse Test Cup	2	*Positive	1	50.00%
Henry Schein OneStep+ DOA	5	*Positive	3	60.00%
Siemens Atellica	1	*Positive	1	100.00%
Confirm Biosciences T Cup	1	*Positive	1	100.00%
Healgen Sci Drug Test Cup	1	*Positive	1	100.00%
Alere iCup Drugs of Abuse	14	Negative	4	28.57%
ASC One Screen Drug Test Card	1	Negative	1	100.00%
Amedica Biotech AmediCheck Instant Test Cup	4	Negative	1	25.00%
Alere USD UScreen Cup	8	Negative	2	25.00%
Alere iScreen Drugs of Abuse Card	1	Negative	1	100.00%
Noble Medical Insta-Screen Test Cup	3	Negative	3	100.00%
Branan Med Corp ToxCup	5	Negative	1	20.00%
McKesson Drugs of Abuse Test Cup	2	Negative	1	50.00%
Henry Schein OneStep+ DOA	5	Negative	2	40.00%
UDS-6		* Denotes Correct Response		
AG: All Methods	115			
		*Negative	112	97.39%
		Positive	3	2.61%
AlereTriage	3	*Negative	3	100.00%
Siemens Dim/AR/ES/RxL/Xpand/EXL	1	*Negative	1	100.00%
E-Screen/Alere	25	*Negative	25	100.00%
Express Diag. Express Check	1	*Negative	1	100.00%
Alere iCup Drugs of Abuse	14	*Negative	12	85.71%
MedTox Profile II	1	*Negative	1	100.00%
Beckman AU	1	*Negative	1	100.00%
CLIAwaived Inc Rapid Drug Cup	3	*Negative	3	100.00%
MedTox Profile V	2	*Negative	2	100.00%
Jant Accutest Multi Drug Screen	1	*Negative	1	100.00%
Amedica Biotech AmediCheck Instant Test Cup	4	*Negative	4	100.00%
Nova Century Scientific Nova Screen DOA Cassette	1	*Negative	1	100.00%
American Screening Corp Reveal Multi-Drug Test Cup	4	*Negative	4	100.00%
Mindray BS-200	2	*Negative	2	100.00%
Medica Easy RA	1	*Negative	1	100.00%
Roche Hitachi 717	1	*Negative	1	100.00%
Alere USD UScreen Cup	8	*Negative	8	100.00%
Alere iScreen Drugs of Abuse Card	1	*Negative	1	100.00%
Noble Medical Insta-Screen Test Cup	3	*Negative	3	100.00%
Hemosure First Sign UDS Dipstick	14	*Negative	14	100.00%
Siemens Viva/ProE System	2	*Negative	2	100.00%
TFS DRI Assay	1	*Negative	1	100.00%
UCP Biosciences Home Drug Screening Test Cup	1	*Negative	1	100.00%
Abbott Architect c and ci Series	1	*Negative	1	100.00%
Consult Diag Drugs of Abuse Cups	2	*Negative	2	100.00%
Branan Med Corp ToxCup	5	*Negative	5	100.00%
American Screening Discover Plus Cup	1	*Negative	1	100.00%
McKesson Drugs of Abuse Test Cup	2	*Negative	2	100.00%
Henry Schein OneStep+ DOA	5	*Negative	5	100.00%
Siemens Atellica	1	*Negative	1	100.00%
Confirm Biosciences T Cup	1	*Negative	1	100.00%
Healgen Sci Drug Test Cup	1	*Negative	1	100.00%
Alere iCup Drugs of Abuse	14	Positive	2	14.29%
ASC One Screen Drug Test Card	1	Positive	1	100.00%
		Methamphetamine		
UDS-5		* Denotes Correct Response		
AG: All Methods	148			
		*Negative	143	96.62%
		Positive	5	3.38%
AlereTriage	4	*Negative	4	100.00%
E-Screen/Alere	48	*Negative	47	97.92%
Alere iCup Drugs of Abuse	18	*Negative	17	94.44%
Instant Technologies i Screen	1	*Negative	1	100.00%
Integrated E-Z Split Key Cup	1	*Negative	1	100.00%
MedTox Profile II	2	*Negative	2	100.00%
Noble Med Split Specimen Cup	1	*Negative	1	100.00%
CLIAwaived Inc Rapid Drug Cup	3	*Negative	3	100.00%
ASC One Screen Drug Test Card	1	*Negative	1	100.00%
MedTox Profile V	2	*Negative	2	100.00%
Jant Accutest Multi Drug Screen	1	*Negative	1	100.00%



Amedica Biotech AmediCheck Instant Test Cup	4	*Negative	4	100.00%
Nova Century Scientific Nova Screen DOA Cassette	1	*Negative	1	100.00%
American Screening Corp Reveal Multi-Drug Test Cup	4	*Negative	4	100.00%
Alere USD UScreen Cup	8	*Negative	7	87.50%
Alere iScreen Drugs of Abuse Card	1	*Negative	1	100.00%
Noble Medical Insta-Screen Test Cup	3	*Negative	2	66.67%
Hemosure First Sign UDS Dipstick	14	*Negative	14	100.00%
Premier Biotech Bio-Cup	12	*Negative	12	100.00%
UCP Biosciences Home Drug Screening Test Cup	1	*Negative	1	100.00%
Consult Diag Drugs of Abuse Cups	2	*Negative	2	100.00%
Branan Med Corp ToxCup	5	*Negative	5	100.00%
American Screening Discover Plus Cup	1	*Negative	1	100.00%
McKesson Drugs of Abuse Test Cup	2	*Negative	2	100.00%
Henry Schein OneStep+ DOA	5	*Negative	5	100.00%
Confirm Biosciences T Cup	1	*Negative	1	100.00%
Healgen Sci Drug Test Cup	1	*Negative	1	100.00%
E-Screen/Alere	48	Positive	1	2.08%
Express Diag. Express Check	1	Positive	1	100.00%
Alere iCup Drugs of Abuse	18	Positive	1	5.56%
Alere USD UScreen Cup	8	Positive	1	12.50%
Noble Medical Insta-Screen Test Cup	3	Positive	1	33.33%
UDS-6		* Denotes Correct Response		
AG: All Methods	148			
		*Negative	147	99.32%
		Positive	1	0.68%
AlereTriage	4	*Negative	3	75.00%
E-Screen/Alere	48	*Negative	48	100.00%
Express Diag. Express Check	1	*Negative	1	100.00%
Alere iCup Drugs of Abuse	18	*Negative	18	100.00%
Instant Technologies i Screen	1	*Negative	1	100.00%
Integrated E-Z Split Key Cup	1	*Negative	1	100.00%
MedTox Profile II	2	*Negative	2	100.00%
Noble Med Split Specimen Cup	1	*Negative	1	100.00%
CLIAwaived Inc Rapid Drug Cup	3	*Negative	3	100.00%
ASC One Screen Drug Test Card	1	*Negative	1	100.00%
MedTox Profile V	2	*Negative	2	100.00%
Jant Accutest Multi Drug Screen	1	*Negative	1	100.00%
Amedica Biotech AmediCheck Instant Test Cup	4	*Negative	4	100.00%
Nova Century Scientific Nova Screen DOA Cassette	1	*Negative	1	100.00%
American Screening Corp Reveal Multi-Drug Test Cup	4	*Negative	4	100.00%
Alere USD UScreen Cup	8	*Negative	8	100.00%
Alere iScreen Drugs of Abuse Card	1	*Negative	1	100.00%
Noble Medical Insta-Screen Test Cup	3	*Negative	3	100.00%
Hemosure First Sign UDS Dipstick	14	*Negative	14	100.00%
Premier Biotech Bio-Cup	12	*Negative	12	100.00%
UCP Biosciences Home Drug Screening Test Cup	1	*Negative	1	100.00%
Consult Diag Drugs of Abuse Cups	2	*Negative	2	100.00%
Branan Med Corp ToxCup	5	*Negative	5	100.00%
American Screening Discover Plus Cup	1	*Negative	1	100.00%
McKesson Drugs of Abuse Test Cup	2	*Negative	2	100.00%
Henry Schein OneStep+ DOA	5	*Negative	5	100.00%
Confirm Biosciences T Cup	1	*Negative	1	100.00%
Healgen Sci Drug Test Cup	1	*Negative	1	100.00%
AlereTriage	4	Positive	1	25.00%
		Methaqualone		
UDS-5		* Denotes Correct Response		
AG: All Methods	1			
		*Negative	1	100.00%
Beckman AU	1	*Negative	1	100.00%
UDS-6		* Denotes Correct Response		
AG: All Methods	1			
		*Negative	1	100.00%
Beckman AU	1	*Negative	1	100.00%
		Opiates		
UDS-5		* Denotes Correct Response		
AG: All Methods	178			
		*Positive	158	88.76%
		Negative	20	11.24%
AlereTriage	4	*Positive	4	100.00%



Siemens Dim/AR/ES/RxL/Xpand/EXL	1	*Positive	1	100.00%
E-Screen/Alere	49	*Positive	34	69.39%
Express Diag. Express Check	1	*Positive	1	100.00%
Alere iCup Drugs of Abuse	19	*Positive	19	100.00%
Instant Technologies i Screen	2	*Positive	1	50.00%
Integrated E-Z Split Key Cup	1	*Positive	1	100.00%
MedTox Profile II	2	*Positive	2	100.00%
Noble Med Split Specimen Cup	1	*Positive	1	100.00%
Beckman AU	1	*Positive	1	100.00%
Innovacon Multi-Clin	5	*Positive	5	100.00%
CLIAwaived Inc Rapid Drug Cup	3	*Positive	3	100.00%
ASC One Screen Drug Test Card	1	*Positive	1	100.00%
MedTox Profile V	2	*Positive	2	100.00%
Jant Accutest Multi Drug Screen	1	*Positive	1	100.00%
Amedica Biotech AmediCheck Instant Test Cup	4	*Positive	4	100.00%
Nova Century Scientific Nova Screen DOA Cassette	1	*Positive	1	100.00%
American Screening Corp Reveal Multi-Drug Test Cup	4	*Positive	2	50.00%
Mindray BS-200	2	*Positive	2	100.00%
Medica Easy RA	1	*Positive	1	100.00%
Alere iCassette	6	*Positive	6	100.00%
Roche Hitachi 717	1	*Positive	1	100.00%
Alere USD UScreen Cup	8	*Positive	8	100.00%
TFS Indiko/Indiko Plus	3	*Positive	3	100.00%
Alere iScreen Drugs of Abuse Card	1	*Positive	1	100.00%
Noble Medical Insta-Screen Test Cup	3	*Positive	3	100.00%
Hemosure First Sign UDS Dipstick	14	*Positive	14	100.00%
Siemens Viva/ProE System	4	*Positive	4	100.00%
TFS DRI Assay	1	*Positive	1	100.00%
Premier Biotech Bio-Cup	12	*Positive	12	100.00%
UCP Biosciences Home Drug Screening Test Cup	1	*Positive	1	100.00%
Abbott Architect c and ci Series	1	*Positive	1	100.00%
Consult Diag Drugs of Abuse Cups	2	*Positive	2	100.00%
Branan Med Corp ToxCup	5	*Positive	4	80.00%
American Screening Discover Plus Cup	1	*Positive	1	100.00%
McKesson Drugs of Abuse Test Cup	2	*Positive	2	100.00%
Henry Schein OneStep+ DOA	5	*Positive	4	80.00%
Siemens Atellica	1	*Positive	1	100.00%
Confirm Biosciences T Cup	1	*Positive	1	100.00%
Healgen Sci Drug Test Cup	1	*Positive	1	100.00%
E-Screen/Alere	49	Negative	15	30.61%
Instant Technologies i Screen	2	Negative	1	50.00%
American Screening Corp Reveal Multi-Drug Test Cup	4	Negative	2	50.00%
Branan Med Corp ToxCup	5	Negative	1	20.00%
Henry Schein OneStep+ DOA	5	Negative	1	20.00%
UDS-6		* Denotes Correct Response		
AG: All Methods	178			
		*Positive	154	86.52%
		Negative	24	13.48%
AlereTriage	4	*Positive	4	100.00%
Siemens Dim/AR/ES/RxL/Xpand/EXL	1	*Positive	1	100.00%
E-Screen/Alere	49	*Positive	30	61.22%
Express Diag. Express Check	1	*Positive	1	100.00%
Alere iCup Drugs of Abuse	19	*Positive	18	94.74%
Instant Technologies i Screen	2	*Positive	2	100.00%
Integrated E-Z Split Key Cup	1	*Positive	1	100.00%
MedTox Profile II	2	*Positive	2	100.00%
Noble Med Split Specimen Cup	1	*Positive	1	100.00%
Beckman AU	1	*Positive	1	100.00%
Innovacon Multi-Clin	5	*Positive	5	100.00%
CLIAwaived Inc Rapid Drug Cup	3	*Positive	3	100.00%
ASC One Screen Drug Test Card	1	*Positive	1	100.00%
MedTox Profile V	2	*Positive	2	100.00%
Jant Accutest Multi Drug Screen	1	*Positive	1	100.00%
Amedica Biotech AmediCheck Instant Test Cup	4	*Positive	3	75.00%
Nova Century Scientific Nova Screen DOA Cassette	1	*Positive	1	100.00%
American Screening Corp Reveal Multi-Drug Test Cup	4	*Positive	2	50.00%
Mindray BS-200	2	*Positive	2	100.00%
Medica Easy RA	1	*Positive	1	100.00%
Alere iCassette	6	*Positive	6	100.00%



Roche Hitachi 717	1	*Positive	1	100.00%
Alere USD UScreen Cup	8	*Positive	8	100.00%
TFS Indiko/Indiko Plus	3	*Positive	3	100.00%
Alere iScreen Drugs of Abuse Card	1	*Positive	1	100.00%
Noble Medical Insta-Screen Test Cup	3	*Positive	3	100.00%
Hemosure First Sign UDS Dipstick	14	*Positive	14	100.00%
Siemens Viva/ProE System	4	*Positive	4	100.00%
TFS DRI Assay	1	*Positive	1	100.00%
Premier Biotech Bio-Cup	12	*Positive	12	100.00%
UCP Biosciences Home Drug Screening Test Cup	1	*Positive	1	100.00%
Abbott Architect c and ci Series	1	*Positive	1	100.00%
Consult Diag Drugs of Abuse Cups	2	*Positive	2	100.00%
Branan Med Corp ToxCup	5	*Positive	5	100.00%
American Screening Discover Plus Cup	1	*Positive	1	100.00%
McKesson Drugs of Abuse Test Cup	2	*Positive	2	100.00%
Henry Schein OneStep+ DOA	5	*Positive	4	80.00%
Siemens Atellica	1	*Positive	1	100.00%
Confirm Biosciences T Cup	1	*Positive	1	100.00%
Healgen Sci Drug Test Cup	1	*Positive	1	100.00%
E-Screen/Alere	49	Negative	19	38.78%
Alere iCup Drugs of Abuse	19	Negative	1	5.26%
Amedica Biotech AmediCheck Instant Test Cup	4	Negative	1	25.00%
American Screening Corp Reveal Multi-Drug Test Cup	4	Negative	2	50.00%
Henry Schein OneStep+ DOA	5	Negative	1	20.00%
		Oxycodone		
UDS-5		* Denotes Correct Response		
AG: All Methods	111			
		*Negative	107	96.40%
		Positive	4	3.60%
E-Screen/Alere	17	*Negative	17	100.00%
Alere iCup Drugs of Abuse	14	*Negative	14	100.00%
MedTox Profile II	1	*Negative	1	100.00%
Beckman AU	1	*Negative	1	100.00%
Innovacon Multi-Clin	5	*Negative	5	100.00%
CLIAwaived Inc Rapid Drug Cup	3	*Negative	3	100.00%
ASC One Screen Drug Test Card	1	*Negative	1	100.00%
MedTox Profile V	2	*Negative	2	100.00%
Jant Accutest Multi Drug Screen	1	*Negative	1	100.00%
Amedica Biotech AmediCheck Instant Test Cup	4	*Negative	1	25.00%
Nova Century Scientific Nova Screen DOA Cassette	1	*Negative	1	100.00%
American Screening Corp Reveal Multi-Drug Test Cup	4	*Negative	4	100.00%
Mindray BS-200	2	*Negative	2	100.00%
Medica Easy RA	1	*Negative	1	100.00%
Alere iCassette	3	*Negative	3	100.00%
Roche Hitachi 717	1	*Negative	1	100.00%
Alere USD UScreen Cup	8	*Negative	8	100.00%
TFS Indiko/Indiko Plus	2	*Negative	2	100.00%
Alere iScreen Drugs of Abuse Card	2	*Negative	2	100.00%
Noble Medical Insta-Screen Test Cup	3	*Negative	3	100.00%
Siemens Viva/ProE System	2	*Negative	2	100.00%
TFS DRI Assay	1	*Negative	1	100.00%
Premier Biotech Bio-Cup	12	*Negative	12	100.00%
UCP Biosciences Home Drug Screening Test Cup	1	*Negative	1	100.00%
Abbott Architect c and ci Series	1	*Negative	1	100.00%
Consult Diag Drugs of Abuse Cups	2	*Negative	2	100.00%
Branan Med Corp ToxCup	5	*Negative	4	80.00%
American Screening Discover Plus Cup	1	*Negative	1	100.00%
McKesson Drugs of Abuse Test Cup	2	*Negative	2	100.00%
Henry Schein OneStep+ DOA	5	*Negative	5	100.00%
Siemens Atellica	1	*Negative	1	100.00%
Confirm Biosciences T Cup	1	*Negative	1	100.00%
Healgen Sci Drug Test Cup	1	*Negative	1	100.00%
Amedica Biotech AmediCheck Instant Test Cup	4	Positive	3	75.00%
Branan Med Corp ToxCup	5	Positive	1	20.00%
UDS-6		* Denotes Correct Response		
AG: All Methods	111			
		*Negative	109	98.20%
		Positive	2	1.80%
E-Screen/Alere	17	*Negative	17	100.00%



Alere iCup Drugs of Abuse	14	*Negative	14	100.00%
MedTox Profile II	1	*Negative	1	100.00%
Beckman AU	1	*Negative	1	100.00%
Innovacon Multi-Clin	5	*Negative	5	100.00%
CLIAwaived Inc Rapid Drug Cup	3	*Negative	3	100.00%
ASC One Screen Drug Test Card	1	*Negative	1	100.00%
MedTox Profile V	2	*Negative	2	100.00%
Jant Accutest Multi Drug Screen	1	*Negative	1	100.00%
Amedica Biotech AmediCheck Instant Test Cup	4	*Negative	2	50.00%
Nova Century Scientific Nova Screen DOA Cassette	1	*Negative	1	100.00%
American Screening Corp Reveal Multi-Drug Test Cup	4	*Negative	4	100.00%
Mindray BS-200	2	*Negative	2	100.00%
Medica Easy RA	1	*Negative	1	100.00%
Alere iCassette	3	*Negative	3	100.00%
Roche Hitachi 717	1	*Negative	1	100.00%
Alere USD UScreen Cup	8	*Negative	8	100.00%
TFS Indiko/Indiko Plus	2	*Negative	2	100.00%
Alere iScreen Drugs of Abuse Card	2	*Negative	2	100.00%
Noble Medical Insta-Screen Test Cup	3	*Negative	3	100.00%
Siemens Viva/ProE System	2	*Negative	2	100.00%
TFS DRI Assay	1	*Negative	1	100.00%
Premier Biotech Bio-Cup	12	*Negative	12	100.00%
UCP Biosciences Home Drug Screening Test Cup	1	*Negative	1	100.00%
Abbott Architect c and ci Series	1	*Negative	1	100.00%
Consult Diag Drugs of Abuse Cups	2	*Negative	2	100.00%
Branan Med Corp ToxCup	5	*Negative	5	100.00%
American Screening Discover Plus Cup	1	*Negative	1	100.00%
McKesson Drugs of Abuse Test Cup	2	*Negative	2	100.00%
Henry Schein OneStep+ DOA	5	*Negative	5	100.00%
Siemens Atellica	1	*Negative	1	100.00%
Confirm Biosciences T Cup	1	*Negative	1	100.00%
Healgen Sci Drug Test Cup	1	*Negative	1	100.00%
Amedica Biotech AmediCheck Instant Test Cup	4	Positive	2	50.00%
		PCP (Phencyclidine)		
UDS-5		* Denotes Correct Response		
AG: All Methods	152			
		*Negative	149	98.03%
		Positive	3	1.97%
AlereTriage	3	*Negative	3	100.00%
E-Screen/Alere	49	*Negative	47	95.92%
Express Diag. Express Check	1	*Negative	1	100.00%
Alere iCup Drugs of Abuse	12	*Negative	12	100.00%
Instant Technologies i Screen	2	*Negative	2	100.00%
Integrated E-Z Split Key Cup	1	*Negative	1	100.00%
MedTox Profile II	3	*Negative	3	100.00%
Beckman AU	1	*Negative	1	100.00%
Innovacon Multi-Clin	5	*Negative	5	100.00%
CLIAwaived Inc Rapid Drug Cup	3	*Negative	3	100.00%
ASC One Screen Drug Test Card	1	*Negative	1	100.00%
MedTox Profile V	2	*Negative	2	100.00%
Jant Accutest Multi Drug Screen	1	*Negative	1	100.00%
Amedica Biotech AmediCheck Instant Test Cup	1	*Negative	1	100.00%
Nova Century Scientific Nova Screen DOA Cassette	1	*Negative	1	100.00%
American Screening Corp Reveal Multi-Drug Test Cup	4	*Negative	4	100.00%
Mindray BS-200	1	*Negative	1	100.00%
Alere iCassette	6	*Negative	6	100.00%
Roche Hitachi 717	1	*Negative	1	100.00%
Alere USD UScreen Cup	6	*Negative	6	100.00%
Alere iScreen Drugs of Abuse Card	1	*Negative	1	100.00%
Noble Medical Insta-Screen Test Cup	3	*Negative	3	100.00%
Hemosure First Sign UDS Dipstick	14	*Negative	14	100.00%
TFS DRI Assay	1	*Negative	1	100.00%
Premier Biotech Bio-Cup	12	*Negative	12	100.00%
UCP Biosciences Home Drug Screening Test Cup	1	*Negative	1	100.00%
Consult Diag Drugs of Abuse Cups	2	*Negative	2	100.00%
Branan Med Corp ToxCup	5	*Negative	4	80.00%
McKesson Drugs of Abuse Test Cup	2	*Negative	2	100.00%
Henry Schein OneStep+ DOA	5	*Negative	5	100.00%
Siemens Atellica	1	*Negative	1	100.00%



Healgen Sci Drug Test Cup	1	*Negative	1	100.00%
E-Screen/Alere	49	Positive	2	4.08%
Branan Med Corp ToxCup	5	Positive	1	20.00%
UDS-6		* Denotes Correct Response		
AG: All Methods	152			
		*Negative	152	100.00%
AlereTriage	3	*Negative	3	100.00%
E-Screen/Alere	49	*Negative	49	100.00%
Express Diag. Express Check	1	*Negative	1	100.00%
Alere iCup Drugs of Abuse	12	*Negative	12	100.00%
Instant Technologies i Screen	2	*Negative	2	100.00%
Integrated E-Z Split Key Cup	1	*Negative	1	100.00%
MedTox Profile II	3	*Negative	3	100.00%
Beckman AU	1	*Negative	1	100.00%
Innovacon Multi-Clin	5	*Negative	5	100.00%
CLIAwaived Inc Rapid Drug Cup	3	*Negative	3	100.00%
ASC One Screen Drug Test Card	1	*Negative	1	100.00%
MedTox Profile V	2	*Negative	2	100.00%
Jant Accutest Multi Drug Screen	1	*Negative	1	100.00%
Amedica Biotech AmediCheck Instant Test Cup	1	*Negative	1	100.00%
Nova Century Scientific Nova Screen DOA Cassette	1	*Negative	1	100.00%
American Screening Corp Reveal Multi-Drug Test Cup	4	*Negative	4	100.00%
Mindray BS-200	1	*Negative	1	100.00%
Alere iCassette	6	*Negative	6	100.00%
Roche Hitachi 717	1	*Negative	1	100.00%
Alere USD UScreen Cup	6	*Negative	6	100.00%
Alere iScreen Drugs of Abuse Card	1	*Negative	1	100.00%
Noble Medical Insta-Screen Test Cup	3	*Negative	3	100.00%
Hemosure First Sign UDS Dipstick	14	*Negative	14	100.00%
TFS DRI Assay	1	*Negative	1	100.00%
Premier Biotech Bio-Cup	12	*Negative	12	100.00%
UCP Biosciences Home Drug Screening Test Cup	1	*Negative	1	100.00%
Consult Diag Drugs of Abuse Cups	2	*Negative	2	100.00%
Branan Med Corp ToxCup	5	*Negative	5	100.00%
McKesson Drugs of Abuse Test Cup	2	*Negative	2	100.00%
Henry Schein OneStep+ DOA	5	*Negative	5	100.00%
Siemens Atellica	1	*Negative	1	100.00%
Healgen Sci Drug Test Cup	1	*Negative	1	100.00%
		Propoxyphene		
UDS-6		* Denotes Correct Response		
AG: All Methods	5			
		*Positive	4	80.00%
		Negative	1	20.00%
UDS-6		* Denotes Correct Response		
AG: All Methods	10			
		*Positive	4	40.00%
		Negative	6	60.00%
UDS-6		* Denotes Correct Response		
AG: All Methods	1			
		*Positive	1	100.00%
UDS-6		* Denotes Correct Response		
AG: All Methods	1			
		*Positive	1	100.00%
UDS-6		* Denotes Correct Response		
AG: All Methods	5			
		*Positive	4	80.00%
		Negative	1	20.00%
UDS-6		* Denotes Correct Response		
AG: All Methods	2			
		*Positive	2	100.00%
UDS-6		* Denotes Correct Response		
AG: All Methods	3			
		*Positive	3	100.00%
UDS-6		* Denotes Correct Response		
AG: All Methods	2			
		*Positive	2	100.00%
UDS-6		* Denotes Correct Response		
AG: All Methods	5			



		*Positive	5	100.00%
UDS-6		* Denotes Correct Response		
AG: All Methods	1			
		*Positive	1	100.00%
UDS-6		* Denotes Correct Response		
AG: All Methods	5			
		*Positive	2	40.00%
		Negative	3	60.00%
UDS-6		* Denotes Correct Response		
AG: All Methods	1			
		*Positive	1	100.00%
UDS-5		* Denotes Correct Response		
AG: All Methods	40			
		*Negative	39	97.50%
		Positive	1	2.50%
E-Screen/Alere	4	*Negative	4	100.00%
Alere iCup Drugs of Abuse	10	*Negative	10	100.00%
MedTox Profile II	1	*Negative	1	100.00%
Beckman AU	1	*Negative	1	100.00%
Innovacon Multi-Clin	5	*Negative	5	100.00%
MedTox Profile V	2	*Negative	2	100.00%
Alere iCassette	3	*Negative	3	100.00%
Consult Diag Drugs of Abuse Cups	2	*Negative	2	100.00%
Branan Med Corp ToxCup	5	*Negative	5	100.00%
Henry Schein OneStep+ DOA	5	*Negative	5	100.00%
Siemens Atellica	1	*Negative	1	100.00%
McKesson Drugs of Abuse Test Cup	1	Positive	1	100.00%
		Tricyclic Antidepressants		
UDS-5		* Denotes Correct Response		
AG: All Methods	53			
		*Negative	52	98.11%
		Positive	1	1.89%
AlereTriage	4	*Negative	4	100.00%
Alere iCup Drugs of Abuse	11	*Negative	11	100.00%
Innovacon Multi-Clin	4	*Negative	4	100.00%
MedTox Profile V	2	*Negative	2	100.00%
Jant Accutest Multi Drug Screen	1	*Negative	1	100.00%
Amedica Biotech AmediCheck Instant Test Cup	1	*Negative	1	100.00%
Nova Century Scientific Nova Screen DOA Cassette	1	*Negative	1	100.00%
Alere iCassette	3	*Negative	3	100.00%
Alere USD UScreen Cup	5	*Negative	5	100.00%
Alere iScreen Drugs of Abuse Card	1	*Negative	1	100.00%
Noble Medical Insta-Screen Test Cup	3	*Negative	3	100.00%
UCP Biosciences Home Drug Screening Test Cup	1	*Negative	1	100.00%
Consult Diag Drugs of Abuse Cups	2	*Negative	2	100.00%
Branan Med Corp ToxCup	5	*Negative	4	80.00%
McKesson Drugs of Abuse Test Cup	2	*Negative	2	100.00%
Henry Schein OneStep+ DOA	5	*Negative	5	100.00%
Confirm Biosciences T Cup	1	*Negative	1	100.00%
Healgen Sci Drug Test Cup	1	*Negative	1	100.00%
Branan Med Corp ToxCup	5	Positive	1	20.00%
UDS-6		* Denotes Correct Response		
AG: All Methods	53			
		*Negative	53	100.00%
AlereTriage	4	*Negative	4	100.00%
Alere iCup Drugs of Abuse	11	*Negative	11	100.00%
Innovacon Multi-Clin	4	*Negative	4	100.00%
MedTox Profile V	2	*Negative	2	100.00%
Jant Accutest Multi Drug Screen	1	*Negative	1	100.00%
Amedica Biotech AmediCheck Instant Test Cup	1	*Negative	1	100.00%
Nova Century Scientific Nova Screen DOA Cassette	1	*Negative	1	100.00%
Alere iCassette	3	*Negative	3	100.00%
Alere USD UScreen Cup	5	*Negative	5	100.00%
Alere iScreen Drugs of Abuse Card	1	*Negative	1	100.00%
Noble Medical Insta-Screen Test Cup	3	*Negative	3	100.00%
UCP Biosciences Home Drug Screening Test Cup	1	*Negative	1	100.00%
Consult Diag Drugs of Abuse Cups	2	*Negative	2	100.00%
Branan Med Corp ToxCup	5	*Negative	5	100.00%
McKesson Drugs of Abuse Test Cup	2	*Negative	2	100.00%



Henry Schein OneStep+ DOA	5	*Negative	5	100.00%
Confirm Biosciences T Cup	1	*Negative	1	100.00%
Healgen Sci Drug Test Cup	1	*Negative	1	100.00%
		Buprenorphine		
UDS-5		* Denotes Correct Response		
AG: All Methods	58			
		*Negative	58	100.00%
Alere iCup Drugs of Abuse	9	*Negative	9	100.00%
CLIAwaived Inc Rapid Drug Cup	3	*Negative	3	100.00%
ASC One Screen Drug Test Card	1	*Negative	1	100.00%
MedTox Profile V	1	*Negative	1	100.00%
Mindray BS-200	2	*Negative	2	100.00%
Medica Easy RA	1	*Negative	1	100.00%
Alere USD UScreen Cup	2	*Negative	2	100.00%
TFS Indiko/Indiko Plus	1	*Negative	1	100.00%
Alere iScreen Drugs of Abuse Card	2	*Negative	2	100.00%
Noble Medical Insta-Screen Test Cup	3	*Negative	3	100.00%
Siemens Viva/ProE System	4	*Negative	4	100.00%
TFS DRI Assay	1	*Negative	1	100.00%
Premier Biotech Bio-Cup	12	*Negative	12	100.00%
Consult Diag Drugs of Abuse Cups	2	*Negative	2	100.00%
Branan Med Corp ToxCup	5	*Negative	5	100.00%
American Screening Discover Plus Cup	1	*Negative	1	100.00%
McKesson Drugs of Abuse Test Cup	1	*Negative	1	100.00%
Henry Schein OneStep+ DOA	5	*Negative	5	100.00%
Siemens Atellica	1	*Negative	1	100.00%
Confirm Biosciences T Cup	1	*Negative	1	100.00%
UDS-6		* Denotes Correct Response		
AG: All Methods	58			
		*Negative	57	98.28%
		Positive	1	1.72%
Alere iCup Drugs of Abuse	9	*Negative	8	88.89%
CLIAwaived Inc Rapid Drug Cup	3	*Negative	3	100.00%
ASC One Screen Drug Test Card	1	*Negative	1	100.00%
MedTox Profile V	1	*Negative	1	100.00%
Mindray BS-200	2	*Negative	2	100.00%
Medica Easy RA	1	*Negative	1	100.00%
Alere USD UScreen Cup	2	*Negative	2	100.00%
TFS Indiko/Indiko Plus	1	*Negative	1	100.00%
Alere iScreen Drugs of Abuse Card	2	*Negative	2	100.00%
Noble Medical Insta-Screen Test Cup	3	*Negative	3	100.00%
Siemens Viva/ProE System	4	*Negative	4	100.00%
TFS DRI Assay	1	*Negative	1	100.00%
Premier Biotech Bio-Cup	12	*Negative	12	100.00%
Consult Diag Drugs of Abuse Cups	2	*Negative	2	100.00%
Branan Med Corp ToxCup	5	*Negative	5	100.00%
American Screening Discover Plus Cup	1	*Negative	1	100.00%
McKesson Drugs of Abuse Test Cup	1	*Negative	1	100.00%
Henry Schein OneStep+ DOA	5	*Negative	5	100.00%
Siemens Atellica	1	*Negative	1	100.00%
Confirm Biosciences T Cup	1	*Negative	1	100.00%
Alere iCup Drugs of Abuse	9	Positive	1	11.11%
		6-Monoacetylmorphine (6 AM)		
UDS-5		* Denotes Correct Response		
AG: All Methods	5			
		*Negative	5	100.00%
TFS Indiko/Indiko Plus	2	*Negative	2	100.00%
Siemens Viva/ProE System	3	*Negative	3	100.00%
UDS-6		* Denotes Correct Response		
AG: All Methods	5			
		*Negative	5	100.00%
TFS Indiko/Indiko Plus	2	*Negative	2	100.00%
Siemens Viva/ProE System	3	*Negative	3	100.00%
		EDDP		
UDS-5		* Denotes Correct Response		
AG: All Methods	16			
		*Negative	16	100.00%
AlereTriage	1	*Negative	1	100.00%



Siemens Viva/ProE System	3	*Negative	3	100.00%
Premier Biotech Bio-Cup	12	*Negative	12	100.00%
UDS-6		* Denotes Correct Response		
AG: All Methods	16			
		*Positive	16	100.00%
AlereTriage	1	*Positive	1	100.00%
Siemens Viva/ProE System	3	*Positive	3	100.00%
Premier Biotech Bio-Cup	12	*Positive	12	100.00%
		Hydrocodone		
UDS-5		* Denotes Correct Response		
AG: All Methods	1			
		*Negative	1	100.00%
Abbott Architect c and ci Series	1	*Negative	1	100.00%
UDS-6		* Denotes Correct Response		
AG: All Methods	1			
		*Negative	1	100.00%
Abbott Architect c and ci Series	1	*Negative	1	100.00%
		Tramadol		
UDS-5		* Denotes Correct Response		
AG: All Methods	1			
		*Positive	1	100.00%
Abbott Architect c and ci Series	1	*Positive	1	100.00%
UDS-6		* Denotes Correct Response		
AG: All Methods	1			
		*Negative	1	100.00%
Abbott Architect c and ci Series	1	*Negative	1	100.00%
705 - Urine Microalbumin & Creatinine	Labs	Response	Number	Percent
		Urine Microalbumin (Semi-Quant)		
UM-5		* Denotes Correct Response		
PG: Siemens Clinitek 50	6			
		*150 mg/L	1	16.67%
		*30 mg/L	1	16.67%
		*80 mg/L	4	66.67%
UM-6		* Denotes Correct Response		
PG: Siemens Clinitek 50	6			
		*150 mg/L	1	16.67%
		*80 mg/L	5	83.33%
UM-5		* Denotes Correct Response		
PG: Siemens Clinitek Status	29			
		*150 mg/L	2	6.90%
		*30 mg/L	2	6.90%
		*80 mg/L	25	86.21%
UM-6		* Denotes Correct Response		
PG: Siemens Clinitek Status	29			
		*150 mg/L	12	41.38%
		*80 mg/L	17	58.62%
UM-5		* Denotes Correct Response		
SG: Select Siemens Clinitek/Visual	37			
		*150 mg/L	3	8.11%
		*30 mg/L	3	8.11%
		*80 mg/L	31	83.78%
UM-6		* Denotes Correct Response		
SG: Select Siemens Clinitek/Visual	37			
		*150 mg/L	14	37.84%
		*80 mg/L	23	62.16%
		Urine Creatinine (Semi-Quant)		
UM-5		* Denotes Correct Response		
PG: Siemens Clinitek 50	6			
		*100 mg/dL	6	100.00%
UM-6		* Denotes Correct Response		
PG: Siemens Clinitek 50	6			
		*100 mg/dL	5	83.33%
		*50 mg/dL	1	16.67%
UM-5		* Denotes Correct Response		
PG: Siemens Clinitek Status	26			
		*100 mg/dL	12	46.15%
		*200 mg/dL	3	11.54%



		*50 mg/dL	9	34.62%
		300 mg/dL	2	7.69%
UM-6		* Denotes Correct Response		
PG: Siemens Clinitek Status	26			
		*100 mg/dL	11	42.31%
		*200 mg/dL	4	15.38%
		*50 mg/dL	10	38.46%
		300 mg/dL	1	3.85%
UM-5		* Denotes Correct Response		
SG: Select Siemens Clinitek/Visual	34			
		*100 mg/dL	19	55.88%
		*200 mg/dL	3	8.82%
		*50 mg/dL	10	29.41%
		300 mg/dL	2	5.88%
UM-6		* Denotes Correct Response		
SG: Select Siemens Clinitek/Visual	34			
		*100 mg/dL	17	50.00%
		*200 mg/dL	4	11.76%
		*50 mg/dL	12	35.29%
		300 mg/dL	1	2.94%
756 - HSV Antibodies	Labs	Response	Number	Percent
		HSV-2 Antibody		
VA-5		* Denotes Correct Response		
AG: All Methods	1			
		*Negative	1	100.00%
Biokit HSV-II	1	*Negative	1	100.00%
VA-6		* Denotes Correct Response		
AG: All Methods	1			
		*Positive	1	100.00%
Biokit HSV-II	1	*Positive	1	100.00%
770 - Viral Antigen Detection;792 - Respiratory Viral Antigen Detection (Waived)	Labs	Response	Number	Percent
		Influenza A - Aq Detect		
VAD-11		* Denotes Correct Response		
AG: All Methods	589			
		*Negative	575	97.62%
		Positive	14	2.38%
Binax Now Flu A	2	*Negative	2	100.00%
Alere BinaxNOW Influenza A+B Card	38	*Negative	36	94.74%
Sekisui Diagnostics OSOM Influenza A & B	23	*Negative	23	100.00%
Quidel QuickVue A+B	42	*Negative	42	100.00%
Quidel Sofia/Sofia 2 Influenza A+B	201	*Negative	198	98.51%
Alere Influenza A & B	5	*Negative	4	80.00%
BD Veritor Flu A+B	65	*Negative	64	98.46%
Henry Schein OneStep + Influenza A&B (distinguish)	5	*Negative	5	100.00%
Medline Influenza A and B	3	*Negative	3	100.00%
Cepheid GeneXpert - Xpert Flu	29	*Negative	28	96.55%
Alere ID NOW Influenza A & B	26	*Negative	25	96.15%
Roche Cobas Liat Influenza A/B	6	*Negative	6	100.00%
Sekisui Diagnostics OSOM Ultra Influenza A & B	85	*Negative	81	95.29%
Consult Diag Influenza A & B	33	*Negative	32	96.97%
Lifesign Status Influenza A & B	26	*Negative	26	100.00%
Alere BinaxNOW Influenza A+B Card	38	Positive	2	5.26%
Quidel Sofia/Sofia 2 Influenza A+B	201	Positive	3	1.49%
Alere Influenza A & B	5	Positive	1	20.00%
BD Veritor Flu A+B	65	Positive	1	1.54%
Cepheid GeneXpert - Xpert Flu	29	Positive	1	3.45%
Alere ID NOW Influenza A & B	26	Positive	1	3.85%
Sekisui Diagnostics OSOM Ultra Influenza A & B	85	Positive	4	4.71%
Consult Diag Influenza A & B	33	Positive	1	3.03%
VAD-12		* Denotes Correct Response		
AG: All Methods	584			
		*Negative	572	97.95%
		Positive	12	2.05%
Binax Now Flu A	2	*Negative	2	100.00%
Alere BinaxNOW Influenza A+B Card	34	*Negative	33	97.06%
Sekisui Diagnostics OSOM Influenza A & B	23	*Negative	22	95.65%
Quidel QuickVue A+B	42	*Negative	42	100.00%



Quidel Sofia/Sofia 2 Influenza A+B	201	*Negative	198	98.51%
Alere Influenza A & B	5	*Negative	4	80.00%
BD Veritor Flu A+B	65	*Negative	64	98.46%
Henry Schein OneStep + Influenza A&B (distinguish)	5	*Negative	5	100.00%
Medline Influenza A and B	3	*Negative	3	100.00%
Cepheid GeneXpert - Xpert Flu	29	*Negative	29	100.00%
Alere ID NOW Influenza A & B	25	*Negative	25	100.00%
Roche Cobas Liat Influenza A/B	6	*Negative	6	100.00%
Sekisui Diagnostics OSOM Ultra Influenza A & B	85	*Negative	82	96.47%
Consult Diag Influenza A & B	33	*Negative	32	96.97%
Lifesign Status Influenza A & B	26	*Negative	25	96.15%
Alere BinaxNOW Influenza A+B Card	34	Positive	1	2.94%
Sekisui Diagnostics OSOM Influenza A & B	23	Positive	1	4.35%
Quidel Sofia/Sofia 2 Influenza A+B	201	Positive	3	1.49%
Alere Influenza A & B	5	Positive	1	20.00%
BD Veritor Flu A+B	65	Positive	1	1.54%
Sekisui Diagnostics OSOM Ultra Influenza A & B	85	Positive	3	3.53%
Consult Diag Influenza A & B	33	Positive	1	3.03%
Lifesign Status Influenza A & B	26	Positive	1	3.85%
VAD-13		* Denotes Correct Response		
AG: All Methods	46			
		*Negative	45	97.83%
		Positive	1	2.17%
Alere BinaxNOW Influenza A+B Card	1	*Negative	1	100.00%
Sekisui Diagnostics OSOM Influenza A & B	8	*Negative	8	100.00%
Quidel QuickVue A+B	12	*Negative	11	91.67%
Quidel Sofia/Sofia 2 Influenza A+B	6	*Negative	6	100.00%
BD Veritor Flu A+B	5	*Negative	5	100.00%
Henry Schein OneStep + Influenza A&B (distinguish)	2	*Negative	2	100.00%
Cepheid GeneXpert - Xpert Flu	8	*Negative	8	100.00%
Consult Diag Influenza A & B	4	*Negative	4	100.00%
Quidel QuickVue A+B	12	Positive	1	8.33%
VAD-14		* Denotes Correct Response		
AG: All Methods	51			
		*Positive	51	100.00%
Alere BinaxNOW Influenza A+B Card	5	*Positive	5	100.00%
Sekisui Diagnostics OSOM Influenza A & B	8	*Positive	8	100.00%
Quidel QuickVue A+B	12	*Positive	12	100.00%
Quidel Sofia/Sofia 2 Influenza A+B	6	*Positive	6	100.00%
BD Veritor Flu A+B	5	*Positive	5	100.00%
Henry Schein OneStep + Influenza A&B (distinguish)	2	*Positive	2	100.00%
Cepheid GeneXpert - Xpert Flu	8	*Positive	8	100.00%
Alere ID NOW Influenza A & B	1	*Positive	1	100.00%
Consult Diag Influenza A & B	4	*Positive	4	100.00%
VAD-15		* Denotes Correct Response		
AG: All Methods	51			
		*Positive	51	100.00%
Alere BinaxNOW Influenza A+B Card	5	*Positive	5	100.00%
Sekisui Diagnostics OSOM Influenza A & B	8	*Positive	8	100.00%
Quidel QuickVue A+B	12	*Positive	12	100.00%
Quidel Sofia/Sofia 2 Influenza A+B	6	*Positive	6	100.00%
BD Veritor Flu A+B	5	*Positive	5	100.00%
Henry Schein OneStep + Influenza A&B (distinguish)	2	*Positive	2	100.00%
Cepheid GeneXpert - Xpert Flu	8	*Positive	8	100.00%
Alere ID NOW Influenza A & B	1	*Positive	1	100.00%
Consult Diag Influenza A & B	4	*Positive	4	100.00%
		Influenza B - Ag Detect		
VAD-11		* Denotes Correct Response		
AG: All Methods	588			
		*Negative	573	97.45%
		Positive	15	2.55%
Binax Now Flu B	2	*Negative	2	100.00%
Alere BinaxNOW Influenza A+B Card	38	*Negative	36	94.74%
Sekisui Diagnostics OSOM Influenza A & B	23	*Negative	23	100.00%
Quidel QuickVue A+B	42	*Negative	42	100.00%
Quidel Sofia/Sofia 2 Influenza A+B	200	*Negative	196	98.00%
Alere Influenza A & B	5	*Negative	5	100.00%
BD Veritor Flu A+B	65	*Negative	64	98.46%
Henry Schein OneStep + Influenza A&B (distinguish)	5	*Negative	5	100.00%



Medline Influenza A and B	3	*Negative	3	100.00%
Cepheid GeneXpert - Xpert Flu	29	*Negative	28	96.55%
Alere ID NOW Influenza A & B	26	*Negative	25	96.15%
Roche Cobas Liat Influenza A/B	6	*Negative	6	100.00%
Sekisui Diagnostics OSOM Ultra Influenza A & B	85	*Negative	80	94.12%
Consult Diag Influenza A & B	33	*Negative	32	96.97%
Lifesign Status Influenza A & B	26	*Negative	26	100.00%
Alere BinaxNOW Influenza A+B Card	38	Positive	2	5.26%
Quidel Sofia/Sofia 2 Influenza A+B	200	Positive	4	2.00%
BD Veritor Flu A+B	65	Positive	1	1.54%
Cepheid GeneXpert - Xpert Flu	29	Positive	1	3.45%
Alere ID NOW Influenza A & B	26	Positive	1	3.85%
Sekisui Diagnostics OSOM Ultra Influenza A & B	85	Positive	5	5.88%
Consult Diag Influenza A & B	33	Positive	1	3.03%
VAD-12		* Denotes Correct Response		
AG: All Methods	588			
		*Positive	573	97.45%
		Negative	15	2.55%
Binax Now Flu B	2	*Positive	2	100.00%
Alere BinaxNOW Influenza A+B Card	38	*Positive	36	94.74%
Sekisui Diagnostics OSOM Influenza A & B	23	*Positive	22	95.65%
Quidel QuickVue A+B	42	*Positive	41	97.62%
Quidel Sofia/Sofia 2 Influenza A+B	200	*Positive	195	97.50%
Alere Influenza A & B	5	*Positive	5	100.00%
BD Veritor Flu A+B	65	*Positive	64	98.46%
Henry Schein OneStep + Influenza A&B (distinguish)	5	*Positive	5	100.00%
Medline Influenza A and B	3	*Positive	3	100.00%
Cepheid GeneXpert - Xpert Flu	29	*Positive	29	100.00%
Alere ID NOW Influenza A & B	26	*Positive	25	96.15%
Roche Cobas Liat Influenza A/B	6	*Positive	6	100.00%
Sekisui Diagnostics OSOM Ultra Influenza A & B	85	*Positive	81	95.29%
Consult Diag Influenza A & B	33	*Positive	33	100.00%
Lifesign Status Influenza A & B	26	*Positive	26	100.00%
Alere BinaxNOW Influenza A+B Card	38	Negative	2	5.26%
Sekisui Diagnostics OSOM Influenza A & B	23	Negative	1	4.35%
Quidel QuickVue A+B	42	Negative	1	2.38%
Quidel Sofia/Sofia 2 Influenza A+B	200	Negative	5	2.50%
BD Veritor Flu A+B	65	Negative	1	1.54%
Alere ID NOW Influenza A & B	26	Negative	1	3.85%
Sekisui Diagnostics OSOM Ultra Influenza A & B	85	Negative	4	4.71%
VAD-13		* Denotes Correct Response		
AG: All Methods	51			
		*Positive	48	94.12%
		Negative	3	5.88%
Alere BinaxNOW Influenza A+B Card	5	*Positive	5	100.00%
Sekisui Diagnostics OSOM Influenza A & B	8	*Positive	8	100.00%
Quidel QuickVue A+B	12	*Positive	9	75.00%
Quidel Sofia/Sofia 2 Influenza A+B	6	*Positive	6	100.00%
BD Veritor Flu A+B	5	*Positive	5	100.00%
Henry Schein OneStep + Influenza A&B (distinguish)	2	*Positive	2	100.00%
Cepheid GeneXpert - Xpert Flu	8	*Positive	8	100.00%
Alere ID NOW Influenza A & B	1	*Positive	1	100.00%
Consult Diag Influenza A & B	4	*Positive	4	100.00%
Quidel QuickVue A+B	12	Negative	3	25.00%
VAD-14		* Denotes Correct Response		
AG: All Methods	46			
		*Negative	46	100.00%
Alere BinaxNOW Influenza A+B Card	1	*Negative	1	100.00%
Sekisui Diagnostics OSOM Influenza A & B	8	*Negative	8	100.00%
Quidel QuickVue A+B	12	*Negative	12	100.00%
Quidel Sofia/Sofia 2 Influenza A+B	6	*Negative	6	100.00%
BD Veritor Flu A+B	5	*Negative	5	100.00%
Henry Schein OneStep + Influenza A&B (distinguish)	2	*Negative	2	100.00%
Cepheid GeneXpert - Xpert Flu	8	*Negative	8	100.00%
Consult Diag Influenza A & B	4	*Negative	4	100.00%
VAD-15		* Denotes Correct Response		
AG: All Methods	46			
		*Negative	46	100.00%
Alere BinaxNOW Influenza A+B Card	1	*Negative	1	100.00%



Sekisui Diagnostics OSOM Influenza A & B	8	*Negative	8	100.00%
Quidel QuickVue A+B	12	*Negative	12	100.00%
Quidel Sofia/Sofia 2 Influenza A+B	6	*Negative	6	100.00%
BD Veritor Flu A+B	5	*Negative	5	100.00%
Henry Schein OneStep + Influenza A&B (distinguish)	2	*Negative	2	100.00%
Cepheid GeneXpert - Xpert Flu	8	*Negative	8	100.00%
Consult Diag Influenza A & B	4	*Negative	4	100.00%
		Influenza A/B - Ag Detect		
VAD-11		* Denotes Correct Response		
AG: All Methods	11			
		*Negative	11	100.00%
Alere BinaxNOW Influenza A+B Card	1	*Negative	1	100.00%
Quidel QuickVue	8	*Negative	8	100.00%
Henry Schein OneStep + Influenza A/B	2	*Negative	2	100.00%
VAD-12		* Denotes Correct Response		
AG: All Methods	11			
		*Positive	11	100.00%
Alere BinaxNOW Influenza A+B Card	1	*Positive	1	100.00%
Quidel QuickVue	8	*Positive	8	100.00%
Henry Schein OneStep + Influenza A/B	2	*Positive	2	100.00%
		RSV - Ag Detect		
VAD-11		* Denotes Correct Response		
AG: All Methods	298			
		*Negative	291	97.65%
		Positive	7	2.35%
BD Directigen RSV	1	*Negative	1	100.00%
BD EZ RSV	1	*Negative	1	100.00%
Alere BinaxNOW RSV	148	*Negative	144	97.30%
Quidel QuickVue RSV	35	*Negative	34	97.14%
Wampole Clearview RSV	4	*Negative	4	100.00%
Quidel QuickVue RSV 10	1	*Negative	1	100.00%
Quidel Sofia/Sofia 2 RSV	46	*Negative	45	97.83%
BD Veritor RSV	33	*Negative	32	96.97%
Cepheid GeneXpert - Xpert RSV	21	*Negative	21	100.00%
Alere ID NOW RSV	8	*Negative	8	100.00%
Alere BinaxNOW RSV	148	Positive	4	2.70%
Quidel QuickVue RSV	35	Positive	1	2.86%
Quidel Sofia/Sofia 2 RSV	46	Positive	1	2.17%
BD Veritor RSV	33	Positive	1	3.03%
VAD-12		* Denotes Correct Response		
AG: All Methods	298			
		*Negative	293	98.32%
		Positive	5	1.68%
BD Directigen RSV	1	*Negative	1	100.00%
BD EZ RSV	1	*Negative	1	100.00%
Alere BinaxNOW RSV	148	*Negative	146	98.65%
Quidel QuickVue RSV	35	*Negative	33	94.29%
Wampole Clearview RSV	4	*Negative	4	100.00%
Quidel QuickVue RSV 10	1	*Negative	1	100.00%
Quidel Sofia/Sofia 2 RSV	46	*Negative	46	100.00%
BD Veritor RSV	33	*Negative	32	96.97%
Cepheid GeneXpert - Xpert RSV	21	*Negative	21	100.00%
Alere ID NOW RSV	8	*Negative	8	100.00%
Alere BinaxNOW RSV	148	Positive	2	1.35%
Quidel QuickVue RSV	35	Positive	2	5.71%
BD Veritor RSV	33	Positive	1	3.03%
VAD-13		* Denotes Correct Response		
AG: All Methods	26			
		*Positive	26	100.00%
BD EZ RSV	1	*Positive	1	100.00%
Alere BinaxNOW RSV	10	*Positive	10	100.00%
Quidel QuickVue RSV	4	*Positive	4	100.00%
Wampole Clearview RSV	3	*Positive	3	100.00%
Quidel QuickVue RSV 10	1	*Positive	1	100.00%
Quidel Sofia/Sofia 2 RSV	2	*Positive	2	100.00%
Cepheid GeneXpert - Xpert RSV	4	*Positive	4	100.00%
Alere ID NOW RSV	1	*Positive	1	100.00%
VAD-14		* Denotes Correct Response		



AG: All Methods	26			
		*Negative	26	100.00%
BD EZ RSV	1	*Negative	1	100.00%
Alere BinaxNOW RSV	10	*Negative	10	100.00%
Quidel QuickVue RSV	4	*Negative	4	100.00%
Wampole Clearview RSV	3	*Negative	3	100.00%
Quidel QuickVue RSV 10	1	*Negative	1	100.00%
Quidel Sofia/Sofia 2 RSV	2	*Negative	2	100.00%
Cepheid GeneXpert - Xpert RSV	4	*Negative	4	100.00%
Alere ID NOW RSV	1	*Negative	1	100.00%
VAD-15		* Denotes Correct Response		
AG: All Methods	26			
		*Positive	25	96.15%
		Negative	1	3.85%
BD EZ RSV	1	*Positive	1	100.00%
Alere BinaxNOW RSV	10	*Positive	10	100.00%
Quidel QuickVue RSV	4	*Positive	3	75.00%
Wampole Clearview RSV	3	*Positive	3	100.00%
Quidel QuickVue RSV 10	1	*Positive	1	100.00%
Quidel Sofia/Sofia 2 RSV	2	*Positive	2	100.00%
Cepheid GeneXpert - Xpert RSV	4	*Positive	4	100.00%
Alere ID NOW RSV	1	*Positive	1	100.00%
Quidel QuickVue RSV	4	Negative	1	25.00%
		Adenovirus		
VAD-11		* Denotes Correct Response		
AG: All Methods	10			
		*Negative	10	100.00%
SAS Adeno Test	10	*Negative	10	100.00%
VAD-12		* Denotes Correct Response		
AG: All Methods	10			
		*Positive	10	100.00%
SAS Adeno Test	10	*Positive	10	100.00%
VAD-13		* Denotes Correct Response		
AG: All Methods	10			
		*Negative	10	100.00%
SAS Adeno Test	10	*Negative	10	100.00%
VAD-14		* Denotes Correct Response		
AG: All Methods	10			
		*Positive	10	100.00%
SAS Adeno Test	10	*Positive	10	100.00%
VAD-15		* Denotes Correct Response		
AG: All Methods	10			
		*Negative	10	100.00%
SAS Adeno Test	10	*Negative	10	100.00%
932 - Hepatitis C Antibody (Waived);767 - Hepatitis Markers	Labs	Response	Number	Percent
		HCV (Waived)		
VM-11		* Denotes Correct Response		
AG: All Methods	4			
		*Negative	4	100.00%
OraSure OraQuick HCV Rapid Antibody Test	4	*Negative	4	100.00%
VM-12		* Denotes Correct Response		
AG: All Methods	4			
		*Positive	4	100.00%
OraSure OraQuick HCV Rapid Antibody Test	4	*Positive	4	100.00%
		HBsAg		
VM-11		* Denotes Correct Response		
AG: All Methods	2			
		*Negative	2	100.00%
Roche Cobas e Series	1	*Negative	1	100.00%
Abbott Architect ci4100	1	*Negative	1	100.00%
VM-12		* Denotes Correct Response		
AG: All Methods	1			
		*Positive	1	100.00%
Roche Cobas e Series	1	*Positive	1	100.00%
767 - Hepatitis Markers	Labs	Response	Number	Percent
		HBsAg		
VM-11		* Denotes Correct Response		



AG: All Methods	2			
		*Negative	2	100.00%
Roche Cobas e Series	1	*Negative	1	100.00%
Abbott Architect ci4100	1	*Negative	1	100.00%
VM-12		* Denotes Correct Response		
AG: All Methods	1			
		*Positive	1	100.00%
Roche Cobas e Series	1	*Positive	1	100.00%
VM-13		* Denotes Correct Response		
AG: All Methods	1			
		*Positive	1	100.00%
Roche Cobas e Series	1	*Positive	1	100.00%
VM-14		* Denotes Correct Response		
AG: All Methods	2			
		*Negative	2	100.00%
Roche Cobas e Series	1	*Negative	1	100.00%
Abbott Architect ci4100	1	*Negative	1	100.00%
VM-15		* Denotes Correct Response		
AG: All Methods	2			
		*Negative	2	100.00%
Roche Cobas e Series	1	*Negative	1	100.00%
Abbott Architect ci4100	1	*Negative	1	100.00%
		HCV		
VM-11		* Denotes Correct Response		
AG: All Methods	1			
		*Negative	1	100.00%
Abbott Architect ci4100	1	*Negative	1	100.00%
VM-12		* Denotes Correct Response		
AG: All Methods	1			
		*Positive	1	100.00%
Abbott Architect ci4100	1	*Positive	1	100.00%
VM-13		* Denotes Correct Response		
AG: All Methods	1			
		*Negative	1	100.00%
Abbott Architect ci4100	1	*Negative	1	100.00%
VM-14		* Denotes Correct Response		
AG: All Methods	1			
		*Negative	1	100.00%
Abbott Architect ci4100	1	*Negative	1	100.00%
VM-15		* Denotes Correct Response		
AG: All Methods	1			
		*Positive	1	100.00%
Abbott Architect ci4100	1	*Positive	1	100.00%
		Anti-HBc IgM		
VM-11		* Denotes Correct Response		
AG: All Methods	1			
		*Negative	1	100.00%
Abbott Architect ci4100	1	*Negative	1	100.00%
VM-12		* Denotes Correct Response		
AG: All Methods	1			
		*Negative	1	100.00%
Abbott Architect ci4100	1	*Negative	1	100.00%
VM-13		* Denotes Correct Response		
AG: All Methods	1			
		*Negative	1	100.00%
Abbott Architect ci4100	1	*Negative	1	100.00%
VM-14		* Denotes Correct Response		
AG: All Methods	1			
		*Negative	1	100.00%
Abbott Architect ci4100	1	*Negative	1	100.00%
VM-15		* Denotes Correct Response		
AG: All Methods	1			
		*Negative	1	100.00%
Abbott Architect ci4100	1	*Negative	1	100.00%
		HAV IgM		
VM-11		* Denotes Correct Response		
AG: All Methods	1			
		*Negative	1	100.00%



Abbott Architect ci4100	1	*Negative	1	100.00%
VM-12		* Denotes Correct Response		
AG: All Methods	1			
		*Negative	1	100.00%
Abbott Architect ci4100	1	*Negative	1	100.00%
VM-13		* Denotes Correct Response		
AG: All Methods	1			
		*Positive	1	100.00%
Abbott Architect ci4100	1	*Positive	1	100.00%
VM-14		* Denotes Correct Response		
AG: All Methods	1			
		*Negative	1	100.00%
Abbott Architect ci4100	1	*Negative	1	100.00%
VM-15		* Denotes Correct Response		
AG: All Methods	1			
		*Negative	1	100.00%
Abbott Architect ci4100	1	*Negative	1	100.00%
		Anti-HBc Total (International)		
VM-11		* Denotes Correct Response		
AG: All Methods	2			
		*Negative	2	100.00%
Intl OCD Vitros ECi/ECiQ	2	*Negative	2	100.00%
VM-12		* Denotes Correct Response		
AG: All Methods	2			
		*Positive	2	100.00%
Intl OCD Vitros ECi/ECiQ	2	*Positive	2	100.00%
VM-13		* Denotes Correct Response		
AG: All Methods	2			
		*Negative	2	100.00%
Intl OCD Vitros ECi/ECiQ	2	*Negative	2	100.00%
VM-14		* Denotes Correct Response		
AG: All Methods	2			
		*Negative	2	100.00%
Intl OCD Vitros ECi/ECiQ	2	*Negative	2	100.00%
VM-15		* Denotes Correct Response		
AG: All Methods	2			
		*Negative	2	100.00%
Intl OCD Vitros ECi/ECiQ	2	*Negative	2	100.00%
		HBsAb (International)		
VM-11		* Denotes Correct Response		
AG: All Methods	2			
		*Positive	2	100.00%
Intl OCD Vitros ECi/ECiQ	2	*Positive	2	100.00%
VM-12		* Denotes Correct Response		
AG: All Methods	2			
		*Negative	2	100.00%
Intl OCD Vitros ECi/ECiQ	2	*Negative	2	100.00%
VM-13		* Denotes Correct Response		
AG: All Methods	2			
		*Negative	2	100.00%
Intl OCD Vitros ECi/ECiQ	2	*Negative	2	100.00%
VM-14		* Denotes Correct Response		
AG: All Methods	2			
		*Negative	2	100.00%
Intl OCD Vitros ECi/ECiQ	2	*Negative	2	100.00%
VM-15		* Denotes Correct Response		
AG: All Methods	2			
		*Negative	2	100.00%
Intl OCD Vitros ECi/ECiQ	2	*Negative	2	100.00%
		HBsAg (International)		
VM-11		* Denotes Correct Response		
AG: All Methods	2			
		*Negative	2	100.00%
Intl OCD Vitros ECi/ECiQ	2	*Negative	2	100.00%
VM-12		* Denotes Correct Response		
AG: All Methods	2			
		*Positive	2	100.00%



Intl OCD Vitros ECi/ECiQ	2	*Positive	2	100.00%
VM-13		* Denotes Correct Response		
AG: All Methods	2			
		*Positive	2	100.00%
Intl OCD Vitros ECi/ECiQ	2	*Positive	2	100.00%
VM-14		* Denotes Correct Response		
AG: All Methods	2			
		*Negative	2	100.00%
Intl OCD Vitros ECi/ECiQ	2	*Negative	2	100.00%
VM-15		* Denotes Correct Response		
AG: All Methods	2			
		*Negative	1	50.00%
		Positive	1	50.00%
Intl OCD Vitros ECi/ECiQ	2	*Negative	1	50.00%
Intl OCD Vitros ECi/ECiQ		Positive	1	50.00%
		HCV (International)		
VM-11		* Denotes Correct Response		
AG: All Methods	2			
		*Negative	2	100.00%
Intl OCD Vitros ECi/ECiQ	2	*Negative	2	100.00%
VM-12		* Denotes Correct Response		
AG: All Methods	2			
		*Positive	2	100.00%
Intl OCD Vitros ECi/ECiQ	2	*Positive	2	100.00%
VM-13		* Denotes Correct Response		
AG: All Methods	2			
		*Negative	2	100.00%
Intl OCD Vitros ECi/ECiQ	2	*Negative	2	100.00%
VM-14		* Denotes Correct Response		
AG: All Methods	2			
		*Negative	2	100.00%
Intl OCD Vitros ECi/ECiQ	2	*Negative	2	100.00%
VM-15		* Denotes Correct Response		
AG: All Methods	2			
		*Positive	2	100.00%
Intl OCD Vitros ECi/ECiQ	2	*Positive	2	100.00%
771 - Vaginal Pathogens I	Labs	Response	Number	Percent
		Candida sp. - Non Culture		
VP-11		* Denotes Correct Response		
AG: All Methods	49			
		*Negative	49	100.00%
Affirm VP	48	*Negative	48	100.00%
BD Max	1	*Negative	1	100.00%
VP-12		* Denotes Correct Response		
AG: All Methods	49			
		*Negative	48	97.96%
		Positive	1	2.04%
Affirm VP	48	*Negative	48	100.00%
BD Max	1	Positive	1	100.00%
VP-13		* Denotes Correct Response		
AG: All Methods	49			
		*Negative	49	100.00%
Affirm VP	48	*Negative	48	100.00%
BD Max	1	*Negative	1	100.00%
VP-14		* Denotes Correct Response		
AG: All Methods	49			
		*Positive	49	100.00%
Affirm VP	48	*Positive	48	100.00%
BD Max	1	*Positive	1	100.00%
VP-15		* Denotes Correct Response		
AG: All Methods	49			
		*Negative	49	100.00%
Affirm VP	48	*Negative	48	100.00%
BD Max	1	*Negative	1	100.00%
		G. vaginalis - Non Culture		
VP-11		* Denotes Correct Response		
AG: All Methods	49			
		*Negative	49	100.00%



American Academy of Family Physicians
11400 Tomahawk Creek Parkway
Leawood, KS 66211-2672
1-800-274-7911
www.aafp.org/pt pt@aafp.org

Proficiency Testing Event Summary

AAFP PT 2019-C

Affirm VP	48	*Negative	48	100.00%
BD Max	1	*Negative	1	100.00%
VP-12		* Denotes Correct Response		
AG: All Methods	49			
		*Negative	49	100.00%
Affirm VP	48	*Negative	48	100.00%
BD Max	1	*Negative	1	100.00%
VP-13		* Denotes Correct Response		
AG: All Methods	49			
		*Negative	49	100.00%
Affirm VP	48	*Negative	48	100.00%
BD Max	1	*Negative	1	100.00%
VP-14		* Denotes Correct Response		
AG: All Methods	49			
		*Negative	49	100.00%
Affirm VP	48	*Negative	48	100.00%
BD Max	1	*Negative	1	100.00%
VP-15		* Denotes Correct Response		
AG: All Methods	49			
		*Positive	49	100.00%
Affirm VP	48	*Positive	48	100.00%
BD Max	1	*Positive	1	100.00%
		T. vaginalis - Non Culture Tech		
VP-11		* Denotes Correct Response		
AG: All Methods	49			
		*Positive	49	100.00%
Affirm VP	48	*Positive	48	100.00%
BD Max	1	*Positive	1	100.00%
VP-12		* Denotes Correct Response		
AG: All Methods	49			
		*Positive	49	100.00%
Affirm VP	48	*Positive	48	100.00%
BD Max	1	*Positive	1	100.00%
VP-13		* Denotes Correct Response		
AG: All Methods	49			
		*Negative	49	100.00%
Affirm VP	48	*Negative	48	100.00%
BD Max	1	*Negative	1	100.00%
VP-14		* Denotes Correct Response		
AG: All Methods	49			
		*Negative	49	100.00%
Affirm VP	48	*Negative	48	100.00%
BD Max	1	*Negative	1	100.00%
VP-15		* Denotes Correct Response		
AG: All Methods	49			
		*Negative	49	100.00%
Affirm VP	48	*Negative	48	100.00%
BD Max	1	*Negative	1	100.00%